




SpyGlass™ Retrieval Basket for use with *SpyGlass™ DS II Direct Visualization System* *featuring Sight Shield Technology*

Boston
Scientific

The **SpyGlass DS II Direct Visualization System** with **Sight Shield Technology** is the complete direct imaging system designed for pancreaticobiliary stone and stricture management of the bile and pancreatic ducts during endoscopic retrograde cholangiopancreatography (ERCP).



5-15%

5-15% of CBD stone cases cannot be treated effectively using standard ERCP techniques¹

34%

34% of patients were found to have residual biliary stones NOT seen in occlusion cholangiogram but were detected with SpyGlass DS.

The SpyGlass Retrieval Basket is used with the SpyGlass DS System to capture and remove residual biliary and pancreatic stones. Without their removal, residual stones can possibly lead to patient discomfort and follow-up ERCP's – impacting procedure time and patient experience.

Made from nitinol wires, the SpyGlass Retrieval Basket has an open diameter of 15mm. The basket also features a Mediglide™ coating to help decrease friction and reduce the force required to deploy and retract the basket.*


Results with the SpyGlass DS II Direct Visualization System include a demonstrated 89% stone clearance after a single ERCP session², which may allow patients to avoid repeat endoscopic interventions or more invasive procedures to treat difficult biliary and pancreatic stones.



Ordering Information

SpyGlass™ DS System


Order Number	Product Description	Packaging
M00546650	SpyGlass™ DS Digital Controller	Each
M00546610	SpyScope™ DS II Access & Delivery Catheter	Each



Capturing and removing stone using
SpyGlass Retrieval Basket.

SpyGlass Retrieval Basket

Order Number	Product Description	Cable Diameter (in/mm)	Jaw Outer Diameter (mm)	Working Length (cm)	Required Endoscope Working Channel (mm)
M00546550	SpyGlass Retrieval Basket	0.039 / 1.0	N/A	286	1.2



*Data on file

1. Parsi et al. Endoscopic management of difficult common bile duct stones. World J Gastroenterol 2013; 19(2): 165-173
2. Shah et al. Performance of a full disposable, digital, single-operator cholangiopancreatroscope. Endoscopy 2017; 49; 651-658.

All photographs taken by Boston Scientific.
All trademarks are the property of their respective owners.
CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.
Indications, contraindications, warnings, and instructions for use can be found in the
product labelling supplied with each device or at www.IFU-BSCI.com.