

Endoscopy Product Catalogue

**Boston
Scientific**
Advancing science for life™

Produktkatalog – Catalogue Produits
Catálogo de Productos – Catalogo dei Prodotti

2020



► RX Biliary System™

► Biliary

► SpyGlass™

► EUS

► Stents

► Dilation

► Guidewires



► Hemostasis

► Tissue Resection

► Biopsy

► Enteral Feeding

► Pulmonary

► Infection Prevention

► Radiofrequency
Ablation

Updated 02/2020



Contents



Alphabetical
Index

Product Code
Index

RX Biliary System™

Cannulas	Rapid Exchange XL Cannula	39
	RX ERCP Cannulas	39
Sphincterotomes	Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	40
	Dreamtome™ RX Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	40
	Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	41
	Jagtome™ Revolution RX Cannulating Sphincterotome (Autotome™ Preloaded with Jagwire Revolution High Performance Guidewire)	41
	Autotome™ RX Cannulating Sphincterotomes	42
	RX Needle-Knife XL Triple-Lumen Needle-Knife	42
Stone Retrieval	Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Distal Injection	43
	Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	43
	Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Distal Injection	44
	Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	44
	Trapezoid™ RX Wireguided Retrieval Baskets	44
Dilation	Hurricane™ RX Biliary Balloon Dilatation Catheters	45
	CRE™ RX Biliary Balloon Dilatation Catheter	45
Cytology	RX Cytology Brush	45
Accessories	RX Accessories	46
	Encore™ 26 Inflation Device	46

Biliary

Cannulas	Tandem™ XL Triple-Lumen ERCP Cannula	47
	Contour™ ERCP Cannulas	47
Sphincterotomes	Stonetome™ Stone Removal Devices	47
	TRUEtome™ Jag Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	48
	TRUEtome™ Dreamwire™ Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	48
	TRUEtome™ Hydra Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	48
	TRUEtome™ Jagwire Revolution Cannulating Sphincterotomes with Pre-loaded Jagwire™ Revolution Guidewire	48
	TRUEtome™ Cannulating Sphincterotomes	49
	Ultratome XL™ Triple-Lumen Sphincterotomes	50
	Ultratome™ Double-Lumen Sphincterotomes	50
	Microknife™ XL Triple-Lumen Needle-Knife	50
Stone Retrieval	Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Distal Injection	51
	Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	51
	Extractor™ Pro DL Double-Lumen Retrieval Balloon Catheters	51
Dilation	MaxForce™ Biliary Balloon Dilatation Catheter	52
Accessories	Biliary Accessories	52

SpyGlass™ Direct Visualisation Systems

SpyGlass™ DS II	SpyGlass™ DS Capital Equipment	53
	SpyGlass™ DS Disposable Equipment	53
	Autolith™ Touch Biliary EHL System	53
	Service Codes	54

Endoscopic Ultrasound (EUS)

Fine Needle Aspiration (FNA)	Expect™ Standard Handle Needles	55
	Expect™ Slimline Handle Needles	55
Fine Needle Biopsy (FNB)	Acquire™ EUS-FNB Devices	56
Stent and Delivery System	Hot AXIOS™ Stent and Electrocautery Enhanced Delivery System	56
	AXIOS™ Stent and Delivery System	56

Stents

RX Biliary Metal Stents	WallFlex™ Biliary RX Fully Covered Stent System RMV	57
	WallFlex™ Biliary RX Partially Covered Stent System (Permalume™ Covering)	57
	WallFlex™ Biliary RX Uncovered Stent System	58
	Wallstent™ RX Biliary Endoprosthesis Stent System (Permalume™ Covering)	58
	Wallstent™ RX Biliary Endoprosthesis Stent System (Uncovered)	59
Biliary Metal Stents	Wallstent™ Endoscopic Biliary Endoprosthesis Stent System (Partially Covered)	59
	Wallstent™ Endoscopic Biliary Endoprosthesis Stent System (Uncovered)	59
RX Biliary Plastic Stents	Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System	60
	Advanix™ Stents – Singles	61–62
	NaviFlex™ RX Delivery System	62
	RX Biliary Plastic Stents with RX Delivery System	63
Biliary Plastic Stents	Flexima™ Biliary Stents with Delivery System	64
	Flexima™ Nasobiliary Catheter	64
	Percuflex™ Duodenal Bend Biliary Stents	65
	Biliary Introducer	65
	C-flex™ Pigtail Biliary Stents	65
Pancreatic Plastic Stents	Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Pigtail	66
	Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Straight	66
	Advanix™ Pancreatic Single Stents – Pigtail	67–68
	Advanix™ Pancreatic Single Stents – Straight	69–70
	NaviFlex™ RX Pancreatic Delivery System & NaviFlex™ Pushers	71
Esophageal Stents	WallFlex™ Fully Covered RMV Esophageal Stent System	72
	WallFlex™ Fully Covered RMV Esophageal Longer Loop Stent System	72
	WallFlex™ Partially Covered Esophageal Stent System	72
	Ultraflex™ Covered Esophageal NG Stent System – Proximal Release	73
	Ultraflex™ Covered Esophageal NG Stent System – Distal Release	73
	Ultraflex™ Uncovered Esophageal NG Stent System – Proximal Release	73
	Ultraflex™ Uncovered Esophageal NG Stent System – Distal Release	73
Colonic and Duodenal Stents	WallFlex™ Colonic Stent System with Anchor Lock Delivery System	74
	WallFlex™ Colonic Soft Stent System with Anchor Lock Delivery System	74
	WallFlex™ Duodenal Stent System with Anchor Lock Delivery System	75
	WallFlex™ Duodenal Soft Stent System with Anchor Lock Delivery System	75
	Wallstent™ Enteral Endoprosthesis Colonic/Duodenal Stent with UNISTEP™ Plus Delivery System	76
	Ultraflex Precision™ Colonic Stent System – Proximal Release	76

Dilation

CRE™ (Controlled Radial Expansion) Balloon Dilatation Catheters	CRE™ PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilatation Catheters (including DASE indication)	77
	CRE™ PRO Wireguided Esophageal/Pyloric/Colonic/Biliary Balloon Dilatation Catheters (including DASE indication)	77
	CRE™ Fixed Wire Balloon Dilatation Catheters	78
MaxForce TTS™ Balloon Dilators	MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter	79
	Rigiflex™ II Achalasia Balloon Dilators	79
Inflation Systems	Alliance™ II Inflation System	80
Dilation Accessories	Inflation/Deflation Devices and Accessories	80

Guidewires

High Performance	Dreamwire™ High Performance Guidewires	81
	Hydra Jagwire™ High Performance Guidewires	81
	Jagwire™ High Performance Guidewires	82
	Jagwire™ High Performance Guidewires – Extendable Jagwire	82
	Jagwire™ Revolution High Performance Guidewire	82
	NovaGold™ High Performance Guidewires	83
	NaviPro™ Hydrophilic Guidewires	83
Standard	Amplatz Super Stiff™ Guidewire	83
	Wallstent™ Super Stiff Guidewire	83

Hemostasis

Clipping Devices	Resolution 360™ Clips	85
	Resolution™ Clips	85
	Speedband Superview Super 7™ Multiple Band Ligators	85
Bipolar Electrohemostasis Catheters	Injection Gold Probe™ Electrohemostasis/Injection Catheters	86
	Gold Probe™ Electrohemostasis Catheters	86
Needles	Interject™ Sclerotherapy Needles	87
	Interject™ Clear Sclerotherapy Needles	87
Electrosurgical Bipolar Cable Adaptor	Electrosurgical Bipolar Cable Adaptor	88

Tissue Resection

EMR Devices for Upper GI	Captivator™ EMR	89
Snares	Rotatable Single-Use Snares	89
	Sensation™ Single-Use Short Throw Snares	90
	Captiflex™ Single-Use Snares	90
	Captivator™ Single-Use Snares	91
	Captivator™ II Single-Use Snares	91
	Captivator™ COLD Single-Use Snare	91
	Profile™ Single-Use Snares	92
Hot Biopsy Forceps	Radial Jaw™ 4 Hot Biopsy Forceps	92
Accessories	Active Cords	92
Retrieval Devices	RescueNet™	93
	Rescue™ Forceps	93
Polyp Traps	TrapEase™ Polyp Traps	93
Lifting Agents	Orise™ Gel	93

Biopsy

Forceps	Radial Jaw™ 4 Standard Capacity Biopsy Forceps	95
	Radial Jaw™ 4 Large Capacity Biopsy Forceps	95
	Radial Jaw™ 4 Jumbo Biopsy Forceps	96
	Radial Jaw™ 4 Gastropediatric Biopsy Forceps	96
	Multibite™ Multiple Sample Biopsy Forceps	96

EndoVive™ Enteral Access Devices

Initial Placement	Standard PEG Kits	97
	Standard PEG Kits with ENFit™	97
	PEG Accessories	97
	PEG Accessories with ENFit™	97
	One-Step Button™ Low Profile PEG Kits	98
	One-Step Button™ Low Profile PEG Kits with ENFit™	98
	One-Step Button™ and Low Profile Replacement Button Accessories	99
	One-Step Button™ and Low Profile Replacement Button Accessories with ENFit™	99
Replacements	Replacement Gastrostomy Feeding Tubes with ENFit™	100
	Low Profile Replacement Buttons	100
Replacement Feeding Tubes	Through-The-PEG (TTP) Jejunal Feeding Tube	101
	Replacement Gastrostomy Feeding Tubes	101

Pulmonary

Alair™ Bronchial Thermoplasty (BT) System	Alair™ Bronchial Thermoplasty (BT) System	103
	Alair™ BT Controller Accessory Kits ATS 201	103
Endoscopic Bronchial Ultrasound	Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needles – Olympus Scope Compatible	103
	Acquire™ Pulmonary Endobronchial Ultrasound Fine Biopsy Needles– Olympus Scope Compatible	104
Airway Stents	Ultraflex™ Tracheobronchial Covered Stent System – Distal Release	104
	Ultraflex™ Tracheobronchial Uncovered Stent System – Proximal Release	105
	Ultraflex™ Tracheobronchial Uncovered Stent System – Distal Release	105
Biopsy Forceps	Radial Jaw™ 4 Pulmonary Biopsy Forceps	106
Mini-Forceps	CoreDx™ Pulmonary Mini-Forceps	106
Transbronchial Aspiration Needles	eXcelon™ Transbronchial Aspiration Needles	106
Cytology Brushes	Cellebri™ Cytology Brushes – Bronchoscope	107
Airway Retrieval Baskets	Zero Tip™ Airway Retrieval Baskets	107
Dilation Balloons	CRE™ Pulmonary Balloon Dilatation Catheter	107
Guidewires	Jagwire™ Pulmonary Guidewire	108
Inflation System	Alliance™ II Inflation System	108

Infection Prevention

Irrigation & Insufflation	Hydra™ Irrigation System	109
Cleaning Brushes	HedgeHog™ Cleaning Brushes	109
Valves	Orca™ Single Use Air/Water and Suction Valves	109
Scope Transport	CinchPad™ Procedure & Endoscope Transport Pad	110
Bite Blocks	Blox™ Bite Blocks	110
Sponges	Koala™ Sponge	110
Procedure Gowns	EndoArmor™ Procedural Gowns	110
Lubricating Gels	EndoGlide™ Lubricating Gels	110

Radiofrequency Ablation

Ablation Catheter	Habib™ EndoHPB Bipolar Radiofrequency Catheter	111
--------------------------	--	-----

Alphabetical Index

Acquire™ EUS-FNB Devices	56
Acquire™ Pulmonary Endobronchial Ultrasound Fine Biopsy Needles – Olympus Scope Compatible	104
Active Cords	92
Advanix™ Pancreatic Single Stents – Pigtail	67 – 68
Advanix™ Pancreatic Single Stents – Straight	69 – 70
Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Pigtail	66
Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Straight	66
Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System	60
Advanix™ Stents – Singles	61 – 62
Alair™ Bronchial Thermoplasty (BT) System	103
Alair™ BT Controller Accessory Kits ATS 201	103
Alliance™ II Inflation System	44, 80, 108
Amplatz Super Stiff™ Guidewire	83
Autolith™ Touch Biliary EHL System	53
Autotome™ RX Cannulating Sphincterotomes	42
AXIOS™ Stent and Delivery System	56
Biliary Accessories	52
Biliary Introducer	65
Blox™ Bite Blocks	110
C-flex™ Pigtail Biliary Stents	65
Captiflex™ Single-Use Snares	90
Captivator™ COLD Single-Use Snare	91
Captivator™ EMR	89
Captivator™ Single-Use Snares	91
Captivator™ II Single-Use Snares	91
Cellebri™ Cytology Brushes – Bronchoscope	107
CinchPad™ Procedure & Endoscope Transport Pad	110
Contour™ ERCP Cannulas	47
CoreDx™ Pulmonary Mini-Forceps	106
CRE™ Fixed Wire Balloon Dilatation Catheters	78
CRE™ PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilatation Catheters (including DASE indication)	77
CRE™ PRO Wireguided/Esophageal/Pyloric/Colonic Biliary Balloon Dilatation Catheters (including DASE indication)	77
CRE™ Pulmonary Balloon Dilatation Catheter	107
CRE™ RX Biliary Balloon Dilatation Catheter	45
Dreamtome™ RX Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	40
Dreamwire™ High Performance Guidewires	81
Electrosurgical Bipolar Cable Adaptor	88
Encore™ 26 Inflation Device	46
EndoArmor™ Procedural Gowns	110
EndoGlide™ Lubricating Gels	110
eXcelon™ Transbronchial Aspiration Needles	106
Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needles – Olympus Scope Compatible	103
Expect™ Slimline Handle Needles	55
Expect™ Standard Handle Needles	55
Extractor™ Pro DL Double-Lumen Retrieval Balloon Catheters	51
Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Distal Injection	43
Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	43
Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Distal Injection	44
Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	44
Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Distal Injection	51



Alphabetical Index (CONTINUED)

Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	51
Flexima™ Biliary Stents with Delivery System	64
Flexima™ Nasobiliary Catheter	64
Gold Probe™ Electrohemostasis Catheters	86
Habib™ EndoHPB Bipolar Radiofrequency Catheter	111
HedgeHog™ Cleaning Brushes	109
Hot AXIOS™ Stent and Electrocautery Enhanced Delivery System	56
Hurricane™ RX Biliary Dilation Balloons	45
Hydra™ Irrigation System	109
Hydra Jagwire™ High Performance Guidewires	81
Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	40
Inflation/Deflation Devices and Accessories	80
Injection Gold Probe™ Electrohemostasis/Injection Catheters	86
Interject™ Clear Sclerotherapy Needles	87
Interject™ Sclerotherapy Needles	87
Jagtome™ Revolution RX Cannulating Sphincterotome (Autotome™ Preloaded with Jagwire Revolution High Performance Guidewire)	41
Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	41
Jagwire™ High Performance Guidewires	82
Jagwire™ High Performance Guidewires – Extendable Jagwire	82
Jagwire™ Pulmonary Guidewire	108
Jagwire™ Revolution High Performance Guidewire	82
Koala™ Sponge	110
Low Profile Replacement Button	100
MaxForce™ Biliary Balloon Dilatation Catheter	52
MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter	79
Microknife™ XL Triple-Lumen Needle-Knife	50
Multibite™ Multiple Sample Biopsy Forceps	96
NaviFlex™ RX Delivery System	62
NaviFlex™ RX Pancreatic Delivery System & NaviFlex™ Pushers	71
NaviPro™ Hydrophilic Guidewires	83
NovaGold™ High Performance Guidewires	83
One-Step Button™ and Low Profile Replacement Button Accessories	99
One-Step Button™ and Low Profile Replacement Button Accessories with ENFit™	99
One-Step Button™ Low Profile PEG Kits	98
One-Step Button™ Low Profile PEG Kits with ENFit™	98
Orca™ Single Use Air/Water and Suction Valves	109
Orise™ Gel	93
PEG Accessories	97
PEG Accessories with ENFit™	97
Percuflex™ Duodenal Bend Biliary Stents	65
Profile™ Single-Use Snares	92
Radial Jaw™ 4 Gastropediatric Biopsy Forceps	96
Radial Jaw™ 4 Hot Biopsy Forceps	92
Radial Jaw™ 4 Jumbo Biopsy Forceps	96
Radial Jaw™ 4 Large Capacity Biopsy Forceps	95
Radial Jaw™ 4 Pulmonary Biopsy Forceps	106
Radial Jaw™ 4 Standard Capacity Biopsy Forceps	95
Rapid Exchange XL Cannula	39
Replacement Gastrostomy Feeding Tubes	101
Replacement Gastrostomy Feeding Tubes with ENFit™	100



Alphabetical Index (CONTINUED)

RescueNet™	93
Rescue™ Forceps	93
Resolution 360™ Clips	85
Resolution™ Clips	85
Rigiflex™ II Achalasia Balloon Dilators	79
Rotatable Single-Use Snares	89
RX Accessories	46
RX Biliary Plastic Stents with RX Delivery System	63
RX Cytology Brush	45
RX ERCP Cannulas	39
RX Needle-Knife XL Triple-Lumen Needle-Knife	42
Sensation™ Single-Use Short Throw Snares	90
Speedband Superview Super 7™ Multiple Band Ligators	85
SpyGlass™ DS Capital Equipment	53
SpyGlass™ DS Disposable Equipment	53
SpyGlass™ Service Codes	54
Standard PEG Kits	97
Standard PEG Kits with ENFit™	97
Stonetome™ Stone Removal Devices	47
Tandem™ XL Triple-Lumen ERCP Cannula	47
Through-The-PEG (TTP) Jejunal Feeding Tube	101
TrapEase™ Polyp Traps	93
Trapezoid™ RX Wireguided Retrieval Baskets	44
TRUEtome™ Cannulating Sphincterotomes	49
TRUEtome™ Dreamwire™ Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	48
TRUEtome™ Hydra Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	48
TRUEtome™ Jag Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	48
TRUEtome™ Jagwire Revolution Cannulating Sphincterotomes with Pre-loaded Jagwire™ Revolution Guidewire	48
Ultraflex™ Covered Esophageal NG Stent System – Distal Release	73
Ultraflex™ Covered Esophageal NG Stent System – Proximal Release	73
Ultraflex Precision™ Colonic Stent System – Proximal Release	76
Ultraflex™ Tracheobronchial Covered Stent System – Distal Release	104
Ultraflex™ Tracheobronchial Uncovered Stent System – Distal Release	105
Ultraflex™ Tracheobronchial Uncovered Stent System – Proximal Release	105
Ultraflex™ Uncovered Esophageal NG Stent System – Distal Release	73
Ultraflex™ Uncovered Esophageal NG Stent System – Proximal Release	73
Ultratome™ Double-Lumen Sphincterotomes	50
Ultratome XL™ Triple-Lumen Sphincterotomes	50
WallFlex™ Biliary RX Fully Covered Stent System RMV	57
WallFlex™ Biliary RX Partially Covered Stent System	57
WallFlex™ Biliary RX Uncovered Stent System	58
WallFlex™ Colonic Soft Stent System with Anchor Lock Delivery System	74
WallFlex™ Colonic Stent System with Anchor Lock Delivery System	74
WallFlex™ Duodenal Soft Stent System with Anchor Lock Delivery System	75
WallFlex™ Duodenal Stent System with Anchor Lock Delivery System	75
WallFlex™ Fully Covered RMV Esophageal Longer Loop Stent System	72
WallFlex™ Fully Covered RMV Esophageal Stent System	72
WallFlex™ Partially Covered Esophageal Stent System	72
Wallstent™ Enteral Endoprosthesis Colonic/Duodenal Stent with UNISTEP™ Plus Delivery System	76
Wallstent™ Endoscopic Biliary Endoprosthesis Stent System (Partially Covered)	59



Alphabetical Index (CONTINUED)

Wallstent™ Endoscopic Biliary Endoprosthesis Stent System (Uncovered)	59
Wallstent™ RX Biliary Endoprosthesis Stent System (Permalume™ Covering)	58
Wallstent™ RX Biliary Endoprosthesis Stent System (Uncovered)	59
Wallstent™ Super Stiff Guidewire	83
Zero Tip™ Airway Retrieval Baskets	107

RX Biliary System™	Biliary	SpyGlass™	EUS	Stents	Dilation	Guidewires	Hemostasis	Tissue Resection	Biopsy	Enteral Feeding	Pulmonary	Infection Prevention	RFA
--------------------	---------	-----------	-----	--------	----------	------------	------------	------------------	--------	-----------------	-----------	----------------------	-----



Product Code Index

Order Number	Description	Page Number
	Rescue™ Forceps	
DGA-274-5	Rescue Alligator Jaw Grasping Forceps	93
DGA-275-5	Rescue Alligator Jaw Grasping Forceps (Long)	93
DGC-383-5	Rescue Combo Grasping Forceps	93
	RescueNet™	
DGN-538-5	RescueNet Retrieval Device	93
	Rescue™ Forceps	
DGR-276-5	Rescue Rat Tooth Jaw Grasping Forceps	93
DGR-277-5	Rescue Rat Tooth Jaw Grasping Forceps (Pediatric)	93
	Wallstent™ Super Stiff Guidewire	
H965180011	Wallstent Super Stiff Guidewire	83
	Wallstent™ Endoscopic Biliary Endoprosthesis Stent System (Uncovered)	
H965430100	Wallstent Biliary	59
H965430200	Wallstent Biliary	59
H965430300	Wallstent Biliary	59
H965430400	Wallstent Biliary	59
H965430500	Wallstent Biliary	59
H965430600	Wallstent Biliary	59
	Wallstent™ Endoscopic Biliary Endoprosthesis Stent System (Partially Covered)	
H965430700	Wallstent Biliary w/Permalume	59
H965430800	Wallstent Biliary w/Permalume	59
H965430900	Wallstent Biliary w/Permalume	59
H965431000	Wallstent Biliary w/Permalume	59
H965431100	Wallstent Biliary w/Permalume	59
H965431200	Wallstent Biliary w/Permalume	59
	Alair™ Bronchial Thermoplasty (BT) System	
M005ATS20000	Alair RF Controller – monopolar (return electrode not included)	103
	Alair™ BT Controller Accessory Kits ATS 201	
M005ATS20120	British Type G	103
M005ATS20130	EU Hybrid Type E/F	103
M005ATS20140	Australian Type I	103
M005ATS20150	Danish Type K	103
M005ATS20160	Brazil IEC 60906-1	103
M005ATS20180	Swiss Type J	103
	Alair™ Bronchial Thermoplasty (BT) System	
M005ATS25020	Alair BT Catheter – Global	103
	Habib™ EndoHPB Bipolar Radiofrequency Catheter	
M00500070	Habib Endo HPB Bipolar Radiofrequency Catheter	111
5100	Adapter Cable Bipolar, 28 mm pins	111
5420	Adapter Cable Bipolar, Internal 4 mm pin	111
5700	Adapter Cable Bipolar, 22 mm spacing	111
	One-Step Button and Low Profile Replacement Button Accessories with ENFit™	
M00501061	EndoVive Right Angle Feeding Set Continuous Adapter w/ENFit, 18Fr	99
M00501071	EndoVive Right Angle Feeding Set Continuous Adapter w/ENFit, 24Fr	99
M00501081	EndoVive Right Angle Feeding Set Continuous Adapter w/ENFit, 28Fr	99
M00501091	EndoVive Bolus Feeding Set Straight Adapter w/ENFit, 18Fr	99
M00501101	EndoVive Bolus Feeding Set Straight Adapter w/ENFit, 24Fr	99
M00501111	EndoVive Bolus Feeding Set Straight Adapter w/ENFit, 28Fr	99



Product Code Index (CONTINUED)

Order Number	Description	Page Number
One-Step Button and Low Profile Replacement Button Accessories with ENFit™ (CONTINUED)		
M00501121	EndoVive Decompression Tube w/ENFit, 18Fr	99
M00501131	EndoVive Decompression Tube w/ENFit, 18Fr	99
M00501141	EndoVive Decompression Tube w/ENFit, 18Fr	99
M00501151	EndoVive Decompression Tube w/ENFit, 18Fr	99
M00501161	EndoVive Decompression Tube w/ENFit, 18Fr	99
M00501181	EndoVive Decompression Tube w/ENFit, 24Fr	99
M00501191	EndoVive Decompression Tube w/ENFit, 24Fr	99
M00501201	EndoVive Decompression Tube w/ENFit, 24Fr	99
M00501211	EndoVive Decompression Tube w/ENFit, 24Fr	99
PEG Accessories with ENFit™		
M00501261	EndoVive Y-Port w/ENFit 20Fr	97
M00501271	EndoVive Y-Port w/ENFit 24Fr	97
Standard PEG Kits with ENFit™		
M00509081	EndoVive StdPEG Pull w/ENFit 20Fr	97
M00509101	EndoVive StdPEG Pull w/ENFit 24Fr	97
M00509091	EndoVive StdPEG Push w/ENFit 20Fr	97
M00509111	EndoVive StdPEG Push w/ENFit 24Fr	97
PEG Accessories with ENFit™		
M00509241	ENFit Feed Port Connector	97
M00509221	ENFit Transition Connector Enteral Distal End	97
One-Step Button™ Low Profile PEG Kits with ENFit™		
M00509370	Low Profile PEG Kit with ENFit Push, 24Fr	98
M00509380	Low Profile PEG Kit with ENFit Push, 24Fr	98
M00509390	Low Profile PEG Kit with ENFit Push, 24Fr	98
M00509400	Low Profile PEG Kit with ENFit Push, 24Fr	98
M00509760	Low Profile PEG Kit with ENFit Pull, 18Fr	98
M00509770	Low Profile PEG Kit with ENFit Pull, 18Fr	98
M00509780	Low Profile PEG Kit with ENFit Pull, 18Fr	98
M00509790	Low Profile PEG Kit with ENFit Pull, 18Fr	98
M00509650	Low Profile PEG Kit with ENFit Pull, 24Fr	98
M00509660	Low Profile PEG Kit with ENFit Pull, 24Fr	98
M00509670	Low Profile PEG Kit with ENFit Pull, 24Fr	98
M00510230	Low Profile PEG Kit with ENFit Pull, 24Fr	98
M00509550	Low Profile PEG Kit with ENFit Push, 18Fr	98
M00509560	Low Profile PEG Kit with ENFit Push, 18Fr	98
M00509570	Low Profile PEG Kit with ENFit Push, 18Fr	98
Replacement Gastrostomy Feeding Tubes with ENFit™		
M00509881	12Fr (4.00mm) Repl. G-Tube, Straight	100
M00509891	14Fr (4.67mm) Repl. G-Tube, Straight	100
M00509901	16Fr (5.33mm) Repl. G-Tube, Straight	100
M00509911	18Fr (6.00mm) Repl. G-Tube, Straight	100
M00509921	20Fr (6.67mm) Repl. G-Tube, Straight	100
M00509931	22Fr (7.33mm) Repl. G-Tube, Straight	100
M00509941	24Fr (8.00mm) Repl. G-Tube, Straight	100
M00509961	12Fr (4.00mm) Repl. G-Tube, Right Angle	100
M00509971	14Fr (4.67mm) Repl. G-Tube, Right Angle	100
M00509981	16Fr (5.33mm) Repl. G-Tube, Right Angle	100



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Replacement Gastrostomy Feeding Tubes with ENFit™ (CONTINUED)		
M00509991	18Fr (6.00mm) Repl. G-Tube, Right Angle	100
M00510001	20Fr (6.67mm) Repl. G-Tube, Right Angle	100
M00510011	22Fr (7.33mm) Repl. G-Tube, Right Angle	100
M00510021	24Fr (8.00mm) Repl. G-Tube, Right Angle	100
Multibite™ Multiple Sample Biopsy Forceps		
M00510101	Multibite Multiple Sample Biopsy Forceps	96
M00510121	Multibite Multiple Sample Biopsy Forceps	96
Trapezoid™ RX Wireguided Retrieval Baskets		
M00510860	Trapezoid RX Wireguided Retrieval Basket	44
M00510870	Trapezoid RX Wireguided Retrieval Basket	44
M00510880	Trapezoid RX Wireguided Retrieval Basket	44
M00510890	Trapezoid RX Wireguided Retrieval Basket	44
Zero Tip™ Airway Retrieval Baskets		
M00513200	Zero Tip Airway Retrieval Basket	107
M00513210	Zero Tip Airway Retrieval Basket	107
Radial Jaw™ 4 Large Capacity Biopsy Forceps		
M00513301	Radial Jaw 4 Large Capacity	95
M00513302	Radial Jaw 4 Large Capacity	95
M00513303	Radial Jaw 4 Large Capacity	95
M00513311	Radial Jaw 4 Large Capacity w/Needle	95
M00513312	Radial Jaw 4 Large Capacity w/Needle	95
M00513313	Radial Jaw 4 Large Capacity w/Needle	95
M00513321	Radial Jaw 4 Large Capacity	95
M00513322	Radial Jaw 4 Large Capacity	95
M00513323	Radial Jaw 4 Large Capacity	95
M00513331	Radial Jaw 4 Large Capacity w/Needle	95
M00513332	Radial Jaw 4 Large Capacity w/Needle	95
M00513333	Radial Jaw 4 Large Capacity w/Needle	95
Radial Jaw™ 4 Jumbo Biopsy Forceps		
M00513343	Radial Jaw 4 Jumbo	96
M00513353	Radial Jaw 4 Jumbo w/Needle	96
M00513361	Radial Jaw 4 Jumbo	96
M00513362	Radial Jaw 4 Jumbo	96
M00513363	Radial Jaw 4 Jumbo	96
M00513371	Radial Jaw 4 Jumbo w/Needle	96
M00513372	Radial Jaw 4 Jumbo w/Needle	96
M00513373	Radial Jaw 4 Jumbo w/Needle	96
Radial Jaw™ 4 Standard Capacity Biopsy Forceps		
M00513381	Radial Jaw 4 Standard Capacity	95
M00513382	Radial Jaw 4 Standard Capacity	95
M00513383	Radial Jaw 4 Standard Capacity	95
M00513391	Radial Jaw 4 Standard Capacity w/Needle	95
M00513392	Radial Jaw 4 Standard Capacity w/Needle	95
M00513393	Radial Jaw 4 Standard Capacity w/Needle	95
M00513401	Radial Jaw 4 Standard Capacity	95
M00513402	Radial Jaw 4 Standard Capacity	95
M00513403	Radial Jaw 4 Standard Capacity	95



Product Code Index (CONTINUED)

Order Number	Description	Page Number
	Radial Jaw™ 4 Standard Capacity Biopsy Forceps (CONTINUED)	
M00513411	Radial Jaw 4 Standard Capacity w/Needle	95
M00513412	Radial Jaw 4 Standard Capacity w/Needle	95
M00513413	Radial Jaw 4 Standard Capacity w/Needle	95
	Radial Jaw™ 4 Gastropediatric Biopsy Forceps	
M00513441	Radial Jaw 4 Gastropediatric	96
M00513443	Radial Jaw 4 Gastropediatric	96
M00513451	Radial Jaw 4 Gastropediatric w/Needle	96
M00513453	Radial Jaw 4 Gastropediatric w/Needle	96
	Ultraflex™ Uncovered Esophageal NG Stent System – Distal Release	
M00513700	Ultraflex Esophageal NG Stent System	73
M00513710	Ultraflex Esophageal NG Stent System	73
M00513720	Ultraflex Esophageal NG Stent System	73
	Ultraflex™ Covered Esophageal NG Stent System – Distal Release	
M00513730	Ultraflex Esophageal NG Stent System	73
M00513740	Ultraflex Esophageal NG Stent System	73
M00513750	Ultraflex Esophageal NG Stent System	73
	Ultraflex™ Uncovered Esophageal NG Stent System – Proximal Release	
M00513800	Ultraflex Esophageal NG Stent System	73
M00513810	Ultraflex Esophageal NG Stent System	73
M00513820	Ultraflex Esophageal NG Stent System	73
M00513830	Ultraflex Esophageal NG Stent System	73
	Ultraflex™ Covered Esophageal NG Stent System – Proximal Release	
M00513840	Ultraflex Esophageal NG Stent System	73
M00513850	Ultraflex Esophageal NG Stent System	73
M00513860	Ultraflex Esophageal NG Stent System	73
	Ultraflex™ Covered Esophageal NG Stent System – Distal Release	
M00514200	Ultraflex Esophageal NG Stent System	73
M00514210	Ultraflex Esophageal NG Stent System	73
	Ultraflex™ Uncovered Esophageal NG Stent System – Proximal Release	
M00514220	Ultraflex Esophageal NG Stent System	73
M00514230	Ultraflex Esophageal NG Stent System	73
	Ultraflex™ Covered Esophageal NG Stent System – Proximal Release	
M00514240	Ultraflex Esophageal NG Stent System	73
M00514250	Ultraflex Esophageal NG Stent System	73
	Radial Jaw™ 4 Hot Biopsy Forceps	
M00515031	Radial Jaw 4 Hot Biopsy Forceps	92
M00515032	Radial Jaw 4 Hot Biopsy Forceps	92
M00515033	Radial Jaw 4 Hot Biopsy Forceps	92
	Jagwire™ Pulmonary Guidewire	
M00515171	Jagwire Pulmonary Guidewire	108
	Radial Jaw™ 4 Pulmonary Biopsy Forceps	
M00515181	Radial Jaw 4 Pulmonary Standard Capacity	106
M00515182	Radial Jaw 4 Pulmonary Standard Capacity	106
M00515191	Radial Jaw 4 Pulmonary Standard Capacity w/Needle	106
M00515192	Radial Jaw 4 Pulmonary Standard Capacity w/Needle	106
M00515201	Radial Jaw 4 Pulmonary Large Capacity	106
M00515202	Radial Jaw 4 Pulmonary Large Capacity	106



Product Code Index (CONTINUED)

Order Number	Description	Page Number
	CoreDx™ Pulmonary Mini-Forceps	
M00515220	Pulmonary Mini-Forceps	106
	Celebrity™ Cytology Brushes – Bronchoscope	
M00516001	Celebrity Cytology Brush	107
M00516011	Celebrity Cytology Brush	107
M00516071	Celebrity Cytology Brush	107
M00516151	Celebrity Cytology Brush	107
	WallFlex™ Fully Covered RMV Esophageal Stent System	
M00516210	WallFlex Fully Covered RMV Esophageal Stent System	72
M00516220	WallFlex Fully Covered RMV Esophageal Stent System	72
M00516230	WallFlex Fully Covered RMV Esophageal Stent System	72
M00516240	WallFlex Fully Covered RMV Esophageal Stent System	72
M00516250	WallFlex Fully Covered RMV Esophageal Stent System	72
M00516260	WallFlex Fully Covered RMV Esophageal Stent System	72
	WallFlex™ Fully Covered RMV Esophageal Longer Loop Stent System	
M00516270	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	72
M00516280	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	72
M00516290	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	72
M00516300	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	72
M00516310	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	72
M00516320	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	72
	WallFlex™ Partially Covered Esophageal Stent System	
M00516900	WallFlex Partially Covered Esophageal Stent System	72
M00516910	WallFlex Partially Covered Esophageal Stent System	72
M00516920	WallFlex Partially Covered Esophageal Stent System	72
M00516930	WallFlex Partially Covered Esophageal Stent System	72
M00516940	WallFlex Partially Covered Esophageal Stent System	72
M00516950	WallFlex Partially Covered Esophageal Stent System	72
	Interject™ Clear Sclerotherapy Needles	
M00518111	Interject Sclerotherapy Needle	87
	Interject™ Sclerotherapy Needles	
M00518151	Interject Sclerotherapy Needle	87
M00518161	Interject Sclerotherapy Needle	87
M00518251	Interject Sclerotherapy Needle	87
M00518261	Interject Sclerotherapy Needle	87
	Interject™ Clear Sclerotherapy Needles	
M00518301	Interject Sclerotherapy Needle	87
M00518311	Interject Sclerotherapy Needle	87
	Interject™ Sclerotherapy Needles	
M00518351	Interject Sclerotherapy Needle	87
M00518361	Interject Sclerotherapy Needle	87
	Orise™ Gel	
M00519201	ORISE Gel Syringe Twin Pack Kit	93
M00519210	ORISE Gel Syringe Twin Pack Kit with Interject Needle 23ga Clear	93
M00519211	ORISE Gel Syringe Twin Pack Kit with Interject Needle 23ga Clear	93
	Resolution 360™ Clips	
M00521230	Resolution 360 Clip	85
M00521231	Resolution 360 Clip	85


[Contents](#)
[Alphabetical Index](#)
[Product Code Index](#)

Product Code Index (CONTINUED)

Order Number	Description	Page Number
Resolution 360™ Clips (CONTINUED)		
M00521232	Resolution 360 Clip	85
M00521233	Resolution 360 Clip	85
M00521240	Resolution 360 Clip	85
M00521241	Resolution 360 Clip	85
M00521242	Resolution 360 Clip	85
M00521243	Resolution 360 Clip	85
Resolution™ Clips		
M00522600	Resolution Clip	85
M00522601	Resolution Clip	85
M00522602	Resolution Clip	85
M00522610	Resolution Clip	85
M00522611	Resolution Clip	85
M00522612	Resolution Clip	85
Ultratome™ Double-Lumen Sphincterotomes		
M00530800	Ultratome Short Nose	50
M00530810	Ultratome Long Nose	50
M00530820	Ultratome Short Nose	50
M00530830	Ultratome Long Nose	50
Contour™ ERCP Cannulas		
M00530850	Contour ERCP Cannula	47
M00530860	Contour ERCP Cannula	47
M00530870	Contour ERCP Cannula	47
M00530880	Contour ERCP Cannula	47
M00530890	Contour ERCP Cannula	47
C-flex™ Pigtail Biliary Stents		
M00532010	C-flex Pigtail Stent	65
M00532020	C-flex Pigtail Stent	65
M00532030	C-flex Pigtail Stent	65
M00532040	C-flex Pigtail Stent	65
M00532050	C-flex Pigtail Stent	65
M00532060	C-flex Pigtail Stent	65
Advanix™ Stents – Singles		
M00532160	Advanix Stent	62
M00532180	Advanix Stent	62
M00532190	Advanix Stent	62
M00532210	Advanix Stent	62
M00532220	Advanix Stent	62
M00532230	Advanix Stent	62
M00532240	Advanix Stent	62
M00532260	Advanix Stent	62
M00532270	Advanix Stent	62
M00532290	Advanix Stent	62
M00532300	Advanix Stent	62
M00532310	Advanix Stent	62
Microknife™ XL Triple-Lumen Needle-Knife		
M00532810	Microknife XL Triple-Lumen Needle-Knife	50



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Advanix™ Stents – Singles		
M00532870	Advanix Stent	61
M00532880	Advanix Stent	61
M00532890	Advanix Stent	61
M00532900	Advanix Stent	61
M00532910	Advanix Stent	61
M00532920	Advanix Stent	61
M00532930	Advanix Stent	61
M00532940	Advanix Stent	61
M00532950	Advanix Stent	61
M00532960	Advanix Stent	61
M00532970	Advanix Stent	61
M00532980	Advanix Stent	61
M00532990	Advanix Stent	61
M00533000	Advanix Stent	61
M00533010	Advanix Stent	61
M00533020	Advanix Stent	61
M00533030	Advanix Stent	61
M00533040	Advanix Stent	61
M00533230	Advanix Stent	61
M00533240	Advanix Stent	61
M00533250	Advanix Stent	61
M00533260	Advanix Stent	61
M00533270	Advanix Stent	61
M00533280	Advanix Stent	61
M00533290	Advanix Stent	61
M00533300	Advanix Stent	61
M00533310	Advanix Stent	61
M00533320	Advanix Stent	61
M00533330	Advanix Stent	61
M00533340	Advanix Stent	61
M00533350	Advanix Stent	61
M00533360	Advanix Stent	61
M00533370	Advanix Stent	61
M00533380	Advanix Stent	61
M00533390	Advanix Stent	61
M00533400	Advanix Stent	61
NaviFlex™ RX Delivery System		
M00533560	RX Delivery System	62
M00533570	RX Delivery System	62
M00533580	RX Delivery System	62
Percuflex™ Duodenal Bend Biliary Stents		
M00533610	Percuflex Duodenal Blend Stent	65
M00533620	Percuflex Duodenal Blend Stent	65
M00533630	Percuflex Duodenal Blend Stent	65
M00533640	Percuflex Duodenal Blend Stent	65
M00533650	Percuflex Duodenal Blend Stent	65
M00533660	Percuflex Duodenal Blend Stent	65



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Percuflex™ Duodenal Bend Biliary Stents (CONTINUED)		
M00533670	Percuflex Duodenal Bend Stent	65
M00533680	Percuflex Duodenal Bend Stent	65
M00533690	Percuflex Duodenal Bend Stent	65
M00533700	Percuflex Duodenal Bend Stent	65
Biliary Introducer		
M00533910	Push catheter for use with 7Fr (2.33mm) O.D. stent	65
M00533920	Push and guide catheter for use with 10Fr (3.33mm) O.D. stent	65
Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System		
M00534200	Advanix Pre-loaded RX Stent	60
M00534210	Advanix Pre-loaded RX Stent	60
M00534220	Advanix Pre-loaded RX Stent	60
M00534230	Advanix Pre-loaded RX Stent	60
M00534240	Advanix Pre-loaded RX Stent	60
M00534250	Advanix Pre-loaded RX Stent	60
M00534260	Advanix Pre-loaded RX Stent	60
M00534270	Advanix Pre-loaded RX Stent	60
M00534280	Advanix Pre-loaded RX Stent	60
M00534290	Advanix Pre-loaded RX Stent	60
M00534300	Advanix Pre-loaded RX Stent	60
M00534310	Advanix Pre-loaded RX Stent	60
M00534320	Advanix Pre-loaded RX Stent	60
M00534330	Advanix Pre-loaded RX Stent	60
M00534340	Advanix Pre-loaded RX Stent	60
M00534350	Advanix Pre-loaded RX Stent	60
M00534360	Advanix Pre-loaded RX Stent	60
M00534370	Advanix Pre-loaded RX Stent	60
M00534560	Advanix Pre-loaded RX Stent	60
M00534570	Advanix Pre-loaded RX Stent	60
M00534580	Advanix Pre-loaded RX Stent	60
M00534590	Advanix Pre-loaded RX Stent	60
M00534600	Advanix Pre-loaded RX Stent	60
M00534610	Advanix Pre-loaded RX Stent	60
M00534620	Advanix Pre-loaded RX Stent	60
M00534630	Advanix Pre-loaded RX Stent	60
M00534640	Advanix Pre-loaded RX Stent	60
M00534650	Advanix Pre-loaded RX Stent	60
M00534660	Advanix Pre-loaded RX Stent	60
M00534670	Advanix Pre-loaded RX Stent	60
M00534680	Advanix Pre-loaded RX Stent	60
M00534690	Advanix Pre-loaded RX Stent	60
M00534700	Advanix Pre-loaded RX Stent	60
M00534710	Advanix Pre-loaded RX Stent	60
M00534720	Advanix Pre-loaded RX Stent	60
M00534730	Advanix Pre-loaded RX Stent	60
NaviFlex™ RX Pancreatic Delivery System & NaviFlex™ Pushers		
M00535000	Pancreatic Stent RX Pusher	71
M00535010	Pancreatic Stent RX Pusher	71



Product Code Index (CONTINUED)

Order Number	Description	Page Number
	NaviFlex™ RX Pancreatic Delivery System & NaviFlex™ Pushers CONTINUED	
M00535020	Pancreatic Stent Pusher	71
M00535030	Pancreatic Stent Pusher	71
M00535060	NaviFlex RX Pancreatic Stent Delivery System	71
M00535080	NaviFlex RX Pancreatic Stent Delivery System	71
	Stonetome™ Stone Removal Devices	
M00535110	Stonetome (above cut wire)	47
M00535130	Stonetome (above cut wire)	47
M00535150	Stonetome (below cut wire)	47
M00535170	Stonetome (below cut wire)	47
M00535190	Stonetome (below cut wire)	47
M00535210	Stonetome (below cut wire)	47
	Tandem™ XL Triple-Lumen ERCP Cannula	
M00535700	Tandem XL Triple-Lumen ERCP Cannula	47
	Ultratome XL™ Triple-Lumen Sphincterotomes	
M00535900	Ultratome XL Short Nose	50
M00535910	Ultratome XL Long Nose	50
M00535920	Ultratome XL Short Nose	50
M00535930	Ultratome XL Long Nose	50
	Advanix™ Pancreatic Single Stents – Straight	
M00536000	Advanix Pancreatic Stent	69
M00536010	Advanix Pancreatic Stent	69
M00536020	Advanix Pancreatic Stent	69
M00536030	Advanix Pancreatic Stent	69
M00536040	Advanix Pancreatic Stent	69
M00536050	Advanix Pancreatic Stent	69
M00536060	Advanix Pancreatic Stent	69
M00536070	Advanix Pancreatic Stent	69
M00536080	Advanix Pancreatic Stent	69
M00536090	Advanix Pancreatic Stent	69
M00536100	Advanix Pancreatic Stent	69
M00536110	Advanix Pancreatic Stent	69
M00536120	Advanix Pancreatic Stent	69
M00536130	Advanix Pancreatic Stent	69
M00536140	Advanix Pancreatic Stent	69
M00536150	Advanix Pancreatic Stent	69
M00536160	Advanix Pancreatic Stent	69
M00536170	Advanix Pancreatic Stent	69
M00536180	Advanix Pancreatic Stent	69
M00536190	Advanix Pancreatic Stent	69
M00536200	Advanix Pancreatic Stent	69
M00536210	Advanix Pancreatic Stent	69
M00536220	Advanix Pancreatic Stent	69
M00536230	Advanix Pancreatic Stent	69
M00536240	Advanix Pancreatic Stent	69
M00536250	Advanix Pancreatic Stent	69
M00536260	Advanix Pancreatic Stent	69
M00536270	Advanix Pancreatic Stent	69



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Advanix™ Pancreatic Single Stents – Straight (CONTINUED)		
M00536280	Advanix Pancreatic Stent	69
M00536290	Advanix Pancreatic Stent	69
M00536300	Advanix Pancreatic Stent	69
M00536310	Advanix Pancreatic Stent	69
M00536320	Advanix Pancreatic Stent	69
M00536330	Advanix Pancreatic Stent	69
M00536340	Advanix Pancreatic Stent	69
M00536350	Advanix Pancreatic Stent	69
M00536360	Advanix Pancreatic Stent	70
M00536370	Advanix Pancreatic Stent	70
M00536380	Advanix Pancreatic Stent	70
M00536390	Advanix Pancreatic Stent	70
M00536400	Advanix Pancreatic Stent	70
M00536410	Advanix Pancreatic Stent	70
M00536420	Advanix Pancreatic Stent	70
M00536430	Advanix Pancreatic Stent	70
M00536440	Advanix Pancreatic Stent	70
M00536450	Advanix Pancreatic Stent	70
M00536460	Advanix Pancreatic Stent	70
M00536470	Advanix Pancreatic Stent	70
M00536480	Advanix Pancreatic Stent	70
M00536490	Advanix Pancreatic Stent	70
M00536500	Advanix Pancreatic Stent	70
M00536510	Advanix Pancreatic Stent	70
M00536520	Advanix Pancreatic Stent	70
M00536530	Advanix Pancreatic Stent	70
M00536540	Advanix Pancreatic Stent	70
M00536550	Advanix Pancreatic Stent	70
M00536560	Advanix Pancreatic Stent	70
M00536570	Advanix Pancreatic Stent	70
M00536580	Advanix Pancreatic Stent	70
M00536590	Advanix Pancreatic Stent	70
M00536600	Advanix Pancreatic Stent	70
M00536610	Advanix Pancreatic Stent	70
M00536620	Advanix Pancreatic Stent	70
M00536630	Advanix Pancreatic Stent	70
M00536640	Advanix Pancreatic Stent	70
M00536650	Advanix Pancreatic Stent	70
M00536660	Advanix Pancreatic Stent	70
M00536670	Advanix Pancreatic Stent	70
M00536680	Advanix Pancreatic Stent	70
M00536690	Advanix Pancreatic Stent	70
M00536700	Advanix Pancreatic Stent	70
M00536710	Advanix Pancreatic Stent	70
M00536720	Advanix Pancreatic Stent	70
M00536730	Advanix Pancreatic Stent	70



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Advanix™ Pancreatic Single Stents – Pigtail		
M00536740	Advanix Pancreatic Stent	67
M00536750	Advanix Pancreatic Stent	67
M00536760	Advanix Pancreatic Stent	67
M00536770	Advanix Pancreatic Stent	67
M00536780	Advanix Pancreatic Stent	67
M00536790	Advanix Pancreatic Stent	67
M00536800	Advanix Pancreatic Stent	67
M00536810	Advanix Pancreatic Stent	67
M00536820	Advanix Pancreatic Stent	67
M00536830	Advanix Pancreatic Stent	67
M00536840	Advanix Pancreatic Stent	67
M00536850	Advanix Pancreatic Stent	67
M00536860	Advanix Pancreatic Stent	67
M00536870	Advanix Pancreatic Stent	67
M00536880	Advanix Pancreatic Stent	67
M00536890	Advanix Pancreatic Stent	67
M00536900	Advanix Pancreatic Stent	67
M00536910	Advanix Pancreatic Stent	67
M00536920	Advanix Pancreatic Stent	67
M00536930	Advanix Pancreatic Stent	67
M00536940	Advanix Pancreatic Stent	67
M00536950	Advanix Pancreatic Stent	67
M00536960	Advanix Pancreatic Stent	67
M00536970	Advanix Pancreatic Stent	67
M00536980	Advanix Pancreatic Stent	67
M00536990	Advanix Pancreatic Stent	67
M00537000	Advanix Pancreatic Stent	67
M00537010	Advanix Pancreatic Stent	67
M00537020	Advanix Pancreatic Stent	67
M00537030	Advanix Pancreatic Stent	67
M00537040	Advanix Pancreatic Stent	67
M00537050	Advanix Pancreatic Stent	67
M00537060	Advanix Pancreatic Stent	67
M00537070	Advanix Pancreatic Stent	67
M00537080	Advanix Pancreatic Stent	67
M00537090	Advanix Pancreatic Stent	67
M00537100	Advanix Pancreatic Stent	68
M00537110	Advanix Pancreatic Stent	68
M00537120	Advanix Pancreatic Stent	68
M00537130	Advanix Pancreatic Stent	68
M00537140	Advanix Pancreatic Stent	68
M00537150	Advanix Pancreatic Stent	68
M00537160	Advanix Pancreatic Stent	68
M00537170	Advanix Pancreatic Stent	68
M00537180	Advanix Pancreatic Stent	68
M00537190	Advanix Pancreatic Stent	68
M00537200	Advanix Pancreatic Stent	68


[Contents](#)

Product Code Index (CONTINUED)

Order Number	Description	Page Number
Advanix™ Pancreatic Single Stents – Pigtail (CONTINUED)		
M00537210	Advanix Pancreatic Stent	68
M00537220	Advanix Pancreatic Stent	68
M00537230	Advanix Pancreatic Stent	68
M00537240	Advanix Pancreatic Stent	68
M00537250	Advanix Pancreatic Stent	68
M00537260	Advanix Pancreatic Stent	68
M00537270	Advanix Pancreatic Stent	68
M00537280	Advanix Pancreatic Stent	68
M00537290	Advanix Pancreatic Stent	68
M00537300	Advanix Pancreatic Stent	68
M00537310	Advanix Pancreatic Stent	68
M00537320	Advanix Pancreatic Stent	68
M00537330	Advanix Pancreatic Stent	68
M00537340	Advanix Pancreatic Stent	68
Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Pigtail		
M00537350	Advanix Pancreatic Stent Kit	66
M00537360	Advanix Pancreatic Stent Kit	66
M00537370	Advanix Pancreatic Stent Kit	66
M00537380	Advanix Pancreatic Stent Kit	66
M00537390	Advanix Pancreatic Stent Kit	66
M00537400	Advanix Pancreatic Stent Kit	66
M00537410	Advanix Pancreatic Stent Kit	66
M00537420	Advanix Pancreatic Stent Kit	66
M00537430	Advanix Pancreatic Stent Kit	66
M00537440	Advanix Pancreatic Stent Kit	66
M00537450	Advanix Pancreatic Stent Kit	66
Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Straight		
M00537460	Advanix Pancreatic Stent Kit	66
M00537470	Advanix Pancreatic Stent Kit	66
M00537480	Advanix Pancreatic Stent Kit	66
M00537490	Advanix Pancreatic Stent Kit	66
M00537500	Advanix Pancreatic Stent Kit	66
M00537510	Advanix Pancreatic Stent Kit	66
M00537520	Advanix Pancreatic Stent Kit	66
M00537530	Advanix Pancreatic Stent Kit	66
M00537540	Advanix Pancreatic Stent Kit	66
M00537550	Advanix Pancreatic Stent Kit	66
M00537560	Advanix Pancreatic Stent Kit	66
Flexima™ Biliary Stents with Delivery System		
M00539200	Flexima Stent with Delivery System	64
M00539210	Flexima Stent with Delivery System	64
M00539220	Flexima Stent with Delivery System	64
M00539230	Flexima Stent with Delivery System	64
M00539240	Flexima Stent with Delivery System	64
M00539250	Flexima Stent with Delivery System	64
M00539260	Flexima Stent with Delivery System	64
M00539270	Flexima Stent with Delivery System	64



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Flexima™ Biliary Stents with Delivery System (CONTINUED)		
M00539280	Flexima Stent with Delivery System	64
M00539290	Flexima Stent with Delivery System	64
M00539300	Flexima Stent with Delivery System	64
M00539310	Flexima Stent with Delivery System	64
M00539320	Flexima Stent with Delivery System	64
M00539330	Flexima Stent with Delivery System	64
M00539340	Flexima Stent with Delivery System	64
M00539350	Flexima Stent with Delivery System	64
M00539360	Flexima Stent with Delivery System	64
M00539370	Flexima Stent with Delivery System	64
M00539380	Flexima Stent with Delivery System	64
M00539390	Flexima Stent with Delivery System	64
Flexima™ Nasobiliary Catheter		
M00540120	Flexima Nasal Catheter	64
M00540130	Flexima Nasal Catheter	64
M00540140	Flexima Nasal Catheter	64
M00540160	Flexima Nasal Catheter	64
M00540170	Flexima Nasal Catheter	64
M00540180	Flexima Nasal Catheter	64
Speedband Superview Super 7™ Multiple Band Ligators		
M00542250	Super 7 Band Ligator	85
M00542251	Super 7 Band Ligator	85
M00542253	Super 7 Band Ligator	85
RX Cytology Brush		
M00545000	RX Cytology Brush	45
Autotome™ RX Cannulating Sphincterotomes		
M00545150	Autotome RX 49	42
M00545160	Autotome RX 49	42
M00545170	Autotome RX 44	42
M00545180	Autotome RX 44	42
M00545190	Autotome RX 39	42
M00545200	Autotome RX 39	42
RX Accessories		
M00545211	RX Locking Device – Olympus / FujiFilm Brand Endoscope Compatible	46
M00545221	RX Locking Device – Pentax Brand Endoscope Compatible	46
M00545261	RX Locking Device and Biopsy Cap System – Olympus / FujiFilm G5 Series Brand Endoscope Compatible	46
M00545271	RX Locking Device and Biopsy Cap System – Pentax Brand Endoscope Compatible	46
M00545281	RX Locking Device and Biopsy Cap System – FujiFilm Type 3 Series Brand Endoscope Compatible	46
RX Biliary Plastic Stents with RX Delivery System		
M00545550	RX Biliary Plastic Stent with RX Delivery System	63
M00545560	RX Biliary Plastic Stent with RX Delivery System	63
M00545570	RX Biliary Plastic Stent with RX Delivery System	63
M00545580	RX Biliary Plastic Stent with RX Delivery System	63
M00545590	RX Biliary Plastic Stent with RX Delivery System	63
M00545600	RX Biliary Plastic Stent with RX Delivery System	63
M00545610	RX Biliary Plastic Stent with RX Delivery System	63
M00545620	RX Biliary Plastic Stent with RX Delivery System	63



Product Code Index (CONTINUED)

Order Number	Description	Page Number
RX Biliary Plastic Stents with RX Delivery System (CONTINUED)		
M00545630	RX Biliary Plastic Stent with RX Delivery System	63
M00545640	RX Biliary Plastic Stent with RX Delivery System	63
M00545650	RX Biliary Plastic Stent with RX Delivery System	63
M00545660	RX Biliary Plastic Stent with RX Delivery System	63
M00545670	RX Biliary Plastic Stent with RX Delivery System	63
M00545680	RX Biliary Plastic Stent with RX Delivery System	63
M00545690	RX Biliary Plastic Stent with RX Delivery System	63
M00545700	RX Biliary Plastic Stent with RX Delivery System	63
M00545710	RX Biliary Plastic Stent with RX Delivery System	63
M00545720	RX Biliary Plastic Stent with RX Delivery System	63
M00545730	RX Biliary Plastic Stent with RX Delivery System	63
M00545740	RX Biliary Plastic Stent with RX Delivery System	63
RX ERCP Cannulas		
M00545760	RX ERCP Cannula Standard Tip	39
M00545770	RX ERCP Cannula Tapered Tip	39
M00545780	RX ERCP Cannula Ball Tip	39
Rapid Exchange XL Cannula		
M00545820	Rapid Exchange XL Cannula	39
RX Needle-Knife XL Triple-Lumen Needle-Knife		
M00545840	RX Needle-Knife XL Triple-Lumen Needle-Knife	42
Hurricane™ RX Biliary Dilation Balloons		
M00545890	Hurricane RX Dilation Balloon	45
M00545900	Hurricane RX Dilation Balloon	45
M00545910	Hurricane RX Dilation Balloon	45
M00545920	Hurricane RX Dilation Balloon	45
M00545930	Hurricane RX Dilation Balloon	45
M00545940	Hurricane RX Dilation Balloon	45
M00545950	Hurricane RX Dilation Balloon	45
M00545960	Hurricane RX Dilation Balloon	45
SpyGlass™ DS II		
M00546470	SpyBite™ Max Biopsy Forceps	53
M00546550	SpyGlass™ Retrieval Basket	53
M00546560	SpyGlass™ Retrieval Snare	53
M00546610	SpyScope™ DS II Access and Delivery Catheter	53
M00546650	SpyGlass DS II Digital Controller	53
Autolith™ Touch Biliary EHL System		
M00546620	1.9Fr, 375cm Biliary EHL Probe	53
M005466800	Autolith Touch EHL Generator	53
M00546750	Autolith Touch Extender Cable	53
M00546760	Autolith Touch Foot Pedal	53
Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Distal Injection		
M00547000	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	43
M00547010	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	43
M00547020	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	43
Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Proximal Injection		
M00547030	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	43
M00547040	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	43



Product Code Index (CONTINUED)

Order Number	Description	Page Number
	Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Proximal Injection (CONTINUED)	
M00547050	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	43
	Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Distal Injection	
M00547100	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	51
M00547110	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	51
M00547120	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	51
	Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	
M00547130	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	51
M00547140	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	51
M00547150	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	51
	Extractor™ Pro DL Double-Lumen Retrieval Balloon Catheters	
M00547200	Extractor Pro DL Double-Lumen Retrieval Balloon Catheter	51
M00547210	Extractor Pro DL Double-Lumen Retrieval Balloon Catheter	51
M00547220	Extractor Pro DL Double-Lumen Retrieval Balloon Catheter	51
	Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Distal Injection	
M00547300	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	44
M00547310	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	44
M00547320	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	44
	Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	
M00547330	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	44
M00547340	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	44
M00547350	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	44
	Expect™ Standard Handle Needles	
M00550000	Expect Standard Handle Needle	55
M00550001	Expect Standard Handle Needle	55
M00550010	Expect Standard Handle Needle	55
M00550011	Expect Standard Handle Needle	55
M00550020	Expect Standard Handle Needle	55
M00550021	Expect Standard Handle Needle	55
M00550040	Expect Standard Handle Needle	55
M00550041	Expect Standard Handle Needle	55
	Amplatz Super Stiff™ Guidewire	
M00550090	Amplatz Super Stiff Guidewire	83
	CRE™ Pulmonary Balloon Dilatation Catheter	
M00550300	CRE Pulmonary Balloon Dilatation Catheter	107
M00550310	CRE Pulmonary Balloon Dilatation Catheter	107
M00550320	CRE Pulmonary Balloon Dilatation Catheter	107
M00550330	CRE Pulmonary Balloon Dilatation Catheter	107
M00550340	CRE Pulmonary Balloon Dilatation Catheter	107
M00550350	CRE Pulmonary Balloon Dilatation Catheter	107
	Alliance™ II Inflation System	
M00550601	Alliance II Single-Use Syringe/Gauge Assembly	80, 108
M00550620	Alliance II Integrated Inflation/Lithotripsy Device	44, 80, 108
	Inflation/Deflation Devices and Accessories	
M00550630	CRE™ SteriFlate Disposable Inflation Device	80
	NovaGold™ High Performance Guidewires	
M00552001	NovaGold High Performance Guidewires	83
M00552011	NovaGold High Performance Guidewires	83


[Contents](#)

Product Code Index (CONTINUED)

Order Number	Description	Page Number
Acquire™ Pulmonary Endobronchial Ultrasound Fine Biopsy Needles – Olympus Scope Compatible		
M00552351	22ga Acquire Pulmonary – Olympus scope compatible	104
M00552401	25ga Acquire Pulmonary – Olympus scope compatible	104
Inflation/Deflation Devices and Accessories		
M00553200	Pneumatic Inflator	80
AXIOS™ Stent and Delivery System		
M00553440	AXIOS Stent and Delivery System	56
M00553450	AXIOS Stent and Delivery System	56
Hot AXIOS™ Stent and Electrocautery Enhanced Delivery System		
M00553520	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	56
M00553530	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	56
M00553540	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	56
M00553550	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	56
M00553560	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	56
Rigiflex™ II Achalasia Balloon Dilators		
M00554500	Rigiflex II Achalasia Balloon Dilator	79
M00554510	Rigiflex II Achalasia Balloon Dilator	79
M00554520	Rigiflex II Achalasia Balloon Dilator	79
Expect™ Slimline Handle Needles		
M00555500	Expect Slimline Handle Needle	55
M00555501	Expect Slimline Handle Needle	55
M00555510	Expect Slimline Handle Needle	55
M00555511	Expect Slimline Handle Needle	55
M00555520	Expect Slimline Handle Needle	55
M00555521	Expect Slimline Handle Needle	55
M00555530	Expect Slimline Handle Needle	55
M00555531	Expect Slimline Handle Needle	55
Acquire™ EUS-FNB Devices		
M00555540	22ga Acquire EUS-FNB Needle	56
M00555541	22ga Acquire EUS-FNB Needle	56
M00555560	25ga Acquire EUS-FNB Needle	56
M00555561	25ga Acquire EUS-FNB Needle	56
M00555580	19ga Acquire FNB Needle, Flexible	56
WallFlex™ Colonic Soft Stent System with Anchor Lock Delivery System		
M00555700	WallFlex Colonic Soft Stent	74
M00555710	WallFlex Colonic Soft Stent	74
M00555720	WallFlex Colonic Soft Stent	74
M00555730	WallFlex Colonic Soft Stent	74
M00555740	WallFlex Colonic Soft Stent	74
M00555750	WallFlex Colonic Soft Stent	74
M00555760	WallFlex Colonic Soft Stent	74
M00555770	WallFlex Colonic Soft Stent	74
M00555780	WallFlex Colonic Soft Stent	74
WallFlex™ Duodenal Soft Stent System with Anchor Lock Delivery System		
M00555800	WallFlex Duodenal Soft Stent	75
M00555810	WallFlex Duodenal Soft Stent	75
M00555820	WallFlex Duodenal Soft Stent	75
M00555830	WallFlex Duodenal Soft Stent	75



Product Code Index (CONTINUED)

Order Number	Description	Page Number
WallFlex™ Duodenal Soft Stent System with Anchor Lock Delivery System (CONTINUED)		
M00555840	WallFlex Duodenal Soft Stent	75
M00555850	WallFlex Duodenal Soft Stent	75
M00555860	WallFlex Duodenal Soft Stent	75
M00555870	WallFlex Duodenal Soft Stent	75
M00555880	WallFlex Duodenal Soft Stent	75
Hydra Jagwire™ High Performance Guidewires		
M00556001	Hydra Jagwire Standard	81
M00556011	Hydra Jagwire Standard	81
M00556021	Hydra Jagwire Stiff Shaft	81
M00556031	Hydra Jagwire Stiff Shaft	81
M00556041	Hydra Jagwire Standard	81
M00556051	Hydra Jagwire Standard	81
M00556061	Hydra Jagwire Stiff Shaft	81
M00556071	Hydra Jagwire Stiff Shaft	81
Dreamwire™ High Performance Guidewires		
M00556101	Dreamwire Standard	81
M00556111	Dreamwire Standard	81
M00556121	Dreamwire Stiff Shaft	81
M00556131	Dreamwire Stiff Shaft	81
M00556141	Dreamwire Standard	81
M00556151	Dreamwire Standard	81
M00556161	Dreamwire Stiff Shaft	81
M00556171	Dreamwire Stiff Shaft	81
NaviPro™ Hydrophilic Guidewires		
M00556181	NaviPro Standard	83
M00556191	NaviPro Standard	83
M00556201	NaviPro Standard	83
M00556211	NaviPro Standard	83
M00556221	NaviPro Standard	83
M00556231	NaviPro Standard	83
M00556241	NaviPro Stiff Shaft	83
M00556251	NaviPro Stiff Shaft	83
Jagwire™ High Performance Guidewires		
M00556461	Jagwire Standard	82
M00556481	Jagwire Stiff Shaft	82
M00556491	Jagwire Stiff Shaft	82
M00556511	Jagwire Standard	82
M00556521	Jagwire Stiff Shaft	82
M00556531	Jagwire Stiff Shaft	82
M00556561	Jagwire Standard	82
M00556581	Jagwire Standard	82
M00556591	Jagwire Standard	82
M00556601	Jagwire Stiff Shaft	82
M00556611	Jagwire Stiff Shaft	82
M00556621	Jagwire Stiff Shaft	82
Jagwire™ High Performance Guidewires – Extendable Jagwire		
M00556841	Jagwire Standard	82



Contents



Alphabetical Index

Product Code Index

Product Code Index (CONTINUED)

Order Number	Description	Page Number
Jagwire™ High Performance Guidewires – Extendable Jagwire (CONTINUED)		
M00556851	Jagwire Standard	82
M00556861	Jagwire Stiff Shaft	82
M00556871	Jagwire Stiff Shaft	82
M00556902	Jagwire Jag Tail	82
Jagwire™ Revolution High Performance Guidewire		
M00557001	Jagwire Revolution High Performance Guidewire	82
M00557011	Jagwire Revolution High Performance Guidewire	82
M00557021	Jagwire Revolution High Performance Guidewire	82
M00557031	Jagwire Revolution High Performance Guidewire	82
M00557041	Jagwire Revolution High Performance Guidewire	82
Ultraflex™ Precision Colonic Stent System – Proximal Release		
M00557330	Ultraflex Precision Colonic Stent	76
M00557360	Ultraflex Precision Colonic Stent	76
M00557380	Ultraflex Precision Colonic Stent	76
Jagwire™ High Performance Guidewires		
M00557641	Jagwire Standard	82
M00557651	Jagwire Standard	82
M00557661	Jagwire Stiff Shaft	82
M00557671	Jagwire Stiff Shaft	82
Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needles – Olympus Scope Compatible		
M00558220	Expect Pulmonary 22ga – Olympus Scope Compatible	103
M00558221	Expect Pulmonary 22ga – Olympus Scope Compatible	103
M00558250	Expect Pulmonary 25ga – Olympus Scope Compatible	103
M00558251	Expect Pulmonary 25ga – Olympus Scope Compatible	103
M00558731	Expect Pulmonary Needle Adaptor – Olympus Scope Compatible	103
CRE™ Fixed Wire Balloon Dilators		
M00558330	CRE Fixed Wire Balloon Dilator	78
M00558340	CRE Fixed Wire Balloon Dilator	78
M00558350	CRE Fixed Wire Balloon Dilator	78
M00558360	CRE Fixed Wire Balloon Dilator	78
M00558361	CRE Fixed Wire Balloon Dilator	78
M00558370	CRE Fixed Wire Balloon Dilator	78
M00558371	CRE Fixed Wire Balloon Dilator	78
M00558380	CRE Fixed Wire Balloon Dilator	78
M00558381	CRE Fixed Wire Balloon Dilator	78
CRE™ PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilatation Catheters (including DASE indication)		
M00558600	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	77
M00558610	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	77
M00558620	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	77
M00558630	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	77
M00558640	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	77
M00558650	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	77
CRE™ PRO Wireguided Esophageal/Pyloric/Colonic/Biliary Balloon Dilatation Catheters (including DASE indication)		
M00558660	CRE PRO Wireguided Colonic/Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	77
M00558670	CRE PRO Wireguided Colonic/Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	77
M00558680	CRE PRO Wireguided Colonic/Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	77
M00558690	CRE PRO Wireguided Colonic/Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	77



Product Code Index (CONTINUED)

Order Number	Description	Page Number
CRE™ PRO Wireguided Esophageal/Pyloric/Colonic/Biliary Balloon Dilatation Catheters (including DASE indication) (CONTINUED)		
M00558700	CRE PRO Wireguided Colonic/Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	77
M00558710	CRE PRO Wireguided Colonic/Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	77
CRE™ RX Biliary Balloon Dilatation Catheter		
M00558920	CRE RX Biliary Balloon Dilatation Catheter	45
M00558930	CRE RX Biliary Balloon Dilatation Catheter	45
M00558940	CRE RX Biliary Balloon Dilatation Catheter	45
M00558860	CRE RX Biliary Balloon Dilatation Catheter	45
M00558870	CRE RX Biliary Balloon Dilatation Catheter	45
M00558880	CRE RX Biliary Balloon Dilatation Catheter	45
M00558890	CRE RX Biliary Balloon Dilatation Catheter	45
M00558900	CRE RX Biliary Balloon Dilatation Catheter	45
M00558910	CRE RX Biliary Balloon Dilatation Catheter	45
Gold Probe™ Electrohemostasis Catheters		
M00560070	Gold Probe Catheter	86
M00560071	Gold Probe Catheter	86
M00560100	Gold Probe Catheter	86
M00560101	Gold Probe Catheter	86
Injection Gold Probe™ Electrohemostasis/Injection Catheters		
M00560150	Injection Gold Probe Catheter	86
M00560160	Injection Gold Probe Catheter	86
Gold Probe™ Electrohemostasis Catheters		
M00560220	Gold Probe Catheter	86
Sensation™ Single-Use Short Throw Snares		
M00560311	Medium Oval	90
M00560321	Small Oval	90
M00560331	Medium Crescent	90
Captivator™ COLD Single-Use Snare		
M00561101	10mm Rounded–Stiff	91
M00561102	10mm Rounded–Stiff	91
M00561103	10mm Rounded–Stiff	91
Captivator™ II Single-Use Snares		
M00561191	25 mm Rounded	91
M00561192	25 mm Rounded	91
M00561193	25 mm Rounded	91
M00561221	10 mm Rounded	91
M00561222	10 mm Rounded	91
M00561223	10 mm Rounded	91
M00561231	15 mm Rounded	91
M00561232	15 mm Rounded	91
M00561233	15 mm Rounded	91
M00561241	20 mm Rounded	91
M00561242	20 mm Rounded	91
M00561243	20 mm Rounded	91
Accessories		
M00561250	Active Cord, Fits Olympus Brand Generator	52, 92
M00561270	Active Cord, Fits Endostat™, Endostat™ II, Endostat™ III, Bovie®, Valleylab®, Aspen Labs and ERBE® Brand Generators	52, 92



Product Code Index (CONTINUED)

Order Number	Description	Page Number
	Captivator™ II Single-Use Snares	
M00561291	33 mm Rounded	91
	Captivator™ Single-Use Snares	
M00561311	Large Oval	91
	Electrosurgical Bipolar Cable Adaptors	
M00561440	Bipolar Cable Adaptor	88
M00561450	Bipolar Cable Adaptor	88
M00561460	Bipolar Cable Adaptor	88
M00561470	Bipolar Cable Adaptor	88
M00561480	Bipolar Cable Adaptor	88
M00561490	Bipolar Cable Adaptor	88
M00561500	Bipolar Cable Adaptor	88
	Captivator™ EMR	
M00561600	Captivator EMR for Standard Gastrosopes	89
M00561601	Captivator EMR for Standard Gastrosopes	89
M00561610	Captivator EMR for Large Gastrosopes	89
M00561611	Captivator EMR for Large Gastrosopes	89
	Rotatable Single-Use Snares	
M00561821	Small Oval	89
M00561831	Medium Oval	89
	Captivator™ Single-Use Snares	
M00562301	Small Oval	91
M00562321	Medium Oval	91
M00562341	Medium Hexagonal	91
M00562371	Medium Crescent	91
M00562391	Large Oval (thin wire)	91
	Captiflex™ Single-Use Snares	
M00562401	Medium Oval	90
M00562402	Medium Oval	90
M00562421	Small Oval	90
M00562422	Small Oval	90
	Captivator™ Single-Use Snares	
M00562451	Small Hexagonal	91
	Captiflex™ Single-Use Snares	
M00562471	Extra Small Oval	90
	Profile™ Single-Use Snares	
M00562531	Extra Small Oval	92
M00562551	Small Oval	92
M00562571	Medium Oval	92
	Sensation™ Single-Use Short Throw Snares	
M00562651	Large Oval	90
M00562652	Large Oval	90
M00562671	Medium Oval	90
M00562672	Medium Oval	90
M00562673	Medium Oval	90
M00562691	Small Oval	90
M00562692	Small Oval	90
M00562693	Small Oval	90


[Contents](#)
[◀ 30 ▶](#)
[Alphabetical Index](#)
[Product Code Index](#)

RX Biliary System™

Biliary

SpyGlass™

EUS

Stents

Dilation

Guidewires

Hemostasis

Tissue Resection

Biopsy

Enteral Feeding

Pulmonary

Infection Prevention

RFA

Product Code Index (CONTINUED)

Order Number	Description	Page Number
eXcelon™ Transbronchial Aspiration Needles		
M00564101	eXcelon Transbronchial Aspiration Needle	106
M00564111	eXcelon Transbronchial Aspiration Needle	106
M00564121	eXcelon Transbronchial Aspiration Needle	106
Ultraflex™ Tracheobronchial Uncovered Stent System – Distal Release		
M00564500	Ultraflex Tracheobronchial Stent System	105
M00564510	Ultraflex Tracheobronchial Stent System	105
M00564520	Ultraflex Tracheobronchial Stent System	105
M00564530	Ultraflex Tracheobronchial Stent System	105
M00564540	Ultraflex Tracheobronchial Stent System	105
M00564560	Ultraflex Tracheobronchial Stent System	105
M00564570	Ultraflex Tracheobronchial Stent System	105
M00564590	Ultraflex Tracheobronchial Stent System	105
M00564600	Ultraflex Tracheobronchial Stent System	105
Ultraflex™ Tracheobronchial Uncovered Stent System – Proximal Release		
M00564640	Ultraflex Tracheobronchial Stent System	105
M00564650	Ultraflex Tracheobronchial Stent System	105
M00564660	Ultraflex Tracheobronchial Stent System	105
M00564670	Ultraflex Tracheobronchial Stent System	105
M00564680	Ultraflex Tracheobronchial Stent System	105
M00564690	Ultraflex Tracheobronchial Stent System	105
M00564700	Ultraflex Tracheobronchial Stent System	105
M00564710	Ultraflex Tracheal NG Stent System	105
M00564720	Ultraflex Tracheal NG Stent System	105
Ultraflex™ Tracheobronchial Covered Stent System – Distal Release		
M00564740	Ultraflex Tracheobronchial Stent System	104
M00564750	Ultraflex Tracheobronchial Stent System	104
M00564760	Ultraflex Tracheobronchial Stent System	104
M00564770	Ultraflex Tracheobronchial Stent System	104
M00564780	Ultraflex Tracheobronchial Stent System	104
M00564790	Ultraflex Tracheobronchial Stent System	104
M00564800	Ultraflex Tracheobronchial Stent System	104
M00564810	Ultraflex Tracheobronchial Stent System	104
M00564820	Ultraflex Tracheal NG Stent System	104
M00564830	Ultraflex Tracheal NG Stent System	104
M00564840	Ultraflex Tracheal NG Stent System	104
M00564850	Ultraflex Tracheal NG Stent System	104
M00564860	Ultraflex Tracheal NG Stent System	104
M00564870	Ultraflex Tracheal NG Stent System	104
M00564880	Ultraflex Tracheal NG Stent System	104
M00564890	Ultraflex Tracheal NG Stent System	104
M00564900	Ultraflex Tracheal NG Stent System	104
M00564910	Ultraflex Tracheal NG Stent System	104
WallFlex™ Duodenal Stent System with Anchor Lock Delivery System		
M00565010	WallFlex Duodenal Stent	75
M00565020	WallFlex Duodenal Stent	75
M00565030	WallFlex Duodenal Stent	75



Product Code Index (CONTINUED)

Order Number	Description	Page Number
WallFlex™ Colonic Stent System with Anchor Lock Delivery System		
M00565040	WallFlex Colonic Stent	74
M00565050	WallFlex Colonic Stent	74
M00565060	WallFlex Colonic Stent	74
M00565070	WallFlex Colonic Stent	74
M00565080	WallFlex Colonic Stent	74
M00565090	WallFlex Colonic Stent	74
M00565100	WallFlex Colonic Stent	74
M00565110	WallFlex Colonic Stent	74
M00565120	WallFlex Colonic Stent	74
M00565130	WallFlex Colonic Stent	74
M00565140	WallFlex Colonic Stent	74
M00565150	WallFlex Colonic Stent	74
Wallstent™ Enteral Endoprosthesis Colonic/Duodenal Stent with UNISTEP™ Plus Delivery System		
M00565560	Wallstent Enteral Endoprosthesis Colonic/Duodenal Stent	76
M00565570	Wallstent Enteral Endoprosthesis Colonic/Duodenal Stent	76
M00565580	Wallstent Enteral Endoprosthesis Colonic/Duodenal Stent	76
M00565590	Wallstent Enteral Endoprosthesis Colonic/Duodenal Stent	76
M00565620	Wallstent Enteral Endoprosthesis Colonic/Duodenal Stent	76
M00565630	Wallstent Enteral Endoprosthesis Colonic/Duodenal Stent	76
M00565640	Wallstent Enteral Endoprosthesis Colonic/Duodenal Stent	76
M00565650	Wallstent Enteral Endoprosthesis Colonic/Duodenal Stent	76
Accessories		
M00566001	Rapid Refill Continuous Injection System	46, 52
Through-The-PEG (TTP) Jejunal Feeding Tubes		
M00566320	8.5Fr (2.83mm) TTP Jejunal Feeding Tube	101
M00566330	8.5Fr (2.83mm) TTP Jejunal Feeding Tube	101
M00566340	12Fr (4.00mm) TTP Jejunal Feeding Tube	101
M00566350	12Fr (4.00mm) TTP Jejunal Feeding Tube	101
Accessories		
M00566670	Encore™ 26 Inflation Device	46, 52
MaxForce™ Biliary Balloon Dilatation Catheter		
M00567330	MaxForce Biliary Balloon Dilatation Catheter	52
M00567340	MaxForce Biliary Balloon Dilatation Catheter	52
M00567370	MaxForce Biliary Balloon Dilatation Catheter	52
M00567380	MaxForce Biliary Balloon Dilatation Catheter	52
M00567410	MaxForce Biliary Balloon Dilatation Catheter	52
M00567420	MaxForce Biliary Balloon Dilatation Catheter	52
Low Profile Replacements Button		
M00568170	28Fr (9.33mm) Low Profile Replacement Button	100
M00568180	28Fr (9.33mm) Low Profile Replacement Button	100
M00568190	28Fr (9.33mm) Low Profile Replacement Button	100
M00568260	18Fr (6.00mm) Low Profile Replacement Button	100
M00568270	18Fr (6.00mm) Low Profile Replacement Button	100
M00568280	18Fr (6.00mm) Low Profile Replacement Button	100
M00568290	18Fr (6.00mm) Low Profile Replacement Button	100
M00568300	18Fr (6.00mm) Low Profile Replacement Button	100
M00568320	24Fr (8.00mm) Low Profile Replacement Button	100



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Low Profile Replacements Button (CONTINUED)		
M00568330	24Fr (8.00mm) Low Profile Replacement Button	100
M00568340	24Fr (8.00mm) Low Profile Replacement Button	100
M00568350	24Fr (8.00mm) Low Profile Replacement Button	100
Standard PEG Kits		
M00568381	20Fr Standard PEG Kit, Pull	97
M00568391	24Fr Standard PEG Kit, Pull	97
M00568401	20Fr Standard PEG Kit, Push	97
M00568411	24Fr Standard PEG Kit, Push	97
One-Step Button™ Low Profile PEG Kits		
M00568450	18Fr (6.00mm) Low Profile PEG Kit, Pull	98
M00568460	18Fr (6.00mm) Low Profile PEG Kit, Pull	98
M00568470	18Fr (6.00mm) Low Profile PEG Kit, Pull	98
M00568480	18Fr (6.00mm) Low Profile PEG Kit, Pull	98
M00568500	24Fr (8.00mm) Low Profile PEG Kit, Pull	98
M00568510	24Fr (8.00mm) Low Profile PEG Kit, Pull	98
M00568520	24Fr (8.00mm) Low Profile PEG Kit, Pull	98
M00568530	24Fr (8.00mm) Low Profile PEG Kit, Pull	98
M00568560	18Fr (6.00mm) Low Profile PEG Kit, Push	98
M00568570	18Fr (6.00mm) Low Profile PEG Kit, Push	98
M00568580	18Fr (6.00mm) Low Profile PEG Kit, Push	98
M00568600	24Fr (8.00mm) Low Profile PEG Kit, Push	98
M00568610	24Fr (8.00mm) Low Profile PEG Kit, Push	98
M00568620	24Fr (8.00mm) Low Profile PEG Kit, Push	98
M00568630	24Fr (8.00mm) Low Profile PEG Kit, Push	98
Wallstent™ RX Biliary Endoprosthesis Stent System (Uncovered)		
M00569610	Wallstent RX Biliary	59
M00569620	Wallstent RX Biliary	59
M00569630	Wallstent RX Biliary	59
M00569640	Wallstent RX Biliary	59
M00569650	Wallstent RX Biliary	59
M00569660	Wallstent RX Biliary	59
Wallstent™ RX Biliary Endoprosthesis Stent System (Permalume™ Covering)		
M00569670	Wallstent RX w/Permalume	58
M00569680	Wallstent RX w/Permalume	58
M00569690	Wallstent RX w/Permalume	58
M00569700	Wallstent RX w/Permalume	58
M00569710	Wallstent RX w/Permalume	58
M00569720	Wallstent RX w/Permalume	58
Wallstent™ RX Biliary Endoprosthesis Stent System (Uncovered)		
M00569800	Wallstent RX Biliary	59
M00569810	Wallstent RX Biliary	59
WallFlex™ Biliary RX Fully Covered Stent System RMV (Permalume™ Covering, includes benign indication)		
M00570450	WallFlex RX Biliary w/Permalume	57
M00570460	WallFlex RX Biliary w/Permalume	57
M00570470	WallFlex RX Biliary w/Permalume	57
M00570480	WallFlex RX Biliary w/Permalume	57
M00570490	WallFlex RX Biliary w/Permalume	57



Product Code Index (CONTINUED)

Order Number	Description	Page Number
WallFlex™ Biliary RX Uncovered Stent System		
M00570600	WallFlex RX Biliary	57
M00570610	WallFlex RX Biliary	57
M00570620	WallFlex RX Biliary	57
M00570630	WallFlex RX Biliary	57
M00570640	WallFlex RX Biliary	57
M00570650	WallFlex RX Biliary	57
M00570660	WallFlex RX Biliary	57
WallFlex™ Biliary RX Partially Covered Stent System (Permalume™ Covering)		
M00570700	WallFlex RX Biliary w/Permalume	57
M00570710	WallFlex RX Biliary w/Permalume	57
M00570720	WallFlex RX Biliary w/Permalume	57
M00570730	WallFlex RX Biliary w/Permalume	57
M00570740	WallFlex RX Biliary w/Permalume	57
WallFlex™ Biliary RX Uncovered Stent System		
M00570890	WallFlex RX Biliary	58
Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire		
M00573000	Jagtome RX 49	41
M00573010	Jagtome RX 49	41
M00573020	Jagtome RX 49	41
M00573030	Jagtome RX 49	41
M00573040	Jagtome RX 44	41
M00573050	Jagtome RX 44	41
M00573060	Jagtome RX 44	41
M00573070	Jagtome RX 44	41
M00573080	Jagtome RX 39	41
M00573090	Jagtome RX 39	41
M00573100	Jagtome RX 39	41
M00573110	Jagtome RX 39	41
M00573500	Jagtome RX 49	41
M00573510	Jagtome RX 49	41
M00573520	Jagtome RX 49	41
M00573530	Jagtome RX 49	41
M00573540	Jagtome RX 44	41
M00573550	Jagtome RX 44	41
M00573560	Jagtome RX 44	41
M00573570	Jagtome RX 44	41
WallFlex™ Biliary RX Fully Covered Stent System RMV (Permalume™ Covering, includes benign indication)		
M00574210	WallFlex RX Biliary w/Permalume	57
M00574230	WallFlex RX Biliary w/Permalume	57
M00574240	WallFlex RX Biliary w/Permalume	57
WallFlex™ Biliary RX Partially Covered Stent System (Permalume™ Covering)		
M00576740	WallFlex RX Biliary w/Permalume	57
M00576750	WallFlex RX Biliary w/Permalume	57
M00576760	WallFlex RX Biliary w/Permalume	57
M00576770	WallFlex RX Biliary w/Permalume	57
WallFlex™ Biliary RX Uncovered Stent System		
M00576780	WallFlex RX Biliary	58


[Contents](#)

Product Code Index (CONTINUED)

Order Number	Description	Page Number
	WallFlex™ Biliary RX Uncovered Stent System (CONTINUED)	
M00576790	WallFlex RX Biliary	58
	WallFlex™ Biliary RX Fully Covered System RMV (Permalume™ Covering, includes benign indication)	
M00576800	WallFlex RX Biliary w/Permalume	57
M00576810	WallFlex RX Biliary w/Permalume	57
M00576820	WallFlex RX Biliary w/Permalume	57
M00576830	WallFlex RX Biliary w/Permalume	57
	One-Step Button™ and Low Profile Replacement Button Accessories	
M00580171	18Fr (6.00mm) Bolus Straight Feeding Adapter	99
M00580181	18Fr (6.00mm) Continuous Right Angle Feeding Adapter	99
M00580231	24Fr (8.00mm) Bolus Straight Feeding Adapter	99
M00580241	24Fr (8.00mm) Continuous Right Angle Feeding Adapter	99
M00580271	28Fr (9.33mm) Bolus Straight Feeding Adapter	99
M00580281	28Fr (9.33mm) Continuous Right Angle Feeding Adapter	99
M00580401	18Fr (6.00mm) Decompression Tube	99
M00580411	18Fr (6.00mm) Decompression Tube	99
M00580421	18Fr (6.00mm) Decompression Tube	99
M00580431	18Fr (6.00mm) Decompression Tube	99
M00580441	18Fr (6.00mm) Decompression Tube	99
M00580461	24Fr (8.00mm) Decompression Tube	99
M00580471	24Fr (8.00mm) Decompression Tube	99
M00580481	24Fr (8.00mm) Decompression Tube	99
M00580491	24Fr (8.00mm) Decompression Tube	99
	PEG Accessories	
M00580641	Y-Port 20Fr	97
M00580651	Y-Port 24Fr	97
M00580731	20Fr (6.67mm) Round Bolsters	97
M00580781	20Fr (6.67mm) and 24F (8.00mm) C-Clamps	97
	Replacement Gastrostomy Feeding Tubes	
M00582040	12Fr (4.00mm) Replacement G-Tube, Straight	101
M00582050	14Fr (4.67mm) Replacement G-Tube, Straight	101
M00582060	16Fr (5.33mm) Replacement G-Tube, Straight	101
M00582070	18Fr (6.00mm) Replacement G-Tube, Straight	101
M00582080	20Fr (6.67mm) Replacement G-Tube, Straight	101
M00582090	24Fr (8.00mm) Replacement G-Tube, Straight	101
M00582100	28Fr (9.33mm) Replacement G-Tube, Straight	101
M00582110	12Fr (4.00mm) Replacement G-Tube, Right Angle	101
M00582120	14Fr (4.67mm) Replacement G-Tube, Right Angle	101
M00582130	16Fr (5.33mm) Replacement G-Tube, Right Angle	101
M00582140	18Fr (6.00mm) Replacement G-Tube, Right Angle	101
M00582150	20Fr (6.67mm) Replacement G-Tube, Right Angle	101
M00582160	24Fr (8.00mm) Replacement G-Tube, Right Angle	101
M00582170	28Fr (9.33mm) Replacement G-Tube, Right Angle	101
M00582860	22Fr (7.33mm) Replacement G-Tube, Straight	101
M00582870	22Fr (7.33mm) Replacement G-Tube, Right Angle	101
	Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	
M00583000	Hydratome RX 49	40
M00583010	Hydratome RX 49	40



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire (CONTINUED)		
M00583020	Hydratome RX 49	40
M00583030	Hydratome RX 49	40
M00583040	Hydratome RX 44	40
M00583050	Hydratome RX 44	40
M00583060	Hydratome RX 44	40
M00583070	Hydratome RX 44	40
TRUEtome™ Jag Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire		
M00583080	TRUEtome Jag 39 Sphincterotome	48
M00583090	TRUEtome Jag 44 Sphincterotome	48
M00583100	TRUEtome Jag 44 Sphincterotome	48
M00583110	TRUEtome Jag 44 Sphincterotome	48
M00583120	TRUEtome Jag 39 Sphincterotome	48
M00583130	TRUEtome Jag 44 Sphincterotome	48
M00583140	TRUEtome Jag 44 Sphincterotome	48
M00583150	TRUEtome Jag 44 Sphincterotome	48
TRUEtome™ Hydra Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire		
M00583160	TRUEtome Hydra 44 Sphincterotome	48
M00583170	TRUEtome Hydra 44 Sphincterotome	48
TRUEtome™ Dreamwire™ Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire		
M00583190	TRUEtome Dreamwire 44 Sphincterotome	48
M00583210	TRUEtome Dreamwire 44 Sphincterotome	48
TRUEtome™ Jagwire Revolution Cannulating Sphincterotomes with Pre-loaded Jagwire™ Revolution Guidewire		
M00583260	TRUEtome 39 Sphincterotome	48
M00583270	TRUEtome 44 Sphincterotome	48
M00583280	TRUEtome 39 Sphincterotome	48
M00583290	TRUEtome 44 Sphincterotome	48
Dreamtome™ RX Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire		
M00584000	Dreamtome RX 49 Sphincterotome	40
M00584010	Dreamtome RX 49 Sphincterotome	40
M00584020	Dreamtome RX 49 Sphincterotome	40
M00584030	Dreamtome RX 49 Sphincterotome	40
M00584040	Dreamtome RX 44 Sphincterotome	40
M00584050	Dreamtome RX 44 Sphincterotome	40
M00584060	Dreamtome RX 44 Sphincterotome	40
M00584070	Dreamtome RX 44 Sphincterotome	40
TRUEtome™ Cannulating Sphincterotomes		
M00584120	TRUEtome 39 Sphincterotome	49
M00584121	TRUEtome 39 Sphincterotome	49
M00584122	TRUEtome 39 Sphincterotome	49
M00584150	TRUEtome 49 Sphincterotome	49
M00584151	TRUEtome 49 Sphincterotome	49
M00584152	TRUEtome 49 Sphincterotome	49
M00584160	TRUEtome 49 Sphincterotome	49
M00584161	TRUEtome 49 Sphincterotome	49
M00584162	TRUEtome 49 Sphincterotome	49
M00584170	TRUEtome 44 Sphincterotome	49
M00584171	TRUEtome 44 Sphincterotome	49



Product Code Index (CONTINUED)



Order Number	Description	Page Number
TRUEtome™ Cannulating Sphincterotomes (CONTINUED)		
M00584172	TRUEtome 44 Sphincterotome	49
M00584180	TRUEtome 44 Sphincterotome	49
M00584181	TRUEtome 44 Sphincterotome	49
M00584182	TRUEtome 44 Sphincterotome	49
M00584190	TRUEtome 39 Sphincterotome	49
M00584191	TRUEtome 39 Sphincterotome	49
M00584192	TRUEtome 39 Sphincterotome	49
Jagtome™ Revolution RX Cannulating Sphincterotome (Autotome™ Preloaded with Jagwire Revolution High Performance Guidewire)		
M00584240	Jagtome Revolution RX 39 Sphincterotome	41
M00584250	Jagtome Revolution RX 39 Sphincterotome	41
M00584260	Jagtome Revolution RX 39 Sphincterotome	41
M00584270	Jagtome Revolution RX 39 Sphincterotome	41
MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter		
M00589060	MaxForce TTS Balloon Dilatation Catheter	79
M00589080	MaxForce TTS Balloon Dilatation Catheter	79
M00589100	MaxForce TTS Balloon Dilatation Catheter	79
M00589120	MaxForce TTS Balloon Dilatation Catheter	79
M00589140	MaxForce TTS Balloon Dilatation Catheter	79
M00589150	MaxForce TTS Balloon Dilatation Catheter	79
M00589160	MaxForce TTS Balloon Dilatation Catheter	79
M00589180	MaxForce TTS Balloon Dilatation Catheter	79
SpyGlass™ Service Codes		
SRVM005ATLT1005	Autolith™ Touch Digital Essentialcare 1-Year Contract	54
SRVM005ATLT1006	Autolith Touch Digital Essentialcare 2-Years Contract	54
SRVM005ATLT1007	Autolith Touch Digital Essentialcare 3-Years Contract	54
SRVM005ATLT1008	Autolith Touch Digital Essentialcare 4-Years Contract	54
SRVM005SPYDIG1001	Point Of Sale (POS) SpyGlass Digital Essentialcare 1-Year Contract	54
SRVM005SPYDIG1002	Point Of Sale (POS) SpyGlass Digital Essentialcare 2-Years Contract	54
SRVM005SPYDIG1003	Point Of Sale (POS) SpyGlass Digital Essentialcare 3-Years Contract	54
SRVM005SPYDIG1004	Point Of Sale (POS) SpyGlass Digital Essentialcare 4-Years Contract	54
SRVM005SPYDIG1005	Post Sale SpyGlass Digital Essentialcare 1-Year Contract	54
SRVM005SPYDIG1006	Post Sale SpyGlass Digital Essentialcare 2-Years Contract	54
SRVM005SPYDIG1007	Post Sale SpyGlass Digital Essentialcare 3-Years Contract	54
SRVM005SPYDIG1008	Post Sale SpyGlass Digital Essentialcare 4-Years Contract	54
SRVM005SPYDIG1009	Point Of Sale (POS) SpyGlass Digital Evercare 1-Year Contract	54
SRVM005SPYDIG1010	Point Of Sale (POS) SpyGlass Digital Evercare 2-Years Contract	54
SRVM005SPYDIG1011	Point Of Sale (POS) SpyGlass Digital Evercare 3-Years Contract	54
SRVM005SPYDIG1012	Point Of Sale (POS) SpyGlass Digital Evercare 4-Years Contract	54
SRVM005SPYDIG1013	Post Sale SpyGlass Digital Evercare 1-Year Contract	54
SRVM005SPYDIG1014	Post Sale SpyGlass Digital Evercare 2-Years Contract	54
SRVM005SPYDIG1015	Post Sale SpyGlass Digital Evercare 3-Years Contract	54
SRVM005SPYDIG1016	Post Sale SpyGlass Digital Evercare 4-Years Contract	54
HedgeHog™ Cleaning Brushes		
SBB-382-50	HedgeHog Double Bundle	109
SBD-228-50	HedgeHog Double-End Channel Brush	109
SBD-289-50	HedgeHog Dual-End Channel/Valve Brush	109
SBD-291-50	HedgeHog Sweeper Brush	109



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Blox™ Bite Blocks		
SBP-236-40	Pediatric Blox Bite Block 38Fr	110
HedgeHog™ Cleaning Brushes		
SBS-227-50	HedgeHog Single-End Channel Brush	109
Blox™ Bite Blocks		
SBT-114-100	Blox Bite Block 54Fr with Dental Rim	110
SBT-321-100	Blox Bite Block 54Fr without Dental Rim	110
SBT-546-100	Blox Bite Block 60Fr with Dental Rim	110
Hydra™ Irrigation System		
SCT-350-10	Hydra Water Bottle Cap	109
SCT-466-10	Hydra Water Bottle Cap	109
SCT-609-10	Hydra Water Bottle Cap with CO2	109
SCT-611-10	Hydra Water Bottle Cap with CO2	109
Koala™ Sponge		
SES-240-50	Koala Cleaning Sponge Imbibed w/ Pure	110
EndoArmor™ Procedural Gowns		
SGE-368-100	EndoArmor GI Procedure Gowns	110
SGE-389-100	EndoArmor+ Procedure Gowns	110
Hydra™ Irrigation System		
SIT-351-10	Hydra Water Bottle Cap with CO2	109
SIT-355-15	Hydra Irrigation Tubing	109
SIT-361-100	Hydra Disposable Auxiliary Water Jet Connectors	109
SIT-467-10	Hydra Water Bottle Cap with CO2	109
SIT-470-100	Hydra Disposable Auxiliary Water Jet Connectors	109
SIT-575-15	Hydra Irrigation Tubing	109
SIT-577-15	Hydra Irrigation Tubing	109
SIT-611-10	Hydra Water Bottle Cap with CO2	109
EndoGlide™ Lubricating Gels		
SLT-394-100	1.1 oz EndoGlide+ Lubricating Gel	110
SLT-395-50	2 oz EndoGlide+ Lubricating Gel	110
CinchPad™ Procedure & Endoscope Transport Pad		
SPX-017-25	CinchPad Endoscope Transport Pad	110
Koala™ Sponge		
SSB-210-50	Koala Cleaning Sponge	110
TrapEase™ Polyp Traps		
STE-297-50	TrapEase Single Trap	93
STE-298-25	TrapEase Quad Trap	93
Orca™ Air/Water and Suction Valves		
SUV-617-50	Orca Air/Water and Suction Valves	109






CANNULAS



Rapid Exchange XL Cannula ⓧ

Order Number	Description	Shaft O.D. (F) / (mm)	Length (cm)	Recommended Guidewire (in) / (mm)
M00545820	Rapid Exchange XL Cannula	7-5.5 / 2.33-1.84	195	.035 / 0.89

Packaged one per box




RX ERCP Cannulas ⓧ

Order Number	Description	Shaft O.D. (F) / (mm)	Length (cm)	Recommended Guidewire (in) / (mm)
M00545760	RX ERCP Cannula Standard Tip	8.5-5 / 2.83-1.67	210	.035 / 0.89
M00545770	RX ERCP Cannula Tapered Tip	8.5-5 / 2.83-1.67	210	.035 / 0.89
M00545780	RX ERCP Cannula Ball Tip	8.5-5 / 2.83-1.67	210	.035 / 0.89

Packaged one per box

SPHINCTEROTOMES



Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire




Order Number	Description	Cut Wire (mm)	Shaft O.D. (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	Catheter Length (cm)	PRE-LOADED GUIDEWIRE** Length (cm)	O.D. (in) / (mm)
M00583000	Hydratome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	260	.035 / 0.89
M00583010	Hydratome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	260	.035 / 0.89
M00583020	Hydratome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	450	.035 / 0.89
M00583030	Hydratome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	450	.035 / 0.89
M00583040*	Hydratome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	260	.035 / 0.89
M00583050*	Hydratome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	260	.035 / 0.89
M00583060*	Hydratome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	450	.035 / 0.89
M00583070*	Hydratome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	450	.035 / 0.89

Packaged one per box

* Includes handle brake

** Pre-loaded with .035in (0.89mm) Hydra Jagwire™ Guidewire



Dreamtome™ RX Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire




Order Number	Description	Cut Wire (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	PRE-LOADED GUIDEWIRE** Length (cm)	O.D. (in) / (mm)
M00584000	Dreamtome RX 49	20	5	4.9 / 1.63	260	.035 / 0.89
M00584010	Dreamtome RX 49	30	5	4.9 / 1.63	260	.035 / 0.89
M00584020	Dreamtome RX 49	20	5	4.9 / 1.63	450	.035 / 0.89
M00584030	Dreamtome RX 49	30	5	4.9 / 1.63	450	.035 / 0.89
M00584040*	Dreamtome RX 44	20	5	4.4 / 1.47	260	.035 / 0.89
M00584050*	Dreamtome RX 44	30	5	4.4 / 1.47	260	.035 / 0.89
M00584060*	Dreamtome RX 44	20	5	4.4 / 1.47	450	.035 / 0.89
M00584070*	Dreamtome RX 44	30	5	4.4 / 1.47	450	.035 / 0.89

Packaged one per box

* Includes handle brake

** Pre-loaded with Straight-tip guidewires only

SPHINCTEROTOMES (CONTINUED)




Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire

Order Number	Description	Cut Wire (mm)	Shaft O.D. (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	Catheter Length (cm)	PRE-LOADED GUIDEWIRE™		
							Shaft Stiffness	Length (cm)	O.D. (in) / (mm)
M00573000	Jagtome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Standard	260	.035 / 0.89
M00573020	Jagtome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Standard	260	.035 / 0.89
M00573010	Jagtome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Standard	450	.035 / 0.89
M00573030	Jagtome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Standard	450	.035 / 0.89
M00573500	Jagtome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Stiff	260	.035 / 0.89
M00573520	Jagtome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Stiff	260	.035 / 0.89
M00573510	Jagtome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Stiff	450	.035 / 0.89
M00573530	Jagtome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Stiff	450	.035 / 0.89
M00573040*	Jagtome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	260	.035 / 0.89
M00573060*	Jagtome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	260	.035 / 0.89
M00573050*	Jagtome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	450	.035 / 0.89
M00573070*	Jagtome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	450	.035 / 0.89
M00573540*	Jagtome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Stiff	260	.035 / 0.89
M00573560*	Jagtome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Stiff	260	.035 / 0.89
M00573550*	Jagtome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Stiff	450	.035 / 0.89
M00573570*	Jagtome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Stiff	450	.035 / 0.89
M00573080*	Jagtome RX 39	20	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Stiff	260	.025 / 0.64
M00573100*	Jagtome RX 39	30	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Stiff	260	.025 / 0.64
M00573090*	Jagtome RX 39	20	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Standard	450	.025 / 0.64
M00573110*	Jagtome RX 39	30	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Standard	450	.025 / 0.64

Packaged one per box

* Includes handle brake


** Pre-loaded with the high performance Jagwire™ Guidewire (Jagtome RX 49 and RX 44 contain the Stiff Shaft Jagwire and Jagtome RX 39 the Standard Jagwire)



NEW Jagtome™ Revolution RX Cannulating Sphincterotome (Autotome™ Preloaded with Jagwire Revolution High Performance Guidewire)

Order Number	Description	Cut Wire (mm)	Tip Length (mm)	Distal Tip O.D. (F)	Preloaded Guidewire (in) / (cm)	Guidewire	Tip Style
M00584240	Jagtome Revolution RX 39 Sphincterotome	20	5	3.9	.025 / 260	Jagwire Revolution Guidewire	Straight
M00584250	Jagtome Revolution RX 39 Sphincterotome	20	5	3.9	.025 / 450	Jagwire Revolution Guidewire	Straight
M00584260	Jagtome Revolution RX 39 Sphincterotome	30	5	3.9	.025 / 260	Jagwire Revolution Guidewire	Straight
M00584270	Jagtome Revolution RX 39 Sphincterotome	30	5	3.9	.025 / 450	Jagwire Revolution Guidewire	Straight

SPHINCTEROTOMES (CONTINUED)



Autotome™ RX Cannulating Sphincterotomes

Order Number	Description	Cut Wire (mm)	Shaft O.D. (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	Catheter Length (cm)	Recommended Guidewire (in) / (mm)
M00545150	Autotome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	.035 / 0.89
M00545160	Autotome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	.035 / 0.89
M00545170*	Autotome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	.035 / 0.89
M00545180*	Autotome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	.035 / 0.89
M00545190*	Autotome RX 39	20	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	.025 / 0.64
M00545200*	Autotome RX 39	30	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	.025 / 0.64

Packaged one per box

Recommended Guidewires: For Autotome RX 44 and 49, recommended guidewire is .035in (0.89mm) Dreamwire™ High Performance Guidewire

For Autotome RX 39, recommended guidewire is .025in (0.69mm) Jagwire™ Guidewire

* Includes handle brake


RX Needle-Knife XL Triple-Lumen Needle-Knife

Order Number	Description	Cut Wire (mm)	Shaft O.D. (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)
M00545840	RX Needle-Knife	5	7-5.5 / 2.3-1.83	5	5.5 / 1.8

Packaged one per box

Recommended Guidewire: .035 in (0.89mm) Dreamwire™ High Performance Guidewire

STONE RETRIEVAL




Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Distal Injection



Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547000	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Above	7-6 / 2.33-2.00	200
M00547010	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Above	7-6 / 2.33-2.00	200
M00547020	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Above	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Proximal Injection




Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547030	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Below	7-6 / 2.33-2.00	200
M00547040	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Below	7-6 / 2.33-2.00	200
M00547050	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Below	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

STONE RETRIEVAL (CONTINUED)




Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Distal Injection

Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547300	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Above	7-6 / 2.33-2.00	200
M00547310	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Above	7-6 / 2.33-2.00	200
M00547320	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Above	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire




Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Proximal Injection

Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547330	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Below	7-6 / 2.33-2.00	200
M00547340	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Below	7-6 / 2.33-2.00	200
M00547350	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Below	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Trapezoid™ RX Wireguided Retrieval Baskets

Order Number	Description	Open Basket Diameter (cm)	Required Working Channel (mm)
M00510860	Trapezoid RX Wireguided Retrieval Basket	1.5	3.2
M00510870	Trapezoid RX Wireguided Retrieval Basket	2.0	3.2
M00510880	Trapezoid RX Wireguided Retrieval Basket	2.5	3.2
M00510890*	Trapezoid RX Wireguided Retrieval Basket	3.0	3.2
M00550620	Alliance™ II Inflation Handle**		


Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

* The Trapezoid 3 x 6cm basket is designed for crushing calculus larger than 1.5cm in diameter

** Alliance™ II recommended for use as Mechanical Lithotripter. Alliance II is reusable

DILATION



Hurricane™ RX Biliary Dilation Balloons



Order Number	Description	Balloon O.D. (mm)	Balloon Length (cm)	Catheter O.D. (F) / (mm)	Maximum Inflation Pressure (ATM) / (kPa)	Minimum Channel Size (mm)	Recommended Guidewire (in) / (mm)
M00545890	Hurricane RX Dilation Balloon	4	2	5.8 / 1.93	12 / 1216	3.2	.035 / 0.89
M00545900	Hurricane RX Dilation Balloon	4	4	5.8 / 1.93	12 / 1216	3.2	.035 / 0.89
M00545910	Hurricane RX Dilation Balloon	6	2	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545920	Hurricane RX Dilation Balloon	6	4	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545930	Hurricane RX Dilation Balloon	8	2	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545940	Hurricane RX Dilation Balloon	8	4	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545950	Hurricane RX Dilation Balloon	10	2	5.8 / 1.93	8 / 810	3.2	.035 / 0.89
M00545960	Hurricane RX Dilation Balloon	10	4	5.8 / 1.93	8 / 810	3.2	.035 / 0.89

Packaged one per box


Recommended Inflation Device: Encore™ 26 Inflation Device, M00566670

CRE™ RX Biliary Balloon Dilatation Catheter



Order Number	Inflated Balloon OD (mm)	Balloon Dilation Length (cm)	Inflation Pressure (ATM)	Inflation Pressure (kPa)	Catheter Size (Fr)	Minimum Working Channel (mm)	Working Length (cm)	
M00558920	8-9-10	24-27-30	3.0	3-5.5-9	304-608-1013	7.5	3.7	180
M00558930	10-11-12	30-33-36	3.0	3-5-8	304-557-912	7.5	3.7	180
M00558940	12-13.5-15	36-40.5-45	3.0	3-4.5-8	304-507-811	7.5	3.7	180
M00558860	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	7.5	3.7	180
M00558870	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	7.5	3.7	180
M00558880	10-11-12	30-33-36	5.5	3-5-8	304-507-811	7.5	3.7	180
M00558890	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	7.5	3.7	180
M00558900	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	7.5	3.7	180
M00558910	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	7.5	3.7	180

CYTOLOGY



RX Cytology Brush



Order Number	Description	Brush Width (mm)	Catheter O.D. (F) / (mm)	Working Length (cm)	Recommended Guidewire (in) / (mm)
M00545000	RX Cytology Brush	2.1	8 / 2.67	200	.035 / 0.89

Packaged one per box


ACCESSORIES

RX Accessories



Order Number	Description	Units
M00545211†	RX Locking Device – Olympus / FujiFilm Brand Endoscope Compatible	box 10
M00545221†	RX Locking Device – Pentax Brand Endoscope Compatible	box 10
M00545261	RX Locking Device and Biopsy Cap System – Olympus / FujiFilm G5 Series Brand Endoscope Compatible	box 10
M00545271	RX Locking Device and Biopsy Cap System – Pentax Brand Endoscope Compatible	box 10
M00545281	RX Locking Device and Biopsy Cap System – FujiFilm Type 3 Series Brand Endoscope Compatible	box 10
M00566001	Rapid Refill Continuous Injection System	box 10

† Product not available in the European Economic Area (EEA), except in Greece. Not for use or distribution in the EEA, except in Greece.
Please check availability of the product with your local sales representative or your local customer service.





Encore™ 26 Inflation Device




Order Number	Description	Units
M00566670	Encore 26 Inflation Device	box 1

RX Biliary Metal and Plastic Stents are located in the Stents Section on pages 57 – 59 and 60 – 63

Biliary Guidewires are located in the Guidewires Section on page 81



CANNULAS



Injection Lumens
Guidewire Lumen


Tapered .035in
(0.89mm)

Tandem™ XL Triple-Lumen ERCP Cannula



Order Number	Description	Shaft O.D. (F) / (mm)	Distal O.D. (F) / (mm)	Length (cm)	Tip	Recommended Guidewire (in) / (mm)
M00535700	Tandem XL Triple-Lumen ERCP Cannula	7.0-5.5 / 2.33-1.83	5.5 / 1.83	200	Tapered	.035 / 0.89

Packaged one per box



Standard .035in
(0.89mm)

Tapered .035in
(0.89mm)

Ultra Tapered .025in
(0.64mm)

5-4-3 .018in
(0.46mm)

Ball .035in
(0.89mm)


Contour™ ERCP Cannulas



Order Number	Description	Shaft O.D. (F) / (mm)	Length (cm)	Tip	Recommended Guidewire (in) / (mm)
M00530850	Contour ERCP Cannula	5 / 1.67	210	Standard	.035 / 0.89
M00530860	Contour ERCP Cannula	5 / 1.67	210	Tapered	.035 / 0.89
M00530870	Contour ERCP Cannula	5 / 1.67	210	Ultra Tapered	.025 / 0.64
M00530880	Contour ERCP Cannula	5 / 1.67	210	5-4-3	.018 / 0.46
M00530890	Contour ERCP Cannula	5 / 1.67	210	Ball Tip	.035 / 0.89

Packaged one per box

SPHINCTEROTOMES



Above Cut Wire

Below Cut Wire

Stonetome™ Stone Removal Devices




Order Number	Description	Cut Wire (mm)	Balloon Diameter (mm)	Tip Length (mm)	Catheter Length (cm)
M00535150	Stonetome (below cut wire)	20	11.5	5	200
M00535170	Stonetome (below cut wire)	20	11.5	20	200
M00535190	Stonetome (below cut wire)	30	11.5	5	200
M00535210	Stonetome (below cut wire)	30	11.5	20	200
M00535110	Stonetome (above cut wire)	20	11.5	20	200
M00535130	Stonetome (above cut wire)	30	11.5	20	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

SPHINCTEROTOMES (CONTINUED)



TRUEtome™ Jag Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire



Order Number	Description	Cut Wire (mm)	Tip O.D. (F) / (mm)	Tip Length (mm)	Pre-loaded Guidewire (in / cm)	Units
M00583080	TRUEtome 39 Sphincterotome	20	3.9 / 1.30	5	.025 / 450 Jagwire Guidewire	box 1
M00583090	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.025 / 450 Jagwire Guidewire	box 1
M00583100	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.035 / 450 Jagwire Guidewire	box 1
M00583110	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.035 / Stiff 450 Jagwire Guidewire	box 1
M00583120	TRUEtome 39 Sphincterotome	30	3.9 / 1.30	5	.025 / 450 Jagwire Guidewire	box 1
M00583130	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.025 / 450 Jagwire Guidewire	box 1
M00583140	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.035 / 450 Jagwire Guidewire	box 1
M00583150	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.035 / Stiff 450 Jagwire Guidewire	box 1

TRUEtome™ Dreamwire™ Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire



Order Number	Description	Cut Wire (mm)	Tip O.D. (F) / (mm)	Tip Length (mm)	Pre-loaded Guidewire (in / cm)	Units
M00583190	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.035 / 450 Dreamwire Guidewire	box 1
M00583210	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.035 / 450 Dreamwire Guidewire	box 1

TRUEtome™ Hydra Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire




Order Number	Description	Cut Wire (mm)	Tip O.D. (F) / (mm)	Tip Length (mm)	Pre-loaded Guidewire (in / cm)	Units
M00583160	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.035 / 450 Hydra Jagwire Guidewire	box 1
M00583170	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.035 / 450 Hydra Jagwire Guidewire	box 1

NEW TRUEtome™ Jagwire Revolution Cannulating Sphincterotomes with Pre-loaded Jagwire™ Revolution Guidewire



Order Number	Description	Cut Wire (mm)	Tip O.D. (F) / (mm)	Tip Length (mm)	Pre-loaded Guidewire (in / cm)	Units
M00583260	TRUEtome 39 Sphincterotome	20	3.9 / 1.30	5	.025 / 450 Jagwire Revolution Guidewire	box 1
M00583270	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.025 / 450 Jagwire Revolution Guidewire	box 1
M00583280	TRUEtome 39 Sphincterotome	30	3.9 / 1.30	5	.025 / 450 Jagwire Revolution Guidewire	box 1
M00583290	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.025 / 450 Jagwire Revolution Guidewire	box 1

SPHINCTEROTOMES (CONTINUED)




TRUEtome™ Cannulating Sphincterotomes



Order Number	Description	Cut Wire (mm)	Distal O.D. (F) / (mm)	Tip Length (mm)	Catheter Length (cm)	Units
M00584190	TRUEtome 39	20	3.9 / 1.30	5	200	box 1
M00584120	TRUEtome 39	30	3.9 / 1.30	5	200	box 1
M00584191	TRUEtome 39	20	3.9 / 1.30	5	200	box 10
M00584121	TRUEtome 39	30	3.9 / 1.30	5	200	box 10
M00584192	TRUEtome 39	20	3.9 / 1.30	5	200	box 20
M00584122	TRUEtome 39	30	3.9 / 1.30	5	200	box 20
M00584170	TRUEtome 44	20	4.4 / 1.47	5	200	box 1
M00584180	TRUEtome 44	30	4.4 / 1.47	5	200	box 1
M00584171	TRUEtome 44	20	4.4 / 1.47	5	200	box 10
M00584181	TRUEtome 44	30	4.4 / 1.47	5	200	box 10
M00584172	TRUEtome 44	20	4.4 / 1.47	5	200	box 20
M00584182	TRUEtome 44	30	4.4 / 1.47	5	200	box 20
M00584150	TRUEtome 49	20	4.9 / 1.63	5	200	box 1
M00584160	TRUEtome 49	30	4.9 / 1.63	5	200	box 1
M00584151	TRUEtome 49	20	4.9 / 1.63	5	200	box 10
M00584161	TRUEtome 49	30	4.9 / 1.63	5	200	box 10
M00584152	TRUEtome 49	20	4.9 / 1.63	5	200	box 20
M00584162	TRUEtome 49	30	4.9 / 1.63	5	200	box 20

Recommended Guidewire: .025in (0.64 mm) for TRUEtome 39; .035in (0.89 mm) for TRUEtome 44 and 49

SPHINCTEROTOMES (CONTINUED)




Ultratome XL™ Triple-Lumen Sphincterotomes



Order Number	Description	Cut Wire (mm)	Distal O.D. (F) / (mm)	Tip Length (mm)	Catheter Length (cm)
M00535900	Ultratome XL Short Nose	20	5.5 / 1.83	5	200
M00535910	Ultratome XL Long Nose	20	5.5 / 1.83	20	200
M00535920	Ultratome XL Short Nose	30	5.5 / 1.83	5	200
M00535930	Ultratome XL Long Nose	30	5.5 / 1.83	20	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire




Ultratome™ Double-Lumen Sphincterotomes



Order Number	Description	Cut Wire (mm)	Distal O.D. (F) / (mm)	Tip Length (mm)	Catheter Length (cm)
M00530800	Ultratome Short Nose	20	5.5 / 1.83	5	200
M00530810	Ultratome Long Nose	20	5.5 / 1.83	20	200
M00530820	Ultratome Short Nose	30	5.5 / 1.83	5	200
M00530830	Ultratome Long Nose	30	5.5 / 1.83	20	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Microknife™ XL Triple-Lumen Needle-Knife




Order Number	Description	Shaft O.D. (F) / (mm)	Catheter Length (cm)
M00532810	Microknife XL Triple-Lumen Needle-Knife	7-5.5 / 2.33-1.83	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

STONE RETRIEVAL




Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Distal Injection

Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547100	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Above	7-6 / 2.33-2.00	200
M00547110	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Above	7-6 / 2.33-2.00	200
M00547120	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Above	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire




Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Proximal Injection

Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547130	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Below	7-6 / 2.33-2.00	200
M00547140	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Below	7-6 / 2.33-2.00	200
M00547150	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Below	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire




Extractor™ Pro DL Double-Lumen Retrieval Balloon Catheters

Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547200	Extractor Pro DL Double-Lumen Retrieval Balloon Catheter	9.0-12.0	N/A	7-6 / 2.33-2.00	200
M00547210	Extractor Pro DL Double-Lumen Retrieval Balloon Catheter	12.0-15.0	N/A	7-6 / 2.33-2.00	200
M00547220	Extractor Pro DL Double-Lumen Retrieval Balloon Catheter	15.0-18.0	N/A	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

DILATION



MaxForce™ Biliary Balloon Dilatation Catheter



Order Number	Description	Inflated O.D. (mm)	Balloon Length (cm)	Catheter Size (F) / (mm)	Required Working Channel (mm)
M00567330	MaxForce Balloon Dilatation Catheter	4	2	5 / 1.67	2.8
M00567340	MaxForce Balloon Dilatation Catheter	4	4	5 / 1.67	2.8
M00567370	MaxForce Balloon Dilatation Catheter	6	2	5 / 1.67	2.8
M00567380	MaxForce Balloon Dilatation Catheter	6	4	5 / 1.67	2.8
M00567410	MaxForce Balloon Dilatation Catheter	8	2	5 / 1.67	2.8
M00567420	MaxForce Balloon Dilatation Catheter	8	3	5 / 1.67	2.8

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

Recommended Inflation Device: Encore™ 26 Inflation Device, M00566670

ACCESSORIES

Biliary Accessories


Order Number	Description	Units
M00561250	Active Cord, Fits Olympus Brand Generator	box 1
M00561270	Active Cord, Fits Endostat™, Endostat™ II, Endostat™ III, Bovie®, Valleylab®, Aspen Labs and ERBE® Brand Generators	box 1
M00566670	⊗ Encore 26 Inflation Device	box 1
M00566001	⊗ Rapid Refill Continuous Injection System	box 10

Biliary Metal and Plastic Stents are located in the Stents Section on pages 59 and 64 – 65

Biliary Guidewires are located in the Guidewires Section on page 81



SPYGLASS™ DS II



SpyGlass™ DS Capital Equipment

Order Number	Description	Units
M00546650	SpyGlass DS II Digital Controller	box 1

SpyGlass™ DS Disposable Equipment



Order Number	Description	Units
M00546610	SpyScope™ DS II Access and Delivery Catheter	box 1
M00546550	SpyGlass™ Retrieval Basket	
M00546560	SpyGlass™ Retrieval Snare	
M00546470	SpyBite™ Max Biopsy Forceps	box 1
SIT-575-15	Hydra Irrigation Tubing	box 15

NEW Autolith™ Touch Biliary EHL System

Autolith™ Touch System

Order Number	Description
M005466800	Autolith Touch EHL Generator

Biliary EHL Probe

Order Number	Description
M00546620	1.9Fr, 375cm Biliary EHL Probe


Autolith™ Touch System Accessories

Order Number	Description
M00546750	Autolith Touch Extender Cable
M00546760	Autolith Touch Foot Pedal

SPYGLASS™

Service Codes

1-Year Contract	2-Years Contract	3-Years Contract	4-Years Contract
Point Of Sale (POS) Spyglass Digital Evercare			
SRVM005SPYDIG1009	SRVM005SPYDIG1010	SRVM005SPYDIG1011	SRVM005SPYDIG1012
Post Sale Spyglass Digital Evercare			
SRVM005SPYDIG1013	SRVM005SPYDIG1014	SRVM005SPYDIG1015	SRVM005SPYDIG1016
POS Spyglass Digital Essentialcare			
SRVM005SPYDIG1001	SRVM005SPYDIG1002	SRVM005SPYDIG1003	SRVM005SPYDIG1004
Post Sale Spyglass Digital Essentialcare			
SRVM005SPYDIG1005	SRVM005SPYDIG1006	SRVM005SPYDIG1007	SRVM005SPYDIG1008
Autolith Touch Digital Essentialcare			
SRVM005ATLT1005	SRVM005ATLT1006	SRVM005ATLT1007	SRVM005ATLT1008



RX Biliary System™	Biliary	SpyGlass™	EUS	Stents	Dilation	Guidewires	Hemostasis	Tissue Resection	Biopsy	Enteral Feeding	Pulmonary	Infection Prevention	RFA
--------------------	---------	-----------	-----	--------	----------	------------	------------	------------------	--------	-----------------	-----------	----------------------	-----



Contents




Alphabetical Index

Product Code Index

Endoscopic Ultrasound (EUS)

FINE NEEDLE ASPIRATION (FNA)




Expect™ Standard Handle Needles*



Order Number	Description	Needle Size (gauge) / (mm)	Minimum Working Channel (mm)	Sheath Diameter (mm)	Units
M00550000	Expect Standard Handle Needle	19 / 1.10	2.8	1.83	box 1
M00550001	Expect Standard Handle Needle	19 / 1.10	2.8	1.83	box 5
M00550040	Expect Standard Handle Needle	19 Flex / 1.14	2.8	1.73	box 1
M00550041	Expect Standard Handle Needle	19 Flex / 1.14	2.8	1.73	box 5
M00550010	Expect Standard Handle Needle	22 / 0.72	2.4	1.65	box 1
M00550011	Expect Standard Handle Needle	22 / 0.72	2.4	1.65	box 5
M00550020	Expect Standard Handle Needle	25 / 0.52	2.4	1.52	box 1
M00550021	Expect Standard Handle Needle	25 / 0.52	2.4	1.52	box 5

* Includes a 20cc syringe and one-way stopcock
Working length: 137.5cm to 141.5cm, adjustable
Needle length: 0cm to 8cm, adjustable



Expect™ Slimline Handle Needles*



Order Number	Description	Needle Size (gauge) / (mm)	Minimum Working Channel (mm)	Sheath Diameter (mm)	Units
M00555500	Expect Slimline Handle Needle	19 / 1.10	2.8	1.83	box 1
M00555501	Expect Slimline Handle Needle	19 / 1.10	2.8	1.83	box 5
M00555530	Expect Slimline Handle Needle	19 Flex / 1.14	2.8	1.73	box 1
M00555531	Expect Slimline Handle Needle	19 Flex / 1.14	2.8	1.73	box 5
M00555510	Expect Slimline Handle Needle	22 / 0.72	2.4	1.65	box 1
M00555511	Expect Slimline Handle Needle	22 / 0.72	2.4	1.65	box 5
M00555520	Expect Slimline Handle Needle	25 / 0.52	2.4	1.52	box 1
M00555521	Expect Slimline Handle Needle	25 / 0.52	2.4	1.52	box 5

* Includes a 20cc syringe and one-way stopcock
Working length: 137.5cm to 141.5cm, adjustable
Needle length: 0cm to 8cm, adjustable



Contents




Alphabetical Index

Product Code Index

Endoscopic Ultrasound (EUS)

FINE NEEDLE BIOPSY (FNB)



Acquire™ EUS-FNB Devices



Order Number	Description	Needle Size (gauge) / (mm)	Minimum Working Channel (mm)	Sheath Diameter (mm)	Units
M00555580	19ga Acquire FNB Needle, Flexible	19 / 1.14	2.8	1.73	box 1
M00555540	22ga Acquire EUS-FNB Needle	22 / 0.72	2.4	1.65	box 1
M00555541	22ga Acquire EUS-FNB Needle	22 / 0.72	2.4	1.65	box 5
M00555560	25ga Acquire EUS-FNB Needle	25 / 0.52	2.4	1.52	box 1
M00555561	25ga Acquire EUS-FNB Needle	25 / 0.52	2.4	1.52	box 5

STENT AND DELIVERY SYSTEM

Hot AXIOS™ Stent and Electrocautery Enhanced Delivery System



Order Number	Description	Flange Diameter (mm)	Lumen Diameter (mm)	Saddle Length (mm)	Catheter Outer Diameter (F)	Catheter Working Length (cm)	Catheter Total Length (cm)	Units
M00553520	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	14	6	8	9	138	146	box 1
M00553530	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	17	8	8	9	138	146	box 1
M00553540	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	21	10	10	10.8	138	146	box 1
M00553550	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	24	15	10	10.8	138	146	box 1
M00553560	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	29	20	10	10.8	138	146	box 1

AXIOS™ Stent and Delivery System




Order Number	Description	Flange Diameter (mm)	Lumen Diameter (mm)	Saddle Length (mm)	Catheter Outer Diameter (F)	Catheter Working Length (cm)	Catheter Total Length (cm)	Units
M00553440	AXIOS Stent and Delivery System	21	10	10	10.8	138	146	box 1
M00553450	AXIOS Stent and Delivery System	24	15	10	10.8	138	146	box 1



Contents



RX BILIARY METAL STENTS




WallFlex™ Biliary RX Fully Covered Stent System RMV (Permalume™ Covering, includes benign indication)



Order Number	Description	Diameter (mm)	Length (mm)	Catheter Diameter (F) / (mm)	Guidewire Diameter (in) / (mm)
M00574210†	WallFlex RX Biliary w/Permalume	6	40	8.5 / 2.83	.035 / 0.89
M00574230†	WallFlex RX Biliary w/Permalume	6	60	8.5 / 2.83	.035 / 0.89
M00574240†	WallFlex RX Biliary w/Permalume	6	80	8.5 / 2.83	.035 / 0.89
M00570450†	WallFlex RX Biliary w/Permalume	8	60	8.5 / 2.83	.035 / 0.89
M00570460†	WallFlex RX Biliary w/Permalume	8	80	8.5 / 2.83	.035 / 0.89
M00576800†	WallFlex RX Biliary w/Permalume	8	100	9.0 / 3.00	.035 / 0.89
M00576810†	WallFlex RX Biliary w/Permalume	8	120	9.0 / 3.00	.035 / 0.89
M00570470†	WallFlex RX Biliary w/Permalume	10	40	8.5 / 2.83	.035 / 0.89
M00570480†	WallFlex RX Biliary w/Permalume	10	60	8.5 / 2.83	.035 / 0.89
M00570490†	WallFlex RX Biliary w/Permalume	10	80	8.5 / 2.83	.035 / 0.89
M00576820†	WallFlex RX Biliary w/Permalume	10	100	9.0 / 3.00	.035 / 0.89
M00576830†	WallFlex RX Biliary w/Permalume	10	120	9.0 / 3.00	.035 / 0.89

† Product commercialised only in countries where benign indication is approved.



WallFlex™ Biliary RX Partially Covered Stent System (Permalume™ Covering)




Order Number	Description	Diameter (mm)	Length (mm)	Covered Length (mm)	Catheter Diameter (F) / (mm)	Guidewire Diameter (in) / (mm)
M00570700	WallFlex RX Biliary w/Permalume	8	60	48	8.5 / 2.83	.035 / 0.89
M00570710	WallFlex RX Biliary w/Permalume	8	80	68	8.5 / 2.83	.035 / 0.89
M00576740	WallFlex RX Biliary w/Permalume	8	100	88	9.0 / 3.00	.035 / 0.89
M00576750	WallFlex RX Biliary w/Permalume	8	120	108	9.0 / 3.00	.035 / 0.89
M00570720	WallFlex RX Biliary w/Permalume	10	40	28	8.5 / 2.83	.035 / 0.89
M00570730	WallFlex RX Biliary w/Permalume	10	60	48	8.5 / 2.83	.035 / 0.89
M00570740	WallFlex RX Biliary w/Permalume	10	80	68	8.5 / 2.83	.035 / 0.89
M00576760	WallFlex RX Biliary w/Permalume	10	100	88	9.0 / 3.00	.035 / 0.89
M00576770	WallFlex RX Biliary w/Permalume	10	120	108	9.0 / 3.00	.035 / 0.89

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm or .035in (0.89mm) Dreamwire™ High Performance Guidewire, 260cm

RX BILIARY METAL STENTS (CONTINUED)




WallFlex™ Biliary RX Uncovered Stent System



Order Number	Description	Diameter (mm)	Length (mm)	Catheter Diameter (F) / (mm)	Guidewire Diameter (in) / (mm)
M00570600	WallFlex RX Biliary	8	40	8 / 2.67	.035 / 0.89
M00570610	WallFlex RX Biliary	8	60	8 / 2.67	.035 / 0.89
M00570620	WallFlex RX Biliary	8	80	8 / 2.67	.035 / 0.89
M00570630	WallFlex RX Biliary	8	100	8 / 2.67	.035 / 0.89
M00576780	WallFlex RX Biliary	8	120	8 / 2.67	.035 / 0.89
M00570890	WallFlex RX Biliary	10	40	8 / 2.67	.035 / 0.89
M00570640	WallFlex RX Biliary	10	60	8 / 2.67	.035 / 0.89
M00570650	WallFlex RX Biliary	10	80	8 / 2.67	.035 / 0.89
M00570660	WallFlex RX Biliary	10	100	8 / 2.67	.035 / 0.89
M00576790	WallFlex RX Biliary	10	120	8 / 2.67	.035 / 0.89

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm or .035in (0.89mm) Dreamwire™ High Performance Guidewire, 260cm



Wallstent™ RX Biliary Endoprosthesis Stent System (Permalume™ Covering)



Order Number	Description	Covered Length (mm)	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System (F) / (mm)	Working Length (cm)
M00569670	Wallstent RX w/Permalume	30	8	40	8 / 2.67	194
M00569680	Wallstent RX w/Permalume	50	8	60	8 / 2.67	194
M00569690	Wallstent RX w/Permalume	70	8	80	8 / 2.67	194
M00569700	Wallstent RX w/Permalume	30	10	40	8 / 2.67	194
M00569710	Wallstent RX w/Permalume	50	10	60	8 / 2.67	194
M00569720	Wallstent RX w/Permalume	70	10	80	8 / 2.67	194

Packaged one per box

Recommended Guidewires: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm or .035in (0.89mm) Dreamwire™ High Performance Guidewire, 260cm

RX BILIARY METAL STENTS (CONTINUED)

Wallstent™ RX Biliary Endoprosthesis Stent System (Uncovered)




Order Number	Description	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System (F) / (mm)	Working Length (cm)
M00569610	Wallstent RX Biliary	8	40	8 / 2.67	194
M00569620	Wallstent RX Biliary	8	60	8 / 2.67	194
M00569630	Wallstent RX Biliary	8	80	8 / 2.67	194
M00569800	Wallstent RX Biliary	8	100	8 / 2.67	194
M00569640	Wallstent RX Biliary	10	40	8 / 2.67	194
M00569650	Wallstent RX Biliary	10	60	8 / 2.67	194
M00569660	Wallstent RX Biliary	10	80	8 / 2.67	194
M00569810	Wallstent RX Biliary	10	100	8 / 2.67	194

Packaged one per box

Recommended Guidewires: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm or .035in (0.89mm) Dreamwire™ High Performance Guidewire, 260cm

BILIARY METAL STENTS



Wallstent™ Endoscopic Biliary Endoprosthesis Stent System (Partially Covered)



Order Number	Description	Covered Length (mm)	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System (F) / (mm)	Working Length (cm)
H965430700	Wallstent Biliary w/Permalume	30	8	40	8 / 2.67	194
H965430800	Wallstent Biliary w/Permalume	50	8	60	8 / 2.67	194
H965430900	Wallstent Biliary w/Permalume	70	8	80	8 / 2.67	194
H965431000	Wallstent Biliary w/Permalume	30	10	40	8 / 2.67	194
H965431100	Wallstent Biliary w/Permalume	50	10	60	8 / 2.67	194
H965431200	Wallstent Biliary w/Permalume	70	10	80	8 / 2.67	194

Packaged one per box

Recommended Guidewires: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft or .035in (0.89mm) Dreamwire™ High Performance Guidewire

Wallstent™ Endoscopic Biliary Endoprosthesis Stent System (Uncovered)




Order Number	Description	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System (F) / (mm)	Working Length (cm)
H965430100	Wallstent Biliary	8	40	7.5 / 2.5	194
H965430200	Wallstent Biliary	8	60	7.5 / 2.5	194
H965430300	Wallstent Biliary	8	80	7.5 / 2.5	194
H965430400	Wallstent Biliary	10	40	7.5 / 2.5	194
H965430500	Wallstent Biliary	10	60	7.5 / 2.5	194
H965430600	Wallstent Biliary	10	80	7.5 / 2.5	194

Packaged one per box

Recommended Guidewires: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft or .035in (0.89mm) Dreamwire™ High Performance Guidewire

RX BILIARY PLASTIC STENTS



Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System



Order Number	Description	Shape	O.D. (F) / (mm)	Length (cm) Barb to Barb
M00534200	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	5
M00534210	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	7
M00534220	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	9
M00534230	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	12
M00534240	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	15
M00534250	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	18
M00534260	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	5
M00534270	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	7
M00534280	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	9
M00534290	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	12
M00534300	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	15
M00534310	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	18
M00534320	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	5
M00534330	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	7
M00534340	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	9
M00534350	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	12
M00534360	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	15
M00534370	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	18
M00534560	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	5
M00534570	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	7
M00534580	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	9
M00534590	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	12
M00534600	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	15
M00534610	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	18
M00534620	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	5
M00534630	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	7
M00534640	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	9
M00534650	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	12
M00534660	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	15
M00534670	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	18
M00534680	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	5
M00534690	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	7
M00534700	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	9
M00534710	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	12
M00534720	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	15
M00534730	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	18

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

Scope Compatibility: 7F (2.3mm) and 8.5F (2.83mm) stents require a 3.2mm working channel; 10F (3.3mm) stents require a 4.2mm working channel

RX BILIARY PLASTIC STENTS (CONTINUED)



Advanix™ Stents – Singles




Order Number	Description	Shape	O.D. (F) / (mm)	Length (cm) Barb to Barb
M00532870	Advanix Stent	Duodenal Bend	7 / 2.3	5
M00532880	Advanix Stent	Duodenal Bend	7 / 2.3	7
M00532890	Advanix Stent	Duodenal Bend	7 / 2.3	9
M00532900	Advanix Stent	Duodenal Bend	7 / 2.3	12
M00532910	Advanix Stent	Duodenal Bend	7 / 2.3	15
M00532920	Advanix Stent	Duodenal Bend	7 / 2.3	18
M00532930	Advanix Stent	Duodenal Bend	8.5 / 2.83	5
M00532940	Advanix Stent	Duodenal Bend	8.5 / 2.83	7
M00532950	Advanix Stent	Duodenal Bend	8.5 / 2.83	9
M00532960	Advanix Stent	Duodenal Bend	8.5 / 2.83	12
M00532970	Advanix Stent	Duodenal Bend	8.5 / 2.83	15
M00532980	Advanix Stent	Duodenal Bend	8.5 / 2.83	18
M00532990	Advanix Stent	Duodenal Bend	10 / 3.3	5
M00533000	Advanix Stent	Duodenal Bend	10 / 3.3	7
M00533010	Advanix Stent	Duodenal Bend	10 / 3.3	9
M00533020	Advanix Stent	Duodenal Bend	10 / 3.3	12
M00533030	Advanix Stent	Duodenal Bend	10 / 3.3	15
M00533040	Advanix Stent	Duodenal Bend	10 / 3.3	18
M00533230	Advanix Stent	Center Bend	7 / 2.3	5
M00533240	Advanix Stent	Center Bend	7 / 2.3	7
M00533250	Advanix Stent	Center Bend	7 / 2.3	9
M00533260	Advanix Stent	Center Bend	7 / 2.3	12
M00533270	Advanix Stent	Center Bend	7 / 2.3	15
M00533280	Advanix Stent	Center Bend	7 / 2.3	18
M00533290	Advanix Stent	Center Bend	8.5 / 2.83	5
M00533300	Advanix Stent	Center Bend	8.5 / 2.83	7
M00533310	Advanix Stent	Center Bend	8.5 / 2.83	9
M00533320	Advanix Stent	Center Bend	8.5 / 2.83	12
M00533330	Advanix Stent	Center Bend	8.5 / 2.83	15
M00533340	Advanix Stent	Center Bend	8.5 / 2.83	18
M00533350	Advanix Stent	Center Bend	10 / 3.3	5
M00533360	Advanix Stent	Center Bend	10 / 3.3	7
M00533370	Advanix Stent	Center Bend	10 / 3.3	9
M00533380	Advanix Stent	Center Bend	10 / 3.3	12
M00533390	Advanix Stent	Center Bend	10 / 3.3	15
M00533400	Advanix Stent	Center Bend	10 / 3.3	18

Packaged one per box

Recommended Guidewire: .035in (0.89mm)

Scope Compatibility: 7F (2.3mm) and 8.5F (2.83mm) stents require a 3.2mm working channel; 10F (3.3mm) stents require a 4.2mm working channel

RX BILIARY PLASTIC STENTS (CONTINUED)



Advanix™ Stents – Singles (CONTINUED)



Order Number	Description	Shape	O.D. (F) / (mm)	Length (cm) Between Pigtails
M00532160	Advanix Stent	Double Pigtail	7 / 2.3	3
M00532180	Advanix Stent	Double Pigtail	7 / 2.3	5
M00532190	Advanix Stent	Double Pigtail	7 / 2.3	7
M00532210	Advanix Stent	Double Pigtail	7 / 2.3	10
M00532220	Advanix Stent	Double Pigtail	7 / 2.3	12
M00532230	Advanix Stent	Double Pigtail	7 / 2.3	15
M00532240	Advanix Stent	Double Pigtail	10 / 3.3	3
M00532260	Advanix Stent	Double Pigtail	10 / 3.3	5
M00532270	Advanix Stent	Double Pigtail	10 / 3.3	7
M00532290	Advanix Stent	Double Pigtail	10 / 3.3	10
M00532300	Advanix Stent	Double Pigtail	10 / 3.3	12
M00532310	Advanix Stent	Double Pigtail	10 / 3.3	15

Packaged one per box

Recommended Guidewire: .035in (0.89mm)

Scope Compatibility: 7F (2.3mm) and 8.5F (2.83mm) stents require a 3.2mm working channel; 10F (3.3mm) stents require a 4.2mm working channel



NaviFlex™ RX Delivery System



Order Number	Description	Catheter Diameter (mm)	Working Length (cm)
M00533560	RX Delivery System	7 / 2.3	202.5
M00533570	RX Delivery System	8.5 / 2.83	202.5
M00533580	RX Delivery System	10 / 3.3	202.5

Packaged one per box

Compatible with Advanix Biliary Stents.

RX BILIARY PLASTIC STENTS (CONTINUED)

RX Biliary Plastic Stents with RX Delivery System



Order Number	Description	O.D. (F) / (mm)	Length (cm) Barb to Barb
M00545550	RX Biliary Plastic Stent with RX Delivery System	7 / 2.33	5
M00545560	RX Biliary Plastic Stent with RX Delivery System	7 / 2.33	7
M00545570	RX Biliary Plastic Stent with RX Delivery System	7 / 2.33	10
M00545580	RX Biliary Plastic Stent with RX Delivery System	7 / 2.33	12
M00545590	RX Biliary Plastic Stent with RX Delivery System	7 / 2.33	15
M00545650	RX Biliary Plastic Stent with RX Delivery System	8.5 / 2.83	5
M00545660	RX Biliary Plastic Stent with RX Delivery System	8.5 / 2.83	7
M00545670	RX Biliary Plastic Stent with RX Delivery System	8.5 / 2.83	10
M00545680	RX Biliary Plastic Stent with RX Delivery System	8.5 / 2.83	12
M00545690	RX Biliary Plastic Stent with RX Delivery System	8.5 / 2.83	15
M00545600	RX Biliary Plastic Stent with RX Delivery System	10 / 3.33	5
M00545610	RX Biliary Plastic Stent with RX Delivery System	10 / 3.33	7
M00545620	RX Biliary Plastic Stent with RX Delivery System	10 / 3.33	10
M00545630	RX Biliary Plastic Stent with RX Delivery System	10 / 3.33	12
M00545640	RX Biliary Plastic Stent with RX Delivery System	10 / 3.33	15
M00545700	RX Biliary Plastic Stent with RX Delivery System	11.5 / 3.83	5
M00545710	RX Biliary Plastic Stent with RX Delivery System	11.5 / 3.83	7
M00545720	RX Biliary Plastic Stent with RX Delivery System	11.5 / 3.83	10
M00545730	RX Biliary Plastic Stent with RX Delivery System	11.5 / 3.83	12
M00545740	RX Biliary Plastic Stent with RX Delivery System	11.5 / 3.83	15

Packaged one per box

Recommended Guidewire: .035in (0.89mm)

BILIARY PLASTIC STENTS



Flexima™ Biliary Stents with Delivery System



Order Number	Description	O.D. (F) / (mm)	Length (cm) Barb to Barb
M00539200	Flexima Stent with Delivery System	7 / 2.33	5
M00539210	Flexima Stent with Delivery System	7 / 2.33	7
M00539220	Flexima Stent with Delivery System	7 / 2.33	10
M00539230	Flexima Stent with Delivery System	7 / 2.33	12
M00539240	Flexima Stent with Delivery System	7 / 2.33	15
M00539250	Flexima Stent with Delivery System	8.5 / 2.83	5
M00539260	Flexima Stent with Delivery System	8.5 / 2.83	7
M00539270	Flexima Stent with Delivery System	8.5 / 2.83	10
M00539280	Flexima Stent with Delivery System	8.5 / 2.83	12
M00539290	Flexima Stent with Delivery System	8.5 / 2.83	15
M00539300	Flexima Stent with Delivery System	10 / 3.33	5
M00539310	Flexima Stent with Delivery System	10 / 3.33	7
M00539320	Flexima Stent with Delivery System	10 / 3.33	10
M00539330	Flexima Stent with Delivery System	10 / 3.33	12
M00539340	Flexima Stent with Delivery System	10 / 3.33	15
M00539350	Flexima Stent with Delivery System	11.5 / 3.83	5
M00539360	Flexima Stent with Delivery System	11.5 / 3.83	7
M00539370	Flexima Stent with Delivery System	11.5 / 3.83	10
M00539380	Flexima Stent with Delivery System	11.5 / 3.83	12
M00539390	Flexima Stent with Delivery System	11.5 / 3.83	15

Packaged one per box

Recommended Guidewire: .035in (0.89mm)

Flexima™ Nasobiliary Catheter



Order Number	Description	O.D. (F) / (mm)	Catheter Length (cm)	Recommended Guidewire (in) / (mm)
M00540120†	Flexima Nasal Catheter	6 / 2.00	250	.035 / 0.89
M00540130†	Flexima Nasal Catheter	7.5 / 2.50	250	.035 / 0.89
M00540140†	Flexima Nasal Catheter	8.5 / 2.83	250	.035 / 0.89
M00540160†	Flexima Nasal Catheter	6 / 2.00	250	Guidewires included
M00540170†	Flexima Nasal Catheter	7.5 / 2.50	250	Guidewires included
M00540180†	Flexima Nasal Catheter	8.5 / 2.83	250	Guidewires included

Packaged one per box

Recommended Guidewire: .035in (0.89mm)

BILIARY PLASTIC STENTS (CONTINUED)



Percuflex™ Duodenal Bend Biliary Stents



Order Number	Description	O.D. (F) / (mm)	Length (cm) Barb to Barb
M00533610	Percuflex Duodenal Bend Stent	7 / 2.33	5
M00533620	Percuflex Duodenal Bend Stent	7 / 2.33	7
M00533630	Percuflex Duodenal Bend Stent	7 / 2.33	10
M00533640	Percuflex Duodenal Bend Stent	7 / 2.33	12
M00533650	Percuflex Duodenal Bend Stent	7 / 2.33	15
M00533660	Percuflex Duodenal Bend Stent	10 / 3.33	5
M00533670	Percuflex Duodenal Bend Stent	10 / 3.33	7
M00533680	Percuflex Duodenal Bend Stent	10 / 3.33	10
M00533690	Percuflex Duodenal Bend Stent	10 / 3.33	12
M00533700	Percuflex Duodenal Bend Stent	10 / 3.33	15

Packaged one per box

Recommended Guidewire: .035in (0.89mm)

Biliary Introducer*



Order Number	Description
M00533910	Push catheter for use with 7Fr (2.33mm) O.D. stent
M00533920	Push and guide catheter for use with 10Fr (3.33mm) O.D. stent

Packaged one per box

Recommended Guidewire: .035in (0.89mm)

* Stents sold separately.



C-flex™ Pigtail Biliary Stents*



Order Number	Description	O.D. (F) / (mm)	Length (cm) Barb to Barb	Colour Code
M00532010	C-flex Pigtail Stent	7 / 2.33	5	Blue
M00532020	C-flex Pigtail Stent	7 / 2.33	10	Blue
M00532030	C-flex Pigtail Stent	7 / 2.33	15	Blue
M00532040	C-flex Pigtail Stent	10 / 3.33	5	Grey
M00532050	C-flex Pigtail Stent	10 / 3.33	10	Grey
M00532060	C-flex Pigtail Stent	10 / 3.33	15	Grey


Packaged one per box

Recommended Guidewire: .035in (0.89mm)

* For use with Biliary Introducer.

Stents

PANCREATIC PLASTIC STENTS



Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Pigtail

Order Number	Description	Leading Barb	Shape	O.D. (F) / (mm)	Length (cm)
M00537350	Advanix Pancreatic Stent Kit ●	No	Pigtail	4 / 1.35	3
M00537360	Advanix Pancreatic Stent Kit ▲	No	Pigtail	4 / 1.35	4
M00537370	Advanix Pancreatic Stent Kit ▲	No	Pigtail	4 / 1.35	5
M00537380	Advanix Pancreatic Stent Kit ▲	No	Pigtail	4 / 1.35	7
M00537390	Advanix Pancreatic Stent Kit ●	No	Pigtail	5 / 1.67	3
M00537400	Advanix Pancreatic Stent Kit ▲	No	Pigtail	5 / 1.67	4
M00537410	Advanix Pancreatic Stent Kit ▲	No	Pigtail	5 / 1.67	5
M00537420	Advanix Pancreatic Stent Kit ▲	No	Pigtail	5 / 1.67	7
M00537430	Advanix Pancreatic Stent Kit ●	Yes	Pigtail	5 / 1.67	3
M00537440	Advanix Pancreatic Stent Kit ▲	Yes	Pigtail	5 / 1.67	4
M00537450	Advanix Pancreatic Stent Kit ▲	Yes	Pigtail	5 / 1.67	5

* Products with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Straight

Order Number	Description	Leading Barb	Shape	O.D. (F) / (mm)	Length (cm)
M00537460	Advanix Pancreatic Stent Kit ● ¹	No	Straight	4 / 1.35	3
M00537470	Advanix Pancreatic Stent Kit ▲ ¹	No	Straight	4 / 1.35	4
M00537480	Advanix Pancreatic Stent Kit ▲ ¹	No	Straight	4 / 1.35	5
M00537490	Advanix Pancreatic Stent Kit ▲ ¹	No	Straight	4 / 1.35	7
M00537500	Advanix Pancreatic Stent Kit ●	No	Straight	5 / 1.67	3
M00537510	Advanix Pancreatic Stent Kit ▲	No	Straight	5 / 1.67	4
M00537520	Advanix Pancreatic Stent Kit ▲	No	Straight	5 / 1.67	5
M00537530	Advanix Pancreatic Stent Kit ▲	No	Straight	5 / 1.67	7
M00537540	Advanix Pancreatic Stent Kit ●	Yes	Straight	5 / 1.67	3
M00537550	Advanix Pancreatic Stent Kit ▲	Yes	Straight	5 / 1.67	4
M00537560	Advanix Pancreatic Stent Kit ▲	Yes	Straight	5 / 1.67	5

* Products with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker

¹ Stiffer stent material compared to 4F pigtail stents.

² Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.



PANCREATIC PLASTIC STENTS (CONTINUED)



Advanix™ Pancreatic Single Stents – Pigtail



Order Number	Description	Leading Barb	Shape	O.D. (F) / (mm)	Length (cm)
M00536740	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	2
M00536750	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	3
M00536760	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	4
M00536770	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	5
M00536780	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	6
M00536790	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	7
M00536800	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	8
M00536810	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	9
M00536820	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	10
M00536830	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	11
M00536840	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	12
M00536850	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	13
M00536860	Advanix Pancreatic Stent ●	No	Pigtail	4 / 1.35	2
M00536870	Advanix Pancreatic Stent ●	No	Pigtail	4 / 1.35	3
M00536880	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	4
M00536890	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	5
M00536900	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	6
M00536910	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	7
M00536920	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	8
M00536930	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	9
M00536940	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	10
M00536950	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	11
M00536960	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	12
M00536970	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	13
M00536980	Advanix Pancreatic Stent ●	No	Pigtail	5 / 1.67	2
M00536990	Advanix Pancreatic Stent ●	No	Pigtail	5 / 1.67	3
M00537000	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	4
M00537010	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	5
M00537020	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	6
M00537030	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	7
M00537040	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	8
M00537050	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	9
M00537060	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	10
M00537070	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	11
M00537080	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	12
M00537090	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	13

* Products with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker

1 Stiffer stent material compared to 4F pigtail stents.

2 Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.



PANCREATIC PLASTIC STENTS (CONTINUED)



Advanix™ Pancreatic Single Stents – Pigtail (CONTINUED)

Order Number	Description	Leading Barb	Shape	O.D. (F) / (mm)	Length (cm)
M00537100	Advanix Pancreatic Stent ●	Yes	Pigtail	5 / 1.67	2
M00537110	Advanix Pancreatic Stent ●	Yes	Pigtail	5 / 1.67	3
M00537120	Advanix Pancreatic Stent ▲	Yes	Pigtail	5 / 1.67	4
M00537130	Advanix Pancreatic Stent ▲	Yes	Pigtail	5 / 1.67	5
M00537140	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	6
M00537150	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	7
M00537160	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	8
M00537170	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	9
M00537180	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	10
M00537190	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	11
M00537200	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	12
M00537210	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	13
M00537220	Advanix Pancreatic Stent ■	Yes	Pigtail	7 / 2.33	3
M00537230	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	4
M00537240	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	5
M00537250	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	6
M00537260	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	7
M00537270	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	8
M00537280	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	9
M00537290	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	10
M00537300	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	11
M00537310	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	12
M00537320	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	13
M00537330	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	15
M00537340	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	18

* Products with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker

¹ Stiffer stent material compared to 4F pigtail stents.

² Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.




Contents

Alphabetical Index

Product Code Index

PANCREATIC PLASTIC STENTS (CONTINUED)



Advanix™ Pancreatic Single Stents – Straight



Order Number	Description	Leading Barb	Shape	O.D. (F) / (mm)	Length (cm)
M00536000	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	2
M00536010	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	3
M00536020	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	4
M00536030	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	5
M00536040	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	6
M00536050	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	7
M00536060	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	8
M00536070	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	9
M00536080	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	10
M00536090	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	11
M00536100	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	12
M00536110	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	13
M00536120	Advanix Pancreatic Stent ● ¹	No	Straight	4 / 1.35	2
M00536130	Advanix Pancreatic Stent ● ¹	No	Straight	4 / 1.35	3
M00536140	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	4
M00536150	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	5
M00536160	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	6
M00536170	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	7
M00536180	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	8
M00536190	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	9
M00536200	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	10
M00536210	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	11
M00536220	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	12
M00536230	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	13
M00536240	Advanix Pancreatic Stent ●	No	Straight	5 / 1.67	2
M00536250	Advanix Pancreatic Stent ●	No	Straight	5 / 1.67	3
M00536260	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	4
M00536270	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	5
M00536280	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	6
M00536290	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	7
M00536300	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	8
M00536310	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	9
M00536320	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	10
M00536330	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	11
M00536340	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	12
M00536350	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	13

* Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.


■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker

¹ Stiffer stent material compared to 4F pigtail stents.

² Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.



PANCREATIC PLASTIC STENTS (CONTINUED)



Advanix™ Pancreatic Single Stents – Straight (CONTINUED)



Order Number	Description	Leading Barb	Shape	O.D. (F) / (mm)	Length (cm)
M00536360	Advanix Pancreatic Stent ●	Yes	Straight	5 / 1.67	2
M00536370	Advanix Pancreatic Stent ●	Yes	Straight	5 / 1.67	3
M00536380	Advanix Pancreatic Stent ▲	Yes	Straight	5 / 1.67	4
M00536390	Advanix Pancreatic Stent ▲	Yes	Straight	5 / 1.67	5
M00536400	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	6
M00536410	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	7
M00536420	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	8
M00536430	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	9
M00536440	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	10
M00536450	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	11
M00536460	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	12
M00536470	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	13
M00536480	Advanix Pancreatic Stent ■	Yes	Straight	7 / 2.33	3
M00536490	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	4
M00536500	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	5
M00536510	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	6
M00536520	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	7
M00536530	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	8
M00536540	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	9
M00536550	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	10
M00536560	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	11
M00536570	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	12
M00536580	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	13
M00536590	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	15
M00536600	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	18
M00536610	Advanix Pancreatic Stent ■	Yes	Straight	10 / 3.33	3
M00536620	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	4
M00536630	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	5
M00536640	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	6
M00536650	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	7
M00536660	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	8
M00536670	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	9
M00536680	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	10
M00536690	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	11
M00536700	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	12
M00536710	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	13
M00536720	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	15
M00536730	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	18

* Products with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.


■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker

¹ Stiffer stent material compared to 4F pigtail stents.

² Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.



RX BILIARY PLASTIC STENTS



NaviFlex™ RX Pancreatic Delivery System & NaviFlex™ Pushers



Order Number	Description	Diameter (F) / (mm)
M00535000	Pancreatic Stent RX Pusher	3 / 1.00
M00535010	Pancreatic Stent RX Pusher	4 & 5 / 1.35 & 1.67
M00535020	Pancreatic Stent Pusher	3 / 1.00
M00535030	Pancreatic Stent Pusher	4 & 5 / 1.35 & 1.67
M00535060	NaviFlex RX Pancreatic Stent Delivery System	7 / 2.33
M00535080	NaviFlex RX Pancreatic Stent Delivery System	10 / 3.33

ESOPHAGEAL STENTS



WallFlex™ Fully Covered RMV Esophageal Stent System

Order Number	Description	Stent O.D. (mm)	Proximal / Distal Flares O.D. (mm)	Stent Length (cm)	Catheter Diameter (F) / (mm)	Working Length (cm)
M00516210†	WallFlex Fully Covered RMV Esophageal Stent System	18	25 / 23	10.3	18.5 / 6.2	78
M00516220†	WallFlex Fully Covered RMV Esophageal Stent System	18	25 / 23	12.3	18.5 / 6.2	78
M00516230†	WallFlex Fully Covered RMV Esophageal Stent System	18	25 / 23	15.3	18.5 / 6.2	78
M00516240†	WallFlex Fully Covered RMV Esophageal Stent System	23	28 / 28	10.5	18.5 / 6.2	78
M00516250†	WallFlex Fully Covered RMV Esophageal Stent System	23	28 / 28	12.5	18.5 / 6.2	78
M00516260†	WallFlex Fully Covered RMV Esophageal Stent System	23	28 / 28	15.5	18.5 / 6.2	78

Packaged one per box

† Product commercialised only in countries where refractory benign indication is approved.

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621

WallFlex™ Fully Covered RMV Esophageal Longer Loop Stent System

Order Number	Description	Stent O.D. (mm)	Proximal / Distal Flares O.D. (mm)	Stent Length (cm)	Catheter Diameter (F) / (mm)	Working Length (cm)
M00516270	Wallflex Fully Covered RMV Esophageal Longer Loop Stent System	18	25 / 23	10.3	18.5 / 6.2	78
M00516280	Wallflex Fully Covered RMV Esophageal Longer Loop Stent System	18	25 / 23	12.3	18.5 / 6.2	78
M00516290	Wallflex Fully Covered RMV Esophageal Longer Loop Stent System	18	25 / 23	15.3	18.5 / 6.2	78
M00516300	Wallflex Fully Covered RMV Esophageal Longer Loop Stent System	23	28 / 28	10.5	18.5 / 6.2	78
M00516310	Wallflex Fully Covered RMV Esophageal Longer Loop Stent System	23	28 / 28	12.5	18.5 / 6.2	78
M00516320	Wallflex Fully Covered RMV Esophageal Longer Loop Stent System	23	28 / 28	15.5	18.5 / 6.2	78

Packaged one per box

Limited availability in all European countries (except UK and France). Please contact your local Territory Manager.

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621




WallFlex™ Partially Covered Esophageal Stent System

Order Number	Description	Stent O.D. (mm)	Proximal / Distal Flares O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Catheter Diameter (F) / (mm)	System Length (cm)
M00516900	Wallflex Partially Covered Esophageal Stent System	18	23 / 23	10.3	7.1	18.5 / 6.2	120
M00516910	Wallflex Partially Covered Esophageal Stent System	18	23 / 23	12.3	9	18.5 / 6.2	120
M00516920	Wallflex Partially Covered Esophageal Stent System	18	23 / 23	15.3	12.1	18.5 / 6.2	120
M00516930	Wallflex Partially Covered Esophageal Stent System	23	28 / 28	10.5	7.3	18.5 / 6.2	120
M00516940	Wallflex Partially Covered Esophageal Stent System	23	28 / 28	12.5	9.3	18.5 / 6.2	120
M00516950	Wallflex Partially Covered Esophageal Stent System	23	28 / 28	15.5	12.3	18.5 / 6.2	120

Packaged one per box

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621

ESOPHAGEAL STENTS (CONTINUED)


**Ultraflex™ Covered Esophageal NG Stent System****Proximal Release**

Order Number	Description	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)
M00513840	Ultraflex Esophageal NG Stent System	18	23	10	7
M00513850	Ultraflex Esophageal NG Stent System	18	23	12	9
M00513860	Ultraflex Esophageal NG Stent System	18	23	15	12
M00514240	Ultraflex Esophageal NG Stent System	23	28	10	7
M00514250	Ultraflex Esophageal NG Stent System	23	28	12	9

Distal Release

Order Number	Description	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)
M00513730	Ultraflex Esophageal NG Stent System	18	23	10	7
M00513740	Ultraflex Esophageal NG Stent System	18	23	12	9
M00513750	Ultraflex Esophageal NG Stent System	18	23	15	12
M00514200	Ultraflex Esophageal NG Stent System	23	28	10	7
M00514210	Ultraflex Esophageal NG Stent System	23	28	12	9

Packaged one per box

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621**Ultraflex™ Uncovered Esophageal NG Stent System****Proximal Release**

Order Number	Description	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)
M00513800	Ultraflex Esophageal NG Stent System	18	23	7
M00513810	Ultraflex Esophageal NG Stent System	18	23	10
M00513820	Ultraflex Esophageal NG Stent System	18	23	12
M00513830	Ultraflex Esophageal NG Stent System	18	23	15
M00514220	Ultraflex Esophageal NG Stent System	23	28	7
M00514230	Ultraflex Esophageal NG Stent System	23	28	12



Distal Release

Order Number	Description	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)
M00513700	Ultraflex Esophageal NG Stent System	18	23	7
M00513710	Ultraflex Esophageal NG Stent System	18	23	10
M00513720	Ultraflex Esophageal NG Stent System	18	28	15

Packaged one per box

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621

COLONIC AND DUODENAL STENTS



WallFlex™ Colonic Stent System with Anchor Lock Delivery System



Order Number	Description	Stent Body Diameter (mm)	Stent Flare Diameter (mm)	Expanded Stent Length (cm)	Working Length (cm)	Catheter Diameter (F) / (mm)
M00565100	WallFlex Colonic Stent	22	27	6	230	10 / 3.33
M00565110	WallFlex Colonic Stent	22	27	9	230	10 / 3.33
M00565120	WallFlex Colonic Stent	22	27	12	230	10 / 3.33
M00565130	WallFlex Colonic Stent	22	27	6	135	10 / 3.33
M00565140	WallFlex Colonic Stent	22	27	9	135	10 / 3.33
M00565150	WallFlex Colonic Stent	22	27	12	135	10 / 3.33
M00565040	WallFlex Colonic Stent	25	30	6	230	10 / 3.33
M00565050	WallFlex Colonic Stent	25	30	9	230	10 / 3.33
M00565060	WallFlex Colonic Stent	25	30	12	230	10 / 3.33
M00565070	WallFlex Colonic Stent	25	30	6	135	10 / 3.33
M00565080	WallFlex Colonic Stent	25	30	9	135	10 / 3.33
M00565090	WallFlex Colonic Stent	25	30	12	135	10 / 3.33

Packaged one per box

Recommended Guidewire Systems 230cm: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 450cm: M00556161 Straight Dream Tip™, M00556171 Angled Dream Tip™ or .035 in (0.89mm) Jagwire™ Guidewire Stiff Shaft, 450cm: M00556601 Straight Tip, M00556611 Angled Tip or .035in (0.89mm) Super Stiff Wire, 500cm: H954180011



Systems 135cm: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm: M00556121 Straight Tip, M00556131 Angled Tip or .035in (0.89mm) Jagwire™ Guidewire Stiff Shaft, 260cm: M00557661 Straight Tip

NEW WallFlex™ Colonic Soft Stent System with Anchor Lock Delivery System



Order Number	Description	Diameter Flare/Body (mm)	Stent Length (cm)	Working Length (cm)	Delivery System Diameter/ Working Channel (F/mm)	Guidewire Diameter (inch/mm)
M00555700	WallFlex Colonic Soft Stent	23 / 18	6	230	9 / 3.0	0.035 / 0.89
M00555710	WallFlex Colonic Soft Stent	23 / 18	9	230	9 / 3.0	0.035 / 0.89
M00555720	WallFlex Colonic Soft Stent	23 / 18	12	230	9 / 3.0	0.035 / 0.89
M00555730	WallFlex Colonic Soft Stent	25 / 20	6	230	9 / 3.0	0.035 / 0.89
M00555740	WallFlex Colonic Soft Stent	25 / 20	9	230	9 / 3.0	0.035 / 0.89
M00555750	WallFlex Colonic Soft Stent	25 / 20	12	230	9 / 3.0	0.035 / 0.89
M00555760	WallFlex Colonic Soft Stent	27 / 22	6	230	9 / 3.0	0.035 / 0.89
M00555770	WallFlex Colonic Soft Stent	27 / 22	9	230	9 / 3.0	0.035 / 0.89
M00555780	WallFlex Colonic Soft Stent	27 / 22	12	230	9 / 3.0	0.035 / 0.89

COLONIC AND DUODENAL STENTS



WallFlex™ Duodenal Stent System with Anchor Lock Delivery System



Order Number	Description	Stent Body Diameter (mm)	Stent Flare Diameter (mm)	Expanded Stent Length (cm)	Working Length (cm)	Catheter Diameter (F) / (mm)
M00565010	WallFlex Duodenal Stent	22	27	6	230	10 / 3.33
M00565020	WallFlex Duodenal Stent	22	27	9	230	10 / 3.33
M00565030	WallFlex Duodenal Stent	22	27	12	230	10 / 3.33

Packaged one per box


Recommended Guidewires: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 450cm: M00556161 Straight Dream Tip™, M00556171 Angled Dream Tip™ or .035 in (0.89mm) Jagwire™ Guidewire Stiff Shaft, 450cm: M00556601 Straight Tip, M00556611 Angled Tip or .035in (0.89mm) Super Stiff Wire, 500cm: H954180011

NEW WallFlex™ Duodenal Soft Stent System with Anchor Lock Delivery System



Order Number	Description	Diameter Flare/Body (mm)	Stent Length (cm)	Working Length (cm)	Delivery System Diameter/ Working Channel (F/mm)	Guidewire Diameter (inch/mm)
M00555800	WallFlex Duodenal Soft Stent	23 / 18	6	230	9 / 3.0	0.035 / 0.89
M00555810	WallFlex Duodenal Soft Stent	23 / 18	9	230	9 / 3.0	0.035 / 0.89
M00555820	WallFlex Duodenal Soft Stent	23 / 18	12	230	9 / 3.0	0.035 / 0.89
M00555830	WallFlex Duodenal Soft Stent	25 / 20	6	230	9 / 3.0	0.035 / 0.89
M00555840	WallFlex Duodenal Soft Stent	25 / 20	9	230	9 / 3.0	0.035 / 0.89
M00555850	WallFlex Duodenal Soft Stent	25 / 20	12	230	9 / 3.0	0.035 / 0.89
M00555860	WallFlex Duodenal Soft Stent	27 / 22	6	230	9 / 3.0	0.035 / 0.89
M00555870	WallFlex Duodenal Soft Stent	27 / 22	9	230	9 / 3.0	0.035 / 0.89
M00555880	WallFlex Duodenal Soft Stent	27 / 22	12	230	9 / 3.0	0.035 / 0.89

COLONIC AND DUODENAL STENTS (CONTINUED)



**Wallstent™ Enteral Endoprosthesis Colonic/Duodenal Stent
with UNISTEP™ Plus Delivery System**


Order Number	Description	Expanded Stent Diameter (mm)	Expanded Stent Length (cm)	Working Length (cm)	Catheter Diameter (F) / (mm)
M00565560	Wallstent Enteral Endoprosthesis Colonic/Duodenal Stent	20	6	230	10 / 3.33
M00565570	Wallstent Enteral Endoprosthesis Colonic/Duodenal Stent	20	9	230	10 / 3.33
M00565620	Wallstent Enteral Endoprosthesis Colonic/Duodenal Stent	20	6	135	10 / 3.33
M00565630	Wallstent Enteral Endoprosthesis Colonic/Duodenal Stent	20	9	135	10 / 3.33
M00565580	Wallstent Enteral Endoprosthesis Colonic/Duodenal Stent	22	6	230	10 / 3.33
M00565590	Wallstent Enteral Endoprosthesis Colonic/Duodenal Stent	22	9	230	10 / 3.33
M00565640	Wallstent Enteral Endoprosthesis Colonic/Duodenal Stent	22	6	135	10 / 3.33
M00565650	Wallstent Enteral Endoprosthesis Colonic/Duodenal Stent	22	9	135	10 / 3.33

Packaged one per box

Recommended Guidewire Systems 230cm: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 450cm: M00556161 Straight Dream Tip™, M00556171 Angled Dream Tip™ or .035in (0.89mm) Jagwire™ Guidewire Stiff Shaft, 450cm: M00556601 Straight Tip, M00556611 Angled Tip or .035in (0.89mm) Super Stiff Wire, 500cm: H954180011

Recommended Guidewire Systems 135cm: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm: M00556121 Straight Tip, M00556131 Angled Tip or .035in (0.89mm) Jagwire™ Guidewire Stiff Shaft, 260cm: M00557661 Straight Tip


**Ultraflex Precision™ Colonic Stent System
Proximal Release**


Order Number	Description	Stent Body Diameter (mm)	Stent Flare Diameter (mm)	Expanded Stent Length (cm)	Working Length (cm)	Catheter Diameter (F) / (mm)
M00557330	Ultraflex Precision Colonic Stent	25	30	5.7	100	22 / 7.33
M00557360	Ultraflex Precision Colonic Stent	25	30	8.7	100	22 / 7.33
M00557380	Ultraflex Precision Colonic Stent	25	30	11.7	100	22 / 7.33

Packaged one per box

Recommended Guidewire: .038in (0.97mm) Jagwire™ Guidewire, Order No. M00556621



RX Biliary
System™

Biliary

SpyGlass™

EUS

Stents

Dilation

Guidewires

Hemostasis

Tissue
Resection

Biopsy

Enteral
Feeding

Pulmonary

Infection
Prevention

RFA




Contents




Alphabetical
Index

Product Code
Index



CRE™ (CONTROLLED RADIAL EXPANSION) BALLOON DILATATION CATHETERS



CRE™ PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilatation Catheters (including DASE indication)




Order Number	BALLOON INFLATED O.D. (mm)	BALLOON INFLATED O.D. (Fr)	Balloon Length (cm)	INFLATION PRESSURE (ATM)	INFLATION PRESSURE (kPa)	Minimum Working Channel (mm)*	Catheter Diameter (F) / (mm)	Working Length (cm)	Units
M00558600	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8/3.2"	7.5 / 2.50	180	box 1
M00558610	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8/3.2"	7.5 / 2.50	180	box 1
M00558620	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8/3.2"	7.5 / 2.50	180	box 1
M00558630	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8/3.2"	7.5 / 2.50	180	box 1
M00558640	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8/3.7"	7.5 / 2.50	180	box 1
M00558650**	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8/3.7"	7.5 / 2.50	180	box 1

Recommended Inflation Device: Alliance™ II Integrated Inflation System Re-usable Handle M00550620, Sterile Syringe M00550601
 Guidewire included and preloaded in balloon catheter.

* Duodenoscope Channel Size

** Esophageal and biliary applications only

This product contains no detectable latex.



CRE™ PRO Wireguided Esophageal/Pyloric/Colonic/Biliary Balloon Dilatation Catheters (including DASE indication)



Order Number	BALLOON INFLATED O.D. (mm)	BALLOON INFLATED O.D. (Fr)	Balloon Length (cm)	INFLATION PRESSURE (ATM)	INFLATION PRESSURE (kPa)	Minimum Working Channel (mm)*	Catheter Diameter (F) / (mm)	Working Length (cm)	Units
M00558660	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8/3.2"	7.5 / 2.50	240	box 1
M00558670	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8/3.2"	7.5 / 2.50	240	box 1
M00558680	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8/3.2"	7.5 / 2.50	240	box 1
M00558690	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8/3.2"	7.5 / 2.50	240	box 1
M00558700	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8/3.7"	7.5 / 2.50	240	box 1
M00558710**	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8/3.7"	7.5 / 2.50	240	box 1

Recommended Inflation Device: Alliance™ II Integrated Inflation System Re-usable Handle M00550620, Sterile Syringe M00550601
 Guidewire included and preloaded in balloon catheter.


* Duodenoscope Channel Size

** Esophageal, colonic and biliary applications only

This product contains no detectable latex.

CRE™ (CONTROLLED RADIAL EXPANSION) BALLOON DILATATION CATHETERS

(CONTINUED)



CRE™ Fixed Wire Balloon Dilatation Catheters




Order Number	Inflated O.D. (mm)	Balloon Length (cm)	INFLATION PRESSURE (ATM)	INFLATION PRESSURE (kPa)	Minimum Working Channel (mm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Units
M00558330	6-7-8	8	3-6-10	304-608-1013	2.8	6 / 2.00	180	box 1
M00558340	8-9-10	8	3-5.5-9	304-557-912	2.8	6 / 2.00	180	box 1
M00558350	10-11-12	8	3-5-8	304-507-810	2.8	6 / 2.00	180	box 1
M00558360	12-13.5-15	8	3-4.5-8	304-456-810	2.8	6 / 2.00	180	box 1
M00558361	12-13.5-15	8	3-4.5-8	304-456-810	2.8	6 / 2.00	180	box 5
M00558370	15-16.5-18	8	3-4.5-7	304-456-709	2.8	6 / 2.00	180	box 1
M00558371	15-16.5-18	8	3-4.5-7	304-456-709	2.8	6 / 2.00	180	box 5
M00558380	18-19-20	8	3-4.5-6	304-456-608	2.8	6 / 2.00	180	box 1
M00558381	18-19-20	8	3-4.5-6	304-456-608	2.8	6 / 2.00	180	box 5

Recommended Inflation Device: Alliance™ II Integrated Inflation System Re-usable Handle M00550620, Sterile Syringe M00550601

This product contains no detectable latex.

MAXFORCE TTS™ BALLOON DILATORS



MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter



Order Number	Inflated O.D. (mm)	Balloon Length (cm)	INFLATION PRESSURE (ATM)	INFLATION PRESSURE (kPa)	Minimum Working Channel (mm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Units
M00589060	6	6	10	1013	2.8	6 / 2.00	180	box 1
M00589080	8	6	10	1013	2.8	6 / 2.00	180	box 1
M00589100	10	6	7	709	2.8	6 / 2.00	180	box 1
M00589120	12	6	7	709	2.8	6 / 2.00	180	box 1
M00589140	14	6	7	709	2.8	6 / 2.00	180	box 1
M00589150	15	6	5	507	2.8	6 / 2.00	180	box 1
M00589160	16	6	5	507	2.8	6 / 2.00	180	box 1
M00589180	18	6	5	507	2.8	6 / 2.00	180	box 1

Recommended Inflation Device: Alliance™ II Integrated Inflation System Re-usable Handle M00550620, Sterile Syringe M00550601
This product contains no detectable latex.

Rigiflex™ II Achalasia Balloon Dilators*




Order Number	Inflated O.D. (mm)	Balloon Length (cm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Recommended Guidewire (in) / (mm)**	Units
M00554500	30	10	14 / 4.67	90	.038 / 0.97	box 1
M00554510	35	10	14 / 4.67	90	.038 / 0.97	box 1
M00554520	40	10	14 / 4.67	90	.038 / 0.97	box 1

* Product shipped sterile. Recommended use with pneumatic inflator, M00553200

** Guidewire included with Rigiflex™ II Balloon Dilators

This product contains no detectable latex.

INFLATION SYSTEMS



Alliance™ II Inflation System

Order Number	Description	Units
M00550620	Alliance II Integrated Inflation/Lithotripsy Device	box 1
M00550601	Alliance II Syringe/Gauge Assembly	box 5


DILATION ACCESSORIES

Inflation/Deflation Devices and Accessories

Order Number	Description	Units
M00553200	Pneumatic Inflator	box 1
M00550630	CRE™ SteriFlate Disposable Inflation Device	box 1

RX Biliary Balloon Dilators are located in the RX Biliary Section on page 45

Biliary Balloon Dilators are located in the Biliary Section on page 52



RX Biliary System™	Biliary	SpyGlass™	EUS	Stents	Dilation	Guidewires	Hemostasis	Tissue Resection	Biopsy	Enteral Feeding	Pulmonary	Infection Prevention	RFA
--------------------	---------	-----------	-----	--------	----------	------------	------------	------------------	--------	-----------------	-----------	----------------------	-----



Contents



Alphabetical Index

Product Code Index

HIGH PERFORMANCE




Dreamwire™ High Performance Guidewires




Order Number	Description	O.D. (in) / (mm)	Length (cm)	Dream Tip Style	Units
M00556101	Dreamwire Standard	.035 / 0.89	260	Straight	box 2
M00556111	Dreamwire Standard	.035 / 0.89	260	Angled	box 2
M00556121	Dreamwire Stiff Shaft	.035 / 0.89	260	Straight	box 2
M00556131	Dreamwire Stiff Shaft	.035 / 0.89	260	Angled	box 2
M00556141	Dreamwire Standard	.035 / 0.89	450	Straight	box 2
M00556151	Dreamwire Standard	.035 / 0.89	450	Angled	box 2
M00556161	Dreamwire Stiff Shaft	.035 / 0.89	450	Straight	box 2
M00556171	Dreamwire Stiff Shaft	.035 / 0.89	450	Angled	box 2


5cm Jagwire™ Guidewire Tip



10 cm Straight Dream Tip™



10 cm Angled Dream Tip™




Hydra Jagwire™ High Performance Guidewires



Order Number	Description	O.D. (in) / (mm)	Length (cm)	Dream Tip Style	Units
M00556001	Hydra Jagwire Standard	.035 / 0.89	260	Straight	box 2
M00556011	Hydra Jagwire Standard	.035 / 0.89	260	Angled	box 2
M00556021	Hydra Jagwire Stiff Shaft	.035 / 0.89	260	Straight	box 2
M00556031	Hydra Jagwire Stiff Shaft	.035 / 0.89	260	Angled	box 2
M00556041	Hydra Jagwire Standard	.035 / 0.89	450	Straight	box 2
M00556051	Hydra Jagwire Standard	.035 / 0.89	450	Angled	box 2
M00556061	Hydra Jagwire Stiff Shaft	.035 / 0.89	450	Straight	box 2
M00556071	Hydra Jagwire Stiff Shaft	.035 / 0.89	450	Angled	box 2

Hydra Jagwire Guidewire may remain in place during sphincterotomy. See "Instructions for Use" for details.

HIGH PERFORMANCE (CONTINUED)



Jagwire™ High Performance Guidewires

Order Number	Description	O.D. (in) / (mm)	Length (cm)	Tip Style	Units
M00556481	Jagwire Stiff Shaft	.025-.019 / 0.64-0.48	260	Straight	box 2
M00556491	Jagwire Stiff Shaft	.025-.019 / 0.64-0.48	260	Angled	box 2
M00556461	Jagwire Standard	.025-.018 / 0.64-0.46	260	Straight	box 2
M00557641	Jagwire Standard	.035-.027 / 0.89-0.69	260	Straight	box 2
M00557651	Jagwire Standard	.035-.027 / 0.89-0.69	260	Angled	box 2
M00557661	Jagwire Stiff Shaft	.035-.027 / 0.89-0.69	260	Straight	box 2
M00557671	Jagwire Stiff Shaft	.035-.027 / 0.89-0.69	260	Angled	box 2
M00556561	Jagwire Standard	.025-.018 / 0.64-0.46	450	Straight	box 2
M00556511	Jagwire Standard	.025-.018 / 0.64-0.46	450	Angled	box 2
M00556521	Jagwire Stiff Shaft	.025-.019 / 0.64-0.48	450	Straight	box 2
M00556531	Jagwire Stiff Shaft	.025-.019 / 0.64-0.48	450	Angled	box 2
M00556581	Jagwire Standard	.035-.027 / 0.89-0.69	450	Straight	box 2
M00556591	Jagwire Standard	.035-.027 / 0.89-0.69	450	Angled	box 2
M00556601	Jagwire Stiff Shaft	.035-.027 / 0.89-0.69	450	Straight	box 2
M00556611	Jagwire Stiff Shaft	.035-.027 / 0.89-0.69	450	Angled	box 2
M00556621	Jagwire Stiff Shaft	.038 / 0.97	260	Straight	box 2

Jagwire Guidewire may remain in place during sphincterotomy. See "Instructions for Use" for details.

Jagwire™ High Performance Guidewires – Extendable Jagwire

Order Number	Description	O.D. (in) / (mm)	Length (cm)	Tip Style	Units
M00556902	Jagwire Jag Tail	.035 / 0.89	200	Guidewire Extension	box 5
M00556841	Jagwire Standard	.035-.027 / 0.89-0.69	260	Straight	box 2
M00556851	Jagwire Standard	.035-.027 / 0.89-0.69	260	Angled	box 2
M00556861	Jagwire Stiff Shaft	.035-.027 / 0.89-0.69	260	Straight	box 2
M00556871	Jagwire Stiff Shaft	.035-.027 / 0.89-0.69	260	Angled	box 2


Jagwire Guidewire may remain in place during sphincterotomy. See "Instructions for Use" for details.

NEW Jagwire™ Revolution High Performance Guidewire

Order Number	Description	OD (in)	Length (cm)	Tip Style	Units
M00557001	Jagwire Revolution High Performance Guidewire	.025	260	Straight Round	box 2
M00557011	Jagwire Revolution High Performance Guidewire	.025	260	Angled Round	box 2
M00557021	Jagwire Revolution High Performance Guidewire	.025	450	Straight Round	box 2
M00557031	Jagwire Revolution High Performance Guidewire	.025	450	Angled Round	box 2
M00557041	Jagwire Revolution High Performance Guidewire	.025	500	Straight Round	box 2


Guidewires

HIGH PERFORMANCE (CONTINUED)



NovaGold™ High Performance Guidewires

Order Number	Description	Guidewire Diameter (in) / (mm)	Length (cm)	Tip Style	Units
M00552001	NovaGold High Performance Guidewire	.018 / 0.46	260	Straight (Shapeable)	box 2
M00552011	NovaGold High Performance Guidewire	.018 / 0.46	480	Straight (Shapeable)	box 2



NaviPro™ Hydrophilic Guidewires

Order Number	Description	O.D. (in) / (mm)	Length (cm)	Tip Style	Units
M00556191	NaviPro Standard	.018 / 0.46	260	Angled	box 5
M00556181	NaviPro Standard	.018 / 0.46	260	Straight	box 5
M00556211	NaviPro Standard	.025 / 0.64	260	Angled	box 5
M00556201	NaviPro Standard	.025 / 0.64	260	Straight	box 5
M00556231	NaviPro Standard	.035 / 0.89	260	Angled	box 5
M00556221	NaviPro Standard	.035 / 0.89	260	Straight	box 5
M00556251	NaviPro Stiff Shaft	.035 / 0.89	260	Angled	box 5
M00556241	NaviPro Stiff Shaft	.035 / 0.89	260	Straight	box 5

STANDARD

Amplatz Super Stiff™ Guidewire

Order Number	Description	O.D. (in) / (mm)	Length (cm)	Tip Style	Units
M00550090	Amplatz Super Stiff Guidewire	.038 / 0.97	260	Straight	box 1

Wallstent™ Super Stiff Guidewire

Order Number	Description	O.D. (in) / (mm)	Length (cm)	Tip Style	Units
H965180011	Wallstent Super Stiff Guidewire	.035 / 0.89	500	Tapered	box 5



RX Biliary System™	Biliary	SpyGlass™	EUS	Stents	Dilation	Guidewires	Hemostasis	Tissue Resection	Biopsy	Enteral Feeding	Pulmonary	Infection Prevention	RFA
--------------------	---------	-----------	-----	--------	----------	------------	------------	------------------	--------	-----------------	-----------	----------------------	-----




Contents



Alphabetical Index

Product Code Index

CLIPPING DEVICES




Resolution 360™ Clips

Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Units
M00521230	Resolution 360 Clip	235	2.8	11	box 1
M00521231	Resolution 360 Clip	235	2.8	11	box 10
M00521232†	Resolution 360 Clip	235	2.8	11	box 20
M00521233†	Resolution 360 Clip	235	2.8	11	box 40
M00521240	Resolution 360 Clip	155	2.8	11	box 1
M00521241	Resolution 360 Clip	155	2.8	11	box 10
M00521242	Resolution 360 Clip	155	2.8	11	box 20
M00521243	Resolution 360 Clip	155	2.8	11	box 40

* MR Conditional. Non-clinical testing has demonstrated the Resolution 360™ Clip is Magnetic Resonance (MR) Conditional according to ASTM F2503.

Refer to the directions for use for MR conditions of use.

† Please check product availability in your country with your sales representative




Resolution™ Clips

Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Units
M00522600	Resolution Clip	155	2.8	11	box 1
M00522601	Resolution Clip	155	2.8	11	box 10
M00522602	Resolution Clip	155	2.8	11	box 20
M00522610	Resolution Clip	235	2.8	11	box 1
M00522611	Resolution Clip	235	2.8	11	box 10
M00522612	Resolution Clip	235	2.8	11	box 20

* MR Conditional. Non-clinical testing has demonstrated the Resolution 360™ Clip is Magnetic Resonance (MR) Conditional according to ASTM F2503.

Refer to the directions for use for MR conditions of use.




Speedband Superview Super 7™ Multiple Band Ligators

Order Number	Description	Minimum Working Channel (mm)	Scope O.D. (mm)	Bands / Ligating Unit	Units
M00542250	Super 7 Band Ligator	2.8	8.6-11.5	7	box 1
M00542251	Super 7 Band Ligator	2.8	8.6-11.5	7	box 2
M00542253	Super 7 Band Ligator	2.8	8.6-11.5	7	box 4

Irrigation catheter supplied in each kit


BIPOLAR ELECTROHEMOSTASIS CATHETERS



Injection Gold Probe™ Electrohemostasis/Injection Catheters

Order Number	Description	O.D. (F) / (mm)	Gauge	NEEDLE O.D. (mm)	I.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Units
M00560150	Injection Gold Probe Catheter	7 / 2.33	25	0.51	0.24	2.8	210	box 1
M00560160	Injection Gold Probe Catheter	10 / 3.33	25	0.51	0.24	3.7	210	box 1

Type Generator/Connector: Please refer to Electrosurgical Bipolar Cable Adaptors




Gold Probe™ Electrohemostasis Catheters

Order Number	Description	O.D. (F) / (mm)	Minimum Working Channel (mm)	Working Length (cm)	Units
M00560070	Gold Probe Catheter	7 / 2.33	2.8	300	box 1
M00560071	Gold Probe Catheter	7 / 2.33	2.8	300	box 5
M00560100	Gold Probe Catheter	10 / 3.33	3.7	300	box 1
M00560101	Gold Probe Catheter	10 / 3.33	3.7	300	box 5
M00560220	Gold Probe Catheter	7 / 2.33	2.8	350	box 1

Type Generator/Connector: Please refer to Electrosurgical Bipolar Cable Adaptors

NEEDLES




Interject™ Sclerotherapy Needles*



Order Number	Description	Sheath Design	Sheath O.D. (F) / (mm)	Gauge	NEEDLE O.D. (mm)	I.D. (mm)	Maximum Needle Extension Length (mm)	Working Length (cm)	Units
M00518151	Interject Sclerotherapy Needle	Contrast	1.8	23	0.64	0.32	4	200	box 5
M00518161	Interject Sclerotherapy Needle	Contrast	1.8	25	0.51	0.24	4	200	box 5
M00518251	Interject Sclerotherapy Needle	Contrast	1.8	23	0.64	0.32	6	200	box 5
M00518261	Interject Sclerotherapy Needle	Contrast	1.8	25	0.51	0.24	6	200	box 5
M00518351	Interject Sclerotherapy Needle	Contrast	2.3	23	0.64	0.32	4	240	box 5
M00518361	Interject Sclerotherapy Needle	Contrast	2.3	25	0.51	0.24	4	240	box 5

* With Star Catheter Technology




Interject™ Clear Sclerotherapy Needles*



Order Number	Description	Sheath Design	Sheath O.D. (F) / (mm)	Gauge	NEEDLE O.D. (mm)	I.D. (mm)	Maximum Needle Extension Length (mm)	Working Length (cm)	Units
M00518111	Interject Sclerotherapy Needle	Clear	1.8	25	0.51	0.24	4	200	box 5
M00518301	Interject Sclerotherapy Needle	Clear	2.3	23	0.64	0.32	4	240	box 5
M00518311	Interject Sclerotherapy Needle	Clear	2.3	25	0.51	0.24	4	240	box 5

* With Star Catheter Technology

ELECTROSURGICAL BIPOLAR CABLE ADAPTOR




Order Number	Description	Fix Pin Spacing (mm)	Generator Compatibility*	Units
M00561440	Bipolar Cable Adaptor	18	EMC (Coagulasem Dolley)	box 1
M00561450	Bipolar Cable Adaptor	22	ERBE, EMC Dolley, Olympus (ESG100, PSD60)**	box 1
M00561460	Bipolar Cable Adaptor	25	Lamidey (Optima 4)	box 1
M00561470	Bipolar Cable Adaptor	28	Valleylab®, Lamidey (Duo), Bovie (Aaron 3250™), Olympus (UES10, UES30, Surgemaster), Bowa, Soring	box 1
M00561480	Bipolar Cable Adaptor	28 with Third Pin	ConMed® (System 5000), Bovie (Aaron 2250™), Ethicon®	box 1
M00561490	Bipolar Cable Adaptor	32 with Third Pin	ConMed (Sabre 180, Excalibur, System 7550)	box 1
M00561500	Bipolar Cable Adaptor	39	Codman, Malis®	box 1

Bipolar Cable Adaptors are used to connect Gold Probe™ and Injection Gold Probe™ with the above mentioned generators.
* For generators and their associated cable adaptors not listed above, please contact your Boston Scientific Representative.
** Some variations in faceplates may exist within a generator model. We recommend that you confirm pin spacing prior to ordering.



Tissue Resection

EMR DEVICES FOR UPPER GI




Captivator™ EMR



Order Number	Description	Units
M00561600	Captivator EMR for Standard Gastrosopes	box 1
M00561601	Captivator EMR for Standard Gastrosopes	box 5
M00561610	Captivator EMR for Large Gastrosopes	box 1
M00561611	Captivator EMR for Large Gastrosopes	box 5

SNARES



Rotatable Snares*

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00561821	Small Oval	Medium Stiff	13	2.4	240	2.8	box 10
M00561831	Medium Oval	Medium Stiff	20	2.4	240	2.8	box 10

* With Star Catheter Technology.




Contents

Alphabetical Index

Product Code Index


Tissue Resection

SNARES (CONTINUED)



Sensation™ Short Throw Snares

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562691	Small Oval	Flexible	13	2.4	240	2.8	box 5
M00562692	Small Oval	Flexible	13	2.4	240	2.8	box 20
M00562693	Small Oval	Flexible	13	2.4	240	2.8	box 40
M00560321	Small Oval	Medium Stiff	13	2.4	240	2.8	box 5
M00562671	Medium Oval	Flexible	27	2.4	240	2.8	box 5
M00562672	Medium Oval	Flexible	27	2.4	240	2.8	box 20
M00562673	Medium Oval	Flexible	27	2.4	240	2.8	box 40
M00560311	Medium Oval	Medium Stiff	27	2.4	240	2.8	box 5
M00562651	Large Oval	Medium Stiff	30	2.4	240	2.8	box 5
M00562652	Large Oval	Medium Stiff	30	2.4	240	2.8	box 20
M00560331	Medium Crescent	Medium Stiff	27	2.4	240	2.8	box 5




Captiflex™ Snares

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562471	Extra Small Oval	Flexible	11	2.4	240	2.8	box 10
M00562421	Small Oval	Flexible	13	2.4	240	2.8	box 10
M00562422	Small Oval	Flexible	13	2.4	240	2.8	box 40
M00562401	Medium Oval	Flexible	27	2.4	240	2.8	box 10
M00562402	Medium Oval	Flexible	27	2.4	240	2.8	box 40


Tissue Resection

SNARES (CONTINUED)



Captivator™ Snares

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562301	Small Oval	Stiff	13	2.4	240	2.8	box 10
M00562321	Medium Oval	Stiff	27	2.4	240	2.8	box 10
M00562391	Large Oval (thin wire)	Flexible	30	2.4	240	2.8	box 10
M00561311	Large Oval	Medium Stiff	30	2.4	240	2.8	box 10
M00562451	Small Hexagonal	Stiff	13	2.4	240	2.8	box 10
M00562341	Medium Hexagonal	Stiff	27	2.4	240	2.8	box 10
M00562371	Medium Crescent	Stiff	27	2.4	240	2.8	box 10



Captivator™ II Snares

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00561221	10 mm Rounded	Stiff	10	2.4	240	2.8	box 10
M00561222	10 mm Rounded	Stiff	10	2.4	240	2.8	box 20
M00561223	10 mm Rounded	Stiff	10	2.4	240	2.8	box 40
M00561231	15 mm Rounded	Stiff	15	2.4	240	2.8	box 10
M00561232	15 mm Rounded	Stiff	15	2.4	240	2.8	box 20
M00561233	15 mm Rounded	Stiff	15	2.4	240	2.8	box 40
M00561241	20 mm Rounded	Stiff	20	2.4	240	2.8	box 10
M00561242	20 mm Rounded	Stiff	20	2.4	240	2.8	box 20
M00561243	20 mm Rounded	Stiff	20	2.4	240	2.8	box 40
M00561191	25 mm Rounded	Stiff	25	2.4	240	2.8	box 10
M00561192	25 mm Rounded	Stiff	25	2.4	240	2.8	box 20
M00561193	25 mm Rounded	Stiff	25	2.4	240	2.8	box 40
M00561291	33 mm Rounded	Stiff	33	2.4	240	2.8	box 10

NEW Captivator™ COLD Single-Use Snare


Order Number	Description	Loop Width (mm)	Sheath OD (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00561101	10mm Rounded-Stiff	10	2.4	240	2.8	box 10
M00561102	10mm Rounded-Stiff	10	2.4	240	2.8	box 20
M00561103	10mm Rounded-Stiff	10	2.4	240	2.8	box 40



Contents

Tissue Resection


SNARES (CONTINUED)



Profile™ Snares

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562531	Extra Small Oval	Flexible	11	1.9	240	2.0	box 10
M00562551	Small Oval	Flexible	13	1.9	240	2.0	box 10
M00562571	Medium Oval	Flexible	27	1.9	240	2.0	box 10

HOT BIOPSY FORCEPS



Radial Jaw™ 4 Hot Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00515031	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	box 5
M00515032	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	box 20
M00515033	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	box 40

ACCESSORIES

Active Cords

Order Number	Description	Units
M00561250	Active Cord, Fits Olympus® Brand Generator	box 1
M00561270	Active Cord, Fits Endostat™, Endostat™ II, Endostat™ III, Bovie®, Valleylab®, Aspen Labs and ERBE® Brand Generators	box 1



Contents

Tissue Resection

RETRIEVAL DEVICES

NEW RescueNet™



Order Number	Description	Net Size (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
DGN-538-5	RescueNet Retrieval Device	30	2.5	230	2.8	box 5

NEW Rescue™ Forceps



Order Number	Description	Jaw Opening (mm)	Sheath O.D. (mm)	Working Length (cm)	Units
DGA-274-5	Rescue Alligator Jaw Grasping Forceps	8	2.4	230	box 5
DGA-275-5	Rescue Alligator Jaw Grasping Forceps (Long)	9	2.4	230	box 5
DGR-276-5	Rescue Rat Tooth Jaw Grasping Forceps	7	2.4	230	box 5
DGR-277-5	Rescue Rat Tooth Jaw Grasping Forceps (Pediatric)	4.5	1.8	180	box 5
DGC-383-5	Rescue Combo Grasping Forceps	8	2.4	230	box 5

POLYP TRAPS

NEW TrapEase™ Polyp Traps



Order Number	Description	Units
STE-298-25	TrapEase Quad Trap	box 25
STE-297-50	TrapEase Single Trap	box 50

LIFTING AGENTS

NEW Orise™ Gel



Order Number	Description	Gauge	Maximum Needle Extension Length (mm)	Sheath O.D. (mm)	Working Length (cm)	Units
M00519201	ORISE Gel Syringe Twin Pack Kit*	n/a	n/a	n/a	n/a	box 10
M00519211	ORISE Gel Syringe Twin Pack Kit with Interject Needle 23ga Clear**	23	4	2,3	240	box 10
M00519210	ORISE Gel Syringe Twin Pack Kit with Interject Needle 23ga Clear**	23	4	2,3	240	box 10

* Two 10 ml Syringes per tray.

** Two 10 ml Syringes per Kit.



Contents

Alphabetical Index

Product Code Index



RX Biliary System™	Biliary	SpyGlass™	EUS	Stents	Dilation	Guidewires	Hemostasis	Tissue Resection	Biopsy	Enteral Feeding	Pulmonary	Infection Prevention	RFA
--------------------	---------	-----------	-----	--------	----------	------------	------------	------------------	--------	-----------------	-----------	----------------------	-----




Contents



Alphabetical Index

Product Code Index

FORCEPS



Radial Jaw™ 4 Standard Capacity Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00513383	Radial Jaw 4 Standard Capacity	2.2	160	2.8	Yellow	box 5
M00513381	Radial Jaw 4 Standard Capacity	2.2	160	2.8	Yellow	box 20
M00513382	Radial Jaw 4 Standard Capacity	2.2	160	2.8	Yellow	box 40
M00513393	Radial Jaw 4 Standard Capacity w/Needle	2.2	160	2.8	Yellow	box 5
M00513391	Radial Jaw 4 Standard Capacity w/Needle	2.2	160	2.8	Yellow	box 20
M00513392	Radial Jaw 4 Standard Capacity w/Needle	2.2	160	2.8	Yellow	box 40
M00513403	Radial Jaw 4 Standard Capacity	2.2	240	2.8	Orange	box 5
M00513401	Radial Jaw 4 Standard Capacity	2.2	240	2.8	Orange	box 20
M00513402	Radial Jaw 4 Standard Capacity	2.2	240	2.8	Orange	box 40
M00513413	Radial Jaw 4 Standard Capacity w/Needle	2.2	240	2.8	Orange	box 5
M00513411	Radial Jaw 4 Standard Capacity w/Needle	2.2	240	2.8	Orange	box 20
M00513412	Radial Jaw 4 Standard Capacity w/Needle	2.2	240	2.8	Orange	box 40

Radial Jaw™ 4 Large Capacity Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00513303	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	box 5
M00513301	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	box 20
M00513302	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	box 40
M00513313	Radial Jaw 4 Large Capacity w/Needle	2.4	160	2.8	Yellow	box 5
M00513311	Radial Jaw 4 Large Capacity w/Needle	2.4	160	2.8	Yellow	box 20
M00513312	Radial Jaw 4 Large Capacity w/Needle	2.4	160	2.8	Yellow	box 40
M00513323	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	box 5
M00513321	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	box 20
M00513322	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	box 40
M00513333	Radial Jaw 4 Large Capacity w/Needle	2.4	240	2.8	Orange	box 5
M00513331	Radial Jaw 4 Large Capacity w/Needle	2.4	240	2.8	Orange	box 20
M00513332	Radial Jaw 4 Large Capacity w/Needle	2.4	240	2.8	Orange	box 40

FORCEPS (CONTINUED)

Radial Jaw™ 4 Jumbo Biopsy Forceps




Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00513343	Radial Jaw 4 Jumbo	2.8	160	3.2	Yellow	box 5
M00513353	Radial Jaw 4 Jumbo w/Needle	2.8	160	3.2	Yellow	box 5
M00513363	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	box 5
M00513361	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	box 20
M00513362	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	box 40
M00513373	Radial Jaw 4 Jumbo w/Needle	2.8	240	3.2	Orange	box 5
M00513371	Radial Jaw 4 Jumbo w/Needle	2.8	240	3.2	Orange	box 20
M00513372	Radial Jaw 4 Jumbo w/Needle	2.8	240	3.2	Orange	box 40

Radial Jaw™ 4 Gastropediatric Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00513443	Radial Jaw 4 Gastropediatric	1.8	160	2.0	Yellow	box 5
M00513441	Radial Jaw 4 Gastropediatric	1.8	160	2.0	Yellow	box 20
M00513453	Radial Jaw 4 Gastropediatric w/Needle	1.8	160	2.0	Yellow	box 5
M00513451	Radial Jaw 4 Gastropediatric w/Needle	1.8	160	2.0	Yellow	box 20



Multibite™ Multiple Sample Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00510101	Multibite Multiple Sample Biopsy Forceps	2.4	160	2.8	Yellow	box 5
M00510121	Multibite Multiple Sample Biopsy Forceps	2.4	240	2.8	Orange	box 5



PERCUTANEOUS ENDOSCOPIC GASTROSTOMY (PEG) DEVICES

Initial Placement

Standard PEG Kits



Order Number	Description	Units
M00568381	20Fr Standard PEG Kit, Pull	box 2
M00568391	24Fr Standard PEG Kit, Pull	box 2
M00568401	20Fr Standard PEG Kit, Push	box 2
M00568411	24Fr Standard PEG Kit, Push	box 2

Standard PEG Kits with ENFit™



Order Number	Description	Units
M00509081	EndoVive StdPEG Pull w/ENFit 20Fr	box 2
M00509101	EndoVive StdPEG Pull w/ENFit 24Fr	box 2
M00509091	EndoVive StdPEG Push w/ENFit 20Fr	box 2
M00509111	EndoVive StdPEG Push w/ENFit 24Fr	box 2

PEG Accessories


Order Number	Description	Units
M00580641	Y-Port 20Fr	box 10
M00580651	Y-Port 24Fr	box 10
M00580731	20Fr (6.67mm) Round Bolsters	box 10
M00580781	20Fr (6.67mm) and 24F (8.00mm) C-Clamps	box 10

PEG Accessories with ENFit™

Order Number	Description	Units
M00501261	EndoVive Y-Port w/ENFit 20Fr	box 10
M00501271	EndoVive Y-Port w/ENFit 24Fr	box 10
M00509241	ENFit Feed Port Connector	box 10
M00509221	ENFit Transition Connector Enteral Distal End	box 10

PERCUTANEOUS ENDOSCOPIC GASTROSTOMY (PEG) DEVICES (CONTINUED)

Initial Placement (CONTINUED)



One-Step Button™ Low Profile PEG Kits



Order Number	Description	Length (cm)	Units
M00568450	18Fr (6.00mm) Low Profile PEG Kit, Pull	1.7	box 1
M00568460	18Fr (6.00mm) Low Profile PEG Kit, Pull	2.4	box 1
M00568470	18Fr (6.00mm) Low Profile PEG Kit, Pull	3.4	box 1
M00568480	18Fr (6.00mm) Low Profile PEG Kit, Pull	4.4	box 1
M00568500	24Fr (8.00mm) Low Profile PEG Kit, Pull	1.7	box 1
M00568510	24Fr (8.00mm) Low Profile PEG Kit, Pull	2.4	box 1
M00568520	24Fr (8.00mm) Low Profile PEG Kit, Pull	3.4	box 1
M00568530	24Fr (8.00mm) Low Profile PEG Kit, Pull	4.4	box 1
M00568560	18Fr (6.00mm) Low Profile PEG Kit, Push	2.4	box 1
M00568570	18Fr (6.00mm) Low Profile PEG Kit, Push	3.4	box 1
M00568580	18Fr (6.00mm) Low Profile PEG Kit, Push	4.4	box 1
M00568600	24Fr (8.00mm) Low Profile PEG Kit, Push	1.7	box 1
M00568610	24Fr (8.00mm) Low Profile PEG Kit, Push	2.4	box 1
M00568620	24Fr (8.00mm) Low Profile PEG Kit, Push	3.4	box 1
M00568630	24Fr (8.00mm) Low Profile PEG Kit, Push	4.4	box 1

NEW One-Step Button™ Low Profile PEG Kits with ENFit™



Order Number	Description	Length (cm)	Units
M00509760	Low Profile PEG Kit with ENFit Pull, 18Fr	1.7	box 1
M00509770	Low Profile PEG Kit with ENFit Pull, 18Fr	2.4	box 1
M00509780	Low Profile PEG Kit with ENFit Pull, 18Fr	3.4	box 1
M00509790	Low Profile PEG Kit with ENFit Pull, 18Fr	4.4	box 1
M00509650	Low Profile PEG Kit with ENFit Pull, 24Fr	1.7	box 1
M00509660	Low Profile PEG Kit with ENFit Pull, 24Fr	2.4	box 1
M00509670	Low Profile PEG Kit with ENFit Pull, 24Fr	3.4	box 1
M00510230	Low Profile PEG Kit with ENFit Pull, 24Fr	4.4	box 1
M00509550	Low Profile PEG Kit with ENFit Push, 18Fr	2.4	box 1
M00509560	Low Profile PEG Kit with ENFit Push, 18Fr	3.4	box 1
M00509570	Low Profile PEG Kit with ENFit Push, 18Fr	4.4	box 1
M00509370	Low Profile PEG Kit with ENFit Push, 24Fr	1.7	box 1
M00509380	Low Profile PEG Kit with ENFit Push, 24Fr	2.4	box 1
M00509390	Low Profile PEG Kit with ENFit Push, 24Fr	3.4	box 1
M00509400	Low Profile PEG Kit with ENFit Push, 24Fr	4.4	box 1

PERCUTANEOUS ENDOSCOPIC GASTROSTOMY (PEG) DEVICES (CONTINUED)

Initial Placement (CONTINUED)

One-Step Button™ and Low Profile Replacement Button Accessories



Order Number	Description	Length (cm)	Units
M00580401	18Fr (6.00mm) Decompression Tube	1.2	box 10
M00580411	18Fr (6.00mm) Decompression Tube	1.7	box 10
M00580421	18Fr (6.00mm) Decompression Tube	2.4	box 10
M00580431	18Fr (6.00mm) Decompression Tube	3.4	box 10
M00580441	18Fr (6.00mm) Decompression Tube	4.4	box 10
M00580461	24Fr (8.00mm) Decompression Tube	1.7	box 10
M00580471	24Fr (8.00mm) Decompression Tube	2.4	box 10
M00580481	24Fr (8.00mm) Decompression Tube	3.4	box 10
M00580491	24Fr (8.00mm) Decompression Tube	4.4	box 10
M00580171	18Fr (6.00mm) Bolus Straight Feeding Adapter		box 10
M00580181	18Fr (6.00mm) Continuous Right Angle Feeding Adapter		box 10
M00580231	24Fr (8.00mm) Bolus Straight Feeding Adapter		box 10
M00580241	24Fr (8.00mm) Continuous Right Angle Feeding Adapter		box 10
M00580271	28Fr (9.33mm) Bolus Straight Feeding Adapter		box 10
M00580281	28Fr (9.33mm) Continuous Right Angle Feeding Adapter		box 10


NEW One-Step Button and Low Profile Replacement Button Accessories with ENFit™



Order Number	Description	Length (cm)	Units
M00501121	Decompression Tube w/ENFit, 18Fr	1.2	box 10
M00501131	Decompression Tube w/ENFit, 18Fr	1.7	box 10
M00501141	Decompression Tube w/ENFit, 18Fr	2.4	box 10
M00501151	Decompression Tube w/ENFit, 18Fr	3.4	box 10
M00501161	Decompression Tube w/ENFit, 18Fr	4.4	box 10
M00501181	Decompression Tube w/ENFit, 24Fr	1.7	box 10
M00501191	Decompression Tube w/ENFit, 24Fr	2.4	box 10
M00501201	Decompression Tube w/ENFit, 24Fr	3.4	box 10
M00501211	Decompression Tube w/ENFit, 24Fr	4.4	box 10
M00501091	Bolus Feeding Set Straight Adapter w/ENFit, 18Fr		box 10
M00501061	Right Angle Feeding Set Continous Adapter w/ENFit, 18Fr		box 10
M00501101	Bolus Feeding Set Straight Adapter w/ENFit, 24Fr		box 10
M00501071	Right Angle Feeding Set Continous Adapter w/ENFit, 24Fr		box 10
M00501111	Bolus Feeding Set Straight Adapter w/ENFit, 28Fr		box 10
M00501081	Right Angle Feeding Set Continous Adapter w/ENFit, 28Fr		box 10


PERCUTANEOUS ENDOSCOPIC GASTROSTOMY (PEG) DEVICES (CONTINUED)

Replacements



Replacement Gastrostomy Feeding Tubes with ENFit™

Enfit Kit	Enfit Description	Units
M00509881	12Fr (4.00mm) Repl. G-Tube, Straight	box 5
M00509891	14Fr (4.67mm) Repl. G-Tube, Straight	box 5
M00509901	16Fr (5.33mm) Repl. G-Tube, Straight	box 5
M00509911	18Fr (6.00mm) Repl. G-Tube, Straight	box 5
M00509921	20Fr (6.67mm) Repl. G-Tube, Straight	box 5
M00509931	22Fr (7.33mm) Repl. G-Tube, Straight	box 5
M00509941	24Fr (8.00mm) Repl. G-Tube, Straight	box 5
M00509961	12Fr (4.00mm) Repl. G-Tube, Right Angle	box 5
M00509971	14Fr (4.67mm) Repl. G-Tube, Right Angle	box 5
M00509981	16Fr (5.33mm) Repl. G-Tube, Right Angle	box 5
M00509991	18Fr (6.00mm) Repl. G-Tube, Right Angle	box 5
M00510001	20Fr (6.67mm) Repl. G-Tube, Right Angle	box 5
M00510011	22Fr (7.33mm) Repl. G-Tube, Right Angle	box 5
M00510021	24Fr (8.00mm) Repl. G-Tube, Right Angle	box 5




Low Profile Replacement Button

Order Number	Description	Length (cm)	Colour Code	Units
M00568260	18Fr (6.00mm) Low Profile Replacement Button	1.2	Yellow	box 1
M00568270	18Fr (6.00mm) Low Profile Replacement Button	1.7	Yellow	box 1
M00568280	18Fr (6.00mm) Low Profile Replacement Button	2.4	Yellow	box 1
M00568290	18Fr (6.00mm) Low Profile Replacement Button	3.4	Yellow	box 1
M00568300	18Fr (6.00mm) Low Profile Replacement Button	4.4	Yellow	box 1
M00568320	24Fr (8.00mm) Low Profile Replacement Button	1.7	Blue	box 1
M00568330	24Fr (8.00mm) Low Profile Replacement Button	2.4	Blue	box 1
M00568340	24Fr (8.00mm) Low Profile Replacement Button	3.4	Blue	box 1
M00568350	24Fr (8.00mm) Low Profile Replacement Button	4.4	Blue	box 1
M00568170	28Fr (9.33mm) Low Profile Replacement Button	1.5	Green	box 1
M00568180	28Fr (9.33mm) Low Profile Replacement Button	2.8	Green	box 1
M00568190	28Fr (9.33mm) Low Profile Replacement Button	4.3	Green	box 1

For accessories please refer to the One-Step Button™ and Low Profile Button Accessories

JEJUNAL FEEDING TUBES AND REPLACEMENT GASTROSTOMY FEEDING TUBES


Replacement Feeding Tubes



Through-The-PEG (TTP) Jejunal Feeding Tube





Order Number	Description	Tip Style	Units
M00566320	8.5Fr (2.83mm) TTP Jejunal Feeding Tube	Bent Tip	box 1
M00566330	8.5Fr (2.83mm) TTP Jejunal Feeding Tube	Pigtail Tip	box 1
M00566340	12Fr (4.00mm) TTP Jejunal Feeding Tube	Bent Tip	box 1
M00566350	12Fr (4.00mm) TTP Jejunal Feeding Tube	Pigtail Tip	box 1




Replacement Gastrostomy Feeding Tubes



Order Number	Description	Units
M00582040	12Fr (4.00mm) Repl. G-Tube, Straight	box 1
M00582050	14Fr (4.67mm) Repl. G-Tube, Straight	box 1
M00582060	16Fr (5.33mm) Repl. G-Tube, Straight	box 1
M00582070	18Fr (6.00mm) Repl. G-Tube, Straight	box 1
M00582080	20Fr (6.67mm) Repl. G-Tube, Straight	box 1
M00582860	22Fr (7.33mm) Repl. G-Tube, Straight	box 1
M00582090	24Fr (8.00mm) Repl. G-Tube, Straight	box 1
M00582100	28Fr (9.33mm) Repl. G-Tube, Straight	box 1
M00582110	12Fr (4.00mm) Repl. G-Tube, Right Angle	box 1
M00582120	14Fr (4.67mm) Repl. G-Tube, Right Angle	box 1
M00582130	16Fr (5.33mm) Repl. G-Tube, Right Angle	box 1
M00582140	18Fr (6.00mm) Repl. G-Tube, Right Angle	box 1
M00582150	20Fr (6.67mm) Repl. G-Tube, Right Angle	box 1
M00582870	22Fr (7.33mm) Repl. G-Tube, Right Angle	box 1
M00582160	24Fr (8.00mm) Repl. G-Tube, Right Angle	box 1
M00582170	28Fr (9.33mm) Repl. G-Tube, Right Angle	box 1



ALAIR™ BRONCHIAL THERMOPLASTY (BT) SYSTEM



Alair™ Bronchial Thermoplasty (BT) System

Order Number	Model	Description	Active Electrode Length (mm)	Tip Diameter (mm)	Minimum Working Channel (mm)
M005ATS25020	ATS 2-5	Alair BT Catheter – Global	5.0	~1.5	2.0
M005ATS20000	ATS 200	Alair RF Controller – monopolar (return electrode not included)			

Packaged one per box

Alair™ BT Controller Accessory Kits ATS 201*

Order Number	Model	Description
M005ATS20120	ATS 201	British Type G
M005ATS20130	ATS 201	EU Hybrid Type E/F
M005ATS20140	ATS 201	Australian Type I
M005ATS20150	ATS 201	Danish Type K
M005ATS20160	ATS 201	Brazil IEC 60906-1
M005ATS20180	ATS 201	Swiss Type J

Packaged one per box

* Including Power Cord with Molded Plug, Foot Switch, Operating Manual.

ENDOSCOPIC BRONCHIAL ULTRASOUND



Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needles Olympus Scope Compatible



Order Number	Description	Minimum Working Channel (mm)	Sheath O.D. (mm)	Units
M00558220	Expect Pulmonary 22ga – Olympus Scope Compatible	2.0	1.6	box 1
M00558221	Expect Pulmonary 22ga – Olympus Scope Compatible	2.0	1.6	box 5
M00558250	Expect Pulmonary 25ga – Olympus Scope Compatible	2.0	1.4	box 1
M00558251	Expect Pulmonary 25ga – Olympus Scope Compatible	2.0	1.4	box 5
M00558731	Expect Pulmonary Needle Adaptor – Olympus Scope Compatible			box 10

ENDOSCOPIC BRONCHIAL ULTRASOUND (CONTINUED)


Acquire™ Pulmonary Endobronchial Ultrasound Fine Biopsy Needles

Olympus Scope Compatible



Order Number	Description	Needle Size (gauge) / (mm)	Minimum Working Channel (mm)	Sheath Diameter (mm)	Units
M00552351	22ga Acquire Pulmonary – Olympus scope compatible	22/0.72	2.0	1.6	box 5
M00552401	25ga Acquire Pulmonary – Olympus scope compatible	25/0.52	2.0	1.4	box 5

AIRWAY STENTS



Ultraflex™ Tracheobronchial Covered Stent System

Distal Release



Order Number	Description	Stent O.D. (mm)	Expanded Stent Length (mm)	Covered Length (mm)
M00564820	Ultraflex Tracheal NG Stent System	14	80	65
M00564830	Ultraflex Tracheal NG Stent System	16	40	25
M00564840	Ultraflex Tracheal NG Stent System	16	60	45
M00564850	Ultraflex Tracheal NG Stent System	16	80	65
M00564860	Ultraflex Tracheal NG Stent System	18	40	25
M00564870	Ultraflex Tracheal NG Stent System	18	60	45
M00564880	Ultraflex Tracheal NG Stent System	18	80	65
M00564890	Ultraflex Tracheal NG Stent System	20	40	25
M00564900	Ultraflex Tracheal NG Stent System	20	60	45
M00564910	Ultraflex Tracheal NG Stent System	20	80	65
M00564740	Ultraflex Tracheobronchial Stent System	8	40	25
M00564750	Ultraflex Tracheobronchial Stent System	10	30	15
M00564760	Ultraflex Tracheobronchial Stent System	10	40	25
M00564770	Ultraflex Tracheobronchial Stent System	12	30	15
M00564780	Ultraflex Tracheobronchial Stent System	12	40	25
M00564790	Ultraflex Tracheobronchial Stent System	14	30	15
M00564800	Ultraflex Tracheobronchial Stent System	14	40	25
M00564810	Ultraflex Tracheobronchial Stent System	14	60	45

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Jagwire™ Pulmonary Guidewire, M00515171 / .038in (0.97mm) Amplatz Super Stiff™ Guidewire, M00550090

Additional Specifications for Ultraflex Tracheobronchial Stents: Delivery System Working Length is 95cm, Maximum Catheter Shaft O.D. is 16F (5.33mm), Compressed Stent on catheter yields Delivery System O.D. range 15-25F (5.00-8.33mm)

AIRWAY STENTS (CONTINUED)



Ultraflex™ Tracheobronchial Uncovered Stent System

Proximal Release




Order Number	Description	Stent O.D. (mm)	Expanded Stent Length (mm)
M00564710	Ultraflex Tracheobronchial Stent System	14	40
M00564720	Ultraflex Tracheobronchial Stent System	14	60
M00564640	Ultraflex Tracheobronchial Stent System	8	20
M00564650	Ultraflex Tracheobronchial Stent System	8	40
M00564660	Ultraflex Tracheobronchial Stent System	10	20
M00564670	Ultraflex Tracheobronchial Stent System	10	40
M00564680	Ultraflex Tracheobronchial Stent System	12	20
M00564690	Ultraflex Tracheobronchial Stent System	12	40
M00564700	Ultraflex Tracheobronchial Stent System	14	20

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Jagwire™ Pulmonary Guidewire, M00515171 / .038in (0.97mm) Amplatz Super Stiff™ Guidewire, M00550090

Additional Specifications for Ultraflex Tracheobronchial Stents: Delivery System Working Length is 95cm, Maximum Catheter Shaft O.D. is 16F (5.33mm), Compressed Stent on catheter yields Delivery System O.D. range 15-18F (5.00-6.00mm)



Ultraflex™ Tracheobronchial Uncovered Stent System

Distal Release




Order Number	Description	Stent O.D. (mm)	Expanded Stent Length (mm)
M00564530	Ultraflex Tracheobronchial Stent System	16	40
M00564540	Ultraflex Tracheobronchial Stent System	16	60
M00564560	Ultraflex Tracheobronchial Stent System	18	40
M00564570	Ultraflex Tracheobronchial Stent System	18	60
M00564590	Ultraflex Tracheobronchial Stent System	20	40
M00564600	Ultraflex Tracheobronchial Stent System	20	60
M00564500	Ultraflex Tracheobronchial Stent System	10	30
M00564510	Ultraflex Tracheobronchial Stent System	12	30
M00564520	Ultraflex Tracheobronchial Stent System	14	30

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Jagwire™ Pulmonary Guidewire, M00515171 / .038in (0.97mm) Amplatz Super Stiff™ Guidewire, M00550090

Additional Specifications for Ultraflex Tracheobronchial Stents: Delivery System Working Length is 95cm, Maximum Catheter Shaft O.D. is 16F (5.33mm), Compressed Stent on catheter yields Delivery System O.D. range 15-20F (5.00-6.67mm)

BIOPSY FORCEPS



Radial Jaw™ 4 Pulmonary Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00515181	Radial Jaw 4 Pulmonary Standard Capacity	1.8	100	2.0	Blue	box 5
M00515182	Radial Jaw 4 Pulmonary Standard Capacity	1.8	100	2.0	Blue	box 20
M00515191	Radial Jaw 4 Pulmonary Standard Capacity w/Needle	1.8	100	2.0	Blue	box 5
M00515192	Radial Jaw 4 Pulmonary Standard Capacity w/Needle	1.8	100	2.0	Blue	box 20
M00515201	Radial Jaw 4 Pulmonary Large Capacity	2.2	100	2.8	Blue	box 5
M00515202	Radial Jaw 4 Pulmonary Large Capacity	2.2	100	2.8	Blue	box 20


MINI-FORCEPS

CoreDx™ Pulmonary Mini-Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00515220	Pulmonary Mini-Forceps	1.0	130	1.2	box 1

TRANSBRONCHIAL ASPIRATION NEEDLES



eXcelon™ Transbronchial Aspiration Needles



Order Number	Description	Gauge	NEEDLE O.D. (mm)	I.D. (mm)	Needle Length (mm)	Catheter Length (cm)	Sheath O.D. (mm)	Units
M00564101	eXcelon Transbronchial Aspiration Needle	19	1.07	0.69	15	130	1.8	box 5
M00564111	eXcelon Transbronchial Aspiration Needle	20	0.90	0.58	15	130	1.8	box 5
M00564121	eXcelon Transbronchial Aspiration Needle	21	0.81	0.50	15	130	1.8	box 5


CYTOLOGY BRUSHES

Celebrity™ Cytology Brushes – Bronchoscope



Order Number	Description	Required Working Channel (mm)	Bristle O.D. (mm)	Sheath Length (cm)	Units
M00516001	Celebrity Cytology Brush	2.0	1.0	140	box 10
M00516011	Celebrity Cytology Brush	2.0	1.5	140	box 10
M00516071	Celebrity Cytology Brush	2.0	1.9	100	box 10
M00516151	Celebrity Cytology Brush	2.0	1.9	150	box 10

AIRWAY RETRIEVAL BASKETS



Zero Tip™ Airway Retrieval Baskets



Order Number	Description	Sheath O.D. (mm)	Working Length (cm)	Basket Opening (mm)
M00513200	Zero Tip Airway Retrieval Basket	0.8	120	12
M00513210	Zero Tip Airway Retrieval Basket	1.0	120	16

Packaged one per box

DILATION BALLOONS



CRE™ Pulmonary Balloon Dilatation Catheter



Order Number	Description	Balloon Length (cm)	Initial Balloon O.D. (mm)	Intermediate Balloon O.D. (mm)	Maximum Balloon O.D. (mm)	Working Catheter Length (cm)
M00550300	CRE Pulmonary Balloon Dilatation Catheter	5.5	12	13.5	15	110
M00550310	CRE Pulmonary Balloon Dilatation Catheter	5.5	15	16.5	18	110
M00550320	CRE Pulmonary Balloon Dilatation Catheter	5.5	18	19.0	20	110
M00550330	CRE Pulmonary Balloon Dilatation Catheter	3.0	8	9.0	10	110
M00550340	CRE Pulmonary Balloon Dilatation Catheter	3.0	10	11.0	12	110
M00550350	CRE Pulmonary Balloon Dilatation Catheter	3.0	12	13.5	15	110

Packaged one per box


GUIDEWIRES

Jagwire™ Pulmonary Guidewire



Order Number	Description	O.D. (in) / (mm)	Length (cm)	Units
M00515171	Jagwire Pulmonary Guidewire	.035 / 0.89	180	box 2


INFLATION SYSTEM



Alliance™ II Inflation System



Order Number	Description	Units
M00550620	Alliance II Integrated Inflation/Lithotripsy Device	box 1
M00550601	Alliance II Single-Use Syringe/Gauge Assembly	box 5



IRRIGATION & INSUFFLATION

NEW Hydra™ Irrigation System

Order Number	Description	Compatibility	Units
SIT-575-15	Hydra Irrigation Tubing	Olympus™ AFU-100, ERBE™ EIP 2, Medivators™ EndoStratus™ & EndoGator™	box 15
SIT-577-15	Hydra Irrigation Tubing	Olympus OFP2	box 15
SCT-611-10	Hydra Water Bottle Cap with CO2	Fuse	box 10
SCT-466-10	Hydra Water Bottle Cap	Fuse	box 10
SIT-470-100	⊗ Hydra Disposable Auxiliary Water Jet Connector	Fuse	box 100
SCT-350-10	Hydra Water Bottle Cap	Olympus	box 10
SIT-361-100	⊗ Hydra Disposable Auxiliary Water Jet Connectors	Olympus 160, 180 & 190 Series	box 100
SCT-609-10	Hydra Water Bottle Cap with CO2	Air/Water Tubing Length: 16" (41 cm) CO2 Tubing Length: 18" (46 cm) Olympus	box 10
SIT-355-15	Hydra Irrigation Tubing		box 15
SIT-611-10	Hydra Water Bottle Cap with CO2	Fuse	box 10
SIT-351-10	Hydra Water Bottle Cap with CO2	Olympus	box 10
SIT-467-10	Hydra Water Bottle Cap with CO2	Fuse	box 10

CLEANING BRUSHES

NEW HedgeHog™ Cleaning Brushes



Order Number	Description	Units
SBS-227-50	HedgeHog Single-End Channel Brush	box 50
SBD-228-50	HedgeHog Double-End Channel Brush	box 50
SBD-289-50	HedgeHog Dual-End Channel/Valve Brush	box 50
SBD-291-50	HedgeHog Sweeper Brush	box 50
SBB-382-50	HedgeHog Double Bundle	box 50

VALVES

NEW Orca™ Air/Water and Suction Valves



Order Number	Description	Units
SUV-617-50	Orca Single-Use Air/Water and Suction Valves	Olympus™ (non EUS Scopes) box 50

Infection Prevention

SCOPE TRANSPORT

NEW CinchPad™ Procedure & Endoscope Transport Pad



Order Number	Description	Units
SPX-017-25	CinchPad Endoscope Transport Pad	box 25

BITE BLOCKS

NEW Blox™ Bite Blocks



Order Number	Description	Units
SBT-546-100	Blox Bite Block 60Fr with Dental Rim	box 100
SBT-114-100	Blox Bite Block 54Fr with Dental Rim	box 100
SBT-321-100	Blox Bite Block 54Fr without Dental Rim	box 100
SBP-236-40	Pediatric Blox Bite Block 38Fr	box 40

SPONGES

NEW Koala™ Sponge



Order Number	Description	Units
SSB-210-50	Koala Cleaning Sponge	box 50
SES-240-50	Koala Cleaning Sponge Imbibed w/ Pure	box 50

PROCEDURE GOWNS

NEW EndoArmor™ Procedural Gowns



Order Number	Description	Units
SGE-368-100	EndoArmor GI Procedure Gowns	box 100
SGE-389-100	EndoArmor+ Procedure Gowns	box 100


LUBRICATING GELS

NEW EndoGlide™ Lubricating Gels



Order Number	Description	Units
SLT-394-100	1.1 oz EndoGlide+ Lubricating Gel	box 100
SLT-395-50	2 oz EndoGlide+ Lubricating Gel	box 50





RX Biliary System™	Biliary	SpyGlass™	EUS	Stents	Dilation	Guidewires	Hemostasis	Tissue Resection	Biopsy	Enteral Feeding	Pulmonary	Infection Prevention	RFA
--------------------	---------	-----------	-----	--------	----------	------------	------------	------------------	--------	-----------------	-----------	----------------------	-----




Contents



Alphabetical Index

Product Code Index

ABLATION CATHETER AND CABLES



NEW Habib™ EndoHPB Bipolar Radiofrequency Catheter

Order Number	Description	Generator
M00500070	Habib Endo HPB Bipolar Radiofrequency Catheter	N/A
5100	Adapter Cable Bipolar, 28 mm pins	KLS Martin Maxium, Olympus ESG-100
5420	Adapter Cable Bipolar, Internal 4 mm pin	Erbe ICC200, ICC300, ICC350, Erbe Vio 200 or 300 (D/S), Genii GI 4000
5700	Adapter Cable Bipolar, 22 mm spacing	Erbe Vio 3, Erbe ICC200, ICC300, ICC350, Erbe Vio 200 or 300 (D/S)



Find more information on Boston Scientific products
by using your smartphone to scan this code or visit
www.bostonscientific.com/global-endoscopy

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Material not intended for use in France.

**Boston
Scientific**
Advancing science for life™
www.bostonscientific.eu

© 2020 Boston Scientific Corporation
or its affiliates. All rights reserved.
DINEND2436EA



Contents



Alphabetical
Index

Product Code
Index