


EndoArmor™+ Non-Sterile Gown

AAMI Level 3

Technical Document

EndoArmor+ Non-Sterile Gown Design Overview



All images owned by Boston Scientific

ID	Component	AAMI Barrier Performance Specifications	Physical and Mechanical Specifications
A	Body (Front)	Critical Zone: AAMI Level 3 Outside of Critical Zone: AAMI Level 1	Material: Polyethylene film laminated with nonwoven spunbond polypropylene Color: Blue
B	Sleeves	Critical Zone: AAMI Level 3	Weight: Polyethylene film 23g/m ² Nonwoven spunbond polypropylene 22g/m ² Lamination Adhesive 3g/m ² Tensile Strength: ≥ 6.31N (≥ 1.42 lbf) Tear Resistance: ≥ 1.5N (≥ 0.34 lbf)
C	Body (Back)	Non-protective	Material: Nonwoven spunbond polypropylene Color: Blue Weight: 35 g/m ² Tensile Strength: ≥ 6.31N (≥ 1.42 lbf) Tear Resistance: ≥ 1.5N (≥ 0.34 lbf)
D	Sleeve/Sleeve Ultrasonic Welded Seam	AAMI Level 3	Seam Strength: ≥ 6.31N (≥ 1.42 lbf)
E	Sleeve/Body (Front) Ultrasonic Welded Seam	AAMI Level 1	Seam Strength: ≥ 6.31N (≥ 1.42 lbf)
F	Shoulder Ultrasonic Welded Seam	Non-protective	Seam Strength: ≥ 6.31N (≥ 1.42 lbf)
G	Body/Body Sewn Seam	Non-protective	Material: Polyester Thread Color: White Weight: 150 Denier Tensile Strength: ≥ 6.31N (≥ 1.42 lbf)
H	Waist Tie	N/A	Material: Nonwoven spunbond polypropylene Color: White Weight: 35 g/m ²
I	Cuff	N/A	Material: Elastic polyester fabric Color: White
J	Hook and Loop Closure	N/A	Material: Polyester Color: White

EndoArmor + Non-Sterile Gown has not been cleared by FDA. Gown is not recommended for use when an FDA-cleared surgical gown is available. Gown is non-sterile and should not be used in sterile surgical settings.

All trademarks are the property of their respective owners. Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device. **CAUTION:** U.S. Federal law restricts this device to sale by or on the order of a physician.