



Endoscopy

2022 Product Reference Guide



We're all in.

Because you're all in on delivering optimal patient care, we've made it our specialty to support yours. From clinical research and education to devices for gastrointestinal and pulmonary conditions that are changing the standard of care, we work with you to help solve the challenges that matter most to you and your patients.

Visit bostonscientific.com/gastro to learn more about:

- › Products
- › Education & Training
- › Ordering Information
- › Solutions
- › BSC's eCommerce solution



For Single-Use Only.
Do Not Reuse.



Magnetic Resonance
Safe



Magnetic Resonance
Conditional¹

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Interactive Product Guide

Tap or click on product name for information

36	Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device
82	Acquire™ Pulmonary Endobronchial Ultrasound Fine Needle Biopsy Device
12 71	Active Cords
17	Advanix™ Center Bend Preloaded RX Stent Systems
18	Advanix Center Bend - Single Stent
19	Advanix Double Pigtail - Single Stent
17	Advanix Duodenal Bend Preloaded RX Stent Systems
18	Advanix Duodenal Bend - Single Stent
22	Advanix Pancreatic Stent Kit - Pigtail, Leading Barb and No Leading Barb
22	Advanix Pancreatic Stent Kit - Straight, No Leading Barb and Leading Barb
26	Advanix Pancreatic Stent - Pigtail, Leading Barb
25	Advanix Pancreatic Stent - Pigtail, No Leading Barb
24	Advanix Pancreatic Stent - Straight, Leading Barb
23	Advanix Pancreatic Stent - Straight, No Leading Barb
76	Aegis One-Hour Indicator
41	Agile Fully Covered Esophageal Stent System
41	Agile Partially Covered Esophageal Stent System
78	Alair™ Bronchial Thermoplasty System
55 84	Alliance™ II Integrated Inflation / Lithotripsy Device
55 84	Alliance II Single-Use Syringe / Gauge Assembly
44 84	Amplatz Super Stiff™ Guidewire
35	Autolith™ Touch Bipolar Electrohydraulic Lithotripter

35	Autolith Touch Extender Cable
35	Autolith Touch Foot Switch
6	Autotome™ RX Cannulating Sphincterotome
37	AXIOS™ Stent and Delivery System
37	AXIOS Stent and Electrocautery Enhanced Delivery System
35	Biliary EHL Probe
73	Blox™ Endoscopic Bite Block with Elastic Cloth Strap (60 FR)
73	Blox Endoscopic Bite Block with Elastic Cloth Strap (54 FR)
73	Blox Endoscopic Bite Block with Elastic Cloth Strap (54 FR) (w/o Dental Retention)
70	Captiflex™ Single-Use Snares
67	Captivator™ COLD Single-Use Snare
68	Captivator Single-Use Snares
67	Captivator EMR Device
68	Captivator II Single-Use Snares
82	Cellebri™ Single-Use Endoscopic Cytology Brush
75	CinchPad™ Endoscope Transport Pad
83	Compliance BronchKit
76	Compliance EndoKit™
11	Contour™ ERCP Cannulas
82	CoreDx™ Pulmonary Mini-Forceps for Endobronchial Ultrasound
53	CRE™ Fixed Wire Esophageal Balloon Dilatation Catheters
30	CRE PRO Wireguided Balloon Dilatation Catheters

Note: RX devices are short wire devices.

Interactive Product Guide

Tap or click on product name for information

51	CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	63	EndoVive One-Step Button™ Low Profile PEG Kits
1	CRE PRO Wireguided Esophageal / Pyloric / Colonic Biliary Balloon Dilatation Catheters	59	EndoVive One-Step Button Low Profile PEG Kits with ENFit
8 4	CRE Pulmonary Balloon Dilatation Catheter	64	EndoVive Replacement Gastrostomy Feeding Tubes, Straight or Right Angle Bolster
50	CRE RX Biliary Balloon Dilatation Catheters	60	EndoVive Replacement Gastrostomy Feeding Tubes with ENFit, Straight or Right Angle Bolster
55	CRE SteriFlate Disposable Inflation Device	62	EndoVive Safety PEG Kit
52	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	58	EndoVive Safety PEG Kit with ENFit
52	CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	62	EndoVive Standard PEG Kit
61	Distal End ENFit™ Transition Connector	58	EndoVive Standard PEG Kit with ENFit
4	Dreamtome™ RX Cannulating Sphincterotome	66	EndoVive Through-The-PEG (TTP) Jejunal Feeding Tube
1	Dreamwire™ High Performance Guidewire	61	ENFit Transition Connectors
49	Electrosurgical Bipolar Cable Adapter	14	Epic™ Biliary Endoscopic Stent System
32 33	Encore™ 26 Inflation Device	35	EXALT™ Controller
74	EndoArmor+ Non-Sterile Surgical Gown	77	EXALT Model B Single-Use Bronchoscope
74	Endoglide™+ Sterile Bacteriostatic Lubricating Gel	35	EXALT Model D Single-Use Duodenoscope
62	EndoVive™ Initial Placement PEG Accessories	83	eXcelon™ Transbronchial Aspiration Needle
58	EndoVive Initial Placement PEG Accessories with ENFit™	82	Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle / Adaptor
65	EndoVive Low Profile Replacement Button	36	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle
65	EndoVive Low Profile Replacement Button Accessories	28	Extractor Pro RX Retrieval Balloon Catheter – Distal Injection
60	EndoVive Low Profile Replacement Button with ENFit	28	Extractor Pro RX Retrieval Balloon Catheter – Proximal Injection
61	EndoVive Low Profile Replacement Button Accessories with ENFit	28	Extractor Pro RX-S Retrieval Balloon Catheter – Distal Injection
63	EndoVive One-Step Button™ and Low Profile Replacement Button Accessories	28	Extractor Pro RX-S Retrieval Balloon Catheter – Proximal Injection
59	EndoVive One-Step Button and Low Profile Replacement Button Accessories with ENFit™	29	Extractor Pro XL Retrieval Balloon Catheter – Distal Injection

Note: RX devices are short wire devices.

Interactive Product Guide

Tap or click on product name for information

29	Extractor Pro XL Retrieval Balloon Catheter – Proximal Injection
20	Flexima™ Duodenal Bend Biliary Stent with Delivery System
21	Flexima Single Pigtail Nasobiliary Catheter
46	Gold Probe™ Single-Use Electrohemostasis Catheters
31	Habib™ EndoHPB Bipolar Radiofrequency Catheter
73	HedgeHog Double-End Channel Brush
73	HedgeHog Dual-End Channel / Valve Brush
73	HedgeHog Single-End Channel Brush
73	HedgeHog Sweeper Channel Brush
73	HedgeHog Valve Brush
31	Hurricane™ RX Biliary Balloon Dilatation Catheter
74	Hydra™ Disposable Auxiliary Water Jet Connector
74	Hydra Irrigation Tubing
1	Hydra Jagwire™ High Performance Guidewire
74	Hydra Water Bottle Cap
4	Hydratome™ RX Cannulating Sphincterotome
55	Inflation / Deflation Devices and Accessories
46	Injection Gold Probe™ Single-Use Electrohemostasis / Injection Catheters
47	Interject™ Clear Single-Use Injection Therapy Needle Catheters
47	Interject Contrast Single-Use Injection Therapy Needle Catheters
2	Jagtail™ Guidewire Extension
5	Jagtome Revolution RX Cannulating Sphincterotome

5	Jagtome™ RX Cannulating Sphincterotome
2	Jagwire™ High Performance Guidewire
2	Jagwire Revolution High Performance Guidewire
75	Koala™ Cleaning Sponge
75	Koala Cleaning Sponge Plus Pure Enzymatic
37	LumiCoil™ Platinum Fiducial Markers – Figure 8 Configuration
37	LumiCoil Platinum Fiducial Markers – Straight Configuration
32	MaxForce™ Biliary Balloon Dilatation Catheter
54	MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter
8	MicroKnife™ XL Triple-Lumen Needle Knife
57	Multibite™ Single-Use Multiple Sample Biopsy Forceps
19	NaviFlex™ RX Delivery System
27	NaviFlex RX Pancreatic Delivery System & NaviFlex Pushers
3	NaviPro™ Hydrophilic Guidewire
3	NovaGold™ High Performance Guidewire
76	Orca™ Sterile Single-Use Air/Water & Suction Valves
76	OrcaPod™ 3 Orca Sterile Single-Use Air/Water, Suction, and Seal™ Biopsy Valves
76	OrcaPod 4 - Orca Sterile Single-Use Air/Water, Suction, and Seal Biopsy Valves with WaterJet Connector
71	ORISE™ Gel
71	ORISE ProKnife
55	Pneumatic Inflator
69	Profile™ Single-Use Snares

Note: RX devices are short wire devices.

Interactive Product Guide

Tap or click on product name for information

84	Pulmonary Jagwire™ Pulmonary Guidewire
73	Pure™ Enzymatic Detergent - Fragrance Free
73	Pure Enzymatic Detergent - Spearmint Fresh Scent
57	Radial Jaw™ 4 Single-Use Hot Biopsy Forceps
56	Radial Jaw™ 4 Single-Use Jumbo Biopsy Forceps
56	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps
57	Radial Jaw 4 Single-Use Pediatric Biopsy Forceps
81	Radial Jaw 4 Single-Use Pulmonary Biopsy Forceps
56	Radial Jaw 4 Single-Use Standard Capacity Biopsy Forceps
10	Rapid Exchange XL Cannula
12	Rapid Refill Continuous Injection System
72	Rescue™ Alligator Grasping Forceps
72	Rescue Alligator Grasping Forceps, Long
72	Rescue Rat Tooth / Alligator Grasping Forceps
72	Rescue Rat Tooth Grasping Forceps
72	Rescue Rat Tooth Grasping Forceps - Pediatric
45	Resolution™ Clip
45	Resolution 360™ Clip
45	Resolution 360 ULTRA Clip
54	Rigiflex™ II Single-Use Achalasia Balloon Dilatation Catheter
6	RX Cholangiogram Kits
32	RX Cytology Brush Wireguided Cytology Brush

20	RX Duodenal Bend Biliary Stent with RX Delivery System
10	RX ERCP Cannulas
6	RX Needle Knife XL Triple-Lumen Needle Knife
76	Seal™ Single-Use Biopsy Valve
70	Sensation™ Single-Use Short Throw Snares
69	Single-Use Rotatable Snares
48	Speedband Superview Super 7™ Multiple Band Ligators
34	SpyBite™ Max Biopsy Forceps
34	SpyGlass™ DS Direct Visualization System – Capital Equipment
33	SpyGlass™ Discover Balloon Dilation Catheter
33	SpyGlass™ Discover Controller
33	SpyGlass™ Discover Digital Catheter
33	SpyGlass™ Discover Imager™ II IOC Catheter
33	SpyGlass™ Discover Discover Jagwire™ Guidewire
33	SpyGlass™ Discover Retrieval Basket Box
34	SpyGlass Retrieval Basket
34	SpyGlass Retrieval Snare
34	SpyScope™ DS Access & Delivery Catheter
34	SpyScope DS II Access & Delivery Catheter
74	Sterile Lubricating Jelly
9	Stonetome™ Stone Removal Device
75	Suction Tubing

Note: RX devices are short wire devices.

Interactive Product Guide

Tap or click on product name for information

33	Super Sheath CBDE for use with the SpyGlass
11	Tandem™ XL Triple-Lumen ERCP Cannula
72	TrapEase™ Polyp Traps
29	Trapezoid™ RX Wireguided Retrieval Baskets
8	TRUEtome™ Cannulating Sphincterotome
7	TRUEtome Dreamwire™ Cannulating Sphincterotome
7	TRUEtome Hydra Cannulating Sphincterotome
7	TRUEtome Jag Cannulating Sphincterotome
7	TRUEtome Revolution RX Cannulating Sphincterotome
40	Ultraflex™ Precision Single-Use Colonic Stent System
43	Ultraflex™ Esophageal NG Covered Stent System – Distal Release
43	Ultraflex Esophageal NG Covered Stent System – Proximal Release
43	Ultraflex Large Esophageal NG Covered Stent System – Distal Release
44	Ultraflex Esophageal NG Covered Stent System – Proximal Release
80	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release
44	Ultraflex Esophageal NG Uncovered Stent System – Distal Release
44	Ultraflex Esophageal NG Uncovered Stent System – Proximal Release
79	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release
79	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release
10	Ultratome™ Double-Lumen Sphincterotome
9	Ultratome XL Triple-Lumen Sphincterotome

13	WallFlex™ Biliary RX Fully Covered Stent System RMV
13	WallFlex Biliary RX Partially Covered Stent System
13	WallFlex Biliary RX Uncovered Stent System
38	WallFlex Colonic Stent System
38	WallFlex Colonic Stent System with Anchor Lock Delivery System
39	WallFlex Duodenal Stent System
39	WallFlex Duodenal Stent System with Anchor Lock Delivery System
42	WallFlex Fully Covered Esophageal Stent System
42	WallFlex Partially Covered Esophageal Stent System
16	WALLSTENT™ Biliary Endoprosthesis Stents
40	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP™ Plus Delivery System
15	WALLSTENT RX Biliary Endoprosthesis Stent System
$\frac{3}{40}$	WALLSTENT Super Stiff Guidewire
83	Zero Tip™ Airway Retrieval Basket

Note: RX devices are short wire devices.

► Biliary Access



Biliary Access

Guidewires

► Amplatz Super Stiff™ Guidewire	1
► Dreamwire™ High Performance Guidewire	1
► Hydra Jagwire™ High Performance Guidewire	1
► Jagtail™ Guidewire Extension	2
► Jagwire™ High Performance Guidewire	2
► Jagwire Revolution High Performance Guidewire	2
► NaviPro™ Hydrophilic Guidewire	3
► NovaGold™ High Performance Guidewire	3
► WALLSTENT™ Super Stiff Guidewire	3

RX Preloaded Sphincterotomes

► Dreamtome™ RX Cannulating Sphincterotome (Preloaded with Dreamwire High Performance Guidewire) ...	4
► Hydratome™ RX Cannulating Sphincterotome (Preloaded with Hydra Jagwire High Performance Guidewire)	4
► Jagtome Revolution RX Cannulating Sphincterotome (Preloaded with Jagwire Revolution High Performance Guidewire)	5
► Jagtome™ RX Cannulating Sphincterotome (Preloaded with Jagwire High Performance Guidewire)	5

RX Sphincterotomes

► Autotome™ RX Cannulating Sphincterotome	6
► RX Cholangiogram Kits	6
► RX Needle Knife XL Triple-Lumen Needle Knife	6

Preloaded Sphincterotomes

► TRUEtome Dreamwire Cannulating Sphincterotome (Preloaded with Dreamwire High Performance Guidewire) ...	7
--	---

► TRUEtome Hydra Cannulating Sphincterotome (Preloaded with Hydra Jagwire High Performance Guidewire) ...	7
► TRUEtome™ Jag Cannulating Sphincterotome (Preloaded with Jagwire High Performance Guidewire)	7
► TRUEtome Revolution RX Cannulating Sphincterotome (Preloaded with Jagwire Revolution High Performance Guidewire)	7

Sphincterotomes

► MicroKnife™ XL Triple-Lumen Needle Knife	8
► Stonetome™ Stone Removal Device	9
► TRUEtome Cannulating Sphincterotome	8
► Ultratome™ Double-Lumen Sphincterotome	10
► Ultratome XL Triple-Lumen Sphincterotome	9

RX Cannulas

► Rapid Exchange XL Cannula	10
► RX ERCP Cannulas	10

Cannulas

► Contour™ ERCP Cannulas	11
► Tandem™ XL Triple-Lumen ERCP Cannula	11

Accessories

► Active Cords (Fits Boston Scientific Endoscopy Sphincterotome)	12
---	----

Rapid Refill Continuous Injection System


► Rapid Refill Continuous Injection System	12
--	----



GUIDEWIRES

Amplatz Super Stiff™ Guidewire

Order Number	GTIN	Description	Distal Tip O.D. (in)	Working Length (cm)	Tip	Packaging
M00550090	8714729301929	Amplatz Super Stiff Guidewire	.038	260	Straight	Each




Dreamwire™ High Performance Guidewire



Order Number	GTIN	Description	Distal Tip O.D. (in)	Working Length (cm)	Shaft Style	Dream Tip Style	Tip	Packaging
M00556101	08714729777816	Dreamwire Guidewire ST	.035	260	Standard	Straight	Round	Box 2
M00556111	08714729777847	Dreamwire Guidewire ANG	.035	260	Standard	Angled	Round	Box 2
M00556121	08714729777830	Dreamwire Guidewire ST Stiff	.035	260	Stiff	Straight	Round	Box 2
M00556131	08714729777823	Dreamwire Guidewire ANG Stiff	.035	260	Stiff	Angled	Round	Box 2
M00556141	08714729777885	Dreamwire Guidewire ST	.035	450	Standard	Straight	Round	Box 2
M00556151	08714729777878	Dreamwire Guidewire ANG	.035	450	Standard	Angled	Round	Box 2
M00556161	08714729777861	Dreamwire Guidewire ST Stiff	.035	450	Stiff	Straight	Round	Box 2
M00556171	08714729777854	Dreamwire Guidewire ANG Stiff	.035	450	Stiff	Angled	Round	Box 2

ST = Straight, ANG = Angled



Hydra Jagwire™ High Performance Guidewire



Order Number	GTIN	Description	Distal Tip O.D. (in)	Working Length (cm)	Shaft Style	Dream Tip Style	Tip	Packaging
M00556001	08714729648437	Hydra Jagwire Guidewire ST	.035	260	Standard	Straight	Round	Box 2
M00556011	08714729648444	Hydra Jagwire Guidewire ANG	.035	260	Standard	Angled	Round	Box 2
M00556021	08714729648451	Hydra Jagwire Guidewire ST Stiff	.035	260	Stiff	Straight	Round	Box 2
M00556031	08714729648468	Hydra Jagwire Guidewire ANG Stiff	.035	260	Stiff	Angled	Round	Box 2
M00556041	08714729648475	Hydra Jagwire Guidewire ST	.035	450	Standard	Straight	Round	Box 2
M00556051	08714729648482	Hydra Jagwire Guidewire ANG	.035	450	Standard	Angled	Round	Box 2
M00556061	08714729648499	Hydra Jagwire Guidewire ST Stiff	.035	450	Stiff	Straight	Round	Box 2
M00556071	08714729648505	Hydra Jagwire Guidewire ANG Stiff	.035	450	Stiff	Angled	Round	Box 2

Hydra Jagwire Guidewire may remain in place during sphincterotomy. See Instructions for Use for details.

ST = Straight, ANG = Angled



Section Table of Contents



Prev Page

1


Next Page



Home




GUIDEWIRES (Continued)



Jagwire™ High Performance Guidewire								
Order Number	GTIN	Description	Distal Tip O.D. (in)	Working Length (cm)	Shaft Style	Tip Style	Tip	Packaging
M00556511	08714729874492	Jagwire Guidewire ANG	.025	450	Standard	Angled	Round	Box 2
M00556521	08714729874508	Jagwire Guidewire ST Stiff	.025	450	Stiff	Straight	Round	Box 2
M00556531	08714729874515	Jagwire Guidewire ANG Stiff	.025	450	Stiff	Angled	Round	Box 2
M00556561	08714729501374	Jagwire Guidewire ST	.025	450	Standard	Straight	Trim	Box 2
M00556581	08714729339151	Jagwire Guidewire ST	.035	450	Standard	Straight	Trim	Box 2
M00556591	08714729501398	Jagwire Guidewire ANG	.035	450	Standard	Angled	Trim	Box 2
M00556601	08714729172154	Jagwire Guidewire ST Stiff	.035	450	Stiff	Straight	Trim	Box 2
M00556611	08714729172161	Jagwire Guidewire ANG Stiff	.035	450	Stiff	Angled	Trim	Box 2
M00556621	08714729172178	Jagwire Guidewire ST	.038	260	Standard	Straight	Trim	Box 2
M00556631	08714729172185	Jagwire Guidewire ST Stiff	.035	260	Stiff	Straight	Trim	Box 2
M00556641	08714729172192	Jagwire Guidewire ST	.035	260	Standard	Straight	Trim	Box 2
M00556651	08714729172208	Jagwire Guidewire ANG	.035	260	Standard	Angled	Trim	Box 2
M00556661	08714729172215	Jagwire Guidewire ANG Stiff	.035	260	Stiff	Angled	Trim	Box 2
M00556461	08714729415183	Jagwire Guidewire ST	.025	260	Standard	Straight	Round	Box 2
M00556481	08714729415190	Jagwire Guidewire ST Stiff	.025	260	Stiff	Straight	Round	Box 2
M00556491	08714729415206	Jagwire Guidewire ANG Stiff	.025	260	Stiff	Angled	Round	Box 2
M00556841	08714729268680	Extendable Jagwire Guidewire ST	.035	260	Standard	Straight	Round	Box 2
M00556851	08714729268697	Extendable Jagwire Guidewire ANG	.035	260	Standard	Angled	Round	Box 2
M00556861	08714729268703	Extendable Jagwire Guidewire ST Stiff	.035	260	Stiff	Straight	Round	Box 2
M00556871	08714729268710	Extendable Jagwire Guidewire ANG Stiff	.035	260	Stiff	Angled	Round	Box 2
M00556902	08714729268734	Jagtail™ Guidewire Extension	.035	200	—	Guidewire extension	—	Box 5


Jagwire Guidewire may remain in place during sphincterotomy. See Instructions for Use for details.
ST = Straight, ANG = Angled



Jagwire Revolution High Performance Guidewire							
Order Number	GTIN	Description	Distal Tip O.D. (in)	Working Length (cm)	Tip Style	Tip	Packaging
M00557001	08714729934974	Jagwire Revolution High Performance Guidewire- ST	.025	260	Straight	Round	Box 2
M00557011	08714729934981	Jagwire Revolution High Performance Guidewire- ANG	.025	260	Angled	Round	Box 2
M00557021	08714729934998	Jagwire Revolution High Performance Guidewire- ST	.025	450	Straight	Round	Box 2
M00557031	08714729935001	Jagwire Revolution High Performance Guidewire- ANG	.025	450	Angled	Round	Box 2
M00557041	08714729935018	Jagwire Revolution High Performance Guidewire- ST	.025	500	Straight	Round	Box 2

Jagwire Guidewire may remain in place during sphincterotomy. See Instructions for Use for details.
ST = Straight, ANG = Angled


GUIDEWIRES (Continued)



NaviPro™ Hydrophilic Guidewire



Order Number	GTIN	Description	Guidewire Diameter (in)	Length (cm)	Shaft Style	Tip Style	Packaging
M00556191	08714729839873	NaviPro Hydrophilic Guidewire	0.018	260	Standard	Angled	Box 5
M00556181	08714729839866	NaviPro Hydrophilic Guidewire	0.018	260	Standard	Straight	Box 5
M00556211	08714729839897	NaviPro Hydrophilic Guidewire	0.025	260	Standard	Angled	Box 5
M00556201	08714729839880	NaviPro Hydrophilic Guidewire	0.025	260	Standard	Straight	Box 5
M00556231	08714729839910	NaviPro Hydrophilic Guidewire	0.035	260	Standard	Angled	Box 5
M00556221	08714729839903	NaviPro Hydrophilic Guidewire	0.035	260	Standard	Straight	Box 5
M00556251	08714729839934	NaviPro Hydrophilic Guidewire	0.035	260	Stiff	Angled	Box 5
M00556241	08714729839927	NaviPro Hydrophilic Guidewire	0.035	260	Stiff	Straight	Box 5



NovaGold™ High Performance Guidewire



Order Number	GTIN	Description	Distal Tip O.D. (in)	Working Length (cm)	Tip Style	Tip Length (cm)	Packaging
M00552001	08714729854418	NovaGold High Performance Guidewire	.018	260	Straight (shapeable)	6	Box 2
M00552011	08714729854425	NovaGold High Performance Guidewire	.018	480	Straight (shapeable)	6	Box 2

WALLSTENT™ Super Stiff Guidewire



Order Number	GTIN	Description	Distal Tip O.D. (in)	Working Length (cm)	Tip Style	Packaging
H965180010	8714729203063	WALLSTENT Super Stiff Guidewire	.035	500	Tapered	Each
H965180011	8714729203070	WALLSTENT Super Stiff Guidewire	.035	500	Tapered	Box 5

RX PRELOADED SPHINCTEROTOMES




Dreamtome™ RX Cannulating Sphincterotome (Preloaded with Dreamwire™ High Performance Guidewire)



Order Number	GTIN	Description	Tip Length (mm)	Cut Wire Length (mm)	Preloaded Guidewire (in / cm)	Guidewire Tip Style	Packaging
M00584000	08714729778226	Dreamtome RX 49 Sphincterotome	5	20	.035 / 260 Dreamwire Guidewire	Straight	Each
M00584010	08714729778233	Dreamtome RX 49 Sphincterotome	5	30	.035 / 260 Dreamwire Guidewire	Straight	Each
M00584020	08714729778240	Dreamtome RX 49 Sphincterotome	5	20	.035 / 450 Dreamwire Guidewire	Straight	Each
M00584030	08714729778257	Dreamtome RX 49 Sphincterotome	5	30	.035 / 450 Dreamwire Guidewire	Straight	Each
M00584040	08714729778264	Dreamtome RX 44 Sphincterotome	5	20	.035 / 260 Dreamwire Guidewire	Straight	Each
M00584050	08714729778271	Dreamtome RX 44 Sphincterotome	5	30	.035 / 260 Dreamwire Guidewire	Straight	Each
M00584060	08714729778288	Dreamtome RX 44 Sphincterotome	5	20	.035 / 450 Dreamwire Guidewire	Straight	Each
M00584070	08714729778295	Dreamtome RX 44 Sphincterotome	5	30	.035 / 450 Dreamwire Guidewire	Straight	Each

*Preloaded with straight-tip guidewires only.
Active Cords sold separately (See page 12).



Hydratome™ RX Cannulating Sphincterotome (Preloaded with Hydra Jagwire™ High Performance Guidewire)



Order Number	GTIN	Description	Cut Wire (mm)	Tip Length (mm)	Distal Tip O.D. (Fr)	Preloaded Guidewire (cm)	Guidewire Tip Style	Packaging
M00583000	08714729755159	Hydratome RX 49 Sphincterotome	20	5	4.9	260 Hydra Jagwire Guidewire	Straight	Each
M00583010	08714729755166	Hydratome RX 49 Sphincterotome	30	5	4.9	260 Hydra Jagwire Guidewire	Straight	Each
M00583020	08714729755173	Hydratome RX 49 Sphincterotome	20	5	4.9	450 Hydra Jagwire Guidewire	Straight	Each
M00583030	08714729755180	Hydratome RX 49 Sphincterotome	30	5	4.9	450 Hydra Jagwire Guidewire	Straight	Each
M00583040	08714729755197	Hydratome RX 44 Sphincterotome	20	5	4.4	260 Hydra Jagwire Guidewire	Straight	Each
M00583050	08714729755203	Hydratome RX 44 Sphincterotome	30	5	4.4	260 Hydra Jagwire Guidewire	Straight	Each
M00583060	08714729755210	Hydratome RX 44 Sphincterotome	20	5	4.4	450 Hydra Jagwire Guidewire	Straight	Each
M00583070	08714729755227	Hydratome RX 44 Sphincterotome	30	5	4.4	450 Hydra Jagwire Guidewire	Straight	Each

Active Cords sold separately (See page 12).



Section Table of Contents



Prev Page

4



Next Page



Home





RX PRELOADED SPHINCTEROTOMES (Continued)



Jagtome™ RX Cannulating Sphincterotome (Preloaded with Jagwire™ High Performance Guidewire)								
Order Number	GTIN	Description	Cut Wire (mm)	Tip Length (mm)	Distal Tip O.D. (Fr)	Preloaded Guidewire (in / cm)	Guidewire Tip Style	Packaging
M00573000	08714729776710	Jagtome RX 49 Sphincterotome	20	5	4.9	.035 / 260 Jagwire Guidewire	Straight	Each
M00573010	08714729776727	Jagtome RX 49 Sphincterotome	20	5	4.9	.035 / 450 Jagwire Guidewire	Straight	Each
M00573020	08714729776734	Jagtome RX 49 Sphincterotome	30	5	4.9	.035 / 260 Jagwire Guidewire	Straight	Each
M00573030	08714729776741	Jagtome RX 49 Sphincterotome	30	5	4.9	.035 / 450 Jagwire Guidewire	Straight	Each
M00573040	08714729776659	Jagtome RX 44 Sphincterotome	20	5	4.4	.035 / 260 Jagwire Guidewire	Straight	Each
M00573050	08714729776666	Jagtome RX 44 Sphincterotome	20	5	4.4	.035 / 450 Jagwire Guidewire	Straight	Each
M00573060	08714729776673	Jagtome RX 44 Sphincterotome	30	5	4.4	.035 / 260 Jagwire Guidewire	Straight	Each
M00573070	08714729776543	Jagtome RX 44 Sphincterotome	30	5	4.4	.035 / 450 Jagwire Guidewire	Straight	Each
M00573080	08714729776642	Jagtome RX 39 Sphincterotome	20	5	3.9	.025 / 260 Jagwire Guidewire	Straight	Each
M00573090	08714729776550	Jagtome RX 39 Sphincterotome	20	5	3.9	.025 / 450 Jagwire Guidewire	Straight	Each
M00573100	08714729776536	Jagtome RX 39 Sphincterotome	30	5	3.9	.025 / 260 Jagwire Guidewire	Straight	Each
M00573110	08714729776529	Jagtome RX 39 Sphincterotome	30	5	3.9	.025 / 450 Jagwire Guidewire	Straight	Each



Active Cords sold separately ► (See page 12).



Jagtome Revolution RX Cannulating Sphincterotome (Preloaded with Jagwire Revolution High Performance Guidewire)								
Order Number	GTIN	Description	Cut Wire (mm)	Tip Length (mm)	Distal Tip O.D. (Fr)	Preloaded Guidewire (in / cm)	Guidewire Tip Style	Packaging
M00584240	08714729953234	Jagtome Revolution RX 39 Sphincterotome	20	5	3.9	.025/260 Jagwire Revolution Guidewire	Straight	Each
M00584250	08714729953241	Jagtome Revolution RX 39 Sphincterotome	20	5	3.9	.025/450 Jagwire Revolution Guidewire	Straight	Each
M00584260	08714729953258	Jagtome Revolution RX 39 Sphincterotome	30	5	3.9	.025/260 Jagwire Revolution Guidewire	Straight	Each
M00584270	08714729953265	Jagtome Revolution RX 39 Sphincterotome	30	5	3.9	.025/450 Jagwire Revolution Guidewire	Straight	Each

Active Cords sold separately ► (See page 12).

RX SPHINCTEROTOMES



Autotome™ RX Cannulating Sphincterotome							
Order Number	GTIN	Description	Cut Wire (mm)	Tip Length (mm)	Distal Tip O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00545150	08714729415138	Autotome RX 49 Cannulating Sphincterotome	20	5	4.9	.035	Each
M00545160	08714729415152	Autotome RX 49 Cannulating Sphincterotome	30	5	4.9	.035	Each
M00545170	08714729444749	Autotome RX 44 Cannulating Sphincterotome	20	5	4.4	.035	Each
M00545180	08714729444756	Autotome RX 44 Cannulating Sphincterotome	30	5	4.4	.035	Each
M00545190	08714729444763	Autotome RX 39 Cannulating Sphincterotome	20	5	3.9	.025	Each
M00545200	08714729444770	Autotome RX 39 Cannulating Sphincterotome	30	5	3.9	.025	Each

Active Cords sold separately (See page 12).

RX Cholangiogram Kits*				
Order Number	Description	Sphincterotome	Guidewire	Accessory
M005XXXX1	Cholangiogram Kit	Autotome or Ultratome™ XL Sphincterotome	Jagwire™ Guidewire	N / A
M005XXXX1	Cholangiogram Kit-Preloaded	Autotome Sphincterotome	Jagwire Guidewire	RX Locking Device

*To order RX Cholangiogram Kits please contact your local Boston Scientific Endoscopy Sales Representative.

RX Needle Knife XL Triple-Lumen Needle Knife							
Order Number	GTIN	Description	Length (cm)	Cut Wire (mm)	Distal Tip O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00545840	08714729283799	RX Needle Knife	200	5	5.5	.035	Each

PRELOADED SPHINCTEROTOMES

TRUEtome™ Jag Cannulating Sphincterotome (Preloaded with Jagwire™ High Performance Guidewire) 							
Order Number	GTIN	Description	Tip Length (mm)	Cut Wire Length (mm)	Preloaded Guidewire (in / cm)	Tip O.D. (Fr)	Packaging
M00583080	08714729858621	TRUEtome 39 Sphincterotome	5	20	.025 / 450 Jagwire Guidewire	3.9	Each
M00583090	08714729858638	TRUEtome 44 Sphincterotome	5	20	.025 / 450 Jagwire Guidewire	4.4	Each
M00583100	08714729858645	TRUEtome 44 Sphincterotome	5	20	.035 / 450 Jagwire Guidewire	4.4	Each
M00583110	08714729858652	TRUEtome 44 Sphincterotome	5	20	.035 / Stiff 450 Jagwire Guidewire	4.4	Each
M00583120	08714729858669	TRUEtome 39 Sphincterotome	5	30	.025 / 450 Jagwire Guidewire	3.9	Each
M00583130	08714729858676	TRUEtome 44 Sphincterotome	5	30	.025 / 450 Jagwire Guidewire	4.4	Each
M00583140	08714729858683	TRUEtome 44 Sphincterotome	5	30	.035 / 450 Jagwire Guidewire	4.4	Each
M00583150	08714729858690	TRUEtome 44 Sphincterotome	5	30	.035 / Stiff 450 Jagwire Guidewire	4.4	Each
TRUEtome Hydra Cannulating Sphincterotome (Preloaded with Hydra Jagwire™ High Performance Guidewire) 							
M00583160	08714729858706	TRUEtome 44 Sphincterotome	5	20	.035 / 450 Hydra Jagwire Guidewire	4.4	Each
M00583170	08714729858713	TRUEtome 44 Sphincterotome	5	30	.035 / 450 Hydra Jagwire Guidewire	4.4	Each
TRUEtome Dreamwire™ Cannulating Sphincterotome (Preloaded with Dreamwire High Performance Guidewire) 							
M00583180	08714729858720	TRUEtome 44 Sphincterotome	5	20	.035 / 260 Dreamwire Guidewire	4.4	Each
M00583190	08714729858737	TRUEtome 44 Sphincterotome	5	20	.035 / 450 Dreamwire Guidewire	4.4	Each
M00583200	08714729858744	TRUEtome 44 Sphincterotome	5	30	.035 / 260 Dreamwire Guidewire	4.4	Each
M00583210	08714729858751	TRUEtome 44 Sphincterotome	5	30	.035 / 450 Dreamwire Guidewire	4.4	Each

*Preloaded with straight-tip guidewires only.
Active Cords sold separately  (See page 12).

TRUEtome Revolution RX Cannulating Sphincterotome (Preloaded with Jagwire Revolution High Performance Guidewire) 							
Order Number	GTIN	Description	Tip Length (mm)	Cut Wire Length (mm)	Preloaded Guidewire (in / cm)	Tip O.D. (Fr)	Packaging
M00583260	08714729953272	TRUEtome 39 Sphincterotome	5	20	.025/450 Jagwire Revolution Guidewire	3.9	Each
M00583270	08714729953289	TRUEtome 44 Sphincterotome	5	20	.025/450 Jagwire Revolution Guidewire	4.4	Each
M00583280	08714729953296	TRUEtome 39 Sphincterotome	5	30	.025/450 Jagwire Revolution Guidewire	3.9	Each
M00583290	08714729953302	TRUEtome 44 Sphincterotome	5	30	.025/450 Jagwire Revolution Guidewire	4.4	Each


*Preloaded with straight-tip guidewires only.
Active Cords sold separately  (See page 12).

SPHINCTEROTOMES



TRUEtome™ Cannulating Sphincterotome							
Order Number	GTIN	Description	Cut Wire Length (mm)	Tip Length (mm)	Distal Tip O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00584190	08714729827658	TRUEtome 39 Sphincterotome	20	5	3.9	.025	Each
M00584120	08714729827771	TRUEtome 39 Sphincterotome	30	5	3.9	.025	Each
M00584170	08714729827634	TRUEtome 44 Sphincterotome	20	5	4.4	.035	Each
M00584180	08714729827641	TRUEtome 44 Sphincterotome	30	5	4.4	.035	Each
M00584150	08714729827610	TRUEtome 49 Sphincterotome	20	5	4.9	.035	Each
M00584160	08714729827627	TRUEtome 49 Sphincterotome	30	5	4.9	.035	Each


Active Cords sold separately (See page 12).



MicroKnife™ XL Triple-Lumen Needle Knife						
Order Number	GTIN	Description	Shaft O.D. (Fr)	Catheter Length (cm)	Recommended Guidewire (in)	Packaging
M00532810	08714729172017	MicroKnife XL Triple-Lumen Needle Knife	7.0 - 5.5	200	.035	Each

Active Cords sold separately (See page 12).

SPHINCTEROTOMES (Continued)




Stonetome™ Stone Removal Device



Order Number	GTIN	Description	Balloon Diameter (mm)	Catheter Length (cm)	Tip Length (mm)	Cut Wire (mm)	Distal Tip O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00535150	08714729146667	Stonetome Device (below cut wire)	11.5	200	5	20	5.5	.035	Each
M00535170	08714729146629	Stonetome Device (below cut wire)	11.5	200	20	20	5.5	.035	Each
M00535190	08714729146612	Stonetome Device (below cut wire)	11.5	200	5	30	5.5	.035	Each
M00535210	08714729146605	Stonetome Device (below cut wire)	11.5	200	20	30	5.5	.035	Each
M00535110	08714729146599	Stonetome Device (above cut wire)	11.5	200	20	20	5.5	.035	Each
M00535130	08714729146575	Stonetome Device (above cut wire)	11.5	200	20	30	5.5	.035	Each

Active Cords sold separately ► (See page 12).




Ultratome™ XL Triple-Lumen Sphincterotome



Order Number	GTIN	Description	Wire Length (mm)	Shaft O.D. (Fr)	Catheter Length (cm)	Tip Length (mm)	Recommended Guidewire (in)	Packaging
M00535900	08714729103257	Ultratome XL Short Nose Sphincterotome	20	5.5	200	5	.035	Each
M00535910	08714729103301	Ultratome XL Long Nose Sphincterotome	20	5.5	200	20	.035	Each
M00535920	08714729103318	Ultratome XL Short Nose Sphincterotome	30	5.5	200	5	.035	Each
M00535930	08714729103349	Ultratome XL Long Nose Sphincterotome	30	5.5	200	20	.035	Each

Active Cords sold separately ► (See page 12).

SPHINCTEROTOMES (Continued)




Ultratome™ Double-Lumen Sphincterotome



Order Number	GTIN	Description	Wire Length (mm)	Distal Tip O.D.		Catheter Length (cm)	Tip Length (mm)	Recommended Guidewire (in)	Packaging
				(mm)	(Fr)				
M00530800	08714729163664	Ultratome Short Nose Sphincterotome	20	1.8	5.5	200	5	.035	Each
M00530810	08714729183143	Ultratome Long Nose Sphincterotome	20	1.8	5.5	200	20	.035	Each
M00530820	08714729193418	Ultratome Short Nose Sphincterotome	30	1.8	5.5	200	5	.035	Each
M00530830	08714729193425	Ultratome Long Nose Sphincterotome	30	1.8	5.5	200	20	.035	Each

Active Cords sold separately (See page 12).


RX CANNULAS



Rapid Exchange XL Cannula



Order Number	GTIN	Description	Working Length (cm)	Tip Style	Recommended Guidewire (in)	Packaging
M00545820	08714729301905	Rapid Exchange XL Cannula	195	Tapered	.035	Each




RX ERCP Cannulas



Order Number	GTIN	Description	Working Length (cm)	Recommended Guidewire (in)	Packaging
M00545760	08714729268642	RX ERCP Cannula Standard Tip	210	.035	Each
M00545770	08714729268659	RX ERCP Cannula Tapered Tip	210	.035	Each
M00545780	08714729268673	RX ERCP Cannula Ball Tip	210	.035	Each


CANNULAS



Tandem™ XL Triple-Lumen ERCP Cannula



Order Number	GTIN	Distal Tip O.D.		Catheter Size (Fr)	Working Length (cm)	Tip Style	Recommended Guidewire (in)	Packaging
		(Fr)	(mm)					
M00535700	08714729146636	5.5	1.8	7.0 – 5.5	200	Tapered	.035	Each



Contour™ ERCP Cannulas



Order Number	GTIN	Description	Distal Tip O.D.		Working Length (cm)	Tip Style	Recommended Guidewire (in)	Packaging
			(Fr)	(mm)				
M00530850	08714729099789	Contour ERCP Cannula	5	1.7	210	Standard	.035	Each
M00530860	08714729099819	Contour ERCP Cannula	5	1.7	210	Tapered	.035	Each
M00530870	08714729099826	Contour ERCP Cannula	5	1.7	210	Ultra Tapered	.025	Each
M00530880	08714729099833	Contour ERCP Cannula	5	1.7	210	5-4-3	.018	Each
M00530890	08714729099840	Contour ERCP Cannula	5	1.7	210	Ball Tip	.035	Each

ACCESSORIES


RX Locking Devices and Biopsy Cap Systems



Order Number	GTIN	Description	Packaging
M00545211	08714729440604	RX Locking Device Olympus™ / FujiFilm™ Brand Endoscope Compatible	Box 10
M00545261	08714729296614	RX Locking Device and Biopsy Cap System – Olympus / FujiFilm G5 Series Brand Endoscope Compatible	Box 10
M00545271	08714729296621	RX Locking Device and Biopsy Cap System – Pentax™ Brand Endoscope Compatible	Box 10
M00545281	08714729296638	RX Locking Device and Biopsy Cap System – FujiFilm Type 3 Series Brand Endoscope Compatible	Box 10

Active Cords (Fit Boston Scientific Endoscopy Sphincterotome)

Order Number	GTIN	Description	Packaging
M00561250	08714729019060	Active Cord, Fits Olympus Brand Generator	Each
M00561270	08714729053361	Active Cord, Fits Endostat™, Endostat II and Endostat III Generators, Bovie™, Valleylab™, Aspen Labs™, and ERBE™ Brand Generators	Each



Rapid Refill Continuous Injection System



Order Number	GTIN	Description	Packaging
M00566001	08714729285151	Rapid Refill Continuous Injection System	Box 10

Notes and Future Product Updates

Notes and Future Product Updates



Section Table of Contents



Prev Page

Next Page



Home



Biliary Stricture and Stone Management



Biliary and Pancreatic Stents

Biliary Metal Stents

➤ Epic™ Biliary Endoscopic Stent System.....	14
➤ WallFlex™ Biliary RX Fully Covered Stent System RMV.....	13
➤ WallFlex Biliary RX Partially Covered Stent System.....	13
➤ WallFlex Biliary RX Uncovered Stent System	13
➤ WALLSTENT™ Biliary Endoprosthesis Stents	16
➤ WALLSTENT RX Biliary Endoprosthesis Stent System.....	15

Biliary Plastic Stents

➤ Advanix™ Center Bend Preloaded RX Stent Systems	17
➤ Advanix Center Bend - Single Stent.....	18
➤ Advanix Double Pigtail - Single Stent.....	19
➤ Advanix Duodenal Bend Preloaded RX Stent Systems.....	17
➤ Advanix Duodenal Bend - Single Stent	18
➤ Flexima™ Duodenal Bend Biliary Stent with Delivery System ..20	
➤ Flexima Single Pigtail Nasobiliary Catheter	21
➤ NaviFlex™ RX Delivery System.....	19
➤ RX Duodenal Bend Biliary Stent with RX Delivery System ...20	

Pancreatic Plastic Stents

➤ Advanix Pancreatic Stent Kit - Straight, No Leading Barb and Leading Barb.....	22
➤ Advanix Pancreatic Stent Kit - Pigtail, Leading Barb and No Leading Barb	22
➤ Advanix Pancreatic Stent - Straight, No Leading Barb	23
➤ Advanix Pancreatic Stent - Straight, Leading Barb	24
➤ Advanix Pancreatic Stent - Pigtail, No Leading Barb	25
➤ Advanix Pancreatic Stent - Pigtail, Leading Barb.....	26
➤ NaviFlex RX Pancreatic Delivery System & NaviFlex Pushers....	27

Biliary Stone Management

RX Stone Management

➤ Extractor™ Pro RX Retrieval Balloon Catheter – Distal Injection ...28	
➤ Extractor Pro RX Retrieval Balloon Catheter – Proximal Injection..28	
➤ Extractor Pro RX-S Retrieval Balloon Catheter – Distal Injection ..28	
➤ Extractor Pro RX-S Retrieval Balloon Catheter – Proximal Injection..28	

Stone Management

➤ Extractor Pro XL Retrieval Balloon Catheter – Distal Injection...29	
➤ Extractor Pro XL Retrieval Balloon Catheter – Proximal Injection..29	
➤ Trapezoid™ RX Wireguided Retrieval Baskets	29

CRE™ Wireguided Balloon Dilatation Catheters

➤ CRE PRO Wireguided Balloon Dilatation Catheters	30
➤ CRE RX Biliary Balloon Dilatation Catheters.....	30

Biliary Stricture Management

Pancreaticobiliary Ablation

➤ Habib™ EndoHPB Bipolar Radiofrequency Catheter	31
--	----

RX Cytology

➤ RX Cytology Brush Wireguided Cytology Brush.....	32
--	----

RX Dilatation

➤ Hurricane™ RX Biliary Balloon Dilatation Catheter	31
---	----

Biliary Dilatation

➤ MaxForce™ Biliary Balloon Dilatation Catheter.....	32
--	----

RX Accessories

➤ Encore™ 26 Inflation Device.....	32
------------------------------------	----

Biliary Metal Stents

WALLFLEX™ BILIARY RX STENTS



Fully Covered



Magnetic Resonance Conditional¹

WallFlex Biliary RX Fully Covered Stent System RMV*



Order Number	GTIN	Diameter (mm)	Length (mm)	Covered Length (mm) (PC Only)	Catheter Diameter (Fr / mm)	Guidewire Diameter (in / mm)	Working Length (cm)	Packaging
M00570340	08714729953159	8	60	—	8.5 / 2.83	.035 / .89	194	Each
M00570350	08714729953166	8	80	—	8.5 / 2.83	.035 / .89	194	Each
M00570360	08714729953173	10	40	—	8.5 / 2.83	.035 / .89	194	Each
M00570370	08714729953180	10	60	—	8.5 / 2.83	.035 / .89	194	Each
M00570380	08714729953197	10	80	—	8.5 / 2.83	.035 / .89	194	Each


*Compatible with Boston Scientific .035" guidewires, including Dreamwire™ High Performance Guidewire (See page 1).

Above (WallFlex Biliary RX Fully Covered Stent System) with Permalume™ Covering.

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Partially Covered



Magnetic Resonance Conditional¹

WallFlex Biliary RX Partially Covered Stent System*



Order Number	GTIN	Diameter (mm)	Length (mm)	Covered Length (mm) (PC Only)	Catheter Diameter (Fr / mm)	Guidewire Diameter (in / mm)	Working Length (cm)	Packaging
M00570700	08714729759294	8	60	48	8.5 / 2.83	.035 / .89	194	Each
M00570710	08714729759300	8	80	68	8.5 / 2.83	.035 / .89	194	Each
M00570720	08714729759317	10	40	28	8.5 / 2.83	.035 / .89	194	Each
M00570730	08714729759324	10	60	48	8.5 / 2.83	.035 / .89	194	Each
M00570740	08714729759331	10	80	68	8.5 / 2.83	.035 / .89	194	Each


*Compatible with Boston Scientific .035" guidewires, including Dreamwire High Performance Guidewire (See page 1).

Above (WallFlex Biliary RX Partially Covered Stent System) with Permalume Covering.

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Uncovered



Magnetic Resonance Conditional¹

WallFlex Biliary RX Uncovered Stent System*



Order Number	GTIN	Diameter (mm)	Length (mm)	Covered Length (mm) (PC Only)	Catheter Diameter (Fr / mm)	Guidewire Diameter (in / mm)	Working Length (cm)	Packaging
M00570600	08714729759249	8	40	—	8.0 / 2.67	.035 / .89	194	Each
M00570610	08714729759256	8	60	—	8.0 / 2.67	.035 / .89	194	Each
M00570620	08714729759263	8	80	—	8.0 / 2.67	.035 / .89	194	Each
M00570630	08714729759270	8	100	—	8.0 / 2.67	.035 / .89	194	Each
M00570890	08714729765332	10	40	—	8.0 / 2.67	.035 / .89	194	Each
M00570640	08714729756231	10	60	—	8.0 / 2.67	.035 / .89	194	Each
M00570650	08714729758150	10	80	—	8.0 / 2.67	.035 / .89	194	Each
M00570660	08714729759287	10	100	—	8.0 / 2.67	.035 / .89	194	Each

*Compatible with Boston Scientific .035" guidewires, including Dreamwire High Performance Guidewire (See page 1).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Biliary Metal Stents

EPIC™ BILIARY ENDOSCOPIC STENT SYSTEM



 Magnetic Resonance Conditional¹

Epic Biliary Endoscopic Stent

Order Number	GTIN	Description	Stent Diameter (mm)	Stent Length (mm)	Catheter Diameter (Fr)	Guidewire Diameter (in)	Catheter Working Length (cm)	Packaging
M00572230	08714729896852	Epic Biliary Stent	6	40	6	0.035	220	Each
M00572240	08714729896869	Epic Biliary Stent	6	60	6	0.035	220	Each
M00572250	08714729896876	Epic Biliary Stent	6	80	6	0.035	220	Each
M00572260	08714729896883	Epic Biliary Stent	6	100	6	0.035	220	Each
M00572270	08714729896890	Epic Biliary Stent	8	40	6	0.035	220	Each
M00572280	08714729896906	Epic Biliary Stent	8	60	6	0.035	220	Each
M00572290	08714729896913	Epic Biliary Stent	8	80	6	0.035	220	Each
M00572300	08714729896920	Epic Biliary Stent	8	100	6	0.035	220	Each
M00572310	08714729896937	Epic Biliary Stent	10	40	6	0.035	220	Each
M00572320	08714729896944	Epic Biliary Stent	10	60	6	0.035	220	Each
M00572330	08714729896951	Epic Biliary Stent	10	80	6	0.035	220	Each
M00572340	08714729896968	Epic Biliary Stent	10	100	6	0.035	220	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Biliary Metal Stents

WALLSTENT™ RX BILIARY STENTS



 Magnetic Resonance Conditional¹

WALLSTENT RX Biliary Endoprosthesis Stent System*								
Order Number	GTIN	Description	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System (Fr)	Cover Length (mm)	Working Length (cm)	Packaging
M00569670	08714729296874	WALLSTENT RX Biliary Partially Covered Stent System	8	40	8	30	194	Each
M00569680	08714729296881	WALLSTENT RX Biliary Partially Covered Stent System	8	60	8	50	194	Each
M00569690	08714729296898	WALLSTENT RX Biliary Partially Covered Stent System	8	80	8	70	194	Each
M00569700	08714729296904	WALLSTENT RX Biliary Partially Covered Stent System	10	40	8	30	194	Each
M00569710	08714729296911	WALLSTENT RX Biliary Partially Covered Stent System	10	60	8	50	194	Each
M00569720	08714729296928	WALLSTENT RX Biliary Partially Covered Stent System	10	80	8	70	194	Each

*Recommended Guidewire: .035", 260cm, including Dreamwire™ High Performance Guidewire  (See page 1).

Above (WALLSTENT RX Biliary Endoprosthesis Stent System) with Permalume™ Covering.

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹


WALLSTENT RX Biliary Endoprosthesis Stent System*								
Order Number	GTIN	Description	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System (Fr)	Working Length (cm)	Packaging	
M00569610	08714729296812	WALLSTENT RX Biliary Uncovered Stent System	8	40	8	194	Each	
M00569620	08714729296829	WALLSTENT RX Biliary Uncovered Stent System	8	60	8	194	Each	
M00569630	08714729296836	WALLSTENT RX Biliary Uncovered Stent System	8	80	8	194	Each	
M00569800	08714729456711	WALLSTENT RX Biliary Uncovered Stent System	8	100	8	194	Each	
M00569640	08714729296843	WALLSTENT RX Biliary Uncovered Stent System	10	40	8	194	Each	
M00569650	08714729296850	WALLSTENT RX Biliary Uncovered Stent System	10	60	8	194	Each	
M00569660	08714729296867	WALLSTENT RX Biliary Uncovered Stent System	10	80	8	194	Each	
M00569810	08714729456728	WALLSTENT RX Biliary Uncovered Stent System	10	100	8	194	Each	

*Recommended Guidewire: .035", 260cm, including Dreamwire High Performance Guidewire  (See page 1).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Biliary Metal Stents

WALLSTENT™ STENTS




WALLSTENT Biliary Endoscopic Endoprosthesis Stents*								
Order Number	GTIN	Description	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System (Fr)	Cover Length (mm)	Working Length (cm)	Packaging
H965430700	08714729229858	WALLSTENT Biliary Partially Covered Stent	8	40	8	30	194	Each
H965430800	08714729317906	WALLSTENT Biliary Partially Covered Stent	8	60	8	50	194	Each
H965430900	08714729229872	WALLSTENT Biliary Partially Covered Stent	8	80	8	70	194	Each
H965431000	08714729317920	WALLSTENT Biliary Partially Covered Stent	10	40	8	30	194	Each
H965431100	08714729317937	WALLSTENT Biliary Partially Covered Stent	10	60	8	50	194	Each
H965431200	08714729229902	WALLSTENT Biliary Partially Covered Stent	10	80	8	70	194	Each

*Recommended Guidewire: .035", 450cm, including Dreamwire™ High Performance Guidewire (See page 1).

Above (WALLSTENT Biliary Endoprosthesis) with Permalume™ Covering.

†Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



WALLSTENT Biliary Endoscopic Endoprosthesis Stents*							
Order Number	GTIN	Description	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System Diameter (Fr)	Working Length (cm)	Packaging
H965430100	08714729204633	WALLSTENT Biliary Uncovered Stent	8	40	7.5	194	Each
H965430200	08714729204640	WALLSTENT Biliary Uncovered Stent	8	60	7.5	194	Each
H965430300	08714729204657	WALLSTENT Biliary Uncovered Stent	8	80	7.5	194	Each
H965430400	08714729204664	WALLSTENT Biliary Uncovered Stent	10	40	7.5	194	Each
H965430500	08714729204671	WALLSTENT Biliary Uncovered Stent	10	60	7.5	194	Each
H965430600	08714729204688	WALLSTENT Biliary Uncovered Stent	10	80	7.5	194	Each

*Recommended Guidewire: .035", 450cm, including Dreamwire High Performance Guidewire (See page 1).

†Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Biliary Plastic Stents

ADVANIX™ BILIARY STENTS



MR Magnetic Resonance Safe

Advanix Duodenal Bend Preloaded RX Stent Systems					
Order Number	GTIN	Description	Stent Diameter (Fr)	Length (cm) Barb to Barb	Packaging
M00534200	08714729787297	Duodenal Bend	7	5	Each
M00534210	08714729787303	Duodenal Bend	7	7	Each
M00534220	08714729787310	Duodenal Bend	7	9	Each
M00534230	08714729787327	Duodenal Bend	7	12	Each
M00534240	08714729787334	Duodenal Bend	7	15	Each
M00534250	08714729787341	Duodenal Bend	7	18	Each
M00534260	08714729787358	Duodenal Bend	8.5	5	Each
M00534270	08714729787365	Duodenal Bend	8.5	7	Each
M00534280	08714729787372	Duodenal Bend	8.5	9	Each
M00534290	08714729787389	Duodenal Bend	8.5	12	Each
M00534300	08714729787396	Duodenal Bend	8.5	15	Each
M00534310	08714729787402	Duodenal Bend	8.5	18	Each
M00534320	08714729787419	Duodenal Bend	10	5	Each
M00534330	08714729787426	Duodenal Bend	10	7	Each
M00534340	08714729787433	Duodenal Bend	10	9	Each
M00534350	08714729787440	Duodenal Bend	10	12	Each
M00534360	08714729787457	Duodenal Bend	10	15	Each
M00534370	08714729787464	Duodenal Bend	10	18	Each

Scope Compatibility: 7Fr and 8.5Fr stents require a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.



MR Magnetic Resonance Safe

Advanix Center Bend Preloaded RX Stent Systems					
Order Number	GTIN	Description	Stent Diameter (Fr)	Length (cm) Barb to Barb	Packaging
M00534560	08714729787112	Center Bend	7	5	Each
M00534570	08714729787129	Center Bend	7	7	Each
M00534580	08714729787136	Center Bend	7	9	Each
M00534590	08714729787143	Center Bend	7	12	Each
M00534600	08714729787150	Center Bend	7	15	Each
M00534610	08714729787167	Center Bend	7	18	Each
M00534620	08714729787174	Center Bend	8.5	5	Each
M00534630	08714729787181	Center Bend	8.5	7	Each
M00534640	08714729787198	Center Bend	8.5	9	Each
M00534650	08714729787204	Center Bend	8.5	12	Each
M00534660	08714729787211	Center Bend	8.5	15	Each
M00534670	08714729787228	Center Bend	8.5	18	Each
M00534680	08714729787235	Center Bend	10	5	Each
M00534690	08714729787242	Center Bend	10	7	Each
M00534700	08714729787259	Center Bend	10	9	Each
M00534710	08714729787266	Center Bend	10	12	Each
M00534720	08714729787273	Center Bend	10	15	Each
M00534730	08714729787280	Center Bend	10	18	Each

Scope Compatibility: 7Fr and 8.5Fr stents require a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.

Biliary Plastic Stents

ADVANIX™ BILIARY STENTS (Continued)



MR Magnetic Resonance Safe

Advanix Duodenal Bend - Single Stent						
Order Number	GTIN	Description	Stent Diameter (Fr)	Length (cm) Barb to Barb	Packaging	
M00532870	08714729786788	Duodenal Bend	7	5	Each	
M00532880	08714729786795	Duodenal Bend	7	7	Each	
M00532890	08714729786801	Duodenal Bend	7	9	Each	
M00532900	08714729786818	Duodenal Bend	7	12	Each	
M00532910	08714729786825	Duodenal Bend	7	15	Each	
M00532920	08714729786832	Duodenal Bend	7	18	Each	
M00532930	08714729786849	Duodenal Bend	8.5	5	Each	
M00532940	08714729786856	Duodenal Bend	8.5	7	Each	
M00532950	08714729786863	Duodenal Bend	8.5	9	Each	
M00532960	08714729786870	Duodenal Bend	8.5	12	Each	
M00532970	08714729786887	Duodenal Bend	8.5	15	Each	
M00532980	08714729786894	Duodenal Bend	8.5	18	Each	
M00532990	08714729786900	Duodenal Bend	10	5	Each	
M00533000	08714729786917	Duodenal Bend	10	7	Each	
M00533010	08714729786924	Duodenal Bend	10	9	Each	
M00533020	08714729786931	Duodenal Bend	10	12	Each	
M00533030	08714729786948	Duodenal Bend	10	15	Each	
M00533040	08714729786955	Duodenal Bend	10	18	Each	

Scope Compatibility: 7Fr and 8.5Fr stents require a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.




MR Magnetic Resonance Safe

Advanix Center Bend - Single Stent						
Order Number	GTIN	Description	Stent Diameter (Fr)	Length (cm) Barb to Barb	Packaging	
M00533230	08714729786603	Center Bend	7	5	Each	
M00533240	08714729786610	Center Bend	7	7	Each	
M00533250	08714729786627	Center Bend	7	9	Each	
M00533260	08714729786634	Center Bend	7	12	Each	
M00533270	08714729786641	Center Bend	7	15	Each	
M00533280	08714729786658	Center Bend	7	18	Each	
M00533290	08714729786665	Center Bend	8.5	5	Each	
M00533300	08714729786672	Center Bend	8.5	7	Each	
M00533310	08714729786689	Center Bend	8.5	9	Each	
M00533320	08714729786696	Center Bend	8.5	12	Each	
M00533330	08714729786702	Center Bend	8.5	15	Each	
M00533340	08714729786719	Center Bend	8.5	18	Each	
M00533350	08714729786726	Center Bend	10	5	Each	
M00533360	08714729786733	Center Bend	10	7	Each	
M00533370	08714729786740	Center Bend	10	9	Each	
M00533380	08714729786757	Center Bend	10	12	Each	
M00533390	08714729786764	Center Bend	10	15	Each	
M00533400	08714729786771	Center Bend	10	18	Each	

Scope Compatibility: 7Fr and 8.5Fr stents require a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.

► Biliary Plastic Stents

ADVANIX™ BILIARY STENTS (Continued)



MR Magnetic Resonance Safe

Advanix Double Pigtail - Single Stent					
Order Number	GTIN	Description	Stent Diameter (Fr)	Length (cm) Between Pigtails	Packaging
M00532160	08714729786962	Double Pigtail	7	3	Each
M00532180	08714729786979	Double Pigtail	7	5	Each
M00532190	08714729786986	Double Pigtail	7	7	Each
M00532210	08714729786993	Double Pigtail	7	10	Each
M00532220	08714729787006	Double Pigtail	7	12	Each
M00532230	08714729787013	Double Pigtail	7	15	Each
M00532240	08714729787020	Double Pigtail	10	3	Each
M00532260	08714729787037	Double Pigtail	10	5	Each
M00532270	08714729787044	Double Pigtail	10	7	Each
M00532290	08714729787051	Double Pigtail	10	10	Each
M00532300	08714729787068	Double Pigtail	10	12	Each
M00532310	08714729787075	Double Pigtail	10	15	Each

Scope Compatibility: 7Fr and 8.5Fr stents require a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.

NAVIFLEX™ RX DELIVERY SYSTEM



NaviFlex RX Delivery System*				
Order Number	GTIN	Description	Catheter Diameter (Fr)	Length (cm)
M00533560	08714729787082	RX Delivery System	7	202.5
M00533570	08714729787099	RX Delivery System	8.5	202.5
M00533580	08714729787105	RX Delivery System	10	202.5

*Compatible with Advanix Biliary Stents.

Biliary Plastic Stents

DUODENAL BEND STENTS

MR Magnetic Resonance Safe

RX Duodenal Bend Biliary Stent with RX Delivery System



Order Number	GTIN	Description	O.D. (Fr)	Length (cm) Barb to Barb	Recommended Guidewire (in)	Packaging
M00545550	08714729337850	Duodenal Bend Biliary Stent with RX Delivery System	7	5	.035	Each
M00545560	08714729337867	Duodenal Bend Biliary Stent with RX Delivery System	7	7	.035	Each
M00545570	08714729337874	Duodenal Bend Biliary Stent with RX Delivery System	7	10	.035	Each
M00545580	08714729337881	Duodenal Bend Biliary Stent with RX Delivery System	7	12	.035	Each
M00545590	08714729337898	Duodenal Bend Biliary Stent with RX Delivery System	7	15	.035	Each
M00545650	08714729252221	Duodenal Bend Biliary Stent with RX Delivery System	8.5	5	.035	Each
M00545660	08714729301837	Duodenal Bend Biliary Stent with RX Delivery System	8.5	7	.035	Each
M00545670	08714729301844	Duodenal Bend Biliary Stent with RX Delivery System	8.5	10	.035	Each
M00545680	08714729301851	Duodenal Bend Biliary Stent with RX Delivery System	8.5	12	.035	Each
M00545690	08714729301868	Duodenal Bend Biliary Stent with RX Delivery System	8.5	15	.035	Each
M00545600	08714729337904	Duodenal Bend Biliary Stent with RX Delivery System	10	5	.035	Each
M00545610	08714729337911	Duodenal Bend Biliary Stent with RX Delivery System	10	7	.035	Each
M00545620	08714729337928	Duodenal Bend Biliary Stent with RX Delivery System	10	10	.035	Each
M00545630	08714729301813	Duodenal Bend Biliary Stent with RX Delivery System	10	12	.035	Each
M00545640	08714729301820	Duodenal Bend Biliary Stent with RX Delivery System	10	15	.035	Each
M00545700	08714729301875	Duodenal Bend Biliary Stent with RX Delivery System	11.5	5	.035	Each
M00545710	08714729301882	Duodenal Bend Biliary Stent with RX Delivery System	11.5	7	.035	Each
M00545720	08714729301899	Duodenal Bend Biliary Stent with RX Delivery System	11.5	10	.035	Each
M00545730	08714729337942	Duodenal Bend Biliary Stent with RX Delivery System	11.5	12	.035	Each
M00545740	08714729337959	Duodenal Bend Biliary Stent with RX Delivery System	11.5	15	.035	Each



MR Magnetic Resonance Safe


Flexima™ Duodenal Bend Biliary Stent with Delivery System



Order Number	GTIN	O.D. (Fr)	Length (cm) Barb to Barb	Recommended Guidewire (in)	Packaging
M00539200	08714729337362	7	5	.035	Each
M00539210	08714729162520	7	7	.035	Each
M00539220	08714729162537	7	10	.035	Each
M00539230	08714729162544	7	12	.035	Each
M00539240	08714729162551	7	15	.035	Each
M00539250	08714729162568	8.5	5	.035	Each
M00539260	08714729162575	8.5	7	.035	Each
M00539270	08714729162582	8.5	10	.035	Each
M00539280	08714729162599	8.5	12	.035	Each
M00539290	08714729162605	8.5	15	.035	Each
M00539300	08714729162612	10	5	.035	Each
M00539310	08714729162629	10	7	.035	Each
M00539320	08714729162636	10	10	.035	Each
M00539330	08714729162643	10	12	.035	Each
M00539340	08714729162650	10	15	.035	Each
M00539350	08714729162667	11.5	5	.035	Each
M00539360	08714729162674	11.5	7	.035	Each
M00539370	08714729162681	11.5	10	.035	Each
M00539380	08714729162698	11.5	12	.035	Each
M00539390	08714729162704	11.5	15	.035	Each

Biliary Plastic Stents

PIGTAIL STENTS




MR Magnetic Resonance Safe

Flexima™ Single Pigtail Nasobiliary Catheter*				
Order Number	GTIN	O.D. (Fr)	Working Length (cm)	Packaging
M00540110	08714729337379	5	250	Each
M00540120	08714729337386	6	250	Each
M00540130	08714729337393	7.5	250	Each
M00540140	08714729337409	8.5	250	Each

*Recommended Guidewire .035" Dreamwire™ High Performance Guidewire, Order #5614, [▶](#) (See page 1).

Pancreatic Plastic Stents

ADVANIX™ PANCREATIC STENTS



 Magnetic Resonance Conditional¹

Advanix Pancreatic Stent Kit - Straight, No Leading Barb and Leading Barb



Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00537460	08714729856245	Straight Stent – No Leading Barb	4	3	Each
M00537470	08714729856252	Straight Stent – No Leading Barb	4	4	Each
M00537480	08714729856269	Straight Stent – No Leading Barb	4	5	Each
M00537490	08714729862956	Straight Stent – No Leading Barb	4	7	Each
M00537500	08714729856276	Straight Stent – No Leading Barb	5	3	Each
M00537510	08714729856283	Straight Stent – No Leading Barb	5	4	Each
M00537520	08714729856290	Straight Stent – No Leading Barb	5	5	Each
M00537530	08714729856306	Straight Stent – No Leading Barb	5	7	Each
M00537540	08714729856313	Straight Stent – Leading Barb	5	3	Each
M00537550	08714729856320	Straight Stent – Leading Barb	5	4	Each
M00537560	08714729856337	Straight Stent – Leading Barb	5	5	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

Advanix Pancreatic Stent Kit - Pigtail, Leading Barb and No Leading Barb



Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00537350	08714729856344	Single Pigtail Stent – No Leading Barb	4	3	Each
M00537360	08714729856351	Single Pigtail Stent – No Leading Barb	4	4	Each
M00537370	08714729856368	Single Pigtail Stent – No Leading Barb	4	5	Each
M00537380	08714729856375	Single Pigtail Stent – No Leading Barb	4	7	Each
M00537390	08714729856382	Single Pigtail Stent – No Leading Barb	5	3	Each
M00537400	08714729856399	Single Pigtail Stent – No Leading Barb	5	4	Each
M00537410	08714729856405	Single Pigtail Stent – No Leading Barb	5	5	Each
M00537420	08714729856412	Single Pigtail Stent – No Leading Barb	5	7	Each
M00537430	08714729856429	Single Pigtail Stent – Leading Barb	5	3	Each
M00537440	08714729856436	Single Pigtail Stent – Leading Barb	5	4	Each
M00537450	08714729856443	Single Pigtail Stent – Leading Barb	5	5	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Pancreatic Plastic Stents

ADVANIX™ PANCREATIC STENTS (Continued)



 Magnetic Resonance Conditional¹

Advanix Pancreatic Stent - Straight, No Leading Barb					
Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00536000	08714729857686	Straight – No Leading Barb	3	2	Each
M00536010	08714729857693	Straight – No Leading Barb	3	3	Each
M00536020	08714729857709	Straight – No Leading Barb	3	4	Each
M00536030	08714729857716	Straight – No Leading Barb	3	5	Each
M00536040	08714729857723	Straight – No Leading Barb	3	6	Each
M00536050	08714729857730	Straight – No Leading Barb	3	7	Each
M00536060	08714729857747	Straight – No Leading Barb	3	8	Each
M00536070	08714729857754	Straight – No Leading Barb	3	9	Each
M00536080	08714729857761	Straight – No Leading Barb	3	10	Each
M00536090	08714729857778	Straight – No Leading Barb	3	11	Each
M00536100	08714729857785	Straight – No Leading Barb	3	12	Each
M00536110	08714729857792	Straight – No Leading Barb	3	13	Each
M00536120	08714729857068	Straight – No Leading Barb	4	2	Each
M00536130	08714729857075	Straight – No Leading Barb	4	3	Each
M00536140	08714729857082	Straight – No Leading Barb	4	4	Each
M00536150	08714729857099	Straight – No Leading Barb	4	5	Each
M00536160	08714729857105	Straight – No Leading Barb	4	6	Each
M00536170	08714729857112	Straight – No Leading Barb	4	7	Each
M00536180	08714729857129	Straight – No Leading Barb	4	8	Each
M00536190	08714729857136	Straight – No Leading Barb	4	9	Each
M00536200	08714729857143	Straight – No Leading Barb	4	10	Each
M00536210	08714729857150	Straight – No Leading Barb	4	11	Each
M00536220	08714729857167	Straight – No Leading Barb	4	12	Each
M00536230	08714729857174	Straight – No Leading Barb	4	13	Each
M00536240	08714729857181	Straight – No Leading Barb	5	2	Each
M00536250	08714729857198	Straight – No Leading Barb	5	3	Each
M00536260	08714729857204	Straight – No Leading Barb	5	4	Each
M00536270	08714729857211	Straight – No Leading Barb	5	5	Each
M00536280	08714729857228	Straight – No Leading Barb	5	6	Each
M00536290	08714729857235	Straight – No Leading Barb	5	7	Each
M00536300	08714729857242	Straight – No Leading Barb	5	8	Each
M00536310	08714729857259	Straight – No Leading Barb	5	9	Each
M00536320	08714729857266	Straight – No Leading Barb	5	10	Each
M00536330	08714729857273	Straight – No Leading Barb	5	11	Each
M00536340	08714729857280	Straight – No Leading Barb	5	12	Each
M00536350	08714729857297	Straight – No Leading Barb	5	13	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Pancreatic Plastic Stents

ADVANIX™ PANCREATIC STENTS (Continued)



 Magnetic Resonance Conditional¹

Advanix Pancreatic Stent - Straight, Leading Barb 					
Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00536360	08714729857303	Straight – Leading Barb	5	2	Each
M00536370	08714729857310	Straight – Leading Barb	5	3	Each
M00536380	08714729857327	Straight – Leading Barb	5	4	Each
M00536390	08714729857334	Straight – Leading Barb	5	5	Each
M00536400	08714729857341	Straight – Leading Barb	5	6	Each
M00536410	08714729857358	Straight – Leading Barb	5	7	Each
M00536420	08714729857365	Straight – Leading Barb	5	8	Each
M00536430	08714729857372	Straight – Leading Barb	5	9	Each
M00536440	08714729857389	Straight – Leading Barb	5	10	Each
M00536450	08714729857396	Straight – Leading Barb	5	11	Each
M00536460	08714729857402	Straight – Leading Barb	5	12	Each
M00536470	08714729857419	Straight – Leading Barb	5	13	Each
M00536480	08714729857426	Straight – Leading Barb	7	3	Each
M00536490	08714729857433	Straight – Leading Barb	7	4	Each
M00536500	08714729857440	Straight – Leading Barb	7	5	Each
M00536510	08714729857457	Straight – Leading Barb	7	6	Each
M00536520	08714729857464	Straight – Leading Barb	7	7	Each
M00536530	08714729857471	Straight – Leading Barb	7	8	Each
M00536540	08714729857488	Straight – Leading Barb	7	9	Each
M00536550	08714729857495	Straight – Leading Barb	7	10	Each
M00536560	08714729857501	Straight – Leading Barb	7	11	Each
M00536570	08714729857518	Straight – Leading Barb	7	12	Each
M00536580	08714729857525	Straight – Leading Barb	7	13	Each
M00536590	08714729857532	Straight – Leading Barb	7	15	Each
M00536600	08714729857549	Straight – Leading Barb	7	18	Each
M00536610	08714729857556	Straight – Leading Barb	10	3	Each
M00536620	08714729857563	Straight – Leading Barb	10	4	Each
M00536630	08714729857570	Straight – Leading Barb	10	5	Each
M00536640	08714729857587	Straight – Leading Barb	10	6	Each
M00536650	08714729857594	Straight – Leading Barb	10	7	Each
M00536660	08714729857600	Straight – Leading Barb	10	8	Each
M00536670	08714729857617	Straight – Leading Barb	10	9	Each
M00536680	08714729857624	Straight – Leading Barb	10	10	Each
M00536690	08714729857631	Straight – Leading Barb	10	11	Each
M00536700	08714729857648	Straight – Leading Barb	10	12	Each
M00536710	08714729857655	Straight – Leading Barb	10	13	Each
M00536720	08714729857662	Straight – Leading Barb	10	15	Each
M00536730	08714729857679	Straight – Leading Barb	10	18	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents



Prev Page

24

Next Page



Home



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pancreatic Plastic Stents

ADVANIX™ PANCREATIC STENTS (Continued)



 Magnetic Resonance Conditional¹

Advanix Pancreatic Stent - Pigtail, No Leading Barb 					
Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00536740	08714729856948	Single Pigtail – No Leading Barb	3	2	Each
M00536750	08714729856955	Single Pigtail – No Leading Barb	3	3	Each
M00536760	08714729856962	Single Pigtail – No Leading Barb	3	4	Each
M00536770	08714729856979	Single Pigtail – No Leading Barb	3	5	Each
M00536780	08714729856986	Single Pigtail – No Leading Barb	3	6	Each
M00536790	08714729856993	Single Pigtail – No Leading Barb	3	7	Each
M00536800	08714729857006	Single Pigtail – No Leading Barb	3	8	Each
M00536810	08714729857013	Single Pigtail – No Leading Barb	3	9	Each
M00536820	08714729857020	Single Pigtail – No Leading Barb	3	10	Each
M00536830	08714729857037	Single Pigtail – No Leading Barb	3	11	Each
M00536840	08714729857044	Single Pigtail – No Leading Barb	3	12	Each
M00536850	08714729857051	Single Pigtail – No Leading Barb	3	13	Each
M00536860	08714729856450	Single Pigtail – No Leading Barb	4	2	Each
M00536870	08714729856467	Single Pigtail – No Leading Barb	4	3	Each
M00536880	08714729856474	Single Pigtail – No Leading Barb	4	4	Each
M00536890	08714729856481	Single Pigtail – No Leading Barb	4	5	Each
M00536900	08714729856498	Single Pigtail – No Leading Barb	4	6	Each
M00536910	08714729856504	Single Pigtail – No Leading Barb	4	7	Each
M00536920	08714729856511	Single Pigtail – No Leading Barb	4	8	Each
M00536930	08714729856528	Single Pigtail – No Leading Barb	4	9	Each
M00536940	08714729856535	Single Pigtail – No Leading Barb	4	10	Each
M00536950	08714729856542	Single Pigtail – No Leading Barb	4	11	Each
M00536960	08714729856559	Single Pigtail – No Leading Barb	4	12	Each
M00536970	08714729856566	Single Pigtail – No Leading Barb	4	13	Each
M00536980	08714729856696	Single Pigtail – No Leading Barb	5	2	Each
M00536990	08714729856702	Single Pigtail – No Leading Barb	5	3	Each
M00537000	08714729856719	Single Pigtail – No Leading Barb	5	4	Each
M00537010	08714729856726	Single Pigtail – No Leading Barb	5	5	Each
M00537020	08714729856733	Single Pigtail – No Leading Barb	5	6	Each
M00537030	08714729856740	Single Pigtail – No Leading Barb	5	7	Each
M00537040	08714729856757	Single Pigtail – No Leading Barb	5	8	Each
M00537050	08714729856764	Single Pigtail – No Leading Barb	5	9	Each
M00537060	08714729856771	Single Pigtail – No Leading Barb	5	10	Each
M00537070	08714729856788	Single Pigtail – No Leading Barb	5	11	Each
M00537080	08714729856795	Single Pigtail – No Leading Barb	5	12	Each
M00537090	08714729856801	Single Pigtail – No Leading Barb	5	13	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Pancreatic Plastic Stents

ADVANIX™ PANCREATIC STENTS (Continued)



 Magnetic Resonance
Conditional¹

Advanix Pancreatic Stent - Pigtail, Leading Barb



Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00537100	08714729856573	Single Pigtail – Leading Barb	5	2	Each
M00537110	08714729856580	Single Pigtail – Leading Barb	5	3	Each
M00537120	08714729856597	Single Pigtail – Leading Barb	5	4	Each
M00537130	08714729856603	Single Pigtail – Leading Barb	5	5	Each
M00537140	08714729856610	Single Pigtail – Leading Barb	5	6	Each
M00537150	08714729856627	Single Pigtail – Leading Barb	5	7	Each
M00537160	08714729856634	Single Pigtail – Leading Barb	5	8	Each
M00537170	08714729856641	Single Pigtail – Leading Barb	5	9	Each
M00537180	08714729856658	Single Pigtail – Leading Barb	5	10	Each
M00537190	08714729856665	Single Pigtail – Leading Barb	5	11	Each
M00537200	08714729856672	Single Pigtail – Leading Barb	5	12	Each
M00537210	08714729856689	Single Pigtail – Leading Barb	5	13	Each
M00537220	08714729856818	Single Pigtail – Leading Barb	7	3	Each
M00537230	08714729856825	Single Pigtail – Leading Barb	7	4	Each
M00537240	08714729856832	Single Pigtail – Leading Barb	7	5	Each
M00537250	08714729856849	Single Pigtail – Leading Barb	7	6	Each
M00537260	08714729856856	Single Pigtail – Leading Barb	7	7	Each
M00537270	08714729856863	Single Pigtail – Leading Barb	7	8	Each
M00537280	08714729856870	Single Pigtail – Leading Barb	7	9	Each
M00537290	08714729856887	Single Pigtail – Leading Barb	7	10	Each
M00537300	08714729856894	Single Pigtail – Leading Barb	7	11	Each
M00537310	08714729856900	Single Pigtail – Leading Barb	7	12	Each
M00537320	08714729856917	Single Pigtail – Leading Barb	7	13	Each
M00537330	08714729856924	Single Pigtail – Leading Barb	7	15	Each
M00537340	08714729856931	Single Pigtail – Leading Barb	7	18	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

► Pancreatic Plastic Stents


NAVIFLEX™ RX PANCREATIC DELIVERY SYSTEM & NAVIFLEX PUSHERS



NaviFlex RX Pancreatic Delivery System & NaviFlex Pushers					
Order Number	GTIN	Description	Diameter (Fr)	Length	Packaging
M00535060	08714729856221	NaviFlex RX Pancreatic Stent Delivery System	7	N/A	Each
M00535080	08714729856238	NaviFlex RX Pancreatic Stent Delivery System	10	N/A	Each
M00535000	08714729856207	Pancreatic Stent RX Pusher	3	N/A	Each
M00535010	08714729856214	Pancreatic Stent RX Pusher	4 / 5	N/A	Each
M00535020	08714729856184	Pancreatic Stent Pusher	3	N/A	Each
M00535030	08714729856191	Pancreatic Stent Pusher	4 / 5	N/A	Each

Biliary Stone Management


RX STONE MANAGEMENT



Extractor™ Pro RX-S Retrieval Balloon Catheter – Distal Injection




Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547300	08714729907138	Extractor Pro RX-S Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547310	08714729907152	Extractor Pro RX-S Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547320	08714729907176	Extractor Pro RX-S Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Extractor Pro RX-S Retrieval Balloon Catheter – Proximal Injection




Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547330	08714729907145	Extractor Pro RX-S Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547340	08714729907169	Extractor Pro RX-S Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547350	08714729907183	Extractor Pro RX-S Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Extractor Pro RX Retrieval Balloon Catheter – Distal Injection



Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547000	08714729790266	Extractor Pro RX Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547010	08714729790297	Extractor Pro RX Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547020	08714729790310	Extractor Pro RX Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Extractor Pro RX Retrieval Balloon Catheter – Proximal Injection




Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547030	08714729790280	Extractor Pro RX Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547040	08714729790273	Extractor Pro RX Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547050	08714729790303	Extractor Pro RX Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Biliary Stone Management


STONE MANAGEMENT



Extractor™ Pro XL Retrieval Balloon Catheter – Distal Injection




Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547100	08714729790389	Extractor Pro XL Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547110	08714729790402	Extractor Pro XL Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547120	08714729790648	Extractor Pro XL Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Extractor Pro XL Retrieval Balloon Catheter – Proximal Injection



Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547130	08714729790396	Extractor Pro XL Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547140	08714729790365	Extractor Pro XL Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547150	08714729790372	Extractor Pro XL Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Trapezoid™ RX Wireguided Retrieval Baskets




Order Number	GTIN	Description	Open Basket Diameter (cm)	Recommended Guidewire (in)	Required Working Channel (mm)	Packaging
M00510860	08714729296379	Trapezoid RX Wireguided Retrieval Basket	1.5	.035	3.2	Each
M00510870	08714729296386	Trapezoid RX Wireguided Retrieval Basket	2.0	.035	3.2	Each
M00510880	08714729296393	Trapezoid RX Wireguided Retrieval Basket	2.5	.035	3.2	Each
M00510890*	08714729296409	Trapezoid RX Wireguided Retrieval Basket	3.0	.035	3.2	Each

Recommend Alliance™ II Inflation System for use as Mechanical Lithotripter; Order #5062, [\(See page 84\)](#).

*Refer to Instructions for Use when using this device; recommended for use with Biliary Calculi of >1.5cm.

Biliary Stone Management

CRE™ WIREGUIDED BALLOON DILATATION CATHETERS




CRE PRO Wireguided Balloon Dilatation Catheters

Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558600	08714729797487	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8 / 3.2*	7.5	180	Each
M00558610	08714729797494	8-9-10	24-27-30	5.5	3-5-5-9	304-557-912	2.8 / 3.2*	7.5	180	Each
M00558620	08714729797500	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8 / 3.2*	7.5	180	Each
M00558630	08714729797517	12-13.5-15	36-40.5-45	5.5	3-4-5-8	304-456-811	2.8 / 3.2*	7.5	180	Each
M00558640	08714729797524	15-16.5-18	45-49.5-54	5.5	3-4-5-7	304-456-709	2.8 / 3.7*	7.5	180	Each
M00558650	08714729797531	18-19-20	54-57-60	5.5	3-4-5-6	304-456-608	2.8 / 3.7*	7.5	180	Each
M00558660	08714729797548	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8 / 3.2*	7.5	240	Each
M00558670	08714729797555	8-9-10	24-27-30	5.5	3-5-5-9	304-557-912	2.8 / 3.2*	7.5	240	Each
M00558680	08714729797562	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8 / 3.2*	7.5	240	Each
M00558690	08714729797579	12-13.5-15	36-40.5-45	5.5	3-4-5-8	304-456-811	2.8 / 3.2*	7.5	240	Each
M00558700	08714729797586	15-16.5-18	45-49.5-54	5.5	3-4-5-7	304-456-709	2.8 / 3.7*	7.5	240	Each
M00558710	08714729797593	18-19-20	54-57-60	5.5	3-4-5-6	304-456-608	2.8 / 3.7*	7.5	240	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060,

(See page 84). Guidewire included and preloaded in balloon dilatation catheter.

*Duodenoscope channel size.



CRE RX Biliary Balloon Dilatation Catheters*


Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558860	08714729965411	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	3.7	7.5	180	Each
M00558870	08714729965428	8-9-10	24-27-30	5.5	3-5-5-9	304-557-912	3.7	7.5	180	Each
M00558880	08714729965435	10-11-12	30-33-36	5.5	3-5-8	304-507-811	3.7	7.5	180	Each
M00558890	08714729965442	12-13.5-15	36-40.5-45	5.5	3-4-5-8	304-456-811	3.7	7.5	180	Each
M00558900	08714729965459	15-16.5-18	45-49.5-54	5.5	3-4-5-7	304-456-709	3.7	7.5	180	Each
M00558910	08714729965466	18-19-20	54-57-60	5.5	3-4-5-6	304-456-608	3.7	7.5	180	Each
M00558920	08714729965473	8-9-10	24-27-30	3.0	3-5-5-9	304-608-1013	3.7	7.5	180	Each
M00558930	08714729965480	10-11-12	30-33-36	3.0	3-5-8	304-557-912	3.7	7.5	180	Each
M00558940	08714729965497	12-13.5-15	36-40.5-45	3.0	3-4-5-8	304-507-811	3.7	7.5	180	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060.

*No Pre-loaded Guidewire in the Balloon Dilatation Catheter.

Biliary Stricture Management


PANCREATICOBILIARY ABLATION



Habib™ EndoHPB Bipolar Radiofrequency Catheter

Order Number	GTIN	Description	Diameter (Fr)	Working Length (cm)	Packaging
M00500070	08714729981909	 Habib EndoHPB Bipolar RF Catheter	8	180	Each
5100	08714729982524	Bipolar Adapter Cable, 28mm spacing	NA	200	Each
5420	08714729982548	Bipolar Adapter Cable, Internal 4mm pin	NA	200	Each
5700	08714729982562	Bipolar Adapter Cable, 22mm spacing	NA	200	Each

RX DILATATION



Hurricane™ RX Biliary Balloon Dilatation Catheter




Order Number	GTIN	Balloon Inflated O.D. (mm)	Balloon Length (cm)	Catheter O.D. (Fr)	Catheter Length (cm)	Maximum Inflation Pressure (atm)	Minimum Channel Size (mm)	Recommended Guidewire (in)	Packaging
M00545890	08714729283805	4	2	5.8	180	12	3.2	.035	Each
M00545900	08714729283812	4	4	5.8	180	12	3.2	.035	Each
M00545910	08714729283829	6	2	5.8	180	11	3.2	.035	Each
M00545920	08714729283836	6	4	5.8	180	11	3.2	.035	Each
M00545930	08714729283843	8	2	5.8	180	11	3.2	.035	Each
M00545940	08714729283850	8	4	5.8	180	11	3.2	.035	Each
M00545950	08714729283867	10	2	5.8	180	8	3.2	.035	Each
M00545960	08714729283874	10	4	5.8	180	8	3.2	.035	Each

Recommended inflation device: Encore™ 26 Inflation Device Order #6667,  (See page 32).

► Biliary Stricture Management

RX CYTOLOGY




RX Cytology Brush Wireguided Cytology Brush



Order Number	GTIN	Description	Brush Width (mm)	Catheter O.D. (Fr)	Working Length (cm)	Recommended Guidewire (in)	Packaging
M00545000	08714729268628	RX Cytology Brush	2.1	8	200	.035	Each

BILIARY DILATATION




MaxForce™ Biliary Balloon Dilatation Catheter



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Catheter Size		Working Length (cm)	Required Working Channel (mm)	Recommended Guidewire (in)	Packaging
		(mm)	(Fr)		(mm)	(Fr)				
M00567330	08714729097037	4	12	2	1.67	5	180	2.8	.035	Each
M00567340	08714729097051	4	12	4	1.67	5	180	2.8	.035	Each
M00567370	08714729097068	6	18	2	1.67	5	180	2.8	.035	Each
M00567380	08714729097075	6	18	4	1.67	5	180	2.8	.035	Each
M00567410	08714729097082	8	24	2	1.67	5	180	2.8	.035	Each
M00567420	08714729097099	8	24	3	1.67	5	180	2.8	.035	Each

RX ACCESSORIES



Encore™ 26 Inflation Device

Order Number	GTIN	Description	Packaging
M00566670	08714729755241	Encore 26 Inflation Device	Each

Notes and Future Product Updates

Notes and Future Product Updates



Section Table of Contents



Prev Page

Next Page



Home



► Cholangioscopy



SpyGlass™ Discover Digital Catheter

► SpyGlass Discover Digital Catheter33	► Accessories33
--	-----------------------

SpyGlass™ DS Direct Visualization System

► Capital Equipment34	► Disposable Devices34
-----------------------------	------------------------------

EXALT™ Model D Single-Use Duodenoscope


► Capital Equipment35	► Disposable Devices35
-----------------------------	------------------------------

Autolith™ Touch Biliary EHL System

► Autolith Touch System35	► Biliary EHL Probe35
► Autolith Touch System Accessories.....35	

SpyGlass™ Discover Digital Catheter & Accessories

SPYGLASS DISCOVER DIGITAL CATHETER & ACCESSORIES



Includes
(1) Y-Adaptor



Type BF Applied Part

SpyGlass Discover Digital Catheter & Accessories



Order Number	GTIN	Description	Packaging
M00546780	8714729994183	SpyGlass™ Discover Digital Catheter	Box 1
M00542950	8714729998723	SpyGlass™ Discover Controller	Box 1
M00543000	00191506013516	SpyGlass™ Discover Balloon Dilation Catheter (6.0x40mmx40cm)	Box 1
M00543010	00191506013509	SpyGlass™ Discover Balloon Dilation Catheter (7.0x40mmx40cm)	Box 1
M00543040	00191506013493	SpyGlass™ Discover Balloon Dilation Catheter (8.0x40mmx40cm)	Box 1
M00566670	08714729755241	Encore™ 26 Inflation Device	Each
M00542961	14543527251959	Super Sheath CBDE for use with the SpyGlass	Box 10
M00542990	191506004255	SpyGlass™ Discover Retrieval Basket Box	Box 1
M00542981	191506005238	SpyGlass™ Discover Imager™ II IOC Catheter	Box 5
M00542971	191506004774	SpyGlass™ Discover Jagwire™ Guidewire	Box 2



Section Table of Contents



Prev Page

33

Next Page




Home



SpyGlass™ DS Direct Visualization System

CAPITAL EQUIPMENT & DEVICES



Capital Equipment

Order Number	GTIN	Description	Packaging
M00546650	08714729874348	SpyGlass DS Digital Controller	Please contact your Boston Scientific Sales Representative for ordering information.

For more information on ordering a compatible EHL or Laser Probe for the SpyGlass DS Direct Visualization System, please contact your Boston Scientific Representative.


Disposable Devices



Order Number	GTIN	Description	Packaging
M00546600	08714729863236	SpyScope™ DS Access & Delivery Catheter	Each
M00546610	08714729965404	SpyScope DS II Access & Delivery Catheter	Each
M00546470	08714729996699	SpyBite™ Max Biopsy Forceps	Each
M00546550	08714729965886	SpyGlass Retrieval Basket	Each
M00546560	08714729965893	SpyGlass Retrieval Snare	Each

EXALT™ Model D Singe-Use Duodenoscope

CAPITAL EQUIPMENT & DEVICES



Capital Equipment

Order Number	GTIN	Description	Packaging
M00542430	08714729990901	EXALT™ Controller	Please contact your Boston Scientific Sales Representative for ordering information.

For more information on ordering a compatible EHL or Laser Probe for the SpyGlass DS Direct Visualization System, please contact your Boston Scientific Representative.

Disposable Devices



Order Number	GTIN	Description	Packaging
M00542421	08714729993605	EXALT™ Model D Single-use Duodenoscope	Box 2
M00542420	08714729983514	EXALT™ Model D Single-use Duodenoscope	Each

Autolith™ Touch Biliary EHL System

AUTOLITH TOUCH BILIARY EHL SYSTEM

Biliary EHL Probe

Order Number	Description	Packaging
M00546620	 1.9, 375cm Biliary Electrohydraulic Lithotripter Probe	Each

Autolith Touch System

M005466800	Autolith Touch Bipolar Electrohydraulic Lithotripter	Each
------------	--	------

Autolith Touch System Accessories

M00546750	Autolith Touch Extender Cable	Each
M00546760	Autolith Touch Foot Switch	Each

Notes and Future Product Updates



Section Table of Contents




Prev Page

Next Page



Home





EUS

Expect™ Slimline (SL) Endoscopic Ultrasound Aspiration Needles

- Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle36

Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device

- Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device36



LumiCoil™ Platinum Fiducial Markers

- LumiCoil Platinum Fiducial Markers – Straight Configuration..... 37
- LumiCoil Platinum Fiducial Markers – Figure 8 Configuration 37

AXIOS™ System

- AXIOS Stent and Electrocautery Enhanced Delivery System..... 37
- AXIOS Stent and Delivery System 37

EXPECT™ SLIMLINE (SL) ENDOSCOPIC ULTRASOUND ASPIRATION NEEDLES



Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle*




Order Number	GTIN	Needle Size	Minimum Working Channel (mm)	Sheath Diameter (mm)	Packaging
M00555500	08714729848141	19 gauge (1.10mm)	2.8	1.83	Each
M00555510	08714729848158	22 gauge (0.72mm)	2.4	1.65	Each
M00555520	08714729848165	25 gauge (0.52mm)	2.4	1.52	Each
M00555530	08714729848172	19 gauge Flex (1.14mm)	2.8	1.73	Each
M00555501	08714729848714	19 gauge (1.10mm)	2.8	1.83	Box 5
M00555511	08714729848721	22 gauge (0.72mm)	2.4	1.65	Box 5
M00555521	08714729848738	25 gauge (0.52mm)	2.4	1.52	Box 5
M00555531	08714729848806	19 gauge Flex (1.14mm)	2.8	1.73	Box 5

*Packaging also includes: 20cc syringe and one-way stopcock

*Needle length: 0cm to 8cm, adjustable

*Working length: 137.5cm to 141.5cm, adjustable

ACQUIRE™ ENDOSCOPIC ULTRASOUND FINE NEEDLE BIOPSY (FNB) DEVICE



Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device*




Order Number	GTIN	Needle Size	Minimum Working Channel (mm)	Sheath Diameter (mm)	Packaging
M00555540	08714729931782	22 gauge (0.72mm)	2.4	1.65	Each
M00555560	08714729931805	25 gauge (0.52mm)	2.4	1.52	Each
M00555580	08714729931829	19 gauge Flex (1.14mm)	2.8	1.73	Each
M00555541	08714729931843	22 gauge (0.72mm)	2.4	1.65	Box 5
M00555561	08714729931867	25 gauge (0.52mm)	2.4	1.52	Box 5

*Packaging also includes: 20cc syringe and one-way stopcock

*Needle length: 0cm to 8cm, adjustable


*Working length: 137.5cm to 141.5cm, adjustable

LUMICOIL™ PLATINUM FIDUCIAL MARKERS



LumiCoil Platinum Fiducial Markers							
Order Number	GTIN	Description	Fiducial O.D. (in)	Overall Fiducial Length (mm)	Implanted Fiducial Length (mm)	Differentiator	Packaging
M00553831	08714729972877	Lumicoil Platinum Fiducial Marker	0.018	5	5	Straight	Box 3
M00553841	08714729972891	Lumicoil Platinum Fiducial Marker	0.018	10	5	Figure 8	Box 3


AXIOS™ SYSTEM



AXIOS Stent and Electrocautery Enhanced Delivery System								
Order Number	GTIN	Flange Diameter (mm)	Lumen Diameter (mm)	Saddle Length (mm)	Catheter O.D. (Fr)	Catheter Working Length (cm)	Catheter Total Length (cm)	Packaging
M00553640	08714729904588	21	10	10	10.8	138	144	Each
M00553650	08714729904595	24	15	10	10.8	138	144	Each
M00553660	08714729951179	29	20	10	10.8	138	144	Each
M00553670	8714729979333	24	15	15	10.8	138	144	Each
M00553680	00191506008086	14	6	8	9	138	144	Each
M00553690	00191506008093	17	8	8	9	138	144	Each

AXIOS Stent and Delivery System								
Order Number	GTIN	Flange Diameter (mm)	Lumen Diameter (mm)	Saddle Length (mm)	Catheter O.D. (Fr)	Catheter Working Length (cm)	Catheter Total Length (cm)	Packaging
M00553740	08714729904625	21	10	10	10.8	138	144	Each
M00553750	08714729904632	24	15	10	10.8	138	144	Each

► Enteral and Esophageal Stents



Enteral and Esophageal Stents

Colonic and Duodenal Stents

► WallFlex™ Colonic Stent System	38
► WallFlex Colonic Stent System with Anchor Lock Delivery System	38
► WallFlex Duodenal Stent System	39
► WallFlex Duodenal Stent System with Anchor Lock Delivery System	39
► WALLSTENT™ Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP™ Plus Delivery System	40
► WALLSTENT Super Stiff Guidewire	40

Esophageal Stents

► Agile Fully Covered Esophageal Stent System	41
► Agile Partially Covered Esophageal Stent System	41
► Amplatzer Super Stiff™ Guidewire	44
► Ultraflex Esophageal NG Covered Stent System – Distal Release	43
► Ultraflex Esophageal NG Covered Stent System – Proximal Release	43
► Ultraflex Esophageal NG Covered Stent System – Distal Release	43
► Ultraflex Esophageal NG Covered Stent System – Proximal Release	44
► Ultraflex Esophageal NG Stent Uncovered System – Distal Release	44
► Ultraflex Esophageal NG Uncovered Stent System – Proximal Release	44
► WallFlex Fully Covered Esophageal Stent System	42
► WallFlex Partially Covered Esophageal Stent System	42

WALLFLEX™ COLONIC AND DUODENAL STENTS

Wallflex Soft Colonic Stent System								
Order Number	GTIN	Stent Body Diameter (mm)	Stent Flare Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery System Diameter (Fr)	Packaging
M00555760	08714729902584	22	27	60	230		9	Each
M00555770	08714729902591	22	27	90	230		9	Each
M00555780	08714729902607	22	27	120	230		9	Each



 Magnetic Resonance Conditional¹

WallFlex Colonic Stent System with Anchor Lock Delivery System*								
Order Number	GTIN	Stent Body Diameter (mm)	Stent Flare Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery System Diameter (Fr)	Packaging
M00565040	08714729456513	25	30	60	230	270	10	Each
M00565050	08714729456520	25	30	90	230	270	10	Each
M00565060	08714729456537	25	30	120	230	270	10	Each
M00565070	08714729456544	25	30	60	135	175	10	Each
M00565080	08714729456551	25	30	90	135	175	10	Each
M00565090	08714729456568	25	30	120	135	175	10	Each
M00565100	08714729456575	22	27	60	230	270	10	Each
M00565110	08714729456582	22	27	90	230	270	10	Each
M00565120	08714729456599	22	27	120	230	270	10	Each
M00565130	08714729456605	22	27	60	135	175	10	Each
M00565140	08714729456612	22	27	90	135	175	10	Each
M00565150	08714729456629	22	27	120	135	175	10	Each

*Recommended Guidewire .035" WALLSTENT™ Super Stiff Guidewire, Order #18001,  (See page 40).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

► Enteral Stents

WALLFLEX™ COLONIC AND DUODENAL STENTS (Continued)

WallFlex Soft Duodenal Stent System									
Order Number	GTIN	Stent Body Diameter (mm)	Flare Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery System Diameter (Fr)	Packaging	
M00555860	08714729902492	22	27	60	230		9	Each	
M00555870	08714729902508	22	27	90	230		9	Each	
M00555880	08714729902515	22	27	120	230		9	Each	



 Magnetic Resonance Conditional¹



WallFlex Duodenal Stent System with Anchor Lock Delivery System*									
Order Number	GTIN	Stent Body Diameter (mm)	Flare Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery System Diameter (Fr)	Packaging	
M00565010	08714729456483	22	27	60	230	270	10	Each	
M00565020	08714729456490	22	27	90	230	270	10	Each	
M00565030	08714729456506	22	27	120	230	270	10	Each	

*Recommended Guidewire .035" Hydra Jagwire™ Guidewire Stiff Shaft 450cm Order #5606,  (See page 1).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

► Enteral Stents

WALLSTENT™ COLONIC AND DUODENAL STENTS



 Magnetic Resonance Conditional¹

WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP™ Plus Delivery System*



Order Number	GTIN	Expanded Stent Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery System Diameter (Fr)	Packaging
M00565560	08714729285052	20	60	230	255	10	Each
M00565570	08714729285069	20	90	230	255	10	Each
M00565580	08714729285076	22	60	230	255	10	Each
M00565590	08714729285083	22	90	230	255	10	Each
M00565620	08714729285113	20	60	135	160	10	Each
M00565630	08714729285120	20	90	135	160	10	Each
M00565640	08714729285137	22	60	135	160	10	Each
M00565650	08714729285144	22	90	135	160	10	Each

*Recommended Guidewire .035" WALLSTENT Super Stiff Guidewire, Order #18001, (see below).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

WALLSTENT Super Stiff Guidewire



Order Number	GTIN	Description	Distal Tip O.D. (in)	Length (cm)	Tip Style	Packaging
H965180011	08714729203070	WALLSTENT Super Stiff Guidewire	.035	500	Tapered	Box 5

Esophageal Stents


AGILE™ ESOPHAGEAL STENTS

Agile Partially Covered Esophageal Stent System						
Order Number	GTIN	Stent O.D. (mm)	Proximal / Distal Flares O.D. (mm)	Stent Length (cm)	Catheter Diameter (Fr / mm)	Packaging
M00517710	191506006587	23	28	6.2	18.5 / 6.2	Each
M00517720	191506006594	23	28	10.1	18.5 / 6.2	Each
M00517730	191506006600	23	28	12	18.5 / 6.2	Each
M00517740	191506006617	23	28	15	18.5 / 6.2	Each
M00517200	08714729972907	14	19	6.2	10.5 / 3.5	Each
M00517210	08714729972914	14	19	10.2	10.5 / 3.5	Each
M00517220	08714729972921	14	19	11.9	10.5 / 3.5	Each
M00517230	08714729972938	14	19	14.8	10.5 / 3.5	Each
M00517240	08714729972945	18	23	5.9	10.5 / 3.5	Each
M00517250	08714729972952	18	23	9.7	10.5 / 3.5	Each
M00517260	08714729972969	18	23	11.9	10.5 / 3.5	Each
M00517270	08714729972976	18	23	14.9	10.5 / 3.5	Each

Agile Fully Covered Esophageal Stent System						
Order Number	GTIN	Stent O.D. (mm)	Proximal / Distal Flares O.D. (mm)	Stent Length (cm)	Catheter Diameter (Fr / mm)	Packaging
M00517750	191506006624	23	28	6.2	18.5 / 6.2	Each
M00517760	191506006631	23	28	10.1	18.5 / 6.2	Each
M00517770	191506006648	23	28	12	18.5 / 6.2	Each
M00517780	191506006655	23	28	15	18.5 / 6.2	Each
M00517400	08714729973027	14	19	6.2	10.5 / 3.5	Each
M00517410	08714729973034	14	19	10.2	10.5 / 3.5	Each
M00517420	08714729973041	14	19	11.9	10.5 / 3.5	Each
M00517430	08714729973058	14	19	14.8	10.5 / 3.5	Each
M00517440	08714729973065	18	23	5.9	10.5 / 3.5	Each
M00517450	08714729973072	18	23	9.7	10.5 / 3.5	Each
M00517460	08714729973089	18	23	11.9	10.5 / 3.5	Each
M00517470	08714729973096	18	23	14.9	10.5 / 3.5	Each

Esophageal Stents

WALLFLEX™ ESOPHAGEAL STENTS



 Magnetic Resonance Conditional¹


WallFlex Fully Covered Esophageal Stent System*



Order Number	GTIN	Stent O.D. (mm)	Proximal / Distal Flares O.D. (mm)	Stent Length (cm)	Catheter Diameter (Fr / mm)	Working Length (cm)	Packaging
M00516700	08714729778035	18	25 / 23	10.3	18.5 / 6.2	78	Each
M00516710	08714729778042	18	25 / 23	12.3	18.5 / 6.2	78	Each
M00516720	08714729778059	18	25 / 23	15.3	18.5 / 6.2	78	Each
M00516730	08714729778066	23	28 / 28	10.5	18.5 / 6.2	78	Each
M00516740	08714729778073	23	28 / 28	12.5	18.5 / 6.2	78	Each
M00516750	08714729778080	23	28 / 28	15.5	18.5 / 6.2	78	Each

*Recommended Guidewire .038" Jagwire™ Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

WallFlex Partially Covered Esophageal Stent System*




Order Number	GTIN	Stent O.D. (mm)	Proximal / Distal Flares O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Catheter Diameter (Fr / mm)	Working Length (cm)	Packaging
M00516900	08714729765240	18	23 / 23	10.3	7.1	18.5 / 6.2	78	Each
M00516910	08714729765257	18	23 / 23	12.3	9.1	18.5 / 6.2	78	Each
M00516920	08714729765264	18	23 / 23	15.3	12.1	18.5 / 6.2	78	Each
M00516930	08714729765271	23	28 / 28	10.5	7.3	18.5 / 6.2	78	Each
M00516940	08714729765288	23	28 / 28	12.5	9.3	18.5 / 6.2	78	Each
M00516950	08714729765295	23	28 / 28	15.5	12.3	18.5 / 6.2	78	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Esophageal Stents

ULTRAFLEX™ STENTS



 Magnetic Resonance Conditional¹

Ultraflex Esophageal NG Covered Stent System – Distal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Packaging
M00513730	08714729716150	18	23	10	7	Each
M00513740	08714729716167	18	23	12	9	Each
M00513750	08714729649113	18	23	15	12	Each

*Recommended Guidewire .038" Jagwire™ Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

 Magnetic Resonance Conditional¹


Ultraflex Esophageal NG Stent Covered System – Proximal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Packaging
M00513840	08714729716099	18	23	10	7	Each
M00513850	08714729716075	18	23	12	9	Each
M00513860	08714729716082	18	23	15	12	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

Ultraflex Esophageal NG Stent Covered System – Distal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Packaging
M00514200	08714729121190	23	28	10	7	Each
M00514210	08714729121206	23	28	12	9	Each

*Recommended Guidewire .038" Jagwire™ Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Esophageal Stents


ULTRAFLEX™ STENTS (Continued)

 Magnetic Resonance Conditional¹

Ultraflex Esophageal NG Covered Stent System – Proximal Release*						
Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Packaging
M00514240	08714729716105	23	28	10	7	Each
M00514250	08714729716112	23	28	12	9	Each

*Recommended Guidewire .038" Jagwire™ Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

Ultraflex Esophageal NG Uncovered Stent System – Distal Release*					
Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Packaging
M00513700	08714729716129	18	23	7	Each
M00513710	08714729716136	18	23	10	Each
M00513720	08714729716143	18	23	15	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

 Magnetic Resonance Conditional¹

Ultraflex Esophageal NG Uncovered Stent System – Proximal Release*					
Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Packaging
M00513800	08714729716174	18	23	7	Each
M00513810	08714729716181	18	23	10	Each
M00513830	08714729716204	18	23	15	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

 Magnetic Resonance Conditional¹

Amplatz Super Stiff™ Guidewire						
Order Number	GTIN	Description	O.D. (in)	Length (cm)	Tip Style	Packaging
M00550090	08714729301929	Amplatz Super Stiff Guidewire	.038	260	Straight	Each



Hemostasis

Clipping Devices

- Resolution 360™ ULTRA Clip 45
- Resolution 360™ Clip 45
- Resolution™ Clip 45

Bipolar Electrohemostasis Catheters

- Gold Probe™ Single-Use Electrohemostasis Catheters 46
- Injection Gold Probe™ Single-Use
Electrohemostasis / Injection Catheters 46

Needles

- Interject™ Clear Single-Use Injection Therapy
Needle Catheters 47
- Interject Contrast Single-Use Injection Therapy
Needle Catheters 47


Multiple Band Ligators

- Speedband Superview Super 7™ Multiple Band Ligators 48


Electrosurgical Bipolar Cable Adapter

- Electrosurgical Bipolar Cable Adapter 49

CLIPPING DEVICES




Resolution 360™ ULTRA Clip					
Order Number	GTIN	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Packaging
M00521400	08714729997245	235	2.8	17	Each
M00521401	08714729997252	235	2.8	17	Box 10
M00521402	08714729997269	235	2.8	17	Box 20
M00521403	08714729997276	235	2.8	17	Box 40



 Magnetic Resonance Conditional¹

Resolution 360™ Clip					
Order Number	GTIN	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Packaging
M00521240	08714729875666	155	2.8	11	Each
M00521241	08714729875673	155	2.8	11	Box 10
M00521242	08714729875680	155	2.8	11	Box 20
M00521230	08714729875628	235	2.8	11	Each
M00521231	08714729875635	235	2.8	11	Box 10
M00521232	08714729875642	235	2.8	11	Box 20
M00521233	08714729875659	235	2.8	11	Box 40

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.




 Magnetic Resonance Conditional¹

Resolution™ Clip					
Order Number	GTIN	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Packaging
M00522600	08714729504764	155	2.8	11	Each
M00522601	08714729504771	155	2.8	11	Box 10
M00522602	08714729504788	155	2.8	11	Box 20
M00522610	08714729504795	235	2.8	11	Each
M00522611	08714729504801	235	2.8	11	Box 10
M00522612	08714729504818	235	2.8	11	Box 20

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

BIPOLAR ELECTROHEMOSTASIS CATHETERS




Injection Gold Probe™ Single-Use Electrohemostasis / Injection Catheters



Order Number	GTIN	Description	O.D. (Fr)	Needle Size Gauge	Working Length (cm)	Type Generator / Connector	Packaging
M00560150	08714729056492	Injection Gold Probe Catheter	7	25	210	Standard Plug*	Each
M00560160	08714729056508	Injection Gold Probe Catheter	10	25	210	Standard Plug*	Each

*Compatible with Endostat™ Bipolar-Monopolar Electrosurgical Generators.

For connections to other generators, please refer to the Electrosurgical Bipolar Cable Adapters, (See page 49).



Gold Probe™ Single-Use Electrohemostasis Catheters




Order Number	GTIN	Description	O.D. (Fr)	Working Length (cm)	Type Generator / Connector	Packaging
M00560070	08714729174691	Gold Probe Catheter	7	300	Standard Plug*	Each
M00560100	08714729174684	Gold Probe Catheter	10	300	Standard Plug*	Each
M00560220	08714729099574	Gold Probe Catheter	7	350	Standard Plug*	Each
M00560071	08714729071013	Gold Probe Catheter	7	300	Standard Plug*	Box 5
M00560101	08714729071020	Gold Probe Catheter	10	300	Standard Plug*	Box 5

*Compatible with Endostat Bipolar-Monopolar Electrosurgical Generators.

For connections to other generators, please refer to the Electrosurgical Bipolar Cable Adapters, (See page 49).

NEEDLES




Interject™ Contrast Single-Use Injection Therapy Needle Catheters*



Order Number	GTIN	Sheath Design	Needle Gauge	Maximum Needle Extension Length (mm)	Sheath O.D. (mm)	Working Length (cm)	Packaging
M00518151	08714729296430	Contrast	23	4	1.8	200	Box 5
M00518161	08714729296447	Contrast	25	4	1.8	200	Box 5
M00518251	08714729296478	Contrast	23	6	1.8	200	Box 5
M00518261	08714729296485	Contrast	25	6	1.8	200	Box 5
M00518351	08714729296522	Contrast	23	4	2.3	240	Box 5
M00518361	08714729296539	Contrast	25	4	2.3	240	Box 5

*With Star Catheter Technology.





Interject Clear Single-Use Injection Therapy Needle Catheters*



Order Number	GTIN	Sheath Design	Needle Gauge	Maximum Needle Extension Length (mm)	Sheath O.D. (mm)	Working Length (cm)	Packaging
M00518101	08714729296416	Clear	23	4	1.8	200	Box 5
M00518111	08714729296423	Clear	25	4	1.8	200	Box 5
M00518301	08714729296508	Clear	23	4	2.3	240	Box 5
M00518311	08714729296515	Clear	25	4	2.3	240	Box 5


*With Star Catheter Technology.

MULTIPLE BAND LIGATORS



Speedband Superview Super 7™ Multiple Band Ligators						
Order Number	GTIN	Description	Required Working Channel (mm)	Scope O.D. (mm)	Bands / Ligating Unit	Packaging
M00542250	08714729201953	Speedband Superview Super 7 Band Ligator	2.8	8.6 – 11.5	7	Each
M00542251	08714729201960	Speedband Superview Super 7 Band Ligator	2.8	8.6 – 11.5	7	Box 2
M00542253	08714729504832	Speedband Superview Super 7 Band Ligator	2.8	8.6 – 11.5	7	Box 4

ELECTROSURGICAL BIPOLAR CABLE ADAPTER



Electrosurgical Bipolar Cable Adapter*

Order Number	GTIN	Description	Pin Spacing (mm)	Generator Compatibility**	Packaging
M00561440	08714729797609	Bipolar Cable Adapter	18	EMC (Coagulasem Dolley)	Each
M00561450	08714729797616	Bipolar Cable Adapter	22	ERBE, EMC Dolley	Each
M00561460	08714729797623	Bipolar Cable Adapter	25	Lamidey (Optima 4)	Each
M00561470	08714729797630	Bipolar Cable Adapter	28	Valleylab™, Lamidey (Duo), Bovie™ (Aaron 3250™), Olympus (ESG100, PSD60, UES10, UES30 & Surgemaster), Bowa, Soring	Each
M00561480	08714729797647	Bipolar Cable Adapter	28 / Third Pin	ConMed™ (System 5000), Bovie (Aaron 2250™), Ethicon™	Each
M00561490	08714729797654	Bipolar Cable Adapter	32 / Third Pin	ConMed (Sabre 180, Excalibur, System 7550)	Each
M00561500	08714729797661	Bipolar Cable Adapter	39	Codman, Malis™	Each

*Bipolar Cable Adapters are used to connect Gold Probe™ and Injection Gold Probe™ Catheters with the above mentioned generators.

Bipolar Cable Adapters are not required for the connection of Gold Probe and Injection Gold Probe to Endostat™ II or III Generators.

** For generators and their associated cable adapters not listed above, please contact your Boston Scientific Representative.

Note: Monopolar active cords are in the Resection & Retrieval Section, (See page 71).

Note: Some variations in generator faceplates may exist within a generator model. We recommend that you confirm pin spacing prior to ordering.

Notes and Future Product Updates



Section Table of Contents



Prev Page

Next Page



Home





Dilation

CRE™ Balloon Dilatation Catheters

- CRE Fixed Wire Esophageal Balloon Dilatation Catheters53
- CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters 51
- CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters 51
- CRE RX Biliary Balloon Dilatation Catheters50
- CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters52
- CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters52

MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter

- MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter54


Rigiflex™ II Single-Use Achalasia Balloon Dilatation Catheter

- Rigiflex II Single-Use Achalasia Balloon Dilators54

Inflation Systems and Dilation Accessories

- Alliance™ II Integrated Inflation / Lithotripsy Device55
- Alliance II Single-Use Syringe / Gauge Assembly55
- CRE SteriFlate Disposable Inflation Device55
- Pneumatic Inflator / Deflation Devices and Accessories55

CRE™ RX BILIARY BALLOON DILATATION CATHETERS




CRE RX Biliary Balloon Dilatation Catheters*

Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558860	08714729965411	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	3.7	7.5	180	Each
M00558870	08714729965428	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	3.7	7.5	180	Each
M00558880	08714729965435	10-11-12	30-33-36	5.5	3-5-8	304-507-811	3.7	7.5	180	Each
M00558890	08714729965442	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	3.7	7.5	180	Each
M00558900	08714729965459	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	3.7	7.5	180	Each
M00558910	08714729965466	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	3.7	7.5	180	Each
M00558920	08714729965473	8-9-10	24-27-30	3.0	3-5.5-9	304-608-1013	3.7	7.5	180	Each
M00558930	08714729965480	10-11-12	30-33-36	3.0	3-5-8	304-557-912	3.7	7.5	180	Each
M00558940	08714729965497	12-13.5-15	36-40.5-45	3.0	3-4.5-8	304-507-811	3.7	7.5	180	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060, (See page 84).

*No Pre-loaded Guidewire in the Balloon Dilatation Catheter.

CRE™ PRO WIREGUIDED BALLOON DILATATION CATHETERS



CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558600	08714729797487	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8 / 3.2*	7.5	180	Each
M00558610	08714729797494	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8 / 3.2*	7.5	180	Each
M00558620	08714729797500	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8 / 3.2*	7.5	180	Each
M00558630	08714729797517	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8 / 3.2*	7.5	180	Each
M00558640	08714729797524	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8 / 3.7*	7.5	180	Each
M00558650**	08714729797531	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8 / 3.7*	7.5	180	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060, (See page 84).

Guidewire included and preloaded in balloon dilatation catheter.

*Duodenoscope channel size.

**Esophageal and Biliary application only.

CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558660	08714729797548	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8 / 3.2*	7.5	240	Each
M00558670	08714729797555	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8 / 3.2*	7.5	240	Each
M00558680	08714729797562	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8 / 3.2*	7.5	240	Each
M00558690	08714729797579	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8 / 3.2*	7.5	240	Each
M00558700	08714729797586	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8 / 3.7*	7.5	240	Each
M00558710**	08714729797593	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8 / 3.7*	7.5	240	Each


Recommended Inflation Device: Alliance II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060, (See page 84).

Guidewire included and preloaded in balloon dilatation catheter.

*Duodenoscope channel size.

**Esophageal, Colonic and Biliary application only.

CRE™ WIREGUIDED BALLOON DILATATION CATHETERS



CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558390	08714729339328	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8	7.5	180	Each
M00558400	08714729339335	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8	7.5	180	Each
M00558410	08714729339342	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8	7.5	180	Each
M00558420	08714729339359	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8	7.5	180	Each
M00558430	08714729339366	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8	7.5	180	Each
M00558440*	08714729339373	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8	7.5	180	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060, (See page 84).

Guidewire included and preloaded in balloon dilatation catheter.

*Esophageal application only.

CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters




Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558450	08714729339380	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8	7.5	240	Each
M00558460	08714729339397	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8	7.5	240	Each
M00558470	08714729339403	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8	7.5	240	Each
M00558480	08714729339410	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8	7.5	240	Each
M00558490	08714729339427	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8	7.5	240	Each
M00558500*	08714729339434	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8	7.5	240	Each

Recommended Inflation Device: Alliance II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060, (See page 84).

Guidewire included and preloaded in balloon dilatation catheter.

*Esophageal and Colonic application only.

CRE™ BALLOON DILATATION CATHETERS (Continued)




CRE Fixed Wire Esophageal Balloon Dilatation Catheters

Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558330	08714729201984	6-7-8	18-21-24	8	3-6-10	304-608-1013	2.8	6	180	Each
M00558340	08714729201991	8-9-10	24-27-30	8	3-5.5-9	304-557-912	2.8	6	180	Each
M00558350	08714729202004	10-11-12	30-33-36	8	3-5-8	304-507-811	2.8	6	180	Each
M00558360	08714729195986	12-13.5-15	36-40.5-45	8	3-4.5-8	304-456-811	2.8	6	180	Each
M00558370	08714729195993	15-16.5-18	45-49.5-54	8	3-4.5-7	304-456-709	2.8	6	180	Each
M00558380	08714729202011	18-19-20	54-57-60	8	3-4.5-6	304-456-608	2.8	6	180	Each
M00558361	08714729339298	12-13.5-15	36-40.5-45	8	3-4.5-8	304-456-811	2.8	6	180	Box 5
M00558371	08714729339304	15-16.5-18	45-49.5-54	8	3-4.5-7	304-456-709	2.8	6	180	Box 5
M00558381	08714729339311	18-19-20	54-57-60	8	3-4.5-6	304-456-608	2.8	6	180	Box 5

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060, (See page 84).

MAXFORCE TTS™ ESOPHAGEAL HIGH PERFORMANCE BALLOON DILATATION CATHETER




MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size		Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)		(mm)	(Fr)		
M00589060	08714729099581	6	18	6	10	1013	2.8	2.0	6	180	Each
M00589080	08714729099598	8	24	6	10	1013	2.8	2.0	6	180	Each
M00589100	08714729099895	10	30	6	7	709	2.8	2.0	6	180	Each
M00589120	08714729099901	12	36	6	7	709	2.8	2.0	6	180	Each
M00589140	08714729080602	14	42	6	7	709	2.8	2.0	6	180	Each
M00589150	08714729099918	15	45	6	5	507	2.8	2.0	6	180	Each
M00589160	08714729080619	16	48	6	5	507	2.8	2.0	6	180	Each
M00589180	08714729099925	18	54	6	5	507	2.8	2.0	6	180	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060, (See page 84).

RIGIFLEX™ II SINGLE-USE ACHALASIA BALLOON DILATATION CATHETER



Rigiflex II Single-Use Achalasia Balloon Dilators*




Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Working Length (cm)	Catheter Size		Recommended Guidewire** (in)	Packaging
		(mm)	(Fr)		(psi)	(kPa)		(mm)	(Fr)		
M00554500	08714729719106	30	90	10	20	138	90	4.6	14	.038	Each
M00554510	08714729719113	35	105	10	20	138	90	4.6	14	.038	Each
M00554520	08714729719120	40	120	10	20	138	90	4.6	14	.038	Each

*Product Shipped Sterile. Recommended use with Pneumatic Inflator Order #5320, (See page 55).

**Guidewire included with Rigiflex II Balloon Dilation Catheter.

INFLATION SYSTEMS AND DILATION ACCESSORIES



Alliance™ II Single-Use Syringe / Gauge Assembly / Integrated Inflation / Lithotripsy Device

Order Number	GTIN	Description	Packaging
M00550620	08714729283881	Alliance II Integrated Inflation / Lithotripsy Device	Each
M00550601	08714729129332	⊗ Alliance II Single-Use Syringe / Gauge Assembly	Box 5

Inflation / Deflation Devices and Accessories

Order Number	GTIN	Description	Packaging
M00553200	08714729195450	Pneumatic Inflator	Each
M00550630	00856877002506	CRE™ SteriFlate Disposable Inflation Device	Each




Biopsy

Forceps

- | | |
|--|---|
| ► Multibite™ Single-Use Multiple Sample Biopsy Forceps57 | ► Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps 56 |
| ► Radial Jaw™ 4 Single-Use Hot Biopsy Forceps57 | ► Radial Jaw 4 Single-Use Pediatric Biopsy Forceps57 |
| ► Radial Jaw 4 Single-Use Jumbo Biopsy Forceps 56 | ► Radial Jaw 4 Single-Use Standard Capacity Biopsy Forceps.... 56 |

FORCEPS



Radial Jaw™ 4 Single-Use Jumbo Biopsy Forceps



Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M00513343	08714729761891	Radial Jaw 4 Jumbo Forceps	2.8	160	3.2	Yellow	Box 5
M00513353	08714729761907	Radial Jaw 4 Jumbo Forceps w/ Needle	2.8	160	3.2	Yellow	Box 5
M00513363	08714729761914	Radial Jaw 4 Jumbo Forceps	2.8	240	3.2	Orange	Box 5
M00513373	08714729761938	Radial Jaw 4 Jumbo Forceps w/ Needle	2.8	240	3.2	Orange	Box 5
M00513361	08714729648871	Radial Jaw 4 Jumbo Forceps	2.8	240	3.2	Orange	Box 20
M00513371	08714729648895	Radial Jaw 4 Jumbo Forceps w/ Needle	2.8	240	3.2	Orange	Box 20
M00513362	08714729648888	Radial Jaw 4 Jumbo Forceps	2.8	240	3.2	Orange	Box 40
M00513372	08714729648901	Radial Jaw 4 Jumbo Forceps w/ Needle	2.8	240	3.2	Orange	Box 40

Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps



Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M00513303	08714729761839	Radial Jaw 4 Large Capacity Forceps	2.4	160	2.8	Yellow	Box 5
M00513313	08714729761853	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	160	2.8	Yellow	Box 5
M00513323	08714729761877	Radial Jaw 4 Large Capacity Forceps	2.4	240	2.8	Orange	Box 5
M00513333	08714729755494	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	240	2.8	Orange	Box 5
M00513301	08714729648758	Radial Jaw 4 Large Capacity Forceps	2.4	160	2.8	Yellow	Box 20
M00513311	08714729648772	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	160	2.8	Yellow	Box 20
M00513321	08714729648796	Radial Jaw 4 Large Capacity Forceps	2.4	240	2.8	Orange	Box 20
M00513331	08714729648819	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	240	2.8	Orange	Box 20
M00513302	08714729648765	Radial Jaw 4 Large Capacity Forceps	2.4	160	2.8	Yellow	Box 40
M00513312	08714729648789	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	160	2.8	Yellow	Box 40
M00513322	08714729648802	Radial Jaw 4 Large Capacity Forceps	2.4	240	2.8	Orange	Box 40
M00513332	08714729648826	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	240	2.8	Orange	Box 40

Radial Jaw 4 Single-Use Standard Capacity Biopsy Forceps



Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M00513383	08714729777915	Radial Jaw 4 Standard Capacity Forceps	2.2	160	2.8	Yellow	Box 5
M00513393	08714729777939	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	160	2.8	Yellow	Box 5
M00513403	08714729777922	Radial Jaw 4 Standard Capacity Forceps	2.2	240	2.8	Orange	Box 5
M00513413	08714729777946	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	240	2.8	Orange	Box 5
M00513381	08714729777953	Radial Jaw 4 Standard Capacity Forceps	2.2	160	2.8	Yellow	Box 20
M00513391	08714729777977	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	160	2.8	Yellow	Box 20
M00513401	08714729777960	Radial Jaw 4 Standard Capacity Forceps	2.2	240	2.8	Orange	Box 20
M00513411	08714729777984	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	240	2.8	Orange	Box 20
M00513382	08714729778011	Radial Jaw 4 Standard Capacity Forceps	2.2	160	2.8	Yellow	Box 40
M00513392	08714729777991	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	160	2.8	Yellow	Box 40
M00513402	08714729778028	Radial Jaw 4 Standard Capacity Forceps	2.2	240	2.8	Orange	Box 40
M00513412	08714729778004	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	240	2.8	Orange	Box 40



FORCEPS (Continued)

Radial Jaw™ 4 Single-Use Pediatric Biopsy Forceps



Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M00513443	08714729790211	Radial Jaw 4 Pediatric Biopsy Forceps	1.8	160	2.0	Yellow	Box 5
M00513453	08714729790228	Radial Jaw 4 Pediatric Biopsy Forceps w/ Needle	1.8	160	2.0	Yellow	Box 5

HOT BIOPSY FORCEPS




Radial Jaw 4 Single-Use Hot Biopsy Forceps



Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M00515031	08714729790136	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	Red	Box 5
M00515032	08714729790105	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	Red	Box 20
M00515033	08714729790143	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	Red	Box 40

FORCEPS



Multibite™ Single-Use Multiple Sample Biopsy Forceps



Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M00510101	08714729164975	Multibite Multiple Sample Biopsy Forceps	2.4	160	2.8	Yellow	Box 5
M00510121	08714729164982	Multibite Multiple Sample Biopsy Forceps	2.4	240	2.8	Orange	Box 5

Notes and Future Product Updates

Notes and Future Product Updates



Section Table of Contents



Prev Page

Next Page



Home





EndoVive™ Enteral Access Devices

Percutaneous Endoscopic Gastrostomy (PEG) Devices with ENFit™

Initial Placement Peg Kits with ENFit

- Safety PEG Kit with ENFit58
- Standard PEG Kit with ENFit.....58
- Initial Placement PEG Accessories with ENFit.....58

One-Step Button™ Low Profile Initial Placement PEG Kits and Accessories with ENFit

- One-Step Button Low Profile PEG Kits with ENFit.....59
- One-Step Button and Low Profile Replacement Button Accessories with ENFit.....59

Replacement Gastrostomy Tubes (G-Tube) with ENFit

- Replacement G-Tubes, Straight Bolster, with ENFit60
- Replacement G-Tubes, Right Angle Bolster, with ENFit.....60

Low Profile Replacement Buttons and Accessories with ENFit

- Low Profile Replacement Button with ENFit.....60
- One-Step Button and Low Profile Replacement Button Accessories with ENFit..... 61

ENFit Transition Connectors

- ENFit Transition Connectors 61

Percutaneous Endoscopic Gastrostomy (PEG) Devices

Initial Placement PEG Kits

- Safety PEG Kit.....62
- Standard PEG Kit62
- Initial Placement PEG Accessories62

Initial Placement Low Profile PEG Kits

- One-Step Button Low Profile PEG Kits63
- One-Step Button and Low Profile Replacement Button Accessories63

Replacement Gastrostomy Tubes (G-Tube)

- Replacement G-Tubes, Straight Bolster64
- Replacement G-Tubes, Right Angle Bolster.....64


Low Profile Replacement Buttons and Accessories

- Low Profile Replacement Button65
- One-Step Button and Low Profile Replacement Button Accessories65

Jejunal Feeding Tube Devices


- Through-The-PEG (TTP) Jejunal Feeding Tube66

INITIAL PLACEMENT PEG KITS WITH ENFit™†



Safety PEG* Kit with ENFit				
Order Number	GTIN	Description		Packaging
M00509001	08714729880394	20Fr Safety PEG Kit, Pull with ENFit		Box 2
M00509011	08714729880417	20Fr Safety PEG Kit, Push with ENFit		Box 2
M00509021	08714729880431	24Fr Safety PEG Kit, Pull with ENFit		Box 2
M00509031	08714729880455	24Fr Safety PEG Kit, Push with ENFit		Box 2

*Includes local anesthetic.



Standard PEG* Kit with ENFit				
Order Number	GTIN	Description		Packaging
M00509041	08714729880479	20Fr Standard PEG Kit, Pull with ENFit		Box 2
M00509051	08714729880493	20Fr Standard PEG Kit, Push with ENFit		Box 2
M00509061	08714729880516	24Fr Standard PEG Kit, Pull with ENFit		Box 2
M00509071	08714729880530	24Fr Standard PEG Kit, Push with ENFit		Box 2


*Includes local anesthetic.

Initial Placement PEG Accessories with ENFit				
Order Number	GTIN	Description		Packaging
M00501261	08714729863564	20Fr Y-Ports with ENFit		Box 10
M00501271	08714729863588	24Fr Y-Ports with ENFit		Box 10

†NOTE: Make sure all appropriate service lines in your healthcare system are ready to use ENFit enteral connectors. Please contact your Boston Scientific Sales Representative for support with ENFit enteral connectors.

Recommended ENFit Transition Connectors [▶](#) on Page 61.

ONE-STEP BUTTON™ LOW PROFILE INITIAL PLACEMENT PEG KITS WITH ENFit™†



One-Step Button Low Profile PEG Kits with ENFit 				
Order Number	GTIN	Description	Length (cm)	Packaging
M00509690	08714729897415	18Fr Low Profile PEG Kit, Pull with ENFit	1.2	Each
M00509700	08714729897422	18Fr Low Profile PEG Kit, Pull with ENFit	1.7	Each
M00509710	08714729897439	18Fr Low Profile PEG Kit, Pull with ENFit	2.4	Each
M00509720	08714729897446	18Fr Low Profile PEG Kit, Pull with ENFit	3.4	Each
M00509730	08714729897453	18Fr Low Profile PEG Kit, Pull with ENFit	4.4	Each
M00509740	08714729897491	18Fr Low Profile PEG Kit, Push with ENFit	1.7	Each
M00509750	08714729897507	18Fr Low Profile PEG Kit, Push with ENFit	2.4	Each
M00509580	08714729897378	24Fr Low Profile PEG Kit, Pull with ENFit	1.7	Each
M00509590	08714729897385	24Fr Low Profile PEG Kit, Pull with ENFit	2.4	Each
M00509600	08714729897392	24Fr Low Profile PEG Kit, Pull with ENFit	3.4	Each
M00509610	08714729897408	24Fr Low Profile PEG Kit, Pull with ENFit	4.4	Each
M00509630	08714729897460	24Fr Low Profile PEG Kit, Push with ENFit	2.4	Each
M00509620	08714729897477	24Fr Low Profile PEG Kit, Push with ENFit	3.4	Each
M00509680	08714729897484	24Fr Low Profile PEG Kit, Push with ENFit	4.4	Each

One-Step Button and Low Profile Replacement Button Accessories with ENFit 				
Order Number	GTIN	Description	Length (cm)	Packaging
M00501121	08714729863380	18Fr Decompression Tube with ENFit	1.2	Box 10
M00501131	08714729863403	18Fr Decompression Tube with ENFit	1.7	Box 10
M00501141	08714729863427	18Fr Decompression Tube with ENFit	2.4	Box 10
M00501151	08714729863441	18Fr Decompression Tube with ENFit	3.4	Box 10
M00501161	08714729863465	18Fr Decompression Tube with ENFit	4.4	Box 10
M00501181	08714729863489	24Fr Decompression Tube with ENFit	1.7	Box 10
M00501191	08714729863502	24Fr Decompression Tube with ENFit	2.4	Box 10
M00501201	08714729863526	24Fr Decompression Tube with ENFit	3.4	Box 10
M00501211	08714729863540	24Fr Decompression Tube with ENFit	4.4	Box 10
M00501091	08714729863335	18Fr Bolus Feeding Set Straight Adapter with ENFit	–	Box 10
M00501061	08714729863281	18Fr Right Angle Feeding Set Continuous Adapter with ENFit	–	Box 10
M00501101	08714729863359	24Fr Bolus Feeding Set Straight Adapter with ENFit	–	Box 10
M00501071	08714729863304	24Fr Right Angle Feeding Set Continuous Adapter with ENFit	–	Box 10

†NOTE: Make sure all appropriate service lines in your healthcare system are ready to use ENFit enteral connectors. Please contact your Boston Scientific Sales Representative for support with ENFit enteral connectors.

Recommended ENFit Transition Connectors  on Page 61.

REPLACEMENT GASTROSTOMY TUBES WITH ENFIT™†




Replacement Gastrostomy Feeding Tube with ENFit, Straight Bolster				
Order Number	GTIN	Description		Packaging
M00509881	00813939023862	12Fr Replacement G-Tube, Straight with ENFit		Box 5
M00509891	00813939023879	14Fr Replacement G-Tube, Straight with ENFit		Box 5
M00509901	00813939023886	16Fr Replacement G-Tube, Straight with ENFit		Box 5
M00509911	00813939023893	18Fr Replacement G-Tube, Straight with ENFit		Box 5
M00509921	00813939023909	20Fr Replacement G-Tube, Straight with ENFit		Box 5
M00509931	00813939023916	22Fr Replacement G-Tube, Straight with ENFit		Box 5
M00509941	00813939023923	24Fr Replacement G-Tube, Straight with ENFit		Box 5



Replacement Gastrostomy Feeding Tube with ENFit, Right Angle Bolster				
Order Number	GTIN	Description		Packaging
M00509961	00813939024005	12Fr Replacement G-Tube, Right Angle with ENFit		Box 5
M00509971	00813939024012	14Fr Replacement G-Tube, Right Angle with ENFit		Box 5
M00509981	00813939024029	16Fr Replacement G-Tube, Right Angle with ENFit		Box 5
M00509991	00813939024036	18Fr Replacement G-Tube, Right Angle with ENFit		Box 5
M00510001	00813939024043	20Fr Replacement G-Tube, Right Angle with ENFit		Box 5
M00510011	00813939024050	22Fr Replacement G-Tube, Right Angle with ENFit		Box 5
M00510021	00813939024067	24Fr Replacement G-Tube, Right Angle with ENFit		Box 5

LOW PROFILE REPLACEMENT BUTTONS & ACCESSORIES WITH ENFIT




Low Profile Replacement Button with ENFit				
Order Number	GTIN	Description	Length (cm)	Packaging
M00509830	08714729902614	18Fr Low Profile Replacement Button with ENFit	1.2	Each
M00509840	08714729902621	18Fr Low Profile Replacement Button with ENFit	1.7	Each
M00509850	08714729902638	18Fr Low Profile Replacement Button with ENFit	2.4	Each
M00509860	08714729902645	18Fr Low Profile Replacement Button with ENFit	3.4	Each
M00509870	08714729902652	18Fr Low Profile Replacement Button with ENFit	4.4	Each
M00509410	08714729902669	24Fr Low Profile Replacement Button with ENFit	1.7	Each
M00509420	08714729902676	24Fr Low Profile Replacement Button with ENFit	2.4	Each
M00509430	08714729902683	24Fr Low Profile Replacement Button with ENFit	3.4	Each
M00509440	08714729902690	24Fr Low Profile Replacement Button with ENFit	4.4	Each

†NOTE: Make sure all appropriate service lines in your healthcare system are ready to use ENFit enteral connectors. Please contact your Boston Scientific Sales Representative for support with ENFit enteral connectors.


Recommended ENFit Transition Connectors  on Page 61.

LOW PROFILE REPLACEMENT BUTTON ACCESSORIES WITH ENFit™†



One-Step Button™ and Low Profile Replacement Button Accessories with ENFit				
Order Number	GTIN	Description	Length (cm)	Packaging
M00501121	08714729863380	18Fr Decompression Tube with ENFit	1.2	Box 10
M00501131	08714729863403	18Fr Decompression Tube with ENFit	1.7	Box 10
M00501141	08714729863427	18Fr Decompression Tube with ENFit	2.4	Box 10
M00501151	08714729863441	18Fr Decompression Tube with ENFit	3.4	Box 10
M00501161	08714729863465	18Fr Decompression Tube with ENFit	4.4	Box 10
M00501181	08714729863489	24Fr Decompression Tube with ENFit	1.7	Box 10
M00501191	08714729863502	24Fr Decompression Tube with ENFit	2.4	Box 10
M00501201	08714729863526	24Fr Decompression Tube with ENFit	3.4	Box 10
M00501211	08714729863540	24Fr Decompression Tube with ENFit	4.4	Box 10
M00501091	08714729863335	18Fr Bolus Feeding Set Straight Adapter with ENFit	–	Box 10
M00501061	08714729863281	18Fr Right Angle Feeding Set Continuous Adapter with ENFit	–	Box 10
M00501101	08714729863359	24Fr Bolus Feeding Set Straight Adapter with ENFit	–	Box 10
M00501071	08714729863304	24Fr Right Angle Feeding Set Continuous Adapter with ENFit	–	Box 10


ENFit TRANSITION CONNECTORS



ENFit Transition Connectors			
Order Number	GTIN	Description	Packaging
M00509241	10813939024972	ENFit Feeding Port Transition Connector	Box 10


†NOTE: Make sure all appropriate service lines in your healthcare system are ready to use ENFit enteral connectors. Please contact your Boston Scientific Sales Representative for support with ENFit enteral connectors.

INITIAL PLACEMENT PEG KITS



Safety PEG* Kit				
Order Number	GTIN	Description		Packaging
M00566461	08714729802624	20Fr Safety PEG Kit, Pull		Box 2
M00566471	08714729802631	20Fr Safety PEG Kit, Push		Box 2
M00566481	08714729802648	24Fr Safety PEG Kit, Pull		Box 2
M00566491	08714729802655	24Fr Safety PEG Kit, Push		Box 2

*Includes local anesthetic.





Standard PEG* Kit				
Order Number	GTIN	Description		Packaging
M00568201	08714729285199	20Fr Standard PEG Kit, Pull		Box 2
M00568211	08714729285205	20Fr Standard PEG Kit, Push		Box 2
M00568241	08714729296799	24Fr Standard PEG Kit, Pull		Box 2
M00568251	08714729296805	24Fr Standard PEG Kit, Push		Box 2

*Includes local anesthetic.

Initial Placement PEG Accessories				
Order Number	GTIN	Description		Packaging
M00580641	08714729169116	20Fr Y-Ports		Box 10
M00580651	08714729169123	24Fr Y-Ports		Box 10

INITIAL PLACEMENT LOW PROFILE PEG KITS




One-Step Button™ Low Profile PEG Kits				
Order Number	GTIN	Description	Length (cm)	Packaging
M00563120	08714729284246	18Fr Low Profile PEG Kit, Pull	1.2	Each
M00563130	08714729284253	18Fr Low Profile PEG Kit, Pull	1.7	Each
M00563140	08714729284260	18Fr Low Profile PEG Kit, Pull	2.4	Each
M00563150	08714729284277	18Fr Low Profile PEG Kit, Pull	3.4	Each
M00563160	08714729284284	18Fr Low Profile PEG Kit, Pull	4.4	Each
M00563180	08714729284307	18Fr Low Profile PEG Kit, Push	1.7	Each
M00563190	08714729284314	18Fr Low Profile PEG Kit, Push	2.4	Each
M00563010	08714729284147	24Fr Low Profile PEG Kit, Pull	1.7	Each
M00563020	08714729284154	24Fr Low Profile PEG Kit, Pull	2.4	Each
M00563030	08714729284161	24Fr Low Profile PEG Kit, Pull	3.4	Each
M00563040	08714729268833	24Fr Low Profile PEG Kit, Pull	4.4	Each
M00563080	08714729284208	24Fr Low Profile PEG Kit, Push	2.4	Each
M00563090	08714729284215	24Fr Low Profile PEG Kit, Push	3.4	Each
M00563100	08714729284222	24Fr Low Profile PEG Kit, Push	4.4	Each

One-Step Button and Low Profile Replacement Button Accessories				
Order Number	GTIN	Description	Length (cm)	Packaging
M00580401	08714729166443	18Fr Decompression Tube	1.2	Box 10
M00580411	08714729166450	18Fr Decompression Tube	1.7	Box 10
M00580421	08714729166467	18Fr Decompression Tube	2.4	Box 10
M00580431	08714729166474	18Fr Decompression Tube	3.4	Box 10
M00580441	08714729166481	18Fr Decompression Tube	4.4	Box 10
M00580461	08714729166504	24Fr Decompression Tube	1.7	Box 10
M00580471	08714729166511	24Fr Decompression Tube	2.4	Box 10
M00580481	08714729339649	24Fr Decompression Tube	3.4	Box 10
M00580491	08714729166535	24Fr Decompression Tube	4.4	Box 10
M00580171	08714729142232	18Fr Bolus Feeding Set Straight Adapter	—	Box 10
M00580181	08714729142249	18Fr Right Angle Feeding Set Continuous Adapter	—	Box 10
M00580231	08714729166405	24Fr Bolus Feeding Set Straight Adapter	—	Box 10
M00580241	08714729166412	24Fr Right Angle Feeding Set Continuous Adapter	—	Box 10

REPLACEMENT GASTROSTOMY TUBES




Replacement Gastrostomy Feeding Tubes, Straight Bolster				
Order Number	GTIN	Description		Packaging
M00582040	00813939020243	12Fr Replacement G-Tube, Straight		Each
M00582050	00813939020250	14Fr Replacement G-Tube, Straight		Each
M00582060	00813939020267	16Fr Replacement G-Tube, Straight		Each
M00582070	00813939020274	18Fr Replacement G-Tube, Straight		Each
M00582080	00813939020281	20Fr Replacement G-Tube, Straight		Each
M00582860	00813939020298	22Fr Replacement G-Tube, Straight		Each
M00582090	00813939020304	24Fr Replacement G-Tube, Straight		Each
M00582100	00813939020311	28Fr Replacement G-Tube, Straight		Each



Replacement Gastrostomy Feeding Tubes, Right Angle Bolster				
Order Number	GTIN	Description		Packaging
M00582110	00813939020007	12Fr Replacement G-Tube, Right Angle		Each
M00582120	00813939020014	14Fr Replacement G-Tube, Right Angle		Each
M00582130	00813939020021	16Fr Replacement G-Tube, Right Angle		Each
M00582140	00813939020038	18Fr Replacement G-Tube, Right Angle		Each
M00582150	00813939020045	20Fr Replacement G-Tube, Right Angle		Each
M00582870	00813939020052	22Fr Replacement G-Tube, Right Angle		Each
M00582160	00813939020069	24Fr Replacement G-Tube, Right Angle		Each
M00582170	00813939020076	28Fr Replacement G-Tube, Right Angle		Each


LOW PROFILE REPLACEMENT BUTTONS & ACCESSORIES



Low Profile Replacement Button 				
Order Number	GTIN	Description	Length (cm)	Packaging
M00562760	08714729284017	18Fr Low Profile Replacement Button	1.2	Each
M00562770	08714729142201	18Fr Low Profile Replacement Button	1.7	Each
M00562780	08714729284024	18Fr Low Profile Replacement Button	2.4	Each
M00562790	08714729284031	18Fr Low Profile Replacement Button	3.4	Each
M00562800	08714729284048	18Fr Low Profile Replacement Button	4.4	Each
M00562820	08714729284062	24Fr Low Profile Replacement Button	1.7	Each
M00562830	08714729142218	24Fr Low Profile Replacement Button	2.4	Each
M00562840	08714729284079	24Fr Low Profile Replacement Button	3.4	Each
M00562850	08714729284086	24Fr Low Profile Replacement Button	4.4	Each

One-Step Button™ and Low Profile Replacement Button Accessories 				
Order Number	GTIN	Description	Length (cm)	Packaging
M00580401	08714729166443	18Fr Decompression Tube	1.2	Box 10
M00580411	08714729166450	18Fr Decompression Tube	1.7	Box 10
M00580421	08714729166467	18Fr Decompression Tube	2.4	Box 10
M00580431	08714729166474	18Fr Decompression Tube	3.4	Box 10
M00580441	08714729166481	18Fr Decompression Tube	4.4	Box 10
M00580461	08714729166504	24Fr Decompression Tube	1.7	Box 10
M00580471	08714729166511	24Fr Decompression Tube	2.4	Box 10
M00580481	08714729339649	24Fr Decompression Tube	3.4	Box 10
M00580491	08714729166535	24Fr Decompression Tube	4.4	Box 10
M00580171	08714729142232	18Fr Bolus Feeding Set Straight Adapter	—	Box 10
M00580181	08714729142249	18Fr Right Angle Feeding Set Continuous Adapter	—	Box 10
M00580231	08714729166405	24Fr Bolus Feeding Set Straight Adapter	—	Box 10
M00580241	08714729166412	24Fr Right Angle Feeding Set Continuous Adapter	—	Box 10

JEJUNAL FEEDING TUBE DEVICES



Through-The-PEG (TTP) Jejunal Feeding Tube					
Order Number	GTIN	Description	Tip Style	Packaging	
M00566320	08714729786313	8.5Fr TTP Jejunal Feeding Tube	Bent Tip	Each	
M00566330	08714729789727	8.5Fr TTP Jejunal Feeding Tube	Pigtail Tip	Each	
M00566340	08714729789666	12Fr TTP Jejunal Feeding Tube	Bent Tip	Each	
M00566350	08714729786252	12Fr TTP Jejunal Feeding Tube	Pigtail Tip	Each	

Notes and Future Product Updates



Section Table of Contents



Prev Page

Next Page



Home



► Resection & Retrieval



Resection Devices

EMR for Upper GI

- Captivator™ EMR Device 67

Stiff Snares

- Captivator COLD Single-Use Snares 67
- Captivator II Single-Use Snares 68
- Captivator Single-Use Snares 68

Specialty Snares

- Profile™ Single-Use Snares 69
- Single-Use Rotatable Snares 69

Snares

- Captiflex™ Single-Use Snares 70
- Sensation™ Single-Use Short Throw Snares 70

Active Cords

- Active Cords 71

Submucosal Lifting

- ORISE™ Gel 71

Electrosurgical Knife

- ORISE ProKnife 71

Retrieval Devices

Retrieval Device

- RescueNet™ Retrieval Device 72

Retrieval Grasping Forceps


- Rescue™ Rat Tooth Grasping Forceps 72
- Rescue Rat Tooth Grasping Forceps - Pediatric 72
- Rescue Alligator Grasping Forceps, Long 72
- Rescue Alligator Grasping Forceps 72
- Rescue Rat Tooth / Alligator Grasping Forceps 72

Polyp Traps

- TrapEase™ Single Chamber Polyp Trap 72
- TrapEase Quad Chamber Polyp Trap 72

► Resection Devices

EMR FOR UPPER GI




Captivator™ EMR Device



Order Number	GTIN	Description	Packaging
M00561600	08714729842675	Captivator EMR Device for Standard Gastrosopes	Box 1
M00561610	08714729842699	Captivator EMR Device for Large Gastrosopes	Box 1
M00561601	08714729842682	Captivator EMR Device for Standard Gastrosopes	Box 5
M00561611	08714729842705	Captivator EMR Device for Large Gastrosopes	Box 5

STIFF SNARES



Captivator COLD Single-Use Snare




Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00561101	08714729971269	10mm Rounded – Stiff	10	2.4	2.8	240	Box 10
M00561102	08714729984979	10mm Rounded – Stiff	10	2.4	2.8	240	Box 20
M00561103	08714729984986	10mm Rounded – Stiff	10	2.4	2.8	240	Box 40



Resection Devices


STIFF SNARES (Continued)



Captivator™ II Single-Use Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00561221	08714729855903	10mm Rounded – Stiff	10	2.4	2.8	240	Box 10
M00561222	08714729861300	10mm Rounded – Stiff	10	2.4	2.8	240	Box 20
M00561223	08714729861317	10mm Rounded – Stiff	10	2.4	2.8	240	Box 40
M00561231	08714729855910	15mm Rounded – Stiff	15	2.4	2.8	240	Box 10
M00561232	08714729861331	15mm Rounded – Stiff	15	2.4	2.8	240	Box 20
M00561233	08714729861348	15mm Rounded – Stiff	15	2.4	2.8	240	Box 40
M00561241	08714729855927	20mm Rounded – Stiff	20	2.4	2.8	240	Box 10
M00561242	08714729861362	20mm Rounded – Stiff	20	2.4	2.8	240	Box 20
M00561243	08714729861379	20mm Rounded – Stiff	20	2.4	2.8	240	Box 40
M00561191	08714729855934	25mm Rounded – Stiff	25	2.4	2.8	240	Box 10
M00561192	08714729861270	25mm Rounded – Stiff	25	2.4	2.8	240	Box 20
M00561193	08714729861287	25mm Rounded – Stiff	25	2.4	2.8	240	Box 40
M00561291	08714729645801	Extra Large Rounded – Stiff	33	2.4	2.8	240	Box 10




Captivator Single-Use Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562301	08714729019312	Small Oval – Stiff	13	2.4	2.8	240	Box 10
M00562321	08714729019336	Medium Oval – Stiff	27	2.4	2.8	240	Box 10
M00561311	08714729019251	Large Oval – Medium Stiff	30	2.4	2.8	240	Box 10
M00562391	08714729019404	Large Oval (thin wire) – Flexible	30	2.4	2.8	240	Box 10
M00562451	08714729071068	Small Hexagonal – Stiff	13	2.4	2.8	240	Box 10
M00562341	08714729019350	Medium Hexagonal – Stiff	27	2.4	2.8	240	Box 10
M00562371	08714729019381	Medium Crescent – Stiff	27	2.4	2.8	240	Box 10

► Resection Devices


SPECIALITY SNARES



Profile™ Single-Use Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562531	08714729156628	Extra Small Oval – Flexible	11	1.9	2.0	240	Box 10
M00562551	08714729156659	Small Oval – Flexible	13	1.9	2.0	240	Box 10
M00562571	08714729156635	Medium Oval – Flexible	27	1.9	2.0	240	Box 10



Single-Use Rotatable Snares*




Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00561821	08714729283980	Small Oval – Medium Stiff	13	2.4	2.8	240	Box 10
M00561831	08714729283997	Medium Oval – Medium Stiff	20	2.4	2.8	240	Box 10

*With Star Catheter Technology.

Resection Devices


SNARES



Captiflex™ Single-Use Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562471	08714729019459	Extra Small Oval – Flexible	11	2.4	2.8	240	Box 10
M00562421	08714729019435	Small Oval – Flexible	13	2.4	2.8	240	Box 10
M00562422	08714729501657	Small Oval – Flexible	13	2.4	2.8	240	Box 40
M00562401	08714729019411	Medium Oval – Flexible	27	2.4	2.8	240	Box 10
M00562402	08714729501640	Medium Oval – Flexible	27	2.4	2.8	240	Box 40



Sensation™ Single-Use Short Throw Snares




Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562691	08714729158158	Small Oval – Flexible	13	2.4	2.8	240	Box 5
M00562692	08714729268826	Small Oval – Flexible	13	2.4	2.8	240	Box 20
M00562693	08714729501671	Small Oval – Flexible	13	2.4	2.8	240	Box 40
M00560321	08714729283928	Small Oval – Medium Stiff	13	2.4	2.8	240	Box 5
M00562671	08714729158141	Medium Oval – Flexible	27	2.4	2.8	240	Box 5
M00562672	08714729268819	Medium Oval – Flexible	27	2.4	2.8	240	Box 20
M00562673	08714729501664	Medium Oval – Flexible	27	2.4	2.8	240	Box 40
M00560311	08714729283904	Medium Oval – Medium Stiff	27	2.4	2.8	240	Box 5
M00562651	08714729158110	Large Oval – Medium Stiff	30	2.4	2.8	240	Box 5
M00562652	08714729268802	Large Oval – Medium Stiff	30	2.4	2.8	240	Box 20
M00560331	08714729283935	Medium Crescent – Medium Stiff	27	2.4	2.8	240	Box 5

► Resection Devices

ACTIVE CORDS


Active Cords			
Order Number	GTIN	Description	Packaging
M00561250	08714729019060	Active Cord, Fits Olympus™ Brand Generator	Each
M00561270	08714729053361	Active Cord, Fits Endostat™, Endostat II and Endostat III Generators, Bovie™, Valleylab™, Aspen Labs™, and ERBE™ Brand Generators	Each

SUBMUCOSAL LIFTING



ORISE™ Gel			
Order Number	GTIN	Description	Packaging
M00519221	08714729993841	ORISE Gel Syringe Single Pack (one 10ml syringe per tray)	Box 10
M00519201	08714729974574	ORISE Gel Syringe Twin Pack (two 10mL syringes per tray)	Box 10
M00519210	08714729974581	ORISE Gel Syringe Kit Twin Pack with Interject™ Injection Therapy Needle Catheter 23ga Clear (two 10mL syringes per kit)	Box 1
M00519231	08714729993865	ORISE Gel Syringe Kit Single Pack with Interject Injection Therapy Needle Catheter 23ga Clear (one 10ml syringe per kit)	Box 10
M00519211	08714729974598	ORISE Gel Syringe Kit Twin Pack with Interject Injection Therapy Needle Catheter 23ga Clear (two 10mL syringes per kit)	Box 10


ELECTROSURGICAL KNIFE



ORISE™ ProKnife							
Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00519371	08714729999904	3.0mm ORISE ProKnife	10	2.4	2.8	240	Each
M00519400	08714729995609	3.0mm ORISE ProKnife, Twin Pack ORISE Gel, Interject Injection Therapy Needle Catheter	10	2.4	2.8	240	Each
M00519351	08714729999881	1.5mm ORISE ProKnife	10	2.4	2.8	240	Each
M00519390	08714729995593	2.0mm ORISE ProKnife, Twin Pack ORISE Gel, Interject Injection Therapy Needle Catheter	10	2.4	2.8	240	Each
M00519361	08714729999898	2.0mm ORISE ProKnife	10	2.4	2.8	240	Each
M00519380	08714729995586	1.5mm ORISE ProKnife, Twin Pack ORISE Gel, Interject Injection Therapy Needle Catheter	10	2.4	2.8	240	Each

► Retrieval Devices

RETRIEVAL DEVICE




RescueNet™ Retrieval Device



Order Number	GTIN	Description	Net Size (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
DGN-538-5	10840253107425	RescueNet Retrieval Device	30	2.5	2.8	230	Box 5

RETRIEVAL GRASPING FORCEPS



Rescue™ Retrieval Single-Use Grasping Forceps



Order Number	GTIN	Description	Jaw Opening Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00527651	26948318968637	Rescue Rat Tooth Grasping Forceps	7	2.4	2.8	230	Box 5
M00527751	26948318968644	Rescue Rat Tooth Grasping Forceps - Pediatric	4.5	1.8	2.0	180	Box 5
M00527551	26948318968620	Rescue Alligator Grasping Forceps, Long	9	2.4	2.8	230	Box 5
M00527451	26948318968613	Rescue Alligator Grasping Forceps	8	2.4	2.8	230	Box 5
M00538351	26948318968606	Rescue Rat Tooth / Alligator Grasping Forceps	8	2.4	2.8	230	Box 5


POLYP TRAPS



TrapEase™ Polyp Traps

Order Number	GTIN	Description	Packaging
STE-297-50	00816849010182	TrapEase Single Chamber Polyp Trap	Box 50
STE-298-25	00816849010199	TrapEase Quad Chamber Polyp Trap	Box 25

► Infection Prevention



Infection Prevention

Cleaning Brushes

► HedgeHog™ Valve Brush	73
► HedgeHog Double Bundle.....	73
► HedgeHog Double-End Channel Brush	73
► HedgeHog Dual-End Channel / Valve Brush	73
► HedgeHog Sweeper Channel Brush	73
► HedgeHog Single-End Channel Brush.....	73

Detergent

► Pure™ Enzymatic Detergent - Spearmint Fresh Scent	73
► Pure Enzymatic Detergent - Fragrance Free.....	73

Bite Blocks

► Blox™ Endoscopic Bite Block with Elastic Cloth Strap (54 FR)	73
► Blox Endoscopic Bite Block with Elastic Cloth Strap (54 FR) (w/o Dental Retention).....	73
► Blox Endoscopic Bite Block with Elastic Cloth Strap (60 FR) ..	73

Lubricant

► Endoglide™ + Sterile Bacteriostatic Lubricating Gel	74
► Sterile Lubricating Jelly	74

Gowns

► EndoArmor+ Non-Sterile Surgical Gown.....	74
---	----

Irrigation & Insufflation

► Hydra™ Disposable Auxiliary Water Jet Connector.....	74
► Hydra Water Bottle Cap	74
► Hydra Irrigation Tubing.....	74

Scope Transport

► CinchPad™ Endoscope Transport Pad.....	75
--	----

Sponges

► Koala™ Cleaning Sponge	75
► Koala Cleaning Sponge plus Pure Enzymatic	75

Suction Tubing

► Suction Tubing	75
------------------------	----

Valves

► Orca™ Sterile Single-Use Air / Water & Suction Valves.....	76
► OrcaPod™ 3 and 4	76
► Seal™ Single-Use Biopsy Valve	76

Compliance EndoKit™

► Compliance EndoKit	76
----------------------------	----

Aegis One-Hour Indicator

► Aegis One-Hour Indicator	76
----------------------------------	----


CLEANING BRUSHES

HedgeHog™ Cleaning Brushes 						
Order Number	GTIN	Description	Catheter/Handle Length (cm)	Brush Length (mm)	Brush Width (mm)	Packaging
SBB-119-50	00816849010274	HedgeHog Valve Brush	6	40, 35	10, 5	Box 50
SBB-382-50	00816849016900	HedgeHog Double Bundle	240, 6	13, 40 / 35	5, 10 / 5	Box 50
SBD-228-50	00816849010083	HedgeHog Double-End Channel Brush	240	13	5	Box 50
SBD-289-50	00816849010090	HedgeHog Dual-End Channel / Valve Brush	45mm, 240	40, 13	10, 5	Box 50
SBD-291-50	00816849016443	HedgeHog Sweeper Channel Brush	240	13, 120	13, 5.5	Box 50
SBS-227-50	00816849010076	HedgeHog Single-End Channel Brush	240	13	5	Box 50

DETERGENT

Pure™ Enzymatic Detergents				
Order Number	GTIN	Description	Packaging	
SEE-220-4	10840253105308	Pure Enzymatic Detergent - Spearmint Fresh Scent	Box 4	
SEE-573-4	10840253105292	Pure Enzymatic Detergent - Fragrance Free	Box 4	

BITE BLOCKS



Blox™ Bite Blocks with Adjustable Straps					
Order Number	GTIN	Description			Packaging
SBT-114-100	00816849010168	Blox Endoscopic Bite Block with Elastic Cloth Strap (54 FR)			Box 100
SBT-321-100	00816849016764	Blox Endoscopic Bite Block with Elastic Cloth Strap (54 FR) (w/o Dental Retention)			Box 100
SBT-546-100	00840253104809	Blox Endoscopic Bite Block with Elastic Cloth Strap (60 FR)			Box 100

► Infection Prevention


LUBRICANT

Endoglide™ + 			
Order Number	GTIN	Description	Packaging
SLT-394-100	00816849017167	Endoglide+ Sterile Bacteriostatic Lubricating Gel (32g / 1.1oz)	Box 100
SLT-395-50	00840253103055	Endoglide+ Sterile Bacteriostatic Lubricating Gel (56.7g / 2oz)	Box 50
SLT-612-50	20851809002900	Sterile Lubricating Jelly (56.7g / 2oz)	Box 50

GOWNS

EndoArmor™ Gown 			
Order Number	GTIN	Description	Packaging
M00501972	00191506005924	EndoArmor+ Non-Sterile Surgical Gown	Box 100

IRRIGATION & INSUFFLATION



Hydra™ Irrigation System							
Order Number	GTIN	Description	Scope Compatibility	CO2	Air/Water Tubing Length (in)	CO2 Tubing Length (in)	Packaging
SIT-361-100	00816849016955	 Hydra Disposable Auxiliary Water Jet Connector	Olympus	N/A	N/A	N/A	Box 100
SCT-350-10	00816849016887	Hydra Water Bottle Cap System	Olympus	No	16	N/A	Box 10
SCT-351-10	00816849016894	Hydra Water Bottle Cap System	Olympus	Yes	16	18	Box 10
SIT-355-15	00840253112200	Hydra Irrigation Tubing System (compatible with ERBE, Medivators pumps)	N/A	N/A	N/A	N/A	Box 15
SIT-576-15	00840253108401	Hydra Irrigation Tubing System (compatible with Olympus pumps)	N/A	N/A	N/A	N/A	Box 15

SCOPE TRANSPORT

CinchPad™ Endoscope Transport Pad



Order Number	GTIN	Description	Packaging
SPX-017-25	00816849010045	CinchPad Endoscope Transport Pad (990.6mm x 1054mm)	Box 25

SPONGES



Koala Sponge

Koala™ Cleaning Sponges



Order Number	GTIN	Description	Packaging
SSB-210-50	00816849015613	Koala Cleaning Sponge	Box 50
SES-240-50	00816849015569	Koala Cleaning Sponge plus Pure Enzymatic	Box 50

SUCTION TUBING


Suction Tubing




Order Number	GTIN	Description	Packaging
SVT-089-100	00816849016788	Suction Tubing, NS, 3/16" x 6'	Box 100
SVT-090-100	00816849016795	Suction Tubing, NS, 3/16" x 10'	Box 100
SVT-095-100	00816849016801	Suction Tubing, NS, 3/16" x 12'	Box 100
SVT-143-100	00816849016818	Suction Tubing, NS, 1/4" X 6'	Box 100
SVT-144-60	00816849016825	Suction Tubing, NS 1/4" x 10'	Box 60
SVT-145-50	00816849016832	Suction Tubing, NS, 1/4" x 12'	Box 50

► Infection Prevention


VALVES



Orca™ Sterile Single-Use Air / Water & Suction Valves



Seal™ Single-Use Biopsy Valve



Hydra Disposable Auxiliary Water Jet Connector

Orca™ Sterile Single-Use Air / Water & Suction Valves				
Order Number	GTIN	Description	Scope Compatibility	Packaging
SUV-617-50	00840253111012	Orca Sterile Single-Use Air / Water & Suction Valves	Olympus	Box 50
OrcaPod™				
Order Number	GTIN	Description	Scope Compatibility	Packaging
SUV-628-50	08714729981565	OrcaPod 3 - Orca Sterile Single-Use Air/Water, Suction, and Seal Single-Use Biopsy Valve	Olympus	Box 50
SUV-629-50	08714729981589	OrcaPod 4 - Orca Sterile Single-Use Air/Water, Suction, and Seal Single-Use Biopsy Valve with WaterJet Connector	Olympus	Box 50
Seal™ Single-Use Biopsy Valve				
Order Number	GTIN	Description	Scope Compatibility	Packaging
SBC-501-100	08714729981602	Seal Single-Use Biopsy Valve	Olympus/Fuji	Box 100


COMPLIANCE ENDOKIT™

Compliance EndoKit				
Order Number	GTIN	Description	Packaging	
CEK-721-20	00840253113337	Compliance EndoKit Kit Components: CinchPad™ Endoscope Transport Pad, Gauze, EcoBowl, Koala™ Cleaning Sponge + Pure™ Enzymatic, Biohazard Label, HedgeHog™ DualEnd Channel/Valve Brush, Endoglide+ Sterile Bacteriostatic Lubricating Gel	Box 20	

Please consult your Boston Scientific Territory Manager on options to customize your Compliance EndoKit™ to meet your individual endoscope reprocessing needs.

AEGIS ONE-HOUR INDICATOR

Aegis One-Hour Indicator				
Order Number	GTIN	Description	Packaging	
M00501991	191506007256	Aegis One-Hour Indicator	Box 100	



Pulmonary Endoscopy

EXALT™ Model B Single-Use Bronchoscope

- EXALT™ Model B Single-Use Bronchoscope 77

Alair™ Bronchial Thermoplasty System

- Alair Bronchial Thermoplasty System 78

Airway Stents

- Ultraflex™ Single-Use Partially Covered Tracheobronchial Stent System – Distal Release 80
- Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release 79
- Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release 79

Biopsy Forceps

- Radial Jaw™ 4 Single-Use Pulmonary Biopsy Forceps 81

Endobronchial Ultrasound (EBUS) Transbronchial Aspiration Needle

- Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle / Adaptor 82

Endobronchial Ultrasound Fine Needle Biopsy (FNB) Device

- Acquire™ Pulmonary Fine Needle Biopsy Device 82

Transbronchial Intra-Nodal Forceps Biopsy (IFB) Device

- CoreDx™ Pulmonary Mini Forceps 82

Transbronchial Aspiration Needles

- eXcelon™ Transbronchial Aspiration Needle 83

Cytology Brushes

- Cellebri™ Single-Use Endoscopic Cytology Brush 83

Airway Retrieval Baskets

- Zero Tip™ Airway Retrieval Basket 83

Infection Prevention

- Compliance BronchKit 83

Controlled Radial Expansion Pulmonary Balloon Dilatation Catheter

- CRE™ Pulmonary Balloon Dilatation Catheter 84


Pulmonary Guidewires

- Amplatz Super Stiff™ Guidewire 84
- Pulmonary Jagwire™ Pulmonary Guidewire 84

Inflation System

- Alliance™ II Integrated Inflation / Lithotripsy Device 84
- Alliance II Single-Use Syringe / Gauge Assembly 84

EXALT™ MODEL B SINGLE-USE BRONCHOSCOPE




EXALT Model B Single-Use Bronchoscope

Order Number	GTIN	Description	Field of View	Articulation Angle	Insertion Tube O.D. (mm)	Insertion Portion Maximum Diameter (mm)	Working Length (cm)	Minimum Working Channel Diameter (mm)	Minimum ETT Size Compatibility (mm)	Minimum DLT Size Compatibility
M00543751	00191506007492	EXALT Model B single-use bronchoscope - Large	90°	180° up, 180° down	5.8	6.3	60	2.6	7.0	N/A
M00542731	00191506007461	EXALT Model B single-use bronchoscope - Regular	90°	180° up, 180° down	5.0	5.5	60	2.0	6.0	41Fr
M00542711	00191506007423	EXALT Model B single-use bronchoscope - Slim	90°	180° up, 180° down	3.8	4.3	60	1.0	5.0	35Fr

EXALT Model B Single-Use Bronchoscope

Order Number	Description	Packaging
M00542760	EXALT Monitor (mount included)	Each
M00542780	EXALT Monitor (no mount included)	Each
M00542820	EXALT Monitor mount	Each
M00542751	EXALT Model B Single-Use Bronchoscope – Large	Box 5
M00542731	EXALT Model B Single-Use Bronchoscope – Regular	Box 5
M00542711	EXALT Model B Single-Use Bronchoscope – Slim	Box 5

ALAIR™ BRONCHIAL THERMOPLASTY SYSTEM



Alair Bronchial Thermoplasty System*

Order Number	GTIN	Model Number	Description	Packaging
M005ATS20000	08714729796695	Alair ATS 200	Alair Bronchial Thermoplasty Controller	Each
M005ATS25010	08714729802792	Alair ATS 2-5	⊗ Alair Bronchial Thermoplasty Catheter**	Each
M005ATS20110	08714729802778	Alair ATS 201	Alair RF Controller Accessory Kit	Each


*Service Agreement available with 1, 2, 3, and 4 year terms.

****Note:** Initial stocking order requires a minimum order of 6 catheters (covering the complete treatment of 2 patients).

Brief Statement of Relevant Indications for Use, Contraindications, Warnings, and Adverse Events: The Alair™ Bronchial Thermoplasty System is indicated for the treatment of severe persistent asthma in patients 18 years and older whose asthma is not well controlled with inhaled corticosteroids and long-acting beta-agonists. The Alair System is not for use in patients with an active implantable electronic device or known sensitivity to medications used in bronchoscopy. Previously treated airways of the lung should not be retreated with the Alair System. Patients should be stable and suitable to undergo bronchoscopy. The most common side effect of BT is an expected transient increase in the frequency and worsening of respiratory-related symptoms. ENDO-551804-AA

UPN	Description
SRVM005ATS1001	ENDO POS 1 Year Service Alair
SRVM005ATS1002	ENDO POS 2 Year Service Alair
SRVM005ATS1003	ENDO POS 3 Year Service Alair
SRVM005ATS1004	ENDO POS 4 Year Service Alair

AIRWAY STENTS



 Magnetic Resonance Conditional¹



Ultraflex™ Single-Use Uncovered Tracheobronchial Stent System – Distal Release*

Order Number	GTIN	Description	Expanded Stent O.D. (mm)	Expanded Stent Length (mm)	Compressed Stent Delivery System O.D. (mm)**	Tip Max O.D. (mm)	Packaging
M00576300	08714729837688	Uncovered Distal Release	10	30	7	4.1	Each
M00576310	08714729837695	Uncovered Distal Release	12	30	7	4.1	Each
M00576320	08714729837701	Uncovered Distal Release	14	30	9	4.1	Each
M00576330	08714729837718	Uncovered Distal Release	16	40	9	5.3	Each
M00576340	08714729837725	Uncovered Distal Release	16	60	9	5.3	Each
M00576350	08714729837732	Uncovered Distal Release	18	40	11	5.3	Each
M00576360	08714729837749	Uncovered Distal Release	18	60	11	5.3	Each
M00576370	08714729837756	Uncovered Distal Release	20	40	11	5.3	Each
M00576380	08714729837763	Uncovered Distal Release	20	60	11	5.3	Each

*Recommended Guidewire .035" Pulmonary Jagwire™ Pulmonary Guidewire Order #1517 and Amplatz Super Stiff™ Guidewire, Order #5009,  (See page 84).

Delivery System working length is 95cm.

**Data on File at Boston Scientific.

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹



Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release*

Order Number	GTIN	Description	Expanded Stent O.D. (mm)	Expanded Stent Length (mm)	Compressed Stent Delivery System O.D. (mm)**	Tip Max O.D. (mm)	Packaging
M00576410	08714729837770	Uncovered Proximal Release	8	20	7	4.1	Each
M00576420	08714729837787	Uncovered Proximal Release	8	40	7	4.1	Each
M00576430	08714729837794	Uncovered Proximal Release	10	20	7	4.1	Each
M00576440	08714729837800	Uncovered Proximal Release	10	40	7	4.1	Each
M00576450	08714729837817	Uncovered Proximal Release	12	20	7	4.1	Each
M00576460	08714729837824	Uncovered Proximal Release	12	40	7	4.1	Each
M00576470	08714729837831	Uncovered Proximal Release	14	20	9	4.1	Each
M00576480	08714729837848	Uncovered Proximal Release	14	40	9	4.1	Each
M00576490	08714729837855	Uncovered Proximal Release	14	60	9	4.1	Each


*Recommended Guidewire .035" Pulmonary Jagwire Pulmonary Guidewire Order #1517 and Amplatz Super Stiff™ Guidewire, Order #5009,  (See page 84).

Delivery System working length is 95cm.

**Data on File at Boston Scientific.

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

AIRWAY STENTS (Continued)



 Magnetic Resonance Conditional¹

Ultraflex™ Single-Use Partially Covered Tracheobronchial Stent System – Distal Release*




Order Number	GTIN	Description	Expanded Stent O.D. (mm)	Expanded Stent Length (mm)	Compressed Stent Delivery System O.D. (mm)**	Cover Length (mm)	Tip Max O.D. (mm)	Packaging
M00576520	08714729842323	Partially Covered Distal Release	8	40	7	25	4.1	Each
M00576530	08714729842330	Partially Covered Distal Release	10	30	7	15	4.1	Each
M00576540	08714729842347	Partially Covered Distal Release	10	40	7	25	4.1	Each
M00576550	08714729842354	Partially Covered Distal Release	12	30	7	15	4.1	Each
M00576560	08714729842361	Partially Covered Distal Release	12	40	7	25	4.1	Each
M00576570	08714729842378	Partially Covered Distal Release	14	30	9	15	4.1	Each
M00576580	08714729842385	Partially Covered Distal Release	14	40	9	25	4.1	Each
M00576590	08714729842392	Partially Covered Distal Release	14	60	9	45	4.1	Each
M00576600	08714729842408	Partially Covered Distal Release	14	80	9	65	4.1	Each
M00576610	08714729842415	Partially Covered Distal Release	16	40	9	25	5.3	Each
M00576620	08714729842422	Partially Covered Distal Release	16	60	9	45	5.3	Each
M00576630	08714729842439	Partially Covered Distal Release	16	80	9	65	5.3	Each
M00576640	08714729842446	Partially Covered Distal Release	18	40	11	25	5.3	Each
M00576650	08714729842453	Partially Covered Distal Release	18	60	11	45	5.3	Each
M00576660	08714729842460	Partially Covered Distal Release	18	80	11	65	5.3	Each
M00576670	08714729842477	Partially Covered Distal Release	20	40	11	25	5.3	Each
M00576680	08714729842484	Partially Covered Distal Release	20	60	11	45	5.3	Each
M00576690	08714729842491	Partially Covered Distal Release	20	80	11	65	5.3	Each

*Recommended Guidewire .035" Pulmonary Jagwire™ Pulmonary Guidewire Order #1517 and Amplatz Super Stiff™ Guidewire, Order #5009,  (See page 84).
Delivery System working length is 95cm.

**Data on File at Boston Scientific.


¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

BIOPSY FORCEPS



Radial Jaw™ 4 Single-Use Pulmonary Biopsy Forceps 							
Order Number	GTIN	Description	Differentiator	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Packaging
M00515181	08714729792864	Single-Use Biopsy Forceps	Pulmonary Standard Capacity	1.8	100	2.0	Box 5
M00515182	08714729799337	Single-Use Biopsy Forceps	Pulmonary Standard Capacity	1.8	100	2.0	Box 20
M00515191	08714729792871	Single-Use Biopsy Forceps	Pulmonary Standard Capacity w/ Needle	1.8	100	2.0	Box 5
M00515192	08714729818083	Single-Use Biopsy Forceps	Pulmonary Standard Capacity w/ Needle	1.8	100	2.0	Box 20
M00515201	08714729793281	Single-Use Biopsy Forceps	Pulmonary Large Capacity	2.2	100	2.8	Box 5
M00515202	08714729818076	Single-Use Biopsy Forceps	Pulmonary Large Capacity	2.2	100	2.8	Box 20

EXPECT™ PULMONARY ENDOBRONCHIAL ULTRASOUND TRANSBRONCHIAL ASPIRATION NEEDLE




Expect Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle



Order Number	GTIN	Needle Size	Minimum Working Channel (mm)	Sheath O.D. (mm)	Packaging
M00558220	08714729861409	22Ga – Olympus scope compatible	2.0	1.6	Each
M00558250	08714729861416	25Ga – Olympus scope compatible	2.0	1.4	Each
M00558731	08714729893189	Expect Pulmonary Needle Adaptor – Olympus scope compatible	N/A	N/A	Box 10

ACQUIRE™ PULMONARY ENDOBRONCHIAL ULTRASOUND FINE NEEDLE BIOPSY DEVICE




Acquire Pulmonary Endobronchial Ultrasound Fine Needle Biopsy Device



Order Number	GTIN	Needle Size	Minimum Working Channel (mm)	Sheath O.D. (mm)	Packaging
M00552350	08714729986225	22Ga – Olympus Scope Compatible	2.0	1.6	Each
M00552400	08714729986249	25Ga – Olympus Scope Compatible	2.0	1.4	Each
M00558731	08714729893189	Expect Pulmonary Needle Adaptor – Olympus scope compatible	N/A	N/A	Box 10

COREDx™ PULMONARY MINI-FORCEPS FOR ENDOBRONCHIAL ULTRASOUND




CoreDx Pulmonary Mini-Forceps for Endobronchial Ultrasound



Order Number	GTIN	Description	Minimum Working Channel (mm)	Catheter Length (cm)	Packaging
M00515220	08714729971054	CoreDx Pulmonary Mini-Forceps	1.2	130	Each

TRANSBRONCHIAL ASPIRATION NEEDLES



eXcelon™ Transbronchial Aspiration Needle*



Order Number	GTIN	Description	Needle Gauge	Needle Length (mm)	Catheter Length (cm)	Sheath O.D. (mm)	Packaging
M00564101	08714729456452	eXcelon Transbronchial Aspiration Needle	19	15	130	1.8	Box 5
M00564111	08714729456469	eXcelon Transbronchial Aspiration Needle	20	15	130	1.8	Box 5
M00564121	08714729456476	eXcelon Transbronchial Aspiration Needle	21	15	130	1.8	Box 5

*Packaging also includes a 20ml syringe.


CYTOLOGY BRUSHES

Cellebri™ Single-Use Endoscopic Cytology Brush



Order Number	GTIN	Required Working Channel (mm)	Bristle O.D. (mm)	Sheath Length (cm)	Packaging
M00516001	08714729138983	2.0	1.0	140	Box 10
M00516011	08714729138990	2.0	1.5	140	Box 10
M00516071	08714729139089	2.0	1.9	100	Box 10
M00516151	08714729014003	2.0	1.9	150	Box 10

AIRWAY RETRIEVAL BASKETS



Zero Tip™ Airway Retrieval Basket



Order Number	GTIN	Description	Sheath O.D. (mm)	Working Length (cm)	Basket Opening (mm)	Packaging
M00513200	08714729414995	Zero Tip Airway Retrieval Basket	0.8	120	12	Each
M00513210	08714729415008	Zero Tip Airway Retrieval Basket	1.0	120	16	Each

INFECTION PREVENTION

Compliance BronchKit

Order Number	GTIN	Description	Packaging
CBK-881-20	00840253116031	Compliance BronchKit Kit Components: CinchPad™ Endoscope Transport Pad, Gauze, EcoBowl, Koala Cleaning Sponge + Pure Enzymatic, Biohazard Label, HedgeHog™ Double Bundle Cleaning Brush, 20ml Syringe	Box 20
CBK-882-20	00840253120281	Compliance BronchKit Kit Components: CinchPad™ Endoscope Transport Pad, Gauze, EcoBowl, Koala Cleaning Sponge + Pure Enzymatic, Biohazard Label, HedgeHog™ Double Bundle Cleaning Brush, 20ml Syringe, EndoArmor+ Gown, Suction Tubing	Box 20

CONTROLLED RADIAL EXPANSION PULMONARY BALLOON DILATATION CATHETER



CRE™ Pulmonary Balloon Dilatation Catheter*



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Catheter Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)		
M00550300	08714729456186	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	110	Each
M00550310	08714729456193	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	110	Each
M00550320	08714729456209	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	110	Each
M00550330	08714729456216	8-9-10	24-27-30	3.0	3-5.5-9	304-557-912	110	Each
M00550340	08714729456223	10-11-12	30-33-36	3.0	3-5-8	304-507-811	110	Each
M00550350	08714729456230	12-13.5-15	36-40.5-45	3.0	3-4.5-8	304-456-811	110	Each

*Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060 (see below).

PULMONARY GUIDEWIRES

Pulmonary Jagwire™ Pulmonary Guidewire




Order Number	GTIN	Description	O.D. (in)	Length (cm)	Packaging
M00515171	08714729455813	Pulmonary Jagwire Pulmonary Guidewire	.035	180	Box 2

Amplatz Super Stiff™ Guidewire



Order Number	GTIN	Description	O.D. (in)	Length (cm)	Packaging
M00550090	08714729301929	Amplatz Super Stiff Guidewire	.038	260	Each

INFLATION SYSTEM



Alliance™ II Integrated Inflation / Lithotripsy Device and Alliance II Single-Use Syringe / Gauge Assembly

Order Number	GTIN	Description	Packaging
M00550620	08714729283881	Alliance II Integrated Inflation / Lithotripsy Device	Each
M00550601	08714729129332	Alliance II Single-Use Syringe / Gauge Assembly	Box 5

Notes and Future Product Updates

Notes and Future Product Updates



Section Table of Contents



Prev Page

Next Page



Home





Solutions

- **GuidePoint Reimbursement Resources and Solutions**.....85
- **ADVANTICS™ Innovative Healthcare Solutions** 86
- **GI Specialized Pathology**..... 86
- **Scope Repair**..... 86
- **GI Lab Optimization**..... 86

GUIDEPOINT REIMBURSEMENT RESOURCES AND SOLUTIONS

The Endoscopy Health Economics and Market Access (HEMA) Team offers a suite of reimbursement services and support for healthcare professionals seeking a better understanding of reimbursement for Boston Scientific's endoscopy and pulmonary products and procedures. We can help you develop a base level of knowledge in coding, coverage and payment for these procedures and offer educational presentations and webinars on changes in healthcare policies and Medicare reimbursement and changes to Medicare reimbursement impacting providers and facilities.

About Our Team

Boston Scientific Endoscopy's HEMA Team is comprised of highly skilled professionals specializing in reimbursement, health economics and health policy issues. The team has a broad background with experience working for public and private payers, medical associations, healthcare consulting agencies and academic research organizations. We use clinical and economic outcomes data and collaborate with outside experts to promote appropriate coding, coverage and reimbursement for Boston Scientific endoscopy and pulmonary products and procedures.

Contact Information

For more information about Boston Scientific Endoscopy's GuidePoint Reimbursement Resources and Solutions, please contact our **Health Economics and Reimbursement department at 800-876-9960 x54510 or ENDOREIMBURSEMENT@bsci.com**.

www.bostonscientific.com/reimbursement

GUIDEPOINT
Reimbursement Resources

Health economic and reimbursement information provided by Boston Scientific Corporation is gathered from third-party sources and is subject to change without notice as a result of complex and frequently changing laws, regulations, rules and policies. This information is presented for illustrative purposes only and does not constitute reimbursement or legal advice. Boston Scientific encourages providers to submit accurate and appropriate claims for services. It is always the provider's responsibility to determine medical necessity, the proper site for delivery of any services and to submit appropriate codes, charges, and modifiers for services that are rendered. Boston Scientific recommends that you consult with your payers, reimbursement specialists and/or legal counsel regarding coding, coverage and reimbursement matters. Boston Scientific does not promote the use of its products outside their FDA-approved label.

► Solutions

Boston Scientific's Endoscopy business develops less invasive medical devices for treating gastrointestinal diseases. We are committed to providing clinically proven technologies, services and customer support to help healthcare providers improve patient outcomes while reducing the total cost of patient care.

ADVANTICS™ INNOVATIVE HEALTHCARE SOLUTIONS

ADVANTICS™ Innovative Healthcare Solutions is a new way to collaborate to meet your specific challenges in the changing healthcare landscape. ADVANTICS recognizes that every provider is unique and that only by working together can we find truly effective solutions.



GI SPECIALIZED PATHOLOGY

Diagnostics for GI Diseases. Getting an accurate and timely diagnosis is critical to patient care. GI Specialized Pathology services include GI subspecialty trained pathologists, GI trained histotechnicians to help ensure proper slide orientation, and customer service teams dedicated to improving efficiencies and the diagnosis process for physicians and their patients. Today, of the top diagnostic labs in the U.S., ours is the only pathology lab focused solely on GI diseases.



SCOPE REPAIR

Capital Equipment Management. As an alternative to original equipment manufacturer (OEM) and costly service contracts, scope repair services and certified pre-owned equipment are offered through a partnership with Total Scope, Inc. Established in 1993, Total Scope is focused on providing excellent customer service, including scope maintenance education and training for its customers. The company is certified to the quality standard ISO 13485, for the service, repair and maintenance of medical instrumentation. Total Scope is a woman-owned business.

GI LAB OPTIMIZATION

Endoscopy Lab and Continuous Improvement. GI Lab Optimization services are designed to address specific customer issues related to quality, cost, safety, and staff and patient satisfaction. A mix of Lean Six Sigma, project management and Human Centered Design expertise can be used to tailor solutions to your needs for endoscopy lab optimization.

To learn more visit, [bostonscientific.com/endoscopysolutions](https://www.bostonscientific.com/endoscopysolutions)

Notes and Future Product Updates

Notes and Future Product Updates



Section Table of Contents



Prev Page

Next Page



Home





Product Code Index

▶ Begins on Page87

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
5100.....	Biopolar Adaptor Cable, 28mm spacing	31	M00501201	EndoVive Low Profile Replacement Button Accessories with ENFit – 24Fr Decompression Tube	59, 61
5420.....	Biopolar Adaptor Cable, 4mm Pin*	31	M00501211	EndoVive Low Profile Replacement Button Accessories with ENFit – 24Fr Decompression Tube	59, 61
5700.....	Biopolar Adaptor Cable, 22mm spacing	31	M00501261	EndoVive Initial Placement PEG Accessories – Y-Ports with ENFit.....	58
CBK-881-20	Compliance BronchKit	83	M00501271	EndoVive Initial Placement PEG Accessories – Y-Ports with ENFit.....	58
CBK-882-20	Compliance BronchKit	83	M00501972	EndoArmor+ Non-Sterile Surgical Gown.....	83
CEK-721-20.....	Compliance EndoKit™	76	M00501991	Aegis One-Hour Indicator	85
DGA-275-5.....	Rescue™ Alligator Long Grasping Forceps	72	M00509001	EndoVive Safety PEG Kit, Pull with ENFit.....	58
DGN-538-5	RescueNet™ Retrieval Device	72	M00509011	EndoVive Safety PEG Kit, Push with ENFit	58
DGR-276-5.....	Rescue Rat Tooth Grasping Forceps	72	M00509021	EndoVive Safety PEG Kit, Pull with ENFit.....	58
H965180010	WALLSTENT™ Super Stiff Guidewire.....	3	M00509031	EndoVive Safety PEG Kit, Push with ENFit	58
H965180011	WALLSTENT Super Stiff Guidewire.....	3, 45	M00509041	EndoVive Standard PEG Kit, Pull with ENFit.....	58
H965430100	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	16	M00509051	EndoVive Standard PEG Kit, Push with ENFit.....	58
H965430200	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	16	M00509061	EndoVive Standard PEG Kit, Pull with ENFit.....	58
H965430300	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	16	M00509071	EndoVive Standard PEG Kit, Push with ENFit.....	58
H965430400	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	16	M00509241	ENFit Feeding Port Transition Connector.....	61
H965430500	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	16	M00509580	EndoVive One-Step Button™ Low Profile PEG Kit, Pull with ENFit.....	59
H965430600	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	16	M00509590	EndoVive One-Step Button Low Profile PEG Kit, Pull with ENFit.....	59
H965430700	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	16	M00509600	EndoVive One-Step Button Low Profile PEG Kit, Pull with ENFit.....	59
H965430800	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	16	M00509610	EndoVive One-Step Button Low Profile PEG Kit, Pull with ENFit.....	59
H965430900	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	16	M00509620	EndoVive One-Step Button Low Profile PEG Kit, Push with ENFit	59
H965431000	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	16	M00509630	EndoVive One-Step Button Low Profile PEG Kit, Push with ENFit	59
H965431100	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	16	M00509680	EndoVive One-Step Button Low Profile PEG Kit, Push with ENFit	59
H965431200	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	16	M00509690	EndoVive One-Step Button Low Profile PEG Kit, Pull with ENFit.....	59
M00500070	Habib™ EndoHPB Bipolar RF Catheter.....	31	M00509700	EndoVive One-Step Button Low Profile PEG Kit, Pull with ENFit.....	59
M00501061	EndoVive™ Low Profile Replacement Button Accessories with ENFit™ – 18Fr Right Angle Feeding Set, Continuous Adapter	59, 61	M00509710	EndoVive One-Step Button Low Profile PEG Kit, Pull with ENFit.....	59
M00501071	EndoVive Low Profile Replacement Button Accessories with ENFit – 24Fr Right Angle Feeding Set, Continuous Adapter.....	59, 61	M00509720	EndoVive One-Step Button Low Profile PEG Kit, Pull with ENFit.....	59
M00501091	EndoVive Low Profile Replacement Button Accessories with ENFit – 18Fr Bolus Feeding Set, Straight Adapter....	59, 61	M00509730	EndoVive One-Step Button Low Profile PEG Kit, Pull with ENFit.....	59
M00501101	EndoVive Low Profile Replacement Button Accessories with ENFit – 24Fr Bolus Feeding Set, Straight Adapter	59, 61	M00509740	EndoVive One-Step Button Low Profile PEG Kit, Push with ENFit	59
M00501121	EndoVive Low Profile Replacement Button Accessories with ENFit – 18Fr Decompression Tube.....	59, 61	M00509750	EndoVive One-Step Button Low Profile PEG Kit, Push with ENFit	59
M00501131	EndoVive Low Profile Replacement Button Accessories with ENFit – 18Fr Decompression Tube.....	59, 61	M00509830	EndoVive Low Profile Replacement Button with ENFit	60
M00501141	EndoVive Low Profile Replacement Button Accessories with ENFit – 18Fr Decompression Tube.....	59, 61	M00509840	EndoVive Low Profile Replacement Button with ENFit.....	60
M00501151	EndoVive Low Profile Replacement Button Accessories with ENFit – 18Fr Decompression Tube.....	59, 61			
M00501161	EndoVive Low Profile Replacement Button Accessories with ENFit – 18Fr Decompression Tube.....	59, 61			
M00501181	EndoVive Low Profile Replacement Button Accessories with ENFit – 24Fr Decompression Tube.....	59, 61			
M00501191	EndoVive Low Profile Replacement Button Accessories with ENFit – 24Fr Decompression Tube.....	59, 61			

*For use with the Habib EndoHPB Bipolar Radiofrequency Catheter

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00509850	EndoVive Low Profile Replacement Button with ENFit.....	60	M00513321	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps	56
M00509860	EndoVive™ Low Profile Replacement Button with ENFit™	60	M00513322	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps	56
M00509870	EndoVive Low Profile Replacement Button with ENFit.....	60	M00513323	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps	56
M00509410	EndoVive™ Low Profile Replacement Button with ENFit.....	60	M00513331	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle.....	56
M00509420	EndoVive Low Profile Replacement Button with ENFit.....	60	M00513332	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle.....	56
M00509430	EndoVive Low Profile Replacement Button with ENFit.....	60	M00513333	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle.....	56
M00509440	EndoVive Low Profile Replacement Button with ENFit.....	60	M00513343	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps	56
M00509881	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight.....	60	M00513353	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps w/ Needle..	56
M00509891	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight.....	60	M00513361	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps	56
M00509901	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight.....	60	M00513362	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps	56
M00509911	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight.....	60	M00513363	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps	56
M00509921	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight.....	60	M00513371	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps w/ Needle..	56
M00509931	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight.....	60	M00513372	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps w/ Needle..	56
M00509941	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight.....	60	M00513373	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps w/ Needle..	56
M00509961	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	60	M00513381	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps..	56
M00509971	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	60	M00513382	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps..	56
M00509981	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	60	M00513383	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps..	56
M00509991	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	60	M00513391	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle.....	56
M00510001	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	60	M00513392	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle.....	56
M00510011	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	60	M00513393	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle.....	56
M00510021	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	60	M00513401	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps..	56
M00510101	Multibite™ Single-Use Multiple Sample Biopsy Forceps.....	57	M00513402	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps....	56
M00510121	Multibite Single-Use Multiple Sample Biopsy Forceps.....	57	M00513403	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps..	56
M00510860	Trapezoid™ RX Wireguided Retrieval Basket	29	M00513411	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle.....	56
M00510870	Trapezoid RX Wireguided Retrieval Basket	29	M00513412	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle.....	56
M00510880	Trapezoid RX Wireguided Retrieval Basket	29	M00513413	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle.....	56
M00510890	Trapezoid RX Wireguided Retrieval Basket	29	M00513443	Radial Jaw 4 Single-Use Pediatric Biopsy Forceps	57
M00513200	Zero Tip™ Airway Retrieval Basket	82	M00513453	Radial Jaw 4 Single-Use Pediatric Biopsy Forceps w/ Needle....	57
M00513210	Zero Tip Airway Retrieval Basket	82	M00513700	Ultraflex™ Esophageal NG Uncovered Stent System – Distal Release.....	44
M00513301	Radial Jaw™ 4 Single-Use Large Capacity Biopsy Forceps	56	M00513710	Ultraflex Esophageal NG Uncovered Stent System – Distal Release.....	44
M00513302	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps	56	M00513720	Ultraflex Esophageal NG Uncovered Stent System – Distal Release.....	44
M00513303	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps	56	M00513730	Ultraflex Esophageal NG Covered Stent System – Distal Release.....	43
M00513311	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle.....	56	M00513740	Ultraflex Esophageal NG Covered Stent System – Distal Release.....	43
M00513312	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle.....	56	M00513750	Ultraflex Esophageal NG Covered Stent System – Distal Release.....	43
M00513313	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle.....	56	M00513800	Ultraflex Esophageal NG Uncovered Stent System – Proximal Release	44

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00513810	Ultraflex Esophageal NG Uncovered Stent System – Proximal Release	44	M00517200	AGILE™ Esophageal Stent System	41
M00513830	Ultraflex™ Esophageal NG Uncovered Stent System – Proximal Release	44	M00517210	AGILE Esophageal Stent System	41
M00513840	Ultraflex Esophageal NG Covered Stent System – Proximal Release	43	M00517220	AGILE Esophageal Stent System	41
M00513850	Ultraflex Esophageal NG Covered Stent System – Proximal Release	43	M00517230	AGILE Esophageal Stent System	41
M00513860	Ultraflex Esophageal NG Covered Stent System – Proximal Release	43	M00517240	AGILE Esophageal Stent System	41
M00514200	Ultraflex Esophageal NG Covered Stent System – Distal Release	43	M00517250	AGILE Esophageal Stent System	41
M00514210	Ultraflex Esophageal NG Covered Stent System – Distal Release	43	M00517260	AGILE Esophageal Stent System	41
M00514240	Ultraflex Esophageal NG Covered Stent System – Proximal Release	44	M00517270	AGILE Esophageal Stent System	41
M00514250	Ultraflex Esophageal NG Covered Stent System – Proximal Release	44	M00517400	AGILE Esophageal Stent System	41
M00515031	Radial Jaw™ 4 Single-Use Hot Biopsy Forceps	57	M00517410	AGILE Esophageal Stent System	41
M00515032	Radial Jaw 4 Single-Use Hot Biopsy Forceps	57	M00517420	AGILE Esophageal Stent System	41
M00515033	Radial Jaw 4 Single-Use Hot Biopsy Forceps	57	M00517430	AGILE Esophageal Stent System	41
M00515171	Pulmonary Jagwire™ Pulmonary Guidewire	84	M00517440	AGILE Esophageal Stent System	41
M00515181	Radial Jaw 4 Single-Use Pulmonary Forceps – Standard Capacity	81	M00517450	AGILE Esophageal Stent System	41
M00515182	Radial Jaw 4 Single-Use Pulmonary Forceps – Standard Capacity	81	M00517460	AGILE Esophageal Stent System	41
M00515191	Radial Jaw 4 Single-Use Pulmonary Forceps – Standard Capacity with Needle	81	M00517470	AGILE Esophageal Stent System	41
M00515192	Radial Jaw 4 Single-Use Pulmonary Forceps – Standard Capacity with Needle	81	M00517710	AGILE ESO Esophageal Stent System	41
M00515201	Radial Jaw 4 Single-Use Pulmonary Forceps – Large Capacity	81	M00517720	AGILE ESO Esophageal Stent System	41
M00515202	Radial Jaw 4 Single-Use Pulmonary Forceps – Large Capacity	81	M00517730	AGILE ESO Esophageal Stent System	41
M00515220	CoreDx™ Pulmonary Mini-Forceps for Endobronchial Ultrasound	82	M00517740	AGILE ESO Esophageal Stent System	41
M00516001	Cellebrity™ Single-Use Endoscopic Cytology Brush	83	M00517750	AGILE ESO Esophageal Stent System	41
M00516011	Cellebrity Single-Use Endoscopic Cytology Brush	83	M00517760	AGILE ESO Esophageal Stent System	41
M00516071	Cellebrity Single-Use Endoscopic Cytology Brush	83	M00517770	AGILE ESO Esophageal Stent System	41
M00516151	Cellebrity Single-Use Endoscopic Cytology Brush	83	M00517780	AGILE ESO Esophageal Stent System	41
M00516700	WallFlex™ Fully Covered Esophageal Stent System	42	M00518101	Interject™ Clear Single-Use Injection Therapy Needle Catheters	47
M00516710	WallFlex Fully Covered Esophageal Stent System	42	M00518111	Interject Clear Single-Use Injection Therapy Needle Catheters	47
M00516720	WallFlex Fully Covered Esophageal Stent System	42	M00518151	Interject Contrast Single-Use Injection Therapy Needle Catheters	47
M00516730	WallFlex Fully Covered Esophageal Stent System	42	M00518161	Interject Contrast Single-Use Injection Therapy Needle Catheters	47
M00516740	WallFlex Fully Covered Esophageal Stent System	42	M00518251	Interject Contrast Single-Use Injection Therapy Needle Catheters	47
M00516750	WallFlex Fully Covered Esophageal Stent System	42	M00518261	Interject Contrast Single-Use Injection Therapy Needle Catheters	47
M00516900	WallFlex Partially Covered Esophageal Stent System	42	M00518301	Interject Clear Single-Use Injection Therapy Needle Catheters	47
M00516910	WallFlex Partially Covered Esophageal Stent System	42	M00518311	Interject Clear Single-Use Injection Therapy Needle Catheters	47
M00516920	WallFlex Partially Covered Esophageal Stent System	42	M00518351	Interject Contrast Single-Use Injection Therapy Needle Catheters	47
M00516930	WallFlex Partially Covered Esophageal Stent System	42	M00518361	Interject Contrast Single-Use Injection Therapy Needle Catheters	47
M00516940	WallFlex Partially Covered Esophageal Stent System	42	M00519201	ORISE™ Gel Syringe Twin Pack (two 10mL syringes per tray)	71
M00516950	WallFlex Partially Covered Esophageal Stent System	42	M00519210	ORISE Gel Syringe Twin Pack Kit with Interject Needle 23ga Clear (two 10mL syringes per kit)	71
			M00519211	ORISE Gel Syringe Twin Pack Kit with Interject Needle 23ga Clear (two 10mL syringes per kit)	71
			M00519221	ORISE Gel Syringe Single Pack (one 10mL syringe per tray)	71

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00519231	ORISE Gel Syringe Kit Single Pack with Interject Injection Therapy Needle Catheter 23ga Clear (one 10ml syringe per kit)	71	M00532300	Advanix Double Pigtail - Single Stent	19
M00519351	ORISE™ ProKnife	71	M00532310	Advanix Double Pigtail - Single Stent	19
M00519361	ORISE ProKnife	71	M00532810	MicroKnife™ XL Triple-Lumen Needle Knife	8
M00519371	ORISE ProKnife	71	M00532870	Advanix Duodenal Bend - Single Stent	18
M00519380	ORISE ProKnife	71	M00532880	Advanix Duodenal Bend - Single Stent	18
M00519390	ORISE ProKnife	71	M00532890	Advanix Duodenal Bend - Single Stent	18
M00519400	ORISE ProKnife	71	M00532900	Advanix Duodenal Bend - Single Stent	18
M00521230	Resolution 360™ Clip	45	M00532910	Advanix Duodenal Bend - Single Stent	18
M00521231	Resolution 360 Clip	45	M00532920	Advanix Duodenal Bend - Single Stent	18
M00521232	Resolution 360 Clip	45	M00532930	Advanix Duodenal Bend - Single Stent	18
M00521233	Resolution 360 Clip	45	M00532940	Advanix Duodenal Bend - Single Stent	18
M00521240	Resolution 360 Clip	45	M00532950	Advanix Duodenal Bend - Single Stent	18
M00521241	Resolution 360 Clip	45	M00532960	Advanix Duodenal Bend - Single Stent	18
M00521242	Resolution 360 Clip	45	M00532970	Advanix Duodenal Bend - Single Stent	18
M00521400	Resolution 360 ULTRA Clip	45	M00532980	Advanix Duodenal Bend - Single Stent	18
M00521401	Resolution 360 ULTRA Clip	45	M00532990	Advanix Duodenal Bend - Single Stent	18
M00521402	Resolution 360 ULTRA Clip	45	M00533000	Advanix Duodenal Bend - Single Stent	18
M00521403	Resolution 360 ULTRA Clip	45	M00533010	Advanix Duodenal Bend - Single Stent	18
M00522600	Resolution™ Clip	45	M00533020	Advanix Duodenal Bend - Single Stent	18
M00522601	Resolution Clip	45	M00533030	Advanix Duodenal Bend - Single Stent	18
M00522602	Resolution Clip	45	M00533040	Advanix Duodenal Bend - Single Stent	18
M00522610	Resolution Clip	45	M00533230	Advanix Center Bend - Single Stent	18
M00522611	Resolution Clip	45	M00533240	Advanix Center Bend - Single Stent	18
M00522612	Resolution Clip	45	M00533250	Advanix Center Bend - Single Stent	18
M00527451	Rescue™ Alligator Grasping Forceps	72	M00533260	Advanix Center Bend - Single Stent	18
M00527751	Rescue Rat Tooth PED Grasping Forceps	72	M00533270	Advanix Center Bend - Single Stent	18
M00530800	Ultratome™ Double-Lumen Sphincterotome – Short Nose	10	M00533280	Advanix Center Bend - Single Stent	18
M00530810	Ultratome Double-Lumen Sphincterotome – Long Nose	10	M00533290	Advanix Center Bend - Single Stent	18
M00530820	Ultratome Double-Lumen Sphincterotome – Short Nose	10	M00533300	Advanix Center Bend - Single Stent	18
M00530830	Ultratome Double-Lumen Sphincterotome – Long Nose	10	M00533310	Advanix Center Bend - Single Stent	18
M00530850	Contour™ ERCP Cannulas – Standard	11	M00533320	Advanix Center Bend - Single Stent	18
M00530860	Contour ERCP Cannulas – Tapered	11	M00533330	Advanix Center Bend - Single Stent	18
M00530870	Contour ERCP Cannulas – Ultra Tapered	11	M00533340	Advanix Center Bend - Single Stent	18
M00530880	Contour ERCP Cannulas – 5-4-3	11	M00533350	Advanix Center Bend - Single Stent	18
M00530890	Contour ERCP Cannulas – Ball Tip	11	M00533360	Advanix Center Bend - Single Stent	18
M00532160	Advanix™ Double Pigtail - Single Stent	19	M00533370	Advanix Center Bend - Single Stent	18
M00532180	Advanix Double Pigtail - Single Stent	19	M00533380	Advanix Center Bend - Single Stent	18
M00532190	Advanix Double Pigtail - Single Stent	19	M00533390	Advanix Center Bend - Single Stent	18
M00532210	Advanix Double Pigtail - Single Stent	19	M00533400	Advanix Center Bend - Single Stent	18
M00532220	Advanix Double Pigtail - Single Stent	19	M00533560	NaviFlex™ RX Delivery System	19
M00532230	Advanix Double Pigtail - Single Stent	19	M00533570	NaviFlex RX Delivery System	19
M00532240	Advanix Double Pigtail - Single Stent	19	M00533580	NaviFlex RX Delivery System	19
M00532260	Advanix Double Pigtail - Single Stent	19	M00534200	Advanix Duodenal Bend Preloaded RX Stent Systems	17
M00532270	Advanix Double Pigtail - Single Stent	19	M00534210	Advanix Duodenal Bend Preloaded RX Stent Systems	17
M00532290	Advanix Double Pigtail - Single Stent	19	M00534220	Advanix Duodenal Bend Preloaded RX Stent Systems	17
			M00534230	Advanix Duodenal Bend Preloaded RX Stent Systems	17
			M00534240	Advanix Duodenal Bend Preloaded RX Stent Systems	17

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00534250	Advanix Duodenal Bend Preloaded RX Stent Systems.....	17	M00535920	Ultratome XL Triple-Lumen Sphincterotome – Short Nose	9
M00534260	Advanix Duodenal Bend Preloaded RX Stent Systems.....	17	M00535930	Ultratome XL Triple-Lumen Sphincterotome – Long Nose	9
M00534270	Advanix Duodenal Bend Preloaded RX Stent Systems.....	17	M00536000	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534280	Advanix™ Duodenal Bend Preloaded RX Stent Systems.....	17	M00536010	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534290	Advanix Duodenal Bend Preloaded RX Stent Systems.....	17	M00536020	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534300	Advanix Duodenal Bend Preloaded RX Stent Systems.....	17	M00536030	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534310	Advanix Duodenal Bend Preloaded RX Stent Systems.....	17	M00536040	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534320	Advanix Duodenal Bend Preloaded RX Stent Systems.....	17	M00536050	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534330	Advanix Duodenal Bend Preloaded RX Stent Systems.....	17	M00536060	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534340	Advanix Duodenal Bend Preloaded RX Stent Systems.....	17	M00536070	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534350	Advanix Duodenal Bend Preloaded RX Stent Systems.....	17	M00536080	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534360	Advanix Duodenal Bend Preloaded RX Stent Systems.....	17	M00536090	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534370	Advanix Duodenal Bend Preloaded RX Stent Systems.....	17	M00536100	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534560	Advanix Center Bend Preloaded RX Stent Systems	17	M00536110	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534570	Advanix Center Bend Preloaded RX Stent Systems	17	M00536120	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534580	Advanix Center Bend Preloaded RX Stent Systems	17	M00536130	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534590	Advanix Center Bend Preloaded RX Stent Systems	17	M00536140	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534600	Advanix Center Bend Preloaded RX Stent Systems	17	M00536150	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534610	Advanix Center Bend Preloaded RX Stent Systems	17	M00536160	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534620	Advanix Center Bend Preloaded RX Stent Systems	17	M00536170	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534630	Advanix Center Bend Preloaded RX Stent Systems	17	M00536180	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534640	Advanix Center Bend Preloaded RX Stent Systems	17	M00536190	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534650	Advanix Center Bend Preloaded RX Stent Systems	17	M00536200	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534660	Advanix Center Bend Preloaded RX Stent Systems	17	M00536210	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534670	Advanix Center Bend Preloaded RX Stent Systems	17	M00536220	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534680	Advanix Center Bend Preloaded RX Stent Systems	17	M00536230	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534690	Advanix Center Bend Preloaded RX Stent Systems	17	M00536240	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534700	Advanix Center Bend Preloaded RX Stent Systems	17	M00536250	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534710	Advanix Center Bend Preloaded RX Stent Systems	17	M00536260	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534720	Advanix Center Bend Preloaded RX Stent Systems	17	M00536270	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00534730	Advanix Center Bend Preloaded RX Stent Systems	17	M00536280	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00535000	NaviFlex™ Pancreatic Stent RX Pusher	27	M00536290	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00535010	NaviFlex Pancreatic Stent RX Pusher	27	M00536300	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00535020	NaviFlex Pancreatic Stent Pusher.....	27	M00536310	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00535030	NaviFlex Pancreatic Stent Pusher.....	27	M00536320	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00535060	NaviFlex RX Pancreatic Stent Delivery System	27	M00536330	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00535080	NaviFlex RX Pancreatic Stent Delivery System	27	M00536340	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00535110	Stonetome™ Stone Removal Device	9	M00536350	Advanix Pancreatic Stent – Straight, No Leading Barb.....	23
M00535130	Stonetome Stone Removal Device	9	M00536360	Advanix Pancreatic Stent – Straight, Leading Barb	24
M00535150	Stonetome Stone Removal Device	9	M00536370	Advanix Pancreatic Stent – Straight, Leading Barb	24
M00535170	Stonetome Stone Removal Device	9	M00536380	Advanix Pancreatic Stent – Straight, Leading Barb	24
M00535190	Stonetome Stone Removal Device	9	M00536390	Advanix Pancreatic Stent – Straight, Leading Barb	24
M00535210	Stonetome Stone Removal Device	9	M00536400	Advanix Pancreatic Stent – Straight, Leading Barb	24
M00535700	Tandem™ XL Triple-Lumen ERCP Cannula – Tapered	11	M00536410	Advanix Pancreatic Stent – Straight, Leading Barb	24
M00535900	Ultratome™ XL Triple-Lumen Sphincterotome – Short Nose.....	9	M00536420	Advanix Pancreatic Stent – Straight, Leading Barb	24
M00535910	Ultratome XL Triple-Lumen Sphincterotome – Long Nose	9	M00536430	Advanix Pancreatic Stent – Straight, Leading Barb	24
			M00536440	Advanix Pancreatic Stent – Straight, Leading Barb	24

► Product Code Index

[illegible]

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00537390	Advanix Pancreatic Stent Kit – Pigtail, No Leading Barb.....	22	M00542420	EXALT™ Model D Single-use Duodenoscope	35
M00537400	Advanix™ Pancreatic Stent Kit – Pigtail, No Leading Barb.....	22	M00542421	EXALT Model D Single-use Duodenoscope.....	35
M00537410	Advanix Pancreatic Stent Kit – Pigtail, No Leading Barb.....	22	M00542430	EXALT Controller.....	35
M00537420	Advanix™ Pancreatic Stent Kit – Pigtail, No Leading Barb.....	22	M00542711	EXALT Model B Single-Use Bronchoscope - Slim.....	77
M00537430	Advanix Pancreatic Stent Kit – Pigtail, Leading Barb	22	M00542731	EXALT Model B Single-Use Bronchoscope - Regular	77
M00537440	Advanix Pancreatic Stent Kit – Pigtail, Leading Barb	22	M00542760	EXALT Monitor (mount included)	77
M00537450	Advanix Pancreatic Stent Kit – Pigtail, Leading Barb	22	M00542780	EXALT Monitor (no mount included)	77
M00537460	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	22	M00542820	EXALT Monitor Mount	77
M00537470	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	22	M00542950	SpyGlass™ Discover Controller	33
M00537480	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	22	M00542961	Super Sheath CBDE for use with the SpyGlass.....	33
M00537490	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	22	M00542971	SpyGlass Discover Jagwire™ Guidewire.....	33
M00537500	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	22	M00542981	SpyGlass Discover Imager™ II IOC Catheter.....	33
M00537510	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	22	M00542990	SpyGlass Discover Retrieval Basket Box.....	33
M00537520	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	22	M00543000	SpyGlass Discover Balloon Dilation Catheter	33
M00537530	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	22	M00543010	SpyGlass Discover Balloon Dilation Catheter	33
M00537540	Advanix Pancreatic Stent Kit – Straight, Leading Barb.....	22	M00543040	SpyGlass Discover Balloon Dilation Catheter	33
M00537550	Advanix Pancreatic Stent Kit – Straight, Leading Barb.....	22	M00543751	EXALT Model B Single-Use Bronchoscope - Large	77
M00537560	Advanix Pancreatic Stent Kit – Straight, Leading Barb.....	22	M00545000	RX Cytology Brush Wireguided Cytology Brush	32
M00538351	Rescue Retrieval Single-Use Grasping Forceps	72	M00545150	Autotome™ RX 49 Cannulating Sphincterotome.....	6
M00539200	Flexima™ Duodenal Bend Biliary Stent with Delivery System.....	20	M00545160	Autotome RX 49 Cannulating Sphincterotome.....	6
M00539210	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545170	Autotome RX 44 Cannulating Sphincterotome.....	6
M00539220	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545180	Autotome RX 44 Cannulating Sphincterotome.....	6
M00539230	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545190	Autotome RX 39 Cannulating Sphincterotome.....	6
M00539240	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545200	Autotome RX 39 Cannulating Sphincterotome.....	6
M00539250	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545211	RX Locking Device Olympus™ / FujiFilm™ Brand Endoscope Compatible.....	12
M00539260	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545261	RX Locking Device and Biopsy Cap System – Olympus FujiFilm G5 Series Brand Endoscope Compatible.....	12
M00539270	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545271	RX Locking Device and Biopsy Cap System – Pentax™ Brand Endoscope Compatible.....	12
M00539280	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545281	RX Locking Device and Biopsy Cap System – FujiFilm Type 3 Brand Endoscope Compatible	12
M00539290	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545550	RX Duodenal Bend Biliary Stent with RX Delivery System	20
M00539300	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545560	RX Duodenal Bend Biliary Stent with RX Delivery System	20
M00539310	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545570	RX Duodenal Bend Biliary Stent with RX Delivery System	20
M00539320	Flexima Duodenal Bend Biliary Stent with Delivery System....	20	M00545580	RX Duodenal Bend Biliary Stent with RX Delivery System	20
M00539330	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545590	RX Duodenal Bend Biliary Stent with RX Delivery System	20
M00539340	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545600	RX Duodenal Bend Biliary Stent with RX Delivery System	20
M00539350	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545610	RX Duodenal Bend Biliary Stent with RX Delivery System	20
M00539360	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545620	RX Duodenal Bend Biliary Stent with RX Delivery System	20
M00539370	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545630	RX Duodenal Bend Biliary Stent with RX Delivery System	20
M00539380	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545640	RX Duodenal Bend Biliary Stent with RX Delivery System	20
M00539390	Flexima Duodenal Bend Biliary Stent with Delivery System..	20	M00545650	RX Duodenal Bend Biliary Stent with RX Delivery System	20
M00540110	Flexima Single Pigtail Nasobiliary Catheter.....	21	M00545660	RX Duodenal Bend Biliary Stent with RX Delivery System	20
M00540120	Flexima Single Pigtail Nasobiliary Catheter.....	21	M00545670	RX Duodenal Bend Biliary Stent with RX Delivery System	20
M00540130	Flexima Single Pigtail Nasobiliary Catheter.....	21	M00545680	RX Duodenal Bend Biliary Stent with RX Delivery System	20
M00540140	Flexima Single Pigtail Nasobiliary Catheter.....	21	M00545690	RX Duodenal Bend Biliary Stent with RX Delivery System	20
M00542250	Speedband Superview Super 7™ Multiple Band Ligators.....	48	M00545700	RX Duodenal Bend Biliary Stent with RX Delivery System	20
M00542251	Speedband Superview Super 7™ Multiple Band Ligators.....	48			
M00542253	Speedband Superview Super 7 Multiple Band Ligators.....	48			

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00545710	RX Duodenal Bend Biliary Stent with RX Delivery System	20	M00547340	Extractor Pro RX-S Retrieval Balloon Catheter – Proximal Injection	28
M00545720	RX Duodenal Bend Biliary Stent with RX Delivery System	20	M00547350	Extractor Pro RX-S Retrieval Balloon Catheter – Proximal Injection	28
M00545730	RX Duodenal Bend Biliary Stent with RX Delivery System	20	M00550090	Amplatz Super Stiff™ Guidewires – Straight Tip	1, 44, 84
M00545740	RX Duodenal Bend Biliary Stent with RX Delivery System	20	M00550300	CRE™ Pulmonary Balloon Dilatation Catheter	84
M00545760	RX ERCP Cannula – Standard Tip	10	M00550310	CRE Pulmonary Balloon Dilatation Catheter	84
M00545770	RX ERCP Cannula – Tapered Tip	10	M00550320	CRE Pulmonary Balloon Dilatation Catheter	84
M00545780	RX ERCP Cannula – Ball Tip	10	M00550330	CRE Pulmonary Balloon Dilatation Catheter	84
M00545820	Rapid Exchange XL Cannula	10	M00550340	CRE Pulmonary Balloon Dilatation Catheter	84
M00545840	RX Needle Knife XL Triple-Lumen Needle Knife	6	M00550350	CRE Pulmonary Balloon Dilatation Catheter	84
M00545890	Hurricane™ RX Biliary Balloon Dilatation Catheter	31	M00550601	Alliance™ II Integrated Inflation / Lithotripsy Device and Alliance II Syringe / Gauge Assembly	55, 84
M00545900	Hurricane RX Biliary Balloon Dilatation Catheter	31	M00550620	Alliance II Integrated Inflation / Lithotripsy Device – Inflation Handle	55, 84
M00545910	Hurricane RX Biliary Balloon Dilatation Catheter	31	M00550630	CRE SteriFlate Disposable Inflation Device	55
M00545920	Hurricane RX Biliary Balloon Dilatation Catheter	31	M00552001	NovaGold™ High Performance Guidewire	3
M00545930	Hurricane RX Biliary Balloon Dilatation Catheter	31	M00552011	NovaGold High Performance Guidewire	3
M00545940	Hurricane RX Biliary Balloon Dilatation Catheter	31	M00552350	Acquire™ Pulmonary Endobronchial Ultrasound Fine Needle Biopsy Device	82
M00545950	Hurricane RX Biliary Balloon Dilatation Catheter	31	M00552400	Acquire™ Pulmonary Endobronchial Ultrasound Fine Needle Biopsy Device	2
M00545960	Hurricane RX Biliary Balloon Dilatation Catheter	31	M00553200	Inflation / Deflation Devices and Accessories – Pneumatic Inflator	55
M00546470	SpyBite™ Max Biopsy Forceps	34	M00553640	AXIOS™ Stent and Electrocautery Enhanced Delivery System	37
M00546550	SpyGlass™ Retrieval Basket	34	M00553650	AXIOS Stent and Electrocautery Enhanced Delivery System	37
M00546560	SpyGlass™ Retrieval Snare	34	M00553660	AXIOS Stent and Electrocautery Enhanced Delivery System	37
M00546600	SpyScope™ DS Access & Delivery Catheter	34	M00553670	AXIOS Stent and Electrocautery Enhanced Delivery System	37
M00546610	SpyScope DS II Access & Delivery Catheter	34	M00553680	AXIOS Stent and Electrocautery Enhanced Delivery System	37
M00546620	Biliary Electrohydraulic Lithotripter Probe	35	M00553690	AXIOS Stent and Electrocautery Enhanced Delivery System	37
M00546650	SpyGlass DS Digital Controller	34	M00553740	AXIOS Stent and Delivery System	37
M005466800	Autolith™ Touch Bipolar Electrohydraulic Lithotripter	35	M00553750	AXIOS Stent and Delivery System	37
M00546750	Autolith Touch Extender Cable	35	M00553831	Lumicoil™ Platinum Fiducial Marker – Straight	37
M00546760	Autolith Touch Foot Switch	35	M00553841	Lumicoil Platinum Fiducial Marker – Figure 8	37
M00546780	SpyGlass™ Discover Digital Catheter	33	M00554500	Rigiflex™ II Single-Use Achalasia Balloon Dilators	54
M00547000	Extractor™ Pro RX Retrieval Balloon Catheter – Distal Injection	28	M00554510	Rigiflex II Single-Use Achalasia Balloon Dilators	54
M00547010	Extractor Pro RX Retrieval Balloon Catheter – Distal Injection	28	M00554520	Rigiflex II Single-Use Achalasia Balloon Dilators	54
M00547020	Extractor Pro RX Retrieval Balloon Catheter – Distal Injection	28	M00555500	Expect™ Slimline (SL) Endoscopic Ultrasound Aspiration Needle	36
M00547030	Extractor Pro RX Retrieval Balloon Catheter – Proximal Injection	28	M00555501	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle	36
M00547040	Extractor Pro RX Retrieval Balloon Catheter – Proximal Injection	28	M00555510	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle	36
M00547050	Extractor Pro RX Retrieval Balloon Catheter – Proximal Injection	28	M00555520	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle	36
M00547100	Extractor Pro XL Retrieval Balloon Catheter – Distal Injection	29	M00555530	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle	36
M00547110	Extractor Pro XL Retrieval Balloon Catheter – Distal Injection	29	M00555511	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle	36
M00547120	Extractor Pro XL Retrieval Balloon Catheter – Distal Injection	29	M00555521	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle	36
M00547130	Extractor Pro XL Retrieval Balloon Catheter – Proximal Injection	29	M00555531	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle	36
M00547140	Extractor Pro XL Retrieval Balloon Catheter – Proximal Injection	29	M00555540	Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device	36
M00547150	Extractor Pro XL Retrieval Balloon Catheter – Proximal Injection	29	M00555541	Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device	36
M00547300	Extractor Pro RX-S Retrieval Balloon Catheter – Distal Injection	28	M00555560	Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device	36
M00547310	Extractor Pro RX-S Retrieval Balloon Catheter – Distal Injection	28	M00555561	Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device	36
M00547320	Extractor Pro RX-S Retrieval Balloon Catheter – Distal Injection	28			
M00547330	Extractor Pro RX-S Retrieval Balloon Catheter – Proximal Injection	28			

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M0055580	Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device	36	M0055631	Jagwire High Performance Guidewire – Stiff Shaft – Straight	2
M00555760	WallFlex™ Soft Colonic Stent System	38	M0055641	Jagwire High Performance Guidewire – Straight	2
M00555770	WallFlex Soft Colonic Stent System	38	M0055651	Jagwire High Performance Guidewire – Angled	2
M00555780	WallFlex Soft Colonic Stent System	38	M0055661	Jagwire High Performance Guidewire – Stiff Shaft – Angled	2
M00555860	WallFlex Soft Duodenal Stent System	39	M0055681	Jagwire Extendable High Performance Guidewire – Straight	2
M00555870	WallFlex Soft Duodenal Stent System	39	M00556851	Jagwire Extendable High Performance Guidewire – Angled	2
M00555880	WallFlex Soft Duodenal Stent System	39	M00556861	Jagwire Extendable High Performance Guidewire – Stiff Shaft – Straight	2
M00556001	Hydra Jagwire™ High Performance Guidewire – Straight	1	M00556871	Jagwire Extendable High Performance Guidewire – Stiff Shaft – Angled	2
M00556011	Hydra Jagwire High Performance Guidewire – Angled	1	M00556902	Jagtail™ High Performance Guidewire Extension	2
M00556021	Hydra Jagwire High Performance Guidewire – Stiff Shaft – Straight	1	M00557001	Jagwire Revolution High Performance Guidewire – Straight	2
M00556031	Hydra Jagwire High Performance Guidewire – Stiff Shaft – Angled	1	M00557011	Jagwire Revolution High Performance Guidewire – Angled	2
M00556041	Hydra Jagwire High Performance Guidewire – Straight	1	M00557021	Jagwire Revolution High Performance Guidewire – Straight	2
M00556051	Hydra Jagwire High Performance Guidewire – Angled	1	M00557031	Jagwire Revolution High Performance Guidewire – Angled	2
M00556061	Hydra Jagwire High Performance Guidewire – Stiff Shaft – Straight	1	M00557041	Jagwire Revolution High Performance Guidewire – Straight	2
M00556071	Hydra Jagwire™ High Performance Guidewire – Stiff Shaft – Angled	1	M00558220	Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle	82
M00556101	Dreamwire™ High Performance Guidewire – Standard Shaft	1	M00558250	Expect Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle	82
M00556111	Dreamwire High Performance Guidewire – Standard Shaft	1	M00558330	CRE™ Fixed Wire Esophageal Balloon Dilatation Catheters	53
M00556121	Dreamwire High Performance Guidewire – Stiff Shaft	1	M00558340	CRE Fixed Wire Esophageal Balloon Dilatation Catheters	53
M00556131	Dreamwire High Performance Guidewire – Stiff Shaft	1	M00558350	CRE Fixed Wire Esophageal Balloon Dilatation Catheters	53
M00556141	Dreamwire High Performance Guidewire – Standard Shaft	1	M00558360	CRE Fixed Wire Esophageal Balloon Dilatation Catheters	53
M00556151	Dreamwire High Performance Guidewire – Standard Shaft	1	M00558361	CRE Fixed Wire Esophageal Balloon Dilatation Catheters	53
M00556161	Dreamwire High Performance Guidewire – Stiff Shaft	1	M00558370	CRE Fixed Wire Esophageal Balloon Dilatation Catheters	53
M00556171	Dreamwire High Performance Guidewire – Stiff Shaft	1	M00558371	CRE Fixed Wire Esophageal Balloon Dilatation Catheters	53
M00556181	NaviPro™ Hydrophilic Guidewire	3	M00558380	CRE Fixed Wire Esophageal Balloon Dilatation Catheters	53
M00556191	NaviPro Hydrophilic Guidewire	3	M00558381	CRE Fixed Wire Esophageal Balloon Dilatation Catheters	53
M00556201	NaviPro Hydrophilic Guidewire	3	M00558390	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	52
M00556211	NaviPro Hydrophilic Guidewire	3	M00558400	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	52
M00556221	NaviPro Hydrophilic Guidewire	3	M00558410	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	52
M00556231	NaviPro Hydrophilic Guidewire	3	M00558420	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	52
M00556241	NaviPro Hydrophilic Guidewire	3	M00558430	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	52
M00556251	NaviPro Hydrophilic Guidewire	3	M00558440	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	52
M00556461	Jagwire High Performance Guidewire – Straight	2	M00558450	CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	52
M00556481	Stiff Shaft Jagwire Guidewire – Straight	2	M00558460	CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	52
M00556491	Stiff Shaft Jagwire Guidewire – Angled	2	M00558470	CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	52
M00556511	Jagwire™ High Performance Guidewire – Stiff Shaft – Angled	2	M00558480	CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	52
M00556521	Jagwire High Performance Guidewire – Stiff Shaft – Straight	2	M00558490	CRE Wireguided Esophageal / Pyloric / Colonic Balloon	
M00556531	Jagwire High Performance Guidewire – Stiff Shaft – Angled	2			
M00556561	Jagwire High Performance Guidewire – Straight	2			
M00556581	Jagwire High Performance Guidewire – Straight	2			
M00556591	Jagwire High Performance Guidewire – Angled	2			
M00556601	Jagwire High Performance Guidewire – Stiff Shaft – Straight	2			
M00556611	Jagwire High Performance Guidewire – Stiff Shaft – Angled	2			
M00556621	Jagwire High Performance Guidewire – Straight	2			

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
	Dilatation Catheters	52		Standard Plug	46
M00558500	CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	52	M00560101	Gold Probe Single-Use Electrohemostasis Catheters – Standard Plug	46
M00558600	CRE PRO Wireguided Biliary Balloon Dilators.....	30	M00560150	Injection Gold Probe™ Single-Use Electrohemostasis / Injection Catheters – Standard Plug.....	46
	CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	51	M00560160	Injection Gold Probe Single-Use Electrohemostasis / Injection Catheters – Standard Plug.....	46
M00558610	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	30	M00560220	Gold Probe Single-Use Electrohemostasis Catheters – Standard Plug	46
	CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	51	M00560311	Sensation™ Single-Use Medium Stiff Wire Short Throw Snares – Medium Oval	70
M00558620	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	30	M00560321	Sensation Single-Use Medium Stiff Wire Short Throw Snares – Small Oval	70
	CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	51	M00560331	Sensation Single-Use Medium Stiff Wire Short Throw Snares – Medium Crescent	70
M00558630	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	30	M00561101	Captivator™ COLD Single-Use Snares.....	67
	CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	51	M00561102	Captivator COLD Single-Use Snares.....	67
M00558640	CRE™ PRO Wireguided Biliary Balloon Dilatation Catheters ..	30	M00561103	Captivator COLD Single-Use Snares.....	67
	CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	51	M00561191	Captivator II Single-Use Snares – Rounded Stiff	68
M00558650	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	30	M00561192	Captivator II Single-Use Snares – Rounded Stiff	68
	CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	51	M00561193	Captivator II Single-Use Snares – Rounded Stiff	68
M00558660	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	30	M00561221	Captivator II Single-Use Snares – Rounded Stiff	68
	CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	51	M00561222	Captivator II Single-Use Snares – Rounded Stiff	68
M00558670	CRE™ PRO Wireguided Biliary Balloon Dilatation Catheters ..	30	M00561223	Captivator II Single-Use Snares – Rounded Stiff	68
	CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	51	M00561231	Captivator II Single-Use Snares – Rounded Stiff	68
M00558680	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	30	M00561232	Captivator II Single-Use Snares – Rounded Stiff	68
	CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	51	M00561233	Captivator II Single-Use Snares – Rounded Stiff	68
M00558690	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	30	M00561241	Captivator II Single-Use Snares – Rounded Stiff	68
	CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	51	M00561242	Captivator II Single-Use Snares – Rounded Stiff	68
M00558700	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	30	M00561243	Captivator II Single-Use Snares – Rounded Stiff	68
	CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	51	M00561250	Active Cords – Fits Olympus™ Brand Generator	12, 71
M00558710	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	30	M00561270	Active Cord, Fits Endostat™, Endostat II and Endostat III Generators, Bovie™, Valleylab™, Aspen Labs™, and ERBE™ Brand Generators	12, 71
	CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	51	M00561291	Captivator II Single-Use Snare – Extra Large Rounded	68
M00558731	Expect™ Pulmonary Adaptor.....	82	M00561311	Captivator Single-Use Snares – Large Oval Medium	68
M00558860	CRE RX Biliary Balloon Dilatation Catheters	30, 50	M00561440	Electrosurgical Bipolar Cable Adapter.....	49
M00558870	CRE RX Biliary Balloon Dilatation Catheters	30, 50	M00561450	Electrosurgical Bipolar Cable Adapter.....	49
M00558880	CRE RX Biliary Balloon Dilatation Catheters	30, 50	M00561460	Electrosurgical Bipolar Cable Adapter.....	49
M00558890	CRE RX Biliary Balloon Dilatation Catheters	30, 50	M00561470	Electrosurgical Bipolar Cable Adapter.....	49
M00558900	CRE RX Biliary Balloon Dilatation Catheters	30, 50	M00561480	Electrosurgical Bipolar Cable Adapter.....	49
M00558910	CRE RX Biliary Balloon Dilatation Catheters	30, 50	M00561490	Electrosurgical Bipolar Cable Adapter.....	49
M00558920	CRE RX Biliary Balloon Dilatation Catheters	30, 50	M00561500	Electrosurgical Bipolar Cable Adapter.....	49
M00558930	CRE RX Biliary Balloon Dilatation Catheters	30, 50	M00561600	Captivator EMR Device for Standard Gastrosopes.....	67
M00558940	CRE RX Biliary Balloon Dilatation Catheters	30, 50	M00561601	Captivator EMR Device for Standard Gastrosopes.....	67
M00560070	Gold Probe™ Single-Use Electrohemostasis Catheters – Standard Plug	46	M00561610	Captivator EMR Device for Large Gastrosopes	67
M00560071	Gold Probe Single-Use Electrohemostasis Catheters – Standard Plug	46	M00561611	Captivator EMR Device for Large Gastrosopes	67
M00560100	Gold Probe Single-Use Electrohemostasis Catheters –		M00561821	Single-Use Rotatable Snares – Small Oval.....	69
			M00561831	Single-Use Rotatable Snares – Medium Oval.....	69

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00562301	Captivator Single-Use Snares – Small Oval	68	M00563150	EndoVive One Step Button Low Profile PEG Kit, Pull	63
M00562321	Captivator Single-Use Snares – Medium Oval	68	M00563160	EndoVive One Step Button Low Profile PEG Kit, Pull	63
M00562341	Captivator Single-Use Snares – Medium Hexagonal	68	M00563180	EndoVive One Step Button Low Profile PEG Kit, Push	63
M00562371	Captivator Single-Use Snares – Medium Crescent	68	M00563190	EndoVive One Step Button Low Profile PEG Kit, Push	63
M00562391	Captivator Single-Use Snares – Large Oval	68	M00564101	eXcelon™ Transbronchial Aspiration Needle	83
M00562401	Captiflex™ Single-Use Snares – Medium Oval	70	M00564111	eXcelon Transbronchial Aspiration Needle	83
M00562402	Captiflex Single-Use Snares – Medium Oval	70	M00564121	eXcelon Transbronchial Aspiration Needle	83
M00562421	Captiflex Single-Use Snares – Small Oval	70	M00565010	WallFlex™ Duodenal Stent System with Anchor Lock Delivery System	39
M00562422	Captiflex Single-Use Snares – Small Oval	70	M00565020	WallFlex Duodenal Stent System with Anchor Lock Delivery System	39
M00562451	Captivator Single-Use Snares – Small Hexagonal	68	M00565030	WallFlex Duodenal Stent System with Anchor Lock Delivery System	39
M00562471	Captiflex Single-Use Snares – Extra Small Oval	70	M00565040	WallFlex Colonic Stent System with Anchor Lock Delivery System	38
M00562531	Profile™ Single-Use Snares – Extra Small Oval	69	M00565050	WallFlex Colonic Stent System with Anchor Lock Delivery System	38
M00562551	Profile Single-Use Snares – Small Oval	69	M00565060	WallFlex Colonic Stent System with Anchor Lock Delivery System	38
M00562571	Profile™ Single-Use Snares – Medium Oval	69	M00565070	WallFlex Colonic Stent System with Anchor Lock Delivery System	38
M00562651	Sensation™ Single-Use Short Throw Snares – Large Oval / Medium Stiff	70	M00565080	WallFlex Colonic Stent System with Anchor Lock Delivery System	38
M00562652	Sensation Single-Use Short Throw Snares – Large Oval / Medium Stiff	70	M00565090	WallFlex Colonic Stent System with Anchor Lock Delivery System	38
M00562671	Sensation Single-Use Short Throw Snares – Medium Oval / Flexible	70	M00565100	WallFlex Colonic Stent System with Anchor Lock Delivery System	38
M00562672	Sensation Single-Use Short Throw Snares – Medium Oval / Flexible	70	M00565110	WallFlex Colonic Stent System with Anchor Lock Delivery System	38
M00562673	Sensation Single-Use Short Throw Snares – Medium Oval / Flexible	70	M00565120	WallFlex Colonic Stent System with Anchor Lock Delivery System	38
M00562691	Sensation Single-Use Short Throw Snares – Small Oval / Flexible	70	M00565130	WallFlex Colonic Stent System with Anchor Lock Delivery System	38
M00562692	Sensation Single-Use Short Throw Snares – Small Oval / Flexible	70	M00565140	WallFlex Colonic Stent System with Anchor Lock Delivery System	38
M00562693	Sensation Single-Use Short Throw Snares – Small Oval / Flexible	70	M00565150	WallFlex Colonic Stent System with Anchor Lock Delivery System	38
M00562760	EndoVive™ Low Profile Replacement Button	65	M00565560	WALLSTENT™ Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP™ Plus Delivery System	40
M00562770	EndoVive Low Profile Replacement Button	65	M00565570	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	40
M00562780	EndoVive Low Profile Replacement Button	65	M00565580	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	40
M00562790	EndoVive Low Profile Replacement Button	65	M00565590	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	40
M00562800	EndoVive Low Profile Replacement Button	65	M00565620	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	40
M00562820	EndoVive Low Profile Replacement Button	65	M00565630	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	40
M00562830	EndoVive Low Profile Replacement Button	65	M00565640	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	40
M00562840	EndoVive Low Profile Replacement Button	65	M00565650	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	40
M00562850	EndoVive Low Profile Replacement Button	65			
M00563010	EndoVive One Step Button™ Low Profile PEG Kit, Pull	63			
M00563020	EndoVive One Step Button Low Profile PEG Kit, Pull	63			
M00563030	EndoVive One Step Button Low Profile PEG Kit, Pull	63			
M00563040	EndoVive One Step Button Low Profile PEG Kit, Pull	63			
M00563080	EndoVive One Step Button Low Profile PEG Kit, Push	63			
M00563090	EndoVive One Step Button Low Profile PEG Kit, Push	63			
M00563100	EndoVive One Step Button Low Profile PEG Kit, Push	63			
M00563120	EndoVive One Step Button Low Profile PEG Kit, Pull	63			
M00563130	EndoVive One Step Button Low Profile PEG Kit, Pull	63			
M00563140	EndoVive One Step Button Low Profile PEG Kit, Pull	63			

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00566001	Rapid Refill Continuous Injection System.....	12	M00570650	WallFlex Biliary RX Uncovered Stent System.....	13
M00566320	EndoVive Through-the-PEG (TTP) Jejunal Feeding Tube – Bent Tip.....	66	M00570660	WallFlex Biliary RX Uncovered Stent System.....	13
M00566330	EndoVive Through-the-PEG (TTP) Jejunal Feeding Tube – Pigtail Tip.....	66	M00570700	WallFlex Biliary RX Partially Covered Stent System	13
M00566340	EndoVive Through-the-PEG (TTP) Jejunal Feeding Tube – Bent Tip.....	66	M00570710	WallFlex Biliary RX Partially Covered Stent System	13
M00566350	EndoVive Through-the-PEG (TTP) Jejunal Feeding Tube – Pigtail Tip.....	66	M00570720	WallFlex Biliary RX Partially Covered Stent System	13
M00566461	EndoVive Safety PEG Kit – Pull.....	62	M00570730	WallFlex Biliary RX Partially Covered Stent System	13
M00566471	EndoVive Safety PEG Kit – Push.....	62	M00570740	WallFlex Biliary RX Partially Covered Stent System	13
M00566481	EndoVive Safety PEG Kit – Pull.....	62	M00570890	WallFlex Biliary RX Uncovered Stent System.....	13
M00566491	EndoVive Safety PEG Kit – Push.....	62	M00572230	Epic™ Biliary Stent	14
M00566670	Encore™ 26 Inflation Device	32, 33	M00572240	Epic Biliary Stent	14
M00567330	MaxForce™ Biliary Balloon Dilatation Catheter	32	M00572250	Epic Biliary Stent	14
M00567340	MaxForce Biliary Balloon Dilatation Catheter.....	32	M00572260	Epic Biliary Stent	14
M00567370	MaxForce Biliary Balloon Dilatation Catheter.....	32	M00572270	Epic Biliary Stent	14
M00567380	MaxForce Biliary Balloon Dilatation Catheter.....	32	M00572280	Epic Biliary Stent	14
M00567410	MaxForce Biliary Balloon Dilatation Catheter.....	32	M00572290	Epic Biliary Stent	14
M00567420	MaxForce Biliary Balloon Dilatation Catheter.....	32	M00572300	Epic Biliary Stent	14
M00568201	EndoVive™ Standard PEG Kit – Pull.....	62	M00572310	Epic Biliary Stent	14
M00568211	EndoVive Standard PEG Kit – Push.....	62	M00572320	Epic Biliary Stent	14
M00568241	EndoVive Standard PEG Kit – Pull	62	M00572330	Epic Biliary Stent	14
M00568251	EndoVive Standard PEG Kit – Push.....	62	M00572340	Epic Biliary Stent	14
M00569610	WALLSTENT™ RX Biliary Endoprosthesis Stent System.....	15	M00573000	Jagtome™ RX 49 Sphincterotome (Preloaded with Jagwire™ Guidewire)	5
M00569620	WALLSTENT RX Biliary Endoprosthesis Stent System.....	15	M00573010	Jagtome RX 49 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00569630	WALLSTENT RX Biliary Endoprosthesis Stent System.....	15	M00573020	Jagtome RX 49 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00569640	WALLSTENT RX Biliary Endoprosthesis Stent System.....	15	M00573030	Jagtome RX 49 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00569650	WALLSTENT RX Biliary Endoprosthesis Stent System.....	15	M00573040	Jagtome RX 44 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00569660	WALLSTENT RX Biliary Endoprosthesis Stent System.....	15	M00573050	Jagtome RX 44 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00569670	WALLSTENT RX Biliary Endoprosthesis Stent System.....	15	M00573060	Jagtome RX 44 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00569680	WALLSTENT RX Biliary Endoprosthesis Stent System.....	15	M00573070	Jagtome RX 44 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00569690	WALLSTENT RX Biliary Endoprosthesis Stent System.....	15	M00573080	Jagtome RX 39 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00569700	WALLSTENT RX Biliary Endoprosthesis Stent System.....	15	M00573090	Jagtome RX 39 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00569710	WALLSTENT RX Biliary Endoprosthesis Stent System.....	15	M00573100	Jagtome RX 39 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00569720	WALLSTENT RX Biliary Endoprosthesis Stent System.....	15	M00573110	Jagtome RX 39 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00569800	WALLSTENT RX Biliary Endoprosthesis Stent System.....	15	M00576300	Ultraflex™ Single-Use Uncovered Tracheobronchial Stent System – Distal Release	79
M00569810	WALLSTENT RX Biliary Endoprosthesis Stent System.....	15	M00576310	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release	79
M00570340	WallFlex™ Biliary RX Fully Covered Stent System RMV	13	M00576320	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release	79
M00570350	WallFlex Biliary RX Fully Covered Stent System RMV	13	M00576330	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release	79
M00570360	WallFlex Biliary RX Fully Covered Stent System RMV	13			
M00570370	WallFlex Biliary RX Fully Covered Stent System RMV	13			
M00570380	WallFlex Biliary RX Fully Covered Stent System RMV	13			
M00570600	WallFlex Biliary RX Uncovered Stent System.....	13			
M00570610	WallFlex Biliary RX Uncovered Stent System.....	13			
M00570620	WallFlex Biliary RX Uncovered Stent System.....	13			
M00570630	WallFlex Biliary RX Uncovered Stent System.....	13			
M00570640	WallFlex Biliary RX Uncovered Stent System.....	13			

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00576340	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release	79	M00576650	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80
M00576350	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release	79	M00576660	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80
M00576360	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release	79	M00576670	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80
M00576370	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release	79	M00576680	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80
M00576380	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release	79	M00576690	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80
M00576410	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	79	M00580171	EndoVive™ Low Profile Replacement Button Accessories – 18Fr Bolus Feeding Set, Straight Adapter	63
M00576420	Ultraflex™ Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	79	M00580181	EndoVive Low Profile Replacement Button Accessories – 18Fr Right Angle Feeding Set, Continuous Adapter	63
M00576430	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	79	M00580231	EndoVive Low Profile Replacement Button Accessories – 24Fr Bolus Feeding Set, Straight Adapter	63
M00576440	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	79	M00580241	EndoVive Low Profile Replacement Button Accessories – 24Fr Right Angle Feeding Set, Continuous Adapter	63
M00576450	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	79	M00580401	EndoVive Low Profile Replacement Button Accessories – 18Fr Decompression Tube	63
M00576460	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	79	M00580411	EndoVive Low Profile Replacement Button Accessories – 18Fr Decompression Tube	63
M00576470	Ultraflex™ Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	79	M00580421	EndoVive Low Profile Replacement Button Accessories – 18Fr Decompression Tube	63
M00576480	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	79	M00580431	EndoVive Low Profile Replacement Button Accessories – 18Fr Decompression Tube	63
M00576490	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	79	M00580441	EndoVive Low Profile Replacement Button Accessories – 18Fr Decompression Tube	63
M00576520	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80	M00580461	EndoVive Low Profile Replacement Button Accessories – 24Fr Decompression Tube	63
M00576530	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80	M00580471	EndoVive Low Profile Replacement Button Accessories – 24Fr Decompression Tube	63
M00576540	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80	M00580481	EndoVive Low Profile Replacement Button Accessories – 24Fr Decompression Tube	63
M00576550	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80	M00580491	EndoVive Low Profile Replacement Button Accessories – 24Fr Decompression Tube	63
M00576560	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80	M00580641	EndoVive Initial Placement PEG Accessories – Y-Ports	62
M00576570	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80	M00580651	EndoVive Initial Placement PEG Accessories – Y-Ports	62
M00576580	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80	M00582040	EndoVive Replacement Gastrostomy Feeding Tubes, Straight Bolster	64
M00576590	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80	M00582050	EndoVive Replacement Gastrostomy Feeding Tubes, Straight Bolster	64
M00576600	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80	M00582060	EndoVive Replacement Gastrostomy Feeding Tubes, Straight Bolster	64
M00576610	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80	M00582070	EndoVive Replacement Gastrostomy Feeding Tubes, Straight Bolster	64
M00576620	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80	M00582080	EndoVive Replacement Gastrostomy Feeding Tubes, Straight Bolster	64
M00576630	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80	M00582090	EndoVive Replacement Gastrostomy Feeding Tubes, Straight Bolster	64
M00576640	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	80	M00582100	EndoVive Replacement Gastrostomy Feeding Tubes, Straight Bolster	64
			M00582110	EndoVive Replacement Gastrostomy Feeding Tubes, Right Angle Bolster	64

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00582120	EndoVive Replacement Gastrostomy Feeding Tubes, Right Angle Bolster	64	M00583190	TRUEtome Dreamwire 44 Cannulating Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	7
M00582130	EndoVive Replacement Gastrostomy Feeding Tubes, Right Angle Bolster	64	M00583200	TRUEtome Dreamwire 44 Cannulating Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	7
M00582140	EndoVive Replacement Gastrostomy Feeding Tubes, Right Angle Bolster	64	M00583210	TRUEtome Dreamwire 44 Cannulating Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	7
M00582150	EndoVive Replacement Gastrostomy Feeding Tubes, Right Angle Bolster	64	M00583260	TRUEtome Revolution RX Cannulating Sphincterotome (Preloaded with Jagwire Revolution High Performance Guidewire)	7
M00582160	EndoVive™ Replacement Gastrostomy Feeding Tubes, Right Angle Bolster	64	M00583270	TRUEtome Revolution RX Cannulating Sphincterotome (Preloaded with Jagwire Revolution High Performance Guidewire)	7
M00582170	EndoVive Replacement Gastrostomy Feeding Tubes, Right Angle Bolster	64	M00583280	TRUEtome Revolution RX Cannulating Sphincterotome (Preloaded with Jagwire Revolution High Performance Guidewire)	7
M00582860	EndoVive Replacement Gastrostomy Feeding Tubes, Straight Bolster	64	M00583290	TRUEtome Revolution RX Cannulating Sphincterotome (Preloaded with Jagwire Revolution High Performance Guidewire)	7
M00582870	EndoVive Replacement Gastrostomy Feeding Tubes, Right Angle Bolster	64	M00584000	Dreamtome™ RX 49 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	4
M00583000	Hydratome™ RX 49 Sphincterotome (Preloaded with Hydra Jagwire™ Guidewire)	4	M00584010	Dreamtome RX 49 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	4
M00583010	Hydratome RX 49 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	4	M00584020	Dreamtome RX 49 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	4
M00583020	Hydratome™ RX 49 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	4	M00584030	Dreamtome RX 49 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	4
M00583030	Hydratome RX 49 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	4	M00584040	Dreamtome RX 44 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	4
M00583040	Hydratome RX 44 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	4	M00584050	Dreamtome RX 44 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	4
M00583050	Hydratome RX 44 Sphincterotome (Preloaded with Hydra Jagwire™ Guidewire)	4	M00584060	Dreamtome RX 44 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	4
M00583060	Hydratome RX 44 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	4	M00584070	Dreamtome RX 44 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	4
M00583070	Hydratome RX 44 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	4	M00584120	TRUEtome 39 Sphincterotome	8
M00583080	TRUEtome™ Jag 39 Cannulating Sphincterotome (Preloaded with Jagwire High Performance Guidewire)	7	M00584150	TRUEtome 49 Sphincterotome	8
M00583090	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with Jagwire High Performance Guidewire)	7	M00584160	TRUEtome 49 Sphincterotome	8
M00583100	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with Jagwire High Performance Guidewire)	7	M00584170	TRUEtome 44 Sphincterotome	8
M00583110	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with Jagwire High Performance Guidewire)	7	M00584180	TRUEtome 44 Sphincterotome	8
M00583120	TRUEtome Jag 39 Cannulating Sphincterotome (Preloaded with Jagwire High Performance Guidewire)	7	M00584190	TRUEtome 39 Sphincterotome	8
M00583130	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with Jagwire High Performance Guidewire)	7	M00584240	Jagtome™ Revolution RX 39 Sphincterotome (Preloaded with Jagwire™ Revolution High Performance Guidewire)	5
M00583140	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with Jagwire High Performance Guidewire)	7	M00584250	Jagtome Revolution RX 39 Sphincterotome (Preloaded with Jagwire Revolution High Performance Guidewire)	5
M00583150	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with Jagwire High Performance Guidewire)	7	M00584260	Jagtome Revolution RX 39 Sphincterotome (Preloaded with Jagwire Revolution High Performance Guidewire)	5
M00583160	TRUEtome Hydra 44 Cannulating Sphincterotome (Preloaded with Hydra Jagwire High Performance Guidewire)	7	M00584270	Jagtome Revolution RX 39 Sphincterotome (Preloaded with Jagwire Revolution High Performance Guidewire)	5
M00583170	TRUEtome Hydra 44 Cannulating Sphincterotome (Preloaded with Hydra Jagwire High Performance Guidewire)	7	M00589060	MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter	54
M00583180	TRUEtome Dreamwire 44 Cannulating Sphincterotome (Preloaded with Dreamwire™ High Performance Guidewire)	7	M00589080	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	54
			M00589100	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	54

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00589120	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	54	SCT-351-10	Hydra Water Bottle Cap	74
M00589140	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	54	SEE-220-4	Pure™ Enzymatic Detergent - Spearmint Scent	73
M00589150	MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter	54	SEE-573-4	Pure Enzymatic Detergent - Fragrance Free	73
M00589160	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	54	SES-240-50	Koala Cleaning Sponge plus Pure Enzymatic with Ecobowl	75
M00589180	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	54	SIT-355-15	Hydra Irrigation Tubing (Compatible with ERBE, Medivators Pumps)	74
M005ATS20000	Alair™ Bronchial Thermoplasty System Controller	78	SIT-361-100	Hydra Disposable Auxiliary Water Jet Connector	74
M005ATS25010	Alair Bronchial Thermoplasty System Catheter	78	SIT-576-15	Hydra Irrigation Tubing (Compatible with Olympus Pumps)	74
M005ATS20110	Alair RF Controller Accessory Kit	78	SLT-394-100	Endoglide™+ Sterile Bacteriostatic Lubricating Gel	74
M005XXX1	RX Cholangiogram Kit-Preloaded – Autotome™ Sphincterotome – RX Locking Device	6	SLT-395-50	Endoglide+ Sterile Bacteriostatic Lubricating Gel	74
M005XXX1	RX Cholangiogram Kits – Autotome or Ultratome™ XL Sphincterotome	6	SLT-612-50	Sterile Lubricating Jelly (56.7g/2oz)	74
SBB-119-50	HedgeHog™ Valve Brush	73	SPX-017-25	CinchPad™ Endoscope Transport Pad	75
SBB-382-50	HedgeHog™ Double Bundle (Double-End Brush + Valve Brush)	73	SSB-210-50	Koala™ Cleaning Sponge	75
SBC-501-100	Seal Single-Use Biopsy Valve	76	STE-297-50	TrapEase™ Single Chamber Polyp Trap	72
SBD-228-50	HedgeHog Double-End Channel Brush	73	STE-298-50	TrapEase Quad Chamber Polyp Trap	72
SBD-289-50	HedgeHog Dual-End Channel / Valve Brush	73	SUV-617-50	Orca™ Sterile Single-Use Air/Water & Suction Valves	76
SBD-291-50	HedgeHog Sweeper Channel Brush	73	SUV-628-50	OrcaPod™ 3 - Orca Sterile Single-Use Air/Water, Suction, and Seal Biopsy Valves	76
SBS-227-50	HedgeHog Single-End Channel Brush	73	SUV-629-50	OrcaPod 4 - Orca Sterile Single-Use Air/Water, Suction, and Seal Biopsy Valves with WaterJet Connector	76
SBT-114-100	Blox™ Endoscopic Bite Block with Elastic Cloth Strap (54 FR)	73	SVT-089-100	Suction Tubing, NS, 3/16" x 6'	75
SBT-321-100	Blox Endoscopic Bite Block with Elastic Cloth Strap (54 FR) (w/o Dental Retention)	73	SVT-090-100	Suction Tubing, NS, 3/16" x 10'	75
SBT-546-100	Blox Endoscopic Bite Block with Elastic Cloth Strap (60 FR)	73	SVT-095-100	Suction Tubing, NS, 3/16" x 12'	75
SCT-350-10	Hydra™ Water Bottle Cap	74	SVT-143-100	Suction Tubing, 1/4" X 6'	75
			SVT-144-60	Suction Tubing, NS 1/4" x 10'	75
			SVT-145-50	Suction Tubing, NS, 1/4" x 12'	75

Notes and Future Product Updates



Section Table of Contents



Prev Page

Next Page



Home





Terms & Conditions

► Begins on Page 102

Certain Terms and Conditions

ORDERING INFORMATION

All orders should be placed directly with Boston Scientific Endoscopy Customer Service by mail, telephone, fax or Electronic Data Interchange (EDI).

Boston Scientific Endoscopy
500 Commander Shea Blvd.
Quincy, MA 02171
(888) 272-1001 Phone
(888) 272-9444 Fax
(877) 272-3340 EDI

Please contact EDI if interested in ordering your endoscopy products via the internet.

Please note that the above address is our mailing address only.

Customer Service is available between the hours of 8:30 am and 9:30 pm (Eastern Standard Time) weekdays excluding national holidays. Our toll-free fax number, (888) 272-9444, is available at all times. In addition, technical and emergency service is available 24 hours a day, 7 days a week, call (888) 272-1001 for further instructions.

Standing Orders

Boston Scientific Endoscopy accepts standing orders on all products. This convenience provides you with scheduled shipments using a single purchase order. Contact the Customer Service Database Department at (888) 325-1060 for details.

Terms

Terms are Net 30 Days. All shipments are F.O.B. Quincy, Massachusetts. Prices are subject to change without notice.

Returned Goods

Permission to return merchandise is required for all returns and, except as provided herein, may be withheld in the sole discretion of Boston Scientific Endoscopy. This authorization must be requested by calling or writing:

Customer Service

Boston Scientific Endoscopy

500 Commander Shea Blvd.

Quincy, MA 02171

(800) 862-1284

The request for authorization must include:

- Reason for the return – “Ordered in Error,” “Shipped in Error,” etc.
- Catalog number, size, quantity, lot number
- Purchase order number, Boston Scientific Endoscopy invoice number and date of invoice

All merchandise returned to Boston Scientific Endoscopy must be accompanied by a copy of the packing slip or invoice enclosed with the original shipment. The Returned Goods Authorization (RGA) number must be clearly marked on the shipping container. Merchandise returned without an RGA number will not be accepted or considered as returned product.

Merchandise that will not be accepted for return:

- Expired product or products with less than 6 months shelf life
- Merchandise manufactured as SPECIALS or CUSTOM
- Merchandise in unsaleable condition (opened, damaged, tampered, etc.)
- Merchandise not listed in this Price List or any Supplement
- Merchandise sterilized or used by customer
- Merchandise that is not packaged properly for shipping when returned
- Incomplete kits or boxes (unless defective)
- Unauthorized returns

► Certain Terms and Conditions

Returned Goods (continued)

Allowance Schedule for Authorized Returns:

Allowance Schedule for Capital Equipment:

- Capital Equipment processed in error by Boston Scientific Endoscopy will receive full credit if returned in saleable condition
- Capital Equipment cannot be returned once the customer has taken title of the equipment

Allowance Schedule for Metal and Polyflex™ Stents and Alair™ Bronchial Thermoplasty Products:

- Merchandise processed in error by Boston Scientific Endoscopy will receive full credit if returned in saleable condition
- Merchandise returned within 30 days of invoice date in saleable condition will receive 100% credit
- Merchandise returned from 31 days to 90 days in saleable condition will receive full credit less a 30% restocking charge
- Merchandise returned from 91 days to 150 days in saleable condition will receive full credit less a 50% restocking fee
- Merchandise returned past 150 days will not receive any credit

Allowance Schedule for All Other Products:

- Merchandise processed in error by Boston Scientific Endoscopy will receive full credit if returned in saleable condition
- Other merchandise returned within 30 days of invoice date in saleable condition will receive full credit less a 15% restocking charge
- Merchandise returned past 31 days in saleable condition will receive full credit less a 40% restocking charge provided there is at least 6 months left on the expiration date

Pricing and Specification Changes

As we are continually introducing products and improving our existing products, this Price List does not purport to be complete, and the equipment supplied may not agree in all details with the description and specifications shown in the catalog and / or price list. All product offerings, prices and specifications are subject to change without notice.

Notice

Indications, contraindications and instructions for use, as well as product warranty information, if applicable, can be found in the product labeling supplied with each device. Caution: Federal law restricts this device to sale by or on the order of a physician.

Warning Re: Single-Use Devices

Products marked as “single-use” are for single use only. Do not reuse, reprocess or resterilize single-use products. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and / or lead to device failure which in turn may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and / or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. **BSC assumes no liability with respect to single-use devices reused, reprocessed or resterilized and makes no warranties, express or implied, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose, with respect to any such instrument.**

©2022 Boston Scientific Corporation or its Affiliates. All rights reserved.

INDICATIONS FOR USE

The WallFlex™ Biliary RX Fully Covered Stent System RMV is indicated for use in the palliative treatment of biliary strictures produced by malignant neoplasms, relief of malignant biliary obstruction prior to surgery and for indwell up to 12 months in the treatment of benign biliary strictures secondary to chronic pancreatitis.

LIMITATIONS

The sale, distribution, and use of the device are restricted to prescription use in accordance with 21 CFR §801.109.

Contraindications:

- The WallFlex Biliary RX Fully Covered Stent System RMV should not be placed in strictures that cannot be dilated enough to pass the delivery system, in a perforated duct, or in very small intrahepatic ducts.
- The WallFlex Biliary RX Fully Covered Stent System RMV should not be used in patients for whom endoscopic techniques are contraindicated.

Warnings:

- The safety and effectiveness of the stent for benign stricture treatment has not been established for indwell periods exceeding 12 months.
- The WallFlex Biliary RX Fully Covered Stent System RMV is for single-use only.
- The safety and effectiveness of the WallFlex Biliary RX Fully Covered Stent System RMV for use in the vascular system has not been established.
- Clinical Study was not able to demonstrate that the use of this device is safe and effective in the treatment of benign biliary post liver transplantation anastomotic strictures and of benign biliary post abdominal surgery strictures.
- Testing of overlapped stents has not been conducted.
- The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

PLEASE REFER TO THE LABELING FOR A MORE COMPLETE LIST OF WARNINGS, PRECAUTIONS AND CONTRAINDICATIONS.



Section Table of Contents



Prev Page

103

Next Page



Home



For pricing information contact your Boston Scientific Sales Representative or Customer Service at **1.888.272.1001**

► Customer Ordering Information

Universal Product Number

Radial Jaw™ 4 Single-Use Jumbo Biopsy Forceps

Order Number	GTIN	Description
M005 13343	08714729761891	Radial Jaw 4 Jumbo Forceps
M005 13353	08714729761907	Radial Jaw 4 Jumbo Forceps
M005 13363	08714729761914	Radial Jaw 4 Jumbo Forceps
M005 13373	08714729761938	Radial Jaw 4 Jumbo Forceps
M005 13361	08714729648871	Radial Jaw 4 Jumbo Forceps
M005 13371	08714729648895	Radial Jaw 4 Jumbo Forceps
M005 13362	08714729648888	Radial Jaw 4 Jumbo Forceps


Bold type
is old order
number

For customer convenience, the Boston Scientific Endoscopy ordering system accepts the UPN (Universal Product Number) e.g., **M00513362** found on product labeling as well as Boston Scientific's old 4 digit order numbers e.g., **1336**.

Please use either one of these numbers when placing an order via phone, fax, internet, or Electronic Data Interchange (EDI).

For those customers who scan product today, Boston Scientific is moving to the GS1 GTIN barcode and labeling standards as part of our ongoing commitment to improve supply chain efficiency. For more information go to the website: www.bostonscientific.com/GS1.

Boston Scientific Endoscopy Customer Service Representatives will be happy to help you with any questions. We are eager to make the ordering procedure as simple and efficient as possible.



**Magnetic Resonance
Conditional¹**

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

All trademarks are the property of their respective owners.

**Boston
Scientific**

Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

Ordering Information
1.888.272.1001

©2022 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO – 1201902-AB