



Orca™

Sterile Single-Use Air/Water
and Suction Valves

One Less
Worry



Orca™

Sterile Single-Use Air/Water and Suction Valves

Recessed, matte finish finger grip

Helps provide greater gripping capability and minimizes finger slippage

Outfitted with a new sealing feature

Designed to enhance flow through the endoscope

Large, flat, color-coded finger perch

Designed for comfort and easy valve identification

Snap-fit boot

Designed to secure attachment of valve to cylinder to help prevent leakage

Silicone-free

In accordance with manufacturer's recommendations

Enclosed spring

Designed to prevent entrapment of fluid or debris which could compromise function

Reinforced one-way seal component

To support consistent insufflation

ORDERING INFORMATION

Order number	GTIN	Description	Scope compatibility	Packaging
SUV-617-50	00840253111012	Orca Sterile Single-Use Air/Water and Suction Valves	Olympus™ 160/180/190/1100 Series, EXALT™ Model D Single-Use Duodenoscope	Box 50

Please consult your Boston Scientific representative to identify which Infection Prevention procedure kits include the Orca Sterile Single-Use Air/Water and Suction Valves to support your endoscope reprocessing needs.

Boston Scientific
Advancing science for life™

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labeling supplied with each device. All trademarks are the property of their respective owners. All images owned by Boston Scientific.

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

© 2026 Boston Scientific Corporation or its affiliates. All rights reserved.

ENDO-470733-AD FEB 2026