

# Cold Snare


For endoscopic resection of gastrointestinal polyp tissue

Available in three shapes and two loop sizes, the Cold Snare line offers 360-degree rotation and a thin braided wire.

**Oval** 2 sizes available

**Shield** 2 sizes available

**Drop** 2 sizes available



## Rotatable Cold Snares for a variety of procedural needs

- **Multiple shapes and sizes available**
- Compatible with endoscopes that have a **minimum working channel of 2.8 mm**
- **230 cm** working length

\*Five users performed a simulated tissue removal step on an anatomical model and classified the tissue removal process as easy without tissue tearing. Bench Test results may not necessarily be indicative of clinical performance. The testing was performed by or on behalf of Boston Scientific. Data on file.



## Ordering information

Order number	GTIN	Size (mm)	Shape	Wire material	Box
M00591201	26926839812370	10	Oval	Stainless Steel	10
M00591202	26926839812387	10	Oval	Stainless Steel	20
M00591211	26926839812400	15	Oval	Stainless Steel	10
M00591212	26926839812417	15	Oval	Stainless Steel	20
M00591221	26926839812431	10	Shield	Stainless Steel	10
M00591222	26926839812448	10	Shield	Stainless Steel	20
M00591231	26926839812462	15	Shield	Stainless Steel	10
M00591232	26926839812479	15	Shield	Stainless Steel	20
M00591241	26926839812493	10	Drop	Nitinol	10
M00591242	26926839812509	10	Drop	Nitinol	20
M00591251	26926839812523	15	Drop	Nitinol	10
M00591252	26926839812530	15	Drop	Nitinol	20



Ready to order?  
Scan the QR code.

**Boston Scientific**  
Advancing science for life™

Boston Scientific data on file.

This product is manufactured by Anrei Medical and distributed by Boston Scientific.

Indications for use can be found in the product labeling supplied with each device.

All trademarks are the property of their respective owners. All photographs are owned by Boston Scientific.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

© 2025 Boston Scientific Corporation or its affiliates. All rights reserved.

Boston Scientific Corporation  
300 Boston Scientific Way  
Marlborough, MA 01752-1234  
©2025 Boston Scientific Corporation or its affiliates.  
All rights reserved.  
ENDO-2387209-AA