


Boa®

Single-Use Polypectomy Snare

**Boston
Scientific**
Advancing science for life™



**Control
Through Innovation**

Boa[®]

Polypectomy Snare

FEATURES AND BENEFITS

The OmniLoop™ Advantage

- 3-in-1 variable loop sizing (10, 20 & 30 mm)
- Designed to provide flat loop deployment and retraction while maintaining shape through multiple deployments
- Maximizes loop diameter while minimizing length

Innovative Design

- Ergonomic handle is designed to provide comfort and tactile feel while reducing hand strain
- Integrated depth stop

Teardrop loop configuration

Depth stop

10 mm, 20 mm & 30 mm markings

Ordering Information

Order Number	Description	Loop Width	Wire Type	Sheath O.D. (mm)	Working Length (cm)	Packaging
DSO-013-10	Polypectomy Snare	10mm, 20 mm & 30 mm	Stiff, Braided	2.3 mm	230 cm	Box 10

Design allows for consistent snare actuation in tortuous scope positions

All trademarks are the property of their respective owners.
Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.
CAUTION: The law restricts these devices to sale by or on the order of a physician.
Information for use only in countries with applicable health authority registrations.
Material not intended for use in France.

Boston Scientific
Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752

© 2017 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-453506-AA March 2017