

# EndoVive™

Enteral Access Devices with ENFit™ connection

## EndoVive Standard PEG Kit Percutaneous Endoscopic Gastrostomy Tube ENFit

ENFit Kit	ENFit Description	Units
M00509081	EndoVive StdPEG Pull w/ENFit 20Fr	Box 2
M00509101	EndoVive StdPEG Pull w/ENFit 24Fr	Box 2
M00509091	EndoVive StdPEG Push w/ENFit 20Fr	Box 2
M00509111	EndoVive StdPEG Push w/ENFit 24Fr	Box 2

## EndoVive PEG Accessories

UPN	Description	Units
M00501261	EndoVive Y-Port Adapt w/ENFit 20Fr	Box 10
M00501271	EndoVive Y-Port Adapt w/ENFit 24Fr	Box 10
M00580731	20Fr Round Bolsters	Box 10
M00580781	20Fr and 24Fr C-Clamps	Box 10
M00509241	ENFit Feed Port Connector	Box 10
M00509221	ENFit Transition Connector Enteral Distal End	Box 10

## EndoVive Gastrostomy Tube Replacement Gastrostomy Feeding Tube ENFit

ENFit Kit	ENFit Description	Units
M00509881	12Fr Repl G-Tube, Straight	Box 5
M00509891	14Fr Repl G-Tube, Straight	Box 5
M00509901	16Fr Repl G-Tube, Straight	Box 5
M00509911	18Fr Repl G-Tube, Straight	Box 5
M00509921	20Fr Repl G-Tube, Straight	Box 5
M00509931	22Fr Repl G-Tube, Straight	Box 5
M00509941	24Fr Repl G-Tube, Straight	Box 5
M00509961	12Fr Repl G-Tube, Right Angle	Box 5
M00509971	14Fr Repl G-Tube, Right Angle	Box 5
M00509981	16Fr Repl G-Tube, Right Angle	Box 5
M00509991	18Fr Repl G-Tube, Right Angle	Box 5
M00510001	20Fr Repl G-Tube, Right Angle	Box 5
M00510011	22Fr Repl G-Tube, Right Angle	Box 5
M00510021	24Fr Repl G-Tube, Right Angle	Box 5

**Boston  
Scientific**  
Advancing science for life™

# EndoVive™

Enteral Access Devices with ENFit™ connection



Product Guide

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Information contained herein is for distribution outside the U.S. only. Illustrations for information purposes – not indicative of actual size or clinical outcome.

Product available in the European Economic Area (EEA) only. Please check availability with your local sales representative or customer service.

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. Please check availability with your local sales representative or customer service.

ENDO-430706-AA OCT2016 Produced in Germany by medicalvision.

**Boston  
Scientific**  
Advancing science for life™

[www.bostonscientific.eu](http://www.bostonscientific.eu)

© 2016 Boston Scientific Corporation  
or its affiliates. All rights reserved.



In the interest of patient safety, Boston Scientific will be offering Enteral Feeding products that comply with the ISO 80369-3 standard to reduce the risk of misconnection via an enteral-specific connector design. The connection is referred to as ENFit™, which employs a threaded twist connection exclusive to enteral feeding to help prevent misconnections with other small bore connectors.



### EndoVive Standard PEG Kit – Percutaneous Endoscopic Gastrostomy Tube ENFit

### EndoVive Gastrostomy Tube – Replacement Gastrostomy Feeding Tube ENFit

#### Product Features

##### Silicone Dome-Bolstered PEG Tube

- Facilitates traction removal

##### Less than or equal to 2.00 mm O.D. x 240 cm Retrieval Snare

- Traverses scopes with a 2.0 mm or larger working channel
- Universal-sized snare, intended to eliminate the cost of a specialty snare

##### Double-Tier Tray Packaging

- Kit organisation promotes procedural efficiency

##### External Bolsters

- Two external bolster options: accommodates physician and patient needs: Round; T-Bar

##### Multiple Sizes and Accessories

- Two French sizes for broad patient coverage
- Separate feeding and medication ports
- C-Clamp is provided to occlude tubing when necessary



.....  
Silicone  
Dome-Bolstered  
PEG Tube

KIT CONTENTS – TRAY 1	
(1) Fenestrated Drape	(4) Drainage Sponges
(4) Gauze Pads	(1) ENFit Y port
(1) Scalpel	(1) C-clamp
(1) Curved Hemostat	(1) Round & T-bar Bolsters
(1) Scissors	(1) Seldinger Needle (Push kit)
(1) Angiocath™ Needle (Trocar)	

KIT CONTENTS – TRAY 2
(1) PEG tube Pull or Push
(1) Insertion Wire (Pull kit) or Guidewire (Push kit)
(1) Retrieval Snare



#### Product Features

##### Low Profile “Apple” Shaped Balloon

- Designed for proper retention during use and easy decompression before device removal
- Balloon inflation varying depending on the size: up to 5 ml for 12 to 16Fr balloons up to 20 ml for 18 to 24Fr balloons.

##### Standard Profile Gastrostomy with Optional Right Angle External Bolster

- Designed to minimise the potential tube kinking when it's taped down
- Designed to minimise premature pull out due to shape of the balloon
- Separate medication and feeding ports
- Clear shaft for internal lumen visualisation

KIT CONTENTS
(1) Standard Balloon Replacement with ENFit connection with Straight or Right Angle External Bolster
(2) 4”x 4” Gauze Pads
(1) 6 ml Syringe
(1) Lubricating Jelly



.....  
Right Angle  
Gastrostomy  
Tube

.....  
Straight  
Gastrostomy Tube