

Rescue™

Retrieval Devices

**Boston
Scientific**
Advancing science for life™



Predictable Performance for Unpredictable Circumstances



Retrieval Net

Retrieval Forceps

Ordering Information



RescueNet™ Retrieval Device

Retrieval Net


Retrieval Forceps

Ordering Information

Designed for optimal visibility and control with excellent mesh strength for quick and efficient foreign body removal.


OmniLoop™ Technology

Enables smooth, flat loop deployment and retraction to minimize scope repositioning





ProMesh™ Netting

- Ultra-strong, porous polyester mesh material designed to minimize tearing, fraying or shearing
- Conforms to the foreign object enhancing ease of capture and removal



Optical Blue Loop Sheath

Blue color enhances visibility in a fluid field.



Rescue™ Forceps



Retrieval Net

Retrieval Forceps

Ordering Information

A portfolio of forceps to provide the right solution for your needs.

Rescue Combo Rat Tooth / Alligator Grasping Forceps


- Combination rat tooth and alligator jaw forceps
- Versatile design with a hybrid jaw to maximize both grasping strength and object adherence

Rescue Alligator Grasping Forceps

Rescue Rat Tooth Grasping Forceps

Insulated sheath


Designed to be atraumatic while maintaining passability and flexibility



Rescue™ Retrieval Devices

[Retrieval Net](#)[Retrieval Forceps](#)[Ordering Information](#)

RescueNet™ Retrieval Device



Order Number	Description	Net Size	Sheath O.D.	Minimal Working Channel	Working Length	Packaging
DGN-538-5	RescueNet Retrieval Net	3cm x 5.5cm	2.5mm	2.8mm	230cm	Box 5

Rescue Forceps



Order Number	Description	Jaw Opening	Sheath O.D.	Working Length	Packaging
DGA-274-5	Rescue Alligator Jaw Grasping Forceps	8mm	2.4mm	230cm	Box 5
DGA-275-5	Rescue Alligator Jaw Grasping Forceps (Long)	9mm	2.4mm	230cm	Box 5
DGR-276-5	Rescue Rat Tooth Jaw Grasping Forceps	7mm	2.4mm	230cm	Box 5
DGR-277-5	Rescue Rat Tooth Jaw Grasping Forceps (Pediatric)	4.5mm	1.8mm	180cm	Box 5
DGC-383-5	Rescue Combo Grasping Forceps	8mm	2.4mm	230cm	Box 5

Rescue™ Retrieval Devices

[Retrieval Net](#)[Retrieval Forceps](#)[Ordering Information](#)

**Boston
Scientific**
Advancing science for life™

All trademarks are the property of their respective owners.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Caution: The law restricts these devices to sale by or on the order of a physician. Information for use only in countries with applicable health authority registrations. Material not intended for use in France.

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com/gastro

©2017 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-465201-AA April 2017