



Rescue™ Single-Use Grasping Forceps

Prescriptive Information

Refer to the device directions for use for complete instructions on device use.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Warning

Use same care as for any other surgical device to avoid punctures or cuts to hands, gloves, etc.

For single use only. Do not reuse, reprocess or re-sterilize. Reuse, reprocessing, or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death.

Inspect the package containing the forceps. If any damage is noticed, do not use the forceps as sterility may be compromised.

DO NOT USE if any abnormalities are found or if the forceps are not operating properly.

Intended Use/Indications for Use

Grasping forceps are intended to be used to grasp tissue and/or retrieve foreign bodies, excised tissue, and stents from within the gastrointestinal tract.

Precautions

None Known

Adverse Events

None Known

Contraindications

None known within the scope of the indicated use.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.

2024 Copyright © Boston Scientific Corporation or its affiliates. All rights reserved. (copyright statement only required if not otherwise on material)