



Pulmonary Endoscopy

Alair™ Bronchial Thermoplasty System

Alair Bronchial Thermoplasty System.....67

Airway Stents

Dynamic™ (Y) Stent67

Polyflex™ Single-Use Self-Expanding Silicone Airway Stent.....69

Ultraflex™ Single-Use Uncovered Tracheobronchial Stent System –
Distal Release.....68

Ultraflex Single-Use Uncovered Tracheobronchial Stent System –
Proximal Release68

Biopsy Forceps

Radial Jaw™ 4 Single-Use Pulmonary Biopsy Forceps.....70

Transbronchial Aspiration Needles

eXcelon™ Single-Use Transbronchial Aspiration Needle.....70

Cytology Brushes

Cellebrity™ Single-Use Endoscopic Cytology Brush – Bronchoscope.....71



Airway Retrieval Baskets

Zero Tip™ Single-Use Airway Retrieval Basket71

Controlled Radial Expansion
Pulmonary Balloon Dilators

CRE™ Single-Use Pulmonary Balloon Dilator72

Pulmonary Guidewires


Pulmonary Jagwire™ Pulmonary Guidewire.....72

Inflation System

Alliance™ II Inflation System.....72

Pulmonary

Alair™ Bronchial Thermoplasty System



Alair Bronchial Thermoplasty System*



Order Number	Model Number	Description	Packaging
M005 ATS20000	Alair ATS 200	Alair Bronchial Thermoplasty Controller	Each
M005 ATS25010	Alair ATS 2-5	Alair Bronchial Thermoplasty Catheter**	Each
M005 ATS20110	Alair ATS 201	Alair RF Controller Accessory Kit	Each


*Service Agreement available with 1, 2, 3, and 4 year terms.

****Note:** Initial stocking order requires a minimum order of 6 catheters (covering the complete treatment of 2 patients).

Alair Bronchial Thermoplasty System Service Agreements

UPN	Description
SRVM005 ATS1001	ENDO POS 1 Year Service Alair
SRVM005 ATS1002	ENDO POS 2 Year Service Alair
SRVM005 ATS1003	ENDO POS 3 Year Service Alair
SRVM005 ATS1004	ENDO POS 4 Year Service Alair

Airway Stents




Dynamic™ (Y) Stent



Order Number	Tracheal Width (I) (mm)	Bronchial Width (I) (mm)	Tracheal Length (mm)	Bronchial Lengths (R / L) (mm)	Packaging
M005 70670	11	8	110	25 / 40	Each
M005 70680	13	10	110	25 / 40	Each
M005 70690	15	12	110	25 / 40	Each

Pulmonary

Airway Stents (Continued)



Ultraflex™ Single-Use Uncovered Tracheobronchial Stent System – Distal Release*



Order Number	Description	Expanded Stent O.D. (mm)	Expanded Stent Length (mm)	Compressed Stent O.D. (Fr)**	Tip Max O.D. (mm)	Packaging
M00576300	Uncovered Distal Release	10	30	15	4.1	Each
M00576310	Uncovered Distal Release	12	30	17	4.1	Each
M00576320	Uncovered Distal Release	14	30	18	4.1	Each
M00576330	Uncovered Distal Release	16	40	20	5.3	Each
M00576340	Uncovered Distal Release	16	60	20	5.3	Each
M00576350	Uncovered Distal Release	18	40	20	5.3	Each
M00576360	Uncovered Distal Release	18	60	20	5.3	Each
M00576370	Uncovered Distal Release	20	40	20	5.3	Each
M00576380	Uncovered Distal Release	20	60	20	5.3	Each

*Recommended Guidewire .035" Pulmonary Jagwire™ Pulmonary Guidewire, Order #1517, (see page 72).

Delivery System working length is 180cm.

**Data on File at Boston Scientific.



Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release*



Order Number	Description	Expanded Stent O.D. (mm)	Expanded Stent Length (mm)	Compressed Stent O.D. (Fr)**	Tip Max O.D. (mm)	Packaging
M00576410	Uncovered Proximal Release	8	20	15	4.1	Each
M00576420	Uncovered Proximal Release	8	40	15	4.1	Each
M00576430	Uncovered Proximal Release	10	20	15	4.1	Each
M00576440	Uncovered Proximal Release	10	40	15	4.1	Each
M00576450	Uncovered Proximal Release	12	20	17	4.1	Each
M00576460	Uncovered Proximal Release	12	40	17	4.1	Each
M00576470	Uncovered Proximal Release	14	20	18	4.1	Each
M00576480	Uncovered Proximal Release	14	40	18	4.1	Each
M00576490	Uncovered Proximal Release	14	60	18	4.1	Each


*Recommended Guidewire .035" Pulmonary Jagwire™ Pulmonary Guidewire, Order #1517, (see page 72).

Delivery System working length is 180cm.

**Data on File at Boston Scientific.

Pulmonary

Airway Stents (Continued)




Polyflex™ Single-Use Self-Expanding Silicone Airway Stent



Order Number	Stent ID (mm)	Length (mm)	Delivery System Diameter (mm)	Packaging
M00570000	8	20	7	Each
M00570010	8	30	7	Each
M00570020	10	20	8	Each
M00570030	10	30	8	Each
M00570040	10	40	8	Each
M00570050	10	50	8	Each
M00570060	12	20	9	Each
M00570070	12	30	9	Each
M00570080	12	40	9	Each
M00570090	12	50	9	Each
M00570100	14	20	9	Each
M00570110	14	30	9	Each
M00570120	14	40	9	Each
M00570130	14	50	9	Each
M00570140	14	60	9	Each
M00570150	16	30	10	Each
M00570160	16	40	10	Each
M00570170	16	50	10	Each
M00570180	16	60	10	Each
M00570190	16	70	10	Each
M00570200	18	30	11	Each
M00570210	18	40	11	Each
M00570220	18	50	11	Each
M00570230	18	60	11	Each
M00570240	18	70	11	Each
M00570250	18	80	11	Each
M00570260	20	40	12	Each
M00570270	20	50	12	Each
M00570280	20	60	12	Each
M00570290	20	70	12	Each
M00570300	20	80	12	Each
M00570310	22	50	13	Each
M00570320	22	60	13	Each
M00570330	22	80	13	Each

Pulmonary

Biopsy Forceps




Radial Jaw™ 4 Single-Use Pulmonary Biopsy Forceps



Order Number	Description	Differentiator	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Packaging
M005 15181	Single-Use Biopsy Forceps	Pulmonary Standard Capacity	1.8	100	2.0	Box 5
M005 15182	Single-Use Biopsy Forceps	Pulmonary Standard Capacity	1.8	100	2.0	Box 20
M005 15191	Single-Use Biopsy Forceps	Pulmonary Standard Capacity w/ Needle	1.8	100	2.0	Box 5
M005 15192	Single-Use Biopsy Forceps	Pulmonary Standard Capacity w/ Needle	1.8	100	2.0	Box 20
M005 15201	Single-Use Biopsy Forceps	Pulmonary Large Capacity	2.2	100	2.8	Box 5
M005 15202	Single-Use Biopsy Forceps	Pulmonary Large Capacity	2.2	100	2.8	Box 20

Transbronchial Aspiration Needles



eXcelon™ Single-Use Transbronchial Aspiration Needle*



Order Number	Description	Needle Gauge	Needle Length (mm)	Catheter Length (cm)	Sheath O.D. (mm)	Packaging
M005 64101	eXcelon Transbronchial Aspiration Needle	19	15	130	1.8	Box 5
M005 64111	eXcelon Transbronchial Aspiration Needle	20	15	130	1.8	Box 5
M005 64121	eXcelon Transbronchial Aspiration Needle	21	15	130	1.8	Box 5

*Packaging also includes a 20ml syringe.

Pulmonary


Cytology Brushes

Celebrity™ Single-Use Endoscopic Cytology Brush – Bronchoscope



Order Number	Required Working Channel (mm)	Bristle O.D. (mm)	Sheath Length (cm)	Packaging
M005 1600 1	2.0	1.0	140	Box 10
M005 1601 1	2.0	1.5	140	Box 10
M005 1607 1	2.0	1.9	100	Box 10
M005 1615 1	2.0	1.9	150	Box 10

Airway Retrieval Baskets



Zero Tip™ Single-Use Airway Retrieval Basket



Order Number	Description	Sheath O.D. (mm)	Working Length (cm)	Basket Opening (mm)	Packaging
M005 13200	Zero Tip Airway Retrieval Basket	0.8	120	12	Each
M005 13210	Zero Tip Airway Retrieval Basket	1.0	120	16	Each

Pulmonary

Controlled Radial Expansion Pulmonary Balloon Dilators



CRE™ Single-Use Pulmonary Balloon Dilator*



Order Number	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Catheter Length (cm)	Packaging
	(mm)	(Fr)		(ATM)	(kPa)		
M005 50300	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	75	Each
M005 50310	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	75	Each
M005 50320	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	75	Each
M005 50330	8-9-10	24-27-30	3.0	3-5.5-9	304-557-912	75	Each
M005 50340	10-11-12	30-33-36	3.0	3-5-8	304-507-811	75	Each
M005 50350	12-13.5-15	36-40.5-45	3.0	3-4.5-8	304-456-811	75	Each

*Recommended Inflation Device: Alliance™ II Integrated Inflation System, Handle Order #5062, Sterile Syringe Order #5060 (see below).


Pulmonary Guidewires

Pulmonary Jagwire™ Pulmonary Guidewire



Order Number	Description	O.D. (in)	Length (cm)	Packaging
M005 15171	Pulmonary Jagwire Pulmonary Guidewire	.035	180	Box 2

Inflation System



Alliance II Inflation System

Order Number	Description	Packaging
M005 50620	Inflation Handle	Each
M005 50601	 Single-Use Syringe / Gauge Assembly	Box 5

[illegible]

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.