


A difference you can feel™.

Boston
Scientific



Extractor™ Pro RX

Retrieval Balloon



Defining tomorrow, today
in Endoscopy.

Extractor™ Pro RX


Triple-lumen Retrieval Balloon

Boston Scientific has earned the confidence and allegiance of physicians around the world who have adopted the RX Biliary System™ – the first short-wire delivery system – for everyday use. The entire family of RX Biliary System devices is specifically engineered to increase efficiency and control during ERCP procedures.



Enhanced Balloon Design

- New and improved balloon is designed to improve strength and reliability.
- Each of the three balloon sizes can be inflated to two distinct diameters: 9–12mm; 12–15mm; 15–18mm to better accommodate anatomical variations and help obviate the need to use an additional balloon and change catheters if an alternative balloon diameter is required.
- “Squared shoulder” balloon is designed to facilitate stone removal.



Procedural Efficiency and Control

- The RX Biliary System locking device allows physicians to maintain guidewire position and helps facilitate multiple sweeps of the duct during stone retrieval.
- Option to front load or back load the guidewire.
- Injection-above or injection-below options to support physician preference and facilitate procedural needs.
- Designed for excellent injection flow rate.

Ordering Information

Extractor Pro RX Triple-lumen Retrieval Balloon

Order Number	Balloon O.D. (mm)	Injection	Catheter O.D. (F/mm)
M00547000	9–12	Above	7–6 / 2.33–2
M00547010	12–15	Above	7–6 / 2.33–2
M00547020	15–18	Above	7–6 / 2.33–2
M00547030	9–12	Below	7–6 / 2.33–2
M00547040	12–15	Below	7–6 / 2.33–2
M00547050	15–18	Below	7–6 / 2.33–2

Packaged one per box

Recommended guidewire: .035" (0.89mm) Dreamwire™ High-Performance Guidewire. System length is 200cm.

European Distribution Centre –

The Netherlands

T: +31 45 54 67 700 F: +31 45 54 67 800

European Headquarters –

Paris

T: +33 1 57 66 80 00 F: +33 1 57 66 84 99

Argentina (Freephone)

T: 0800 555 2678 F: +5411 4896 8517

Australia/New Zealand

T: +61 2 8063 8100 F: +61 2 9330 1404

Austria

T: +43 1 60 810 F: +43 1 60 810 60

Belgium (Freephone)

T: 0800 94 494 F: 0800 93 343

Brazil

T: +55 11 5502 8500 F: +55 11 5103 2212

Canada

T: +1 888 359 9691 F: +1 888 575 73965

Chile

T: +562 445 4904 F: +562 445 4915

China – Beijing

T: +86 10 8525 1588 F: +86 10 8525 1566

China – Guangzhou

T: +86 20 8767 9791 F: +86 20 8767 9789

China – Shanghai

T: +86 21 6141 5959 F: +86 21 6141 5900

Colombia

T: +57 1 629 5045 F: +57 1 612 4761

Czech Republic

T: +420 296 331 901 F: +420 296 331 935

Denmark (Freephone)

T: 80 30 80 02 F: 80 30 80 05

Eire (Freephone)

T: 1800 882 969 F: 1800 882 968

Finland

T: +358 20 762 8882 F: +358 20 762 8883

France

T: +33 1 39 30 49 00 F: +33 1 39 30 49 01

Germany (Freephone)

T: 08000 723300 F: 08000 723319

Greece

T: +30 210 9542 300 F: +30 210 9542 310

Hong Kong

T: +852 2960 7100 F: +852 2563 5276

Hungary

T: +36 1 456 30 40 F: +36 1 456 30 41

India – Bangalore

T: +91 80 2212 4928/9 F: +91 80 2207 5153

India – Chennai

T: +91 44 2220 1879 F: +91 44 2220 0716

India – Delhi

T: +91 11 4243 2222 F: +91 11 2610 0808

India – Mumbai

T: +91 22 4030 9165 F: +91 22 4040 9199

Italy

T: +39 010 60 60 1 F: +39 010 60 60 200

Korea

T: +82 2 3476 2121 F: +82 2 3476 1776

Mexico

T: +52 55 5687 63 90 F: +52 55 5687 62 28

Middle East/Gulf/North Africa

T: +961 1 805 410 F: +961 1 805 445

The Netherlands

T: +31 30 602 55 44 F: +31 30 602 55 05

Norway (Freephone)

T: 800 104 04 F: 800 101 90

Poland

T: +48 22 435 14 14 F: +48 22 435 14 10

Portugal

T: +351 1 381 25 40 F: +351 21 381 25 58

South Africa

T: +27 11 840 8600 F: +27 11 463 6077

South East Asia – Malaysia

T: +60 3 2283 3813 F: +60 3 2284 3813

South East Asia – Philippines

T: +63 2 687 6994 F: +63 2 687 3047

South East Asia – Singapore

T: +65 6418 8888 F: +65 6418 8899

South East Asia – Thailand

T: +66 2 6380 100 F: +66 2 6380 400

South East Europe

T: +30 210 9542 300 F: +30 210 9542 310

Spain

T: +34 901 11 12 15 F: +34 91 319 50 03

Sweden

T: +46 42 25 69 00 F: +46 42 25 69 69

Switzerland

T: +41 32 626 57 00 F: +41 32 626 57 01

Taiwan

T: +886 2 2747 7278 F: +886 2 2747 7270

Turkey

T: +90 216 464 36 66 F: +90 216 464 36 67

UK (Freephone)

T: 0844 800 4512 F: 0844 800 4513

Uruguay

T: +59 82 900 6212 F: +59 82 900 6212

Venezuela

T: +58 212 959 6275 F: +58 212 959 5328

Boston Scientific

Defining tomorrow, today.™

www.bostonscientific-international.com