


NovaGold™

High Performance Guidewire

**Boston
Scientific**
Advancing science for life™



DELIVERING DURABILITY AND VISIBILITY IN AN .018" GUIDEWIRE

*Optimizing wireguided cannulation for the
treatment of pancreatico-biliary diseases*

ABOUT

FEATURES


ACCESSORIES

ORDERING INFO



About the NovaGold Guidewire

The NovaGold High Performance Guidewire is intended for use in selective cannulation of the biliary ducts, including the common bile, pancreatic, cystic, right and left hepatic ducts, and to aid in the placement of diagnostic and therapeutic devices during endoscopic procedures.



Available in both short-wire and long-wire lengths to accommodate physician preference.


NovaGold™

High Performance Guidewire

Boston
Scientific

Device Features

1 of 3

A detailed image of the NovaGold guidewire. It features a long, thin, gold-colored wire that is inserted into a white, flexible, multi-stranded catheter. The catheter has several colored bands: a blue band near the top, a green band further down, and a brown band near the bottom. The wire is shown curving through the catheter, demonstrating its flexibility and pushability.

The **Triton™ Alloy Core** provides a durable wire for working with accessory devices. The stiffness of the wire makes it pushable yet flexible for loading and rail support through the scope to facilitate device exchange.

ABOUT

FEATURES

ACCESSORIES

ORDERING INFO



NovaGold™

High Performance Guidewire

Boston
Scientific

Device Features

< 2 of 3 >

The **6cm Atraumatic Tip** provides flexibility needed to access the pancreatico-biliary system.

The **Radiopaque Guidewire Tip** is designed to help confirm location within the duct under fluoroscopy.

ABOUT

FEATURES

ACCESSORIES

ORDERING INFO



NovaGold™


High Performance Guidewire

Boston
Scientific

Device Features

< 3 of 3

The **NovaGold Guidewire's Steerability** enables maneuverability for cannulation and working in tortuous anatomy.



ABOUT

FEATURES

ACCESSORIES

ORDERING INFO



NovaGold™

High Performance Guidewire

Boston
Scientific

Device Accessories

Boston Scientific's wide range of endoscopic access and therapeutic devices provides physicians with options to implement the most effective course of treatment.

Through innovation and continuous improvement, Boston Scientific offers a range of accessories for the diagnosis and treatment of pancreatico-biliary stone and stricture management.

Visit www.bostonscientific.com for treatment options, product information and links to procedural videos.



ABOUT

FEATURES

ACCESSORIES

ORDERING INFO



NovaGold™

High Performance Guidewire

Ordering Information

Product Code	Guidewire Diameter (inches)	Guidewire Length (cm)	Tip	Packaging
M00552001	.018	260	Straight (Shapeable)	Box of 2
M00552011	.018	480	Straight (Shapeable)	Box of 2

**Boston
Scientific**
Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752
www.bostonscientific.com

Ordering Information
1.888.272.1001

©2015 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-207914-AB May 2015

ENDO-428801-AA October 2016

All trademarks are the property of their respective owners.
Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.
CAUTION: Federal (U.S) law restricts this device to sale by or on the order of a physician.
Distributed by Boston Scientific Corporation, Marlborough, Massachusetts; manufactured by NeoMetrics Inc., Plymouth, Minnesota.

ABOUT

FEATURES

ACCESSORIES

ORDERING INFO

