

Radial Jaw[®] 4 Hot Biopsy Forceps

Boston
Scientific

Single-Use Biopsy Forceps



Samples
Like Never Before™



Endoscopy

Radial Jaw® 4 Hot Biopsy Forceps

Single-Use Biopsy Forceps

*New jaw configuration, now available
in our family of market-leading single-use
biopsy forceps.*

NEW JAW CONFIGURATION

Designed to Provide:

- Tissue specimens for excellent sample handling, preparation and accurate histological diagnosis
- Consistent sample retention for precise sample identification
- Enhanced passability through tortuous anatomy

IMPROVED MICRO-MESH TOOTH DESIGN

- Engineered for a clean, precise bite

POSITIONING MARKERS

- Engineered to assist in the insertion and withdrawal of forceps to enhance procedural safety and efficiency



STREAMLINED CATHETER

Designed to Provide:

- The right balance of columnar strength and flexibility for excellent pushability and control during scope passage

FENESTRATION HOLE

- Designed to provide consistent sample retention and removal

SINGLE-USE

- Designed to eliminate the risk of transmitting patient-to-patient disease

NEW DISTAL END TUBE

- Improved visibility
- Prevents inadvertent lodging of the cap in the scope working channel



*Boston Scientific, setting new
standards for endoscopic biopsy.*

ORDERING INFORMATION

Order Number	Product Description	Jaw OD (mm)	Working Length (cm)	Required Working Channel (mm)	Units
M00515031	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	≥2.8	Box 5
M00515032	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	≥2.8	Box 20
M00515033	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	≥2.8	Box 40

Boston Scientific

Defining tomorrow, today.™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientific.com/endoscopy

Ordering Information
1.800.225.3226

© 2012 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-91502-AB October 2012

Radial Jaw is a registered trademark of Boston Scientific Corporation and its affiliates.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.