

Radial Jaw[®] 4

Single-Use Biopsy Forceps

Boston
Scientific

Jumbo

Large Capacity

Samples
Like Never Before

Radial Jaw® 4

Single-Use Biopsy Forceps

Introducing Radial Jaw 4

The next generation of our market-leading single-use biopsy forceps

Featuring...

New Jaw Configuration Designed to Provide:

- Large tissue specimens for excellent sample handling, preparation and accurate histological diagnosis

Jumbo: More than 1.7 X the Jaw Volume of Radial Jaw® 3 Large Capacity Forceps

Large Capacity: More than 1.2 X the Jaw Volume of Radial Jaw 3 Large Capacity Forceps

- Consistent sample retention for precise sample identification
- Enhanced passability through tortuous anatomy

Improved Micro-Mesh Tooth Design

- Engineered for a clean, precise bite

New Streamlined Catheter Designed to Provide:

- The right balance of columnar strength and flexibility for excellent pushability and control during scope passage

New Positioning Markers

- Engineered to assist in the insertion and withdrawal of forceps to enhance procedural safety and efficiency

Single-Use

- Designed to eliminate the risk of transmitting patient-to-patient disease



Boston Scientific, setting new standards for endoscopic biopsy.

ORDERING INFORMATION

Order Number	Description	Working Length (cm)	Required Working Channel (mm)	Units
Radial Jaw® 4 Single-Use Jumbo Biopsy Forceps				
M00513341	Radial Jaw 4 Jumbo	160	3.2	Box 20
M00513351	Radial Jaw 4 Jumbo w/Needle	160	3.2	Box 20
M00513361	Radial Jaw 4 Jumbo	240	3.2	Box 20
M00513371	Radial Jaw 4 Jumbo w/Needle	240	3.2	Box 20
M00513342	Radial Jaw 4 Jumbo	160	3.2	Box 40
M00513352	Radial Jaw 4 Jumbo w/Needle	160	3.2	Box 40
M00513362	Radial Jaw 4 Jumbo	240	3.2	Box 40
M00513372	Radial Jaw 4 Jumbo w/Needle	240	3.2	Box 40
Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps				
M00513301	Radial Jaw 4 Large Capacity	160	2.8	Box 20
M00513311	Radial Jaw 4 Large Capacity w/Needle	160	2.8	Box 20
M00513321	Radial Jaw 4 Large Capacity	240	2.8	Box 20
M00513331	Radial Jaw 4 Large Capacity w/Needle	240	2.8	Box 20
M00513302	Radial Jaw 4 Large Capacity	160	2.8	Box 40
M00513312	Radial Jaw 4 Large Capacity w/Needle	160	2.8	Box 40
M00513322	Radial Jaw 4 Large Capacity	240	2.8	Box 40
M00513332	Radial Jaw 4 Large Capacity w/Needle	240	2.8	Box 40

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.
CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician. Data on file

Boston Scientific

Delivering what's next.™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientific.com

Ordering Information
1.800.225.3226

© 2006 Boston Scientific Corporation
or its affiliates. All rights reserved.

DVG1800 5M 9/06

ENDO-150605-AA September 2013