

# Radial Jaw® 4 Standard Capacity

Boston  
Scientific

Single-Use Biopsy Forceps

Samples  
Like Never Before™



# Radial Jaw® 4 Standard Capacity

Single-Use Biopsy Forceps

*New jaw configuration, now available in our family of market-leading single-use biopsy forceps.*

## JAW CONFIGURATION SIZE

**Designed to Provide:**

- Tissue specimens for excellent sample handling, preparation and accurate histological diagnosis
- Consistent sample retention for precise sample identification
- Enhanced passability through tortuous anatomy

## MICRO-MESH TOOTH DESIGN

- Engineered for a clean, precise bite

## POSITIONING MARKERS

- Engineered to assist in the insertion and withdrawal of forceps to enhance procedural safety and efficiency



## STREAMLINED CATHETER

**Designed to Provide:**

- The right balance of columnar strength and flexibility for excellent pushability and control during scope passage

## FENESTRATION HOLE

- Designed to provide consistent sample retention and removal

## SINGLE-USE

- Designed to eliminate the risk of transmitting patient-to-patient disease

## NEW DISTAL END TUBE

- Improved visibility
- Prevents inadvertent lodging of the tube in the scope working channel



*Boston Scientific, setting new standards for endoscopic biopsy.*

## ORDERING INFORMATION

Order Number	Product Description	Working Length (cm)	Required Working Channel (mm)	Units
M00513383	Radial Jaw 4 Standard Capacity	160	≥2.8	Box 5
M00513393	Radial Jaw 4 Standard Capacity w/Needle	160	≥2.8	Box 5
M00513403	Radial Jaw 4 Standard Capacity	240	≥2.8	Box 5
M00513413	Radial Jaw 4 Standard Capacity w/Needle	240	≥2.8	Box 5
M00513381	Radial Jaw 4 Standard Capacity	160	≥2.8	Box 20
M00513391	Radial Jaw 4 Standard Capacity w/Needle	160	≥2.8	Box 20
M00513401	Radial Jaw 4 Standard Capacity	240	≥2.8	Box 20
M00513411	Radial Jaw 4 Standard Capacity w/Needle	240	≥2.8	Box 20
M00513382	Radial Jaw 4 Standard Capacity	160	≥2.8	Box 40
M00513392	Radial Jaw 4 Standard Capacity w/Needle	160	≥2.8	Box 40
M00513402	Radial Jaw 4 Standard Capacity	240	≥2.8	Box 40
M00513412	Radial Jaw 4 Standard Capacity w/Needle	240	≥2.8	Box 40

Radial Jaw is a registered trademark of Boston Scientific Corporation and its affiliates.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.

# Boston Scientific

*Defining tomorrow, today.™*

Boston Scientific Corporation  
One Boston Scientific Place  
Natick, MA 01760-1537  
www.bostonscientific.com/endoscopy

**Ordering Information**  
**1.800.225.3226**

© 2012 Boston Scientific Corporation  
or its affiliates. All rights reserved.

ENDO-91501-AA July 2012