



Biopsy

Biopsy


Forceps

Multibite™ Single-Use Multiple Sample Biopsy Forceps.....	55
Radial Jaw™ 3 Single-Use Max Capacity Biopsy Forceps	55
Radial Jaw 4 Single-Use Hot Biopsy Forceps	54
Radial Jaw 4 Single-Use Jumbo Biopsy Forceps	53
Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps	53
Radial Jaw 4 Single-Use Pediatric Biopsy Forceps	54
Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps.....	53



Biopsy

Forceps



Radial Jaw™ 4 Single-Use Jumbo Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M00513343	Radial Jaw 4 Jumbo Forceps	2.8	160	3.2	Yellow	Box 5
M00513353	Radial Jaw 4 Jumbo Forceps w/ Needle	2.8	160	3.2	Yellow	Box 5
M00513363	Radial Jaw 4 Jumbo Forceps	2.8	240	3.2	Orange	Box 5
M00513373	Radial Jaw 4 Jumbo Forceps w/ Needle	2.8	240	3.2	Orange	Box 5
M00513361	Radial Jaw 4 Jumbo Forceps	2.8	240	3.2	Orange	Box 20
M00513371	Radial Jaw 4 Jumbo Forceps w/ Needle	2.8	240	3.2	Orange	Box 20
M00513362	Radial Jaw 4 Jumbo Forceps	2.8	240	3.2	Orange	Box 40
M00513372	Radial Jaw 4 Jumbo Forceps w/ Needle	2.8	240	3.2	Orange	Box 40

Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M00513303	Radial Jaw 4 Large Capacity Forceps	2.4	160	2.8	Yellow	Box 5
M00513313	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	160	2.8	Yellow	Box 5
M00513323	Radial Jaw 4 Large Capacity Forceps	2.4	240	2.8	Orange	Box 5
M00513333	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	240	2.8	Orange	Box 5
M00513301	Radial Jaw 4 Large Capacity Forceps	2.4	160	2.8	Yellow	Box 20
M00513311	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	160	2.8	Yellow	Box 20
M00513321	Radial Jaw 4 Large Capacity Forceps	2.4	240	2.8	Orange	Box 20
M00513331	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	240	2.8	Orange	Box 20
M00513302	Radial Jaw 4 Large Capacity Forceps	2.4	160	2.8	Yellow	Box 40
M00513312	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	160	2.8	Yellow	Box 40
M00513322	Radial Jaw 4 Large Capacity Forceps	2.4	240	2.8	Orange	Box 40
M00513332	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	240	2.8	Orange	Box 40

Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M00513383	Radial Jaw 4 Standard Capacity Forceps	2.2	160	2.8	Yellow	Box 5
M00513393	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	160	2.8	Yellow	Box 5
M00513403	Radial Jaw 4 Standard Capacity Forceps	2.2	240	2.8	Orange	Box 5
M00513413	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	240	2.8	Orange	Box 5
M00513381	Radial Jaw 4 Standard Capacity Forceps	2.2	160	2.8	Yellow	Box 20
M00513391	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	160	2.8	Yellow	Box 20
M00513401	Radial Jaw 4 Standard Capacity Forceps	2.2	240	2.8	Orange	Box 20
M00513411	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	240	2.8	Orange	Box 20
M00513382	Radial Jaw 4 Standard Capacity Forceps	2.2	160	2.8	Yellow	Box 40
M00513392	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	160	2.8	Yellow	Box 40
M00513402	Radial Jaw 4 Standard Capacity Forceps	2.2	240	2.8	Orange	Box 40
M00513412	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	240	2.8	Orange	Box 40

Biopsy


Forceps (Continued)

Radial Jaw™ 4 Single-Use Pediatric Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Packaging
M00513443	Radial Jaw 4 Pediatric Biopsy Forceps	1.8	160	≥ 2.0	Box 5
M00513453	Radial Jaw 4 Pediatric Biopsy Forceps w/ Needle	1.8	160	≥ 2.0	Box 5

Hot Biopsy Forceps




Radial Jaw 4 Single-Use Hot Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Packaging
M00515031	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	≥ 2.8	Box 5
M00515032	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	≥ 2.8	Box 20
M00515033	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	≥ 2.8	Box 40

BIOPSY


Forceps (Continued)



Radial Jaw™ 3 Single-Use Max Capacity Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M005 1586 1	Radial Jaw 3 Max Capacity Forceps	3.3	160	3.7	Yellow	Box 5
M005 1587 1	Radial Jaw 3 Max Capacity Forceps w/ Needle	3.3	160	3.7	Yellow	Box 5
M005 1588 1	Radial Jaw 3 Max Capacity Forceps	3.3	240	3.7	Orange	Box 5
M005 1589 1	Radial Jaw 3 Max Capacity Forceps w/ Needle	3.3	240	3.7	Orange	Box 5
M005 1589 2	Radial Jaw 3 Max Capacity Forceps w/ Needle	3.3	240	3.7	Orange	Box 20



Multibite™ Single-Use Multiple Sample Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M005 1010 1	Multibite Multiple Sample Biopsy Forceps	2.4	160	2.8	Yellow	Box 5
M005 1012 1	Multibite Multiple Sample Biopsy Forceps	2.4	240	2.8	Orange	Box 5

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.