

Take Control

The Resolution 360 Clip builds on the strong clinical and economic track record established by Resolution Clip while offering new benefits to providers and patients including controlled rotation in tortuous anatomy designed for more accurate clip placement and shorter, more efficient procedures.

The Resolution 360 Clip has received CE mark and is indicated for:

1. Hemostasis for:
 - Mucosal/sub-mucosal defects < 3 cm,
 - Bleeding ulcers, Arteries < 2 mm, Polyps < 1.5 cm in diameter, Diverticula in the colon,
 - Prophylactic clipping to reduce the risk of delayed bleeding post lesion resection.
2. As a supplementary method, closure of GI tract luminal perforations < 20 mm that can be treated conservatively.
3. Endoscopic marking.
4. Anchoring to affix jejunal feeding tubes to the wall of the small bowel.



Resolution 360™ Clip



Economic and Clinical Value

Features	Potential Benefits
Rotation in either direction	<ul style="list-style-type: none"> • 1 to 1 rotation designed for more accurate clip placement, and shorter, more efficient procedures.
Choice of two rotation mechanisms	<ul style="list-style-type: none"> • Innovative braided design enables the physician or nurse/technician to provide controlled clip placement for potential procedural efficiency and accuracy. • Ability to maintain hand in handle while rotating to preserve clip positioning for deployment.
Device Built on Resolution Platform	<ul style="list-style-type: none"> • Resolution Clip has a reliable clinical history with 11 years of proven success, 6 million clips sold and over 220 clinical articles.
More Clinical Indications, including Prophylactic Indication	<ul style="list-style-type: none"> • More indications than its main competitors*, presenting opportunities to realize efficiencies in inventory carrying cost, lab storage space, and labor associated with managing more product.
Tactile Feedback	<ul style="list-style-type: none"> • Tactile feedback may allow the user to confirm hemostasis before deploying the clips, reducing the likelihood of using a second clip. • Upon deployment, Resolution Clips are designed to lock securely in place for improved retention, which may lead to fewer hospital readmissions and repeat procedures. In a randomized, prospective study on hemoclipping of chronic canine ulcers, "the Resolution Clip was retained significantly longer" than other hemoclips in the study with no delay of ulcer healing. (The results of this canine study may not be representative of performance in humans.)**
Easy Setup/Use	<ul style="list-style-type: none"> • Elimination of plastic sheath allows for fewer setup steps.

Product Codes

Order Number	Description	Unit	Availability
M00521230	Resolution 360 Clip 235cm	Box 1	Available now
M00521231	Resolution 360 Clip 235cm	Box 10	Available now
M00521232	Resolution 360 Clip 235cm	Box 20	Coming Soon
M00521233	Resolution 360 Clip 235cm	Box 40	Coming Soon
M00521240	Resolution 360 Clip 155cm	Box 1	Coming Soon
M00521241	Resolution 360 Clip 155cm	Box 10	Coming Soon
M00521242	Resolution 360 Clip 155cm	Box 20	Coming Soon
M00521243	Resolution 360 Clip 155cm	Box 40	Coming Soon



* Main competitors in EU are Instinct, QuickClip pro and AGS clips.

** Jensen DM, Machiicado GA. Hemoclipping of chronic canine ulcers: a randomized, prospective study of initial deployment success, clip retention rates, and ulcer healing. *Gastrointestinal Endoscopy*, 2009 Nov;70(5):969-75. doi: 10.1016/j.gie.2009.04.052. Epub 2009 Jul 28.

All cited trademarks are the property of their respective owners.

Caution: Law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations.

ENDO-369001-AA January 2016. Printed in The Netherlands by DeBudelse/creativeservices.

Boston Scientific
Advancing science for life™

www.bostonscientific.eu

© 2016 Boston Scientific Corporation or its affiliates. All rights reserved.

CE 0344