

Polypectomy Snares


**Boston
Scientific**
Advancing science for life™



Rotatable Snare
Designed to provide 360° rotation
to help access difficult polyps




Sensation™ Short Throw Snare
Short throw handle is designed for
nurse and technician comfort



Captivator™ II Snare
Rounded shape intended to facilitate polyp capture



Captiflex™ Snare
Incorporates a braided wire design
engineered to produce a flexible snare



Captivator™ Snare
Offers a variety of snare loop configurations to
accommodate a wide range of therapeutic functions

Profile™ Pediatric Snare
Engineered to fit through most 2.0mm
pediatric working channels



Polypectomy Snares

Available in a broad range of loop sizes and shape configurations.

Ordering Information

CAPTIVATOR™ II SINGLE-USE SNARE

Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00561291	Extra Large Rounded	Stiff	33	2.4	240	2.8	Box 10

CAPTIVATOR™ SINGLE-USE SNARES

Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00562301	Small Oval	Stiff	13	2.4	240	2.8	Box 10
M00562321	Medium Oval	Stiff	27	2.4	240	2.8	Box 10
M00561311	Large Oval	Medium Stiff	30	2.4	240	2.8	Box 10
M00562391	Large Oval (thin wire)	Flexible	30	2.4	240	2.8	Box 10
M00562451	Small Hexagonal	Stiff	13	2.4	240	2.8	Box 10
M00562341	Medium Hexagonal	Stiff	27	2.4	240	2.8	Box 10
M00562371	Medium Crescent	Stiff	27	2.4	240	2.8	Box 10

SINGLE-USE ROTATABLE SNARES*

Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00561821	Small Oval	Medium Stiff	13	2.4	195	2.8	Box 10
M00561831	Medium Oval	Medium Stiff	20	2.4	195	2.8	Box 10

* With Star Catheter Technology

SENSATION™ SINGLE-USE SHORT THROW SNARES

Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00562691	Small Oval	Flexible	13	2.4	240	2.8	Box 5
M00562692	Small Oval	Flexible	13	2.4	240	2.8	Box 20
M00562693	Small Oval	Flexible	13	2.4	240	2.8	Box 40
M00560321	Small Oval	Medium Stiff	13	2.4	240	2.8	Box 5
M00562671	Medium Oval	Flexible	27	2.4	240	2.8	Box 5
M00562672	Medium Oval	Flexible	27	2.4	240	2.8	Box 20
M00562673	Medium Oval	Flexible	27	2.4	240	2.8	Box 40
M00560311	Medium Oval	Medium Stiff	27	2.4	240	2.8	Box 5
M00562651	Large Oval	Medium Stiff	30	2.4	240	2.8	Box 5
M00562652	Large Oval	Medium Stiff	30	2.4	240	2.8	Box 20
M00560331	Medium Crescent	Medium Stiff	27	2.4	240	2.8	Box 5

CAPTIFLEX™ SINGLE-USE SNARES

Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00562471	Extra Small Oval	Flexible	11	2.4	240	2.8	Box 10
M00562421	Small Oval	Flexible	13	2.4	240	2.8	Box 10
M00562422	Small Oval	Flexible	13	2.4	240	2.8	Box 40
M00562401	Medium Oval	Flexible	27	2.4	240	2.8	Box 10
M00562402	Medium Oval	Flexible	27	2.4	240	2.8	Box 40

PROFILE™ SINGLE-USE SNARES

Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00562531	Extra Small Oval	Flexible	11	1.9	240	2.0	Box 10
M00562551	Small Oval	Flexible	13	1.9	240	2.0	Box 10
M00562571	Medium Oval	Flexible	27	1.9	240	2.0	Box 10

Distributed by Boston Scientific Corporation, Natick, Massachusetts; manufactured by BSC Coyoil, Costa Rica.

Boston Scientific
Advancing science for life™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientific.com/endo-resources

Ordering Information
1.888.272.1001

©2013 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-121901-AA February 2013

Captiflex, Captivator, Profile and Sensation are registered or unregistered trademarks of Boston Scientific Corporation or its affiliates.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.