

# Polypectomy Snares

**Boston  
Scientific**  
Advancing science for life™



**Rotatable Snare**  
Designed to provide 360° rotation  
to help access difficult polyps



**Sensation™ Short Throw Snare**  
Short throw handle is designed for  
nurse and technician comfort



**Captivator™ II Snare**  
Rounded shape intended to facilitate polyp capture



**Captiflex™ Snare**  
Incorporates a braided wire design  
engineered to produce a flexible snare



**Captivator™ Snare**  
Offers a variety of snare loop configurations to  
accommodate a wide range of therapeutic functions




**Profile™ Pediatric Snare**  
Engineered to fit through most 2.0mm  
pediatric working channels

# Polypectomy Snares

Available in a broad range of loop sizes and shape configurations.

## Ordering Information



CAPTIVATOR™ II SINGLE-USE SNARE							
Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00561291	Extra Large <b>Rounded</b>	Stiff	33	2.4	240	2.8	Box 10

CAPTIVATOR™ SINGLE-USE SNARES							
Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00562301	Small <b>Oval</b>	Stiff	13	2.4	240	2.8	Box 10
M00562321	Medium <b>Oval</b>	Stiff	27	2.4	240	2.8	Box 10
M00561311	Large <b>Oval</b>	Medium Stiff	30	2.4	240	2.8	Box 10
M00562391	Large <b>Oval</b> (thin wire)	Flexible	30	2.4	240	2.8	Box 10
M00562451	Small <b>Hexagonal</b>	Stiff	13	2.4	240	2.8	Box 10
M00562341	Medium <b>Hexagonal</b>	Stiff	27	2.4	240	2.8	Box 10
M00562371	Medium <b>Crescent</b>	Stiff	27	2.4	240	2.8	Box 10

SINGLE-USE ROTATABLE SNARES*							
Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00561821	Small <b>Oval</b>	Medium Stiff	13	2.4	195	2.8	Box 10
M00561831	Medium <b>Oval</b>	Medium Stiff	20	2.4	195	2.8	Box 10

\* With Star Catheter Technology

SENSATION™ SINGLE-USE SHORT THROW SNARES							
Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00562691	Small <b>Oval</b>	Flexible	13	2.4	240	2.8	Box 5
M00562692	Small <b>Oval</b>	Flexible	13	2.4	240	2.8	Box 20
M00562693	Small <b>Oval</b>	Flexible	13	2.4	240	2.8	Box 40
M00560321	Small <b>Oval</b>	Medium Stiff	13	2.4	240	2.8	Box 5
M00562671	Medium <b>Oval</b>	Flexible	27	2.4	240	2.8	Box 5
M00562672	Medium <b>Oval</b>	Flexible	27	2.4	240	2.8	Box 20
M00562673	Medium <b>Oval</b>	Flexible	27	2.4	240	2.8	Box 40
M00560311	Medium <b>Oval</b>	Medium Stiff	27	2.4	240	2.8	Box 5
M00562651	Large <b>Oval</b>	Medium Stiff	30	2.4	240	2.8	Box 5
M00562652	Large <b>Oval</b>	Medium Stiff	30	2.4	240	2.8	Box 20
M00560331	Medium <b>Crescent</b>	Medium Stiff	27	2.4	240	2.8	Box 5

CAPTIFLEX™ SINGLE-USE SNARES							
Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00562471	Extra Small <b>Oval</b>	Flexible	11	2.4	240	2.8	Box 10
M00562421	Small <b>Oval</b>	Flexible	13	2.4	240	2.8	Box 10
M00562422	Small <b>Oval</b>	Flexible	13	2.4	240	2.8	Box 40
M00562401	Medium <b>Oval</b>	Flexible	27	2.4	240	2.8	Box 10
M00562402	Medium <b>Oval</b>	Flexible	27	2.4	240	2.8	Box 40

PROFILE™ SINGLE-USE SNARES							
Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00562531	Extra Small <b>Oval</b>	Flexible	11	1.9	240	2.0	Box 10
M00562551	Small <b>Oval</b>	Flexible	13	1.9	240	2.0	Box 10
M00562571	Medium <b>Oval</b>	Flexible	27	1.9	240	2.0	Box 10

Distributed by Boston Scientific Corporation, Natick, Massachusetts; manufactured by BSC Coyal, Costa Rica.

**Boston Scientific**  
Advancing science for life™

Boston Scientific Corporation  
One Boston Scientific Place  
Natick, MA 01760-1537  
[www.bostonscientific.com/endo-resources](http://www.bostonscientific.com/endo-resources)

**Ordering Information**  
**1.888.272.1001**

©2013 Boston Scientific Corporation or its affiliates. All rights reserved.

ENDO-121901-AA February 2013

Captiflex, Captivator, Profile and Sensation are registered or unregistered trademarks of Boston Scientific Corporation or its affiliates.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

**CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.**