

Polypectomy Snares

Boston
Scientific
Advancing science for life™

UPDATED SNARE PRODUCT LABELS AND SIZING INFORMATION

SIZE NAMING CONVENTION		
Legacy Size Description	New Size Description	Loop Width
Mini Micro	Extra Small	11mm
Micro	Small	13mm
Standard	Medium	27mm
Jumbo	Large	30mm
NA	Extra Large	33mm

Single-Use Polypectomy Snare Ordering Information

SENSATION™		SIZE DESCRIPTOR
		LARGE OVAL
		MEDIUM OVAL
		SMALL OVAL
		MEDIUM CRESCENT

SENSATION SHORT THROW

Order Number	Old Descriptor	NEW Descriptor	Loop Width (mm)	Minimum Working Channel (mm)	Quantity
M00562691	Micro Oval	Small Oval - Flexible	13	2.8	Box 5
M00562692	Micro Oval	Small Oval - Flexible	13	2.8	Box 20
M00562693	Micro Oval	Small Oval - Flexible	13	2.8	Box 40
M00560321	Micro Oval	Small Oval - Medium Stiff	13	2.8	Box 5
M00562671	Standard Oval	Medium Oval - Flexible	27	2.8	Box 5
M00562672	Standard Oval	Medium Oval - Flexible	27	2.8	Box 20
M00562673	Standard Oval	Medium Oval - Flexible	27	2.8	Box 40
M00560311	Standard Oval	Medium Oval - Medium Stiff	27	2.8	Box 5
M00562651	Jumbo Oval	Large Oval - Medium Stiff	30	2.8	Box 5
M00562652	Jumbo Oval	Large Oval - Medium Stiff	30	2.8	Box 20
M00560331	Crescent	Medium Crescent - Medium Stiff	27	2.8	Box 5

CAPTIFLEX™		SIZE DESCRIPTOR
		MEDIUM OVAL
		SMALL OVAL
		EXTRA SMALL OVAL

CAPTIFLEX

Order Number	Old Descriptor	NEW Descriptor	Loop Width (mm)	Minimum Working Channel (mm)	Quantity
M00562471	Mini-Micro Oval	Extra Small Oval - Flexible	11	2.8	Box 10
M00562421	Micro Oval	Small Oval - Flexible	13	2.8	Box 10
M00562422	Micro Oval	Small Oval - Flexible	13	2.8	Box 40
M00562401	Standard Oval	Medium Oval - Flexible	27	2.8	Box 10
M00562402	Standard Oval	Medium Oval - Flexible	27	2.8	Box 40

PROFILE™		SIZE DESCRIPTOR
		MEDIUM OVAL
		SMALL OVAL
		EXTRA SMALL OVAL

PROFILE

Order Number	Old Descriptor	NEW Descriptor	Loop Width (mm)	Minimum Working Channel (mm)	Quantity
M00562531	Mini Micro Oval	Extra Small Oval - Flexible	11	2.0	Box 10
M00562551	Micro Oval	Small Oval - Flexible	13	2.0	Box 10
M00562571	Wide Oval	Medium Oval - Flexible	27	2.0	Box 10

NEW INFORMATION
Snare Size loop width indicated

EXTRA LARGE 33mm

CAPTIVATOR™ II
EXTRA LARGE ROUNDED - STIFF

Single-Use Polypectomy Snare. Loop para polipectomia do cólon. Polipectomia para uso único. Einweg-Polypektomieschlinge. Cappo monouso per polipectomia del colon. Polipectomia para Lina Única a Utilização. Polipektomiszynge für Polipektomie für Lina Única a Utilização. Polipektomická smyčka. Jednorazowa petla do polipektomii. 일회용 용종절제술 제거기. Tek Kullu. 一次性用内窥镜术器械. 一次性用内窥镜术器械. 一次性用内窥镜术器械.

Contents (10)

WC 2.8mm Minimum Required Working Channel

240cm

33mm

REF Catalog No. 6129

Use By 2010-07

Label Specification Part Number: 90593289-01A

CAPTIVATOR™		SIZE DESCRIPTOR
		MEDIUM CRESCENT
		LARGE OVAL
		MEDIUM OVAL
		LARGE THIN-WIRE OVAL
		SMALL OVAL
		MEDIUM HEXAGONAL
		SMALL HEXAGONAL
		SMALL OVAL
		MEDIUM HEXAGONAL

CAPTIVATOR

Order Number	Old Descriptor	NEW Descriptor	Loop Width (mm)	Minimum Working Channel (mm)	Quantity
M00562301	Micro Oval	Small Oval - Stiff	13	2.8	Box 10
M00562321	Standard Oval	Medium Oval - Stiff	27	2.8	Box 10
M00561311	Jumbo Oval	Large Oval - Medium Stiff	30	2.8	Box 10
M00562391	Thin-Wire Jumbo Oval	Large Oval (Thin-Wire) - Flexible	30	2.8	Box 10
M00562451	Micro-Hex	Small Hexagonal - Stiff	13	2.8	Box 10
M00562341	Hexagonal	Medium Hexagonal - Stiff	27	2.8	Box 10
M00562371	Crescent	Medium Crescent - Stiff	27	2.8	Box 10

CAPTIVATOR™ II		SIZE DESCRIPTOR
		EXTRA LARGE ROUND
		EXTRA LARGE ROUND

CAPTIVATOR II

Order Number	Old Descriptor	NEW Descriptor	Loop Width (mm)	Minimum Working Channel (mm)	Quantity
M00561291	Rounded Shape	Extra Large Rounded - Stiff	33	2.8	Box 10

SINGLE USE ROTATABLE SNARES		SIZE DESCRIPTOR
		STAR CATHETER TECHNOLOGY
		STAR CATHETER TECHNOLOGY

ROTATABLE SNARES

Order Number	Product Description	NEW Descriptor	Loop Width Diameter (mm)	Sheath O.D. (mm)	Quantity
M00561821	Single-Use Polypectomy Snare	Small Oval Medium Stiff	13	2.8	Box 10
M00561831	Single-Use Polypectomy Snare	Medium Oval Medium Stiff	20	2.8	Box 10