



WallFlex™
Biliary RX Stents

A proven
portfolio.



Choose the biliary metal stent that's right for the patient and their unique stricture with the comprehensive Boston Scientific portfolio.



500 publications



15+ years

WallFlex Biliary RX Fully Covered Stent System

A recent addition to the WallFlex Biliary suite, these fully covered stents feature a long length, designed to help clinicians reach strictures that were previously difficult to access.




WallFlex Biliary RX Partially Covered Stent System

- Studies have shown high patency rates with a low migration rate¹
- Demonstrated migration rate as low as 0%¹



WallFlex Biliary RX Uncovered Stent System

A well-trusted tool for biliary malignancy treatment



Expand your options for treating patients.

Learn more about Boston Scientific biliary stents at [bostonscientific.com/WallFlex](https://www.bostonscientific.com/WallFlex).

	WallFlex Biliary RX Fully Covered Stent System	WallFlex Biliary RX Partially Covered Stent System		WallFlex Biliary RX Uncovered Stent System
Diameters (mm)	10	8	10	10
Lengths (mm)	100	60, 80, 100	40, 60, 80, 100	40, 60, 80, 100
Malignant biliary neoplasm	✓		✓	✓
Bridge to surgery	✓		✓	✓


Ordering information

Fully Covered Stent System RMV with Permalume™ covering




Order Number	Diameter (mm)	Length (mm)	Covered Length (mm) (PC Only)	Catheter Diameter (Fr/mm)	Guidewire Diameter (in/mm)
M00570470	10	40	-	8.5 (2.83)	.035" (.89)
M00570480	10	60	-	8.5 (2.83)	.035" (.89)
M00570490	10	80	-	8.5 (2.83)	.035" (.89)

Partially Covered Stents with Permalume covering



M00570720	10	40	28	8.5 (2.83)	.035" (.89)
M00570730	10	60	48	8.5 (2.83)	.035" (.89)
M00570740	10	80	68	8.5 (2.83)	.035" (.89)
M00570760	10	100	88	8.5 (2.83)	.035" (.89)

Uncovered Stents



M00570890	10	40	-	8.0 (2.67)	.035" (.89)
M00570640	10	60	-	8.0 (2.67)	.035" (.89)
M00570650	10	80	-	8.0 (2.67)	.035" (.89)
M00570660	10	100	-	8.0 (2.67)	.035" (.89)

Indications for use

WallFlex Biliary RX Stent System, WallFlex Biliary PLUS RX Stent System: indicated for use in the palliative treatment of biliary strictures produced by malignant neoplasms and relief of malignant biliary obstruction prior to surgery.

WallFlex Biliary RX Stent System RMV: indicated for use in the palliative treatment of biliary strictures produced by malignant neoplasms, relief of malignant biliary obstruction prior to surgery, and for indwell up to 12 months in the treatment of benign biliary strictures secondary to chronic pancreatitis.

Warning: The safety and effectiveness of this device for use in vascular system has not been established.

Sources

- Yokota Y, Fukasawa M, Takano S, et al. Partially covered metal stents have longer patency than uncovered and fully covered metal stents in the management of distal malignant biliary obstruction: A retrospective study. *BMC Gastroenterol.* 2017;17:105.

The information presented in this document by Boston Scientific is intended for healthcare professionals in India. [The information presented in this document is intended to provide balanced, scientific, and evidence based answers to unsolicited medical questions.] The information should not be, however, treated as comprehensive and does not intend to provide diagnosis, treatment or any medical or health advice. For full information about the product, please refer to the appropriate product labeling. Responsibility for patient care resides with the healthcare professional on the basis of his or her professional licence, experience and knowledge of the patient. Healthcare professionals must rely on their judgment when deciding which treatments and procedures to use with patients.

All cited trademarks are the property of their respective owners. ©2026 Boston Scientific Corporation or its affiliates. All rights reserved.