

Ordering Information

CAPTIVATOR™ II SINGLE-USE POLYPECTOMY SNARES

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00561221	10mm Rounded	Stiff	10	2.4	240	2.8	Box 10
M00561222	10mm Rounded	Stiff	10	2.4	240	2.8	Box 20
M00561223	10mm Rounded	Stiff	10	2.4	240	2.8	Box 40
M00561231	15mm Rounded	Stiff	15	2.4	240	2.8	Box 10
M00561232	15mm Rounded	Stiff	15	2.4	240	2.8	Box 20
M00561233	15mm Rounded	Stiff	15	2.4	240	2.8	Box 40
M00561241	20mm Rounded	Stiff	20	2.4	240	2.8	Box 10
M00561242	20mm Rounded	Stiff	20	2.4	240	2.8	Box 20
M00561243	20mm Rounded	Stiff	20	2.4	240	2.8	Box 40
M00561191	25mm Rounded	Stiff	25	2.4	240	2.8	Box 10
M00561192	25mm Rounded	Stiff	25	2.4	240	2.8	Box 20
M00561193	25mm Rounded	Stiff	25	2.4	240	2.8	Box 40
M00561291	33mm Rounded	Stiff	33	2.4	240	2.8	Box 10

CAPTIVATOR™ SINGLE-USE POLYPECTOMY SNARES

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00562301	Small Oval	Stiff	13	2.4	240	2.8	Box 10
M00562321	Medium Oval	Stiff	27	2.4	240	2.8	Box 10
M00562391	Large Oval (thin wire)	Flexible	30	2.4	240	2.8	Box 10
M00561311	Large Oval	Medium Stiff	30	2.4	240	2.8	Box 10
M00562451	Small Hexagonal	Stiff	13	2.4	240	2.8	Box 10
M00562341	Medium Hexagonal	Stiff	27	2.4	240	2.8	Box 10
M00562371	Medium Crescent	Stiff	27	2.4	240	2.8	Box 10

ROTATABLE SINGLE-USE POLYPECTOMY SNARES*

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00561821	Small Oval	Medium Stiff	13	2.4	240	2.8	Box 10
M00561831	Medium Oval	Medium Stiff	20	2.4	240	2.8	Box 10

* With Star Catheter Technology

SENSATION™ SINGLE-USE POLYPECTOMY SNARES - SHORT THROW SNARES

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00562691	Small Oval	Flexible	13	2.4	240	2.8	Box 5
M00562692	Small Oval	Flexible	13	2.4	240	2.8	Box 20
M00562693	Small Oval	Flexible	13	2.4	240	2.8	Box 40
M00560321	Small Oval	Medium Stiff	13	2.4	240	2.8	Box 5
M00562671	Medium Oval	Flexible	27	2.4	240	2.8	Box 5
M00562672	Medium Oval	Flexible	27	2.4	240	2.8	Box 20
M00562673	Medium Oval	Flexible	27	2.4	240	2.8	Box 40
M00560311	Medium Oval	Medium Stiff	27	2.4	240	2.8	Box 5
M00562651	Large Oval	Medium Stiff	30	2.4	240	2.8	Box 5
M00562652	Large Oval	Medium Stiff	30	2.4	240	2.8	Box 20
M00560331	Medium Crescent	Medium Stiff	27	2.4	240	2.8	Box 5

CAPTIFLEX™ SINGLE-USE POLYPECTOMY SNARES

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00562471	Extra Small Oval	Flexible	11	2.4	240	2.8	Box 10
M00562421	Small Oval	Flexible	13	2.4	240	2.8	Box 10
M00562422	Small Oval	Flexible	13	2.4	240	2.8	Box 40
M00562401	Medium Oval	Flexible	27	2.4	240	2.8	Box 10
M00562402	Medium Oval	Flexible	27	2.4	240	2.8	Box 40

RADIAL JAW™ 4 HOT BIOPSY FORCEPS

Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00515031	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	Box 5
M00515032	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	Box 20
M00515033	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	Box 40

⊗ These products are for single use only. Do not reuse.

Boston Scientific

Advancing science for life™

www.bostonscientific.eu

© 2024 Boston Scientific Corporation or its affiliates. All rights reserved.

ENDO-436716-AC

CE 2797


All trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.

Tissue Resection Devices

Boston Scientific



Polypectomy Snares



Available in a broad range of loop sizes and shape configurations. All of them approved for both hot and cold snaring.

Captivator™ II Snare

- Rounded shape to facilitate polyps and mucosal resection. Short distance from scope tip designed for visualisation and tip control

Captivator™ Snare

- Offers a variety of snare loop configurations
- The therapeutic functions indicated for Captivator are: removal or cauterisation of diminutive polyps, sessile polyps, pedunculated polyps and tissue from within the gastrointestinal tract

Rotatable Snare

- Designed to provide 360° rotation to help access difficult polyps

Sensation™ Short Throw Handle Snare

- Short throw handle is designed for nurse and technician comfort

Captiflex™ Snare

- Incorporates a braided wire design engineered to produce a flexible snare

Radial Jaw™ 4 Hot Biopsy Forceps

Hot forceps available in our family of market-leading single-use biopsy forceps.

Jaw Configuration

Designed to provide:

- Tissue specimens for excellent sample handling, preparation and accurate histological diagnosis
- Enhanced passability through tortuous anatomy

Micro-Mesh Tooth Design

- Engineered for a clean, precise bite

Streamlined Catheter

- Designed to provide the right balance of columnar strength and flexibility for excellent pushability and control during scope passage

Positioning Markers

- Engineered to assist in the insertion and withdrawal of forceps to enhance procedural safety and efficiency

