



FARAPULSE™

Pulsed Field Ablation System



FARAWAVE™
Pulsed Field Ablation Catheter


FARAWAVE™ NAV
Pulsed Field Ablation Catheter



FARAWAVE™ and FARAWAVE™ NAV Pulsed Field Ablation (PFA) Catheters
Ordering Information

Reference Number	
FARAWAVE PFA Catheter M004PF41M401	31 mm (distal diameter)
FARAWAVE PFA Catheter M004PF41M402	35 mm (distal diameter)
FARAWAVE NAV PFA Catheter M004PF41M411	31 mm (distal diameter)
FARAWAVE NAV PFA Catheter M004PF41M412	35 mm (distal diameter)
Catheter Specifications	
Catheter type	Over-the-wire
Shaft diameter	12F
Usable length	115 cm
Number of electrodes	20
Compatibility	
Guidewire	0.035" diameter, 180 cm minimum length
Sheath	FARADrive (REF M004PF21M402)
Accessories	
FARASTAR catheter connection cable for FARAWAVE	REF M004PF41M434
PFA Rx Gen 2 Connection Cable for FARAWAVE NAV	REF M004PF41M444
Merit inqWire	REF IQ35F180J1O5RS

FARAWAVE and FARAWAVE NAV PFA Catheters are sterile and single-use only






FARADrive™ Steerable Sheath Ordering Information

Reference Number	
M004PF21M402	Clear shaft
Sheath Specifications	
Usable length	74 cm
Sheath overall length	91 cm
Inner diameter	13F
Deflection	Unidirectional, 35 mm reach
Radiopaque	Distal end
Dilator Specifications	
Total length	97 cm
Usable length	94 cm
Compatibility	
Guidewires	Maximum diameter 0.035"
Pulsed Field Ablation Catheter	FARAWAVE (REF M004PF41M401, M004PF41M402) FARAWAVE NAV (REF M004PF41M411, REF M004PF41M412)
VersaCross™ integration	VersaCross Connect™ Access Solution for FARADrive (180 cm Pigtail, VXAK0041) (180 cm J-Tip, VXAK0045)

FARADrive Sheath and dilator are sterile and single-use only





FARASTAR™ Pulsed Field Ablation Generator Ordering Information

Reference Number	
M004PF61M401	FARASTAR Pulsed Field Ablation Generator
M004PF61M407	FARASTAR Recording System Module
Shipping Information – FARASTAR Pulsed Field Ablation Generator	
Weight	250 kg
Dimensions	L 120 x H 200 x W 80 cm
Shipping Information – FARASTAR Recording System Module	
Weight	552 lbs (250 kg)
Dimensions	L 120 x H 70 x W 80 cm
Technical Specifications	
Input voltage	120V / 60Hz
Output voltage	Proprietary FARAPULSE biphasic, bipolar waveform
Weight	243 lbs (110 kg)
Dimensions – FARASTAR Pulsed Field Ablation Generator	L 79 x W 61 x H 117 cm (160 cm with monitor)
Dimensions – FARASTAR Recording System Module	L 41 x W 19 x H 46 cm
User Interface	Articulating touchscreen
Accessories	
M004PF61M404	FARASTAR Stimulation Module Cable
M004PF61M405	FARASTAR EGM Cable
M004PF61M409	FARASTAR Stimulation Module Male Cable
M004PF61M410	FARASTAR Stimulation Module Female Cable
M004PF61M411	FARASTAR Stimulation Module Y-Cable Long
M004PF61M412	FARASTAR Stimulation Module Y-Cable Short
M004PF61M413	FARASTAR Recording System Module Catheter Pin Cable
M004PF61M415	FARASTAR Recording System Module ECG Trunk Cable
M004PF61M416	FARASTAR Recording System Module ECG Output Module
M004PF61M417	FARASTAR Recording System Module EGM Input Module

All dimensions are nominal unless otherwise stated.

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Instructions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

© 2024 Boston Scientific Corporation or its affiliates. All rights reserved. All trademarks are the property of their respective owners.

EP-1779710-AB

**Boston
Scientific**
Advancing science for life™

Cardiology
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

Medical Professionals:
1.800.CARDIAC (227.3422)
Customer Service:
1.888.272.1001