


## PERIVAC™

Pericardiocentesis Kits

### PRE-ASSEMBLED KIT FOR PERICARDIAL DRAINAGE PROCEDURES



Note: Plastic Fenestrated Drape Towel (not shown)  
PeriVac Kits do not include hypodermic needles

### PROVIDES THE ESSENTIALS FOR PERICARDIAL DRAINAGE/ASPIRATION



Featuring the Flexima™ Pigtail Drainage Catheter

8.3F (2.76mm)  
Catheter size with  
sideholes for drainage

# PERIVAC™

Pericardiocentesis Kits

## Clinical Benefits

- **PeriVac™ Kits** provide the essentials for pericardial drainage/aspiration
- **Pre-assembled**, “Ready to Use” kit, designed for safety and speed
- **Convenient and “Ready to Use”**, when initiating the procedure in emergencies
- **New atraumatic Flexima™ 8.3F (2.76 mm) pigtail or straight drainage catheters** with sideholes, designed for proper drainage and prevention of occlusions
- **Graduated 1000ml drainage bag** for easy measurement of collected fluid

## Technical Information

### All Kits include:

- 1 PeriVac™ straight or pigtail catheter
- 1 Guidewire with 3 mm “J” tip, 0.89 cm x 80 cm
- 6 Gauze sponges, 10 cm x 10 cm
- 1 Three-Way Stopcock
- 1 Fluid drainage bag with numbered increments, 1000 ml
- 1 Needle with stylet, 18 G x 15 cm
- 1 Dilator, 50 cm length / 8F (2.66 mm)
- 2 Luer Lock™ Syringes, 10 ml and 20 ml
- 1 Blade with handle, #11
- 2 Prep applicators
- 1 Silk Suture with curved needle, 3-0
- 2 Needles, 21 G x 2.5 - 3.8 cm and 25 G x 2.5 cm
- Fenestrated Drape, 40 cm<sup>2</sup>, not shown



## Ordering Information

### PeriVac™ Pericardiocentesis Kits

Model Number	Description	Quantity
M004 <b>4304</b> 1	PeriVac™ with 8.3F (2.76 mm) Pigtail Catheter	Box of 3 Kits
M004 <b>4314</b> 1	PeriVac™ with 8.3F (2.76 mm) Straight Catheter	Box of 3 Kits
M004 <b>400</b> 1	Spare 1000 ml Drainage Bag (individually sterilized)	Box of 3 Bags

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Information contained herein for distribution outside the USA, France and Japan only.

EP-361601-AB Printed in Germany by medicalvision.

**Boston  
Scientific**  
Advancing science for life™

www.bostonscientific.eu  
© 2018 Boston Scientific Corporation  
or its affiliates. All rights reserved.  
DINP2244EB