

CHILLI II™ COOLED ABLATION CATHETER



CHILLI II Cooled Ablation Catheter Product Information

Product Information

Tip: 7F/4mm Electrode Spacing: 2.5mm

Model Number	Type	Shaft Size	Curve Style
M004 9031 0	Thermocouple	7.5F	Standard
M004 9031K2 0	Thermocouple	7.5F	Large
M004 9031N4 0	Thermocouple	7.5F	Asymmetric 4

Cables and Accessories

Model Number	Description
M004 352104 0	CircuCool™ Pump Tubing Kit


For use with MAESTRO 3000™ and MAESTRO 4000™ Controller

M004 681 0	Cable: MAESTRO 3000™ Pod to CHILLI II Cooled Ablation Catheter (10ft)
------------	---

For use with EPT-1000XP™ Controller

M004 62062 0	Cable: CHILLI II Cooled Ablation Catheter to TTT for EPT-1000XP Controllers (8ft)
--------------	---

*Note: All CHILLI II Catheters are 110 cm in length.
Chilli II Catheters require the use of 1 Valleylab™ Ground Pad (Model M004 354 1).*



CHILLI II™ Cooled Ablation System Product Information

Product Information

CHILLI II Cooled Ablation System

Model Number	Description
M004 8005 0	CircuCool™ Fluid Circulation Pump

Technical Information

CircuCool™ Pump

Model Number	Description
M004 352104 0	CircuCool Pump Tubing Kit
M004 354 1	Valleylab™ Ground Pad

For use with MAESTRO 3000™ and MAESTRO 4000™ Generator

M004 681 0	Cable: MAESTRO 3000 or MAESTRO 4000 Pod to Chilli II Cooled Ablation Catheter (10ft)
------------	--

For use with EPT-1000XP™ Generator

M004 62062 0	Cable: CHILLI II Cooled Ablation Catheter to TTT for EPT-1000XP Controllers (8ft)
--------------	---

Cables and Accessories

Description	Specification
AC Power Requirements	100-240 VAC

Environmental Conditions:

Ambient Temperature	18-35° C
Relative humidity	15-90% noncondensing relative humidity
Pump flow rate	0.6ml/second
Pressure limit	Less than 300psi
Dimensions	11in W x 11.8in D x 9in H
Weight	9.6lbs

Note: All CHILLI II Catheters are 110 cm in length. CHILLI II Catheters require the use of 1 Valleylab™ Ground Pad (Model M004 354 1).



CircuCool™ Pump



Tubing Kit

Boston Scientific
Advancing science for life™

Rhythm Management
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

Medical Professionals:
1.800.CARDIAC (227.3422)

Customer Service: 1-888-272-1001

© 2015 Boston Scientific Corporation
or its affiliates. All rights reserved.

EP-143806-AC NOV2015

The Chilli II Catheter is indicated for ventricular tachycardia. BSC currently has no FDA-approved ablation catheters for the treatment of atrial fibrillation.

CHILLI II, Circucool, MAESTRO 3000, EPT-1000XP are registered or unregistered trademarks of Boston Scientific Corporation or its affiliates. All other trademarks are the property of their respective owners.