

Dynamic XT™, Dynamic Tip™ and CS Assist™

Steerable Diagnostic Catheters

The most widely used steerable catheter for the coronary sinus.




Smooth, precise positioning with unparalleled tip control.



Technical Information

Description	Specifications
Shaft Diameter	6F
Usable Length	110cm and 65cm
Electrode Material	Platinum
Curve Size	Large 4.0
Curve Direction	Unidirectional
Configuration	Bipolar/Quadripolar/ Octapolar/Decapolar

Large 4.0 Curve



Referenced catheter configurations are illustrative representations only and may not reflect actual performance.

Ordering Information

Dynamic XT™ Dynamic XT Steerable Catheter 6F

Catalog number	Electrodes	Curve style	Electrode Spacing	Length	Cables
M004 201112 0	4	Large 4.0	2 mm	110 cm	M004 200088P 0
M004 201103 0	4	Large 4.0	2,5,2 mm	110 cm	M004 200088P 0
M004 201104 0	4	Large 4.0	5 mm	110 cm	M004 200088P 0
M004 201110 0	4	Large 4.0	10 mm	110 cm	M004 200088P 0
M004 201114 0	6	Large 4.0	1 mm	110 cm	M004 200089P 0
M004 201113 0	6	Large 4.0	5 mm	110 cm	M004 200089P 0
M004 201109 0	6	Large 4.0	10 mm	110 cm	M004 200089P 0
M004 201106 0	8	Large 4.0	2 mm	110 cm	M004 200774P 0
M004 201105 0	8	Large 4.0	2,5,2 mm	110 cm	M004 200774P 0
M004 201108 0	8	Large 4.0	2,10,2 mm	110 cm	M004 200774P 0
M004 201107 0	8	Large 4.0	5 mm	110 cm	M004 200774P 0
M004 201101 0	10	Large 4.0	2,5,2 mm	110 cm	M004 560004A 0
M004 201102 0	10	Large 4.0	2,6,2 mm	110 cm	M004 560004A 0
M004 6DYNXT009 0	4	Large 4.0	2,5,2 mm	110 cm	M004 200088P 0
M004 6DYNXT011 0	4	Large 4.0	5 mm	110 cm	M004 200088P 0

Ordering Information

Dynamic Tip™ Steerable Catheter 6F

Catalog number	Electrodes	Curve style	Electrode Spacing	Length	Cables
M004 6DYNTP001 0	10	Large 4.0	2,5,2 mm	110 cm	M004 560004A 0

Ordering Information

CS Assist™ Steerable Diagnostic Catheter 6F

Catalog number	Electrodes	Curve style	Electrode Spacing	Length	Cables
M004 6DYNXT014 0	2	Large 4.0	5 mm	65 cm	M004 560002A 0