

EP·XT™

Steerable Diagnostic Catheter

Built on the trusted Dynamic™ XT Steerable Catheter platform.




Ergonomic handle. Precise tip control.



Technical Information

Description	Specifications
Shaft Diameter	6F
Usable Length	110cm
Electrode Material	Platinum
Curve Size	Large 4.0
Curve Direction	Unidirectional
Configuration	Quadripolar/ Octapolar/ Decapolar

Large 4.0 Curve



Referenced catheter configurations are illustrative representations only and may not reflect actual performance.

Ordering Information

EP·XT™ Steerable Diagnostic Catheter

Catalog number	Electrodes	Curve style	Electrode Spacing	Length	Cables
M004 200770 0	4	Large 4.0	2,5,2mm	110 cm	M004 200088P 0
M004 200769 0	4	Large 4.0	5 mm	110 cm	M004 200088P 0
M004 200795 0	8	Large 4.0	2 mm	110 cm	M004 200774P 0
M004 200794 0	8	Large 4.0	2,5,2 mm	110 cm	M004 200774P 0
M004 200797 0	8	Large 4.0	2,10,2 mm	110 cm	M004 200774P 0
M004 200796 0	8	Large 4.0	5 mm	110 cm	M004 200774P 0
M004 201007 0	10	Large 4.0	2,5,2 mm	110 cm	M004 560004A 0

**Boston
Scientific**
Advancing science for life™

www.bostonscientific.eu

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Information not intended for distribution in France.

EP-403503-AA JUN2016 Printed in Germany by medicalvision.

© 2016 Boston Scientific Corporation
or its affiliates. All rights reserved.
DINEP2253EA