



Endoscopy Product Catalogue

2026

Produktkatalog
Catalogue Produits
Catálogo de Productos
Catalogo dei Prodotti

Biliary	RX Biliary System™	Biliary Lumen Patency	Endoscopic Ultrasound (EUS)	Guidewires
Biopsy	Dilatation	Imaging Systems	Pulmonary	Closure
Haemostasis	Tissue Resection	Esophageal and Enteral Metal Stents	Infection Prevention	Endobariatrics



[Contents](#)

[Alphabetical
Index](#)

[Product Code
Index](#)

Contents

Biliary

Cannulas	Tandem™ XL Triple-Lumen ERCP Cannula	38
	Contour™ ERCP Cannulas	38
Sphincterotomes	TRUEtome™ Jag Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	39
	TRUEtome™ Dreamwire™ Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	39
	TRUEtome™ Hydra Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	39
	TRUEtome™ Jagwire Revolution Cannulating Sphincterotomes with Pre-loaded Jagwire™ Revolution Guidewire	39
	TRUEtome™ Cannulating Sphincterotomes	40
	Microknife™ XL Triple-Lumen Needle-Knife	40
	Ultratome XL™ Triple-Lumen Sphincterotomes	41
Stone Retrieval	Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Distal Injection	41
	Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	41
Accessories	Biliary Accessories	42

RX Biliary System™

Cannulas	Rapid Exchange XL Cannula	44
	RX ERCP Cannulas	44
Sphincterotomes	Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	45
	Dreamtome™ RX Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	45
	Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	46
	Jagtome™ Revolution RX Cannulating Sphincterotome (Autotome™ Pre-loaded with Jagwire Revolution High Performance Guidewire)	46
	Autotome™ Pro RX Cannulating Sphincterotomes	47
	Autotome™ RX Preloaded Sphincterotomes	47
	Autotome™ RX Cannulating Sphincterotomes	47
	RX Needle-Knife XL Triple-Lumen Needle-Knife	48
Stone Retrieval	Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Distal Injection	48
	Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	48
	Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Distal Injection	49
	Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	49
	Trapezoid™ RX Wireguided Retrieval Baskets	49
Dilatation	Hurricane™ RX Biliary Balloon Dilatation Catheters	50
	CRE™ RX Biliary Balloon Dilatation Catheter	50
Cytology	RX Cytology Brush	50
Locking Devices	AutoCap™ RX	51
	RX Locking Devices	51
	Encore™ 26 Inflation Device	51



Biliary Lumen Patency

RX Biliary Plastic Stents	Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System	53
	Advanix™ Stents – Singles	54-55
	NaviFlex™ RX Delivery System	55
	RX Biliary Plastic Stents with RX Delivery System	56
Biliary Plastic Stents	Flexima™ Biliary Stents with Delivery System	57
Pancreatic Plastic Stents	Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Pigtail	58
	Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Straight	58
	Advanix™ Pancreatic Single Stents – Pigtail	59-60
	Advanix™ Pancreatic Single Stents – Straight	61-62
	NaviFlex™ RX Pancreatic Delivery System & NaviFlex™ Pushers	63
RX Biliary Metal Stents	WallFlex™ Biliary RX Fully Covered Stent System RMV (Permalume™ Covering, includes benign indication)	63
	WallFlex™ Biliary RX Partially Covered Stent System (Permalume™ Covering)	64
	WallFlex™ Biliary RX Uncovered Stent System	64
	WallFlex™ Biliary RX Fully Covered Stent System	64
	WallFlex™ Biliary PLUS RX Fully Covered Stent System	65
Radiofrequency Ablation Catheter and Cables	Habib™ EndoHPB Bipolar Radiofrequency Catheter	66

Endoscopic Ultrasound (EUS)

Fine Needle Aspiration (FNA)	Expect™ Slimline Handle Needles	68
Fine Needle Biopsy (FNB)	Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Devices	68
	Acquire™ S Endoscopic Ultrasound Fine Needle Biopsy (FNB) Devices	69
Lumen Apposing Metal Stent System	Hot AXIOS™ Stent and Electrocautery Enhanced Delivery System	69

Guidewires

High Performance	Dreamwire™ High Performance Guidewires	71
	Hydra Jagwire™ High Performance Guidewires	71
	Jagwire™ High Performance Guidewires	72
	Jagwire™ Revolution High Performance Guidewire	72
	NaviPro™ Hydrophilic Guidewires	73
Standard	Amplatz Super Stiff™ Guidewire	73

Biopsy

Forceps	Radial Jaw™ 4 Standard Capacity Biopsy Forceps	75
	Radial Jaw™ 4 Large Capacity Biopsy Forceps	75
	Radial Jaw™ 4 Jumbo Biopsy Forceps	76
	Radial Jaw™ 4 Gastropediatric Biopsy Forceps	76

Contents

Dilatation

CRE™ (Controlled Radial Expansion) Balloon Dilatation Catheters	CRE™ PRO Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic/Biliary) (including DASE indication)	78
	CRE™ Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)	78
	CRE™ PRO Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Biliary) (including DASE indication)	79
	CRE™ Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)	79
	CRE™ Fixed Wire Balloon Dilatation Catheters	80
Balloon Dilators	Rigiflex™ II Achalasia Balloon Dilators	80
Inflation Systems	Alliance™ II Inflation System	81
Dilatation Accessories	Inflation/Deflation Devices and Accessories	81

Imaging Systems

Imaging Systems	EXALT™ Model D Single-Use Duodenoscope	83
	EXALT™ Model D Service Codes	83
Direct Visualisation System	SpyGlass™ DS Capital Equipment	84
	SpyGlass™ DS Disposable Equipment	84
	SpyGlass™ DS Service Codes	84
	Scivita Medical™ Capital Equipment	84
	Scivita Medical™ Disposable Equipment	85
	Scivita Medical™ Service Codes	85
	SpyGlass™ Discover Digital Catheter	86
	SpyGlass™ Discover Service Codes	86
	Autolith™ Touch Biliary EHL System	86
	Autolith™ Service Codes	86

Pulmonary

Broncho Consumables	Scivita Medical™ Capital Equipment	88
	Scivita Medical™ Disposable Equipment	88
	Scivita Medical™ Service Codes	88
Airway Stents	Ultraflex™ Tracheobronchial Covered Stent System – Distal Release	90
	Ultraflex™ Tracheobronchial Uncovered Stent System – Proximal Release	90
	Ultraflex™ Tracheobronchial Uncovered Stent System – Distal Release	91
Biopsy Forceps	Radial Jaw™ 4 Pulmonary Biopsy Forceps	91
Mini-Forceps	CoreDx™ Pulmonary Mini-Forceps	91
Cytology Brushes	Cellebriety™ Cytology Brushes – Bronchoscope	92
Airway Retrieval Baskets	Zero Tip™ Airway Retrieval Baskets	92
Dilatation Balloons	CRE™ Pulmonary Balloon Dilatation Catheter	93
Guidewires	Jagwire™ Pulmonary Guidewire	93
Inflation System	Alliance™ II Inflation System	93



Contents

Closure

Clipping Devices	Mantis™ Clip	95
Endoscopic Suturing	OverStitch™ Systems	95
	OverStitch™ Accessories	95
	X-Tack™	95
	X-Tack™ Accessories	95
Endoluminal Vacuum Therapy	Eso-SPONGE®	96
	Endo-SPONGE®	96
	EVT Vacuum Sources	96

Haemostasis

Clipping Devices	Resolution™ Clips	98
	Resolution 360™ Clips	98
	Resolution 360™ ULTRA Clips	99
	Speedband Superview Super 7™ Multiple Band Ligators	99
Bipolar Electrohaemostasis Catheters	Injection Gold Probe™ Electrohaemostasis/Injection Catheters	99
	Gold Probe™ Electrohaemostasis Catheters	100
Needles	Interject™ Sclerotherapy Needles	100
	Interject™ Clear Sclerotherapy Needles	101
Electrosurgical Bipolar Cable Adaptors	Electrosurgical Bipolar Cable Adaptor	101

Tissue Resection

EMR Devices for Upper GI	Captivator™ EMR	103
Electrosurgical Knife	ORISE™ ProKnife	103
Snares	Rotatable Single-Use Polypectomy Snares	103
	Sensation™ Single-Use Short Throw Snares	104
	Captiflex™ Single-Use Snares	104
	Captivator™ Single-Use Snares	105
	Captivator™ II Single-Use Snares	105
	Captivator™ COLD Single-Use Snare	106
Hot Biopsy Forceps	Radial Jaw™ 4 Hot Biopsy Forceps	106
Accessories	Active Cords	106
Retrieval Devices	RescueNet™	107
	Rescue™ Forceps	107
Polyp Traps	TrapEase™ Polyp Traps	107



Esophageal and Enteral Metal Stents

Esophageal Stents	WallFlex™ Fully Covered RMV Esophageal Longer Loop Stent System	109
	WallFlex™ Fully Covered RMV Esophageal Stent	109
	WallFlex™ Partially Covered Esophageal Stent System	110
	WallFlex™ Fully Covered Esophageal Stent	110
	Ultraflex™ Covered Esophageal NG Stent System – Proximal Release	111
	Ultraflex™ Covered Esophageal NG Stent System – Distal Release	111
	Ultraflex™ Uncovered Esophageal NG Stent System – Proximal Release	111
	Ultraflex™ Uncovered Esophageal NG Stent System – Distal Release	111
	Agile™ Fully Covered Esophageal Stent System – Removable	112
	Agile™ Partially Covered Esophageal Stent System	112
	Agile™ OTW Fully Covered Esophageal Stent System – Removable	112
	Agile™ OTW Partially Covered Esophageal Stent System	112
Colonic and Duodenal Stents	WallFlex™ Colonic Stent System with Anchor Lock Delivery System	113
	WallFlex™ Colonic Soft Stent System with Anchor Lock Delivery System	113
	WallFlex™ Duodenal Stent System with Anchor Lock Delivery System	114
	WallFlex™ Duodenal Soft Stent System with Anchor Lock Delivery System	114

Infection Prevention

Irrigation & Insufflation	Hydra™ Irrigation System	116
Cleaning Brushes	HedgeHog™ Cleaning Brushes	116
Valves	Orca™	116
	OrcaPod™	116
	Seal™	116
	Porpoise™	116
Scope Transport	CinchPad™ Procedure & Endoscope Transport Pad	117
Bite Blocks	Blox™ Bite Blocks	117
Sponges	Koala™ Sponge	117
	Compliance Endokits	117
Lubricating Gels	EndoGlide™ Lubricating Gels	117

Endobariatrics

Endoscopic Suturing	OverStitch™ Systems	119
Intragastric Balloon	BIB™ System	119
	Orbera365™	119
Overstitch Accessories	Overtube™	119
	Ensizor™	119
Intragastric Balloon Accessories	Extraction Kit	119

Alphabetical Index

Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Devices	68
Acquire™ S Endoscopic Ultrasound Fine Needle Biopsy (FNB) Devices	69
Acquire™ Pulmonary Endobronchial Ultrasound Fine Biopsy Needles – Olympus Scope Compatible	89
Active Cords	106
Advanix™ Pancreatic Single Stents – Pigtail	59-60
Advanix™ Pancreatic Single Stents – Straight	61-62
Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Pigtail	58
Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Straight	58
Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System	53
Advanix™ Stents – Singles	54-55
Agile™ Fully Covered Esophageal Stent System – Removable	112
Agile™ OTW Fully Covered Esophageal Stent System – Removable	112
Agile™ OTW Partially Covered Esophageal Stent System	112
Agile™ Partially Covered Esophageal Stent System	112
Alliance™ II Inflation System	81, 93
Amplatz Super Stiff™ Guidewire	73
AutoCap™ RX	51
Autolith™ Service Codes	86
Autolith™ Touch Biliary EHL System	86
Autotome™ Pro RX Cannulating Sphincterotomes	47
Autotome™ Pro RX Preloaded Sphincterotomes	47
Autotome™ RX Cannulating Sphincterotomes	47
BIB™ System	119
Biliary Accessories	42
Blox™ Bite Blocks	117
Captiflex™ Single-Use Snares	104
Captivator™ COLD Single-Use Snare	106
Captivator™ EMR	103
Captivator™ Single-Use Snares	105
Captivator™ II Single-Use Snares	105
Celebrity™ Cytology Brushes – Bronchoscope	92
CinchPad™ Procedure & Endoscope Transport Pad	117
Compliance Endokits	117
Contour™ ERCP Cannulas	38
CoreDx™ Pulmonary Mini-Forceps	91
CRE™ Fixed Wire Balloon Dilatation Catheters	80
CRE™ PRO Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Biliary) (including DASE indication)	79
CRE™ PRO Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic/Biliary) (including DASE indication)	78
CRE™ Pulmonary Balloon Dilatation Catheters	93
CRE™ RX Biliary Balloon Dilatation Catheters	50
CRE™ Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)	79
CRE™ Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)	78
Dreamtome™ RX Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	45
Dreamwire™ High Performance Guidewires	71
Electrosurgical Bipolar Cable Adaptor	101
Encore™ 26 Inflation Device	51
EndoGlide™ Lubricating Gels	117
Endo-SPONGE®	96
Ensizor™	119
Eso-SPONGE®	96



Alphabetical Index

EVT Vacuum Sources	96
EXALT™ Model D Service Codes	83
EXALT™ Model D Single-Use Duodenoscope	83
Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needles – Olympus Scope Compatible	89
Expect™ Slimline Handle Needles	68
Extraction Kit	119
Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Distal Injection	48
Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	48
Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Distal Injection	49
Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	49
Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Distal Injection	41
Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	41
Flexima™ Biliary Stents with Delivery System	57
Gold Probe™ Electrohaemostasis Catheters	100
Habib™ EndoHPB Bipolar Radiofrequency Catheter	66
HedgeHog™ Cleaning Brushes	116
Hot AXIOS™ Stent and Electrocautery Enhanced Delivery System	69
Hurricane™ RX Biliary Balloon Dilatation Catheters	50
Hydra™ Irrigation System	116
Hydra Jagwire™ High Performance Guidewires	71
Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	45
Inflation/Deflation Devices and Accessories	81
Injection Gold Probe™ Electrohaemostasis/Injection Catheters	99
Interject™ Clear Sclerotherapy Needles	101
Interject™ Sclerotherapy Needles	100
Jagtome™ Revolution RX Cannulating Sphincterotome (Autotome™ Preloaded with Jagwire Revolution High Performance Guidewire)	46
Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	46
Jagwire™ High Performance Guidewires	72
Jagwire™ Pulmonary Guidewire	93
Jagwire™ Revolution High Performance Guidewire	72
Koala™ Sponge	117
Mantis™ Clip	95
Microknife™ XL Triple-Lumen Needle-Knife	40
NaviFlex™ RX Delivery System	55
NaviFlex™ RX Pancreatic Delivery System & NaviFlex™ Pushers	63
NaviPro™ Hydrophilic Guidewires	73
Orbera365™	119
Orca™	116
OrcaPod™	116
ORISE™ ProKnife	103
OverStitch™ Accessories	95
OverStitch™ Systems	95
Overtube™	119
Porpoise™	116
Radial Jaw™ 4 Gastropediatric Biopsy Forceps	76
Radial Jaw™ 4 Hot Biopsy Forceps	106
Radial Jaw™ 4 Jumbo Biopsy Forceps	76
Radial Jaw™ 4 Large Capacity Biopsy Forceps	75
Radial Jaw™ 4 Pulmonary Biopsy Forceps	91
Radial Jaw™ 4 Standard Capacity Biopsy Forceps	75



Alphabetical Index

Rapid Exchange XL Cannula	44
RescueNet™	107
Rescue™ Forceps	107
Resolution 360™ Clips	98
Resolution 360™ ULTRA Clips	99
Resolution™ Clips	98
Rigiflex™ II Achalasia Balloon Dilators	80
Rotatable Single-Use Polypectomy Snares	103
RX Locking Devices	51
RX Biliary Plastic Stents with RX Delivery System	56
RX Cytology Brush	50
RX ERCP Cannulas	44
RX Needle-Knife XL Triple-Lumen Needle-Knife	48
Scivita Medical™ Capital Equipment	84, 88
Scivita Medical™ Disposable Equipment	85, 88
Scivita Medical™ Service Codes	85, 88
Seal™	116
Sensation™ Single-Use Short Throw Snares	104
Speedband Superview Super 7™ Multiple Band Ligators	99
SpyGlass™ Discover Digital Catheter	86
SpyGlass™ Discover Service Codes	86
SpyGlass™ DS Capital Equipment	84
SpyGlass™ DS Disposable Equipment	84
SpyGlass™ DS Service Codes	84
Tandem™ XL Triple-Lumen ERCP Cannula	38
TrapEase™ Polyp Traps	107
Trapezoid™ RX Wireguided Retrieval Baskets	49
TRUEtome™ Cannulating Sphincterotomes	40
TRUEtome™ Dreamwire™ Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	39
TRUEtome™ Hydra Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	39
TRUEtome™ Jagwire Revolution Cannulating Sphincterotomes with Pre-loaded Jagwire™ Revolution Guidewire	39
TRUEtome™ Jagwire Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	39
Ultraflex™ Covered Esophageal NG Stent System – Distal Release	111
Ultraflex™ Covered Esophageal NG Stent System – Proximal Release	111
Ultraflex™ Tracheobronchial Covered Stent System – Distal Release	90
Ultraflex™ Tracheobronchial Uncovered Stent System – Distal Release	91
Ultraflex™ Tracheobronchial Uncovered Stent System – Proximal Release	90
Ultraflex™ Uncovered Esophageal NG Stent System – Distal Release	111
Ultraflex™ Uncovered Esophageal NG Stent System – Proximal Release	111
Ultratome XL™ Triple-Lumen Sphincterotomes	41
WallFlex™ Biliary RX Fully Covered Stent System RMV (Permalume™ Covering, includes benign indication)	63
WallFlex™ Biliary RX Partially Covered Stent System (Permalume™ Covering)	64
WallFlex™ Biliary RX Uncovered Stent System	64
WallFlex™ Biliary RX Fully Covered Stent System	65
WallFlex™ Colonic Soft Stent System with Anchor Lock Delivery System	113
WallFlex™ Colonic Stent System with Anchor Lock Delivery System	113
WallFlex™ Duodenal Soft Stent System with Anchor Lock Delivery System	114
WallFlex™ Duodenal Stent System with Anchor Lock Delivery System	114
WallFlex™ Fully Covered Esophageal Stent	110
WallFlex™ Fully Covered RMV Esophageal Longer Loop Stent System	109



Alphabetical Index

WallFlex™ Fully Covered RMV Esophageal Stent	109
WallFlex™ Partially Covered Esophageal Stent System	110
X-Tack™	95
X-Tack™ Accessories	95
Zero Tip™ Airway Retrieval Baskets	92



Product Code Index

Order Number	Description	Page Number
Endo-SPONGE®		
5526510	Endo-SPONGE	96
5526520	Endo-SPONGE	96
5526530	Endo-SPONGE	96
Eso-SPONGE®		
5526540	Eso-SPONGE with 15 mm Overtube	96
5526550	Eso-SPONGE with 13 mm Overtube	96
EVT Vacuum Sources		
5526650	Endo-SPONGE® Pump	96
5526651	Endo-SPONGE® Pump Carrying Bag	96
5526653	Disposable Cannister 300 ml with tubing and Y-connector	96
5526661	Redyrob Trans Plus	96
19116	MV1 Mobile Programmable Low-Vacuum-System	96
18032	Collection Bottle MV1	96
18022	Bacterial Filter MV1	96
Extraction Kit		
B-40815	Extraction Kit	119
Compliance Endokits		
CIK-115-20	OP4 I-KIT 1.1 YellowSponge	117
CIK-122-20	OP3 I-KIT 1.1 YellowSponge	117
RescueNet™		
DGN-538-5	RescueNet Retrieval Device	107
Ensizor™		
ES26165-3	Ensizor Endoscopic Scissors Gastric Length	119
ES26135-3	Ensizor Endoscopic Scissors Colonic Length	119
Habib™ EndoHPB Bipolar Radiofrequency Catheter		
M00500070	Habib Endo HPB Bipolar Radiofrequency Catheter	66
5100	Adapter Cable Bipolar, 28 mm pins	66
5420	Adapter Cable Bipolar, Internal 4 mm pin	66
5700	Adapter Cable Bipolar, 22 mm pins	66
Koala™ Sponge		
M00501762	Koala Cleaning Sponge	117
SES-240-50	Koala Cleaning Sponge Imbibed w/ Pure	117
Porpoise™		
M00501891	Porpoise Air/Water Cleaning Adapter (For use with Olympus 160, 180 and 190 gastrointestinal endoscopes)	116
EndoGlide™ Lubricating Gels		
M00501902	EndoGlide+ 2.0 oz Lubricant Gel	117
M00501905	EndoGlide+ 1.1 oz Lubricant Gel	117
Electrosurgical Bipolar Cable Adaptors		
M00505200	Bipolar Cable Adaptor	101
M00505210	Bipolar Cable Adaptor	101
M00505220	Bipolar Cable Adaptor	101
OverStitch™ Systems		
M00505400	OverStitch Endoscopic Suturing System	95, 119
M00505420	OverStitch NXT Endoscopic Suture System	95, 119
M00505441	OverStitch 2-0 Polypropylene Suture	95, 119
M00505470	OverStitch Tissue Helix	95, 119
M00505480	OverStitch NXT Tissue Helix Pro	95, 119
M00505540	OverStitch Suture Cinch	95, 119



Product Code Index

Order Number	Description	Page Number
	X-Tack™	
M00505500	X-Tack Endoscopic HeliX Tacking System (Gastric)	95
M00505510	X-Tack Endoscopic HeliX Tacking System (Colonic)	95
	Overtube™	
M00505560	OverTube Endoscopic Access System	95, 119
	BIB™ System	
M00505680	BIB 6 Month Balloon	119
	Orbera365™	
M00505690	Orbera365 12 Month Balloon	119
AS2290718	Endoline Viper	119
AS2640718	Endoline Punc	119
	Orca™	
M00508001	Orca Fuji Air/Water and Suction Valves	116
	Trapezoid™ RX Wireguided Retrieval Baskets	
M00510860	Trapezoid RX Wireguided Retrieval Basket	49
M00510870	Trapezoid RX Wireguided Retrieval Basket	49
M00510880	Trapezoid RX Wireguided Retrieval Basket	49
M00510890	Trapezoid RX Wireguided Retrieval Basket	49
	Zero Tip™ Airway Retrieval Baskets	
M00513200	Zero Tip Airway Retrieval Basket	92
M00513210	Zero Tip Airway Retrieval Basket	92
	Radial Jaw™ 4 Large Capacity Biopsy Forceps	
M00513301	Radial Jaw 4 Large Capacity	75
M00513302	Radial Jaw 4 Large Capacity	75
M00513303	Radial Jaw 4 Large Capacity	75
M00513311	Radial Jaw 4 Large Capacity w/Needle	75
M00513312	Radial Jaw 4 Large Capacity w/Needle	75
M00513313	Radial Jaw 4 Large Capacity w/Needle	75
M00513321	Radial Jaw 4 Large Capacity	75
M00513322	Radial Jaw 4 Large Capacity	75
M00513323	Radial Jaw 4 Large Capacity	75
M00513331	Radial Jaw 4 Large Capacity w/Needle	75
M00513332	Radial Jaw 4 Large Capacity w/Needle	75
M00513333	Radial Jaw 4 Large Capacity w/Needle	75
	Radial Jaw™ 4 Jumbo Biopsy Forceps	
M00513343	Radial Jaw 4 Jumbo	76
M00513353	Radial Jaw 4 Jumbo w/Needle	76
M00513361	Radial Jaw 4 Jumbo	76
M00513362	Radial Jaw 4 Jumbo	76
M00513363	Radial Jaw 4 Jumbo	76
M00513371	Radial Jaw 4 Jumbo w/Needle	76
M00513372	Radial Jaw 4 Jumbo w/Needle	76
M00513373	Radial Jaw 4 Jumbo w/Needle	76
	Radial Jaw™ 4 Standard Capacity Biopsy Forceps	
M00513381	Radial Jaw 4 Standard Capacity	75
M00513382	Radial Jaw 4 Standard Capacity	75
M00513383	Radial Jaw 4 Standard Capacity	75
M00513391	Radial Jaw 4 Standard Capacity w/Needle	75



Product Code Index

Order Number	Description	Page Number
	Radial Jaw™ 4 Standard Capacity Biopsy Forceps (CONTINUED)	
M00513392	Radial Jaw 4 Standard Capacity w/Needle	75
M00513393	Radial Jaw 4 Standard Capacity w/Needle	75
M00513401	Radial Jaw 4 Standard Capacity	75
M00513402	Radial Jaw 4 Standard Capacity	75
M00513403	Radial Jaw 4 Standard Capacity	75
M00513411	Radial Jaw 4 Standard Capacity w/Needle	75
M00513412	Radial Jaw 4 Standard Capacity w/Needle	75
M00513413	Radial Jaw 4 Standard Capacity w/Needle	75
	Radial Jaw™ 4 Gastropediatric Biopsy Forceps	
M00513441	Radial Jaw 4 Gastropediatric	76
M00513443	Radial Jaw 4 Gastropediatric	76
M00513451	Radial Jaw 4 Gastropediatric w/Needle	76
M00513453	Radial Jaw 4 Gastropediatric w/Needle	76
	Ultraflex™ Uncovered Esophageal NG Stent System – Distal Release	
M00513700	Ultraflex Esophageal NG Stent System	111
M00513710	Ultraflex Esophageal NG Stent System	111
M00513720	Ultraflex Esophageal NG Stent System	111
	Ultraflex™ Covered Esophageal NG Stent System – Distal Release	
M00513730	Ultraflex Esophageal NG Stent System	111
M00513740	Ultraflex Esophageal NG Stent System	111
M00513750	Ultraflex Esophageal NG Stent System	111
M00514200	Ultraflex Esophageal NG Stent System	111
M00514210	Ultraflex Esophageal NG Stent System	111
	Ultraflex™ Uncovered Esophageal NG Stent System – Proximal Release	
M00513800	Ultraflex Esophageal NG Stent System	111
M00513810	Ultraflex Esophageal NG Stent System	111
M00513820	Ultraflex Esophageal NG Stent System	111
M00513830	Ultraflex Esophageal NG Stent System	111
M00514220	Ultraflex Esophageal NG Stent System	111
M00514230	Ultraflex Esophageal NG Stent System	111
	Ultraflex™ Covered Esophageal NG Stent System – Proximal Release	
M00513840	Ultraflex Esophageal NG Stent System	111
M00513850	Ultraflex Esophageal NG Stent System	111
M00513860	Ultraflex Esophageal NG Stent System	111
M00514240	Ultraflex Esophageal NG Stent System	111
M00514250	Ultraflex Esophageal NG Stent System	111
	Radial Jaw™ 4 Hot Biopsy Forceps	
M00515031	Radial Jaw 4 Hot Biopsy Forceps	106
M00515032	Radial Jaw 4 Hot Biopsy Forceps	106
M00515033	Radial Jaw 4 Hot Biopsy Forceps	106
	Jagwire™ Pulmonary Guidewire	
M00515171	Jagwire Pulmonary Guidewire	93
	Radial Jaw™ 4 Pulmonary Biopsy Forceps	
M00515181	Radial Jaw 4 Pulmonary Standard Capacity	91
M00515182	Radial Jaw 4 Pulmonary Standard Capacity	91
M00515191	Radial Jaw 4 Pulmonary Standard Capacity w/Needle	91
M00515192	Radial Jaw 4 Pulmonary Standard Capacity w/Needle	91



Product Code Index

Order Number	Description	Page Number
	Radial Jaw™ 4 Pulmonary Biopsy Forceps (CONTINUED)	
M00515201	Radial Jaw 4 Pulmonary Large Capacity	91
M00515202	Radial Jaw 4 Pulmonary Large Capacity	91
	CoreDx™ Pulmonary Mini-Forceps	
M00515220	Pulmonary Mini-Forceps	91
	Cellebriety™ Cytology Brushes – Bronchoscope	
M00516001	Cellebriety Cytology Brush	92
M00516011	Cellebriety Cytology Brush	92
M00516071	Cellebriety Cytology Brush	92
M00516151	Cellebriety Cytology Brush	92
	WallFlex™ Fully Covered RMV Esophageal Stent	
M00516210	WallFlex Fully Covered RMV Esophageal Stent	109
M00516220	WallFlex Fully Covered RMV Esophageal Stent	109
M00516230	WallFlex Fully Covered RMV Esophageal Stent	109
M00516240	WallFlex Fully Covered RMV Esophageal Stent	109
M00516250	WallFlex Fully Covered RMV Esophageal Stent	109
M00516260	WallFlex Fully Covered RMV Esophageal Stent	109
	WallFlex™ Fully Covered RMV Esophageal Longer Loop Stent System	
M00516270	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	109
M00516280	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	109
M00516290	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	109
M00516300	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	109
M00516310	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	109
M00516320	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	109
	WallFlex™ Fully Covered Esophageal Stent	
M00516700	WallFlex Fully Covered Esophageal Stent	110
M00516710	WallFlex Fully Covered Esophageal Stent	110
M00516720	WallFlex Fully Covered Esophageal Stent	110
M00516730	WallFlex Fully Covered Esophageal Stent	110
M00516740	WallFlex Fully Covered Esophageal Stent	110
M00516750	WallFlex Fully Covered Esophageal Stent	110
	WallFlex™ Partially Covered Esophageal Stent System	
M00516900	WallFlex Partially Covered Esophageal Stent System	110
M00516910	WallFlex Partially Covered Esophageal Stent System	110
M00516920	WallFlex Partially Covered Esophageal Stent System	110
M00516930	WallFlex Partially Covered Esophageal Stent System	110
M00516940	WallFlex Partially Covered Esophageal Stent System	110
M00516950	WallFlex Partially Covered Esophageal Stent System	110
	Agile™ Partially Covered Esophageal Stent System	
M00517200	Agile Partially Covered Esophageal Stent System	112
M00517210	Agile Partially Covered Esophageal Stent System	112
M00517220	Agile Partially Covered Esophageal Stent System	112
M00517230	Agile Partially Covered Esophageal Stent System	112
M00517240	Agile Partially Covered Esophageal Stent System	112
M00517250	Agile Partially Covered Esophageal Stent System	112
M00517260	Agile Partially Covered Esophageal Stent System	112
M00517270	Agile Partially Covered Esophageal Stent System	112



Product Code Index

Order Number	Description	Page Number
Agile™ Fully Covered Esophageal Stent System – Removable		
M00517580	Agile Fully Covered Esophageal Stent System – Removable	112
M00517590	Agile Fully Covered Esophageal Stent System – Removable	112
M00517600	Agile Fully Covered Esophageal Stent System – Removable	112
M00517610	Agile Fully Covered Esophageal Stent System – Removable	112
M00517620	Agile Fully Covered Esophageal Stent System – Removable	112
M00517630	Agile Fully Covered Esophageal Stent System – Removable	112
M00517640	Agile Fully Covered Esophageal Stent System – Removable	112
M00517650	Agile Fully Covered Esophageal Stent System – Removable	112
Agile™ OTW Partially Covered Esophageal Stent System		
M00517710	Agile OTW Partially Covered Esophageal Stent System	112
M00517720	Agile OTW Partially Covered Esophageal Stent System	112
M00517730	Agile OTW Partially Covered Esophageal Stent System	112
M00517740	Agile OTW Partially Covered Esophageal Stent System	112
Agile™ OTW Fully Covered Esophageal Stent System – Removable		
M00517790	Agile OTW Fully Covered Esophageal Stent System – Removable	112
M00517800	Agile OTW Fully Covered Esophageal Stent System – Removable	112
M00517810	Agile OTW Fully Covered Esophageal Stent System – Removable	112
M00517820	Agile OTW Fully Covered Esophageal Stent System – Removable	112
Interject™ Sclerotherapy Needles		
M00518151	Interject Sclerotherapy Needle	100
M00518161	Interject Sclerotherapy Needle	100
M00518251	Interject Sclerotherapy Needle	100
M00518261	Interject Sclerotherapy Needle	100
M00518351	Interject Sclerotherapy Needle	100
M00518361	Interject Sclerotherapy Needle	100
Interject™ Clear Sclerotherapy Needles		
M00518111	Interject Sclerotherapy Needle	101
M00518301	Interject Sclerotherapy Needle	101
M00518311	Interject Sclerotherapy Needle	101
ORISE™ ProKnife		
M00519351	1.5mm ORISE ProKnife	103
M00519361	2.0mm ORISE ProKnife	103
M00519371	3.0mm ORISE ProKnife	103
Resolution 360™ Clips		
M00521230	Resolution 360 Clip	98
M00521231	Resolution 360 Clip	98
M00521232	Resolution 360 Clip	98
M00521233	Resolution 360 Clip	98
M00521240	Resolution 360 Clip	98
M00521241	Resolution 360 Clip	98
M00521242	Resolution 360 Clip	98
M00521243	Resolution 360 Clip	98
Resolution 360™ ULTRA Clips		
M00521400	Resolution 360 ULTRA Clip	99
M00521401	Resolution 360 ULTRA Clip	99
M00521402	Resolution 360 ULTRA Clip	99
M00521403	Resolution 360 ULTRA Clip	99



Product Code Index

Order Number	Description	Page Number
Mantis™ Clips		
M00521420	Mantis Clip	95
M00521421	Mantis Clip	95
M00521422	Mantis Clip	95
M00521423	Mantis Clip	95
Resolution™ Clips		
M00522600	Resolution Clip	98
M00522601	Resolution Clip	98
M00522602	Resolution Clip	98
M00522610	Resolution Clip	98
M00522611	Resolution Clip	98
M00522612	Resolution Clip	98
Rescue™ Forceps		
M00527451	Rescue Alligator Grasping Forceps	107
M00527551	Rescue Alligator Grasping Forceps Long	107
M00527651	Rescue Rat Tooth Grasping Forceps	107
M00527751	Rescue Rat Tooth Grasping Forceps – Pediatric	107
M00538351	Rescue Rat Tooth / Alligator Grasping Forceps	107
Contour™ ERCP Cannulas		
M00530850	Contour ERCP Cannula	38
M00530860	Contour ERCP Cannula	38
M00530870	Contour ERCP Cannula	38
M00530880	Contour ERCP Cannula	38
M00530890	Contour ERCP Cannula	38
Advanix™ Stents – Singles		
M00532160	Advanix Stent	55
M00532180	Advanix Stent	55
M00532190	Advanix Stent	55
M00532210	Advanix Stent	55
M00532220	Advanix Stent	55
M00532230	Advanix Stent	55
M00532240	Advanix Stent	55
M00532260	Advanix Stent	55
M00532270	Advanix Stent	55
M00532290	Advanix Stent	55
M00532300	Advanix Stent	55
M00532310	Advanix Stent	55
Microknife™ XL Triple-Lumen Needle-Knife		
M00532810	Microknife XL Triple-Lumen Needle-Knife	40
Advanix™ Stents – Singles		
M00532870	Advanix Stent	54
M00532880	Advanix Stent	54
M00532890	Advanix Stent	54
M00532900	Advanix Stent	54
M00532910	Advanix Stent	54
M00532920	Advanix Stent	54
M00532930	Advanix Stent	54
M00532940	Advanix Stent	54



Product Code Index

Order Number	Description	Page Number
Advanix™ Stents – Singles (CONTINUED)		
M00532950	Advanix Stent	54
M00532960	Advanix Stent	54
M00532970	Advanix Stent	54
M00532980	Advanix Stent	54
M00532990	Advanix Stent	54
M00533000	Advanix Stent	54
M00533010	Advanix Stent	54
M00533020	Advanix Stent	54
M00533030	Advanix Stent	54
M00533040	Advanix Stent	54
M00533230	Advanix Stent	54
M00533240	Advanix Stent	54
M00533250	Advanix Stent	54
M00533260	Advanix Stent	54
M00533270	Advanix Stent	54
M00533280	Advanix Stent	54
M00533290	Advanix Stent	54
M00533300	Advanix Stent	54
M00533310	Advanix Stent	54
M00533320	Advanix Stent	54
M00533330	Advanix Stent	54
M00533340	Advanix Stent	54
M00533350	Advanix Stent	54
M00533360	Advanix Stent	54
M00533370	Advanix Stent	54
M00533380	Advanix Stent	54
M00533390	Advanix Stent	54
M00533400	Advanix Stent	54
NaviFlex™ RX Delivery System		
M00533560	RX Delivery System	55
M00533570	RX Delivery System	55
M00533580	RX Delivery System	55
Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System		
M00534200	Advanix Pre-loaded RX Stent	53
M00534210	Advanix Pre-loaded RX Stent	53
M00534220	Advanix Pre-loaded RX Stent	53
M00534230	Advanix Pre-loaded RX Stent	53
M00534240	Advanix Pre-loaded RX Stent	53
M00534250	Advanix Pre-loaded RX Stent	53
M00534260	Advanix Pre-loaded RX Stent	53
M00534270	Advanix Pre-loaded RX Stent	53
M00534280	Advanix Pre-loaded RX Stent	53
M00534290	Advanix Pre-loaded RX Stent	53
M00534300	Advanix Pre-loaded RX Stent	53
M00534310	Advanix Pre-loaded RX Stent	53
M00534320	Advanix Pre-loaded RX Stent	53
M00534330	Advanix Pre-loaded RX Stent	53



Product Code Index

Order Number	Description	Page Number
	Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System (CONTINUED)	
M00534340	Advanix Pre-loaded RX Stent	53
M00534350	Advanix Pre-loaded RX Stent	53
M00534360	Advanix Pre-loaded RX Stent	53
M00534370	Advanix Pre-loaded RX Stent	53
M00534560	Advanix Pre-loaded RX Stent	53
M00534570	Advanix Pre-loaded RX Stent	53
M00534580	Advanix Pre-loaded RX Stent	53
M00534590	Advanix Pre-loaded RX Stent	53
M00534600	Advanix Pre-loaded RX Stent	53
M00534610	Advanix Pre-loaded RX Stent	53
M00534620	Advanix Pre-loaded RX Stent	53
M00534630	Advanix Pre-loaded RX Stent	53
M00534640	Advanix Pre-loaded RX Stent	53
M00534650	Advanix Pre-loaded RX Stent	53
M00534660	Advanix Pre-loaded RX Stent	53
M00534670	Advanix Pre-loaded RX Stent	53
M00534680	Advanix Pre-loaded RX Stent	53
M00534690	Advanix Pre-loaded RX Stent	53
M00534700	Advanix Pre-loaded RX Stent	53
M00534710	Advanix Pre-loaded RX Stent	53
M00534720	Advanix Pre-loaded RX Stent	53
M00534730	Advanix Pre-loaded RX Stent	53
	NaviFlex™ RX Pancreatic Delivery System & NaviFlex™ Pushers	
M00535000	Pancreatic Stent RX Pusher	63
M00535010	Pancreatic Stent RX Pusher	63
M00535020	Pancreatic Stent Pusher	63
M00535030	Pancreatic Stent Pusher	63
M00535060	NaviFlex RX Pancreatic Stent Delivery System	63
M00535080	NaviFlex RX Pancreatic Stent Delivery System	63
	Tandem™ XL Triple-Lumen ERCP Cannula	
M00535700	Tandem XL Triple-Lumen ERCP Cannula	38
	Ultratome XL™ Triple-Lumen Sphincterotomes	
M00535900	Ultratome XL Short Nose	41
M00535910	Ultratome XL Long Nose	41
M00535920	Ultratome XL Short Nose	41
M00535930	Ultratome XL Long Nose	41
	Advanix™ Pancreatic Single Stents – Straight	
M00536000	Advanix Pancreatic Stent	61
M00536010	Advanix Pancreatic Stent	61
M00536020	Advanix Pancreatic Stent	61
M00536030	Advanix Pancreatic Stent	61
M00536040	Advanix Pancreatic Stent	61
M00536050	Advanix Pancreatic Stent	61
M00536060	Advanix Pancreatic Stent	61
M00536070	Advanix Pancreatic Stent	61
M00536080	Advanix Pancreatic Stent	61
M00536090	Advanix Pancreatic Stent	61



Product Code Index

Order Number	Description	Page Number
Advanix™ Pancreatic Single Stents – Straight (CONTINUED)		
M00536100	Advanix Pancreatic Stent	61
M00536110	Advanix Pancreatic Stent	61
M00536120	Advanix Pancreatic Stent	61
M00536130	Advanix Pancreatic Stent	61
M00536140	Advanix Pancreatic Stent	61
M00536150	Advanix Pancreatic Stent	61
M00536160	Advanix Pancreatic Stent	61
M00536170	Advanix Pancreatic Stent	61
M00536180	Advanix Pancreatic Stent	61
M00536190	Advanix Pancreatic Stent	61
M00536200	Advanix Pancreatic Stent	61
M00536210	Advanix Pancreatic Stent	61
M00536220	Advanix Pancreatic Stent	61
M00536230	Advanix Pancreatic Stent	61
M00536240	Advanix Pancreatic Stent	61
M00536250	Advanix Pancreatic Stent	61
M00536260	Advanix Pancreatic Stent	61
M00536270	Advanix Pancreatic Stent	61
M00536280	Advanix Pancreatic Stent	61
M00536290	Advanix Pancreatic Stent	61
M00536300	Advanix Pancreatic Stent	61
M00536310	Advanix Pancreatic Stent	61
M00536320	Advanix Pancreatic Stent	61
M00536330	Advanix Pancreatic Stent	61
M00536340	Advanix Pancreatic Stent	61
M00536350	Advanix Pancreatic Stent	61
M00536360	Advanix Pancreatic Stent	62
M00536370	Advanix Pancreatic Stent	62
M00536380	Advanix Pancreatic Stent	62
M00536390	Advanix Pancreatic Stent	62
M00536400	Advanix Pancreatic Stent	62
M00536410	Advanix Pancreatic Stent	62
M00536420	Advanix Pancreatic Stent	62
M00536430	Advanix Pancreatic Stent	62
M00536440	Advanix Pancreatic Stent	62
M00536450	Advanix Pancreatic Stent	62
M00536460	Advanix Pancreatic Stent	62
M00536470	Advanix Pancreatic Stent	62
M00536480	Advanix Pancreatic Stent	62
M00536490	Advanix Pancreatic Stent	62
M00536500	Advanix Pancreatic Stent	62
M00536510	Advanix Pancreatic Stent	62
M00536520	Advanix Pancreatic Stent	62
M00536530	Advanix Pancreatic Stent	62
M00536540	Advanix Pancreatic Stent	62
M00536550	Advanix Pancreatic Stent	62
M00536560	Advanix Pancreatic Stent	62



Product Code Index

Order Number	Description	Page Number
M00536570	Advanix Pancreatic Stent	62
M00536580	Advanix Pancreatic Stent	62
M00536590	Advanix Pancreatic Stent	62
M00536600	Advanix Pancreatic Stent	62
M00536610	Advanix Pancreatic Stent	62
M00536620	Advanix Pancreatic Stent	62
M00536630	Advanix Pancreatic Stent	62
M00536640	Advanix Pancreatic Stent	62
M00536650	Advanix Pancreatic Stent	62
M00536660	Advanix Pancreatic Stent	62
M00536670	Advanix Pancreatic Stent	62
M00536680	Advanix Pancreatic Stent	62
M00536690	Advanix Pancreatic Stent	62
M00536700	Advanix Pancreatic Stent	62
M00536710	Advanix Pancreatic Stent	62
M00536720	Advanix Pancreatic Stent	62
M00536730	Advanix Pancreatic Stent	62
Advanix™ Pancreatic Single Stents – Pigtail		
M00536740	Advanix Pancreatic Stent	59
M00536750	Advanix Pancreatic Stent	59
M00536760	Advanix Pancreatic Stent	59
M00536770	Advanix Pancreatic Stent	59
M00536780	Advanix Pancreatic Stent	59
M00536790	Advanix Pancreatic Stent	59
M00536800	Advanix Pancreatic Stent	59
M00536810	Advanix Pancreatic Stent	59
M00536820	Advanix Pancreatic Stent	59
M00536830	Advanix Pancreatic Stent	59
M00536840	Advanix Pancreatic Stent	59
M00536850	Advanix Pancreatic Stent	59
M00536860	Advanix Pancreatic Stent	59
M00536870	Advanix Pancreatic Stent	59
M00536880	Advanix Pancreatic Stent	59
M00536890	Advanix Pancreatic Stent	59
M00536900	Advanix Pancreatic Stent	59
M00536910	Advanix Pancreatic Stent	59
M00536920	Advanix Pancreatic Stent	59
M00536930	Advanix Pancreatic Stent	59
M00536940	Advanix Pancreatic Stent	59
M00536950	Advanix Pancreatic Stent	59
M00536960	Advanix Pancreatic Stent	59
M00536970	Advanix Pancreatic Stent	59
M00536980	Advanix Pancreatic Stent	59
M00536990	Advanix Pancreatic Stent	59
M00537000	Advanix Pancreatic Stent	59
M00537010	Advanix Pancreatic Stent	59
M00537020	Advanix Pancreatic Stent	59
M00537030	Advanix Pancreatic Stent	59



Product Code Index

Order Number	Description	Page Number
	Advanix™ Pancreatic Single Stents – Pigtail (CONTINUED)	
M00537040	Advanix Pancreatic Stent	59
M00537050	Advanix Pancreatic Stent	59
M00537060	Advanix Pancreatic Stent	59
M00537070	Advanix Pancreatic Stent	59
M00537080	Advanix Pancreatic Stent	59
M00537090	Advanix Pancreatic Stent	59
M00537100	Advanix Pancreatic Stent	60
M00537110	Advanix Pancreatic Stent	60
M00537120	Advanix Pancreatic Stent	60
M00537130	Advanix Pancreatic Stent	60
M00537140	Advanix Pancreatic Stent	60
M00537150	Advanix Pancreatic Stent	60
M00537160	Advanix Pancreatic Stent	60
M00537170	Advanix Pancreatic Stent	60
M00537180	Advanix Pancreatic Stent	60
M00537190	Advanix Pancreatic Stent	60
M00537200	Advanix Pancreatic Stent	60
M00537210	Advanix Pancreatic Stent	60
M00537220	Advanix Pancreatic Stent	60
M00537230	Advanix Pancreatic Stent	60
M00537240	Advanix Pancreatic Stent	60
M00537250	Advanix Pancreatic Stent	60
M00537260	Advanix Pancreatic Stent	60
M00537270	Advanix Pancreatic Stent	60
M00537280	Advanix Pancreatic Stent	60
M00537290	Advanix Pancreatic Stent	60
M00537300	Advanix Pancreatic Stent	60
M00537310	Advanix Pancreatic Stent	60
M00537320	Advanix Pancreatic Stent	60
M00537330	Advanix Pancreatic Stent	60
M00537340	Advanix Pancreatic Stent	60
	Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Pigtail	
M00537350	Advanix Pancreatic Stent Kit	58
M00537360	Advanix Pancreatic Stent Kit	58
M00537370	Advanix Pancreatic Stent Kit	58
M00537380	Advanix Pancreatic Stent Kit	58
M00537390	Advanix Pancreatic Stent Kit	58
M00537400	Advanix Pancreatic Stent Kit	58
M00537410	Advanix Pancreatic Stent Kit	58
M00537420	Advanix Pancreatic Stent Kit	58
M00537430	Advanix Pancreatic Stent Kit	58
M00537440	Advanix Pancreatic Stent Kit	58
M00537450	Advanix Pancreatic Stent Kit	58
	Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Straight	
M00537460	Advanix Pancreatic Stent Kit	58
M00537470	Advanix Pancreatic Stent Kit	58
M00537480	Advanix Pancreatic Stent Kit	58



Product Code Index

Order Number	Description	Page Number
	Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Straight (CONTINUED)	
M00537490	Advanix Pancreatic Stent Kit	58
M00537500	Advanix Pancreatic Stent Kit	58
M00537510	Advanix Pancreatic Stent Kit	58
M00537520	Advanix Pancreatic Stent Kit	58
M00537530	Advanix Pancreatic Stent Kit	58
M00537540	Advanix Pancreatic Stent Kit	58
M00537550	Advanix Pancreatic Stent Kit	58
M00537560	Advanix Pancreatic Stent Kit	58
	Flexima™ Biliary Stents with Delivery System	
M00539200	Flexima Stent with Delivery System	57
M00539210	Flexima Stent with Delivery System	57
M00539220	Flexima Stent with Delivery System	57
M00539230	Flexima Stent with Delivery System	57
M00539240	Flexima Stent with Delivery System	57
M00539250	Flexima Stent with Delivery System	57
M00539260	Flexima Stent with Delivery System	57
M00539270	Flexima Stent with Delivery System	57
M00539280	Flexima Stent with Delivery System	57
M00539290	Flexima Stent with Delivery System	57
M00539300	Flexima Stent with Delivery System	57
M00539310	Flexima Stent with Delivery System	57
M00539320	Flexima Stent with Delivery System	57
M00539330	Flexima Stent with Delivery System	57
M00539340	Flexima Stent with Delivery System	57
M00539350	Flexima Stent with Delivery System	57
M00539360	Flexima Stent with Delivery System	57
M00539370	Flexima Stent with Delivery System	57
M00539380	Flexima Stent with Delivery System	57
M00539390	Flexima Stent with Delivery System	57
	Speedband Superview Super 7™ Multiple Band Ligators	
M00542250	Super 7 Band Ligator	99
M00542251	Super 7 Band Ligator	99
M00542253	Super 7 Band Ligator	99
	EXALT™ Model D Single-Use Duodenoscope	
M0054242CE0	EXALT Model D Single-Use Duodenoscope	83
M0054242CE1	EXALT Model D Single-Use Duodenoscope	83
M00542430	EXALT Controller	83
M00542460	EXALT Controller	83
	SpyGlass™ Discover Digital Catheter	
M00542950	SpyGlass Discover Controller	86
	RX Cytology Brush	
M00545000	RX Cytology Brush	50
	Autocap™ RX	
M00545101	AutoCap RX Olympus Compatible (TJF-160, TJF-180, TJF-190 and JF-260V Series Duodenoscopes)	51
M00545131	AutoCap RX EXALT™ Model D Compatible Integrated Biopsy	51



Product Code Index

Order Number	Description	Page Number
Autotome™ RX Cannulating Sphincterotomes		
M00545150	Autotome RX 49	47
M00545160	Autotome RX 49	47
M00545170	Autotome RX 44	47
M00545180	Autotome RX 44	47
M00545190	Autotome RX 39	47
M00545200	Autotome RX 39	47
RX Locking Devices		
M00545211	RX Locking Device – Olympus / FujiFilm Brand Endoscope Compatible	51
M00545261	RX Locking Device and Biopsy Cap System – Olympus / FujiFilm G5 Series Brand Endoscope Compatible	51
M00545271	RX Locking Device and Biopsy Cap System – Pentax Brand Endoscope Compatible	51
M00545281	RX Locking Device and Biopsy Cap System – FujiFilm Type 3 Series Brand Endoscope Compatible	51
RX Biliary Plastic Stents with RX Delivery System		
M00545550	RX Biliary Plastic Stent with RX Delivery System	56
M00545560	RX Biliary Plastic Stent with RX Delivery System	56
M00545570	RX Biliary Plastic Stent with RX Delivery System	56
M00545580	RX Biliary Plastic Stent with RX Delivery System	56
M00545590	RX Biliary Plastic Stent with RX Delivery System	56
M00545600	RX Biliary Plastic Stent with RX Delivery System	56
M00545610	RX Biliary Plastic Stent with RX Delivery System	56
M00545620	RX Biliary Plastic Stent with RX Delivery System	56
M00545630	RX Biliary Plastic Stent with RX Delivery System	56
M00545640	RX Biliary Plastic Stent with RX Delivery System	56
M00545650	RX Biliary Plastic Stent with RX Delivery System	56
M00545660	RX Biliary Plastic Stent with RX Delivery System	56
M00545670	RX Biliary Plastic Stent with RX Delivery System	56
M00545680	RX Biliary Plastic Stent with RX Delivery System	56
M00545690	RX Biliary Plastic Stent with RX Delivery System	56
M00545700	RX Biliary Plastic Stent with RX Delivery System	56
M00545710	RX Biliary Plastic Stent with RX Delivery System	56
M00545720	RX Biliary Plastic Stent with RX Delivery System	56
M00545730	RX Biliary Plastic Stent with RX Delivery System	56
M00545740	RX Biliary Plastic Stent with RX Delivery System	56
RX ERCP Cannulas		
M00545760	RX ERCP Cannula Standard Tip	44
M00545770	RX ERCP Cannula Tapered Tip	44
M00545780	RX ERCP Cannula Ball Tip	44
Rapid Exchange XL Cannula		
M00545820	Rapid Exchange XL Cannula	44
RX Needle-Knife XL Triple-Lumen Needle-Knife		
M00545840	RX Needle-Knife XL Triple-Lumen Needle-Knife	48
Hurricane™ RX Biliary Dilatation Balloons		
M00545890	Hurricane RX Dilatation Balloon	50
M00545900	Hurricane RX Dilatation Balloon	50
M00545910	Hurricane RX Dilatation Balloon	50
M00545920	Hurricane RX Dilatation Balloon	50
M00545930	Hurricane RX Dilatation Balloon	50
M00545940	Hurricane RX Dilatation Balloon	50



Product Code Index

Order Number	Description	Page Number
	Hurricane™ RX Biliary Dilatation Balloons (CONTINUED)	
M00545950	Hurricane RX Dilatation Balloon	50
M00545960	Hurricane RX Dilatation Balloon	50
	SpyGlass™ DS Direct Visualisation System	
M00546470	SpyBite™ Max Biopsy Forceps	84
M00546550	SpyGlass™ Retrieval Basket	84
M00546560	SpyGlass™ Retrieval Snare	84
M00546610	SpyScope™ DS II Access and Delivery Catheter	84
M00546650	SpyGlass DS Digital Controller	84
	Autolith™ Touch Biliary EHL System	
M00546620	1.9Fr, 375cm Biliary EHL Probe	86
M00546680	Autolith Touch EHL Generator	86
M00546750	Autolith Touch Extender Cable	86
M00546760	Autolith Touch Foot Pedal	86
	SpyGlass™ Discover Digital Catheter	
M00546780	SpyGlass Discover Digital Catheter	86
	Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Distal Injection	
M00547000	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	48
M00547010	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	48
M00547020	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	48
	Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	
M00547030	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	48
M00547040	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	48
M00547050	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	48
	Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Distal Injection	
M00547100	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	41
M00547110	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	41
M00547120	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	41
	Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	
M00547130	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	41
M00547140	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	41
M00547150	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	41
	Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Distal Injection	
M00547300	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	49
M00547310	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	49
M00547320	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	49
	Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	
M00547330	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	49
M00547340	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	49
M00547350	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	49
	Amplatz Super Stiff™ Guidewire	
M00550090	Amplatz Super Stiff Guidewire	73
	CRE™ Pulmonary Balloon Dilatation Catheter	
M00550300	CRE Pulmonary Balloon Dilatation Catheter	93
M00550310	CRE Pulmonary Balloon Dilatation Catheter	93
M00550320	CRE Pulmonary Balloon Dilatation Catheter	93
M00550330	CRE Pulmonary Balloon Dilatation Catheter	93
M00550340	CRE Pulmonary Balloon Dilatation Catheter	93
M00550350	CRE Pulmonary Balloon Dilatation Catheter	93



Product Code Index

Order Number	Description	Page Number
Alliance™ II Inflation System		
M00550601	Alliance II Single-Use Syringe/Gauge Assembly	81, 93
M00550620	Alliance II Integrated Inflation/Lithotripsy Device	81, 93
Inflation/Deflation Devices and Accessories		
M00550630	CRE™ SteriFlate Disposable Inflation Device	81
M00553200	Pneumatic Inflator	81
Acquire™ Pulmonary Endobronchial Ultrasound Fine Biopsy Needles – Olympus Scope Compatible		
M00552350	22ga Acquire Pulmonary – Olympus scope compatible	89
M00552351	22ga Acquire Pulmonary – Olympus scope compatible	89
M00552400	22ga Acquire Pulmonary – Olympus scope compatible	89
M00552401	25ga Acquire Pulmonary – Olympus scope compatible	89
Hot AXIOS™ Stent and Electrocautery Enhanced Delivery System		
M00553520	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	69
M00553530	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	69
M00553540	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	69
M00553550	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	69
M00553560	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	69
M00553570	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	69
Rigiflex™ II Achalasia Balloon Dilators		
M00554500	Rigiflex II Achalasia Balloon Dilator	80
M00554510	Rigiflex II Achalasia Balloon Dilator	80
M00554520	Rigiflex II Achalasia Balloon Dilator	80
Expect™ Slimline Handle Needles		
M00555500	Expect Slimline Handle Needle	68
M00555501	Expect Slimline Handle Needle	68
M00555510	Expect Slimline Handle Needle	68
M00555511	Expect Slimline Handle Needle	68
M00555520	Expect Slimline Handle Needle	68
M00555521	Expect Slimline Handle Needle	68
M00555530	Expect Slimline Handle Needle	68
M00555531	Expect Slimline Handle Needle	68
Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Devices		
M00555540	22ga Acquire EUS-FNB Needle	68
M00555541	22ga Acquire EUS-FNB Needle	68
M00555560	25ga Acquire EUS-FNB Needle	68
M00555561	25ga Acquire EUS-FNB Needle	68
M00555580	19ga Acquire FNB Needle, Flexible	68
Acquire™ S Endoscopic Ultrasound Fine Needle Biopsy (FNB) Devices		
M00555640	22ga Acquire S EUS-FNB Needle	69
M00555641	22ga Acquire S EUS-FNB Needle	69
M00555660	25ga Acquire S EUS-FNB Needle	69
M00555661	25ga Acquire S EUS-FNB Needle	69
WallFlex™ Colonic Soft Stent System with Anchor Lock Delivery System		
M00555700	WallFlex Colonic Soft Stent	113
M00555710	WallFlex Colonic Soft Stent	113
M00555720	WallFlex Colonic Soft Stent	113
M00555730	WallFlex Colonic Soft Stent	113
M00555740	WallFlex Colonic Soft Stent	113



Product Code Index

Order Number	Description	Page Number
	WallFlex™ Colonic Soft Stent System with Anchor Lock Delivery System (CONTINUED)	
M00555750	WallFlex Colonic Soft Stent	113
M00555760	WallFlex Colonic Soft Stent	113
M00555770	WallFlex Colonic Soft Stent	113
M00555780	WallFlex Colonic Soft Stent	113
	WallFlex™ Duodenal Soft Stent System with Anchor Lock Delivery System	
M00555800	WallFlex Duodenal Soft Stent	114
M00555810	WallFlex Duodenal Soft Stent	114
M00555820	WallFlex Duodenal Soft Stent	114
M00555830	WallFlex Duodenal Soft Stent	114
M00555840	WallFlex Duodenal Soft Stent	114
M00555850	WallFlex Duodenal Soft Stent	114
M00555860	WallFlex Duodenal Soft Stent	114
M00555870	WallFlex Duodenal Soft Stent	114
M00555880	WallFlex Duodenal Soft Stent	114
	Hydra Jagwire™ High Performance Guidewires	
M00556001	Hydra Jagwire Standard	71
M00556011	Hydra Jagwire Standard	71
M00556021	Hydra Jagwire Stiff Shaft	71
M00556031	Hydra Jagwire Stiff Shaft	71
M00556041	Hydra Jagwire Standard	71
M00556051	Hydra Jagwire Standard	71
M00556061	Hydra Jagwire Stiff Shaft	71
M00556071	Hydra Jagwire Stiff Shaft	71
	Dreamwire™ High Performance Guidewires	
M00556101	Dreamwire Standard	71
M00556111	Dreamwire Standard	71
M00556121	Dreamwire Stiff Shaft	71
M00556131	Dreamwire Stiff Shaft	71
M00556141	Dreamwire Standard	71
M00556151	Dreamwire Standard	71
M00556161	Dreamwire Stiff Shaft	71
M00556171	Dreamwire Stiff Shaft	71
	NaviPro™ Hydrophilic Guidewires	
M00556181	NaviPro Standard	73
M00556191	NaviPro Standard	73
M00556201	NaviPro Standard	73
M00556211	NaviPro Standard	73
M00556221	NaviPro Standard	73
M00556231	NaviPro Standard	73
M00556241	NaviPro Stiff Shaft	73
M00556251	NaviPro Stiff Shaft	73
	Jagwire™ High Performance Guidewires	
M00556461	Jagwire Guidewire ST	72
M00556481	Jagwire Guidewire ST Stiff	72
M00556491	Jagwire Guidewire ANG Stiff	72
M00556511	Jagwire Guidewire ANG	72
M00556521	Jagwire Guidewire ST Stiff	72



Product Code Index

Order Number	Description	Page Number
Jagwire™ High Performance Guidewires (CONTINUED)		
M00556531	Jagwire Guidewire ANG Stiff	72
M00556561	Jagwire Guidewire ST	72
M00556581	Jagwire Guidewire ST	72
M00556591	Jagwire Guidewire ANG	72
M00556601	Jagwire Guidewire ST Stiff	72
M00556611	Jagwire Guidewire ANG Stiff	72
M00556621	Jagwire Guidewire ST	72
M00557661	Jagwire Guidewire ST Stiff	72
M00557651	Jagwire Guidewire ANG	72
M00556661	Jagwire Guidewire ANG Stiff	72
M00556841	Extendable Jagwire Guidewire ST	72
M00556851	Extendable Jagwire Guidewire ANG	72
M00556861	Extendable Jagwire Guidewire ST Stiff	72
M00556871	Extendable Jagwire Guidewire ANG Stiff	72
M00556902	Jagtail™ Guidewire Extension	72
M00557641	Jagwire Guidewire ST	72
Jagwire™ Revolution High Performance Guidewire		
M00557001	Jagwire Revolution High Performance Guidewire	72
M00557011	Jagwire Revolution High Performance Guidewire	72
M00557021	Jagwire Revolution High Performance Guidewire	72
M00557031	Jagwire Revolution High Performance Guidewire	72
M00557041	Jagwire Revolution High Performance Guidewire	72
Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needles – Olympus Scope Compatible		
M00558220	Expect Pulmonary 22ga – Olympus Scope Compatible	89
M00558221	Expect Pulmonary 22ga – Olympus Scope Compatible	89
M00558250	Expect Pulmonary 25ga – Olympus Scope Compatible	89
M00558251	Expect Pulmonary 25ga – Olympus Scope Compatible	89
M00558731	Expect Pulmonary Needle Adaptor – Olympus Scope Compatible	89
CRE™ Fixed Wire Balloon Dilators		
M00558330	CRE Fixed Wire Balloon Dilator	80
M00558340	CRE Fixed Wire Balloon Dilator	80
M00558350	CRE Fixed Wire Balloon Dilator	80
M00558360	CRE Fixed Wire Balloon Dilator	80
M00558361	CRE Fixed Wire Balloon Dilator	80
M00558370	CRE Fixed Wire Balloon Dilator	80
M00558371	CRE Fixed Wire Balloon Dilator	80
M00558380	CRE Fixed Wire Balloon Dilator	80
M00558381	CRE Fixed Wire Balloon Dilator	80
CRE™ Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)		
M00558390	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)	79
M00558400	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)	79
M00558410	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)	79
M00558420	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)	79
M00558430	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)	79
M00558440	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)	79



Product Code Index

Order Number	Description	Page Number
CRE™ Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)		
M00558450	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)	78
M00558460	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)	78
M00558470	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)	78
M00558480	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)	78
M00558490	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)	78
M00558500	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)	78
CRE™ PRO Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Biliary) (including DASE indication)		
M00558600	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	79
M00558610	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	79
M00558620	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	79
M00558630	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	79
M00558640	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	79
M00558650	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	79
CRE™ PRO Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic/Biliary) (including DASE indication)		
M00558660	CRE PRO Wireguided Esophageal/Pyloric/Colonic/Biliary Balloon Dilators (including DASE indication)	78
M00558670	CRE PRO Wireguided Esophageal/Pyloric/Colonic/Biliary Balloon Dilators (including DASE indication)	78
M00558680	CRE PRO Wireguided Esophageal/Pyloric/Colonic/Biliary Balloon Dilators (including DASE indication)	78
M00558690	CRE PRO Wireguided Esophageal/Pyloric/Colonic/Biliary Balloon Dilators (including DASE indication)	78
M00558700	CRE PRO Wireguided Esophageal/Pyloric/Colonic/Biliary Balloon Dilators (including DASE indication)	78
M00558710	CRE PRO Wireguided Esophageal/Pyloric/Colonic/Biliary Balloon Dilators (including DASE indication)	78
CRE™ RX Biliary Balloon Dilatation Catheter		
M00558920	CRE RX Biliary Balloon Dilatation Catheter	50
M00558930	CRE RX Biliary Balloon Dilatation Catheter	50
M00558940	CRE RX Biliary Balloon Dilatation Catheter	50
M00558860	CRE RX Biliary Balloon Dilatation Catheter	50
M00558870	CRE RX Biliary Balloon Dilatation Catheter	50
M00558880	CRE RX Biliary Balloon Dilatation Catheter	50
M00558890	CRE RX Biliary Balloon Dilatation Catheter	50
M00558900	CRE RX Biliary Balloon Dilatation Catheter	50
M00558910	CRE RX Biliary Balloon Dilatation Catheter	50
Gold Probe™ Electrohaemostasis Catheters		
M00560070	Gold Probe Catheter	100
M00560100	Gold Probe Catheter	100
M00560220	Gold Probe Catheter	100
Injection Gold Probe™ Electrohaemostasis/Injection Catheters		
M00560150	Injection Gold Probe Catheter	99
M00560160	Injection Gold Probe Catheter	99
Sensation™ Single-Use Short Throw Snares		
M00560311	Medium Oval	104
M00560321	Small Oval	104
M00560331	Medium Crescent	104
Captivator™ COLD Single-Use Snare		
M00561101	10mm Rounded–Stiff	106
M00561102	10mm Rounded–Stiff	106
M00561103	10mm Rounded–Stiff	106



Product Code Index

Order Number	Description	Page Number
Captivator™ II Single-Use Snares		
M00561191	25 mm Rounded	105
M00561192	25 mm Rounded	105
M00561193	25 mm Rounded	105
M00561221	10 mm Rounded	105
M00561222	10 mm Rounded	105
M00561223	10 mm Rounded	105
M00561231	15 mm Rounded	105
M00561232	15 mm Rounded	105
M00561233	15 mm Rounded	105
M00561241	20 mm Rounded	105
M00561242	20 mm Rounded	105
M00561243	20 mm Rounded	105
M00561291	33 mm Rounded	105
Accessories		
M00561250	Active Cord, Fits Olympus Brand Generator	42, 106
M00561270	Active Cord, Fits Endostat™, Endostat™ II, Endostat™ III, Bovie®, Valleylab®, Aspen Labs and ERBE® Brand Generators	42, 106
Captivator™ EMR		
M00561600	Captivator EMR for Standard Gastrosopes	103
M00561601	Captivator EMR for Standard Gastrosopes	103
M00561610	Captivator EMR for Large Gastrosopes	103
M00561611	Captivator EMR for Large Gastrosopes	103
Rotatable Single-Use Polypectomy Snares		
M00561821	Small Oval	103
M00561831	Medium Oval	103
Captivator™ Single-Use Snares		
M00561311	Large Oval	105
M00562301	Small Oval	105
M00562321	Medium Oval	105
M00562341	Medium Hexagonal	105
M00562371	Medium Crescent	105
M00562391	Large Oval (thin wire)	105
M00562451	Small Hexagonal	105
Captiflex™ Single-Use Snares		
M00562401	Medium Oval	104
M00562402	Medium Oval	104
M00562421	Small Oval	104
M00562422	Small Oval	104
M00562471	Extra Small Oval	104
Sensation™ Single-Use Short Throw Snares		
M00562651	Large Oval	104
M00562652	Large Oval	104
M00562671	Medium Oval	104
M00562672	Medium Oval	104
M00562673	Medium Oval	104
M00562691	Small Oval	104
M00562692	Small Oval	104
M00562693	Small Oval	104



Product Code Index

Order Number	Description	Page Number
WallFlex™ Duodenal Stent System with Anchor Lock Delivery System		
M00565010	WallFlex Duodenal Stent	114
M00565020	WallFlex Duodenal Stent	114
M00565030	WallFlex Duodenal Stent	114
WallFlex™ Colonic Stent System with Anchor Lock Delivery System		
M00565040	WallFlex Colonic Stent	113
M00565050	WallFlex Colonic Stent	113
M00565060	WallFlex Colonic Stent	113
M00565070	WallFlex Colonic Stent	113
M00565080	WallFlex Colonic Stent	113
M00565090	WallFlex Colonic Stent	113
M00565100	WallFlex Colonic Stent	113
M00565110	WallFlex Colonic Stent	113
M00565120	WallFlex Colonic Stent	113
M00565130	WallFlex Colonic Stent	113
M00565140	WallFlex Colonic Stent	113
M00565150	WallFlex Colonic Stent	113
Accessories		
M00566001	Rapid Refill Continuous Injection System	42, 51
Encore™ 26 Inflation Device		
M00566670	Encore 26 Inflation Device	42, 51
WallFlex™ Biliary RX Fully Covered Stent System RMV (Permalume™ Covering, includes benign indication)		
M00570450	WallFlex RX Biliary (Permalume Covering, includes benign indication)	63
M00570460	WallFlex RX Biliary (Permalume Covering, includes benign indication)	63
M00570470	WallFlex RX Biliary (Permalume Covering, includes benign indication)	63
M00570480	WallFlex RX Biliary (Permalume Covering, includes benign indication)	63
M00570490	WallFlex RX Biliary (Permalume Covering, includes benign indication)	63
WallFlex™ Biliary RX Fully Covered Stent System		
M00570500	WallFlex Biliary Stent System (RX Fully Covered)	65
M00570510	WallFlex Biliary Stent System (RX Fully Covered)	65
M00576700	WallFlex Biliary Stent System (RX Fully Covered)	65
M00576710	WallFlex Biliary Stent System (RX Fully Covered)	65
M00570520	WallFlex Biliary Stent System (RX Fully Covered)	65
M00570530	WallFlex Biliary Stent System (RX Fully Covered)	65
M00570540	WallFlex Biliary Stent System (RX Fully Covered)	65
M00576720	WallFlex Biliary Stent System (RX Fully Covered)	65
M00576730	WallFlex Biliary Stent System (RX Fully Covered)	65
WallFlex™ Biliary RX Uncovered Stent System		
M00570600	WallFlex RX Biliary	64
M00570610	WallFlex RX Biliary	64
M00570620	WallFlex RX Biliary	64
M00570630	WallFlex RX Biliary	64
M00570640	WallFlex RX Biliary	64
M00570650	WallFlex RX Biliary	64
M00570660	WallFlex RX Biliary	64
M00570890	WallFlex RX Biliary	64
WallFlex™ Biliary RX Partially Covered Stent System (Permalume™ Covering)		
M00570700	WallFlex RX Biliary (Permalume Covering)	64
M00570710	WallFlex RX Biliary (Permalume Covering)	64



Product Code Index

Order Number	Description	Page Number
	WallFlex™ Biliary RX Partially Covered Stent System (Permalume™ Covering) (CONTINUED)	
M00570720	WallFlex RX Biliary (Permalume Covering)	64
M00570730	WallFlex RX Biliary (Permalume Covering)	64
M00570740	WallFlex RX Biliary (Permalume Covering)	64
	Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	
M00573040	Jagtome RX 44	46
M00573050	Jagtome RX 44	46
M00573060	Jagtome RX 44	46
M00573070	Jagtome RX 44	46
M00573080	Jagtome RX 39	46
M00573090	Jagtome RX 39	46
M00573100	Jagtome RX 39	46
M00573110	Jagtome RX 39	46
	WallFlex™ Biliary RX Fully Covered Stent System RMV (Permalume™ Covering, includes benign indication)	
M00574210	WallFlex RX Biliary (Permalume Covering, includes benign indication)	63
M00574230	WallFlex RX Biliary (Permalume Covering, includes benign indication)	63
M00574240	WallFlex RX Biliary (Permalume Covering, includes benign indication)	63
	Ultraflex™ Tracheobronchial Uncovered Stent System – Distal Release	
M00576300	Ultraflex Tracheobronchial Stent System	91
M00576310	Ultraflex Tracheobronchial Stent System	91
M00576320	Ultraflex Tracheobronchial Stent System	91
M00576330	Ultraflex Tracheobronchial Stent System	91
M00576340	Ultraflex Tracheobronchial Stent System	91
M00576350	Ultraflex Tracheobronchial Stent System	91
M00576360	Ultraflex Tracheobronchial Stent System	91
M00576370	Ultraflex Tracheobronchial Stent System	91
M00576380	Ultraflex Tracheobronchial Stent System	91
	Ultraflex™ Tracheobronchial Uncovered Stent System – Proximal Release	
M00576410	Ultraflex Tracheobronchial Stent System	90
M00576420	Ultraflex Tracheobronchial Stent System	90
M00576430	Ultraflex Tracheobronchial Stent System	90
M00576440	Ultraflex Tracheobronchial Stent System	90
M00576450	Ultraflex Tracheobronchial Stent System	90
M00576460	Ultraflex Tracheobronchial Stent System	90
M00576470	Ultraflex Tracheobronchial Stent System	90
M00576480	Ultraflex Tracheobronchial Stent System	90
M00576490	Ultraflex Tracheobronchial Stent System	90
	WallFlex™ Biliary Stent System (RX Fully Covered)	
M00576700	WallFlex Biliary Stent System (RX Fully Covered)	65
M00576710	WallFlex Biliary Stent System (RX Fully Covered)	65
M00576720	WallFlex Biliary Stent System (RX Fully Covered)	65
M00576730	WallFlex Biliary Stent System (RX Fully Covered)	65
	WallFlex™ Biliary RX Partially Covered Stent System (Permalume™ Covering)	
M00576740	WallFlex RX Biliary (Permalume Covering)	64
M00576750	WallFlex RX Biliary (Permalume Covering)	64
M00576760	WallFlex RX Biliary (Permalume Covering)	64
M00576770	WallFlex RX Biliary (Permalume Covering)	64



Product Code Index

Order Number	Description	Page Number
	WallFlex™ Biliary RX Uncovered Stent System	
M00576780	WallFlex RX Biliary	64
M00576790	WallFlex RX Biliary	64
	WallFlex™ Biliary RX Fully Covered System RMV (Permalume™ Covering, includes benign indication)	
M00576800	WallFlex RX Biliary (Permalume Covering, includes benign indication)	63
M00576810	WallFlex RX Biliary (Permalume Covering, includes benign indication)	63
M00576820	WallFlex RX Biliary (Permalume Covering, includes benign indication)	63
M00576830	WallFlex RX Biliary (Permalume Covering, includes benign indication)	63
	Ultraflex™ Tracheobronchial Covered Stent System – Distal Release	
M00576520	Ultraflex Tracheobronchial Stent System	90
M00576530	Ultraflex Tracheobronchial Stent System	90
M00576540	Ultraflex Tracheobronchial Stent System	90
M00576550	Ultraflex Tracheobronchial Stent System	90
M00576560	Ultraflex Tracheobronchial Stent System	90
M00576570	Ultraflex Tracheobronchial Stent System	90
M00576580	Ultraflex Tracheobronchial Stent System	90
M00576590	Ultraflex Tracheobronchial Stent System	90
M00576600	Ultraflex Tracheal NG Stent System	90
M00576610	Ultraflex Tracheal NG Stent System	90
M00576620	Ultraflex Tracheal NG Stent System	90
M00576630	Ultraflex Tracheal NG Stent System	90
M00576640	Ultraflex Tracheal NG Stent System	90
M00576650	Ultraflex Tracheal NG Stent System	90
M00576660	Ultraflex Tracheal NG Stent System	90
M00576670	Ultraflex Tracheal NG Stent System	90
M00576680	Ultraflex Tracheal NG Stent System	90
M00576690	Ultraflex Tracheal NG Stent System	90
	Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	
M00583000	Hydratome RX 49	45
M00583010	Hydratome RX 49	45
M00583020	Hydratome RX 49	45
M00583030	Hydratome RX 49	45
M00583040	Hydratome RX 44	45
M00583050	Hydratome RX 44	45
M00583060	Hydratome RX 44	45
M00583070	Hydratome RX 44	45
	TRUEtome™ Jag Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	
M00583080	TRUEtome Jag 39 Sphincterotome	39
M00583090	TRUEtome Jag 44 Sphincterotome	39
M00583100	TRUEtome Jag 44 Sphincterotome	39
M00583110	TRUEtome Jag 44 Sphincterotome	39
M00583120	TRUEtome Jag 39 Sphincterotome	39
M00583130	TRUEtome Jag 44 Sphincterotome	39
M00583140	TRUEtome Jag 44 Sphincterotome	39
M00583150	TRUEtome Jag 44 Sphincterotome	39
	TRUEtome™ Hydra Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	
M00583160	TRUEtome Hydra 44 Sphincterotome	39
M00583170	TRUEtome Hydra 44 Sphincterotome	39



Product Code Index

Order Number	Description	Page Number
TRUEtome™ Dreamwire™ Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire		
M00583190	TRUEtome Dreamwire 44 Sphincterotome	39
M00583210	TRUEtome Dreamwire 44 Sphincterotome	39
TRUEtome™ Jagwire Revolution Cannulating Sphincterotomes with Pre-loaded Jagwire™ Revolution Guidewire		
M00583260	TRUEtome 39 Sphincterotome	39
M00583270	TRUEtome 44 Sphincterotome	39
M00583280	TRUEtome 39 Sphincterotome	39
M00583290	TRUEtome 44 Sphincterotome	39
Dreamtome™ RX Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire		
M00584000	Dreamtome RX 49 Sphincterotome	45
M00584010	Dreamtome RX 49 Sphincterotome	45
M00584020	Dreamtome RX 49 Sphincterotome	45
M00584030	Dreamtome RX 49 Sphincterotome	45
M00584040	Dreamtome RX 44 Sphincterotome	45
M00584050	Dreamtome RX 44 Sphincterotome	45
M00584060	Dreamtome RX 44 Sphincterotome	45
M00584070	Dreamtome RX 44 Sphincterotome	45
TRUEtome™ Cannulating Sphincterotomes		
M00584120	TRUEtome 39 Sphincterotome	40
M00584121	TRUEtome 39 Sphincterotome	40
M00584122	TRUEtome 39 Sphincterotome	40
M00584150	TRUEtome 49 Sphincterotome	40
M00584151	TRUEtome 49 Sphincterotome	40
M00584152	TRUEtome 49 Sphincterotome	40
M00584160	TRUEtome 49 Sphincterotome	40
M00584161	TRUEtome 49 Sphincterotome	40
M00584162	TRUEtome 49 Sphincterotome	40
M00584170	TRUEtome 44 Sphincterotome	40
M00584171	TRUEtome 44 Sphincterotome	40
M00584172	TRUEtome 44 Sphincterotome	40
M00584180	TRUEtome 44 Sphincterotome	40
M00584181	TRUEtome 44 Sphincterotome	40
M00584182	TRUEtome 44 Sphincterotome	40
M00584190	TRUEtome 39 Sphincterotome	40
M00584191	TRUEtome 39 Sphincterotome	40
M00584192	TRUEtome 39 Sphincterotome	40
Jagtome™ Revolution RX Cannulating Sphincterotome (Autotome™ Preloaded with Jagwire Revolution High Performance Guidewire)		
M00584240	Jagtome Revolution RX 39 Sphincterotome	46
M00584250	Jagtome Revolution RX 39 Sphincterotome	46
M00584260	Jagtome Revolution RX 39 Sphincterotome	46
M00584270	Jagtome Revolution RX 39 Sphincterotome	46
Autotome™ Pro RX Cannulating Sphincterotomes		
M00584380	Autotome Pro RX 44	47
M00584390	Autotome Pro RX 39	47
Autotome™ Pro RX Preloaded Sphincterotomes		
M00584400	Dreamtome Pro RX 44	47
M00584410	Dreamtome Pro RX 44	47
M00584420	Hydratome Pro RX 44	47



Product Code Index

Order Number	Description	Page Number
	Autotome™ Pro RX Preloaded Sphincterotomes (CONTINUED)	
M00584430	Hydratome Pro RX 44	47
M00584440	Jagtome Pro RX 44	47
M00584450	Jagtome Pro RX 44	47
M00584460	Jagtome Pro RX 39	47
M00584470	Jagtome Revolution Pro RX 39	47
M00584480	Jagtome Revolution Pro RX 39	47
	Seal™	
SBC-501-100	Seal Single Use Biopsy Valve	116
	HedgeHog™ Cleaning Brushes	
SBD-228-50	HedgeHog Double-End Channel Brush	116
SBD-289-50	HedgeHog Dual-End Channel/Valve Brush	116
SBD-291-50	HedgeHog Sweeper Brush	116
SBS-227-50	HedgeHog Single-End Channel Brush	116
SBB-382-50	HedgeHog Double Bundle	116
	Blox™ Bite Blocks	
SBT-114-100	Blox Bite Block 54Fr with Dental Rim	117
SBT-321-100	Blox Bite Block 54Fr without Dental Rim	117
SBT-546-100	Blox Bite Block 60Fr with Dental Rim	117
	Scivita Medical™ Disposable Equipment	
SBV-1A-B	Single-use Broncho Videoscope	85, 88
SBV-1B-B	Single-use Broncho Videoscope	85, 88
SBV-1C-B	Single-use Broncho Videoscope	85, 88
	Hydra™ Irrigation System	
SCT-609-10	Hydra Water Bottle Cap with CO2	116
SIT-361-100	Hydra Disposable Auxiliary Water Jet Connectors	116
SIT-575-15	Hydra Irrigation Tubing	116
SIT-577-15	Hydra Irrigation Tubing	116
	CinchPad™ Procedure & Endoscope Transport Pad	
SPX-017-25	CinchPad Endoscope Transport Pad	117
	Scivita Medical™ Capital Equipment	
HDVS-S300D	Full HD Visualization Endoscopic Image Processor / Tablet	84
HDVS-S200D	Endoscopic Image Processing Device / Controller	84
	Scivita Medical™ Disposable Equipment	
SCV-P-01G	Single-use Video Cholangioscope 3.70mm OD / 2.0mm WC	85
SCV-P-02G	Single-use Video Cholangioscope 2.95mm OD / 1.2mm WC	85
	Scivita Medical™ Service Codes	
SRVM005HDVS200EVC3	HDVS-S200D 1YR EverCare POST	85, 88
SRVM005HDVS200EVC4	HDVS-S200D 1YR EverCare POS	85, 88
SRVM005HDVS200EVC5	HDVS-S200D 2YR EverCare POST	85, 88
SRVM005HDVS200EVC6	HDVS-S200D 2YR EverCare POS	85, 88
SRVM005HDVS200EVC1	HDVS-S200D 3YR EverCare POST	85, 88
SRVM005HDVS200EVC2	HDVS-S200D 3YR EverCare POS	85, 88
SRVM005HDVS200EVC7	HDVS-S200D 4YR EverCare POST	85, 88
SRVM005HDVS200EVC8	HDVS-S200D 4YR EverCare POS	85, 88
SRVM005HDVS300EVC3	HDVS-S300D 1YR EverCare POST	85, 88
SRVM005HDVS300EVC4	HDVS-S300D 1YR EverCare POS	85, 88
SRVM005HDVS300EVC5	HDVS-S300D 2YR EverCare POST	85, 88



Product Code Index



Order Number	Description	Page Number
Scivita Medical™ Service Codes (CONTINUED)		
SRVM005HDVS300EVC6	HDVS-S300D 2YR EverCare POS	85, 88
SRVM005HDVS300EVC1	HDVS-S300D 3YR EverCare POST	85, 88
SRVM005HDVS300EVC2	HDVS-S300D 3YR EverCare POS	85, 88
SRVM005HDVS300EVC7	HDVS-S300D 4YR EverCare POST	85, 88
SRVM005HDVS300EVC8	HDVS-S300D 4YR EverCare POS	85, 88
SRVM005HDVSS200EC1	EssentialCare S200 -1 year POS	85, 88
SRVM005HDVSS200EC2	EssentialCare S200 -1 year POST	85, 88
SRVM005HDVSS200EC3	EssentialCare S200 -2 years POS	85, 88
SRVM005HDVSS200EC4	EssentialCare S200 -2 years POST	85, 88
SRVM005HDVSS200EC5	EssentialCare S200 -3 years POS	85, 88
SRVM005HDVSS200EC6	EssentialCare S200 -3 years POST	85, 88
SRVM005HDVSS200EC7	EssentialCare S200 -4 years POS	85, 88
SRVM005HDVSS200EC8	EssentialCare S200 -4 years POST	85, 88
SRVM005HDVSS300EC1	EssentialCare for S300 -1 year POS	85, 88
SRVM005HDVSS300EC2	EssentialCare for S300 -1 year POST	85, 88
SRVM005HDVSS300EC3	EssentialCare for S300 -2 years POS	85, 88
SRVM005HDVSS300EC4	EssentialCare for S300 -2 years POST	85, 88
SRVM005HDVSS300EC5	EssentialCare for S300 -3 years POS	85, 88
SRVM005HDVSS300EC6	EssentialCare for S300 -3 years POST	85, 88
SRVM005HDVSS300EC7	EssentialCare for S300 -4 years POS	85, 88
SRVM005HDVSS300EC8	EssentialCare for S300 -4 years POST	85, 88
SRVM005HDVSS200STK	HDVS-S200 Safety Testing (STK)	85, 88
SRVM005HDVSS200EST	HDVS-S200 Electrical Safety Test (EST)	85, 88
SRVM005HDVSS300STK	HDVS-S300D Safety Testing (STK)	85, 88
SRVM005HDVSS300EST	HDVS-S300D Electrical Safety Test (EST)	85, 88
SRVM005HDVFSFEFFEE	Scivita FSE Flat Fee	85, 88
Autolith™ Service Codes		
SRVM005ATLT1001	Autolith EssentialCare POS 1 Year	86
SRVM005ATLT1005	Autolith EssentialCare POST 1 Year	86
SRVM005ATLT1002	Autolith EssentialCare POS 2 Years	86
SRVM005ATLT1006	Autolith EssentialCare POST 2 Years	86
SRVM005ATLT1003	Autolith EssentialCare POS 3 Years	86
SRVM005ATLT1007	Autolith EssentialCare POST 3 Years	86
SRVM005ATLT1004	Autolith EssentialCare POS 4 Years	86
SRVM005ATLT1008	Autolith EssentialCare POST 4 Years	86
SRVM005ATLTECP1001	Autolith EssentialCare PLUS POS 3 Years	86
SRVM005ATLTECP1002	Autolith EssentialCare PLUS POS 4 Years	86
SRVM005ATLTECP1003	Autolith EssentialCare PLUS POS 5 Years	86
SRVM005ATLTECP1004	Autolith EssentialCare PLUS POST 3 Years	86
SRVM005ATLTECP1005	Autolith EssentialCare PLUS POST 4 Years	86
SRVM005ATLTECP1006	Autolith EssentialCare PLUS POST 5 Years	86
EXALT™ Model D Service Codes		
SRVM005EXALTDTC01	EXALT Model D TotalCare POS 1 Year	83
SRVM005EXALTDTC02	EXALT Model D TotalCare POS 2 Years	83
SRVM005EXALTDTC03	EXALT Model D TotalCare POS 3 Years	83
SRVM005EXALTDTC04	EXALT Model D TotalCare POS 4 Years	83
SRVM005EXALTDTC05	EXALT Model D TotalCare POS 5 Years	83
SRVM005EXALTDTC06	EXALT Model D TotalCare POST 1 Year	83



Product Code Index


Order Number	Description	Page Number
EXALT™ Model D Service Codes (CONTINUED)		
SRVM005EXALDTC07	EXALT Model D TotalCare POST 2 Years	83
SRVM005EXALDTC08	EXALT Model D TotalCare POST 3 Years	83
SRVM005EXALDTC09	EXALT Model D TotalCare POST 4 Years	83
SRVM005EXALDTC10	EXALT Model D TotalCare POST 5 Years	83
SpyGlass™ DS Service Codes		
SRVM005SPYDIG1001	SpyGlass Digital EssentialCare POS 1 Year	84
SRVM005SPYDIG1002	SpyGlass Digital EssentialCare POS 2 Year	84
SRVM005SPYDIG1003	SpyGlass Digital EssentialCare POS 3 Year	84
SRVM005SPYDIG1004	SpyGlass Digital EssentialCare POS 4 Year	84
SRVM005SPYDIG1005	SpyGlass Digital EssentialCare POST 1 Year	84
SRVM005SPYDIG1006	SpyGlass Digital EssentialCare POST 2 Year	84
SRVM005SPYDIG1007	SpyGlass Digital EssentialCare POST 3 Year	84
SRVM005SPYDIG1008	SpyGlass Digital EssentialCare POST 4 Year	84
SRVM005SPYDIG1009	SpyGlass Digital EverCare POS 1 Year	84
SRVM005SPYDIG1010	SpyGlass Digital EverCare POS 2 Year	84
SRVM005SPYDIG1011	SpyGlass Digital EverCare POS 3 Year	84
SRVM005SPYDIG1012	SpyGlass Digital EverCare POS 4 Year	84
SRVM005SPYDIG1013	SpyGlass Digital EverCare POST 1 Year	84
SRVM005SPYDIG1014	SpyGlass Digital EverCare POST 2 Year	84
SRVM005SPYDIG1015	SpyGlass Digital EverCare POST 3 Year	84
SRVM005SPYDIG1016	SpyGlass Digital EverCare POST 4 Year	84
SRVM005SPYDSEVCS01	SpyGlass DS EverCare Select 3 Years POS	84
SRVM005SPYDSEVCS02	SpyGlass DS EverCare Select 3 Years POST	84
SRVM005SPYDSEVCS03	SpyGlass DS EverCare Select 4 Years POS	84
SRVM005SPYDSEVCS04	SpyGlass DS EverCare Select 4 Years POST	84
SpyGlass™ Discover Service Codes		
SRVM005SPYDISEC01	SpyGlass Digital EssentialCare POS 1 Year	86
SRVM005SPYDISEC02	SpyGlass Digital EssentialCare POS 2 Year	86
SRVM005SPYDISEC03	SpyGlass Digital EssentialCare POS 3 Year	86
SRVM005SPYDISEC04	SpyGlass Digital EssentialCare POS 4 Year	86
SRVM005SPYDISEC05	SpyGlass Digital EssentialCare POST 1 Year	86
SRVM005SPYDISEC06	SpyGlass Digital EssentialCare POST 2 Year	86
SRVM005SPYDISEC07	SpyGlass Digital EssentialCare POST 3 Year	86
SRVM005SPYDISEC08	SpyGlass Digital EssentialCare POST 4 Year	86
TrapEase™ Polyp Traps		
STE-297-50	TrapEase Single Trap	107
STE-298-25	TrapEase Quad Trap	107
Orca™		
SUV-617-50	Orca Sterile, Single Use Air/Water and Suction Valves	116
OrcaPod™		
SUV-628-50	OrcaPod Sterile, Single Use Air/Water and Suction Valves with Seal Biopsy Valve	116
SUV-629-50	OrcaPod Sterile, Single Use Air/Water and Suction Valves, Seal Biopsy Valve and Hydra Water Jet Connector	116





Biliary	RX Biliary™ System™	Biliary Lumen Patency	Endoscopic Ultrasound (EUS)	Guidewires	Biopsy	Dilation	Imaging Systems	Pulmonary	Closure	Haemostasis	Tissue Resection	Esophageal and Enteral Metal Stents	Infection Prevention	Endobariatrics
---------	---------------------	-----------------------	-----------------------------	------------	--------	----------	-----------------	-----------	---------	-------------	------------------	-------------------------------------	----------------------	----------------

CANNULAS




Tandem™ XL Triple-Lumen ERCP Cannula



Order Number	Description	Shaft OD (F) / (mm)	Distal OD (F) / (mm)	Length (cm)	Tip	Recommended Guidewire (in) / (mm)
M00535700	Tandem XL Triple-Lumen ERCP Cannula	7.0-5.5 / 2.33-1.83	5.5 / 1.83	200	Tapered	.035 / 0.89

Packaged one per box




Contour™ ERCP Cannulas



Order Number	Description	Shaft OD (F) / (mm)	Length (cm)	Tip	Recommended Guidewire (in) / (mm)
M00530850	Contour ERCP Cannula	5 / 1.67	210	Standard	.035 / 0.89
M00530860	Contour ERCP Cannula	5 / 1.67	210	Tapered	.035 / 0.89
M00530870	Contour ERCP Cannula	5 / 1.67	210	Ultra Tapered	.025 / 0.64
M00530880	Contour ERCP Cannula	5 / 1.67	210	5-4-3	.018 / 0.46
M00530890	Contour ERCP Cannula	5 / 1.67	210	Ball Tip	.035 / 0.89

Packaged one per box

SPHINCTEROTOMES



TRUEtome™ Jagwire Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire



Order Number	Description	Cut Wire (mm)	Tip OD (F) / (mm)	Tip Length (mm)	Pre-loaded Guidewire (in / cm)	Units
M00583080	TRUEtome 39 Sphincterotome	20	3.9 / 1.30	5	.025 / 450 Jagwire Guidewire	box 1
M00583090	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.025 / 450 Jagwire Guidewire	box 1
M00583100	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.035 / 450 Jagwire Guidewire	box 1
M00583110	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.035 / Stiff 450 Jagwire Guidewire	box 1
M00583120	TRUEtome 39 Sphincterotome	30	3.9 / 1.30	5	.025 / 450 Jagwire Guidewire	box 1
M00583130	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.025 / 450 Jagwire Guidewire	box 1
M00583140	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.035 / 450 Jagwire Guidewire	box 1
M00583150	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.035 / Stiff 450 Jagwire Guidewire	box 1

TRUEtome™ Dreamwire™ Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire



Order Number	Description	Cut Wire (mm)	Tip OD (F) / (mm)	Tip Length (mm)	Pre-loaded Guidewire (in / cm)	Units
M00583190	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.035 / 450 Dreamwire Guidewire	box 1
M00583210	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.035 / 450 Dreamwire Guidewire	box 1

TRUEtome™ Hydra Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire




Order Number	Description	Cut Wire (mm)	Tip OD (F) / (mm)	Tip Length (mm)	Pre-loaded Guidewire (in / cm)	Units
M00583160	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.035 / 450 Hydra Jagwire Guidewire	box 1
M00583170	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.035 / 450 Hydra Jagwire Guidewire	box 1

TRUEtome™ Jagwire Revolution Cannulating Sphincterotomes with Pre-loaded Jagwire™ Revolution Guidewire



Order Number	Description	Cut Wire (mm)	Tip OD (F) / (mm)	Tip Length (mm)	Pre-loaded Guidewire (in / cm)	Units
M00583260	TRUEtome 39 Sphincterotome	20	3.9 / 1.30	5	.025 / 450 Jagwire Revolution Guidewire	box 1
M00583270	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.025 / 450 Jagwire Revolution Guidewire	box 1
M00583280	TRUEtome 39 Sphincterotome	30	3.9 / 1.30	5	.025 / 450 Jagwire Revolution Guidewire	box 1
M00583290	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.025 / 450 Jagwire Revolution Guidewire	box 1


SPHINCTEROTOMES (CONTINUED)



TRUEtome™ Cannulating Sphincterotomes

Order Number	Description	Cut Wire (mm)	Distal OD (F) / (mm)	Tip Length (mm)	Catheter Length (cm)	Units
M00584190	TRUEtome 39	20	3.9 / 1.30	5	200	box 1
M00584120	TRUEtome 39	30	3.9 / 1.30	5	200	box 1
M00584191	TRUEtome 39	20	3.9 / 1.30	5	200	box 10
M00584121	TRUEtome 39	30	3.9 / 1.30	5	200	box 10
M00584192	TRUEtome 39	20	3.9 / 1.30	5	200	box 20
M00584122	TRUEtome 39	30	3.9 / 1.30	5	200	box 20
M00584170	TRUEtome 44	20	4.4 / 1.47	5	200	box 1
M00584180	TRUEtome 44	30	4.4 / 1.47	5	200	box 1
M00584171	TRUEtome 44	20	4.4 / 1.47	5	200	box 10
M00584181	TRUEtome 44	30	4.4 / 1.47	5	200	box 10
M00584172	TRUEtome 44	20	4.4 / 1.47	5	200	box 20
M00584182	TRUEtome 44	30	4.4 / 1.47	5	200	box 20
M00584150	TRUEtome 49	20	4.9 / 1.63	5	200	box 1
M00584160	TRUEtome 49	30	4.9 / 1.63	5	200	box 1
M00584151	TRUEtome 49	20	4.9 / 1.63	5	200	box 10
M00584161	TRUEtome 49	30	4.9 / 1.63	5	200	box 10
M00584152	TRUEtome 49	20	4.9 / 1.63	5	200	box 20
M00584162	TRUEtome 49	30	4.9 / 1.63	5	200	box 20

Recommended Guidewire: .025in (0.64 mm) for TRUEtome 39; .035in (0.89 mm) for TRUEtome 44 and 49




Microknife™ XL Triple-Lumen Needle-Knife

Order Number	Description	Shaft OD (F) / (mm)	Catheter Length (cm)
M00532810	Microknife XL Triple-Lumen Needle-Knife	7-5.5 / 2.33-1.83	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

SPHINCTEROTOMES (CONTINUED)




Ultratome XL™ Triple-Lumen Sphincterotomes

Order Number	Description	Cut Wire (mm)	Distal OD (F) / (mm)	Tip Length (mm)	Catheter Length (cm)
M00535900	Ultratome XL Short Nose	20	5.5 / 1.83	5	200
M00535910	Ultratome XL Long Nose	20	5.5 / 1.83	20	200
M00535920	Ultratome XL Short Nose	30	5.5 / 1.83	5	200
M00535930	Ultratome XL Long Nose	30	5.5 / 1.83	20	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

STONE RETRIEVAL




Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Distal Injection

Order Number	Description	Balloon OD (mm)	Injection	Catheter OD (F) / (mm)	System Length (cm)
M00547100	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Above	7-6 / 2.33-2.00	200
M00547110	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Above	7-6 / 2.33-2.00	200
M00547120	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Above	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Proximal Injection

Order Number	Description	Balloon OD (mm)	Injection	Catheter OD (F) / (mm)	System Length (cm)
M00547130	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Below	7-6 / 2.33-2.00	200
M00547140	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Below	7-6 / 2.33-2.00	200
M00547150	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Below	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire


ACCESSORIES

Biliary Accessories

Order Number	Description	Units
M00561250	Active Cord, Fits Olympus Brand Generator	box 1
M00561270	Active Cord, Fits Endostat™, Endostat™ II, Endostat™ III, Bovie®, Valleylab®, Aspen Labs and ERBE® Brand Generators	box 1
M00566670	⊗ Encore 26 Inflation Device	box 1
M00566001	⊗ Rapid Refill Continuous Injection System	box 10


Biliary Plastic Stents are located in the Stents Section on pages 53–57

Biliary Guidewires are located in the Guidewires Section on pages 71 – 73



Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics


CANNULAS



Rapid Exchange XL Cannula

Order Number	Description	Shaft OD (F) / (mm)	Length (cm)	Recommended Guidewire (in) / (mm)
M00545820	Rapid Exchange XL Cannula	7-5.5 / 2.33-1.84	195	.035 / 0.89

Packaged one per box




RX ERCP Cannulas

Order Number	Description	Shaft OD (F) / (mm)	Length (cm)	Recommended Guidewire (in) / (mm)
M00545760	RX ERCP Cannula Standard Tip	8.5-5 / 2.83-1.67	210	.035 / 0.89
M00545770	RX ERCP Cannula Tapered Tip	8.5-5 / 2.83-1.67	210	.035 / 0.89
M00545780	RX ERCP Cannula Ball Tip	8.5-5 / 2.83-1.67	210	.035 / 0.89

Packaged one per box

SPHINCTEROTOMES



Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire




Order Number	Description	Cut Wire (mm)	Shaft OD (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	Catheter Length (cm)	PRE-LOADED GUIDEWIRE** Length (cm)	OD (in) / (mm)
M00583000	Hydratome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	260	.035 / 0.89
M00583010	Hydratome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	260	.035 / 0.89
M00583020	Hydratome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	450	.035 / 0.89
M00583030	Hydratome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	450	.035 / 0.89
M00583040*	Hydratome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	260	.035 / 0.89
M00583050*	Hydratome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	260	.035 / 0.89
M00583060*	Hydratome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	450	.035 / 0.89
M00583070*	Hydratome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	450	.035 / 0.89

Packaged one per box

* Includes handle brake

** Pre-loaded with .035in (0.89mm) Hydra Jagwire™ Guidewire



Dreamtome™ RX Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire




Order Number	Description	Cut Wire (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	PRE-LOADED GUIDEWIRE** Length (cm)	OD (in) / (mm)
M00584000	Dreamtome RX 49	20	5	4.9 / 1.63	260	.035 / 0.89
M00584010	Dreamtome RX 49	30	5	4.9 / 1.63	260	.035 / 0.89
M00584020	Dreamtome RX 49	20	5	4.9 / 1.63	450	.035 / 0.89
M00584030	Dreamtome RX 49	30	5	4.9 / 1.63	450	.035 / 0.89
M00584040*	Dreamtome RX 44	20	5	4.4 / 1.47	260	.035 / 0.89
M00584050*	Dreamtome RX 44	30	5	4.4 / 1.47	260	.035 / 0.89
M00584060*	Dreamtome RX 44	20	5	4.4 / 1.47	450	.035 / 0.89
M00584070*	Dreamtome RX 44	30	5	4.4 / 1.47	450	.035 / 0.89

Packaged one per box

* Includes handle brake

** Pre-loaded with Straight-tip guidewires only

SPHINCTEROTOMES (CONTINUED)




Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire

Order Number	Description	Cut Wire (mm)	Shaft OD (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	Catheter Length (cm)	PRE-LOADED GUIDEWIRE**		
							Shaft Stiffness	Length (cm)	OD (in) / (mm)
M00573040*	Jagtome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	260	.035 / 0.89
M00573060*	Jagtome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	260	.035 / 0.89
M00573050*	Jagtome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	450	.035 / 0.89
M00573070*	Jagtome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	450	.035 / 0.89
M00573080*	Jagtome RX 39	20	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Stiff	260	.025 / 0.64
M00573100*	Jagtome RX 39	30	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Stiff	260	.025 / 0.64
M00573090*	Jagtome RX 39	20	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Standard	450	.025 / 0.64
M00573110*	Jagtome RX 39	30	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Standard	450	.025 / 0.64

Packaged one per box

* Includes handle brake


** Pre-loaded with the high performance Jagwire™ Guidewire (Jagtome RX 49 and RX 44 contain the Stiff Shaft Jagwire and Jagtome RX 39 the Standard Jagwire)



Jagtome™ Revolution RX Cannulating Sphincterotome (Autotome™ Preloaded with Jagwire Revolution High Performance Guidewire)

Order Number	Description	Cut Wire (mm)	Tip Length (mm)	Distal Tip OD (F)	Preloaded Guidewire (in) / (cm)	Guidewire	Tip Style
M00584240	Jagtome Revolution RX 39 Sphincterotome	20	5	3.9	.025 / 260	Jagwire Revolution Guidewire	Straight
M00584250	Jagtome Revolution RX 39 Sphincterotome	20	5	3.9	.025 / 450	Jagwire Revolution Guidewire	Straight
M00584260	Jagtome Revolution RX 39 Sphincterotome	30	5	3.9	.025 / 260	Jagwire Revolution Guidewire	Straight
M00584270	Jagtome Revolution RX 39 Sphincterotome	30	5	3.9	.025 / 450	Jagwire Revolution Guidewire	Straight

SPHINCTEROTOMES (CONTINUED)





Autotome™ Pro RX Cannulating Sphincterotomes

Order Number	Description	Cut Wire (mm)	Tip Length (mm)	Distal Tip OD (Fr)	Recommended Guidewire (in / mm)	Units
M00584380	Autotome Pro RX 44	25	5	4.4	.025 / .035	box 1
M00584390	Autotome Pro RX 39	25	5	3.9	.025	box 1

Autotome™ Pro RX Preloaded Sphincterotomes

Order Number	Description	Distal Tip OD (Fr)	Cut Wire (mm)	Tip Length (mm)	Guidewire Diameter (in)	Guidewire Length (cm)	Guidewire	Units
M00584400	Dreamtome™ Pro RX 44 - 25mm/260cm	4.4	25	5	.025/.035	260	Dreamwire	box 1
M00584410	Dreamtome Pro RX 44 - 25mm/450cm	4.4	25	5	.025/.035	450	Dreamwire	box 1
M00584420	Hydratome™ Pro RX 44 - 25mm/260cm	4.4	25	5	.025/.035	260	Hydra Jagwire	box 1
M00584430	Hydratome Pro RX 44 - 25mm/450cm	4.4	25	5	.025/.035	450	Hydra Jagwire	box 1
M00584440	Jagtome™ Pro RX 44 - 25-260-035	4.4	25	5	.025/.035	260	Jagwire	box 1
M00584450	Jagtome Pro RX 44 - 25-450-035	4.4	25	5	.025/.035	450	Jagwire	box 1
M00584460	Jagtome Pro RX 39 - 25-260-025	3.9	25	5	.025	260	Jagwire	box 1
M00584470	Jagtome™ Revolution Pro RX 39 260cm	3.9	25	5	.025	260	Jagwire Revolution	box 1
M00584480	Jagtome Revolution Pro RX 39 450cm	3.9	25	5	.025	450	Jagwire Revolution	box 1



Autotome™ RX Cannulating Sphincterotomes

Order Number	Description	Cut Wire (mm)	Shaft OD (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	Catheter Length (cm)	Recommended Guidewire (in) / (mm)
M00545150	Autotome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	.035 / 0.89
M00545160	Autotome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	.035 / 0.89
M00545170*	Autotome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	.035 / 0.89
M00545180*	Autotome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	.035 / 0.89
M00545190*	Autotome RX 39	20	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	.025 / 0.64
M00545200*	Autotome RX 39	30	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	.025 / 0.64

Packaged one per box

Recommended Guidewires: For Autotome RX 44 and 49, recommended guidewire is .035in (0.89mm) Dreamwire™ High Performance Guidewire

For Autotome RX 39, recommended guidewire is .025in (0.69mm) Jagwire™ Guidewire

* Includes handle brake

SPHINCTEROTOMES (CONTINUED)

RX Needle-Knife XL Triple-Lumen Needle-Knife




Order Number	Description	Cut Wire (mm)	Shaft OD (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)
M00545840	RX Needle-Knife	5	7-5.5 / 2.3-1.83	5	5.5 / 1.8

Packaged one per box

Recommended Guidewire: .035 in (0.89mm) Dreamwire™ High Performance Guidewire

STONE RETRIEVAL




Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Distal Injection



Order Number	Description	Balloon OD (mm)	Injection	Catheter OD (F) / (mm)	System Length (cm)
M00547000	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Above	7-6 / 2.33-2.00	200
M00547010	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Above	7-6 / 2.33-2.00	200
M00547020	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Above	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Proximal Injection




Order Number	Description	Balloon OD (mm)	Injection	Catheter OD (F) / (mm)	System Length (cm)
M00547030	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Below	7-6 / 2.33-2.00	200
M00547040	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Below	7-6 / 2.33-2.00	200
M00547050	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Below	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

STONE RETRIEVAL (CONTINUED)




Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Distal Injection

Order Number	Description	Balloon OD (mm)	Injection	Catheter OD (F) / (mm)	System Length (cm)
M00547300	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Above	7-6 / 2.33-2.00	200
M00547310	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Above	7-6 / 2.33-2.00	200
M00547320	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Above	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire




Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Proximal Injection

Order Number	Description	Balloon OD (mm)	Injection	Catheter OD (F) / (mm)	System Length (cm)
M00547330	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Below	7-6 / 2.33-2.00	200
M00547340	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Below	7-6 / 2.33-2.00	200
M00547350	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Below	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Trapezoid™ RX Wireguided Retrieval Baskets

Order Number	Description	Open Basket Diameter (cm)	Required Working Channel (mm)
M00510860	Trapezoid RX Wireguided Retrieval Basket	1.5	3.2
M00510870	Trapezoid RX Wireguided Retrieval Basket	2.0	3.2
M00510880	Trapezoid RX Wireguided Retrieval Basket	2.5	3.2
M00510890*	Trapezoid RX Wireguided Retrieval Basket	3.0	3.2
M00550620	Alliance™ II Inflation Handle**		


Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

* The Trapezoid 3 x 6cm basket is designed for crushing calculus larger than 1.5cm in diameter

** Alliance™ II recommended for use as Mechanical Lithotripter. Alliance II is reusable

DILATATION



Hurricane™ RX Biliary Balloon Dilatation Catheters



Order Number	Description	Balloon OD (mm)	Balloon Length (cm)	Catheter OD (F) / (mm)	Maximum Inflation Pressure (ATM) / (kPa)	Minimum Channel Size (mm)	Recommended Guidewire (in) / (mm)
M00545890	Hurricane RX Dilatation Balloon	4	2	5.8 / 1.93	12 / 1216	3.2	.035 / 0.89
M00545900	Hurricane RX Dilatation Balloon	4	4	5.8 / 1.93	12 / 1216	3.2	.035 / 0.89
M00545910	Hurricane RX Dilatation Balloon	6	2	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545920	Hurricane RX Dilatation Balloon	6	4	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545930	Hurricane RX Dilatation Balloon	8	2	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545940	Hurricane RX Dilatation Balloon	8	4	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545950	Hurricane RX Dilatation Balloon	10	2	5.8 / 1.93	8 / 810	3.2	.035 / 0.89
M00545960	Hurricane RX Dilatation Balloon	10	4	5.8 / 1.93	8 / 810	3.2	.035 / 0.89

Packaged one per box


Recommended Inflation Device: Encore™ 26 Inflation Device, M00566670

CRE™ RX Biliary Balloon Dilatation Catheters



Order Number	INFLATED BALLOON OD (mm)	INFLATED BALLOON OD (Fr)	Balloon Dilatation Length (cm)	INFLATION PRESSURE (ATM)	INFLATION PRESSURE (kPa)	Catheter Size (Fr)	Minimum Working Channel (mm)	Working Length (cm)
M00558920	8-9-10	24-27-30	3.0	3-5.5-9	304-608-1013	7.5	3.7	180
M00558930	10-11-12	30-33-36	3.0	3-5-8	304-557-912	7.5	3.7	180
M00558940	12-13.5-15	36-40.5-45	3.0	3-4.5-8	304-507-811	7.5	3.7	180
M00558860	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	7.5	3.7	180
M00558870	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	7.5	3.7	180
M00558880	10-11-12	30-33-36	5.5	3-5-8	304-507-811	7.5	3.7	180
M00558890	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	7.5	3.7	180
M00558900	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	7.5	3.7	180
M00558910	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	7.5	3.7	180

CYTOLOGY



RX Cytology Brush



Order Number	Description	Brush Width (mm)	Catheter OD (F) / (mm)	Working Length (cm)	Recommended Guidewire (in) / (mm)
M00545000	RX Cytology Brush	2.1	8 / 2.67	200	.035 / 0.89

Packaged one per box

LOCKING DEVICES

AutoCap™ RX




Order Number	Description	Units
M00545101	AutoCap RX Olympus Compatible (TJF-160, TJF-180, TJF-190 and JF-260V Series Duodenoscopes)	box 10
M00545131	AutoCap RX EXALT™ Model D Compatible Integrated Biopsy	box 10

RX Locking Devices



Order Number	Description	Units
M00545211†	RX Locking Device – Olympus / FujiFilm Brand Endoscope Compatible	box 10
M00545261	RX Locking Device and Biopsy Cap System – Olympus / FujiFilm G5 Series Brand Endoscope Compatible	box 10
M00545271	RX Locking Device and Biopsy Cap System – Pentax Brand Endoscope Compatible	box 10
M00545281	RX Locking Device and Biopsy Cap System – FujiFilm Type 3 Series Brand Endoscope Compatible	box 10
M00566001	Rapid Refill Continuous Injection System	box 10

† Product not available in the European Economic Area (EEA), except in Greece. Not for use or distribution in the EEA, except in Greece. Please check availability of the product with your local sales representative or your local customer service.



Encore™ 26 Inflation Device




Order Number	Description	Units
M00566670	Encore 26 Inflation Device	box 1

RX Biliary Plastic and Metal Stents are located in the Stents Section on pages 53–56 and 63–65

Biliary Guidewires are located in the Guidewires Section on pages 71–73


Biliary Lumen Patency



Biliary	RX Biliary System™	Biliary Lumen Patency	Endoscopic Ultrasound (EUS)	Guidewires	Biopsy	Dilation	Imaging Systems	Pulmonary	Closure	Haemostasis	Tissue Resection	Esophageal and Enteral Metal Stents	Infection Prevention	Endobariatrics
---------	--------------------	-----------------------	-----------------------------	------------	--------	----------	-----------------	-----------	---------	-------------	------------------	-------------------------------------	----------------------	----------------

Biliary Lumen Patency

RX BILIARY PLASTIC STENTS



Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System



Order Number	Description	Shape	OD (F) / (mm)	Length (cm) Barb to Barb
M00534200	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	5
M00534210	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	7
M00534220	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	9
M00534230	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	12
M00534240	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	15
M00534250	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	18
M00534260	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	5
M00534270	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	7
M00534280	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	9
M00534290	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	12
M00534300	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	15
M00534310	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	18
M00534320	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	5
M00534330	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	7
M00534340	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	9
M00534350	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	12
M00534360	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	15
M00534370	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	18
M00534560	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	5
M00534570	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	7
M00534580	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	9
M00534590	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	12
M00534600	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	15
M00534610	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	18
M00534620	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	5
M00534630	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	7
M00534640	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	9
M00534650	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	12
M00534660	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	15
M00534670	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	18
M00534680	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	5
M00534690	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	7
M00534700	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	9
M00534710	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	12
M00534720	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	15
M00534730	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	18

Packaged one per box


Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

Scope Compatibility: 7F (2.3mm) and 8.5F (2.83mm) stents require a 3.2mm working channel; 10F (3.3mm) stents require a 4.2mm working channel



Biliary Lumen Patency

RX BILIARY PLASTIC STENTS (CONTINUED)



Advanix™ Stents – Singles



Order Number	Description	Shape	OD (F) / (mm)	Length (cm) Barb to Barb
M00532870	Advanix Stent	Duodenal Bend	7 / 2.3	5
M00532880	Advanix Stent	Duodenal Bend	7 / 2.3	7
M00532890	Advanix Stent	Duodenal Bend	7 / 2.3	9
M00532900	Advanix Stent	Duodenal Bend	7 / 2.3	12
M00532910	Advanix Stent	Duodenal Bend	7 / 2.3	15
M00532920	Advanix Stent	Duodenal Bend	7 / 2.3	18
M00532930	Advanix Stent	Duodenal Bend	8.5 / 2.83	5
M00532940	Advanix Stent	Duodenal Bend	8.5 / 2.83	7
M00532950	Advanix Stent	Duodenal Bend	8.5 / 2.83	9
M00532960	Advanix Stent	Duodenal Bend	8.5 / 2.83	12
M00532970	Advanix Stent	Duodenal Bend	8.5 / 2.83	15
M00532980	Advanix Stent	Duodenal Bend	8.5 / 2.83	18
M00532990	Advanix Stent	Duodenal Bend	10 / 3.3	5
M00533000	Advanix Stent	Duodenal Bend	10 / 3.3	7
M00533010	Advanix Stent	Duodenal Bend	10 / 3.3	9
M00533020	Advanix Stent	Duodenal Bend	10 / 3.3	12
M00533030	Advanix Stent	Duodenal Bend	10 / 3.3	15
M00533040	Advanix Stent	Duodenal Bend	10 / 3.3	18
M00533230	Advanix Stent	Center Bend	7 / 2.3	5
M00533240	Advanix Stent	Center Bend	7 / 2.3	7
M00533250	Advanix Stent	Center Bend	7 / 2.3	9
M00533260	Advanix Stent	Center Bend	7 / 2.3	12
M00533270	Advanix Stent	Center Bend	7 / 2.3	15
M00533280	Advanix Stent	Center Bend	7 / 2.3	18
M00533290	Advanix Stent	Center Bend	8.5 / 2.83	5
M00533300	Advanix Stent	Center Bend	8.5 / 2.83	7
M00533310	Advanix Stent	Center Bend	8.5 / 2.83	9
M00533320	Advanix Stent	Center Bend	8.5 / 2.83	12
M00533330	Advanix Stent	Center Bend	8.5 / 2.83	15
M00533340	Advanix Stent	Center Bend	8.5 / 2.83	18
M00533350	Advanix Stent	Center Bend	10 / 3.3	5
M00533360	Advanix Stent	Center Bend	10 / 3.3	7
M00533370	Advanix Stent	Center Bend	10 / 3.3	9
M00533380	Advanix Stent	Center Bend	10 / 3.3	12
M00533390	Advanix Stent	Center Bend	10 / 3.3	15
M00533400	Advanix Stent	Center Bend	10 / 3.3	18

Packaged one per box


Recommended Guidewire: .035in (0.89mm)

Scope Compatibility: 7F (2.3mm) and 8.5F (2.83mm) stents require a 3.2mm working channel; 10F (3.3mm) stents require a 4.2mm working channel



Biliary Lumen Patency

RX BILIARY PLASTIC STENTS (CONTINUED)



Advanix™ Stents – Singles (CONTINUED)



Order Number	Description	Shape	OD (F) / (mm)	Length (cm) Between Pigtails
M00532160	Advanix Stent	Double Pigtail	7 / 2.3	3
M00532180	Advanix Stent	Double Pigtail	7 / 2.3	5
M00532190	Advanix Stent	Double Pigtail	7 / 2.3	7
M00532210	Advanix Stent	Double Pigtail	7 / 2.3	10
M00532220	Advanix Stent	Double Pigtail	7 / 2.3	12
M00532230	Advanix Stent	Double Pigtail	7 / 2.3	15
M00532240	Advanix Stent	Double Pigtail	10 / 3.3	3
M00532260	Advanix Stent	Double Pigtail	10 / 3.3	5
M00532270	Advanix Stent	Double Pigtail	10 / 3.3	7
M00532290	Advanix Stent	Double Pigtail	10 / 3.3	10
M00532300	Advanix Stent	Double Pigtail	10 / 3.3	12
M00532310	Advanix Stent	Double Pigtail	10 / 3.3	15

Packaged one per box

Recommended Guidewire: .035in (0.89mm)

Scope Compatibility: 7F (2.3mm) and 8.5F (2.83mm) stents require a 3.2mm working channel; 10F (3.3mm) stents require a 4.2mm working channel



NaviFlex™ RX Delivery System



Order Number	Description	Catheter Diameter (mm)	Working Length (cm)
M00533560	RX Delivery System	7 / 2.3	202.5
M00533570	RX Delivery System	8.5 / 2.83	202.5
M00533580	RX Delivery System	10 / 3.3	202.5

Packaged one per box

Compatible with Advanix Biliary Stents.



Biliary Lumen Patency

RX BILIARY PLASTIC STENTS (CONTINUED)

RX Biliary Plastic Stents with RX Delivery System



Order Number	Description	OD (F) / (mm)	Length (cm) Barb to Barb
M00545550	RX Biliary Plastic Stent with RX Delivery System	7 / 2.33	5
M00545560	RX Biliary Plastic Stent with RX Delivery System	7 / 2.33	7
M00545570	RX Biliary Plastic Stent with RX Delivery System	7 / 2.33	10
M00545580	RX Biliary Plastic Stent with RX Delivery System	7 / 2.33	12
M00545590	RX Biliary Plastic Stent with RX Delivery System	7 / 2.33	15
M00545650	RX Biliary Plastic Stent with RX Delivery System	8.5 / 2.83	5
M00545660	RX Biliary Plastic Stent with RX Delivery System	8.5 / 2.83	7
M00545670	RX Biliary Plastic Stent with RX Delivery System	8.5 / 2.83	10
M00545680	RX Biliary Plastic Stent with RX Delivery System	8.5 / 2.83	12
M00545690	RX Biliary Plastic Stent with RX Delivery System	8.5 / 2.83	15
M00545600	RX Biliary Plastic Stent with RX Delivery System	10 / 3.33	5
M00545610	RX Biliary Plastic Stent with RX Delivery System	10 / 3.33	7
M00545620	RX Biliary Plastic Stent with RX Delivery System	10 / 3.33	10
M00545630	RX Biliary Plastic Stent with RX Delivery System	10 / 3.33	12
M00545640	RX Biliary Plastic Stent with RX Delivery System	10 / 3.33	15
M00545700	RX Biliary Plastic Stent with RX Delivery System	11.5 / 3.83	5
M00545710	RX Biliary Plastic Stent with RX Delivery System	11.5 / 3.83	7
M00545720	RX Biliary Plastic Stent with RX Delivery System	11.5 / 3.83	10
M00545730	RX Biliary Plastic Stent with RX Delivery System	11.5 / 3.83	12
M00545740	RX Biliary Plastic Stent with RX Delivery System	11.5 / 3.83	15

Packaged one per box

Recommended Guidewire: .035in (0.89mm)



Biliary Lumen Patency

BILIARY PLASTIC STENTS



Flexima™ Biliary Stents with Delivery System




Order Number	Description	OD (F) / (mm)	Length (cm) Barb to Barb
M00539200	Flexima Stent with Delivery System	7 / 2.33	5
M00539210	Flexima Stent with Delivery System	7 / 2.33	7
M00539220	Flexima Stent with Delivery System	7 / 2.33	10
M00539230	Flexima Stent with Delivery System	7 / 2.33	12
M00539240	Flexima Stent with Delivery System	7 / 2.33	15
M00539250	Flexima Stent with Delivery System	8.5 / 2.83	5
M00539260	Flexima Stent with Delivery System	8.5 / 2.83	7
M00539270	Flexima Stent with Delivery System	8.5 / 2.83	10
M00539280	Flexima Stent with Delivery System	8.5 / 2.83	12
M00539290	Flexima Stent with Delivery System	8.5 / 2.83	15
M00539300	Flexima Stent with Delivery System	10 / 3.33	5
M00539310	Flexima Stent with Delivery System	10 / 3.33	7
M00539320	Flexima Stent with Delivery System	10 / 3.33	10
M00539330	Flexima Stent with Delivery System	10 / 3.33	12
M00539340	Flexima Stent with Delivery System	10 / 3.33	15
M00539350	Flexima Stent with Delivery System	11.5 / 3.83	5
M00539360	Flexima Stent with Delivery System	11.5 / 3.83	7
M00539370	Flexima Stent with Delivery System	11.5 / 3.83	10
M00539380	Flexima Stent with Delivery System	11.5 / 3.83	12
M00539390	Flexima Stent with Delivery System	11.5 / 3.83	15

Packaged one per box

Recommended Guidewire: .035in (0.89mm)



PANCREATIC PLASTIC STENTS



Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Pigtail



Order Number	Description	Leading Barb	Shape	OD (F) / (mm)	Length (cm)
M00537350	Advanix Pancreatic Stent Kit ●	No	Pigtail	4 / 1.35	3
M00537360	Advanix Pancreatic Stent Kit ▲	No	Pigtail	4 / 1.35	4
M00537370	Advanix Pancreatic Stent Kit ▲	No	Pigtail	4 / 1.35	5
M00537380	Advanix Pancreatic Stent Kit ▲	No	Pigtail	4 / 1.35	7
M00537390	Advanix Pancreatic Stent Kit ●	No	Pigtail	5 / 1.67	3
M00537400	Advanix Pancreatic Stent Kit ▲	No	Pigtail	5 / 1.67	4
M00537410	Advanix Pancreatic Stent Kit ▲	No	Pigtail	5 / 1.67	5
M00537420	Advanix Pancreatic Stent Kit ▲	No	Pigtail	5 / 1.67	7
M00537430	Advanix Pancreatic Stent Kit ●	Yes	Pigtail	5 / 1.67	3
M00537440	Advanix Pancreatic Stent Kit ▲	Yes	Pigtail	5 / 1.67	4
M00537450	Advanix Pancreatic Stent Kit ▲	Yes	Pigtail	5 / 1.67	5

* Products with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Straight



Order Number	Description	Leading Barb	Shape	OD (F) / (mm)	Length (cm)
M00537460	Advanix Pancreatic Stent Kit ● ¹	No	Straight	4 / 1.35	3
M00537470	Advanix Pancreatic Stent Kit ▲ ¹	No	Straight	4 / 1.35	4
M00537480	Advanix Pancreatic Stent Kit ▲ ¹	No	Straight	4 / 1.35	5
M00537490	Advanix Pancreatic Stent Kit ▲ ¹	No	Straight	4 / 1.35	7
M00537500	Advanix Pancreatic Stent Kit ●	No	Straight	5 / 1.67	3
M00537510	Advanix Pancreatic Stent Kit ▲	No	Straight	5 / 1.67	4
M00537520	Advanix Pancreatic Stent Kit ▲	No	Straight	5 / 1.67	5
M00537530	Advanix Pancreatic Stent Kit ▲	No	Straight	5 / 1.67	7
M00537540	Advanix Pancreatic Stent Kit ●	Yes	Straight	5 / 1.67	3
M00537550	Advanix Pancreatic Stent Kit ▲	Yes	Straight	5 / 1.67	4
M00537560	Advanix Pancreatic Stent Kit ▲	Yes	Straight	5 / 1.67	5

* Products with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker

¹ Stiffer stent material compared to 4F pigtail stents.

² Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.



PANCREATIC PLASTIC STENTS (CONTINUED)



Advanix™ Pancreatic Single Stents – Pigtail

Order Number	Description	Leading Barb	Shape	OD (F) / (mm)	Length (cm)
M00536740	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	2
M00536750	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	3
M00536760	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	4
M00536770	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	5
M00536780	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	6
M00536790	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	7
M00536800	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	8
M00536810	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	9
M00536820	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	10
M00536830	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	11
M00536840	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	12
M00536850	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	13
M00536860	Advanix Pancreatic Stent ●	No	Pigtail	4 / 1.35	2
M00536870	Advanix Pancreatic Stent ●	No	Pigtail	4 / 1.35	3
M00536880	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	4
M00536890	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	5
M00536900	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	6
M00536910	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	7
M00536920	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	8
M00536930	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	9
M00536940	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	10
M00536950	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	11
M00536960	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	12
M00536970	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	13
M00536980	Advanix Pancreatic Stent ●	No	Pigtail	5 / 1.67	2
M00536990	Advanix Pancreatic Stent ●	No	Pigtail	5 / 1.67	3
M00537000	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	4
M00537010	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	5
M00537020	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	6
M00537030	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	7
M00537040	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	8
M00537050	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	9
M00537060	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	10
M00537070	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	11
M00537080	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	12
M00537090	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	13

* Products with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker

1 Stiffer stent material compared to 4F pigtail stents.

2 Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.



Biliary Lumen Patency

PANCREATIC PLASTIC STENTS (CONTINUED)



Advanix™ Pancreatic Single Stents – Pigtail (CONTINUED)

Order Number	Description	Leading Barb	Shape	OD (F) / (mm)	Length (cm)
M00537100	Advanix Pancreatic Stent ●	Yes	Pigtail	5 / 1.67	2
M00537110	Advanix Pancreatic Stent ●	Yes	Pigtail	5 / 1.67	3
M00537120	Advanix Pancreatic Stent ▲	Yes	Pigtail	5 / 1.67	4
M00537130	Advanix Pancreatic Stent ▲	Yes	Pigtail	5 / 1.67	5
M00537140	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	6
M00537150	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	7
M00537160	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	8
M00537170	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	9
M00537180	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	10
M00537190	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	11
M00537200	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	12
M00537210	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	13
M00537220	Advanix Pancreatic Stent ■	Yes	Pigtail	7 / 2.33	3
M00537230	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	4
M00537240	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	5
M00537250	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	6
M00537260	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	7
M00537270	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	8
M00537280	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	9
M00537290	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	10
M00537300	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	11
M00537310	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	12
M00537320	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	13
M00537330	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	15
M00537340	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	18

* Products with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker

¹ Stiffer stent material compared to 4F pigtail stents.

² Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.



PANCREATIC PLASTIC STENTS (CONTINUED)



Advanix™ Pancreatic Single Stents – Straight



Order Number	Description	Leading Barb	Shape	OD (F) / (mm)	Length (cm)
M00536000	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	2
M00536010	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	3
M00536020	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	4
M00536030	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	5
M00536040	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	6
M00536050	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	7
M00536060	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	8
M00536070	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	9
M00536080	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	10
M00536090	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	11
M00536100	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	12
M00536110	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	13
M00536120	Advanix Pancreatic Stent ● ¹	No	Straight	4 / 1.35	2
M00536130	Advanix Pancreatic Stent ● ¹	No	Straight	4 / 1.35	3
M00536140	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	4
M00536150	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	5
M00536160	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	6
M00536170	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	7
M00536180	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	8
M00536190	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	9
M00536200	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	10
M00536210	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	11
M00536220	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	12
M00536230	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	13
M00536240	Advanix Pancreatic Stent ●	No	Straight	5 / 1.67	2
M00536250	Advanix Pancreatic Stent ●	No	Straight	5 / 1.67	3
M00536260	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	4
M00536270	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	5
M00536280	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	6
M00536290	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	7
M00536300	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	8
M00536310	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	9
M00536320	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	10
M00536330	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	11
M00536340	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	12
M00536350	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	13

* Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker


¹ Stiffer stent material compared to 4F pigtail stents.

² Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.



Biliary Lumen Patency

PANCREATIC PLASTIC STENTS (CONTINUED)



Advanix™ Pancreatic Single Stents – Straight (CONTINUED)



Order Number	Description	Leading Barb	Shape	OD (F) / (mm)	Length (cm)
M00536360	Advanix Pancreatic Stent ●	Yes	Straight	5 / 1.67	2
M00536370	Advanix Pancreatic Stent ●	Yes	Straight	5 / 1.67	3
M00536380	Advanix Pancreatic Stent ▲	Yes	Straight	5 / 1.67	4
M00536390	Advanix Pancreatic Stent ▲	Yes	Straight	5 / 1.67	5
M00536400	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	6
M00536410	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	7
M00536420	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	8
M00536430	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	9
M00536440	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	10
M00536450	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	11
M00536460	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	12
M00536470	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	13
M00536480	Advanix Pancreatic Stent ■	Yes	Straight	7 / 2.33	3
M00536490	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	4
M00536500	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	5
M00536510	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	6
M00536520	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	7
M00536530	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	8
M00536540	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	9
M00536550	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	10
M00536560	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	11
M00536570	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	12
M00536580	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	13
M00536590	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	15
M00536600	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	18
M00536610	Advanix Pancreatic Stent ■	Yes	Straight	10 / 3.33	3
M00536620	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	4
M00536630	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	5
M00536640	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	6
M00536650	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	7
M00536660	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	8
M00536670	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	9
M00536680	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	10
M00536690	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	11
M00536700	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	12
M00536710	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	13
M00536720	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	15
M00536730	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	18

* Products with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker


¹ Stiffer stent material compared to 4F pigtail stents.

² Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.



Biliary Lumen Patency

PANCREATIC PLASTIC STENTS (CONTINUED)




NaviFlex™ RX Pancreatic Delivery System & NaviFlex™ Pushers



Order Number	Description	Diameter (F) / (mm)
M00535000	Pancreatic Stent RX Pusher	3 / 1.00
M00535010	Pancreatic Stent RX Pusher	4 & 5 / 1.35 & 1.67
M00535020	Pancreatic Stent Pusher	3 / 1.00
M00535030	Pancreatic Stent Pusher	4 & 5 / 1.35 & 1.67
M00535060	NaviFlex RX Pancreatic Stent Delivery System	7 / 2.33
M00535080	NaviFlex RX Pancreatic Stent Delivery System	10 / 3.33

RX BILIARY METAL STENTS



WallFlex™ Biliary RX Fully Covered Stent System RMV (Permalume™ Covering, includes benign indication)




Order Number	Description	Diameter (mm)	Length (mm)	Catheter Diameter (F) / (mm)	Guidewire Diameter (in) / (mm)
M00574210†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	6	40	8.5 / 2.83	.035 / 0.89
M00574230†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	6	60	8.5 / 2.83	.035 / 0.89
M00574240†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	6	80	8.5 / 2.83	.035 / 0.89
M00570450†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	8	60	8.5 / 2.83	.035 / 0.89
M00570460†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	8	80	8.5 / 2.83	.035 / 0.89
M00576800†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	8	100	9.0 / 3.00	.035 / 0.89
M00576810†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	8	120	9.0 / 3.00	.035 / 0.89
M00570470†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	10	40	8.5 / 2.83	.035 / 0.89
M00570480†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	10	60	8.5 / 2.83	.035 / 0.89
M00570490†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	10	80	8.5 / 2.83	.035 / 0.89
M00576820†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	10	100	9.0 / 3.00	.035 / 0.89
M00576830†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	10	120	9.0 / 3.00	.035 / 0.89

† Product commercialised only in countries where benign indication is approved.



Biliary Lumen Patency

RX BILIARY METAL STENTS (CONTINUED)





WallFlex™ Biliary RX Partially Covered Stent System (Permalume™ Covering)



Order Number	Description	Diameter (mm)	Length (mm)	Covered Length (mm)	Catheter Diameter (F) / (mm)	Guidewire Diameter (in) / (mm)
M00570700	WallFlex RX Biliary (Permalume Covering)	8	60	48	8.5 / 2.83	.035 / 0.89
M00570710	WallFlex RX Biliary (Permalume Covering)	8	80	68	8.5 / 2.83	.035 / 0.89
M00576740	WallFlex RX Biliary (Permalume Covering)	8	100	88	9.0 / 3.00	.035 / 0.89
M00576750	WallFlex RX Biliary (Permalume Covering)	8	120	108	9.0 / 3.00	.035 / 0.89
M00570720	WallFlex RX Biliary (Permalume Covering)	10	40	28	8.5 / 2.83	.035 / 0.89
M00570730	WallFlex RX Biliary (Permalume Covering)	10	60	48	8.5 / 2.83	.035 / 0.89
M00570740	WallFlex RX Biliary (Permalume Covering)	10	80	68	8.5 / 2.83	.035 / 0.89
M00576760	WallFlex RX Biliary (Permalume Covering)	10	100	88	9.0 / 3.00	.035 / 0.89
M00576770	WallFlex RX Biliary (Permalume Covering)	10	120	108	9.0 / 3.00	.035 / 0.89

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm or .035in (0.89mm) Dreamwire™ High Performance Guidewire, 260cm



WallFlex™ Biliary RX Uncovered Stent System



Order Number	Description	Diameter (mm)	Length (mm)	Catheter Diameter (F) / (mm)	Guidewire Diameter (in) / (mm)
M00570600	WallFlex RX Biliary	8	40	8 / 2.67	.035 / 0.89
M00570610	WallFlex RX Biliary	8	60	8 / 2.67	.035 / 0.89
M00570620	WallFlex RX Biliary	8	80	8 / 2.67	.035 / 0.89
M00570630	WallFlex RX Biliary	8	100	8 / 2.67	.035 / 0.89
M00576780	WallFlex RX Biliary	8	120	8 / 2.67	.035 / 0.89
M00570890	WallFlex RX Biliary	10	40	8 / 2.67	.035 / 0.89
M00570640	WallFlex RX Biliary	10	60	8 / 2.67	.035 / 0.89
M00570650	WallFlex RX Biliary	10	80	8 / 2.67	.035 / 0.89
M00570660	WallFlex RX Biliary	10	100	8 / 2.67	.035 / 0.89
M00576790	WallFlex RX Biliary	10	120	8 / 2.67	.035 / 0.89

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm or .035in (0.89mm) Dreamwire™ High Performance Guidewire, 260cm

Biliary Lumen Patency

RX BILIARY METAL STENTS (CONTINUED)

WallFlex™ Biliary RX Fully Covered Stent System



Order Number	Description	Delivery System Working Length (cm)	Stent Diameter (mm)	Stent Length (mm)
M00570500	WallFlex Biliary Stent System (RX Fully Covered)	194	8	60
M00570510	WallFlex Biliary Stent System (RX Fully Covered)	194	8	80
M00576700	WallFlex Biliary Stent System (RX Fully Covered)	194	8	100
M00576710	WallFlex Biliary Stent System (RX Fully Covered)	194	8	120
M00570520	WallFlex Biliary Stent System (RX Fully Covered)	194	10	40
M00570530	WallFlex Biliary Stent System (RX Fully Covered)	194	10	60
M00570540	WallFlex Biliary Stent System (RX Fully Covered)	194	10	80
M00576720	WallFlex Biliary Stent System (RX Fully Covered)	194	10	100
M00576730	WallFlex Biliary Stent System (RX Fully Covered)	194	10	120

Packaged one per box

Limited availability in all European countries (except UK and France). Please contact your local Territory Manager.

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621

WallFlex™ Biliary PLUS RX Fully Covered Stent System




Order Number	Description	Delivery System Working Length (cm)	Stent Diameter (mm)	Stent Length (mm)
M005747500	WallFlex™ Biliary PLUS RX Fully Covered Stent System	194	8	60
M005747510	WallFlex™ Biliary PLUS RX Fully Covered Stent System	194	8	80
M005747520	WallFlex™ Biliary PLUS RX Fully Covered Stent System	194	8	100
M005747530	WallFlex™ Biliary PLUS RX Fully Covered Stent System	194	10	40
M005747540	WallFlex™ Biliary PLUS RX Fully Covered Stent System	194	10	60
M005747550	WallFlex™ Biliary PLUS RX Fully Covered Stent System	194	10	80
M005747560	WallFlex™ Biliary PLUS RX Fully Covered Stent System	194	10	100




RADIOFREQUENCY ABLATION CATHETER AND CABLES

Habib™ EndoHPB Bipolar Radiofrequency Catheter




Order Number	Description
M00500070	Habib Endo HPB Bipolar Radiofrequency Catheter



Order Number	Description
5100	Adapter Cable Bipolar, 28 mm pins



Order Number	Description
5700	Adapter Cable Bipolar, 22 mm pins



Order Number	Description
5420	Adapter Cable Bipolar, Internal 4 mm pin



Biliary
RX Biliary™ System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics

FINE NEEDLE ASPIRATION (FNA)



Expect™ Slimline Handle Needles*



Order Number	Description	Needle Size (gauge) / (mm)	Minimum Working Channel (mm)	Sheath Diameter (mm)	Units
M00555500	Expect Slimline Handle Needle	19 / 1.10	2.8	1.83	box 1
M00555501	Expect Slimline Handle Needle	19 / 1.10	2.8	1.83	box 5
M00555530	Expect Slimline Handle Needle	19 Flex / 1.14	2.8	1.73	box 1
M00555531	Expect Slimline Handle Needle	19 Flex / 1.14	2.8	1.73	box 5
M00555510	Expect Slimline Handle Needle	22 / 0.72	2.4	1.65	box 1
M00555511	Expect Slimline Handle Needle	22 / 0.72	2.4	1.65	box 5
M00555520	Expect Slimline Handle Needle	25 / 0.52	2.4	1.52	box 1
M00555521	Expect Slimline Handle Needle	25 / 0.52	2.4	1.52	box 5

* Includes a 20cc syringe and one-way stopcock
Working length: 137.5cm to 141.5cm, adjustable
Needle length: 0cm to 8cm, adjustable

FINE NEEDLE BIOPSY (FNB)



Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Devices



Order Number	Description	Needle Size (gauge) / (mm)	Minimum Working Channel (mm)	Sheath Diameter (mm)	Units
M00555580	19ga Acquire FNB Needle, Flexible	19 / 1.14	2.8	1.73	box 1
M00555540	22ga Acquire EUS-FNB Needle	22 / 0.72	2.4	1.65	box 1
M00555541	22ga Acquire EUS-FNB Needle	22 / 0.72	2.4	1.65	box 5
M00555560	25ga Acquire EUS-FNB Needle	25 / 0.52	2.4	1.52	box 1
M00555561	25ga Acquire EUS-FNB Needle	25 / 0.52	2.4	1.52	box 5

FINE NEEDLE BIOPSY (FNB)



Acquire™ S Endoscopic Ultrasound Fine Needle Biopsy (FNB) Devices



Order Number	Description	Needle Size (gauge) / (mm)	Minimum Working Channel (mm)	Sheath Diameter (mm)	Units
M00555640	22ga Acquire S EUS-FNB Needle	22 / 0.72	2.4	1.65	box 1
M00555660	25ga Acquire S EUS-FNB Needle	25 / 0.52	2.4	1.52	box 1
M00555641	22ga Acquire S EUS-FNB Needle	22 / 0.72	2.4	1.65	box 5
M00555661	25ga Acquire S EUS-FNB Needle	25 / 0.52	2.4	1.52	box 5


LUMEN APPOSING METAL STENT SYSTEM

Hot AXIOS™ Stent and Electrocautery Enhanced Delivery System



Order Number	Description	Flange Diameter (mm)	Lumen Diameter (mm)	Saddle Length (mm)	Catheter Outer Diameter (F)	Catheter Working Length (cm)	Catheter Total Length (cm)	Units
M00553520	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	14	6	8	9	138	146	box 1
M00553530	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	17	8	8	9	138	146	box 1
M00553540	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	21	10	10	10.8	138	146	box 1
M00553550	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	24	15	10	10.8	138	146	box 1
M00553560	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	29	20	10	10.8	138	146	box 1
M00553570*	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	24	15	15	10.8	138	144	box 1

*The 15x15 HOT AXIOS Stent and Electrocautery Enhanced Delivery System is indicated for use to facilitate transgastric or transduodenal endoscopic drainage of a pancreatic pseudocyst or a walled-off necrosis with ≥70 % fluid content.



Guidewires

Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics

HIGH PERFORMANCE




Dreamwire™ High Performance Guidewires




Order Number	Description	OD (in) / (mm)	Length (cm)	Dream Tip Style	Units
M00556101	Dreamwire Standard	.035 / 0.89	260	Straight	box 2
M00556111	Dreamwire Standard	.035 / 0.89	260	Angled	box 2
M00556121	Dreamwire Stiff Shaft	.035 / 0.89	260	Straight	box 2
M00556131	Dreamwire Stiff Shaft	.035 / 0.89	260	Angled	box 2
M00556141	Dreamwire Standard	.035 / 0.89	450	Straight	box 2
M00556151	Dreamwire Standard	.035 / 0.89	450	Angled	box 2
M00556161	Dreamwire Stiff Shaft	.035 / 0.89	450	Straight	box 2
M00556171	Dreamwire Stiff Shaft	.035 / 0.89	450	Angled	box 2


5cm Jagwire™ Guidewire Tip



10 cm Straight Dream Tip™



10 cm Angled Dream Tip™




Hydra Jagwire™ High Performance Guidewires



Order Number	Description	OD (in) / (mm)	Length (cm)	Dream Tip Style	Units
M00556001	Hydra Jagwire Standard	.035 / 0.89	260	Straight	box 2
M00556011	Hydra Jagwire Standard	.035 / 0.89	260	Angled	box 2
M00556021	Hydra Jagwire Stiff Shaft	.035 / 0.89	260	Straight	box 2
M00556031	Hydra Jagwire Stiff Shaft	.035 / 0.89	260	Angled	box 2
M00556041	Hydra Jagwire Standard	.035 / 0.89	450	Straight	box 2
M00556051	Hydra Jagwire Standard	.035 / 0.89	450	Angled	box 2
M00556061	Hydra Jagwire Stiff Shaft	.035 / 0.89	450	Straight	box 2
M00556071	Hydra Jagwire Stiff Shaft	.035 / 0.89	450	Angled	box 2

Hydra Jagwire Guidewire may remain in place during sphincterotomy. See "Instructions for Use" for details.

HIGH PERFORMANCE (CONTINUED)



Jagwire™ High Performance Guidewires



Order Number	Description	Distal Tip OD (in)	Working Length (cm)	Shaft Style	Tip Style	Units
M00556511	Jagwire Guidewire ANG	.025	450	Standard	Angled	box 2
M00556521	Jagwire Guidewire ST Stiff	.025	450	Stiff	Straight	box 2
M00556531	Jagwire Guidewire ANG Stiff	.025	450	Stiff	Angled	box 2
M00556561	Jagwire Guidewire ST	.025	450	Standard	Straight	box 2
M00556581	Jagwire Guidewire ST	.035	450	Standard	Straight	box 2
M00556591	Jagwire Guidewire ANG	.035	450	Standard	Angled	box 2
M00556601	Jagwire Guidewire ST Stiff	.035	450	Stiff	Straight	box 2
M00556611	Jagwire Guidewire ANG Stiff	.035	450	Stiff	Angled	box 2
M00556621	Jagwire Guidewire ST	.038	260	Standard	Straight	box 2
M00557661	Jagwire Guidewire ST Stiff	.035	260	Stiff	Straight	box 2
M00557641	Jagwire Guidewire ST	.035	260	Standard	Straight	box 2
M00557651	Jagwire Guidewire ANG	.035	260	Standard	Angled	box 2
M00557671	Jagwire Guidewire ANG Stiff	.035	260	Stiff	Angled	box 2
M00556461	Jagwire Guidewire ST	.025	260	Standard	Straight	box 2
M00556481	Jagwire Guidewire ST Stiff	.025	260	Stiff	Straight	box 2
M00556491	Jagwire Guidewire ANG Stiff	.025	260	Stiff	Angled	box 2
M00556841	Extendable Jagwire Guidewire ST	.035	260	Standard	Straight	box 2
M00556851	Extendable Jagwire Guidewire ANG	.035	260	Standard	Angled	box 2
M00556861	Extendable Jagwire Guidewire ST Stiff	.035	260	Stiff	Straight	box 2
M00556871	Extendable Jagwire Guidewire ANG Stiff	.035	260	Stiff	Angled	box 2
M00556902	Jagtail™ Guidewire Extension	.035	200	-	Guidewire Extension	box 5

Jagwire Guidewire may remain in place during sphincterotomy. See "Instructions for Use" for details. ST = Straight, ANG = Angled


Jagwire™ Revolution High Performance Guidewire



Order Number	Description	OD (in) / (mm)	Length (cm)	Tip Style	Units
M00557001	Jagwire Revolution High Performance Guidewire	.025 / 0.64	260	Straight Round	box 2
M00557011	Jagwire Revolution High Performance Guidewire	.025 / 0.64	260	Angled Round	box 2
M00557021	Jagwire Revolution High Performance Guidewire	.025 / 0.64	450	Straight Round	box 2
M00557031	Jagwire Revolution High Performance Guidewire	.025 / 0.64	450	Angled Round	box 2
M00557041	Jagwire Revolution High Performance Guidewire	.025 / 0.64	500	Straight Round	box 2

Guidewires

HIGH PERFORMANCE (CONTINUED)



NaviPro™ Hydrophilic Guidewires



Order Number	Description	OD (in) / (mm)	Length (cm)	Tip Style	Units
M00556191	NaviPro Standard	.018 / 0.46	260	Angled	box 5
M00556181	NaviPro Standard	.018 / 0.46	260	Straight	box 5
M00556211	NaviPro Standard	.025 / 0.64	260	Angled	box 5
M00556201	NaviPro Standard	.025 / 0.64	260	Straight	box 5
M00556231	NaviPro Standard	.035 / 0.89	260	Angled	box 5
M00556221	NaviPro Standard	.035 / 0.89	260	Straight	box 5
M00556251	NaviPro Stiff Shaft	.035 / 0.89	260	Angled	box 5
M00556241	NaviPro Stiff Shaft	.035 / 0.89	260	Straight	box 5

STANDARD

Amplatz Super Stiff™ Guidewire




Order Number	Description	OD (in) / (mm)	Length (cm)	Tip Style	Units
M00550090	Amplatz Super Stiff Guidewire	.038 / 0.97	260	Straight	box 1





Biliary
RX Biliary™ System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics

FORCEPS



Radial Jaw™ 4 Standard Capacity Biopsy Forceps



Order Number	Description	Jaw OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00513383	Radial Jaw 4 Standard Capacity	2.2	160	2.8	Yellow	box 5
M00513381	Radial Jaw 4 Standard Capacity	2.2	160	2.8	Yellow	box 20
M00513382	Radial Jaw 4 Standard Capacity	2.2	160	2.8	Yellow	box 40
M00513393	Radial Jaw 4 Standard Capacity w/Needle	2.2	160	2.8	Yellow	box 5
M00513391	Radial Jaw 4 Standard Capacity w/Needle	2.2	160	2.8	Yellow	box 20
M00513392	Radial Jaw 4 Standard Capacity w/Needle	2.2	160	2.8	Yellow	box 40
M00513403	Radial Jaw 4 Standard Capacity	2.2	240	2.8	Orange	box 5
M00513401	Radial Jaw 4 Standard Capacity	2.2	240	2.8	Orange	box 20
M00513402	Radial Jaw 4 Standard Capacity	2.2	240	2.8	Orange	box 40
M00513413	Radial Jaw 4 Standard Capacity w/Needle	2.2	240	2.8	Orange	box 5
M00513411	Radial Jaw 4 Standard Capacity w/Needle	2.2	240	2.8	Orange	box 20
M00513412	Radial Jaw 4 Standard Capacity w/Needle	2.2	240	2.8	Orange	box 40

Radial Jaw™ 4 Large Capacity Biopsy Forceps



Order Number	Description	Jaw OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00513303	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	box 5
M00513301	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	box 20
M00513302	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	box 40
M00513313	Radial Jaw 4 Large Capacity w/Needle	2.4	160	2.8	Yellow	box 5
M00513311	Radial Jaw 4 Large Capacity w/Needle	2.4	160	2.8	Yellow	box 20
M00513312	Radial Jaw 4 Large Capacity w/Needle	2.4	160	2.8	Yellow	box 40
M00513323	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	box 5
M00513321	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	box 20
M00513322	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	box 40
M00513333	Radial Jaw 4 Large Capacity w/Needle	2.4	240	2.8	Orange	box 5
M00513331	Radial Jaw 4 Large Capacity w/Needle	2.4	240	2.8	Orange	box 20
M00513332	Radial Jaw 4 Large Capacity w/Needle	2.4	240	2.8	Orange	box 40

FORCEPS (CONTINUED)

Radial Jaw™ 4 Jumbo Biopsy Forceps




Order Number	Description	Jaw OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00513343	Radial Jaw 4 Jumbo	2.8	160	3.2	Yellow	box 5
M00513353	Radial Jaw 4 Jumbo w/Needle	2.8	160	3.2	Yellow	box 5
M00513363	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	box 5
M00513361	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	box 20
M00513362	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	box 40
M00513373	Radial Jaw 4 Jumbo w/Needle	2.8	240	3.2	Orange	box 5
M00513371	Radial Jaw 4 Jumbo w/Needle	2.8	240	3.2	Orange	box 20
M00513372	Radial Jaw 4 Jumbo w/Needle	2.8	240	3.2	Orange	box 40

Radial Jaw™ 4 Gastropediatric Biopsy Forceps




Order Number	Description	Jaw OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00513443	Radial Jaw 4 Gastropediatric	1.8	160	2.0	Yellow	box 5
M00513441	Radial Jaw 4 Gastropediatric	1.8	160	2.0	Yellow	box 20
M00513453	Radial Jaw 4 Gastropediatric w/Needle	1.8	160	2.0	Yellow	box 5
M00513451	Radial Jaw 4 Gastropediatric w/Needle	1.8	160	2.0	Yellow	box 20


Dilatation



Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics



CRE™ (CONTROLLED RADIAL EXPANSION) BALLOON DILATATION CATHETERS



CRE™ PRO Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic/Biliary) (including DASE indication)



Order Number	BALLOON INFLATED OD (mm) (Fr)		Balloon Length (cm)	INFLATION PRESSURE (ATM) (kPa)		Minimum Working Channel (mm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Units
M00558660	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8	7.5 / 2.50	240	box 1
M00558670	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8	7.5 / 2.50	240	box 1
M00558680	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8	7.5 / 2.50	240	box 1
M00558690	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8	7.5 / 2.50	240	box 1
M00558700	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8	7.5 / 2.50	240	box 1
M00558710	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8	7.5 / 2.50	240	box 1


Recommended Inflation Device: Alliance™ II Integrated Inflation/Lithotripsy Device (Re-usable Handle) M00550620, Alliance™ II Syringe/Gauge Assembly M00550601.
This product contains no detectable latex.

CRE™ Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)




Order Number	BALLOON INFLATED OD (mm) (Fr)		Balloon Length (cm)	INFLATION PRESSURE (ATM) (kPa)		Minimum Working Channel (mm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Units
M00558450	6-7-8	18-21-24	5.5	3-6-10	304-557-912	2.8	7.5 / 2.50	240	box 1
M00558460	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8	7.5 / 2.50	240	box 1
M00558470	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8	7.5 / 2.50	240	box 1
M00558480	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8	7.5 / 2.50	240	box 1
M00558490	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8	7.5 / 2.50	240	box 1
M00558500	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8	7.5 / 2.50	240	box 1

Recommended Inflation Device: Alliance™ II Integrated Inflation/Lithotripsy Device (Re-usable Handle) M00550620, Alliance™ II Syringe/Gauge Assembly M00550601.
This product contains no detectable latex.



CRE™ (CONTROLLED RADIAL EXPANSION) BALLOON DILATATION CATHETERS (CONTINUED)



CRE™ PRO Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Biliary) (including DASE indication)



Order Number	BALLOON INFLATED OD (mm) (Fr)		Balloon Length (cm)	INFLATION PRESSURE (ATM) (kPa)		Minimum Working Channel (mm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Units
M00558600	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8	7.5 / 2.50	180	box 1
M00558610	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8	7.5 / 2.50	180	box 1
M00558620	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8	7.5 / 2.50	180	box 1
M00558630	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8	7.5 / 2.50	180	box 1
M00558640	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8	7.5 / 2.50	180	box 1
M00558650	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8	7.5 / 2.50	180	box 1


Recommended Inflation Device: Alliance™ II Integrated Inflation/Lithotripsy Device (Re-usable Handle) M00550620, Alliance™ II Syringe/Gauge Assembly M00550601.
This product contains no detectable latex.

CRE™ Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)




Order Number	BALLOON INFLATED OD (mm) (Fr)		Balloon Length (cm)	INFLATION PRESSURE (ATM) (kPa)		Minimum Working Channel (mm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Units
M00558390	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8	7.5 / 2.50	180	box 1
M00558400	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8	7.5 / 2.50	180	box 1
M00558410	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8	7.5 / 2.50	180	box 1
M00558420	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8	7.5 / 2.50	180	box 1
M00558430	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8	7.5 / 2.50	180	box 1
M00558440	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8	7.5 / 2.50	180	box 1

Recommended Inflation Device: Alliance™ II Integrated Inflation/Lithotripsy Device (Re-usable Handle) M00550620, Alliance™ II Syringe/Gauge Assembly M00550601.
This product contains no detectable latex.



CRE™ (CONTROLLED RADIAL EXPANSION) BALLOON DILATATION CATHETERS (CONTINUED)



CRE™ Fixed Wire Balloon Dilatation Catheters

Order Number	Inflated OD (mm)	Balloon Length (cm)	INFLATION PRESSURE (ATM)	INFLATION PRESSURE (kPa)	Minimum Working Channel (mm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Units
M00558330	6-7-8	8	3-6-10	304-608-1013	2.8	6 / 2.00	180	box 1
M00558340	8-9-10	8	3-5.5-9	304-557-912	2.8	6 / 2.00	180	box 1
M00558350	10-11-12	8	3-5-8	304-507-810	2.8	6 / 2.00	180	box 1
M00558360	12-13.5-15	8	3-4.5-8	304-456-810	2.8	6 / 2.00	180	box 1
M00558361	12-13.5-15	8	3-4.5-8	304-456-810	2.8	6 / 2.00	180	box 5
M00558370	15-16.5-18	8	3-4.5-7	304-456-709	2.8	6 / 2.00	180	box 1
M00558371	15-16.5-18	8	3-4.5-7	304-456-709	2.8	6 / 2.00	180	box 5
M00558380	18-19-20	8	3-4.5-6	304-456-608	2.8	6 / 2.00	180	box 1
M00558381	18-19-20	8	3-4.5-6	304-456-608	2.8	6 / 2.00	180	box 5

Recommended Inflation Device: Alliance™ II Integrated Inflation/Lithotripsy Device (Re-usable Handle) M00550620, Alliance™ II Syringe/Gauge Assembly M00550601.
This product contains no detectable latex.

BALLOON DILATATORS


Rigiflex™ II Achalasia Balloon Dilators*

Order Number	Inflated OD (mm)	Balloon Length (cm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Recommended Guidewire (in) / (mm)**	Units
M00554500	30	10	14 / 4.67	90	.038 / 0.97	box 1
M00554510	35	10	14 / 4.67	90	.038 / 0.97	box 1
M00554520	40	10	14 / 4.67	90	.038 / 0.97	box 1


* Product shipped sterile. Recommended use with pneumatic inflator, M00553200

** Guidewire included with Rigiflex™ II Balloon Dilators

This product contains no detectable latex.



INFLATION SYSTEMS



Alliance™ II Inflation System


Order Number	Description	Units
M00550620	Alliance II Integrated Inflation/Lithotripsy Device	box 1
M00550601	Alliance II Syringe/Gauge Assembly	box 5

DILATATION ACCESSORIES

Inflation/Deflation Devices and Accessories

Order Number	Description	Units
M00553200	Pneumatic Inflator	box 1
M00550630	CRE™ SteriFlate Disposable Inflation Device	box 1

RX Biliary Balloon Dilators are located in the RX Biliary Section on page 50



Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics

IMAGING SYSTEMS

EXALT™ Model D Single-Use Duodenoscope



Order Number	Item	Description	Languages	Units
M0054242CE0	EXALT Model D Single-Use Duodenoscope	Sterile Single-Use Duodenoscope for use with EXALT Controller	N/A	box 1
M0054242CE1	EXALT Model D Single-Use Duodenoscope	Sterile Single-Use Duodenoscope for use with EXALT Controller	N/A	box 2
M00542430	EXALT Controller	100–240 VAC Controller for use with Single-Use Duodenoscope (including detachable power cable, DVI and HDMI cables)	English, Spanish, French, German, Italian, Dutch and Portuguese	box 1
M00542460	EXALT Controller	100–240 VAC Controller for use with Single-Use Duodenoscope (including detachable power cable, DVI and HDMI cables)	English, Danish, Greek, Portuguese-Brazil, Swedish, Hungarian, Czech, Polish, Norwegian, Korean, Turkish, Finnish, Romanian and Slovak	box 1

EXALT™ Model D Service Codes


EXALT Model D TotalCare POS

Order Number	Description
SRVM005EXALTDTC01	EXALT Model D TotalCare POS 1 Year
SRVM005EXALTDTC02	EXALT Model D TotalCare POS 2 Year
SRVM005EXALTDTC03	EXALT Model D TotalCare POS 3 Year
SRVM005EXALTDTC04	EXALT Model D TotalCare POS 4 Year
SRVM005EXALTDTC05	EXALT Model D TotalCare POS 5 Year

EXALT Model D TotalCare Post Sale

Order Number	Description
SRVM005EXALTDTC06	EXALT Model D TotalCare POST 1 Year
SRVM005EXALTDTC07	EXALT Model D TotalCare POST 2 Year
SRVM005EXALTDTC08	EXALT Model D TotalCare POST 3 Year
SRVM005EXALTDTC09	EXALT Model D TotalCare POST 4 Year
SRVM005EXALTDTC10	EXALT Model D TotalCare POST 5 Year

DIRECT VISUALISATION SYSTEM



SpyGlass™ DS Capital Equipment

Order Number	Description	Units
M00546650	SpyGlass DS Digital Controller	box 1

SpyGlass™ DS Disposable Equipment



Order Number	Description	Units
M00546610	SpyScope™ DS II Access and Delivery Catheter	box 1
M00546550	SpyGlass™ Retrieval Basket	box 1
M00546560	SpyGlass™ Retrieval Snare	box 1
M00546470	SpyBite™ Max Biopsy Forceps	box 1
SIT-575-15	Hydra Irrigation Tubing	box 15

SpyGlass™ DS Service Codes

SpyGlass DS Digital EverCare POS

Order Number	Description
SRVM005SPYDIG1009	SpyGlass Digital EverCare POS 1 Year
SRVM005SPYDIG1010	SpyGlass Digital EverCare POS 2 Years
SRVM005SPYDIG1011	SpyGlass Digital EverCare POS 3 Years
SRVM005SPYDIG1012	SpyGlass Digital EverCare POS 4 Years

SpyGlass DS Digital EverCare Post Sale

Order Number	Description
SRVM005SPYDIG1013	SpyGlass Digital EverCare POST 1 Year
SRVM005SPYDIG1014	SpyGlass Digital EverCare POST 2 Years
SRVM005SPYDIG1015	SpyGlass Digital EverCare POST 3 Years
SRVM005SPYDIG1016	SpyGlass Digital EverCare POST 4 Years

SpyGlass DS Digital EssentialCare POS

Order Number	Description
SRVM005SPYDIG1001	SpyGlass Digital EssentialCare POS 1 Year
SRVM005SPYDIG1002	SpyGlass Digital EssentialCare POS 2 Years
SRVM005SPYDIG1003	SpyGlass Digital EssentialCare POS 3 Years
SRVM005SPYDIG1004	SpyGlass Digital EssentialCare POS 4 Years

SpyGlass DS Digital EssentialCare Post Sale

Order Number	Description
SRVM005SPYDIG1005	SpyGlass Digital EssentialCare POST 1 Year
SRVM005SPYDIG1006	SpyGlass Digital EssentialCare POST 2 Years
SRVM005SPYDIG1007	SpyGlass Digital EssentialCare POST 3 Years
SRVM005SPYDIG1008	SpyGlass Digital EssentialCare POST 4 Years

SpyGlass DS Digital EverCare Select POS

Order Number	Description
SRVM005SPYDSEVCS01	SpyGlass DS EverCare Select 3 Years POS
SRVM005SPYDSEVCS03	SpyGlass DS EverCare Select 4 Years POS

SpyGlass DS Digital EverCare Select Post Sale

Order Number	Description
SRVM005SPYDSEVCS02	SpyGlass DS EverCare Select 3 Years POST
SRVM005SPYDSEVCS04	SpyGlass DS EverCare Select 4 Years POST

Scivita Medical™ Capital Equipment



Order Number	Description	Units
HDVS-S300D	Full HD Visualization Endoscopic Image Processor / Tablet	box 1
HDVS-S200D	Endoscopic Image Processing Device / Controller	box 1

DIRECT VISUALISATION SYSTEM (CONTINUED)

Scivita Medical™ Disposable Equipment



Order Number	Description	Outer Diameter (mm)	Working Channel (mm)	Units
SCV-P-01G	Single-use Video Cholangioscope	3.70	2.0	box 1
SCV-P-02G	Single-use Video Cholangioscope	2.95	1.2	box 1
SBV-1A-B	Single-use Broncho Videoscope	2.8	1.2	box 1
SBV-1B-B	Single-use Broncho Videoscope	4.2	2.2	box 1
SBV-1C-B	Single-use Broncho Videoscope	5.6	2.8	box 1

Scivita Medical™ Service Codes

Scivita Medical EverCare for S200

Order Number	Description
SRVM005HDVS200EVC3	HDVS-S200D 1YR EverCare POST
SRVM005HDVS200EVC4	HDVS-S200D 1YR EverCare POS
SRVM005HDVS200EVC5	HDVS-S200D 2YR EverCare POST
SRVM005HDVS200EVC6	HDVS-S200D 2YR EverCare POS
SRVM005HDVS200EVC1	HDVS-S200D 3YR EverCare POST
SRVM005HDVS200EVC2	HDVS-S200D 3YR EverCare POS
SRVM005HDVS200EVC7	HDVS-S200D 4YR EverCare POST
SRVM005HDVS200EVC8	HDVS-S200D 4YR EverCare POS

Scivita Medical EverCare for S300

Order Number	Description
SRVM005HDVS300EVC3	HDVS-S300D 1YR EverCare POST
SRVM005HDVS300EVC4	HDVS-S300D 1YR EverCare POS
SRVM005HDVS300EVC5	HDVS-S300D 2YR EverCare POST
SRVM005HDVS300EVC6	HDVS-S300D 2YR EverCare POS
SRVM005HDVS300EVC1	HDVS-S300D 3YR EverCare POST
SRVM005HDVS300EVC2	HDVS-S300D 3YR EverCare POS
SRVM005HDVS300EVC7	HDVS-S300D 4YR EverCare POST
SRVM005HDVS300EVC8	HDVS-S300D 4YR EverCare POS

Scivita Medical EssentialCare for S200

Order Number	Description
SRVM005HDVSS200EC1	EssentialCare S200 -1 year POS
SRVM005HDVSS200EC2	EssentialCare S200 -1 year POST
SRVM005HDVSS200EC3	EssentialCare S200 -2 years POS
SRVM005HDVSS200EC4	EssentialCare S200 -2 years POST
SRVM005HDVSS200EC5	EssentialCare S200 -3 years POS
SRVM005HDVSS200EC6	EssentialCare S200 -3 years POST
SRVM005HDVSS200EC7	EssentialCare S200 -4 years POS
SRVM005HDVSS200EC8	EssentialCare S200 -4 years POST

Scivita Medical EssentialCare for S300

Order Number	Description
SRVM005HDVSS300EC1	EssentialCare for S300 -1 year POS
SRVM005HDVSS300EC2	EssentialCare for S300 -1 year POST
SRVM005HDVSS300EC3	EssentialCare for S300 -2 years POS
SRVM005HDVSS300EC4	EssentialCare for S300 -2 years POST
SRVM005HDVSS300EC5	EssentialCare for S300 -3 years POS
SRVM005HDVSS300EC6	EssentialCare for S300 -3 years POST
SRVM005HDVSS300EC7	EssentialCare for S300 -4 years POS
SRVM005HDVSS300EC8	EssentialCare for S300 -4 years POST

Scivita Medical (Electrical) Safety Testing (EST / STK)

Order Number	Description
SRVM005HDVSS200STK	HDVS-S200 Safety Testing (STK)
SRVM005HDVSS200EST	HDVS-S200 Electrical Safety Test (EST)
SRVM005HDVSS300STK	HDVS-S300D Safety Testing (STK)
SRVM005HDVSS300EST	HDVS-S300D Electrical Safety Test (EST)
SRVM005HDVSFSEFFEE	Scivita FSE Flat Fee

DIRECT VISUALISATION SYSTEM (CONTINUED)

SpyGlass™ Discover Digital Catheter



Order Number	Description	Units
M00546780	SpyGlass Discover Digital Catheter	box 1
M00542950	SpyGlass Discover Controller	box 1

SpyGlass™ Discover Service Codes

SpyGlass Discover EssentialCare POS

Order Number	Description
SRVM005SPYDISEC01	SpyGlass Digital EssentialCare POS 1 Year
SRVM005SPYDISEC02	SpyGlass Digital EssentialCare POS 2 Years
SRVM005SPYDISEC03	SpyGlass Digital EssentialCare POS 3 Years
SRVM005SPYDISEC04	SpyGlass Digital EssentialCare POS 4 Years

SpyGlass Discover EssentialCare Post Sale

Order Number	Description
SRVM005SPYDISEC05	SpyGlass Digital EssentialCare POST 1 Year
SRVM005SPYDISEC06	SpyGlass Digital EssentialCare POST 2 Years
SRVM005SPYDISEC07	SpyGlass Digital EssentialCare POST 3 Years
SRVM005SPYDISEC08	SpyGlass Digital EssentialCare POST 4 Years

Autolith™ Touch Biliary EHL System

Order Number	Description	Units
M00546680*	Autolith Touch EHL Generator	box 1
M00546620	1.9 Fr, 375 cm Biliary EHL Probe	box 1
M00546750	Autolith Touch Extender Cable	box 1
M00546760	Autolith Touch Foot Pedal	box 1

*Contains Extender Cable and Foot Pedal

Autolith™ Service Codes

Autolith EssentialCare POS





Order Number	Description
SRVM005ATLT1001	Autolith EssentialCare POS 1 Year
SRVM005ATLT1005	Autolith EssentialCare POST 1 Year
SRVM005ATLT1002	Autolith EssentialCare POS 2 Years
SRVM005ATLT1006	Autolith EssentialCare POST 2 Years
SRVM005ATLT1003	Autolith EssentialCare POS 3 Years
SRVM005ATLT1007	Autolith EssentialCare POST 3 Years
SRVM005ATLT1004	Autolith EssentialCare POS 4 Years
SRVM005ATLT1008	Autolith EssentialCare POST 4 Years

Autolith EssentialCare PLUS POS

Order Number	Description
SRVM005ATLTECP1001	Autolith EssentialCare PLUS POS 3 Years
SRVM005ATLTECP1002	Autolith EssentialCare PLUS POS 4 Years
SRVM005ATLTECP1003	Autolith EssentialCare PLUS POS 5 Years

Autolith EssentialCare PLUS Post Sale

Order Number	Description
SRVM005ATLTECP1004	Autolith EssentialCare PLUS POST 3 Years
SRVM005ATLTECP1005	Autolith EssentialCare PLUS POST 4 Years
SRVM005ATLTECP1006	Autolith EssentialCare PLUS POST 5 Years



Biliary	RX Biliary™ System™	Biliary Lumen Patency	Endoscopic Ultrasound (EUS)	Guidewires	Biopsy	Dilation	Imaging Systems	Pulmonary	Closure	Haemostasis	Tissue Resection	Esophageal and Enteral Metal Stents	Infection Prevention	Endobariatrics
---------	---------------------	-----------------------	-----------------------------	------------	--------	----------	-----------------	-----------	---------	-------------	------------------	-------------------------------------	----------------------	----------------

BRONCHO CONSUMABLES

Scivita Medical™ Capital Equipment



Order Number	Description	Units
HDVS-S300D	Full HD Visualization Endoscopic Image Processor / Tablet	box 1
HDVS-S200D	Endoscopic Image Processing Device / Controller	box 1

Scivita Medical™ Disposable Equipment



Order Number	Description	Size (mm)	Units
SBV-1A-B	Single-use Broncho Videoscope	2.7	box 1
SBV-1B-B	Single-use Broncho Videoscope	4	box 1
SBV-1C-B	Single-use Broncho Videoscope	5.6	box 1

Scivita Medical™ Service Codes

Scivita Medical EverCare for S200

Order Number	Description
SRVM005HDVS200EVC3	HDVS-S200D 1YR EverCare POST
SRVM005HDVS200EVC4	HDVS-S200D 1YR EverCare POS
SRVM005HDVS200EVC5	HDVS-S200D 2YR EverCare POST
SRVM005HDVS200EVC6	HDVS-S200D 2YR EverCare POS
SRVM005HDVS200EVC1	HDVS-S200D 3YR EverCare POST
SRVM005HDVS200EVC2	HDVS-S200D 3YR EverCare POS
SRVM005HDVS200EVC7	HDVS-S200D 4YR EverCare POST
SRVM005HDVS200EVC8	HDVS-S200D 4YR EverCare POS

Scivita Medical EverCare for S300

Order Number	Description
SRVM005HDVS300EVC3	HDVS-S300D 1YR EverCare POST
SRVM005HDVS300EVC4	HDVS-S300D 1YR EverCare POS
SRVM005HDVS300EVC5	HDVS-S300D 2YR EverCare POST
SRVM005HDVS300EVC6	HDVS-S300D 2YR EverCare POS
SRVM005HDVS300EVC1	HDVS-S300D 3YR EverCare POST
SRVM005HDVS300EVC2	HDVS-S300D 3YR EverCare POS
SRVM005HDVS300EVC7	HDVS-S300D 4YR EverCare POST
SRVM005HDVS300EVC8	HDVS-S300D 4YR EverCare POS

Scivita Medical EssentialCare for S200

Order Number	Description
SRVM005HDVSS200EC1	EssentialCare S200 -1 year POS
SRVM005HDVSS200EC2	EssentialCare S200 -1 year POST
SRVM005HDVSS200EC3	EssentialCare S200 -2 years POS
SRVM005HDVSS200EC4	EssentialCare S200 -2 years POST
SRVM005HDVSS200EC5	EssentialCare S200 -3 years POS
SRVM005HDVSS200EC6	EssentialCare S200 -3 years POST
SRVM005HDVSS200EC7	EssentialCare S200 -4 years POS
SRVM005HDVSS200EC8	EssentialCare S200 -4 years POST

Scivita Medical EssentialCare for S300

Order Number	Description
SRVM005HDVSS300EC1	EssentialCare for S300 -1 year POS
SRVM005HDVSS300EC2	EssentialCare for S300 -1 year POST
SRVM005HDVSS300EC3	EssentialCare for S300 -2 years POS
SRVM005HDVSS300EC4	EssentialCare for S300 -2 years POST
SRVM005HDVSS300EC5	EssentialCare for S300 -3 years POS
SRVM005HDVSS300EC6	EssentialCare for S300 -3 years POST
SRVM005HDVSS300EC7	EssentialCare for S300 -4 years POS
SRVM005HDVSS300EC8	EssentialCare for S300 -4 years POST

Scivita Medical (Electrical) Safety Testing (EST / STK)

Order Number	Description
SRVM005HDVSS200STK	HDVS-S200 Safety Testing (STK)
SRVM005HDVSS200EST	HDVS-S200 Electrical Safety Test (EST)
SRVM005HDVSS300STK	HDVS-S300D Safety Testing (STK)
SRVM005HDVSS300EST	HDVS-S300D Electrical Safety Test (EST)
SRVM005HDVSFSEFFEE	Scivita FSE Flat Fee

ENDOSCOPIC BRONCHIAL ULTRASOUND



Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needles Olympus Scope Compatible




Order Number	Description	Minimum Working Channel (mm)	Sheath OD (mm)	Units
M00558220	Expect Pulmonary 22ga – Olympus Scope Compatible	2.0	1.6	box 1
M00558221	Expect Pulmonary 22ga – Olympus Scope Compatible	2.0	1.6	box 5
M00558250	Expect Pulmonary 25ga – Olympus Scope Compatible	2.0	1.4	box 1
M00558251	Expect Pulmonary 25ga – Olympus Scope Compatible	2.0	1.4	box 5
M00558731	Expect Pulmonary Needle Adaptor – Olympus Scope Compatible			box 10

Acquire™ Pulmonary Endobronchial Ultrasound Fine Biopsy Needles Olympus Scope Compatible



Order Number	Description	Needle Size (gauge) / (mm)	Minimum Working Channel (mm)	Sheath Diameter (mm)	Units
M00552350	22ga Acquire Pulmonary – Olympus scope compatible	22/0.72	2.0	1.6	box 1
M00552351	22ga Acquire Pulmonary – Olympus scope compatible	22/0.72	2.0	1.6	box 5
M00552400	25ga Acquire Pulmonary – Olympus scope compatible	25/0.52	2.0	1.4	box 1
M00552401	25ga Acquire Pulmonary – Olympus scope compatible	25/0.52	2.0	1.4	box 5

AIRWAY STENTS



Ultraflex™ Tracheobronchial Covered Stent System Distal Release



Order Number	Description	Stent OD (mm)	Expanded Stent Length (mm)	Covered Length (mm)
M00576600	Ultraflex Tracheal NG Stent System	14	80	65
M00576610	Ultraflex Tracheal NG Stent System	16	40	25
M00576620	Ultraflex Tracheal NG Stent System	16	60	45
M00576630	Ultraflex Tracheal NG Stent System	16	80	65
M00576640	Ultraflex Tracheal NG Stent System	18	40	25
M00576650	Ultraflex Tracheal NG Stent System	18	60	45
M00576660	Ultraflex Tracheal NG Stent System	18	80	65
M00576670	Ultraflex Tracheal NG Stent System	20	40	25
M00576680	Ultraflex Tracheal NG Stent System	20	60	45
M00576690	Ultraflex Tracheal NG Stent System	20	80	65
M00576520	Ultraflex Tracheobronchial Stent System	8	40	25
M00576530	Ultraflex Tracheobronchial Stent System	10	30	15
M00576540	Ultraflex Tracheobronchial Stent System	10	40	25
M00576550	Ultraflex Tracheobronchial Stent System	12	30	15
M00576560	Ultraflex Tracheobronchial Stent System	12	40	25
M00576570	Ultraflex Tracheobronchial Stent System	14	30	15
M00576580	Ultraflex Tracheobronchial Stent System	14	40	25
M00576590	Ultraflex Tracheobronchial Stent System	14	60	45

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Jagwire™ Pulmonary Guidewire, M00515171 / .038in (0.97mm) Amplatz Super Stiff™ Guidewire, M00550090

Additional Specifications for Ultraflex Tracheobronchial Stents: Delivery System Working Length is 95cm, Maximum Catheter Shaft OD is 16F (5.33mm), Compressed Stent on catheter yields Delivery System OD range 15-25F (5.00-8.33mm)



Ultraflex™ Tracheobronchial Uncovered Stent System Proximal Release




Order Number	Description	Stent OD (mm)	Expanded Stent Length (mm)
M00576480	Ultraflex Tracheobronchial Stent System	14	40
M00576490	Ultraflex Tracheobronchial Stent System	14	60
M00576410	Ultraflex Tracheobronchial Stent System	8	20
M00576420	Ultraflex Tracheobronchial Stent System	8	40
M00576430	Ultraflex Tracheobronchial Stent System	10	20
M00576440	Ultraflex Tracheobronchial Stent System	10	40
M00576450	Ultraflex Tracheobronchial Stent System	12	20
M00576460	Ultraflex Tracheobronchial Stent System	12	40
M00576470	Ultraflex Tracheobronchial Stent System	14	20

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Jagwire™ Pulmonary Guidewire, M00515171 / .038in (0.97mm) Amplatz Super Stiff™ Guidewire, M00550090

Additional Specifications for Ultraflex Tracheobronchial Stents: Delivery System Working Length is 95cm, Maximum Catheter Shaft OD is 16F (5.33mm), Compressed Stent on catheter yields Delivery System OD range 15-18F (5.00-6.00mm)

AIRWAY STENTS (CONTINUED)



Ultraflex™ Tracheobronchial Uncovered Stent System Distal Release




Order Number	Description	Stent OD (mm)	Expanded Stent Length (mm)
M00576330	Ultraflex Tracheobronchial Stent System	16	40
M00576340	Ultraflex Tracheobronchial Stent System	16	60
M00576350	Ultraflex Tracheobronchial Stent System	18	40
M00576360	Ultraflex Tracheobronchial Stent System	18	60
M00576370	Ultraflex Tracheobronchial Stent System	20	40
M00576380	Ultraflex Tracheobronchial Stent System	20	60
M00576300	Ultraflex Tracheobronchial Stent System	10	30
M00576310	Ultraflex Tracheobronchial Stent System	12	30
M00576320	Ultraflex Tracheobronchial Stent System	14	30

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Jagwire™ Pulmonary Guidewire, M00515171 / .038in (0.97mm) Amplatz Super Stiff™ Guidewire, M00550090

Additional Specifications for Ultraflex Tracheobronchial Stents: Delivery System Working Length is 95cm, Maximum Catheter Shaft OD is 16F (5.33mm), Compressed Stent on catheter yields Delivery System OD range 15-20F (5.00-6.67mm)

BIOPSY FORCEPS



Radial Jaw™ 4 Pulmonary Biopsy Forceps



Order Number	Description	Jaw OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00515181	Radial Jaw 4 Pulmonary Standard Capacity	1.8	100	2.0	Blue	box 5
M00515182	Radial Jaw 4 Pulmonary Standard Capacity	1.8	100	2.0	Blue	box 20
M00515191	Radial Jaw 4 Pulmonary Standard Capacity w/Needle	1.8	100	2.0	Blue	box 5
M00515192	Radial Jaw 4 Pulmonary Standard Capacity w/Needle	1.8	100	2.0	Blue	box 20
M00515201	Radial Jaw 4 Pulmonary Large Capacity	2.2	100	2.8	Blue	box 5
M00515202	Radial Jaw 4 Pulmonary Large Capacity	2.2	100	2.8	Blue	box 20

MINI-FORCEPS

CoreDx™ Pulmonary Mini-Forceps



Order Number	Description	Jaw OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00515220	Pulmonary Mini-Forceps	1.0	130	1.2	box 1


CYTOLOGY BRUSHES

Cellebrity™ Cytology Brushes – Bronchoscope



Order Number	Description	Required Working Channel (mm)	Bristle OD (mm)	Sheath Length (cm)	Units
M00516001	Cellebrity Cytology Brush	2.0	1.0	140	box 10
M00516011	Cellebrity Cytology Brush	2.0	1.5	140	box 10
M00516071	Cellebrity Cytology Brush	2.0	1.9	100	box 10
M00516151	Cellebrity Cytology Brush	2.0	1.9	150	box 10

AIRWAY RETRIEVAL BASKETS



Zero Tip™ Airway Retrieval Baskets



Order Number	Description	Sheath OD (mm)	Working Length (cm)	Basket Opening (mm)
M00513200	Zero Tip Airway Retrieval Basket	0.8	120	12
M00513210	Zero Tip Airway Retrieval Basket	1.0	120	16

Packaged one per box

DILATATION BALLOONS



CRE™ Pulmonary Balloon Dilatation Catheter



Order Number	Description	Balloon Length (cm)	Initial Balloon OD (mm)	Intermediate Balloon OD (mm)	Maximum Balloon OD (mm)	Working Catheter Length (cm)
M00550300	CRE Pulmonary Balloon Dilatation Catheter	5.5	12	13.5	15	110
M00550310	CRE Pulmonary Balloon Dilatation Catheter	5.5	15	16.5	18	110
M00550320	CRE Pulmonary Balloon Dilatation Catheter	5.5	18	19.0	20	110
M00550330	CRE Pulmonary Balloon Dilatation Catheter	3.0	8	9.0	10	110
M00550340	CRE Pulmonary Balloon Dilatation Catheter	3.0	10	11.0	12	110
M00550350	CRE Pulmonary Balloon Dilatation Catheter	3.0	12	13.5	15	110

Packaged one per box


GUIDEWIRES

Jagwire™ Pulmonary Guidewire



Order Number	Description	OD (in) / (mm)	Length (cm)	Units
M00515171	Jagwire Pulmonary Guidewire	.035 / 0.89	180	box 2


INFLATION SYSTEM



Alliance™ II Inflation System



Order Number	Description	Units
M00550620	Alliance II Integrated Inflation/Lithotripsy Device	box 1
M00550601	Alliance II Single-Use Syringe/Gauge Assembly	box 5



Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics

CLIPPING DEVICES

Mantis™ Clip



Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Units
M00521420	Mantis Clip	235	2.8	box 1
M00521421	Mantis Clip	235	2.8	box 10
M00521422	Mantis Clip	235	2.8	box 20
M00521423	Mantis Clip	235	2.8	box 40

ENDOSCOPIC SUTURING

OverStitch™ Systems



Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Units
M00505400	OverStitch Endoscopic Suturing System	*	-	box 1
M00505441	OverStitch 2-0 Polypropylene Suture	*	2.8	box 12
M00505540	OverStitch Suture Cinch	160	3.2	box 1
M00505470	OverStitch Tissue Helix	*	2.8	box 1
M00505480	OverStitch NXT Tissue Helix Pro	165	2.8	box 1
M00505420	OverStitch NXT Endoscopic Suture System	*	-	box 1

* Suitable for endoscopes with max 110 cm working length.

OverStitch™ Accessories



Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Units
M00505560	OverTube Endoscopic Access System	-	-	box 1
ES26165-3	Ensizor Endoscopic Scissors Gastric Length	165	2.8	box 3

X-Tack™



Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Units
M00505500	X-Tack Endoscopic HeliX Tacking System (Gastric)	> 230	2.8	box 1
M00505510	X-Tack Endoscopic HeliX Tacking System (Colonic)	> 230	2.8	box 1

X-Tack™ Accessories



Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Units
M00505550	Suture Cinch – Long	≥ 204	2.8	box 1
M00505540	OverStitch Suture Cinch	≥ 141	2.8	box 1
ES26165-3	Ensizor Endoscopic Scissors Gastric Length	165	2.8	box 3
ES26235-3	Ensizor Endoscopic Scissors Colonic Length	235	2.8	box 3

ENDOLUMINAL VACUUM THERAPY

Eso-SPONGE®



Order Number	Description	Inner Diameter (mm)	Outer Diameter (mm)	Units
5526540	Eso-SPONGE with 15 mm Overtube	15	19	box 1
5526550	Eso-SPONGE with 13 mm Overtube	13	17	box 1

Endo-SPONGE®



Order Number	Description	Units
5526510	Endo-SPONGE	box 10
5526520	Endo-SPONGE	box 5
5526530	Endo-SPONGE	box 1

EVT Vacuum Sources




Order Number	Description	Units
5526661	Redyrob Trans Plus	box 10
5526650	Endo-SPONGE® Pump*	box 1
5526651	Endo-SPONGE® Pump Carrying Bag	box 5
5526653	Disposable Cannister 300 ml with tubing and Y-connector	box 16
19116	MV1 Mobile Programmable Low-Vacuum-System*	box 1
18032	Collection Bottle MV1	box 15
18022	Bacterial Filter MV1	box 20

* The "Do Not Reuse" symbol does not apply to this product.



Biliary
RX Biliary™ System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics

CLIPPING DEVICES




Resolution™ Clips



Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Units
M00522600	Resolution Clip	155	2.8	11	box 1
M00522601	Resolution Clip	155	2.8	11	box 10
M00522602	Resolution Clip	155	2.8	11	box 20
M00522610	Resolution Clip	235	2.8	11	box 1
M00522611	Resolution Clip	235	2.8	11	box 10
M00522612	Resolution Clip	235	2.8	11	box 20

* MR Conditional. Non-clinical testing has demonstrated the Resolution 360™ Clip is Magnetic Resonance (MR) Conditional according to ASTM F2503. Refer to the directions for use for MR conditions of use.



Resolution 360™ Clips



Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Units
M00521230	Resolution 360 Clip	235	2.8	11	box 1
M00521231	Resolution 360 Clip	235	2.8	11	box 10
M00521232†	Resolution 360 Clip	235	2.8	11	box 20
M00521233†	Resolution 360 Clip	235	2.8	11	box 40
M00521240	Resolution 360 Clip	155	2.8	11	box 1
M00521241	Resolution 360 Clip	155	2.8	11	box 10
M00521242	Resolution 360 Clip	155	2.8	11	box 20
M00521243	Resolution 360 Clip	155	2.8	11	box 40

* MR Conditional. Non-clinical testing has demonstrated the Resolution 360™ Clip is Magnetic Resonance (MR) Conditional according to ASTM F2503. Refer to the directions for use for MR conditions of use.


† Please check product availability in your country with your sales representative

CLIPPING DEVICES (CONTINUED)

Resolution 360™ ULTRA Clips



Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Units
M00521400	Resolution 360 ULTRA Clip	235	2.8	17	box 1
M00521401	Resolution 360 ULTRA Clip	235	2.8	17	box 10
M00521402	Resolution 360 ULTRA Clip	235	2.8	17	box 20
M00521403	Resolution 360 ULTRA Clip	235	2.8	17	box 40




Speedband Superview Super 7™ Multiple Band Ligators



Order Number	Description	Minimum Working Channel (mm)	Scope OD (mm)	Bands / Ligating Unit	Units
M00542250	Super 7 Band Ligator	2.8	8.6-11.5	7	box 1
M00542251	Super 7 Band Ligator	2.8	8.6-11.5	7	box 2
M00542253	Super 7 Band Ligator	2.8	8.6-11.5	7	box 4

Irrigation catheter supplied in each kit

BIPOLAR ELECTROHAEMOSTASIS CATHETERS




Injection Gold Probe™ Electrohaemostasis/Injection Catheters



Order Number	Description	OD (F) / (mm)	Gauge	NEEDLE OD (mm)	ID (mm)	Minimum Working Channel (mm)	Working Length (cm)	Units
M00560150	Injection Gold Probe Catheter	7 / 2.33	25	0.51	0.24	2.8	210	box 1
M00560160	Injection Gold Probe Catheter	10 / 3.33	25	0.51	0.24	3.7	210	box 1

Type Generator/Connector: Please refer to Electrosurgical Bipolar Cable Adaptors

BIPOLAR ELECTROHAEMOSTASIS CATHETERS (CONTINUED)




Gold Probe™ Electrohaemostasis Catheters



Order Number	Description	OD (F) / (mm)	Minimum Working Channel (mm)	Working Length (cm)	Units
M00560070	Gold Probe Catheter	7 / 2.33	2.8	300	box 1
M00560100	Gold Probe Catheter	10 / 3.33	3.7	300	box 1
M00560220	Gold Probe Catheter	7 / 2.33	2.8	350	box 1

Type Generator/Connector: Please refer to Electrosurgical Bipolar Cable Adaptors

NEEDLES




Interject™ Sclerotherapy Needles*



Order Number	Description	Sheath Design	Sheath OD (F) / (mm)	Gauge	NEEDLE OD (mm)	ID (mm)	Maximum Needle Extension Length (mm)	Working Length (cm)	Units
M00518151	Interject Sclerotherapy Needle	Contrast	1.8	23	0.64	0.32	4	200	box 5
M00518161	Interject Sclerotherapy Needle	Contrast	1.8	25	0.51	0.24	4	200	box 5
M00518251	Interject Sclerotherapy Needle	Contrast	1.8	23	0.64	0.32	6	200	box 5
M00518261	Interject Sclerotherapy Needle	Contrast	1.8	25	0.51	0.24	6	200	box 5
M00518351	Interject Sclerotherapy Needle	Contrast	2.3	23	0.64	0.32	4	240	box 5
M00518361	Interject Sclerotherapy Needle	Contrast	2.3	25	0.51	0.24	4	240	box 5

* With Star Catheter Technology

NEEDLES (CONTINUED)




Interject™ Clear Sclerotherapy Needles*



Order Number	Description	Sheath Design	Sheath OD (F) / (mm)	Gauge	NEEDLE OD (mm)	ID (mm)	Maximum Needle Extension Length (mm)	Working Length (cm)	Units
M00518111	Interject Sclerotherapy Needle	Clear	1.8	25	0.51	0.24	4	200	box 5
M00518301	Interject Sclerotherapy Needle	Clear	2.3	23	0.64	0.32	4	240	box 5
M00518311	Interject Sclerotherapy Needle	Clear	2.3	25	0.51	0.24	4	240	box 5

* With Star Catheter Technology

ELECTROSURGICAL BIPOLAR CABLE ADAPTOR




Order Number	Description	Fix Pin Spacing (mm)	Generator Compatibility*	Units
M00505200	Bipolar Cable Adaptor	22	ERBE, EMC Dolley, Olympus (ESG100, PSD60)**	box 1
M00505210	Bipolar Cable Adaptor	28	Valleylab®, Lamidey (Duo), Bovie (Aaron 3250™), Olympus (UES10, UES30, Surgemaster), Bowa, Soring	box 1
M00505220	Bipolar Cable Adaptor	28 with Third Pin	ConMed® (System 5000), Bovie (Aaron 2250™), Ethicon®	box 1

Bipolar Cable Adaptors are used to connect Gold Probe™ and Injection Gold Probe™ with the above mentioned generators.

* For generators and their associated cable adaptors not listed above, please contact your Boston Scientific Representative.

** Some variations in faceplates may exist within a generator model. We recommend that you confirm pin spacing prior to ordering.


Tissue Resection



Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics

Tissue Resection

EMR DEVICES FOR UPPER GI



Captivator™ EMR



Order Number	Description	Units
M00561600	Captivator EMR for Standard Gastrosopes	box 1
M00561601	Captivator EMR for Standard Gastrosopes	box 5
M00561610	Captivator EMR for Large Gastrosopes	box 1
M00561611	Captivator EMR for Large Gastrosopes	box 5


ELECTROSURGICAL KNIFE

NEW ORISE™ ProKnife



Order Number	Description	Catheter OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00519351	1.5mm ORISE ProKnife	2.21	230	2.8	box 1
M00519361	2.0mm ORISE ProKnife	2.21	230	2.8	box 1
M00519371	3.0mm ORISE ProKnife	2.21	230	2.8	box 1

SNARES



Rotatable Single-Use Polypectomy Snares*




Order Number	Description	Stiffness	Loop Width (mm)	Sheath OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00561821	Small Oval	Medium Stiff	13	2.4	240	2.8	box 10
M00561831	Medium Oval	Medium Stiff	20	2.4	240	2.8	box 10

* With Star Catheter Technology.


Tissue Resection

SNARES (CONTINUED)



Sensation™ Single-Use Short Throw Snares


Order Number	Description	Stiffness	Loop Width (mm)	Sheath OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562691	Small Oval	Flexible	13	2.4	240	2.8	box 5
M00562692	Small Oval	Flexible	13	2.4	240	2.8	box 20
M00562693	Small Oval	Flexible	13	2.4	240	2.8	box 40
M00560321	Small Oval	Medium Stiff	13	2.4	240	2.8	box 5
M00562671	Medium Oval	Flexible	27	2.4	240	2.8	box 5
M00562672	Medium Oval	Flexible	27	2.4	240	2.8	box 20
M00562673	Medium Oval	Flexible	27	2.4	240	2.8	box 40
M00560311	Medium Oval	Medium Stiff	27	2.4	240	2.8	box 5
M00562651	Large Oval	Medium Stiff	30	2.4	240	2.8	box 5
M00562652	Large Oval	Medium Stiff	30	2.4	240	2.8	box 20
M00560331	Medium Crescent	Medium Stiff	27	2.4	240	2.8	box 5



Captiflex™ Single-Use Snares


Order Number	Description	Stiffness	Loop Width (mm)	Sheath OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562471	Extra Small Oval	Flexible	11	2.4	240	2.8	box 10
M00562421	Small Oval	Flexible	13	2.4	240	2.8	box 10
M00562422	Small Oval	Flexible	13	2.4	240	2.8	box 40
M00562401	Medium Oval	Flexible	27	2.4	240	2.8	box 10
M00562402	Medium Oval	Flexible	27	2.4	240	2.8	box 40

SNARES (CONTINUED)



Captivator™ Single-Use Snares

Order Number	Description	Stiffness	Loop Width (mm)	Sheath OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562301	Small Oval	Stiff	13	2.4	240	2.8	box 10
M00562321	Medium Oval	Stiff	27	2.4	240	2.8	box 10
M00562391	Large Oval (thin wire)	Flexible	30	2.4	240	2.8	box 10
M00561311	Large Oval	Medium Stiff	30	2.4	240	2.8	box 10
M00562451	Small Hexagonal	Stiff	13	2.4	240	2.8	box 10
M00562341	Medium Hexagonal	Stiff	27	2.4	240	2.8	box 10
M00562371	Medium Crescent	Stiff	27	2.4	240	2.8	box 10




Captivator™ II Single-Use Snares

Order Number	Description	Stiffness	Loop Width (mm)	Sheath OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00561221	10 mm Rounded	Stiff	10	2.4	240	2.8	box 10
M00561222	10 mm Rounded	Stiff	10	2.4	240	2.8	box 20
M00561223	10 mm Rounded	Stiff	10	2.4	240	2.8	box 40
M00561231	15 mm Rounded	Stiff	15	2.4	240	2.8	box 10
M00561232	15 mm Rounded	Stiff	15	2.4	240	2.8	box 20
M00561233	15 mm Rounded	Stiff	15	2.4	240	2.8	box 40
M00561241	20 mm Rounded	Stiff	20	2.4	240	2.8	box 10
M00561242	20 mm Rounded	Stiff	20	2.4	240	2.8	box 20
M00561243	20 mm Rounded	Stiff	20	2.4	240	2.8	box 40
M00561191	25 mm Rounded	Stiff	25	2.4	240	2.8	box 10
M00561192	25 mm Rounded	Stiff	25	2.4	240	2.8	box 20
M00561193	25 mm Rounded	Stiff	25	2.4	240	2.8	box 40
M00561291	33 mm Rounded	Stiff	33	2.4	240	2.8	box 10

Tissue Resection

SNARES (CONTINUED)




Captivator™ COLD Single-Use Snare



Order Number	Description	Loop Width (mm)	Sheath OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00561101	10mm Rounded-Stiff	10	2.4	240	2.8	box 10
M00561102	10mm Rounded-Stiff	10	2.4	240	2.8	box 20
M00561103	10mm Rounded-Stiff	10	2.4	240	2.8	box 40

HOT BIOPSY FORCEPS



Radial Jaw™ 4 Hot Biopsy Forceps



Order Number	Description	Jaw OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00515031	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	box 5
M00515032	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	box 20
M00515033	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	box 40

ACCESSORIES

Active Cords

Order Number	Description	Units
M00561250	Active Cord, Fits Olympus® Brand Generator	box 1
M00561270	Active Cord, Fits Endostat™, Endostat™ II, Endostat™ III, Bovie®, Valleylab®, Aspen Labs and ERBE® Brand Generators	box 1



Tissue Resection

RETRIEVAL DEVICES

RescueNet™



Order Number	Description	Net Size (mm)	Sheath OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
DGN-538-5	RescueNet Retrieval Device	30	2.5	230	2.8	box 5

Rescue™ Forceps



Order Number	Description	Jaw Opening (mm)	Sheath OD (mm)	Minimum Working Channel (cm)	Working Length (cm)	Units
M00527651	Rescue Rat Tooth Grasping Forceps	7	2.4	2.8	230	box 5
M00527751	Rescue Rat Tooth Grasping Forceps – Pediatric	4.5	1.8	2.0	180	box 5
M00527551	Rescue Alligator Grasping Forceps Long	9	2.4	2.8	230	box 5
M00527451	Rescue Alligator Grasping Forceps	8	2.4	2.8	230	box 5
M00538351	Rescue Rat Tooth / Alligator Grasping Forceps	8	2.4	2.8	230	box 5

POLYP TRAPS


TrapEase™ Polyp Traps



Order Number	Description	Units
STE-298-25	TrapEase Quad Trap	box 25
STE-297-50	TrapEase Single Trap	box 50




Esophageal and Enteral Metal Stents



Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics

Esophageal and Enteral Metal Stents

ESOPHAGEAL STENTS



WallFlex™ Fully Covered RMV Esophageal Longer Loop Stent System

Order Number	Description	Stent OD (mm)	Proximal / Distal Flares OD (mm)	Stent Length (mm)	Catheter Diameter (F)	Working Length (cm)
M00516270	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	18	25 / 23	103	18.5	120
M00516280	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	18	25 / 23	123	18.5	120
M00516290	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	18	25 / 23	153	18.5	120
M00516300	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	23	28 / 28	105	18.5	120
M00516310	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	23	28 / 28	125	18.5	120
M00516320	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	23	28 / 28	155	18.5	120

Packaged one per box

Limited availability in all European countries (except UK and France). Please contact your local Territory Manager.

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621

WallFlex™ Fully Covered RMV Esophageal Stent

Order Number	Description	Stent OD (mm)	Stent Length (mm)	Working Length (cm)
M00516210	WallFlex Fully Covered RMV Esophageal Stent	18	103	78
M00516220	WallFlex Fully Covered RMV Esophageal Stent	18	123	78
M00516230	WallFlex Fully Covered RMV Esophageal Stent	18	153	78
M00516240	WallFlex Fully Covered RMV Esophageal Stent	23	105	78
M00516250	WallFlex Fully Covered RMV Esophageal Stent	23	125	78
M00516260	WallFlex Fully Covered RMV Esophageal Stent	23	155	78

Packaged one per box

Limited availability in all European countries (except UK and France). Please contact your local Territory Manager.

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621



Esophageal and Enteral Metal Stents

ESOPHAGEAL STENTS (CONTINUED)



WallFlex™ Partially Covered Esophageal Stent System



Order Number	Description	Stent OD (mm)	Proximal / Distal Flares OD (mm)	Stent Length (mm)	Covered Length (mm)	Catheter Diameter (F)	Working Length (cm)
M00516900	WallFlex Partially Covered Esophageal Stent System	18	23 / 23	103	71	18.5	120
M00516910	WallFlex Partially Covered Esophageal Stent System	18	23 / 23	123	90	18.5	120
M00516920	WallFlex Partially Covered Esophageal Stent System	18	23 / 23	153	121	18.5	120
M00516930	WallFlex Partially Covered Esophageal Stent System	23	28 / 28	105	73	18.5	120
M00516940	WallFlex Partially Covered Esophageal Stent System	23	28 / 28	125	93	18.5	120
M00516950	WallFlex Partially Covered Esophageal Stent System	23	28 / 28	155	123	18.5	120

Packaged one per box

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621

WallFlex™ Fully Covered Esophageal Stent



Order Number	Description	Stent OD (mm)	Stent Length (mm)	Working Length (cm)
M00516700	WallFlex Fully Covered Esophageal Stent	18	103	78
M00516710	WallFlex Fully Covered Esophageal Stent	18	123	78
M00516720	WallFlex Fully Covered Esophageal Stent	18	153	78
M00516730	WallFlex Fully Covered Esophageal Stent	23	105	78
M00516740	WallFlex Fully Covered Esophageal Stent	23	125	78
M00516750	WallFlex Fully Covered Esophageal Stent	23	155	78


Packaged one per box

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621



Esophageal and Enteral Metal Stents

ESOPHAGEAL STENTS (CONTINUED)



Ultraflex™ Covered Esophageal NG Stent System

Proximal Release




Order Number	Description	Stent OD (mm)	Proximal Flare OD (mm)	Stent Length (mm)	Covered Length (mm)
M00513840	Ultraflex Esophageal NG Stent System	18	23	100	70
M00513850	Ultraflex Esophageal NG Stent System	18	23	120	90
M00513860	Ultraflex Esophageal NG Stent System	18	23	150	120
M00514240	Ultraflex Esophageal NG Stent System	23	28	100	70
M00514250	Ultraflex Esophageal NG Stent System	23	28	120	90

Distal Release

Order Number	Description	Stent OD (mm)	Proximal Flare OD (mm)	Stent Length (mm)	Covered Length (mm)
M00513730	Ultraflex Esophageal NG Stent System	18	23	100	70
M00513740	Ultraflex Esophageal NG Stent System	18	23	120	90
M00513750	Ultraflex Esophageal NG Stent System	18	23	150	120
M00514200	Ultraflex Esophageal NG Stent System	23	28	100	70
M00514210	Ultraflex Esophageal NG Stent System	23	28	120	90

Packaged one per box

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621



Ultraflex™ Uncovered Esophageal NG Stent System

Proximal Release



Order Number	Description	Stent OD (mm)	Proximal Flare OD (mm)	Stent Length (mm)
M00513800	Ultraflex Esophageal NG Stent System	18	23	70
M00513810	Ultraflex Esophageal NG Stent System	18	23	100
M00513820	Ultraflex Esophageal NG Stent System	18	23	120
M00513830	Ultraflex Esophageal NG Stent System	18	23	150
M00514220	Ultraflex Esophageal NG Stent System	23	28	70
M00514230	Ultraflex Esophageal NG Stent System	23	28	120

Distal Release

Order Number	Description	Stent OD (mm)	Proximal Flare OD (mm)	Stent Length (mm)
M00513700	Ultraflex Esophageal NG Stent System	18	23	70
M00513710	Ultraflex Esophageal NG Stent System	18	23	100
M00513720	Ultraflex Esophageal NG Stent System	18	28	150

Packaged one per box

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621



Esophageal and Enteral Metal Stents

ESOPHAGEAL STENTS (CONTINUED)

Agile™ Fully Covered Esophageal Stent System – Removable



Order Number	Description	Stent OD (mm)	Proximal / Distal Flares OD (mm)	Stent Length (mm)	Delivery System Diameter (Fr)	Working Length (cm)
M00517580	Agile Fully Covered Esophageal Stent System – Removable	14	19	60	10.5	180
M00517590	Agile Fully Covered Esophageal Stent System – Removable	14	19	100	10.5	180
M00517600	Agile Fully Covered Esophageal Stent System – Removable	14	19	120	10.5	180
M00517610	Agile Fully Covered Esophageal Stent System – Removable	14	19	150	10.5	180
M00517620	Agile Fully Covered Esophageal Stent System – Removable	18	23	60	10.5	180
M00517630	Agile Fully Covered Esophageal Stent System – Removable	18	23	100	10.5	180
M00517640	Agile Fully Covered Esophageal Stent System – Removable	18	23	120	10.5	180
M00517650	Agile Fully Covered Esophageal Stent System – Removable	18	23	150	10.5	180

Agile™ Partially Covered Esophageal Stent System



Order Number	Description	Stent OD (mm)	Proximal / Distal Flares OD (mm)	Stent Length (mm)	Delivery System Diameter (Fr)	Working Length (cm)
M00517200	Agile Partially Covered Esophageal Stent System	14	19	60	10.5	180
M00517210	Agile Partially Covered Esophageal Stent System	14	19	100	10.5	180
M00517220	Agile Partially Covered Esophageal Stent System	14	19	120	10.5	180
M00517230	Agile Partially Covered Esophageal Stent System	14	19	150	10.5	180
M00517240	Agile Partially Covered Esophageal Stent System	18	23	60	10.5	180
M00517250	Agile Partially Covered Esophageal Stent System	18	23	100	10.5	180
M00517260	Agile Partially Covered Esophageal Stent System	18	23	120	10.5	180
M00517270	Agile Partially Covered Esophageal Stent System	18	23	150	10.5	180

Agile™ OTW Fully Covered Esophageal Stent System – Removable



Order Number	Description	Stent OD (mm)	Proximal / Distal Flares OD (mm)	Stent Length (mm)	Delivery System Diameter (Fr)	Working Length (cm)
M00517790	Agile OTW Fully Covered Esophageal Stent System – Removable	23	28	60	18.5	78
M00517800	Agile OTW Fully Covered Esophageal Stent System – Removable	23	28	100	18.5	78
M00517810	Agile OTW Fully Covered Esophageal Stent System – Removable	23	28	120	18.5	78
M00517820	Agile OTW Fully Covered Esophageal Stent System – Removable	23	28	150	18.5	78



Agile™ OTW Partially Covered Esophageal Stent System



Order Number	Description	Stent OD (mm)	Proximal / Distal Flares OD (mm)	Stent Length (mm)	Delivery System Diameter (Fr)	Working Length (cm)
M00517710	Agile OTW Partially Covered Esophageal Stent System	23	28	60	18.5	78
M00517720	Agile OTW Partially Covered Esophageal Stent System	23	28	100	18.5	78
M00517730	Agile OTW Partially Covered Esophageal Stent System	23	28	120	18.5	78
M00517740	Agile OTW Partially Covered Esophageal Stent System	23	28	150	18.5	78

Esophageal and Enteral Metal Stents

COLONIC AND DUODENAL STENTS



WallFlex™ Colonic Stent System with Anchor Lock Delivery System



Order Number	Description	Diameter Flare / Body (mm)	Stent Length (cm)	Working Length (cm)	Catheter Diameter (F) / (mm)
M00565100	WallFlex Colonic Stent	27 / 22	6	230	10 / 3.33
M00565110	WallFlex Colonic Stent	27 / 22	9	230	10 / 3.33
M00565120	WallFlex Colonic Stent	27 / 22	12	230	10 / 3.33
M00565130	WallFlex Colonic Stent	27 / 22	6	135	10 / 3.33
M00565140	WallFlex Colonic Stent	27 / 22	9	135	10 / 3.33
M00565150	WallFlex Colonic Stent	27 / 22	12	135	10 / 3.33
M00565040	WallFlex Colonic Stent	30 / 25	6	230	10 / 3.33
M00565050	WallFlex Colonic Stent	30 / 25	9	230	10 / 3.33
M00565060	WallFlex Colonic Stent	30 / 25	12	230	10 / 3.33
M00565070	WallFlex Colonic Stent	30 / 25	6	135	10 / 3.33
M00565080	WallFlex Colonic Stent	30 / 25	9	135	10 / 3.33
M00565090	WallFlex Colonic Stent	30 / 25	12	135	10 / 3.33

Packaged one per box

Recommended Guidewire Systems 230cm: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 450cm: M00556161 Straight Dream Tip™, M00556171 Angled Dream Tip™ or .035 in (0.89mm) Jagwire™ Guidewire Stiff Shaft, 450cm: M00556601 Straight Tip, M00556611 Angled Tip or .035in (0.89mm) Super Stiff Wire, 500cm: H954180011

Systems 135cm: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm: M00556121 Straight Tip, M00556131 Angled Tip or .035in (0.89mm) Jagwire™ Guidewire Stiff Shaft, 260cm: M00557661 Straight Tip

WallFlex™ Colonic Soft Stent System with Anchor Lock Delivery System





Order Number	Description	Diameter Flare / Body (mm)	Stent Length (cm)	Working Length (cm)	Delivery System Diameter/ Working Channel (F/mm)	Guidewire Diameter (inch/mm)
M00555700	WallFlex Colonic Soft Stent	23 / 18	6	230	9 / 3.0	0.035 / 0.89
M00555710	WallFlex Colonic Soft Stent	23 / 18	9	230	9 / 3.0	0.035 / 0.89
M00555720	WallFlex Colonic Soft Stent	23 / 18	12	230	9 / 3.0	0.035 / 0.89
M00555730	WallFlex Colonic Soft Stent	25 / 20	6	230	9 / 3.0	0.035 / 0.89
M00555740	WallFlex Colonic Soft Stent	25 / 20	9	230	9 / 3.0	0.035 / 0.89
M00555750	WallFlex Colonic Soft Stent	25 / 20	12	230	9 / 3.0	0.035 / 0.89
M00555760	WallFlex Colonic Soft Stent	27 / 22	6	230	9 / 3.0	0.035 / 0.89
M00555770	WallFlex Colonic Soft Stent	27 / 22	9	230	9 / 3.0	0.035 / 0.89
M00555780	WallFlex Colonic Soft Stent	27 / 22	12	230	9 / 3.0	0.035 / 0.89



Esophageal and Enteral Metal Stents

COLONIC AND DUODENAL STENTS



WallFlex™ Duodenal Stent System with Anchor Lock Delivery System



Order Number	Description	Diameter Flare / Body (mm)	Stent Length (cm)	Working Length (cm)	Catheter Diameter (F) / (mm)
M00565010	WallFlex Duodenal Stent	27 / 22	6	230	10 / 3.33
M00565020	WallFlex Duodenal Stent	27 / 22	9	230	10 / 3.33
M00565030	WallFlex Duodenal Stent	27 / 22	12	230	10 / 3.33

Packaged one per box

Recommended Guidewires: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 450cm: M00556161 Straight Dream Tip™, M00556171 Angled Dream Tip™ or .035 in (0.89mm) Jagwire™ Guidewire Stiff Shaft, 450cm: M00556601 Straight Tip, M00556611 Angled Tip or .035in (0.89mm) Super Stiff Wire, 500cm: H954180011




WallFlex™ Duodenal Soft Stent System with Anchor Lock Delivery System



Order Number	Description	Diameter Flare / Body (mm)	Stent Length (cm)	Working Length (cm)	Delivery System Diameter/ Working Channel (F/mm)	Guidewire Diameter (inch/mm)
M00555800	WallFlex Duodenal Soft Stent	23 / 18	6	230	9 / 3.0	0.035 / 0.89
M00555810	WallFlex Duodenal Soft Stent	23 / 18	9	230	9 / 3.0	0.035 / 0.89
M00555820	WallFlex Duodenal Soft Stent	23 / 18	12	230	9 / 3.0	0.035 / 0.89
M00555830	WallFlex Duodenal Soft Stent	25 / 20	6	230	9 / 3.0	0.035 / 0.89
M00555840	WallFlex Duodenal Soft Stent	25 / 20	9	230	9 / 3.0	0.035 / 0.89
M00555850	WallFlex Duodenal Soft Stent	25 / 20	12	230	9 / 3.0	0.035 / 0.89
M00555860	WallFlex Duodenal Soft Stent	27 / 22	6	230	9 / 3.0	0.035 / 0.89
M00555870	WallFlex Duodenal Soft Stent	27 / 22	9	230	9 / 3.0	0.035 / 0.89
M00555880	WallFlex Duodenal Soft Stent	27 / 22	12	230	9 / 3.0	0.035 / 0.89



Infection Prevention



Biliary	RX Biliary™ System™	Biliary Lumen Patency	Endoscopic Ultrasound (EUS)	Guidewires	Biopsy	Dilation	Imaging Systems	Pulmonary	Closure	Haemostasis	Tissue Resection	Esophageal and Enteral Metal Stents	Infection Prevention	Endobariatrics
---------	---------------------	-----------------------	-----------------------------	------------	--------	----------	-----------------	-----------	---------	-------------	------------------	-------------------------------------	----------------------	----------------

IRRIGATION & INSUFFLATION

Hydra™ Irrigation System

Order Number	Description	Compatibility	Units
SIT-575-15	Hydra Irrigation Tubing	Olympus™ AFU-100, ERBE™ EIP 2, Medivators™ EndoStratus™ & EndoGator™	box 15
SIT-577-15	Hydra Irrigation Tubing	Olympus OFP2	box 15
SIT-361-100	⊗ Hydra Disposable Auxiliary Water Jet Connectors	Olympus 160, 180 & 190 Series	box 100
SCT-609-10	Hydra Water Bottle Cap with CO2	Air/Water Tubing Length: 16" (41 cm) CO2 Tubing Length: 18" (46 cm) Olympus	box 10

CLEANING BRUSHES

HedgeHog™ Cleaning Brushes

Order Number	Description	Units
SBS-227-50	HedgeHog Single-End Channel Brush	box 50
SBD-228-50	HedgeHog Double-End Channel Brush	box 50
SBD-289-50	HedgeHog Dual-End Channel/Valve Brush	box 50
SBD-291-50	HedgeHog Sweeper Brush	box 50
SBB-382-50	HedgeHog Double Bundle	box 50

VALVES

Orca™

Order Number	Description	Units
SUV-617-50	Orca Sterile, Single Use Air/Water and Suction Valves (compatible with Olympus)	box 50
M00508001	Orca Fuji Air/Water and Suction Valves (compatible with Fuji)	box 50

OrcaPod™

Order Number	Description	Units
SUV-628-50	OrcaPod Sterile, Single Use Air/Water and Suction Valves with Seal Biopsy Valve	box 50
SUV-629-50	OrcaPod Sterile, Single Use Air/Water and Suction Valves, Seal Biopsy Valve and Hydra Water Jet Connector	box 50

Seal™

Order Number	Description	Units
SBC-501-100	Seal Single Use Biopsy Valve (compatible with both Olympus and newer Fuji endoscopes)	box 100

Porpoise™

Order Number	Description	Units
M00501891	Porpoise Air/Water Cleaning Adapter (For use with Olympus 160, 180 and 190 gastrointestinal endoscopes)	box 50

SCOPE TRANSPORT

CinchPad™ Procedure & Endoscope Transport Pad

Order Number	Description	Units
SPX-017-25	CinchPad Endoscope Transport Pad	box 25

BITE BLOCKS

Blox™ Bite Blocks

Order Number	Description	Units
SBT-546-100	Blox Bite Block 60Fr with Dental Rim	box 100
SBT-114-100	Blox Bite Block 54Fr with Dental Rim	box 100
SBT-321-100	Blox Bite Block 54Fr without Dental Rim	box 100

SPONGES

Koala™ Sponge

Order Number	Description	Units
M00501762	Koala Cleaning Sponge	box 50
SES-240-50	Koala Cleaning Sponge Imbibed w/ Pure	box 50

Compliance Endokits

Order Number	Description	Units
CIK-115-20	OP4 I-KIT 1.1 YellowSponge	box 20
CIK-122-20	OP3 I-KIT 1.1 YellowSponge	box 20

LUBRICATING GELS

EndoGlide™ Lubricating Gels

Order Number	Description	Units
M00501902	EndoGlide+ 2.0 oz Lubricant Gel	box 50
M00501905	EndoGlide+ 1.1 oz Lubricant Gel	box 100



Biliary
RX Biliary™ System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics

ENDOSCOPIC SUTURING

OverStitch™ Systems

Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Units
M00505400	OverStitch Endoscopic Suture System	*	-	box 1
M00505420	OverStitch NXT Endoscopic Suture System	*	-	box 1
M00505441	OverStitch 2-0 Polypropylene Suture	*	2.8	box 12
M00505540	OverStitch Suture Cinch	160	3.2	box 1
M00505470	OverStitch Tissue Helix	*	2.8	box 1
M00505480	OverStitch NXT Tissue Helix Pro	165	2.8	box 1

* Suitable for endoscopes with max 110 cm working length.

INTRAGASTRIC BALLOON

BIB™ System

Order Number	Description	Units
M00505680	BIB 6 Month Balloon	box 1

Orbera365™

Order Number	Description	Units
M00505690	Orbera365 12 Month Balloon	box 1
AS2290718	Endoline Viper	box 5
AS2640718	Endoline Punc	box 5

OVERSTITCH ACCESSORIES

Overtube™

Order Number	Description	Units
M00505560	OverTube Endoscopic Access System	box 1


Ensizor™

Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Units
ES26165-3	Ensizor Endoscopic Scissors Gastric Length	165	2.8	box 3
ES26235-3	Ensizor Endoscopic Scissors Colonic Length	235	2.8	box 3

INTRAGASTRIC BALLOON ACCESSORIES

Extraction Kit

Order Number	Description	Units
B-40815	Extraction Kit	box 5



**Boston
Scientific**
Advancing science for life™

www.bostonscientific.eu

© 2026 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-2174601-AC

All trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material is not intended for use in France.

Biliary	RX Biliary System™	Biliary Lumen Patency	Endoscopic Ultrasound (EUS)	Guidewires
Biopsy	Dilatation	Imaging Systems	Pulmonary	Closure
Haemostasis	Tissue Resection	Esophageal and Enteral Metal Stents	Infection Prevention	Endobariatrics



Contents

Alphabetical
Index

Product Code
Index