



Endoscopy Product Catalogue

2025

Produktkatalog
Catalogue Produits
Catálogo de Productos
Catalogo dei Prodotti

Biliary	RX Biliary System™	Biliary Lumen Patency	Endoscopic Ultrasound (EUS)	Guidewires
Biopsy	Dilation	Imaging Systems	Pulmonary	Closure
Haemostasis	Tissue Resection	Esophageal and Enteral Metal Stents	Infection Prevention	Endobariatrics



[Contents](#)



[Alphabetical
Index](#)

[Product Code
Index](#)

Biliary

Cannulas	Tandem™ XL Triple-Lumen ERCP Cannula	37
	Contour™ ERCP Cannulas	37
Sphincterotomes	TRUEtome™ Jag Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	38
	TRUEtome™ Dreamwire™ Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	38
	TRUEtome™ Hydra Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	38
	TRUEtome™ Jagwire Revolution Cannulating Sphincterotomes with Pre-loaded Jagwire™ Revolution Guidewire	38
	TRUEtome™ Cannulating Sphincterotomes	39
	Microknife™ XL Triple-Lumen Needle-Knife	39
	Ultratome XL™ Triple-Lumen Sphincterotomes	40
Stone Retrieval	Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Distal Injection	40
	Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	40
Accessories	Biliary Accessories	41

RX Biliary System™

Cannulas	Rapid Exchange XL Cannula	43
	RX ERCP Cannulas	43
Sphincterotomes	Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	44
	Dreamtome™ RX Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	44
	Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	45
	Jagtome™ Revolution RX Cannulating Sphincterotome (Autotome™ Pre-loaded with Jagwire Revolution High Performance Guidewire)	45
	Autotome™ Pro RX Cannulating Sphincterotomes	46
	Autotome™ RX Preloaded Sphincterotomes	46
	Autotome™ RX Cannulating Sphincterotomes	46
	RX Needle-Knife XL Triple-Lumen Needle-Knife	47
Stone Retrieval	Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Distal Injection	47
	Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	47
	Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Distal Injection	48
	Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	48
	Trapezoid™ RX Wireguided Retrieval Baskets	48
Dilation	Hurricane™ RX Biliary Balloon Dilatation Catheters	49
	CRE™ RX Biliary Balloon Dilatation Catheter	49
Cytology	RX Cytology Brush	49
Locking Devices	AutoCap™ RX	50
	RX Locking Devices	50
	Encore™ 26 Inflation Device	50

Biliary Lumen Patency

RX Biliary Plastic Stents	Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System	51
	Advanix™ Stents – Singles	52 – 53
	NaviFlex™ RX Delivery System	53
	RX Biliary Plastic Stents with RX Delivery System	54
Biliary Plastic Stents	Flexima™ Biliary Stents with Delivery System	55
Pancreatic Plastic Stents	Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Pigtail	56
	Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Straight	56
	Advanix™ Pancreatic Single Stents – Pigtail	57 – 58
	Advanix™ Pancreatic Single Stents – Straight	59 – 60
	NaviFlex™ RX Pancreatic Delivery System & NaviFlex™ Pushers	61
RX Biliary Metal Stents	WallFlex™ Biliary RX Fully Covered Stent System RMV (Permalume™ Covering, includes benign indication)	61
	WallFlex™ Biliary RX Partially Covered Stent System (Permalume™ Covering)	62
	WallFlex™ Biliary RX Uncovered Stent System	62
Radiofrequency Ablation Catheter and Cables	Habib™ EndoHPB Bipolar Radiofrequency Catheter	63

Endoscopic Ultrasound (EUS)

Fine Needle Aspiration (FNA)	Expect™ Slimline Handle Needles	65
Fine Needle Biopsy (FNB)	Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Devices	65
	Acquire™ S Endoscopic Ultrasound Fine Needle Biopsy (FNB) Devices	66
Lumen Apposing Metal Stent System	Hot AXIOS™ Stent and Electrocautery Enhanced Delivery System	66

Guidewires

High Performance	Dreamwire™ High Performance Guidewires	67
	Hydra Jagwire™ High Performance Guidewires	67
	Jagwire™ High Performance Guidewires	68
	Jagwire™ Revolution High Performance Guidewire	68
	NaviPro™ Hydrophilic Guidewires	69
Standard	Amplatz Super Stiff™ Guidewire	69

Biopsy

Forceps	Radial Jaw™ 4 Standard Capacity Biopsy Forceps	71
	Radial Jaw™ 4 Large Capacity Biopsy Forceps	71
	Radial Jaw™ 4 Jumbo Biopsy Forceps	72
	Radial Jaw™ 4 Gastropediatric Biopsy Forceps	72

Dilation

CRE™ (Controlled Radial Expansion) Balloon Dilatation Catheters	CRE™ PRO Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic/Biliary) (including DASE indication)	73
	CRE™ Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)	73
	CRE™ PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilatation Catheters (including DASE indication)	74
	CRE™ Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)	74
	CRE™ Fixed Wire Balloon Dilatation Catheters	75
Balloon Dilators	Rigiflex™ II Achalasia Balloon Dilators	75
Inflation Systems	Alliance™ II Inflation System	76
Dilation Accessories	Inflation/Deflation Devices and Accessories	76

Imaging Systems

Imaging Systems	EXALT™ Model D Single-Use Duodenoscope	77
	EXALT™ Model D Service Codes	77
	EXALT™ Model B Single-Use Bronchoscope	78
	EXALT™ Model B Service Codes	78
Direct Visualisation System	SpyGlass™ DS Capital Equipment	79
	SpyGlass™ DS Disposable Equipment	79
	SpyGlass™ DS Service Codes	79
	SpyGlass™ Discover Digital Catheter	80
	SpyGlass™ Discover Service Codes	80
	Autolith™ Touch Biliary EHL System	80
	Autolith™ Service Codes	80

Pulmonary

Imaging Systems	EXALT™ Model B Single-Use Bronchoscope	81
	EXALT™ Model B Service Codes	81
Endoscopic Bronchial Ultrasound	Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needles – Olympus Scope Compatible	82
	Acquire™ Pulmonary Endobronchial Ultrasound Fine Biopsy Needles – Olympus Scope Compatible	82
Airway Stents	Ultraflex™ Tracheobronchial Covered Stent System – Distal Release	83
	Ultraflex™ Tracheobronchial Uncovered Stent System – Proximal Release	83
	Ultraflex™ Tracheobronchial Uncovered Stent System – Distal Release	84
Biopsy Forceps	Radial Jaw™ 4 Pulmonary Biopsy Forceps	84
Mini-Forceps	CoreDx™ Pulmonary Mini-Forceps	84
Transbronchial Aspiration Needles	eXcelon™ Transbronchial Aspiration Needles	85
Cytology Brushes	Cellebriety™ Cytology Brushes – Bronchoscope	85
Airway Retrieval Baskets	Zero Tip™ Airway Retrieval Baskets	85
Dilation Balloons	CRE™ Pulmonary Balloon Dilatation Catheter	86
Guidewires	Jagwire™ Pulmonary Guidewire	86
Inflation System	Alliance™ II Inflation System	86

Contents

Closure

Clipping Devices	X-Tack™	87
	X-Tack™ Accessories	87
Endoscopic Suturing	OverStitch™ Systems	87
	OverStitch™ Accessories	87
	Mantis™	88
Endoluminal Vacuum Therapy	Eso-SPONGE®	88
	Endo-SPONGE®	88
	EVT Vacuum Sources	88

Haemostasis

Clipping Devices	Resolution™ Clips	89
	Resolution 360™ Clips	89
	Resolution 360™ ULTRA Clips	90
	Speedband Superview Super 7™ Multiple Band Ligators	90
Bipolar Electrohaemostasis Catheters	Injection Gold Probe™ Electrohaemostasis/Injection Catheters	90
	Gold Probe™ Electrohaemostasis Catheters	91
Needles	Interject™ Sclerotherapy Needles	91
	Interject™ Clear Sclerotherapy Needles	92
Electrosurgical Bipolar Cable Adaptor	Electrosurgical Bipolar Cable Adaptor	92

Tissue Resection

EMR Devices for Upper GI	Captivator™ EMR	93
Electrosurgical Knife	ORISE™ ProKnife	93
Snares	Rotatable Single-Use Snares	93
	Sensation™ Single-Use Short Throw Snares	94
	Captiflex™ Single-Use Snares	94
	Captivator™ Single-Use Snares	95
	Captivator™ II Single-Use Snares	95
	Captivator™ COLD Single-Use Snare	96
Hot Biopsy Forceps	Radial Jaw™ 4 Hot Biopsy Forceps	96
Accessories	Active Cords	96
Retrieval Devices	RescueNet™	97
	Rescue™ Forceps	97
Polyp Traps	TrapEase™ Polyp Traps	97



Esophageal and Enteral Metal Stents

Esophageal Stents	WallFlex™ Fully Covered RMV Esophageal Longer Loop Stent System	99
	WallFlex™ Fully Covered RMV Esophageal Stent	99
	WallFlex™ Partially Covered Esophageal Stent System	100
	WallFlex™ Fully Covered Esophageal Stent	100
	Ultraflex™ Covered Esophageal NG Stent System – Proximal Release	101
	Ultraflex™ Covered Esophageal NG Stent System – Distal Release	101
	Ultraflex™ Uncovered Esophageal NG Stent System – Proximal Release	101
	Ultraflex™ Uncovered Esophageal NG Stent System – Distal Release	101
	Agile™ Fully Covered Esophageal Stent System – Removable	102
	Agile™ Partially Covered Esophageal Stent System	102
	Agile™ OTW Fully Covered Esophageal Stent System – Removable	102
	Agile™ OTW Partially Covered Esophageal Stent System	102
Colonic and Duodenal Stents	WallFlex™ Colonic Stent System with Anchor Lock Delivery System	103
	WallFlex™ Colonic Soft Stent System with Anchor Lock Delivery System	103
	WallFlex™ Duodenal Stent System with Anchor Lock Delivery System	104
	WallFlex™ Duodenal Soft Stent System with Anchor Lock Delivery System	104

Infection Prevention

Irrigation & Insufflation	Hydra™ Irrigation System	105
Cleaning Brushes	HedgeHog™ Cleaning Brushes	105
Valves	Orca™	105
	OrcaPod™	105
	Seal™	105
	Porpoise™	105
Scope Transport	CinchPad™ Procedure & Endoscope Transport Pad	106
Bite Blocks	Blox™ Bite Blocks	106
Sponges	Koala™ Sponge	106
	Compliance Endokits	106
Lubricating Gels	EndoGlide™ Lubricating Gels	106

Endobariatrics

Endoscopic Suturing	OverStitch™ Systems	107
Intragastric Balloon	BIB™ System	107
	Orbera365™	107
Overstitch Accessories	Overtube™	107
	Ensizor™	107
Intragastric Balloon Accessories	Extraction Kit	107

Alphabetical Index

Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Devices	65
Acquire™ S Endoscopic Ultrasound Fine Needle Biopsy (FNB) Devices	66
Acquire™ Pulmonary Endobronchial Ultrasound Fine Biopsy Needles – Olympus Scope Compatible	82
Active Cords	96
Advanix™ Pancreatic Single Stents – Pigtail	57–58
Advanix™ Pancreatic Single Stents – Straight	59–60
Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Pigtail	56
Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Straight	56
Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System	51
Advanix™ Stents – Singles	52–53
Agile™ Fully Covered Esophageal Stent System – Removable	102
Agile™ OTW Fully Covered Esophageal Stent System – Removable	102
Agile™ OTW Partially Covered Esophageal Stent System	102
Agile™ Partially Covered Esophageal Stent System	102
Alliance™ II Inflation System	76
Amplatz Super Stiff™ Guidewire	69
AutoCap™ RX	50
Autolith™ Service Codes	80
Autolith™ Touch Biliary EHL System	80
Autotome™ Pro RX Cannulating Sphincterotomes	46
Autotome™ Pro RX Preloaded Sphincterotomes	46
Autotome™ RX Cannulating Sphincterotomes	46
BIB™ System	107
Biliary Accessories	41
Blox™ Bite Blocks	106
Captiflex™ Single-Use Snares	94
Captivator™ COLD Single-Use Snare	96
Captivator™ EMR	93
Captivator™ Single-Use Snares	95
Captivator™ II Single-Use Snares	95
Cellebri™ Cytology Brushes – Bronchoscope	85
CinchPad™ Procedure & Endoscope Transport Pad	106
Compliance Endokits	106
Contour™ ERCP Cannulas	37
CoreDx™ Pulmonary Mini-Forceps	84
CRE™ Fixed Wire Balloon Dilatation Catheters	75
CRE™ PRO Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Biliary) (including DASE indication)	74
CRE™ PRO Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic/Biliary) (including DASE indication)	73
CRE™ Pulmonary Balloon Dilatation Catheter	86
CRE™ RX Biliary Balloon Dilatation Catheter	49
CRE™ Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)	74
CRE™ Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)	73
Dreamtome™ RX Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	44
Dreamwire™ High Performance Guidewires	67
Electrosurgical Bipolar Cable Adaptor	92
Encore™ 26 Inflation Device	41, 50
EndoGlide™ Lubricating Gels	106
Endo-SPONGE®	88
Ensizor™	87, 107
Eso-SPONGE®	88

Alphabetical Index

EVT Vacuum Sources	88
EXALT™ Model B Service Codes	78
EXALT™ Model B Single-Use Bronchoscope	78, 81
EXALT™ Model D Service Codes	77
EXALT™ Model D Single-Use Duodenoscope	77
eXcelon™ Transbronchial Aspiration Needles	85
Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needles – Olympus Scope Compatible	82
Expect™ Slimline Handle Needles	65
Extraction Kit	107
Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Distal Injection	47
Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	47
Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Distal Injection	48
Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	48
Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Distal Injection	40
Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	40
Flexima™ Biliary Stents with Delivery System	55
Gold Probe™ Electrohaemostasis Catheters	91
Habib™ EndoHPB Bipolar Radiofrequency Catheter	63
HedgeHog™ Cleaning Brushes	105
Hot AXIOS™ Stent and Electrocautery Enhanced Delivery System	66
Hurricane™ RX Biliary Dilation Balloons	49
Hydra™ Irrigation System	105
Hydra Jagwire™ High Performance Guidewires	67
Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	44
Inflation/Deflation Devices and Accessories	76
Injection Gold Probe™ Electrohaemostasis/Injection Catheters	90
Interject™ Clear Sclerotherapy Needles	92
Interject™ Sclerotherapy Needles	91
Jagtome™ Revolution RX Cannulating Sphincterotome (Autotome™ Preloaded with Jagwire Revolution High Performance Guidewire)	45
Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	45
Jagwire™ High Performance Guidewires	68
Jagwire™ Pulmonary Guidewire	86
Jagwire™ Revolution High Performance Guidewire	68
Koala™ Sponge	106
Mantis™	88
Microknife™ XL Triple-Lumen Needle-Knife	39
NaviFlex™ RX Delivery System	53
NaviFlex™ RX Pancreatic Delivery System & NaviFlex™ Pushers	61
NaviPro™ Hydrophilic Guidewires	69
Orbera365™	107
Orca™	105
Orca Pod™	105
ORISE™ ProKnife	93
OverStitch™ Accessories	87
OverStitch™ Systems	87, 107
Overtube™	107
Porpoise™	105
Radial Jaw™ 4 Gastropediatric Biopsy Forceps	72



Alphabetical Index

Radial Jaw™ 4 Hot Biopsy Forceps	96
Radial Jaw™ 4 Jumbo Biopsy Forceps	72
Radial Jaw™ 4 Large Capacity Biopsy Forceps	71
Radial Jaw™ 4 Pulmonary Biopsy Forceps	84
Radial Jaw™ 4 Standard Capacity Biopsy Forceps	71
Rapid Exchange XL Cannula	43
RescueNet™	97
Rescue™ Forceps	97
Resolution 360™ Clips	89
Resolution 360™ ULTRA Clips	90
Resolution™ Clips	89
Rigiflex™ II Achalasia Balloon Dilators	75
Rotatable Single-Use Snares	93
RX Locking Devices	50
RX Biliary Plastic Stents with RX Delivery System	54
RX Cytology Brush	49
RX ERCP Cannulas	43
RX Needle-Knife XL Triple-Lumen Needle-Knife	47
Seal™	105
Sensation™ Single-Use Short Throw Snares	94
Speedband Superview Super 7™ Multiple Band Ligators	90
SpyGlass™ Discover	80
SpyGlass™ Discover Service Codes	80
SpyGlass™ DS Capital Equipment	79
SpyGlass™ DS Disposable Equipment	79
SpyGlass™ DS Service Codes	79
Tandem™ XL Triple-Lumen ERCP Cannula	37
TrapEase™ Polyp Traps	97
Trapezoid™ RX Wireguided Retrieval Baskets	48
TRUEtome™ Cannulating Sphincterotomes	39
TRUEtome™ Dreamwire™ Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	38
TRUEtome™ Hydra Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	38
TRUEtome™ Jag Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	38
TRUEtome™ Jagwire Revolution Cannulating Sphincterotomes with Pre-loaded Jagwire™ Revolution Guidewire	38
Ultraflex™ Covered Esophageal NG Stent System – Distal Release	101
Ultraflex™ Covered Esophageal NG Stent System – Proximal Release	101
Ultraflex™ Tracheobronchial Covered Stent System – Distal Release	83
Ultraflex™ Tracheobronchial Uncovered Stent System – Distal Release	84
Ultraflex™ Tracheobronchial Uncovered Stent System – Proximal Release	83
Ultraflex™ Uncovered Esophageal NG Stent System – Distal Release	101
Ultraflex™ Uncovered Esophageal NG Stent System – Proximal Release	101
Ultratome XL™ Triple-Lumen Sphincterotomes	40
WallFlex™ Biliary RX Fully Covered Stent System RMV (Permalume™ Covering, includes benign indication)	61
WallFlex™ Biliary RX Partially Covered Stent System (Permalume™ Covering)	62
WallFlex™ Biliary RX Uncovered Stent System	62
WallFlex™ Colonic Soft Stent System with Anchor Lock Delivery System	103
WallFlex™ Colonic Stent System with Anchor Lock Delivery System	103
WallFlex™ Duodenal Soft Stent System with Anchor Lock Delivery System	104



Alphabetical Index

WallFlex™ Duodenal Stent System with Anchor Lock Delivery System	104
WallFlex™ Fully Covered Esophageal Stent	100
WallFlex™ Fully Covered RMV Esophageal Longer Loop Stent System	99
WallFlex™ Fully Covered RMV Esophageal Stent	99
WallFlex™ Partially Covered Esophageal Stent System	100
X-Tack™	87
X-Tack™ Accessories	87
Zero Tip™ Airway Retrieval Baskets	85



Product Code Index

Order Number	Description	Page Number
Endo-SPONGE®		
5526510	Endo-SPONGE	88
5526520	Endo-SPONGE	88
5526530	Endo-SPONGE	88
Eso-SPONGE®		
5526540	Eso-SPONGE with 15 mm Overtube	88
5526550	Eso-SPONGE with 13 mm Overtube	88
EVT Vacuum Sources		
5526604	Redyrob Trans Plus	88
5526650	Endo-SPONGE® Pump	88
5526651	Endo-SPONGE® Pump Carrying Bag	88
5526653	Disposable Cannister 300 ml with tubing and Y-connector	88
19116	MV1 Mobile Programmable Low-Vacuum-System	88
18032	Collection Bottle MV1	88
18022	Bacterial Filter MV1	88
EXALT™ Model B Single-Use Bronchoscope Power Cords		
8500-098BL9	Power Cord (EU) – Plug Type F	78, 81
8520-098BL2	Power Cord (CH) – Plug Type J	78, 81
8540-098BL2	Power Cord (DK) – Plug Type K	78, 81
8550-098BL3	Power Cord (IT) – Plug Type L	78, 81
9650-098BL7	Power Cord (UK) – Plug Type G	78, 81
8560-098BL3	Power Cord (IL) – Plug Type H	78, 81
8580-098BL5	Power Cord (ZA) – Plug Type M	78, 81
Extraction Kit		
B-40815	Extraction Kit	107
Compliance Endokits		
CIK-115-20	OP4 I-KIT 1.1 YellowSponge	106
CIK-122-20	OP3 I-KIT 1.1 YellowSponge	106
RescueNet™		
DGN-538-5	RescueNet Retrieval Device	97
Ensizor™		
ES26165-3	Ensizor Endoscopic Scissors Gastric Length	87, 107
ES26135-3	Ensizor Endoscopic Scissors Colonic Length	87, 107
Habib™ EndoHPB Bipolar Radiofrequency Catheter		
M00500070	Habib Endo HPB Bipolar Radiofrequency Catheter	63
5100	Adapter Cable Bipolar, 28 mm pins	63
5420	Adapter Cable Bipolar, Internal 4 mm pin	63
5700	Adapter Cable Bipolar, 22 mm pins	63
Koala™ Sponge		
M00501762	Koala Cleaning Sponge	106
SES-240-50	Koala Cleaning Sponge Imbibed w/ Pure	106
Porpoise™		
M00501891	Porpoise™ Air/Water Cleaning Adapter (For use with Olympus 160, 180 and 190 gastrointestinal endoscopes)	105
EndoGlide™ Lubricating Gels		
M00501902	EndoGlide+ 2.0 oz Lubricant Gel	106
M00501905	EndoGlide+ 1.1 oz Lubricant Gel	106
Electrosurgical Bipolar Cable Adaptors		
M00505200	Bipolar Cable Adaptor	92



Product Code Index

Order Number	Description	Page Number
	Electrosurgical Bipolar Cable Adaptors (CONTINUED)	
M00505210	Bipolar Cable Adaptor	92
M00505220	Bipolar Cable Adaptor	92
	OverStitch™ Systems	
M00505400	OverStitch Endoscopic Suturing System (OLD UPN: ESS-G02-160)	87, 107
M00505420	OverStitch NXT Endoscopic Suture System (OLD UPN: ESS-G02-NXT)	87, 107
M00505441	OverStitch 2-0 Polypropylene Suture (OLD UPN: PLY-G02-020-APL)	87, 107
M00505470	OverStitch Tissue Helix (OLD UPN: THX-165-028)	87, 107
M00505480	OverStitch NXT Tissue Helix Pro (OLD UPN: NXT-THXP-130)	87, 107
M00505540	OverStitch Suture Cinch (OLD UPN: CNH-G01-000)	87, 107
M00505550	OverStitch Suture Cinch – Long (OLD UPN: CNH-C01-213-L)	87
	X-Tack™ Systems	
M00505500	X-Tack Endoscopic HeliX Tacking System (OLD UPN: X-TACK-160-H)	87
M00505510	X-Tack Endoscopic HeliX Tacking System (OLD UPN: XTACK-235-H)	87
	Overtube™	
M00505560	OverTube Endoscopic Access System (OLD UPN: OVT-027-160)	87, 107
	BIB™ System	
M00505680	BIB 6 Month Balloon (OLD UPN: B-40800)	107
	Orbera365™	
M00505690	Orbera365 12 Month Balloon (OLD UPN: B-50012)	107
AS2290718	Endoline Viper	107
AS2640718	Endoline Punc	107
	Trapezoid™ RX Wireguided Retrieval Baskets	
M00510860	Trapezoid RX Wireguided Retrieval Basket	48
M00510870	Trapezoid RX Wireguided Retrieval Basket	48
M00510880	Trapezoid RX Wireguided Retrieval Basket	48
M00510890	Trapezoid RX Wireguided Retrieval Basket	48
	Zero Tip™ Airway Retrieval Baskets	
M00513200	Zero Tip Airway Retrieval Basket	85
M00513210	Zero Tip Airway Retrieval Basket	85
	Radial Jaw™ 4 Large Capacity Biopsy Forceps	
M00513301	Radial Jaw 4 Large Capacity	71
M00513302	Radial Jaw 4 Large Capacity	71
M00513303	Radial Jaw 4 Large Capacity	71
M00513311	Radial Jaw 4 Large Capacity w/Needle	71
M00513312	Radial Jaw 4 Large Capacity w/Needle	71
M00513313	Radial Jaw 4 Large Capacity w/Needle	71
M00513321	Radial Jaw 4 Large Capacity	71
M00513322	Radial Jaw 4 Large Capacity	71
M00513323	Radial Jaw 4 Large Capacity	71
M00513331	Radial Jaw 4 Large Capacity w/Needle	71
M00513332	Radial Jaw 4 Large Capacity w/Needle	71
M00513333	Radial Jaw 4 Large Capacity w/Needle	71
	Radial Jaw™ 4 Jumbo Biopsy Forceps	
M00513343	Radial Jaw 4 Jumbo	72
M00513353	Radial Jaw 4 Jumbo w/Needle	72
M00513361	Radial Jaw 4 Jumbo	72
M00513362	Radial Jaw 4 Jumbo	72



Product Code Index

Order Number	Description	Page Number
	Radial Jaw™ 4 Jumbo Biopsy Forceps (CONTINUED)	
M00513363	Radial Jaw 4 Jumbo	72
M00513371	Radial Jaw 4 Jumbo w/Needle	72
M00513372	Radial Jaw 4 Jumbo w/Needle	72
M00513373	Radial Jaw 4 Jumbo w/Needle	72
	Radial Jaw™ 4 Standard Capacity Biopsy Forceps	
M00513381	Radial Jaw 4 Standard Capacity	71
M00513382	Radial Jaw 4 Standard Capacity	71
M00513383	Radial Jaw 4 Standard Capacity	71
M00513391	Radial Jaw 4 Standard Capacity w/Needle	71
M00513392	Radial Jaw 4 Standard Capacity w/Needle	71
M00513393	Radial Jaw 4 Standard Capacity w/Needle	71
M00513401	Radial Jaw 4 Standard Capacity	71
M00513402	Radial Jaw 4 Standard Capacity	71
M00513403	Radial Jaw 4 Standard Capacity	71
M00513411	Radial Jaw 4 Standard Capacity w/Needle	71
M00513412	Radial Jaw 4 Standard Capacity w/Needle	71
M00513413	Radial Jaw 4 Standard Capacity w/Needle	71
	Radial Jaw™ 4 Gastropediatric Biopsy Forceps	
M00513441	Radial Jaw 4 Gastropediatric	72
M00513443	Radial Jaw 4 Gastropediatric	72
M00513451	Radial Jaw 4 Gastropediatric w/Needle	72
M00513453	Radial Jaw 4 Gastropediatric w/Needle	72
	Ultraflex™ Uncovered Esophageal NG Stent System – Distal Release	
M00513700	Ultraflex Esophageal NG Stent System	101
M00513710	Ultraflex Esophageal NG Stent System	101
M00513720	Ultraflex Esophageal NG Stent System	101
	Ultraflex™ Covered Esophageal NG Stent System – Distal Release	
M00513730	Ultraflex Esophageal NG Stent System	101
M00513740	Ultraflex Esophageal NG Stent System	101
M00513750	Ultraflex Esophageal NG Stent System	101
M00514200	Ultraflex Esophageal NG Stent System	101
M00514210	Ultraflex Esophageal NG Stent System	101
	Ultraflex™ Uncovered Esophageal NG Stent System – Proximal Release	
M00513800	Ultraflex Esophageal NG Stent System	101
M00513810	Ultraflex Esophageal NG Stent System	101
M00513820	Ultraflex Esophageal NG Stent System	101
M00513830	Ultraflex Esophageal NG Stent System	101
M00514220	Ultraflex Esophageal NG Stent System	101
M00514230	Ultraflex Esophageal NG Stent System	101
	Ultraflex™ Covered Esophageal NG Stent System – Proximal Release	
M00513840	Ultraflex Esophageal NG Stent System	101
M00513850	Ultraflex Esophageal NG Stent System	101
M00513860	Ultraflex Esophageal NG Stent System	101
M00514240	Ultraflex Esophageal NG Stent System	101
M00514250	Ultraflex Esophageal NG Stent System	101
	Radial Jaw™ 4 Hot Biopsy Forceps	
M00515031	Radial Jaw 4 Hot Biopsy Forceps	96



Product Code Index

Order Number	Description	Page Number
	Radial Jaw™ 4 Hot Biopsy Forceps (CONTINUED)	
M00515032	Radial Jaw 4 Hot Biopsy Forceps	96
M00515033	Radial Jaw 4 Hot Biopsy Forceps	96
	Jagwire™ Pulmonary Guidewire	
M00515171	Jagwire Pulmonary Guidewire	86
	Radial Jaw™ 4 Pulmonary Biopsy Forceps	
M00515181	Radial Jaw 4 Pulmonary Standard Capacity	84
M00515182	Radial Jaw 4 Pulmonary Standard Capacity	84
M00515191	Radial Jaw 4 Pulmonary Standard Capacity w/Needle	84
M00515192	Radial Jaw 4 Pulmonary Standard Capacity w/Needle	84
M00515201	Radial Jaw 4 Pulmonary Large Capacity	84
M00515202	Radial Jaw 4 Pulmonary Large Capacity	84
	CoreDx™ Pulmonary Mini-Forceps	
M00515220	Pulmonary Mini-Forceps	84
	Celebrity™ Cytology Brushes – Bronchoscope	
M00516001	Celebrity Cytology Brush	85
M00516011	Celebrity Cytology Brush	85
M00516071	Celebrity Cytology Brush	85
M00516151	Celebrity Cytology Brush	85
	WallFlex™ Fully Covered RMV Esophageal Stent	
M00516210	WallFlex Fully Covered RMV Esophageal Stent	99
M00516220	WallFlex Fully Covered RMV Esophageal Stent	99
M00516230	WallFlex Fully Covered RMV Esophageal Stent	99
M00516240	WallFlex Fully Covered RMV Esophageal Stent	99
M00516250	WallFlex Fully Covered RMV Esophageal Stent	99
M00516260	WallFlex Fully Covered RMV Esophageal Stent	99
	WallFlex™ Fully Covered RMV Esophageal Longer Loop Stent System	
M00516270	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	99
M00516280	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	99
M00516290	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	99
M00516300	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	99
M00516310	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	99
M00516320	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	99
	WallFlex™ Fully Covered Esophageal Stent	
M00516700	WallFlex Fully Covered Esophageal Stent	100
M00516710	WallFlex Fully Covered Esophageal Stent	100
M00516720	WallFlex Fully Covered Esophageal Stent	100
M00516730	WallFlex Fully Covered Esophageal Stent	100
M00516740	WallFlex Fully Covered Esophageal Stent	100
M00516750	WallFlex Fully Covered Esophageal Stent	100
	WallFlex™ Partially Covered Esophageal Stent System	
M00516900	WallFlex Partially Covered Esophageal Stent System	100
M00516910	WallFlex Partially Covered Esophageal Stent System	100
M00516920	WallFlex Partially Covered Esophageal Stent System	100
M00516930	WallFlex Partially Covered Esophageal Stent System	100
M00516940	WallFlex Partially Covered Esophageal Stent System	100
M00516950	WallFlex Partially Covered Esophageal Stent System	100



Product Code Index

Order Number	Description	Page Number
Agile™ Partially Covered Esophageal Stent System		
M00517200	Agile Partially Covered Esophageal Stent System	102
M00517210	Agile Partially Covered Esophageal Stent System	102
M00517220	Agile Partially Covered Esophageal Stent System	102
M00517230	Agile Partially Covered Esophageal Stent System	102
M00517240	Agile Partially Covered Esophageal Stent System	102
M00517250	Agile Partially Covered Esophageal Stent System	102
M00517260	Agile Partially Covered Esophageal Stent System	102
M00517270	Agile Partially Covered Esophageal Stent System	102
Agile™ Fully Covered Esophageal Stent System – Removable		
M00517580	Agile Fully Covered Esophageal Stent System – Removable	102
M00517590	Agile Fully Covered Esophageal Stent System – Removable	102
M00517600	Agile Fully Covered Esophageal Stent System – Removable	102
M00517610	Agile Fully Covered Esophageal Stent System – Removable	102
M00517620	Agile Fully Covered Esophageal Stent System – Removable	102
M00517630	Agile Fully Covered Esophageal Stent System – Removable	102
M00517640	Agile Fully Covered Esophageal Stent System – Removable	102
M00517650	Agile Fully Covered Esophageal Stent System – Removable	102
Agile™ OTW Partially Covered Esophageal Stent System		
M00517710	Agile OTW Partially Covered Esophageal Stent System	102
M00517720	Agile OTW Partially Covered Esophageal Stent System	102
M00517730	Agile OTW Partially Covered Esophageal Stent System	102
M00517740	Agile OTW Partially Covered Esophageal Stent System	102
Agile™ OTW Fully Covered Esophageal Stent System – Removable		
M00517790	Agile OTW Fully Covered Esophageal Stent System – Removable	102
M00517800	Agile OTW Fully Covered Esophageal Stent System – Removable	102
M00517810	Agile OTW Fully Covered Esophageal Stent System – Removable	102
M00517820	Agile OTW Fully Covered Esophageal Stent System – Removable	102
Interject™ Sclerotherapy Needles		
M00518151	Interject Sclerotherapy Needle	91
M00518161	Interject Sclerotherapy Needle	91
M00518251	Interject Sclerotherapy Needle	91
M00518261	Interject Sclerotherapy Needle	91
M00518351	Interject Sclerotherapy Needle	91
M00518361	Interject Sclerotherapy Needle	91
Interject™ Clear Sclerotherapy Needles		
M00518111	Interject Sclerotherapy Needle	92
M00518301	Interject Sclerotherapy Needle	92
M00518311	Interject Sclerotherapy Needle	92
ORISE™ ProKnife		
M00519351	1.5mm ORISE ProKnife	93
M00519361	2.0mm ORISE ProKnife	93
M00519371	3.0mm ORISE ProKnife	93
Resolution 360™ Clips		
M00521230	Resolution 360 Clip	89
M00521231	Resolution 360 Clip	89
M00521232	Resolution 360 Clip	89
M00521233	Resolution 360 Clip	89



Product Code Index

Order Number	Description	Page Number
	Resolution 360™ Clips (CONTINUED)	
M00521240	Resolution 360 Clip	89
M00521241	Resolution 360 Clip	89
M00521242	Resolution 360 Clip	89
M00521243	Resolution 360 Clip	89
	Resolution 360™ ULTRA Clips	
M00521400	Resolution 360 ULTRA Clip	90
M00521401	Resolution 360 ULTRA Clip	90
M00521402	Resolution 360 ULTRA Clip	90
M00521403	Resolution 360 ULTRA Clip	90
	Mantis™	
M00521420	Mantis	88
M00521421	Mantis	88
M00521422	Mantis	88
M00521423	Mantis	88
	Resolution™ Clips	
M00522600	Resolution Clip	89
M00522601	Resolution Clip	89
M00522602	Resolution Clip	89
M00522610	Resolution Clip	89
M00522611	Resolution Clip	89
M00522612	Resolution Clip	89
	Rescue™ Forceps	
M00527451	Rescue Alligator Grasping Forceps	97
M00527551	Rescue Alligator Grasping Forceps Long	97
M00527651	Rescue Rat Tooth Grasping Forceps	97
M00527751	Rescue Rat Tooth Grasping Forceps – Pediatric	97
M00538351	Rescue Rat Tooth / Alligator Grasping Forceps	97
	Contour™ ERCP Cannulas	
M00530850	Contour ERCP Cannula	37
M00530860	Contour ERCP Cannula	37
M00530870	Contour ERCP Cannula	37
M00530880	Contour ERCP Cannula	37
M00530890	Contour ERCP Cannula	37
	Advanix™ Stents – Singles	
M00532160	Advanix Stent	53
M00532180	Advanix Stent	53
M00532190	Advanix Stent	53
M00532210	Advanix Stent	53
M00532220	Advanix Stent	53
M00532230	Advanix Stent	53
M00532240	Advanix Stent	53
M00532260	Advanix Stent	53
M00532270	Advanix Stent	53
M00532290	Advanix Stent	53
M00532300	Advanix Stent	53
M00532310	Advanix Stent	53



Product Code Index

Order Number	Description	Page Number
M00532810	Microknife™ XL Triple-Lumen Needle-Knife Microknife XL Triple-Lumen Needle-Knife	39
	Advanix™ Stents – Singles	
M00532870	Advanix Stent	52
M00532880	Advanix Stent	52
M00532890	Advanix Stent	52
M00532900	Advanix Stent	52
M00532910	Advanix Stent	52
M00532920	Advanix Stent	52
M00532930	Advanix Stent	52
M00532940	Advanix Stent	52
M00532950	Advanix Stent	52
M00532960	Advanix Stent	52
M00532970	Advanix Stent	52
M00532980	Advanix Stent	52
M00532990	Advanix Stent	52
M00533000	Advanix Stent	52
M00533010	Advanix Stent	52
M00533020	Advanix Stent	52
M00533030	Advanix Stent	52
M00533040	Advanix Stent	52
M00533230	Advanix Stent	52
M00533240	Advanix Stent	52
M00533250	Advanix Stent	52
M00533260	Advanix Stent	52
M00533270	Advanix Stent	52
M00533280	Advanix Stent	52
M00533290	Advanix Stent	52
M00533300	Advanix Stent	52
M00533310	Advanix Stent	52
M00533320	Advanix Stent	52
M00533330	Advanix Stent	52
M00533340	Advanix Stent	52
M00533350	Advanix Stent	52
M00533360	Advanix Stent	52
M00533370	Advanix Stent	52
M00533380	Advanix Stent	52
M00533390	Advanix Stent	52
M00533400	Advanix Stent	52
	NaviFlex™ RX Delivery System	
M00533560	RX Delivery System	53
M00533570	RX Delivery System	37, 53
M00533580	RX Delivery System	53
	Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System	
M00534200	Advanix Pre-loaded RX Stent	51
M00534210	Advanix Pre-loaded RX Stent	51
M00534220	Advanix Pre-loaded RX Stent	51
M00534230	Advanix Pre-loaded RX Stent	51



Product Code Index

Order Number	Description	Page Number
Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System (CONTINUED)		
M00534240	Advanix Pre-loaded RX Stent	51
M00534250	Advanix Pre-loaded RX Stent	51
M00534260	Advanix Pre-loaded RX Stent	51
M00534270	Advanix Pre-loaded RX Stent	51
M00534280	Advanix Pre-loaded RX Stent	51
M00534290	Advanix Pre-loaded RX Stent	51
M00534300	Advanix Pre-loaded RX Stent	51
M00534310	Advanix Pre-loaded RX Stent	51
M00534320	Advanix Pre-loaded RX Stent	51
M00534330	Advanix Pre-loaded RX Stent	51
M00534340	Advanix Pre-loaded RX Stent	51
M00534350	Advanix Pre-loaded RX Stent	51
M00534360	Advanix Pre-loaded RX Stent	51
M00534370	Advanix Pre-loaded RX Stent	51
M00534560	Advanix Pre-loaded RX Stent	51
M00534570	Advanix Pre-loaded RX Stent	51
M00534580	Advanix Pre-loaded RX Stent	51
M00534590	Advanix Pre-loaded RX Stent	51
M00534600	Advanix Pre-loaded RX Stent	51
M00534610	Advanix Pre-loaded RX Stent	51
M00534620	Advanix Pre-loaded RX Stent	51
M00534630	Advanix Pre-loaded RX Stent	51
M00534640	Advanix Pre-loaded RX Stent	51
M00534650	Advanix Pre-loaded RX Stent	51
M00534660	Advanix Pre-loaded RX Stent	51
M00534670	Advanix Pre-loaded RX Stent	51
M00534680	Advanix Pre-loaded RX Stent	51
M00534690	Advanix Pre-loaded RX Stent	51
M00534700	Advanix Pre-loaded RX Stent	51
M00534710	Advanix Pre-loaded RX Stent	51
M00534720	Advanix Pre-loaded RX Stent	51
M00534730	Advanix Pre-loaded RX Stent	51
NaviFlex™ RX Pancreatic Delivery System & NaviFlex™ Pushers		
M00535000	Pancreatic Stent RX Pusher	61
M00535010	Pancreatic Stent RX Pusher	61
M00535020	Pancreatic Stent Pusher	61
M00535030	Pancreatic Stent Pusher	61
M00535060	NaviFlex RX Pancreatic Stent Delivery System	61
M00535080	NaviFlex RX Pancreatic Stent Delivery System	61
Tandem™ XL Triple-Lumen ERCP Cannula		
M00535700	Tandem XL Triple-Lumen ERCP Cannula	37
Ultratome XL™ Triple-Lumen Sphincterotomes		
M00535900	Ultratome XL Short Nose	40
M00535910	Ultratome XL Long Nose	40
M00535920	Ultratome XL Short Nose	40
M00535930	Ultratome XL Long Nose	40



Product Code Index

Order Number	Description	Page Number
	Advanix™ Pancreatic Single Stents – Straight	
M00536000	Advanix Pancreatic Stent	59
M00536010	Advanix Pancreatic Stent	59
M00536020	Advanix Pancreatic Stent	59
M00536030	Advanix Pancreatic Stent	59
M00536040	Advanix Pancreatic Stent	59
M00536050	Advanix Pancreatic Stent	59
M00536060	Advanix Pancreatic Stent	59
M00536070	Advanix Pancreatic Stent	59
M00536080	Advanix Pancreatic Stent	59
M00536090	Advanix Pancreatic Stent	59
M00536100	Advanix Pancreatic Stent	59
M00536110	Advanix Pancreatic Stent	59
M00536120	Advanix Pancreatic Stent	59
M00536130	Advanix Pancreatic Stent	59
M00536140	Advanix Pancreatic Stent	59
M00536150	Advanix Pancreatic Stent	59
M00536160	Advanix Pancreatic Stent	59
M00536170	Advanix Pancreatic Stent	59
M00536180	Advanix Pancreatic Stent	59
M00536190	Advanix Pancreatic Stent	59
M00536200	Advanix Pancreatic Stent	59
M00536210	Advanix Pancreatic Stent	59
M00536220	Advanix Pancreatic Stent	59
M00536230	Advanix Pancreatic Stent	59
M00536240	Advanix Pancreatic Stent	59
M00536250	Advanix Pancreatic Stent	59
M00536260	Advanix Pancreatic Stent	59
M00536270	Advanix Pancreatic Stent	59
M00536280	Advanix Pancreatic Stent	59
M00536290	Advanix Pancreatic Stent	59
M00536300	Advanix Pancreatic Stent	59
M00536310	Advanix Pancreatic Stent	59
M00536320	Advanix Pancreatic Stent	59
M00536330	Advanix Pancreatic Stent	59
M00536340	Advanix Pancreatic Stent	59
M00536350	Advanix Pancreatic Stent	59
M00536360	Advanix Pancreatic Stent	60
M00536370	Advanix Pancreatic Stent	60
M00536380	Advanix Pancreatic Stent	60
M00536390	Advanix Pancreatic Stent	60
M00536400	Advanix Pancreatic Stent	60
M00536410	Advanix Pancreatic Stent	60
M00536420	Advanix Pancreatic Stent	60
M00536430	Advanix Pancreatic Stent	60
M00536440	Advanix Pancreatic Stent	60
M00536450	Advanix Pancreatic Stent	60
M00536460	Advanix Pancreatic Stent	60



Product Code Index

Order Number	Description	Page Number
Advanix™ Pancreatic Single Stents – Straight (CONTINUED)		
M00536470	Advanix Pancreatic Stent	60
M00536480	Advanix Pancreatic Stent	60
M00536490	Advanix Pancreatic Stent	60
M00536500	Advanix Pancreatic Stent	60
M00536510	Advanix Pancreatic Stent	60
M00536520	Advanix Pancreatic Stent	60
M00536530	Advanix Pancreatic Stent	60
M00536540	Advanix Pancreatic Stent	60
M00536550	Advanix Pancreatic Stent	60
M00536560	Advanix Pancreatic Stent	60
M00536570	Advanix Pancreatic Stent	60
M00536580	Advanix Pancreatic Stent	60
M00536590	Advanix Pancreatic Stent	60
M00536600	Advanix Pancreatic Stent	60
M00536610	Advanix Pancreatic Stent	60
M00536620	Advanix Pancreatic Stent	60
M00536630	Advanix Pancreatic Stent	60
M00536640	Advanix Pancreatic Stent	60
M00536650	Advanix Pancreatic Stent	60
M00536660	Advanix Pancreatic Stent	60
M00536670	Advanix Pancreatic Stent	60
M00536680	Advanix Pancreatic Stent	60
M00536690	Advanix Pancreatic Stent	60
M00536700	Advanix Pancreatic Stent	60
M00536710	Advanix Pancreatic Stent	60
M00536720	Advanix Pancreatic Stent	60
M00536730	Advanix Pancreatic Stent	60
Advanix™ Pancreatic Single Stents – Pigtail		
M00536740	Advanix Pancreatic Stent	57
M00536750	Advanix Pancreatic Stent	57
M00536760	Advanix Pancreatic Stent	57
M00536770	Advanix Pancreatic Stent	57
M00536780	Advanix Pancreatic Stent	57
M00536790	Advanix Pancreatic Stent	57
M00536800	Advanix Pancreatic Stent	57
M00536810	Advanix Pancreatic Stent	57
M00536820	Advanix Pancreatic Stent	57
M00536830	Advanix Pancreatic Stent	57
M00536840	Advanix Pancreatic Stent	57
M00536850	Advanix Pancreatic Stent	57
M00536860	Advanix Pancreatic Stent	57
M00536870	Advanix Pancreatic Stent	57
M00536880	Advanix Pancreatic Stent	57
M00536890	Advanix Pancreatic Stent	57
M00536900	Advanix Pancreatic Stent	57
M00536910	Advanix Pancreatic Stent	57
M00536920	Advanix Pancreatic Stent	57



Product Code Index

Order Number	Description	Page Number
	Advanix™ Pancreatic Single Stents – Pigtail (CONTINUED)	
M00536930	Advanix Pancreatic Stent	57
M00536940	Advanix Pancreatic Stent	57
M00536950	Advanix Pancreatic Stent	57
M00536960	Advanix Pancreatic Stent	57
M00536970	Advanix Pancreatic Stent	57
M00536980	Advanix Pancreatic Stent	57
M00536990	Advanix Pancreatic Stent	57
M00537000	Advanix Pancreatic Stent	57
M00537010	Advanix Pancreatic Stent	57
M00537020	Advanix Pancreatic Stent	57
M00537030	Advanix Pancreatic Stent	57
M00537040	Advanix Pancreatic Stent	57
M00537050	Advanix Pancreatic Stent	57
M00537060	Advanix Pancreatic Stent	57
M00537070	Advanix Pancreatic Stent	57
M00537080	Advanix Pancreatic Stent	57
M00537090	Advanix Pancreatic Stent	57
M00537100	Advanix Pancreatic Stent	58
M00537110	Advanix Pancreatic Stent	58
M00537120	Advanix Pancreatic Stent	58
M00537130	Advanix Pancreatic Stent	58
M00537140	Advanix Pancreatic Stent	58
M00537150	Advanix Pancreatic Stent	58
M00537160	Advanix Pancreatic Stent	58
M00537170	Advanix Pancreatic Stent	58
M00537180	Advanix Pancreatic Stent	58
M00537190	Advanix Pancreatic Stent	58
M00537200	Advanix Pancreatic Stent	58
M00537210	Advanix Pancreatic Stent	58
M00537220	Advanix Pancreatic Stent	58
M00537230	Advanix Pancreatic Stent	58
M00537240	Advanix Pancreatic Stent	58
M00537250	Advanix Pancreatic Stent	58
M00537260	Advanix Pancreatic Stent	58
M00537270	Advanix Pancreatic Stent	58
M00537280	Advanix Pancreatic Stent	58
M00537290	Advanix Pancreatic Stent	58
M00537300	Advanix Pancreatic Stent	58
M00537310	Advanix Pancreatic Stent	58
M00537320	Advanix Pancreatic Stent	58
M00537330	Advanix Pancreatic Stent	58
M00537340	Advanix Pancreatic Stent	58
	Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Pigtail	
M00537350	Advanix Pancreatic Stent Kit	56
M00537360	Advanix Pancreatic Stent Kit	56
M00537370	Advanix Pancreatic Stent Kit	56
M00537380	Advanix Pancreatic Stent Kit	56



Product Code Index

Order Number	Description	Page Number
Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Pigtail (CONTINUED)		
M00537390	Advanix Pancreatic Stent Kit	56
M00537400	Advanix Pancreatic Stent Kit	56
M00537410	Advanix Pancreatic Stent Kit	56
M00537420	Advanix Pancreatic Stent Kit	56
M00537430	Advanix Pancreatic Stent Kit	56
M00537440	Advanix Pancreatic Stent Kit	56
M00537450	Advanix Pancreatic Stent Kit	56
Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Straight		
M00537460	Advanix Pancreatic Stent Kit	56
M00537470	Advanix Pancreatic Stent Kit	56
M00537480	Advanix Pancreatic Stent Kit	56
M00537490	Advanix Pancreatic Stent Kit	56
M00537500	Advanix Pancreatic Stent Kit	56
M00537510	Advanix Pancreatic Stent Kit	56
M00537520	Advanix Pancreatic Stent Kit	56
M00537530	Advanix Pancreatic Stent Kit	56
M00537540	Advanix Pancreatic Stent Kit	56
M00537550	Advanix Pancreatic Stent Kit	56
M00537560	Advanix Pancreatic Stent Kit	56
Flexima™ Biliary Stents with Delivery System		
M00539200	Flexima Stent with Delivery System	55
M00539210	Flexima Stent with Delivery System	55
M00539220	Flexima Stent with Delivery System	55
M00539230	Flexima Stent with Delivery System	55
M00539240	Flexima Stent with Delivery System	55
M00539250	Flexima Stent with Delivery System	55
M00539260	Flexima Stent with Delivery System	55
M00539270	Flexima Stent with Delivery System	55
M00539280	Flexima Stent with Delivery System	55
M00539290	Flexima Stent with Delivery System	55
M00539300	Flexima Stent with Delivery System	55
M00539310	Flexima Stent with Delivery System	55
M00539320	Flexima Stent with Delivery System	55
M00539330	Flexima Stent with Delivery System	55
M00539340	Flexima Stent with Delivery System	55
M00539350	Flexima Stent with Delivery System	55
M00539360	Flexima Stent with Delivery System	55
M00539370	Flexima Stent with Delivery System	55
M00539380	Flexima Stent with Delivery System	55
M00539390	Flexima Stent with Delivery System	55
Speedband Superview Super 7™ Multiple Band Ligators		
M00542250	Super 7 Band Ligator	90
M00542251	Super 7 Band Ligator	90
M00542253	Super 7 Band Ligator	90
EXALT™ Model D Single-Use Duodenoscope		
M0054242CE0	EXALT Model D Single-Use Duodenoscope	77
M0054242CE1	EXALT Model D Single-Use Duodenoscope	77



Product Code Index

Order Number	Description	Page Number
	EXALT™ Model D Single-Use Duodenoscope (CONTINUED)	
M00542430	EXALT Controller	77
M00542460	EXALT Controller	77
	EXALT™ Model B Single-Use Bronchoscope	
M00542711	EXALT Model B Single-Use Bronchoscope – Slim	78, 81
M00542731	EXALT Model B Single-Use Bronchoscope – Regular	78, 81
M00542751	EXALT Model B Single-Use Bronchoscope – Large	78, 81
M00542760	EXALT Monitor (mount included)	78, 81
M00542770	EXALT Monitor (mount included) in alternative languages	78, 81
M00542780	EXALT Monitor (mount not included)	78, 81
M00542790	EXALT Monitor (mount not included) in alternative languages	78, 81
M00542820	EXALT Monitor mount	78, 81
	SpyGlass™ Discover Digital Catheter	
M00542950	SpyGlass Discover Controller	80
	RX Cytology Brush	
M00545000	RX Cytology Brush	49
	Autocap™ RX	
M00545101	AutoCap RX Olympus Compatible (TJF-160, TJF-180, TJF-190 and JF-260V Series Duodenoscopes)	50
M00545131	AutoCap RX EXALT™ Model D Compatible Integrated Biopsy	50
	Autotome™ RX Cannulating Sphincterotomes	
M00545150	Autotome RX 49	46
M00545160	Autotome RX 49	46
M00545170	Autotome RX 44	46
M00545180	Autotome RX 44	46
M00545190	Autotome RX 39	46
M00545200	Autotome RX 39	46
	RX Locking Devices	
M00545211	RX Locking Device – Olympus / FujiFilm Brand Endoscope Compatible	50
M00545261	RX Locking Device and Biopsy Cap System – Olympus / FujiFilm G5 Series Brand Endoscope Compatible	50
M00545271	RX Locking Device and Biopsy Cap System – Pentax Brand Endoscope Compatible	50
M00545281	RX Locking Device and Biopsy Cap System – FujiFilm Type 3 Series Brand Endoscope Compatible	50
	RX Biliary Plastic Stents with RX Delivery System	
M00545550	RX Biliary Plastic Stent with RX Delivery System	54
M00545560	RX Biliary Plastic Stent with RX Delivery System	54
M00545570	RX Biliary Plastic Stent with RX Delivery System	54
M00545580	RX Biliary Plastic Stent with RX Delivery System	54
M00545590	RX Biliary Plastic Stent with RX Delivery System	54
M00545600	RX Biliary Plastic Stent with RX Delivery System	54
M00545610	RX Biliary Plastic Stent with RX Delivery System	54
M00545620	RX Biliary Plastic Stent with RX Delivery System	54
M00545630	RX Biliary Plastic Stent with RX Delivery System	54
M00545640	RX Biliary Plastic Stent with RX Delivery System	54
M00545650	RX Biliary Plastic Stent with RX Delivery System	54
M00545660	RX Biliary Plastic Stent with RX Delivery System	54
M00545670	RX Biliary Plastic Stent with RX Delivery System	54
M00545680	RX Biliary Plastic Stent with RX Delivery System	54
M00545690	RX Biliary Plastic Stent with RX Delivery System	54
M00545700	RX Biliary Plastic Stent with RX Delivery System	54



Product Code Index

Order Number	Description	Page Number
	RX Biliary Plastic Stents with RX Delivery System (CONTINUED)	
M00545710	RX Biliary Plastic Stent with RX Delivery System	54
M00545720	RX Biliary Plastic Stent with RX Delivery System	54
M00545730	RX Biliary Plastic Stent with RX Delivery System	54
M00545740	RX Biliary Plastic Stent with RX Delivery System	54
	RX ERCP Cannulas	
M00545760	RX ERCP Cannula Standard Tip	43
M00545770	RX ERCP Cannula Tapered Tip	43
M00545780	RX ERCP Cannula Ball Tip	43
	Rapid Exchange XL Cannula	
M00545820	Rapid Exchange XL Cannula	43
	RX Needle-Knife XL Triple-Lumen Needle-Knife	
M00545840	RX Needle-Knife XL Triple-Lumen Needle-Knife	47
	Hurricane™ RX Biliary Dilation Balloons	
M00545890	Hurricane RX Dilation Balloon	49
M00545900	Hurricane RX Dilation Balloon	49
M00545910	Hurricane RX Dilation Balloon	49
M00545920	Hurricane RX Dilation Balloon	49
M00545930	Hurricane RX Dilation Balloon	49
M00545940	Hurricane RX Dilation Balloon	49
M00545950	Hurricane RX Dilation Balloon	49
M00545960	Hurricane RX Dilation Balloon	49
	SpyGlass™ DS Direct Visualisation System	
M00546470	SpyBite™ Max Biopsy Forceps	79
M00546550	SpyGlass™ Retrieval Basket	79
M00546560	SpyGlass™ Retrieval Snare	79
M00546610	SpyScope™ DS II Access and Delivery Catheter	79
M00546650	SpyGlass DS Digital Controller	79
	Autolith™ Touch Biliary EHL System	
M00546620	1.9Fr, 375cm Biliary EHL Probe	80
M005466800	Autolith Touch EHL Generator	80
M00546750	Autolith Touch Extender Cable	80
M00546760	Autolith Touch Foot Pedal	80
	SpyGlass™ Discover Digital Catheter	
M00546780	SpyGlass Discover Digital Catheter	80
	Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Distal Injection	
M00547000	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	47
M00547010	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	47
M00547020	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	47
	Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Proximal Injection	
M00547030	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	47
M00547040	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	47
M00547050	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	47
	Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Distal Injection	
M00547100	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	40
M00547110	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	40
M00547120	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	40



Product Code Index

Order Number	Description	Page Number
Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Proximal Injection		
M00547130	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	40
M00547140	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	40
M00547150	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	40
Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Distal Injection		
M00547300	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	48
M00547310	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	48
M00547320	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	48
Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Proximal Injection		
M00547330	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	48
M00547340	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	48
M00547350	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	48
Amplatz Super Stiff™ Guidewire		
M00550090	Amplatz Super Stiff Guidewire	69
CRE™ Pulmonary Balloon Dilatation Catheter		
M00550300	CRE Pulmonary Balloon Dilatation Catheter	86
M00550310	CRE Pulmonary Balloon Dilatation Catheter	86
M00550320	CRE Pulmonary Balloon Dilatation Catheter	86
M00550330	CRE Pulmonary Balloon Dilatation Catheter	86
M00550340	CRE Pulmonary Balloon Dilatation Catheter	86
M00550350	CRE Pulmonary Balloon Dilatation Catheter	86
Alliance™ II Inflation System		
M00550601	Alliance II Single-Use Syringe/Gauge Assembly	76, 86
M00550620	Alliance II Integrated Inflation/Lithotripsy Device	48, 76, 86
Inflation/Deflation Devices and Accessories		
M00550630	CRE™ SteriFlate Disposable Inflation Device	76
M00553200	Pneumatic Inflator	76
Acquire™ Pulmonary Endobronchial Ultrasound Fine Biopsy Needles – Olympus Scope Compatible		
M00552350	22ga Acquire Pulmonary – Olympus scope compatible	82
M00552351	22ga Acquire Pulmonary – Olympus scope compatible	82
M00552400	22ga Acquire Pulmonary – Olympus scope compatible	82
M00552401	25ga Acquire Pulmonary – Olympus scope compatible	82
Hot AXIOS™ Stent and Electrocautery Enhanced Delivery System		
M00553520	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	66
M00553530	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	66
M00553540	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	66
M00553550	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	66
M00553560	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	66
M00553570	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	66
Rigiflex™ II Achalasia Balloon Dilators		
M00554500	Rigiflex II Achalasia Balloon Dilator	75
M00554510	Rigiflex II Achalasia Balloon Dilator	75
M00554520	Rigiflex II Achalasia Balloon Dilator	75
Expect™ Slimline Handle Needles		
M00555500	Expect Slimline Handle Needle	65
M00555501	Expect Slimline Handle Needle	65
M00555510	Expect Slimline Handle Needle	65
M00555511	Expect Slimline Handle Needle	65



Product Code Index

Order Number	Description	Page Number
	Expect™ Slimline Handle Needles (CONTINUED)	
M0055520	Expect Slimline Handle Needle	65
M0055521	Expect Slimline Handle Needle	65
M0055530	Expect Slimline Handle Needle	65
M0055531	Expect Slimline Handle Needle	65
	Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Devices	
M0055540	22ga Acquire EUS-FNB Needle	65
M0055541	22ga Acquire EUS-FNB Needle	65
M0055560	25ga Acquire EUS-FNB Needle	65
M0055561	25ga Acquire EUS-FNB Needle	65
M0055580	19ga Acquire FNB Needle, Flexible	65
	Acquire™ S Endoscopic Ultrasound Fine Needle Biopsy (FNB) Devices	
M00555640	22ga Acquire S EUS-FNB Needle	66
M00555641	22ga Acquire S EUS-FNB Needle	66
M00555660	25ga Acquire S EUS-FNB Needle	66
M00555661	25ga Acquire S EUS-FNB Needle	66
	WallFlex™ Colonic Soft Stent System with Anchor Lock Delivery System	
M00555700	WallFlex Colonic Soft Stent	103
M00555710	WallFlex Colonic Soft Stent	103
M00555720	WallFlex Colonic Soft Stent	103
M00555730	WallFlex Colonic Soft Stent	103
M00555740	WallFlex Colonic Soft Stent	103
M00555750	WallFlex Colonic Soft Stent	103
M00555760	WallFlex Colonic Soft Stent	103
M00555770	WallFlex Colonic Soft Stent	103
M00555780	WallFlex Colonic Soft Stent	103
	WallFlex™ Duodenal Soft Stent System with Anchor Lock Delivery System	
M00555800	WallFlex Duodenal Soft Stent	104
M00555810	WallFlex Duodenal Soft Stent	104
M00555820	WallFlex Duodenal Soft Stent	104
M00555830	WallFlex Duodenal Soft Stent	104
M00555840	WallFlex Duodenal Soft Stent	104
M00555850	WallFlex Duodenal Soft Stent	104
M00555860	WallFlex Duodenal Soft Stent	104
M00555870	WallFlex Duodenal Soft Stent	104
M00555880	WallFlex Duodenal Soft Stent	104
	Hydra Jagwire™ High Performance Guidewires	
M00556001	Hydra Jagwire Standard	67
M00556011	Hydra Jagwire Standard	67
M00556021	Hydra Jagwire Stiff Shaft	67
M00556031	Hydra Jagwire Stiff Shaft	67
M00556041	Hydra Jagwire Standard	67
M00556051	Hydra Jagwire Standard	67
M00556061	Hydra Jagwire Stiff Shaft	67
M00556071	Hydra Jagwire Stiff Shaft	67
	Dreamwire™ High Performance Guidewires	
M00556101	Dreamwire Standard	67
M00556111	Dreamwire Standard	67



Product Code Index

Order Number	Description	Page Number
Dreamwire™ High Performance Guidewires (CONTINUED)		
M00556121	Dreamwire Stiff Shaft	67
M00556131	Dreamwire Stiff Shaft	67
M00556141	Dreamwire Standard	67
M00556151	Dreamwire Standard	67
M00556161	Dreamwire Stiff Shaft	67
M00556171	Dreamwire Stiff Shaft	67
NaviPro™ Hydrophilic Guidewires		
M00556181	NaviPro Standard	69
M00556191	NaviPro Standard	69
M00556201	NaviPro Standard	69
M00556211	NaviPro Standard	69
M00556221	NaviPro Standard	69
M00556231	NaviPro Standard	69
M00556241	NaviPro Stiff Shaft	69
M00556251	NaviPro Stiff Shaft	69
Jagwire™ High Performance Guidewires		
M00556461	Jagwire Guidewire ST	68
M00556481	Jagwire Guidewire ST Stiff	68
M00556491	Jagwire Guidewire ANG Stiff	68
M00556511	Jagwire Guidewire ANG	68
M00556521	Jagwire Guidewire ST Stiff	68
M00556531	Jagwire Guidewire ANG Stiff	68
M00556561	Jagwire Guidewire ST	68
M00556581	Jagwire Guidewire ST	68
M00556591	Jagwire Guidewire ANG	68
M00556601	Jagwire Guidewire ST Stiff	68
M00556611	Jagwire Guidewire ANG Stiff	68
M00556621	Jagwire Guidewire ST	68
M00557661	Jagwire Guidewire ST Stiff	68
M00557651	Jagwire Guidewire ANG	68
M00556661	Jagwire Guidewire ANG Stiff	68
M00556841	Extendable Jagwire Guidewire ST	68
M00556851	Extendable Jagwire Guidewire ANG	68
M00556861	Extendable Jagwire Guidewire ST Stiff	68
M00556871	Extendable Jagwire Guidewire ANG Stiff	68
M00556902	Jagtail™ Guidewire Extension	68
M00557641	Jagwire Guidewire ST	68
Jagwire™ Revolution High Performance Guidewire		
M00557001	Jagwire Revolution High Performance Guidewire	68
M00557011	Jagwire Revolution High Performance Guidewire	68
M00557021	Jagwire Revolution High Performance Guidewire	68
M00557031	Jagwire Revolution High Performance Guidewire	68
M00557041	Jagwire Revolution High Performance Guidewire	68
Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needles – Olympus Scope Compatible		
M00558220	Expect Pulmonary 22ga – Olympus Scope Compatible	82
M00558221	Expect Pulmonary 22ga – Olympus Scope Compatible	82
M00558250	Expect Pulmonary 25ga – Olympus Scope Compatible	82



Product Code Index

Order Number	Description	Page Number
	Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needles – Olympus Scope Compatible (CONTINUED)	
M00558251	Expect Pulmonary 25ga – Olympus Scope Compatible	82
M00558731	Expect Pulmonary Needle Adaptor – Olympus Scope Compatible	82
	CRE™ Fixed Wire Balloon Dilators	
M00558330	CRE Fixed Wire Balloon Dilator	75
M00558340	CRE Fixed Wire Balloon Dilator	75
M00558350	CRE Fixed Wire Balloon Dilator	75
M00558360	CRE Fixed Wire Balloon Dilator	75
M00558361	CRE Fixed Wire Balloon Dilator	75
M00558370	CRE Fixed Wire Balloon Dilator	75
M00558371	CRE Fixed Wire Balloon Dilator	75
M00558380	CRE Fixed Wire Balloon Dilator	75
M00558381	CRE Fixed Wire Balloon Dilator	75
	CRE™ Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)	
M00558390	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)	74
M00558400	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)	74
M00558410	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)	74
M00558420	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)	74
M00558430	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)	74
M00558440	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric)	74
	CRE™ Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)	
M00558450	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)	73
M00558460	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)	73
M00558470	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)	73
M00558480	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)	73
M00558490	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)	73
M00558500	CRE Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)	73
	CRE™ PRO Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Biliary) (including DASE indication)	
M00558600	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	74
M00558610	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	74
M00558620	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	74
M00558630	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	74
M00558640	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	74
M00558650	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	74
	CRE™ PRO Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic/Biliary) (including DASE indication)	
M00558660	CRE PRO Wireguided Esophageal/Pyloric/Colonic/Biliary Balloon Dilators (including DASE indication)	73
M00558670	CRE PRO Wireguided Esophageal/Pyloric/Colonic/Biliary Balloon Dilators (including DASE indication)	73
M00558680	CRE PRO Wireguided Esophageal/Pyloric/Colonic/Biliary Balloon Dilators (including DASE indication)	73
M00558690	CRE PRO Wireguided Esophageal/Pyloric/Colonic/Biliary Balloon Dilators (including DASE indication)	73
M00558700	CRE PRO Wireguided Esophageal/Pyloric/Colonic/Biliary Balloon Dilators (including DASE indication)	73
M00558710	CRE PRO Wireguided Esophageal/Pyloric/Colonic/Biliary Balloon Dilators (including DASE indication)	73
	CRE™ RX Biliary Balloon Dilatation Catheter	
M00558920	CRE RX Biliary Balloon Dilatation Catheter	49
M00558930	CRE RX Biliary Balloon Dilatation Catheter	49
M00558940	CRE RX Biliary Balloon Dilatation Catheter	49
M00558860	CRE RX Biliary Balloon Dilatation Catheter	49
M00558870	CRE RX Biliary Balloon Dilatation Catheter	49
M00558880	CRE RX Biliary Balloon Dilatation Catheter	49



Product Code Index

Order Number	Description	Page Number
	CRE™ RX Biliary Balloon Dilatation Catheter (CONTINUED)	
M00558890	CRE RX Biliary Balloon Dilatation Catheter	49
M00558900	CRE RX Biliary Balloon Dilatation Catheter	49
M00558910	CRE RX Biliary Balloon Dilatation Catheter	49
	Gold Probe™ Electrohaemostasis Catheters	
M00560070	Gold Probe Catheter	91
M00560100	Gold Probe Catheter	91
M00560220	Gold Probe Catheter	91
	Injection Gold Probe™ Electrohaemostasis/Injection Catheters	
M00560150	Injection Gold Probe Catheter	90
M00560160	Injection Gold Probe Catheter	90
	Sensation™ Single-Use Short Throw Snares	
M00560311	Medium Oval	94
M00560321	Small Oval	94
M00560331	Medium Crescent	94
	Captivator™ COLD Single-Use Snare	
M00561101	10mm Rounded–Stiff	96
M00561102	10mm Rounded–Stiff	96
M00561103	10mm Rounded–Stiff	96
	Captivator™ II Single-Use Snares	
M00561191	25 mm Rounded	95
M00561192	25 mm Rounded	95
M00561193	25 mm Rounded	95
M00561221	10 mm Rounded	95
M00561222	10 mm Rounded	95
M00561223	10 mm Rounded	95
M00561231	15 mm Rounded	95
M00561232	15 mm Rounded	95
M00561233	15 mm Rounded	95
M00561241	20 mm Rounded	95
M00561242	20 mm Rounded	95
M00561243	20 mm Rounded	95
M00561291	33 mm Rounded	95
	Accessories	
M00561250	Active Cord, Fits Olympus Brand Generator	41, 96
M00561270	Active Cord, Fits Endostat™, Endostat™ II, Endostat™ III, Bovie®, Valleylab®, Aspen Labs and ERBE® Brand Generators	41, 96
	Captivator™ EMR	
M00561600	Captivator EMR for Standard Gastrosopes	93
M00561601	Captivator EMR for Standard Gastrosopes	93
M00561610	Captivator EMR for Large Gastrosopes	93
M00561611	Captivator EMR for Large Gastrosopes	93
	Rotatable Single-Use Snares	
M00561821	Small Oval	93
M00561831	Medium Oval	93
	Captivator™ Single-Use Snares	
M00561311	Large Oval	95
M00562301	Small Oval	95
M00562321	Medium Oval	95



Product Code Index

Order Number	Description	Page Number
	Captivator™ Single-Use Snares (CONTINUED)	
M00562341	Medium Hexagonal	95
M00562371	Medium Crescent	95
M00562391	Large Oval (thin wire)	95
M00562451	Small Hexagonal	95
	Captiflex™ Single-Use Snares	
M00562401	Medium Oval	94
M00562402	Medium Oval	94
M00562421	Small Oval	94
M00562422	Small Oval	94
M00562471	Extra Small Oval	94
	Sensation™ Single-Use Short Throw Snares	
M00562651	Large Oval	94
M00562652	Large Oval	94
M00562671	Medium Oval	94
M00562672	Medium Oval	94
M00562673	Medium Oval	94
M00562691	Small Oval	94
M00562692	Small Oval	94
M00562693	Small Oval	94
	eXcelon™ Transbronchial Aspiration Needles	
M00564101	eXcelon Transbronchial Aspiration Needle	85
M00564111	eXcelon Transbronchial Aspiration Needle	85
M00564121	eXcelon Transbronchial Aspiration Needle	85
	Ultraflex™ Tracheobronchial Uncovered Stent System – Distal Release	
M00564500	Ultraflex Tracheobronchial Stent System	84
M00564510	Ultraflex Tracheobronchial Stent System	84
M00564520	Ultraflex Tracheobronchial Stent System	84
M00564530	Ultraflex Tracheobronchial Stent System	84
M00564540	Ultraflex Tracheobronchial Stent System	84
M00564560	Ultraflex Tracheobronchial Stent System	84
M00564570	Ultraflex Tracheobronchial Stent System	84
M00564590	Ultraflex Tracheobronchial Stent System	84
M00564600	Ultraflex Tracheobronchial Stent System	84
	Ultraflex™ Tracheobronchial Uncovered Stent System – Proximal Release	
M00564640	Ultraflex Tracheobronchial Stent System	83
M00564650	Ultraflex Tracheobronchial Stent System	83
M00564660	Ultraflex Tracheobronchial Stent System	83
M00564670	Ultraflex Tracheobronchial Stent System	83
M00564680	Ultraflex Tracheobronchial Stent System	83
M00564690	Ultraflex Tracheobronchial Stent System	83
M00564700	Ultraflex Tracheobronchial Stent System	83
M00564710	Ultraflex Tracheobronchial Stent System	83
M00564720	Ultraflex Tracheobronchial Stent System	83
	Ultraflex™ Tracheobronchial Covered Stent System – Distal Release	
M00564740	Ultraflex Tracheobronchial Stent System	83
M00564750	Ultraflex Tracheobronchial Stent System	83
M00564760	Ultraflex Tracheobronchial Stent System	83



Product Code Index

Order Number	Description	Page Number
	Ultraflex™ Tracheobronchial Covered Stent System – Distal Release (CONTINUED)	
M00564770	Ultraflex Tracheobronchial Stent System	83
M00564780	Ultraflex Tracheobronchial Stent System	83
M00564790	Ultraflex Tracheobronchial Stent System	83
M00564800	Ultraflex Tracheobronchial Stent System	83
M00564810	Ultraflex Tracheobronchial Stent System	83
M00564820	Ultraflex Tracheal NG Stent System	83
M00564830	Ultraflex Tracheal NG Stent System	83
M00564840	Ultraflex Tracheal NG Stent System	83
M00564850	Ultraflex Tracheal NG Stent System	83
M00564860	Ultraflex Tracheal NG Stent System	83
M00564870	Ultraflex Tracheal NG Stent System	83
M00564880	Ultraflex Tracheal NG Stent System	83
M00564890	Ultraflex Tracheal NG Stent System	83
M00564900	Ultraflex Tracheal NG Stent System	83
M00564910	Ultraflex Tracheal NG Stent System	83
	WallFlex™ Duodenal Stent System with Anchor Lock Delivery System	
M00565010	WallFlex Duodenal Stent	104
M00565020	WallFlex Duodenal Stent	104
M00565030	WallFlex Duodenal Stent	104
	WallFlex™ Colonic Stent System with Anchor Lock Delivery System	
M00565040	WallFlex Colonic Stent	103
M00565050	WallFlex Colonic Stent	103
M00565060	WallFlex Colonic Stent	103
M00565070	WallFlex Colonic Stent	103
M00565080	WallFlex Colonic Stent	103
M00565090	WallFlex Colonic Stent	103
M00565100	WallFlex Colonic Stent	103
M00565110	WallFlex Colonic Stent	103
M00565120	WallFlex Colonic Stent	103
M00565130	WallFlex Colonic Stent	103
M00565140	WallFlex Colonic Stent	103
M00565150	WallFlex Colonic Stent	103
	Accessories	
M00566001	Rapid Refill Continuous Injection System	50
	Encore™ 26 Inflation Device	
M00566670	Encore 26 Inflation Device	50
	WallFlex™ Biliary RX Fully Covered Stent System RMV (Permalume™ Covering, includes benign indication)	
M00570450	WallFlex RX Biliary (Permalume Covering, includes benign indication)	61
M00570460	WallFlex RX Biliary (Permalume Covering, includes benign indication)	61
M00570470	WallFlex RX Biliary (Permalume Covering, includes benign indication)	61
M00570480	WallFlex RX Biliary (Permalume Covering, includes benign indication)	61
M00570490	WallFlex RX Biliary (Permalume Covering, includes benign indication)	61
	WallFlex™ Biliary RX Uncovered Stent System	
M00570600	WallFlex RX Biliary	62
M00570610	WallFlex RX Biliary	62
M00570620	WallFlex RX Biliary	62
M00570630	WallFlex RX Biliary	62



Product Code Index

Order Number	Description	Page Number
	WallFlex™ Biliary RX Uncovered Stent System (CONTINUED)	
M00570640	WallFlex RX Biliary	62
M00570650	WallFlex RX Biliary	62
M00570660	WallFlex RX Biliary	62
M00570890	WallFlex RX Biliary	62
	WallFlex™ Biliary RX Partially Covered Stent System (Permalume™ Covering)	
M00570700	WallFlex RX Biliary (Permalume Covering)	62
M00570710	WallFlex RX Biliary (Permalume Covering)	62
M00570720	WallFlex RX Biliary (Permalume Covering)	62
M00570730	WallFlex RX Biliary (Permalume Covering)	62
M00570740	WallFlex RX Biliary (Permalume Covering)	62
	Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	
M00573040	Jagtome RX 44	45
M00573050	Jagtome RX 44	45
M00573060	Jagtome RX 44	45
M00573070	Jagtome RX 44	45
M00573080	Jagtome RX 39	45
M00573090	Jagtome RX 39	45
M00573100	Jagtome RX 39	45
M00573110	Jagtome RX 39	45
	WallFlex™ Biliary RX Fully Covered Stent System RMV (Permalume™ Covering, includes benign indication)	
M00574210	WallFlex RX Biliary (Permalume Covering, includes benign indication)	61
M00574230	WallFlex RX Biliary (Permalume Covering, includes benign indication)	61
M00574240	WallFlex RX Biliary (Permalume Covering, includes benign indication)	61
	WallFlex™ Biliary RX Partially Covered Stent System (Permalume™ Covering)	
M00576740	WallFlex RX Biliary (Permalume Covering)	62
M00576750	WallFlex RX Biliary (Permalume Covering)	62
M00576760	WallFlex RX Biliary (Permalume Covering)	62
M00576770	WallFlex RX Biliary (Permalume Covering)	62
	WallFlex™ Biliary RX Uncovered Stent System	
M00576780	WallFlex RX Biliary	62
M00576790	WallFlex RX Biliary	62
	WallFlex™ Biliary RX Fully Covered System RMV (Permalume™ Covering, includes benign indication)	
M00576800	WallFlex RX Biliary (Permalume Covering, includes benign indication)	61
M00576810	WallFlex RX Biliary (Permalume Covering, includes benign indication)	61
M00576820	WallFlex RX Biliary (Permalume Covering, includes benign indication)	61
M00576830	WallFlex RX Biliary (Permalume Covering, includes benign indication)	61
	Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	
M00583000	Hydratome RX 49	44
M00583010	Hydratome RX 49	44
M00583020	Hydratome RX 49	44
M00583030	Hydratome RX 49	44
M00583040	Hydratome RX 44	44
M00583050	Hydratome RX 44	44
M00583060	Hydratome RX 44	44
M00583070	Hydratome RX 44	44
	TRUEtome™ Jag Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	
M00583080	TRUEtome Jag 39 Sphincterotome	38



Product Code Index

Order Number	Description	Page Number
	TRUEtome™ Jag Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire (CONTINUED)	
M00583090	TRUEtome Jag 44 Sphincterotome	38
M00583100	TRUEtome Jag 44 Sphincterotome	38
M00583110	TRUEtome Jag 44 Sphincterotome	38
M00583120	TRUEtome Jag 39 Sphincterotome	38
M00583130	TRUEtome Jag 44 Sphincterotome	38
M00583140	TRUEtome Jag 44 Sphincterotome	38
M00583150	TRUEtome Jag 44 Sphincterotome	38
	TRUEtome™ Hydra Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	
M00583160	TRUEtome Hydra 44 Sphincterotome	38
M00583170	TRUEtome Hydra 44 Sphincterotome	38
	TRUEtome™ Dreamwire™ Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	
M00583190	TRUEtome Dreamwire 44 Sphincterotome	38
M00583210	TRUEtome Dreamwire 44 Sphincterotome	38
	TRUEtome™ Jagwire Revolution Cannulating Sphincterotomes with Pre-loaded Jagwire™ Revolution Guidewire	
M00583260	TRUEtome 39 Sphincterotome	38
M00583270	TRUEtome 44 Sphincterotome	38
M00583280	TRUEtome 39 Sphincterotome	38
M00583290	TRUEtome 44 Sphincterotome	38
	Dreamtome™ RX Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	
M00584000	Dreamtome RX 49 Sphincterotome	44
M00584010	Dreamtome RX 49 Sphincterotome	44
M00584020	Dreamtome RX 49 Sphincterotome	44
M00584030	Dreamtome RX 49 Sphincterotome	44
M00584040	Dreamtome RX 44 Sphincterotome	44
M00584050	Dreamtome RX 44 Sphincterotome	44
M00584060	Dreamtome RX 44 Sphincterotome	44
M00584070	Dreamtome RX 44 Sphincterotome	44
	TRUEtome™ Cannulating Sphincterotomes	
M00584120	TRUEtome 39 Sphincterotome	39
M00584121	TRUEtome 39 Sphincterotome	39
M00584122	TRUEtome 39 Sphincterotome	39
M00584150	TRUEtome 49 Sphincterotome	39
M00584151	TRUEtome 49 Sphincterotome	39
M00584152	TRUEtome 49 Sphincterotome	39
M00584160	TRUEtome 49 Sphincterotome	39
M00584161	TRUEtome 49 Sphincterotome	39
M00584162	TRUEtome 49 Sphincterotome	39
M00584170	TRUEtome 44 Sphincterotome	39
M00584171	TRUEtome 44 Sphincterotome	39
M00584172	TRUEtome 44 Sphincterotome	39
M00584180	TRUEtome 44 Sphincterotome	39
M00584181	TRUEtome 44 Sphincterotome	39
M00584182	TRUEtome 44 Sphincterotome	39
M00584190	TRUEtome 39 Sphincterotome	39
M00584191	TRUEtome 39 Sphincterotome	39
M00584192	TRUEtome 39 Sphincterotome	39



Product Code Index

Order Number	Description	Page Number
	Jagtome™ Revolution RX Cannulating Sphincterotome (Autotome™ Preloaded with Jagwire Revolution High Performance Guidewire)	
M00584240	Jagtome Revolution RX 39 Sphincterotome	45
M00584250	Jagtome Revolution RX 39 Sphincterotome	45
M00584260	Jagtome Revolution RX 39 Sphincterotome	45
M00584270	Jagtome Revolution RX 39 Sphincterotome	45
	Autotome™ Pro RX Cannulating Sphincterotomes	
M00584380	Autotome Pro RX 44	46
M00584390	Autotome Pro RX 39	46
	Autotome™ Pro RX Preloaded Sphincterotomes	
M00584400	Dreamtome Pro RX 44	46
M00584410	Dreamtome Pro RX 44	46
M00584420	Hydratome Pro RX 44	46
M00584430	Hydratome Pro RX 44	46
M00584440	Jagtome Pro RX 44	46
M00584450	Jagtome Pro RX 44	46
M00584460	Jagtome Pro RX 39	46
M00584470	Jagtome Revolution Pro RX 39	46
M00584480	Jagtome Revolution Pro RX 39	46
	Seal™	
SBC-501-100	Seal Single Use Biopsy Valve	105
	HedgeHog™ Cleaning Brushes	
SBD-228-50	HedgeHog Double-End Channel Brush	105
SBD-289-50	HedgeHog Dual-End Channel/Valve Brush	105
SBD-291-50	HedgeHog Sweeper Brush	105
SBS-227-50	HedgeHog Single-End Channel Brush	105
SBB-382-50	HedgeHog Double Bundle	105
	Blox™ Bite Blocks	
SBT-114-100	Blox Bite Block 54Fr with Dental Rim	106
SBT-321-100	Blox Bite Block 54Fr without Dental Rim	106
SBT-546-100	Blox Bite Block 60Fr with Dental Rim	106
	Hydra™ Irrigation System	
SCT-609-10	Hydra Water Bottle Cap with CO2	105
SIT-361-100	Hydra Disposable Auxiliary Water Jet Connectors	105
SIT-575-15	Hydra Irrigation Tubing	79, 105
SIT-577-15	Hydra Irrigation Tubing	105
	CinchPad™ Procedure & Endoscope Transport Pad	
SPX-017-25	CinchPad Endoscope Transport Pad	106
	Autolith™ Service Codes	
SRVM005ATLT1001	Autolith EssentialCare POS 1 Year	80
SRVM005ATLT1002	Autolith EssentialCare POS 2 Years	80
SRVM005ATLT1003	Autolith EssentialCare POS 3 Years	80
SRVM005ATLT1004	Autolith EssentialCare POS 4 Years	80
SRVM005ATLTECP1001	Autolith EssentialCare PLUS POS 3 Years	80
SRVM005ATLTECP1002	Autolith EssentialCare PLUS POS 4 Years	80
SRVM005ATLTECP1003	Autolith EssentialCare PLUS POS 5 Years	80
SRVM005ATLTECP1004	Autolith EssentialCare PLUS POST 3 Years	80
SRVM005ATLTECP1005	Autolith EssentialCare PLUS POST 4 Years	80
SRVM005ATLTECP1006	Autolith EssentialCare PLUS POST 5 Years	80



Product Code Index



Order Number	Description	Page Number
EXALT™ Model D Service Codes		
SRVM005EXALTDTC01	EXALT Model D TotalCare POS 1 Year	77
SRVM005EXALTDTC02	EXALT Model D TotalCare POS 2 Years	77
SRVM005EXALTDTC03	EXALT Model D TotalCare POS 3 Years	77
SRVM005EXALTDTC04	EXALT Model D TotalCare POS 4 Years	77
SRVM005EXALTDTC05	EXALT Model D TotalCare POS 5 Years	77
SRVM005EXALTDTC06	EXALT Model D TotalCare POST 1 Year	77
SRVM005EXALTDTC07	EXALT Model D TotalCare POST 2 Years	77
SRVM005EXALTDTC08	EXALT Model D TotalCare POST 3 Years	77
SRVM005EXALTDTC09	EXALT Model D TotalCare POST 4 Years	77
SRVM005EXALTDTC10	EXALT Model D TotalCare POST 5 Years	77
EXALT™ Model B Service Codes		
SVC20EXALTBANNFEE	EXALT B Annual Rental Fee for Capital	78, 81
SVC20EXALTBDFEE	EXALT B Daily Rental Fee for Capital	78, 81
SVC20EXALTBDMNTFEE	EXALT B Monthly Rental Fee for Capital	78, 81
SRVM005EXBMONEC01	EXALT B Monitor EssentialCare POS 1 Year	78, 81
SRVM005EXBMONEC02	EXALT B Monitor EssentialCare POS 2 Years	78, 81
SRVM005EXBMONEC05	EXALT B Monitor EssentialCare POST 1 Year	78, 81
SRVM005EXBMONEC06	EXALT B Monitor EssentialCare POST 2 Years	78, 81
SRVM005EXBMONECP01	EXALT B Monitor EssCare PLUS POS 1 Year	78, 81
SRVM005EXBMONECP02	EXALT B Monitor EssCare PLUS POS 2 Years	78, 81
SRVM005EXBMONECP05	EXALT B Monitor EssCare PLUS POST 1 Year	78, 81
SRVM005EXBMONECP06	EXALT B Monitor EssCare PLUS POST 2 Years	78, 81
SRVM005EXBMONEST1	EXALT B Monitor Elect Safety Test (EST)	78, 81
SRVM005EXBMONSTK1	EXALT B Monitor Elect Safety Test (STK)	78, 81
SpyGlass™ DS Service Codes		
SRVM005SPYDIG1001	SpyGlass Digital EssentialCare POS 1 Year	79
SRVM005SPYDIG1002	SpyGlass Digital EssentialCare POS 2 Year	79
SRVM005SPYDIG1003	SpyGlass Digital EssentialCare POS 3 Year	79
SRVM005SPYDIG1004	SpyGlass Digital EssentialCare POS 4 Year	79
SRVM005SPYDIG1005	SpyGlass Digital EssentialCare POST 1 Year	79
SRVM005SPYDIG1006	SpyGlass Digital EssentialCare POST 2 Year	79
SRVM005SPYDIG1007	SpyGlass Digital EssentialCare POST 3 Year	79
SRVM005SPYDIG1008	SpyGlass Digital EssentialCare POST 4 Year	79
SRVM005SPYDIG1009	SpyGlass Digital EverCare POS 1 Year	79
SRVM005SPYDIG1010	SpyGlass Digital EverCare POS 2 Year	79
SRVM005SPYDIG1011	SpyGlass Digital EverCare POS 3 Year	79
SRVM005SPYDIG1012	SpyGlass Digital EverCare POS 4 Year	79
SRVM005SPYDIG1013	SpyGlass Digital EverCare POST 1 Year	79
SRVM005SPYDIG1014	SpyGlass Digital EverCare POST 2 Year	79
SRVM005SPYDIG1015	SpyGlass Digital EverCare POST 3 Year	79
SRVM005SPYDIG1016	SpyGlass Digital EverCare POST 4 Year	79
SpyGlass™ Discover Service Codes		
SRVM005SPYDISEC01	SpyGlass Digital EssentialCare POS 1 Year	80
SRVM005SPYDISEC02	SpyGlass Digital EssentialCare POS 2 Year	80
SRVM005SPYDISEC03	SpyGlass Digital EssentialCare POS 3 Year	80
SRVM005SPYDISEC04	SpyGlass Digital EssentialCare POS 4 Year	80
SRVM005SPYDISEC05	SpyGlass Digital EssentialCare POST 1 Year	80




Product Code Index

Order Number	Description	Page Number
	SpyGlass™ Discover Service Codes (CONTINUED)	
SRVM005SPYDISEC06	SpyGlass Digital EssentialCare POST 2 Year	80
SRVM005SPYDISEC07	SpyGlass Digital EssentialCare POST 3 Year	80
SRVM005SPYDISEC08	SpyGlass Digital EssentialCare POST 4 Year	80
	TrapEase™ Polyp Traps	
STE-297-50	TrapEase Single Trap	97
STE-298-25	TrapEase Quad Trap	97
	Orca™	
SUV-617-50	Orca Sterile, Single Use Air/Water and Suction Valves	105
	OrcaPod™	
SUV-628-50	OrcaPod Sterile, Single Use Air/Water and Suction Valves with Seal Biopsy Valve	105
SUV-629-50	OrcaPod Sterile, Single Use Air/Water and Suction Valves, Seal Biopsy Valve and Hydra Water Jet Connector	105





CANNULAS




Tandem™ XL Triple-Lumen ERCP Cannula



Order Number	Description	Shaft OD (F) / (mm)	Distal OD (F) / (mm)	Length (cm)	Tip	Recommended Guidewire (in) / (mm)
M00535700	Tandem XL Triple-Lumen ERCP Cannula	7.0-5.5 / 2.33-1.83	5.5 / 1.83	200	Tapered	.035 / 0.89

Packaged one per box




Contour™ ERCP Cannulas



Order Number	Description	Shaft OD (F) / (mm)	Length (cm)	Tip	Recommended Guidewire (in) / (mm)
M00530850	Contour ERCP Cannula	5 / 1.67	210	Standard	.035 / 0.89
M00530860	Contour ERCP Cannula	5 / 1.67	210	Tapered	.035 / 0.89
M00530870	Contour ERCP Cannula	5 / 1.67	210	Ultra Tapered	.025 / 0.64
M00530880	Contour ERCP Cannula	5 / 1.67	210	5-4-3	.018 / 0.46
M00530890	Contour ERCP Cannula	5 / 1.67	210	Ball Tip	.035 / 0.89

Packaged one per box

SPHINCTEROTOMES



TRUEtome™ Jag Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire



Order Number	Description	Cut Wire (mm)	Tip OD (F) / (mm)	Tip Length (mm)	Pre-loaded Guidewire (in / cm)	Units
M00583080	TRUEtome 39 Sphincterotome	20	3.9 / 1.30	5	.025 / 450 Jagwire Guidewire	box 1
M00583090	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.025 / 450 Jagwire Guidewire	box 1
M00583100	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.035 / 450 Jagwire Guidewire	box 1
M00583110	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.035 / Stiff 450 Jagwire Guidewire	box 1
M00583120	TRUEtome 39 Sphincterotome	30	3.9 / 1.30	5	.025 / 450 Jagwire Guidewire	box 1
M00583130	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.025 / 450 Jagwire Guidewire	box 1
M00583140	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.035 / 450 Jagwire Guidewire	box 1
M00583150	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.035 / Stiff 450 Jagwire Guidewire	box 1

TRUEtome™ Dreamwire™ Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire



Order Number	Description	Cut Wire (mm)	Tip OD (F) / (mm)	Tip Length (mm)	Pre-loaded Guidewire (in / cm)	Units
M00583190	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.035 / 450 Dreamwire Guidewire	box 1
M00583210	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.035 / 450 Dreamwire Guidewire	box 1

TRUEtome™ Hydra Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire




Order Number	Description	Cut Wire (mm)	Tip OD (F) / (mm)	Tip Length (mm)	Pre-loaded Guidewire (in / cm)	Units
M00583160	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.035 / 450 Hydra Jagwire Guidewire	box 1
M00583170	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.035 / 450 Hydra Jagwire Guidewire	box 1

TRUEtome™ Jagwire Revolution Cannulating Sphincterotomes with Pre-loaded Jagwire™ Revolution Guidewire




Order Number	Description	Cut Wire (mm)	Tip OD (F) / (mm)	Tip Length (mm)	Pre-loaded Guidewire (in / cm)	Units
M00583260	TRUEtome 39 Sphincterotome	20	3.9 / 1.30	5	.025 / 450 Jagwire Revolution Guidewire	box 1
M00583270	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.025 / 450 Jagwire Revolution Guidewire	box 1
M00583280	TRUEtome 39 Sphincterotome	30	3.9 / 1.30	5	.025 / 450 Jagwire Revolution Guidewire	box 1
M00583290	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.025 / 450 Jagwire Revolution Guidewire	box 1

SPHINCTEROTOMES (CONTINUED)

**TRUEtome™ Cannulating Sphincterotomes**

Order Number	Description	Cut Wire (mm)	Distal OD (F) / (mm)	Tip Length (mm)	Catheter Length (cm)	Units
M00584190	TRUEtome 39	20	3.9 / 1.30	5	200	box 1
M00584120	TRUEtome 39	30	3.9 / 1.30	5	200	box 1
M00584191	TRUEtome 39	20	3.9 / 1.30	5	200	box 10
M00584121	TRUEtome 39	30	3.9 / 1.30	5	200	box 10
M00584192	TRUEtome 39	20	3.9 / 1.30	5	200	box 20
M00584122	TRUEtome 39	30	3.9 / 1.30	5	200	box 20
M00584170	TRUEtome 44	20	4.4 / 1.47	5	200	box 1
M00584180	TRUEtome 44	30	4.4 / 1.47	5	200	box 1
M00584171	TRUEtome 44	20	4.4 / 1.47	5	200	box 10
M00584181	TRUEtome 44	30	4.4 / 1.47	5	200	box 10
M00584172	TRUEtome 44	20	4.4 / 1.47	5	200	box 20
M00584182	TRUEtome 44	30	4.4 / 1.47	5	200	box 20
M00584150	TRUEtome 49	20	4.9 / 1.63	5	200	box 1
M00584160	TRUEtome 49	30	4.9 / 1.63	5	200	box 1
M00584151	TRUEtome 49	20	4.9 / 1.63	5	200	box 10
M00584161	TRUEtome 49	30	4.9 / 1.63	5	200	box 10
M00584152	TRUEtome 49	20	4.9 / 1.63	5	200	box 20
M00584162	TRUEtome 49	30	4.9 / 1.63	5	200	box 20

Recommended Guidewire: .025in (0.64 mm) for TRUEtome 39; .035in (0.89 mm) for TRUEtome 44 and 49


**Microknife™ XL Triple-Lumen Needle-Knife**

Order Number	Description	Shaft OD (F) / (mm)	Catheter Length (cm)
M00532810	Microknife XL Triple-Lumen Needle-Knife	7-5.5 / 2.33-1.83	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

SPHINCTEROTOMES (CONTINUED)



Ultratome XL™ Triple-Lumen Sphincterotomes




Order Number	Description	Cut Wire (mm)	Distal OD (F) / (mm)	Tip Length (mm)	Catheter Length (cm)
M00535900	Ultratome XL Short Nose	20	5.5 / 1.83	5	200
M00535910	Ultratome XL Long Nose	20	5.5 / 1.83	20	200
M00535920	Ultratome XL Short Nose	30	5.5 / 1.83	5	200
M00535930	Ultratome XL Long Nose	30	5.5 / 1.83	20	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

STONE RETRIEVAL




Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Distal Injection



Order Number	Description	Balloon OD (mm)	Injection	Catheter OD (F) / (mm)	System Length (cm)
M00547100	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Above	7-6 / 2.33-2.00	200
M00547110	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Above	7-6 / 2.33-2.00	200
M00547120	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Above	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Extractor™ Pro XL Triple-Lumen Retrieval Balloon Catheters – Proximal Injection



Order Number	Description	Balloon OD (mm)	Injection	Catheter OD (F) / (mm)	System Length (cm)
M00547130	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Below	7-6 / 2.33-2.00	200
M00547140	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Below	7-6 / 2.33-2.00	200
M00547150	Extractor Pro XL Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Below	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

ACCESSORIES


Biliary Accessories

Order Number	Description	Units
M00561250	Active Cord, Fits Olympus Brand Generator	box 1
M00561270	Active Cord, Fits Endostat™, Endostat™ II, Endostat™ III, Bovie®, Valleylab®, Aspen Labs and ERBE® Brand Generators	box 1
M00566670	⊗ Encore 26 Inflation Device	box 1
M00566001	⊗ Rapid Refill Continuous Injection System	box 10

Biliary Plastic Stents are located in the Stents Section on page 55

Biliary Guidewires are located in the Guidewires Section on pages 67 – 69



RX Biliary System™



Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics



CANNULAS





Rapid Exchange XL Cannula



Order Number	Description	Shaft OD (F) / (mm)	Length (cm)	Recommended Guidewire (in) / (mm)
M00545820	Rapid Exchange XL Cannula	7-5.5 / 2.33-1.84	195	.035 / 0.89

Packaged one per box



RX ERCP Cannulas



Order Number	Description	Shaft OD (F) / (mm)	Length (cm)	Recommended Guidewire (in) / (mm)
M00545760	RX ERCP Cannula Standard Tip	8.5-5 / 2.83-1.67	210	.035 / 0.89
M00545770	RX ERCP Cannula Tapered Tip	8.5-5 / 2.83-1.67	210	.035 / 0.89
M00545780	RX ERCP Cannula Ball Tip	8.5-5 / 2.83-1.67	210	.035 / 0.89

Packaged one per box

SPHINCTEROTOMES



Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire



Order Number	Description	Cut Wire (mm)	Shaft OD (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	Catheter Length (cm)	PRE-LOADED GUIDEWIRE** Length (cm)	OD (in) / (mm)
M00583000	Hydratome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	260	.035 / 0.89
M00583010	Hydratome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	260	.035 / 0.89
M00583020	Hydratome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	450	.035 / 0.89
M00583030	Hydratome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	450	.035 / 0.89
M00583040*	Hydratome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	260	.035 / 0.89
M00583050*	Hydratome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	260	.035 / 0.89
M00583060*	Hydratome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	450	.035 / 0.89
M00583070*	Hydratome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	450	.035 / 0.89

Packaged one per box

* Includes handle brake

** Pre-loaded with .035in (0.89mm) Hydra Jagwire™ Guidewire



Dreamtome™ RX Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire



Order Number	Description	Cut Wire (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	PRE-LOADED GUIDEWIRE** Length (cm)	OD (in) / (mm)
M00584000	Dreamtome RX 49	20	5	4.9 / 1.63	260	.035 / 0.89
M00584010	Dreamtome RX 49	30	5	4.9 / 1.63	260	.035 / 0.89
M00584020	Dreamtome RX 49	20	5	4.9 / 1.63	450	.035 / 0.89
M00584030	Dreamtome RX 49	30	5	4.9 / 1.63	450	.035 / 0.89
M00584040*	Dreamtome RX 44	20	5	4.4 / 1.47	260	.035 / 0.89
M00584050*	Dreamtome RX 44	30	5	4.4 / 1.47	260	.035 / 0.89
M00584060*	Dreamtome RX 44	20	5	4.4 / 1.47	450	.035 / 0.89
M00584070*	Dreamtome RX 44	30	5	4.4 / 1.47	450	.035 / 0.89

Packaged one per box

* Includes handle brake

** Pre-loaded with Straight-tip guidewires only

SPHINCTEROTOMES (CONTINUED)




Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire

Order Number	Description	Cut Wire (mm)	Shaft OD (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	Catheter Length (cm)	PRE-LOADED GUIDEWIRE™		
							Shaft Stiffness	Length (cm)	OD (in) / (mm)
M00573040*	Jagtome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	260	.035 / 0.89
M00573060*	Jagtome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	260	.035 / 0.89
M00573050*	Jagtome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	450	.035 / 0.89
M00573070*	Jagtome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	450	.035 / 0.89
M00573080*	Jagtome RX 39	20	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Stiff	260	.025 / 0.64
M00573100*	Jagtome RX 39	30	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Stiff	260	.025 / 0.64
M00573090*	Jagtome RX 39	20	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Standard	450	.025 / 0.64
M00573110*	Jagtome RX 39	30	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Standard	450	.025 / 0.64

Packaged one per box

* Includes handle brake


** Pre-loaded with the high performance Jagwire™ Guidewire (Jagtome RX 49 and RX 44 contain the Stiff Shaft Jagwire and Jagtome RX 39 the Standard Jagwire)



Jagtome™ Revolution RX Cannulating Sphincterotome (Autotome™ Preloaded with Jagwire Revolution High Performance Guidewire)

Order Number	Description	Cut Wire (mm)	Tip Length (mm)	Distal Tip OD (F)	Preloaded Guidewire (in) / (cm)	Guidewire	Tip Style
M00584240	Jagtome Revolution RX 39 Sphincterotome	20	5	3.9	.025 / 260	Jagwire Revolution Guidewire	Straight
M00584250	Jagtome Revolution RX 39 Sphincterotome	20	5	3.9	.025 / 450	Jagwire Revolution Guidewire	Straight
M00584260	Jagtome Revolution RX 39 Sphincterotome	30	5	3.9	.025 / 260	Jagwire Revolution Guidewire	Straight
M00584270	Jagtome Revolution RX 39 Sphincterotome	30	5	3.9	.025 / 450	Jagwire Revolution Guidewire	Straight

SPHINCTEROTOMES (CONTINUED)




Autotome™ Pro RX Cannulating Sphincterotomes

Order Number	Description	Cut Wire (mm)	Tip Length (mm)	Distal Tip OD (Fr)	Recommended Guidewire (in / mm)	Units
M00584380	Autotome Pro RX 44	25	5	4.4	.025 / .035	box 1
M00584390	Autotome Pro RX 39	25	5	3.9	.025	box 1

Autotome™ Pro RX Preloaded Sphincterotomes

Order Number	Description	Cut Wire (mm)	Tip Length (mm)	Distal Tip OD (Fr)	Preloaded Guidewire Type	Units
M00584400	Dreamtome Pro RX 44	25	5	4.4	0.035" Dreamwire	box 1
M00584410	Dreamtome Pro RX 44	25	5	4.4	0.035" Dreamwire	box 1
M00584420	Hydratome Pro RX 44	25	5	4.4	0.035" Hydra Jagwire	box 1
M00584430	Hydratome Pro RX 44	25	5	4.4	0.035" Hydra Jagwire	box 1
M00584440	Jagtome Pro RX 44	25	5	4.4	0.035" Jagwire	box 1
M00584450	Jagtome Pro RX 44	25	5	4.4	0.035" Jagwire	box 1
M00584460	Jagtome Pro RX 39	25	5	3.9	0.025" Jagwire	box 1
M00584470	Jagtome Revolution Pro RX 39	25	5	3.9	0.025" Jagtome Revolution Guidewire	box 1
M00584480	Jagtome Revolution Pro RX 39	25	5	3.9	0.025" Jagtome Revolution Guidewire	box 1



Autotome™ RX Cannulating Sphincterotomes

Order Number	Description	Cut Wire (mm)	Shaft OD (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	Catheter Length (cm)	Recommended Guidewire (in) / (mm)
M00545150	Autotome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	.035 / 0.89
M00545160	Autotome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	.035 / 0.89
M00545170*	Autotome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	.035 / 0.89
M00545180*	Autotome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	.035 / 0.89
M00545190*	Autotome RX 39	20	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	.025 / 0.64
M00545200*	Autotome RX 39	30	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	.025 / 0.64

Packaged one per box

Recommended Guidewires: For Autotome RX 44 and 49, recommended guidewire is .035in (0.89mm) Dreamwire™ High Performance Guidewire
For Autotome RX 39, recommended guidewire is .025in (0.69mm) Jagwire™ Guidewire

* Includes handle brake

SPHINCTEROTOMES (CONTINUED)

RX Needle-Knife XL Triple-Lumen Needle-Knife




Order Number	Description	Cut Wire (mm)	Shaft OD (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)
M00545840	RX Needle-Knife	5	7-5.5 / 2.3-1.83	5	5.5 / 1.8

Packaged one per box

Recommended Guidewire: .035 in (0.89mm) Dreamwire™ High Performance Guidewire

STONE RETRIEVAL




Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Distal Injection



Order Number	Description	Balloon OD (mm)	Injection	Catheter OD (F) / (mm)	System Length (cm)
M00547000	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Above	7-6 / 2.33-2.00	200
M00547010	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Above	7-6 / 2.33-2.00	200
M00547020	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Above	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Extractor™ Pro RX Triple-Lumen Retrieval Balloon Catheters – Proximal Injection




Order Number	Description	Balloon OD (mm)	Injection	Catheter OD (F) / (mm)	System Length (cm)
M00547030	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Below	7-6 / 2.33-2.00	200
M00547040	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Below	7-6 / 2.33-2.00	200
M00547050	Extractor Pro RX Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Below	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

STONE RETRIEVAL (CONTINUED)




Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Distal Injection

Order Number	Description	Balloon OD (mm)	Injection	Catheter OD (F) / (mm)	System Length (cm)
M00547300	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Above	7-6 / 2.33-2.00	200
M00547310	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Above	7-6 / 2.33-2.00	200
M00547320	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Above	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire




Extractor™ Pro RX-S Triple-Lumen Retrieval Balloon Catheters – Proximal Injection

Order Number	Description	Balloon OD (mm)	Injection	Catheter OD (F) / (mm)	System Length (cm)
M00547330	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	9.0-12.0	Below	7-6 / 2.33-2.00	200
M00547340	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	12.0-15.0	Below	7-6 / 2.33-2.00	200
M00547350	Extractor Pro RX-S Triple-Lumen Retrieval Balloon Catheter	15.0-18.0	Below	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Trapezoid™ RX Wireguided Retrieval Baskets

Order Number	Description	Open Basket Diameter (cm)	Required Working Channel (mm)
M00510860	Trapezoid RX Wireguided Retrieval Basket	1.5	3.2
M00510870	Trapezoid RX Wireguided Retrieval Basket	2.0	3.2
M00510880	Trapezoid RX Wireguided Retrieval Basket	2.5	3.2
M00510890*	Trapezoid RX Wireguided Retrieval Basket	3.0	3.2
M00550620	Alliance™ II Inflation Handle**		


Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

* The Trapezoid 3 x 6cm basket is designed for crushing calculus larger than 1.5cm in diameter

** Alliance™ II recommended for use as Mechanical Lithotripter. Alliance II is reusable

DILATION



Hurricane™ RX Biliary Dilation Balloons



Order Number	Description	Balloon OD (mm)	Balloon Length (cm)	Catheter OD (F) / (mm)	Maximum Inflation Pressure (ATM) / (kPa)	Minimum Channel Size (mm)	Recommended Guidewire (in) / (mm)
M00545890	Hurricane RX Dilation Balloon	4	2	5.8 / 1.93	12 / 1216	3.2	.035 / 0.89
M00545900	Hurricane RX Dilation Balloon	4	4	5.8 / 1.93	12 / 1216	3.2	.035 / 0.89
M00545910	Hurricane RX Dilation Balloon	6	2	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545920	Hurricane RX Dilation Balloon	6	4	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545930	Hurricane RX Dilation Balloon	8	2	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545940	Hurricane RX Dilation Balloon	8	4	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545950	Hurricane RX Dilation Balloon	10	2	5.8 / 1.93	8 / 810	3.2	.035 / 0.89
M00545960	Hurricane RX Dilation Balloon	10	4	5.8 / 1.93	8 / 810	3.2	.035 / 0.89

Packaged one per box


Recommended Inflation Device: Encore™ 26 Inflation Device, M00566670

CRE™ RX Biliary Balloon Dilatation Catheter



Order Number	INFLATED BALLOON OD (mm) (Fr)		Balloon Dilation Length (cm)	INFLATION PRESSURE (ATM) (kPa)		Catheter Size (Fr)	Minimum Working Channel (mm)	Working Length (cm)
M00558920	8-9-10	24-27-30	3.0	3-5.5-9	304-608-1013	7.5	3.7	180
M00558930	10-11-12	30-33-36	3.0	3-5-8	304-557-912	7.5	3.7	180
M00558940	12-13.5-15	36-40.5-45	3.0	3-4.5-8	304-507-811	7.5	3.7	180
M00558860	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	7.5	3.7	180
M00558870	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	7.5	3.7	180
M00558880	10-11-12	30-33-36	5.5	3-5-8	304-507-811	7.5	3.7	180
M00558890	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	7.5	3.7	180
M00558900	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	7.5	3.7	180
M00558910	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	7.5	3.7	180

CYTOLOGY



RX Cytology Brush



Order Number	Description	Brush Width (mm)	Catheter OD (F) / (mm)	Working Length (cm)	Recommended Guidewire (in) / (mm)
M00545000	RX Cytology Brush	2.1	8 / 2.67	200	.035 / 0.89

Packaged one per box

LOCKING DEVICES

AutoCap™ RX




Order Number	Description	Units
M00545101	AutoCap RX Olympus Compatible (TJF-160, TJF-180, TJF-190 and JF-260V Series Duodenoscopes)	box 10
M00545131	AutoCap RX EXALT™ Model D Compatible Integrated Biopsy	box 10

RX Locking Devices



Order Number	Description	Units
M00545211†	RX Locking Device – Olympus / FujiFilm Brand Endoscope Compatible	box 10
M00545261	RX Locking Device and Biopsy Cap System – Olympus / FujiFilm G5 Series Brand Endoscope Compatible	box 10
M00545271	RX Locking Device and Biopsy Cap System – Pentax Brand Endoscope Compatible	box 10
M00545281	RX Locking Device and Biopsy Cap System – FujiFilm Type 3 Series Brand Endoscope Compatible	box 10
M00566001	Rapid Refill Continuous Injection System	box 10

† Product not available in the European Economic Area (EEA), except in Greece. Not for use or distribution in the EEA, except in Greece.
Please check availability of the product with your local sales representative or your local customer service.



Encore™ 26 Inflation Device




Order Number	Description	Units
M00566670	Encore 26 Inflation Device	box 1

RX Biliary Plastic and Metal Stents are located in the Stents Section on pages 51 – 54 and 61 – 62

Biliary Guidewires are located in the Guidewires Section on pages 67 – 69

Biliary Lumen Patency



WALLFLEX™ BILIARY RX
APPROVED FOR
BENIGN
INDICATION
FULLY COVERED STENT

Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics




Contents



Alphabetical Index Product Code Index

Biliary Lumen Patency

RX BILIARY PLASTIC STENTS



Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System



Order Number	Description	Shape	OD (F) / (mm)	Length (cm) Barb to Barb
M00534200	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	5
M00534210	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	7
M00534220	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	9
M00534230	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	12
M00534240	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	15
M00534250	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	18
M00534260	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	5
M00534270	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	7
M00534280	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	9
M00534290	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	12
M00534300	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	15
M00534310	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	18
M00534320	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	5
M00534330	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	7
M00534340	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	9
M00534350	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	12
M00534360	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	15
M00534370	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	18
M00534560	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	5
M00534570	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	7
M00534580	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	9
M00534590	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	12
M00534600	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	15
M00534610	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	18
M00534620	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	5
M00534630	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	7
M00534640	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	9
M00534650	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	12
M00534660	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	15
M00534670	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	18
M00534680	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	5
M00534690	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	7
M00534700	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	9
M00534710	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	12
M00534720	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	15
M00534730	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	18

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

Scope Compatibility: 7F (2.3mm) and 8.5F (2.83mm) stents require a 3.2mm working channel; 10F (3.3mm) stents require a 4.2mm working channel



Contents

Alphabetical
Index

Product Code
Index

Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics

Biliary Lumen Patency

RX BILIARY PLASTIC STENTS (CONTINUED)



Advanix™ Stents – Singles



Order Number	Description	Shape	OD (F) / (mm)	Length (cm) Barb to Barb
M00532870	Advanix Stent	Duodenal Bend	7 / 2.3	5
M00532880	Advanix Stent	Duodenal Bend	7 / 2.3	7
M00532890	Advanix Stent	Duodenal Bend	7 / 2.3	9
M00532900	Advanix Stent	Duodenal Bend	7 / 2.3	12
M00532910	Advanix Stent	Duodenal Bend	7 / 2.3	15
M00532920	Advanix Stent	Duodenal Bend	7 / 2.3	18
M00532930	Advanix Stent	Duodenal Bend	8.5 / 2.83	5
M00532940	Advanix Stent	Duodenal Bend	8.5 / 2.83	7
M00532950	Advanix Stent	Duodenal Bend	8.5 / 2.83	9
M00532960	Advanix Stent	Duodenal Bend	8.5 / 2.83	12
M00532970	Advanix Stent	Duodenal Bend	8.5 / 2.83	15
M00532980	Advanix Stent	Duodenal Bend	8.5 / 2.83	18
M00532990	Advanix Stent	Duodenal Bend	10 / 3.3	5
M00533000	Advanix Stent	Duodenal Bend	10 / 3.3	7
M00533010	Advanix Stent	Duodenal Bend	10 / 3.3	9
M00533020	Advanix Stent	Duodenal Bend	10 / 3.3	12
M00533030	Advanix Stent	Duodenal Bend	10 / 3.3	15
M00533040	Advanix Stent	Duodenal Bend	10 / 3.3	18
M00533230	Advanix Stent	Center Bend	7 / 2.3	5
M00533240	Advanix Stent	Center Bend	7 / 2.3	7
M00533250	Advanix Stent	Center Bend	7 / 2.3	9
M00533260	Advanix Stent	Center Bend	7 / 2.3	12
M00533270	Advanix Stent	Center Bend	7 / 2.3	15
M00533280	Advanix Stent	Center Bend	7 / 2.3	18
M00533290	Advanix Stent	Center Bend	8.5 / 2.83	5
M00533300	Advanix Stent	Center Bend	8.5 / 2.83	7
M00533310	Advanix Stent	Center Bend	8.5 / 2.83	9
M00533320	Advanix Stent	Center Bend	8.5 / 2.83	12
M00533330	Advanix Stent	Center Bend	8.5 / 2.83	15
M00533340	Advanix Stent	Center Bend	8.5 / 2.83	18
M00533350	Advanix Stent	Center Bend	10 / 3.3	5
M00533360	Advanix Stent	Center Bend	10 / 3.3	7
M00533370	Advanix Stent	Center Bend	10 / 3.3	9
M00533380	Advanix Stent	Center Bend	10 / 3.3	12
M00533390	Advanix Stent	Center Bend	10 / 3.3	15
M00533400	Advanix Stent	Center Bend	10 / 3.3	18

Packaged one per box

Recommended Guidewire: .035in (0.89mm)

Scope Compatibility: 7F (2.3mm) and 8.5F (2.83mm) stents require a 3.2mm working channel; 10F (3.3mm) stents require a 4.2mm working channel




Contents

Alphabetical
Index

Product Code
Index

Biliary Lumen Patency

RX BILIARY PLASTIC STENTS (CONTINUED)



Advanix™ Stents – Singles (CONTINUED)



Order Number	Description	Shape	OD (F) / (mm)	Length (cm) Between Pigtails
M00532160	Advanix Stent	Double Pigtail	7 / 2.3	3
M00532180	Advanix Stent	Double Pigtail	7 / 2.3	5
M00532190	Advanix Stent	Double Pigtail	7 / 2.3	7
M00532210	Advanix Stent	Double Pigtail	7 / 2.3	10
M00532220	Advanix Stent	Double Pigtail	7 / 2.3	12
M00532230	Advanix Stent	Double Pigtail	7 / 2.3	15
M00532240	Advanix Stent	Double Pigtail	10 / 3.3	3
M00532260	Advanix Stent	Double Pigtail	10 / 3.3	5
M00532270	Advanix Stent	Double Pigtail	10 / 3.3	7
M00532290	Advanix Stent	Double Pigtail	10 / 3.3	10
M00532300	Advanix Stent	Double Pigtail	10 / 3.3	12
M00532310	Advanix Stent	Double Pigtail	10 / 3.3	15

Packaged one per box

Recommended Guidewire: .035in (0.89mm)

Scope Compatibility: 7F (2.3mm) and 8.5F (2.83mm) stents require a 3.2mm working channel; 10F (3.3mm) stents require a 4.2mm working channel



NaviFlex™ RX Delivery System



Order Number	Description	Catheter Diameter (mm)	Working Length (cm)
M00533560	RX Delivery System	7 / 2.3	202.5
M00533570	RX Delivery System	8.5 / 2.83	202.5
M00533580	RX Delivery System	10 / 3.3	202.5

Packaged one per box

Compatible with Advanix Biliary Stents.



Contents

Alphabetical
Index

Product Code
Index

Biliary Lumen Patency

RX BILIARY PLASTIC STENTS (CONTINUED)

RX Biliary Plastic Stents with RX Delivery System



Order Number	Description	OD (F) / (mm)	Length (cm) Barb to Barb
M00545550	RX Biliary Plastic Stent with RX Delivery System	7 / 2.33	5
M00545560	RX Biliary Plastic Stent with RX Delivery System	7 / 2.33	7
M00545570	RX Biliary Plastic Stent with RX Delivery System	7 / 2.33	10
M00545580	RX Biliary Plastic Stent with RX Delivery System	7 / 2.33	12
M00545590	RX Biliary Plastic Stent with RX Delivery System	7 / 2.33	15
M00545650	RX Biliary Plastic Stent with RX Delivery System	8.5 / 2.83	5
M00545660	RX Biliary Plastic Stent with RX Delivery System	8.5 / 2.83	7
M00545670	RX Biliary Plastic Stent with RX Delivery System	8.5 / 2.83	10
M00545680	RX Biliary Plastic Stent with RX Delivery System	8.5 / 2.83	12
M00545690	RX Biliary Plastic Stent with RX Delivery System	8.5 / 2.83	15
M00545600	RX Biliary Plastic Stent with RX Delivery System	10 / 3.33	5
M00545610	RX Biliary Plastic Stent with RX Delivery System	10 / 3.33	7
M00545620	RX Biliary Plastic Stent with RX Delivery System	10 / 3.33	10
M00545630	RX Biliary Plastic Stent with RX Delivery System	10 / 3.33	12
M00545640	RX Biliary Plastic Stent with RX Delivery System	10 / 3.33	15
M00545700	RX Biliary Plastic Stent with RX Delivery System	11.5 / 3.83	5
M00545710	RX Biliary Plastic Stent with RX Delivery System	11.5 / 3.83	7
M00545720	RX Biliary Plastic Stent with RX Delivery System	11.5 / 3.83	10
M00545730	RX Biliary Plastic Stent with RX Delivery System	11.5 / 3.83	12
M00545740	RX Biliary Plastic Stent with RX Delivery System	11.5 / 3.83	15

Packaged one per box

Recommended Guidewire: .035in (0.89mm)



Contents

Alphabetical Index

Product Code Index

Biliary Lumen Patency

BILIARY PLASTIC STENTS



Flexima™ Biliary Stents with Delivery System



Order Number	Description	OD (F) / (mm)	Length (cm) Barb to Barb
M00539200	Flexima Stent with Delivery System	7 / 2.33	5
M00539210	Flexima Stent with Delivery System	7 / 2.33	7
M00539220	Flexima Stent with Delivery System	7 / 2.33	10
M00539230	Flexima Stent with Delivery System	7 / 2.33	12
M00539240	Flexima Stent with Delivery System	7 / 2.33	15
M00539250	Flexima Stent with Delivery System	8.5 / 2.83	5
M00539260	Flexima Stent with Delivery System	8.5 / 2.83	7
M00539270	Flexima Stent with Delivery System	8.5 / 2.83	10
M00539280	Flexima Stent with Delivery System	8.5 / 2.83	12
M00539290	Flexima Stent with Delivery System	8.5 / 2.83	15
M00539300	Flexima Stent with Delivery System	10 / 3.33	5
M00539310	Flexima Stent with Delivery System	10 / 3.33	7
M00539320	Flexima Stent with Delivery System	10 / 3.33	10
M00539330	Flexima Stent with Delivery System	10 / 3.33	12
M00539340	Flexima Stent with Delivery System	10 / 3.33	15
M00539350	Flexima Stent with Delivery System	11.5 / 3.83	5
M00539360	Flexima Stent with Delivery System	11.5 / 3.83	7
M00539370	Flexima Stent with Delivery System	11.5 / 3.83	10
M00539380	Flexima Stent with Delivery System	11.5 / 3.83	12
M00539390	Flexima Stent with Delivery System	11.5 / 3.83	15

Packaged one per box

Recommended Guidewire: .035in (0.89mm)



Contents

Alphabetical Index

Product Code Index

PANCREATIC PLASTIC STENTS



Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Pigtail



Order Number	Description	Leading Barb	Shape	OD (F) / (mm)	Length (cm)
M00537350	Advanix Pancreatic Stent Kit ●	No	Pigtail	4 / 1.35	3
M00537360	Advanix Pancreatic Stent Kit ▲	No	Pigtail	4 / 1.35	4
M00537370	Advanix Pancreatic Stent Kit ▲	No	Pigtail	4 / 1.35	5
M00537380	Advanix Pancreatic Stent Kit ▲	No	Pigtail	4 / 1.35	7
M00537390	Advanix Pancreatic Stent Kit ●	No	Pigtail	5 / 1.67	3
M00537400	Advanix Pancreatic Stent Kit ▲	No	Pigtail	5 / 1.67	4
M00537410	Advanix Pancreatic Stent Kit ▲	No	Pigtail	5 / 1.67	5
M00537420	Advanix Pancreatic Stent Kit ▲	No	Pigtail	5 / 1.67	7
M00537430	Advanix Pancreatic Stent Kit ●	Yes	Pigtail	5 / 1.67	3
M00537440	Advanix Pancreatic Stent Kit ▲	Yes	Pigtail	5 / 1.67	4
M00537450	Advanix Pancreatic Stent Kit ▲	Yes	Pigtail	5 / 1.67	5

* Products with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Straight



Order Number	Description	Leading Barb	Shape	OD (F) / (mm)	Length (cm)
M00537460	Advanix Pancreatic Stent Kit ● ¹	No	Straight	4 / 1.35	3
M00537470	Advanix Pancreatic Stent Kit ▲ ¹	No	Straight	4 / 1.35	4
M00537480	Advanix Pancreatic Stent Kit ▲ ¹	No	Straight	4 / 1.35	5
M00537490	Advanix Pancreatic Stent Kit ▲ ¹	No	Straight	4 / 1.35	7
M00537500	Advanix Pancreatic Stent Kit ●	No	Straight	5 / 1.67	3
M00537510	Advanix Pancreatic Stent Kit ▲	No	Straight	5 / 1.67	4
M00537520	Advanix Pancreatic Stent Kit ▲	No	Straight	5 / 1.67	5
M00537530	Advanix Pancreatic Stent Kit ▲	No	Straight	5 / 1.67	7
M00537540	Advanix Pancreatic Stent Kit ●	Yes	Straight	5 / 1.67	3
M00537550	Advanix Pancreatic Stent Kit ▲	Yes	Straight	5 / 1.67	4
M00537560	Advanix Pancreatic Stent Kit ▲	Yes	Straight	5 / 1.67	5

* Products with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker

¹ Stiffer stent material compared to 4F pigtail stents.

² Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.



PANCREATIC PLASTIC STENTS (CONTINUED)



Advanix™ Pancreatic Single Stents – Pigtail

Order Number	Description	Leading Barb	Shape	OD (F) / (mm)	Length (cm)
M00536740	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	2
M00536750	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	3
M00536760	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	4
M00536770	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	5
M00536780	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	6
M00536790	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	7
M00536800	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	8
M00536810	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	9
M00536820	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	10
M00536830	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	11
M00536840	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	12
M00536850	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	13
M00536860	Advanix Pancreatic Stent ●	No	Pigtail	4 / 1.35	2
M00536870	Advanix Pancreatic Stent ●	No	Pigtail	4 / 1.35	3
M00536880	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	4
M00536890	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	5
M00536900	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	6
M00536910	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	7
M00536920	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	8
M00536930	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	9
M00536940	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	10
M00536950	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	11
M00536960	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	12
M00536970	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	13
M00536980	Advanix Pancreatic Stent ●	No	Pigtail	5 / 1.67	2
M00536990	Advanix Pancreatic Stent ●	No	Pigtail	5 / 1.67	3
M00537000	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	4
M00537010	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	5
M00537020	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	6
M00537030	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	7
M00537040	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	8
M00537050	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	9
M00537060	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	10
M00537070	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	11
M00537080	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	12
M00537090	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	13

* Products with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker

1 Stiffer stent material compared to 4F pigtail stents.

2 Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.



PANCREATIC PLASTIC STENTS (CONTINUED)



Advanix™ Pancreatic Single Stents – Pigtail (CONTINUED)



Order Number	Description	Leading Barb	Shape	OD (F) / (mm)	Length (cm)
M00537100	Advanix Pancreatic Stent ●	Yes	Pigtail	5 / 1.67	2
M00537110	Advanix Pancreatic Stent ●	Yes	Pigtail	5 / 1.67	3
M00537120	Advanix Pancreatic Stent ▲	Yes	Pigtail	5 / 1.67	4
M00537130	Advanix Pancreatic Stent ▲	Yes	Pigtail	5 / 1.67	5
M00537140	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	6
M00537150	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	7
M00537160	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	8
M00537170	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	9
M00537180	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	10
M00537190	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	11
M00537200	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	12
M00537210	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	13
M00537220	Advanix Pancreatic Stent ■	Yes	Pigtail	7 / 2.33	3
M00537230	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	4
M00537240	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	5
M00537250	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	6
M00537260	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	7
M00537270	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	8
M00537280	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	9
M00537290	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	10
M00537300	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	11
M00537310	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	12
M00537320	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	13
M00537330	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	15
M00537340	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	18

* Products with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker

¹ Stiffer stent material compared to 4F pigtail stents.

² Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.



PANCREATIC PLASTIC STENTS (CONTINUED)



Advanix™ Pancreatic Single Stents – Straight



Order Number	Description	Leading Barb	Shape	OD (F) / (mm)	Length (cm)
M00536000	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	2
M00536010	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	3
M00536020	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	4
M00536030	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	5
M00536040	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	6
M00536050	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	7
M00536060	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	8
M00536070	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	9
M00536080	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	10
M00536090	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	11
M00536100	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	12
M00536110	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	13
M00536120	Advanix Pancreatic Stent ● ¹	No	Straight	4 / 1.35	2
M00536130	Advanix Pancreatic Stent ● ¹	No	Straight	4 / 1.35	3
M00536140	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	4
M00536150	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	5
M00536160	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	6
M00536170	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	7
M00536180	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	8
M00536190	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	9
M00536200	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	10
M00536210	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	11
M00536220	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	12
M00536230	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	13
M00536240	Advanix Pancreatic Stent ●	No	Straight	5 / 1.67	2
M00536250	Advanix Pancreatic Stent ●	No	Straight	5 / 1.67	3
M00536260	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	4
M00536270	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	5
M00536280	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	6
M00536290	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	7
M00536300	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	8
M00536310	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	9
M00536320	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	10
M00536330	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	11
M00536340	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	12
M00536350	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	13

* Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker

¹ Stiffer stent material compared to 4F pigtail stents.

² Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.



PANCREATIC PLASTIC STENTS (CONTINUED)



Advanix™ Pancreatic Single Stents – Straight (CONTINUED)



Order Number	Description	Leading Barb	Shape	OD (F) / (mm)	Length (cm)
M00536360	Advanix Pancreatic Stent ●	Yes	Straight	5 / 1.67	2
M00536370	Advanix Pancreatic Stent ●	Yes	Straight	5 / 1.67	3
M00536380	Advanix Pancreatic Stent ▲	Yes	Straight	5 / 1.67	4
M00536390	Advanix Pancreatic Stent ▲	Yes	Straight	5 / 1.67	5
M00536400	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	6
M00536410	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	7
M00536420	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	8
M00536430	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	9
M00536440	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	10
M00536450	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	11
M00536460	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	12
M00536470	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	13
M00536480	Advanix Pancreatic Stent ■	Yes	Straight	7 / 2.33	3
M00536490	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	4
M00536500	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	5
M00536510	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	6
M00536520	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	7
M00536530	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	8
M00536540	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	9
M00536550	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	10
M00536560	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	11
M00536570	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	12
M00536580	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	13
M00536590	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	15
M00536600	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	18
M00536610	Advanix Pancreatic Stent ■	Yes	Straight	10 / 3.33	3
M00536620	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	4
M00536630	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	5
M00536640	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	6
M00536650	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	7
M00536660	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	8
M00536670	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	9
M00536680	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	10
M00536690	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	11
M00536700	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	12
M00536710	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	13
M00536720	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	15
M00536730	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	18

* Products with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker


¹ Stiffer stent material compared to 4F pigtail stents.

² Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.



Biliary Lumen Patency

PANCREATIC PLASTIC STENTS (CONTINUED)



NaviFlex™ RX Pancreatic Delivery System & NaviFlex™ Pushers



Order Number	Description	Diameter (F) / (mm)
M00535000	Pancreatic Stent RX Pusher	3 / 1.00
M00535010	Pancreatic Stent RX Pusher	4 & 5 / 1.35 & 1.67
M00535020	Pancreatic Stent Pusher	3 / 1.00
M00535030	Pancreatic Stent Pusher	4 & 5 / 1.35 & 1.67
M00535060	NaviFlex RX Pancreatic Stent Delivery System	7 / 2.33
M00535080	NaviFlex RX Pancreatic Stent Delivery System	10 / 3.33

RX BILIARY METAL STENTS





WallFlex™ Biliary RX Fully Covered Stent System RMV (Permalume™ Covering, includes benign indication)



Order Number	Description	Diameter (mm)	Length (mm)	Catheter Diameter (F) / (mm)	Guidewire Diameter (in) / (mm)
M00574210†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	6	40	8.5 / 2.83	.035 / 0.89
M00574230†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	6	60	8.5 / 2.83	.035 / 0.89
M00574240†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	6	80	8.5 / 2.83	.035 / 0.89
M00570450†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	8	60	8.5 / 2.83	.035 / 0.89
M00570460†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	8	80	8.5 / 2.83	.035 / 0.89
M00576800†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	8	100	9.0 / 3.00	.035 / 0.89
M00576810†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	8	120	9.0 / 3.00	.035 / 0.89
M00570470†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	10	40	8.5 / 2.83	.035 / 0.89
M00570480†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	10	60	8.5 / 2.83	.035 / 0.89
M00570490†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	10	80	8.5 / 2.83	.035 / 0.89
M00576820†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	10	100	9.0 / 3.00	.035 / 0.89
M00576830†	WallFlex RX Biliary (Permalume Covering, includes benign indication)	10	120	9.0 / 3.00	.035 / 0.89

† Product commercialised only in countries where benign indication is approved.

RX BILIARY METAL STENTS (CONTINUED)



WallFlex™ Biliary RX Partially Covered Stent System (Permalume™ Covering)



Order Number	Description	Diameter (mm)	Length (mm)	Covered Length (mm)	Catheter Diameter (F) / (mm)	Guidewire Diameter (in) / (mm)
M00570700	WallFlex RX Biliary (Permalume Covering)	8	60	48	8.5 / 2.83	.035 / 0.89
M00570710	WallFlex RX Biliary (Permalume Covering)	8	80	68	8.5 / 2.83	.035 / 0.89
M00576740	WallFlex RX Biliary (Permalume Covering)	8	100	88	9.0 / 3.00	.035 / 0.89
M00576750	WallFlex RX Biliary (Permalume Covering)	8	120	108	9.0 / 3.00	.035 / 0.89
M00570720	WallFlex RX Biliary (Permalume Covering)	10	40	28	8.5 / 2.83	.035 / 0.89
M00570730	WallFlex RX Biliary (Permalume Covering)	10	60	48	8.5 / 2.83	.035 / 0.89
M00570740	WallFlex RX Biliary (Permalume Covering)	10	80	68	8.5 / 2.83	.035 / 0.89
M00576760	WallFlex RX Biliary (Permalume Covering)	10	100	88	9.0 / 3.00	.035 / 0.89
M00576770	WallFlex RX Biliary (Permalume Covering)	10	120	108	9.0 / 3.00	.035 / 0.89

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm or .035in (0.89mm) Dreamwire™ High Performance Guidewire, 260cm



WallFlex™ Biliary RX Uncovered Stent System




Order Number	Description	Diameter (mm)	Length (mm)	Catheter Diameter (F) / (mm)	Guidewire Diameter (in) / (mm)
M00570600	WallFlex RX Biliary	8	40	8 / 2.67	.035 / 0.89
M00570610	WallFlex RX Biliary	8	60	8 / 2.67	.035 / 0.89
M00570620	WallFlex RX Biliary	8	80	8 / 2.67	.035 / 0.89
M00570630	WallFlex RX Biliary	8	100	8 / 2.67	.035 / 0.89
M00576780	WallFlex RX Biliary	8	120	8 / 2.67	.035 / 0.89
M00570890	WallFlex RX Biliary	10	40	8 / 2.67	.035 / 0.89
M00570640	WallFlex RX Biliary	10	60	8 / 2.67	.035 / 0.89
M00570650	WallFlex RX Biliary	10	80	8 / 2.67	.035 / 0.89
M00570660	WallFlex RX Biliary	10	100	8 / 2.67	.035 / 0.89
M00576790	WallFlex RX Biliary	10	120	8 / 2.67	.035 / 0.89

Packaged one per box


Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm or .035in (0.89mm) Dreamwire™ High Performance Guidewire, 260cm

RADIOFREQUENCY ABLATION CATHETER AND CABLES

Habib™ EndoHPB Bipolar Radiofrequency Catheter




Order Number	Description
M00500070	Habib Endo HPB Bipolar Radiofrequency Catheter



Order Number	Description
5100	Adapter Cable Bipolar, 28 mm pins



Order Number	Description
5700	Adapter Cable Bipolar, 22 mm pins



Order Number	Description
5420	Adapter Cable Bipolar, Internal 4 mm pin



Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics

FINE NEEDLE ASPIRATION (FNA)



Expect™ Slimline Handle Needles*



Order Number	Description	Needle Size (gauge) / (mm)	Minimum Working Channel (mm)	Sheath Diameter (mm)	Units
M00555500	Expect Slimline Handle Needle	19 / 1.10	2.8	1.83	box 1
M00555501	Expect Slimline Handle Needle	19 / 1.10	2.8	1.83	box 5
M00555530	Expect Slimline Handle Needle	19 Flex / 1.14	2.8	1.73	box 1
M00555531	Expect Slimline Handle Needle	19 Flex / 1.14	2.8	1.73	box 5
M00555510	Expect Slimline Handle Needle	22 / 0.72	2.4	1.65	box 1
M00555511	Expect Slimline Handle Needle	22 / 0.72	2.4	1.65	box 5
M00555520	Expect Slimline Handle Needle	25 / 0.52	2.4	1.52	box 1
M00555521	Expect Slimline Handle Needle	25 / 0.52	2.4	1.52	box 5

* Includes a 20cc syringe and one-way stopcock
Working length: 137.5cm to 141.5cm, adjustable
Needle length: 0cm to 8cm, adjustable

FINE NEEDLE BIOPSY (FNB)



Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Devices



Order Number	Description	Needle Size (gauge) / (mm)	Minimum Working Channel (mm)	Sheath Diameter (mm)	Units
M00555580	19ga Acquire FNB Needle, Flexible	19 / 1.14	2.8	1.73	box 1
M00555540	22ga Acquire EUS-FNB Needle	22 / 0.72	2.4	1.65	box 1
M00555541	22ga Acquire EUS-FNB Needle	22 / 0.72	2.4	1.65	box 5
M00555560	25ga Acquire EUS-FNB Needle	25 / 0.52	2.4	1.52	box 1
M00555561	25ga Acquire EUS-FNB Needle	25 / 0.52	2.4	1.52	box 5

FINE NEEDLE BIOPSY (FNB)



Acquire™ S Endoscopic Ultrasound Fine Needle Biopsy (FNB) Devices



Order Number	Description	Needle Size (gauge) / (mm)	Minimum Working Channel (mm)	Sheath Diameter (mm)	Units
M00555640	22ga Acquire S EUS-FNB Needle	22 / 0.72	2.4	1.65	box 1
M00555660	25ga Acquire S EUS-FNB Needle	25 / 0.52	2.4	1.52	box 1
M00555641	22ga Acquire S EUS-FNB Needle	22 / 0.72	2.4	1.65	box 5
M00555661	25ga Acquire S EUS-FNB Needle	25 / 0.52	2.4	1.52	box 5


LUMEN APPOSING METAL STENT SYSTEM

Hot AXIOS™ Stent and Electrocautery Enhanced Delivery System



Order Number	Description	Flange Diameter (mm)	Lumen Diameter (mm)	Saddle Length (mm)	Catheter Outer Diameter (F)	Catheter Working Length (cm)	Catheter Total Length (cm)	Units
M00553520	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	14	6	8	9	138	146	box 1
M00553530	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	17	8	8	9	138	146	box 1
M00553540	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	21	10	10	10.8	138	146	box 1
M00553550	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	24	15	10	10.8	138	146	box 1
M00553560	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	29	20	10	10.8	138	146	box 1
M00553570*	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	24	15	15	10.8	138	144	box 1

*The 15x15 HOT AXIOS Stent and Electrocautery Enhanced Delivery System is indicated for use to facilitate transgastric or transduodenal endoscopic drainage of a pancreatic pseudocyst or a walled-off necrosis with ≥70 % fluid content.



Guidewires

Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics



Contents



Alphabetical Index Product Code Index

HIGH PERFORMANCE




Dreamwire™ High Performance Guidewires




Order Number	Description	OD (in) / (mm)	Length (cm)	Dream Tip Style	Units
M00556101	Dreamwire Standard	.035 / 0.89	260	Straight	box 2
M00556111	Dreamwire Standard	.035 / 0.89	260	Angled	box 2
M00556121	Dreamwire Stiff Shaft	.035 / 0.89	260	Straight	box 2
M00556131	Dreamwire Stiff Shaft	.035 / 0.89	260	Angled	box 2
M00556141	Dreamwire Standard	.035 / 0.89	450	Straight	box 2
M00556151	Dreamwire Standard	.035 / 0.89	450	Angled	box 2
M00556161	Dreamwire Stiff Shaft	.035 / 0.89	450	Straight	box 2
M00556171	Dreamwire Stiff Shaft	.035 / 0.89	450	Angled	box 2


5cm Jagwire™ Guidewire Tip



10 cm Straight Dream Tip™



10 cm Angled Dream Tip™



Hydra Jagwire™ High Performance Guidewires




Order Number	Description	OD (in) / (mm)	Length (cm)	Dream Tip Style	Units
M00556001	Hydra Jagwire Standard	.035 / 0.89	260	Straight	box 2
M00556011	Hydra Jagwire Standard	.035 / 0.89	260	Angled	box 2
M00556021	Hydra Jagwire Stiff Shaft	.035 / 0.89	260	Straight	box 2
M00556031	Hydra Jagwire Stiff Shaft	.035 / 0.89	260	Angled	box 2
M00556041	Hydra Jagwire Standard	.035 / 0.89	450	Straight	box 2
M00556051	Hydra Jagwire Standard	.035 / 0.89	450	Angled	box 2
M00556061	Hydra Jagwire Stiff Shaft	.035 / 0.89	450	Straight	box 2
M00556071	Hydra Jagwire Stiff Shaft	.035 / 0.89	450	Angled	box 2

Hydra Jagwire Guidewire may remain in place during sphincterotomy. See "Instructions for Use" for details.

Guidewires

HIGH PERFORMANCE (CONTINUED)



Jagwire™ High Performance Guidewires



Order Number	Description	Distal Tip OD (in)	Working Length (cm)	Shaft Style	Tip Style	Units
M00556511	Jagwire Guidewire ANG	.025	450	Standard	Angled	box 2
M00556521	Jagwire Guidewire ST Stiff	.025	450	Stiff	Straight	box 2
M00556531	Jagwire Guidewire ANG Stiff	.025	450	Stiff	Angled	box 2
M00556561	Jagwire Guidewire ST	.025	450	Standard	Straight	box 2
M00556581	Jagwire Guidewire ST	.035	450	Standard	Straight	box 2
M00556591	Jagwire Guidewire ANG	.035	450	Standard	Angled	box 2
M00556601	Jagwire Guidewire ST Stiff	.035	450	Stiff	Straight	box 2
M00556611	Jagwire Guidewire ANG Stiff	.035	450	Stiff	Angled	box 2
M00556621	Jagwire Guidewire ST	.038	260	Standard	Straight	box 2
M00557661	Jagwire Guidewire ST Stiff	.035	260	Stiff	Straight	box 2
M00557641	Jagwire Guidewire ST	.035	260	Standard	Straight	box 2
M00557651	Jagwire Guidewire ANG	.035	260	Standard	Angled	box 2
M00557671	Jagwire Guidewire ANG Stiff	.035	260	Stiff	Angled	box 2
M00556461	Jagwire Guidewire ST	.025	260	Standard	Straight	box 2
M00556481	Jagwire Guidewire ST Stiff	.025	260	Stiff	Straight	box 2
M00556491	Jagwire Guidewire ANG Stiff	.025	260	Stiff	Angled	box 2
M00556841	Extendable Jagwire Guidewire ST	.035	260	Standard	Straight	box 2
M00556851	Extendable Jagwire Guidewire ANG	.035	260	Standard	Angled	box 2
M00556861	Extendable Jagwire Guidewire ST Stiff	.035	260	Stiff	Straight	box 2
M00556871	Extendable Jagwire Guidewire ANG Stiff	.035	260	Stiff	Angled	box 2
M00556902	Jagtail™ Guidewire Extension	.035	200	-	Guidewire Extension	box 5

Jagwire Guidewire may remain in place during sphincterotomy. See "Instructions for Use" for details. ST = Straight, ANG = Angled

Jagwire™ Revolution High Performance Guidewire




Order Number	Description	OD (in) / (mm)	Length (cm)	Tip Style	Units
M00557001	Jagwire Revolution High Performance Guidewire	.025 / 0.64	260	Straight Round	box 2
M00557011	Jagwire Revolution High Performance Guidewire	.025 / 0.64	260	Angled Round	box 2
M00557021	Jagwire Revolution High Performance Guidewire	.025 / 0.64	450	Straight Round	box 2
M00557031	Jagwire Revolution High Performance Guidewire	.025 / 0.64	450	Angled Round	box 2
M00557041	Jagwire Revolution High Performance Guidewire	.025 / 0.64	500	Straight Round	box 2



Guidewires

HIGH PERFORMANCE (CONTINUED)



NaviPro™ Hydrophilic Guidewires



Order Number	Description	OD (in) / (mm)	Length (cm)	Tip Style	Units
M00556191	NaviPro Standard	.018 / 0.46	260	Angled	box 5
M00556181	NaviPro Standard	.018 / 0.46	260	Straight	box 5
M00556211	NaviPro Standard	.025 / 0.64	260	Angled	box 5
M00556201	NaviPro Standard	.025 / 0.64	260	Straight	box 5
M00556231	NaviPro Standard	.035 / 0.89	260	Angled	box 5
M00556221	NaviPro Standard	.035 / 0.89	260	Straight	box 5
M00556251	NaviPro Stiff Shaft	.035 / 0.89	260	Angled	box 5
M00556241	NaviPro Stiff Shaft	.035 / 0.89	260	Straight	box 5


STANDARD

Amplatz Super Stiff™ Guidewire



Order Number	Description	OD (in) / (mm)	Length (cm)	Tip Style	Units
M00550090	Amplatz Super Stiff Guidewire	.038 / 0.97	260	Straight	box 1






Biopsy

Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics



FORCEPS



Radial Jaw™ 4 Standard Capacity Biopsy Forceps



Order Number	Description	Jaw OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00513383	Radial Jaw 4 Standard Capacity	2.2	160	2.8	Yellow	box 5
M00513381	Radial Jaw 4 Standard Capacity	2.2	160	2.8	Yellow	box 20
M00513382	Radial Jaw 4 Standard Capacity	2.2	160	2.8	Yellow	box 40
M00513393	Radial Jaw 4 Standard Capacity w/Needle	2.2	160	2.8	Yellow	box 5
M00513391	Radial Jaw 4 Standard Capacity w/Needle	2.2	160	2.8	Yellow	box 20
M00513392	Radial Jaw 4 Standard Capacity w/Needle	2.2	160	2.8	Yellow	box 40
M00513403	Radial Jaw 4 Standard Capacity	2.2	240	2.8	Orange	box 5
M00513401	Radial Jaw 4 Standard Capacity	2.2	240	2.8	Orange	box 20
M00513402	Radial Jaw 4 Standard Capacity	2.2	240	2.8	Orange	box 40
M00513413	Radial Jaw 4 Standard Capacity w/Needle	2.2	240	2.8	Orange	box 5
M00513411	Radial Jaw 4 Standard Capacity w/Needle	2.2	240	2.8	Orange	box 20
M00513412	Radial Jaw 4 Standard Capacity w/Needle	2.2	240	2.8	Orange	box 40

Radial Jaw™ 4 Large Capacity Biopsy Forceps



Order Number	Description	Jaw OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00513303	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	box 5
M00513301	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	box 20
M00513302	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	box 40
M00513313	Radial Jaw 4 Large Capacity w/Needle	2.4	160	2.8	Yellow	box 5
M00513311	Radial Jaw 4 Large Capacity w/Needle	2.4	160	2.8	Yellow	box 20
M00513312	Radial Jaw 4 Large Capacity w/Needle	2.4	160	2.8	Yellow	box 40
M00513323	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	box 5
M00513321	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	box 20
M00513322	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	box 40
M00513333	Radial Jaw 4 Large Capacity w/Needle	2.4	240	2.8	Orange	box 5
M00513331	Radial Jaw 4 Large Capacity w/Needle	2.4	240	2.8	Orange	box 20
M00513332	Radial Jaw 4 Large Capacity w/Needle	2.4	240	2.8	Orange	box 40

FORCEPS (CONTINUED)

Radial Jaw™ 4 Jumbo Biopsy Forceps



Order Number	Description	Jaw OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00513343	Radial Jaw 4 Jumbo	2.8	160	3.2	Yellow	box 5
M00513353	Radial Jaw 4 Jumbo w/Needle	2.8	160	3.2	Yellow	box 5
M00513363	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	box 5
M00513361	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	box 20
M00513362	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	box 40
M00513373	Radial Jaw 4 Jumbo w/Needle	2.8	240	3.2	Orange	box 5
M00513371	Radial Jaw 4 Jumbo w/Needle	2.8	240	3.2	Orange	box 20
M00513372	Radial Jaw 4 Jumbo w/Needle	2.8	240	3.2	Orange	box 40

Radial Jaw™ 4 Gastropediatric Biopsy Forceps



Order Number	Description	Jaw OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00513443	Radial Jaw 4 Gastropediatric	1.8	160	2.0	Yellow	box 5
M00513441	Radial Jaw 4 Gastropediatric	1.8	160	2.0	Yellow	box 20
M00513453	Radial Jaw 4 Gastropediatric w/Needle	1.8	160	2.0	Yellow	box 5
M00513451	Radial Jaw 4 Gastropediatric w/Needle	1.8	160	2.0	Yellow	box 20

Dilation



Contents



Alphabetical
Index

Product Code
Index

Biliary

RX Biliary
System™

Biliary Lumen
Patency

Endoscopic
Ultrasound
(EUS)

Guidewires

Biopsy

Dilation

Imaging
Systems

Pulmonary

Closure

Haemostasis


Tissue
Resection

Esophageal
and Enteral
Metal Stents

Infection
Prevention

Endobariatrics

CRE™ (CONTROLLED RADIAL EXPANSION) BALLOON DILATATION CATHETERS



CRE™ PRO Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic/Biliary) (including DASE indication)



Order Number	BALLOON INFLATED OD (mm)	BALLOON INFLATED OD (Fr)	Balloon Length (cm)	INFLATION PRESSURE (ATM)	INFLATION PRESSURE (kPa)	Minimum Working Channel (mm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Units
M00558660	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8	7.5 / 2.50	240	box 1
M00558670	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8	7.5 / 2.50	240	box 1
M00558680	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8	7.5 / 2.50	240	box 1
M00558690	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8	7.5 / 2.50	240	box 1
M00558700	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8	7.5 / 2.50	240	box 1
M00558710	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8	7.5 / 2.50	240	box 1

Recommended Inflation Device: Alliance™ II Integrated Inflation/Lithotripsy Device (Re-usable Handle) M00550620, Alliance™ II Syringe/Gauge Assembly M00550601.
This product contains no detectable latex.


CRE™ Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Colonic)



Order Number	BALLOON INFLATED OD (mm)	BALLOON INFLATED OD (Fr)	Balloon Length (cm)	INFLATION PRESSURE (ATM)	INFLATION PRESSURE (kPa)	Minimum Working Channel (mm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Units
M00558450	6-7-8	18-21-24	5.5	3-6-10	304-557-912	2.8	7.5 / 2.50	240	box 1
M00558460	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8	7.5 / 2.50	240	box 1
M00558470	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8	7.5 / 2.50	240	box 1
M00558480	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8	7.5 / 2.50	240	box 1
M00558490	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8	7.5 / 2.50	240	box 1
M00558500	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8	7.5 / 2.50	240	box 1

Recommended Inflation Device: Alliance™ II Integrated Inflation/Lithotripsy Device (Re-usable Handle) M00550620, Alliance™ II Syringe/Gauge Assembly M00550601.
This product contains no detectable latex.

CRE™ (CONTROLLED RADIAL EXPANSION) BALLOON DILATATION CATHETERS (CONTINUED)



CRE™ PRO Wireguided Balloon Dilatation Catheters (Esophageal/Pyloric/Biliary) (including DASE indication)



Order Number	BALLOON INFLATED OD (mm) (Fr)		Balloon Length (cm)	INFLATION PRESSURE (ATM) (kPa)		Minimum Working Channel (mm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Units
M00558600	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8	7.5 / 2.50	180	box 1
M00558610	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8	7.5 / 2.50	180	box 1
M00558620	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8	7.5 / 2.50	180	box 1
M00558630	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8	7.5 / 2.50	180	box 1
M00558640	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8	7.5 / 2.50	180	box 1
M00558650	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8	7.5 / 2.50	180	box 1

Recommended Inflation Device: Alliance™ II Integrated Inflation/Lithotripsy Device (Re-usable Handle) M00550620, Alliance™ II Syringe/Gauge Assembly M00550601.
This product contains no detectable latex.


CRE™ Wireguided Balloon Dilatation Catheter (Esophageal/Pyloric)



Order Number	BALLOON INFLATED OD (mm) (Fr)		Balloon Length (cm)	INFLATION PRESSURE (ATM) (kPa)		Minimum Working Channel (mm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Units
M00558390	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8	7.5 / 2.50	180	box 1
M00558400	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8	7.5 / 2.50	180	box 1
M00558410	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8	7.5 / 2.50	180	box 1
M00558420	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8	7.5 / 2.50	180	box 1
M00558430	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8	7.5 / 2.50	180	box 1
M00558440	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8	7.5 / 2.50	180	box 1

Recommended Inflation Device: Alliance™ II Integrated Inflation/Lithotripsy Device (Re-usable Handle) M00550620, Alliance™ II Syringe/Gauge Assembly M00550601.
This product contains no detectable latex.

CRE™ (CONTROLLED RADIAL EXPANSION) BALLOON DILATATION CATHETERS (CONTINUED)



CRE™ Fixed Wire Balloon Dilatation Catheters

Order Number	Inflated OD (mm)	Balloon Length (cm)	INFLATION PRESSURE (ATM)	INFLATION PRESSURE (kPa)	Minimum Working Channel (mm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Units
M00558330	6-7-8	8	3-6-10	304-608-1013	2.8	6 / 2.00	180	box 1
M00558340	8-9-10	8	3-5.5-9	304-557-912	2.8	6 / 2.00	180	box 1
M00558350	10-11-12	8	3-5-8	304-507-810	2.8	6 / 2.00	180	box 1
M00558360	12-13.5-15	8	3-4.5-8	304-456-810	2.8	6 / 2.00	180	box 1
M00558361	12-13.5-15	8	3-4.5-8	304-456-810	2.8	6 / 2.00	180	box 5
M00558370	15-16.5-18	8	3-4.5-7	304-456-709	2.8	6 / 2.00	180	box 1
M00558371	15-16.5-18	8	3-4.5-7	304-456-709	2.8	6 / 2.00	180	box 5
M00558380	18-19-20	8	3-4.5-6	304-456-608	2.8	6 / 2.00	180	box 1
M00558381	18-19-20	8	3-4.5-6	304-456-608	2.8	6 / 2.00	180	box 5

Recommended Inflation Device: Alliance™ II Integrated Inflation/Lithotripsy Device (Re-usable Handle) M00550620, Alliance™ II Syringe/Gauge Assembly M00550601.
This product contains no detectable latex.

BALLOON DILATORS

Rigiflex™ II Achalasia Balloon Dilators*


Order Number	Inflated OD (mm)	Balloon Length (cm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Recommended Guidewire (in) / (mm)**	Units
M00554500	30	10	14 / 4.67	90	.038 / 0.97	box 1
M00554510	35	10	14 / 4.67	90	.038 / 0.97	box 1
M00554520	40	10	14 / 4.67	90	.038 / 0.97	box 1

* Product shipped sterile. Recommended use with pneumatic inflator, M00553200

** Guidewire included with Rigiflex™ II Balloon Dilators

This product contains no detectable latex.

INFLATION SYSTEMS



Alliance™ II Inflation System

Order Number	Description	Units
M00550620	Alliance II Integrated Inflation/Lithotripsy Device	box 1
M00550601	Alliance II Syringe/Gauge Assembly	box 5

DILATION ACCESSORIES

Inflation/Deflation Devices and Accessories

Order Number	Description	Units
M00553200	Pneumatic Inflator	box 1
M00550630	CRE™ SteriFlate Disposable Inflation Device	box 1

RX Biliary Balloon Dilators are located in the RX Biliary Section on page 49

Imaging Systems



Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics

IMAGING SYSTEMS

EXALT™ Model D Single-Use Duodenoscope



Order Number	Item	Description	Languages	Units
M0054242CE0	EXALT Model D Single-Use Duodenoscope	Sterile Single-Use Duodenoscope for use with EXALT Controller	N/A	box 1
M0054242CE1	EXALT Model D Single-Use Duodenoscope	Sterile Single-Use Duodenoscope for use with EXALT Controller	N/A	box 2
M00542430	EXALT Controller	100–240 VAC Controller for use with Single-Use Duodenoscope (including detachable power cable, DVI and HDMI cables)	English, Spanish, French, German, Italian, Dutch and Portuguese	box 1
M00542460	EXALT Controller	100–240 VAC Controller for use with Single-Use Duodenoscope (including detachable power cable, DVI and HDMI cables)	English, Danish, Greek, Portuguese-Brazil, Swedish, Hungarian, Czech, Polish, Norwegian, Korean, Turkish, Finnish, Romanian and Slovak	box 1

EXALT™ Model D Service Codes

EXALT Model D TotalCare POS

Order Number	Description
SRVM005EXALTDTC01	EXALT Model D TotalCare POS 1 Year
SRVM005EXALTDTC02	EXALT Model D TotalCare POS 2 Year
SRVM005EXALTDTC03	EXALT Model D TotalCare POS 3 Year
SRVM005EXALTDTC04	EXALT Model D TotalCare POS 4 Year
SRVM005EXALTDTC05	EXALT Model D TotalCare POS 5 Year

EXALT Model D TotalCare Post Sale

Order Number	Description
SRVM005EXALTDTC06	EXALT Model D TotalCare POST 1 Year
SRVM005EXALTDTC07	EXALT Model D TotalCare POST 2 Year
SRVM005EXALTDTC08	EXALT Model D TotalCare POST 3 Year
SRVM005EXALTDTC09	EXALT Model D TotalCare POST 4 Year
SRVM005EXALTDTC10	EXALT Model D TotalCare POST 5 Year

IMAGING SYSTEMS (CONTINUED)

EXALT™ Model B Single-Use Bronchoscope



Order Number	Item	Units
M00542760	EXALT Monitor (mount included)	box 1
M00542780	EXALT Monitor (mount not included)	box 1
M00542770	EXALT Monitor (mount included) in alternative languages	box 1
M00542790	EXALT Monitor (mount not included) in alternative languages	box 1
M00542820	EXALT Monitor mount	box 1
M00542751	EXALT Model B Single-Use Bronchoscope – Large	box 5
M00542731	EXALT Model B Single-Use Bronchoscope – Regular	box 5
M00542711	EXALT Model B Single-Use Bronchoscope – Slim	box 5
8500-098BL9	Power Cord (EU) – Plug Type F	box 1
8520-098BL2	Power Cord (CH) – Plug Type J	box 1
8540-098BL2	Power Cord (DK) – Plug Type K	box 1
8550-098BL3	Power Cord (IT) – Plug Type L	box 1
8560-098BL3	Power Cord (IL) – Plug Type H	box 1
8580-098BL5	Power Cord (ZA) – Plug Type M	box 1
9650-098BL7	Power Cord (UK) – Plug Type G	box 1

EXALT™ Model B Service Codes

EXALT Model B Monitor EssentialCare POS

Order Number	Description
SRVM005EXBMONEC01	EXALT B Monitor EssentialCare POS 1 Year
SRVM005EXBMONEC02	EXALT B Monitor EssentialCare POS 2 Years

EXALT Model B Monitor EssentialCare PLUS POS

Order Number	Description
SRVM005EXBMONECP01	EXALT B Monitor EssCare PLUS POS 1 Year
SRVM005EXBMONECP02	EXALT B Monitor EssCare PLUS POS 2 Years

EXALT Model B Monitor Electric Safety Test

Order Number	Description
SRVM005EXBMONEST1	EXALT B Monitor Elect Safety Test (EST)
SRVM005EXBMONSTK1	EXALT B Monitor Elect Safety Test (STK)

EXALT Model B Monitor EssentialCare Post Sale

Order Number	Description
SRVM005EXBMONEC05	EXALT B Monitor EssentialCare POST 1 Year
SRVM005EXBMONEC06	EXALT B Monitor EssentialCare POST 2 Years


EXALT Model B Monitor EssentialCare PLUS Post Sale

Order Number	Description
SRVM005EXBMONECP05	EXALT B Monitor EssCare PLUS POST 1 Year
SRVM005EXBMONECP06	EXALT B Monitor EssCare PLUS POST 2 Years

EXALT Model B Rental Fees

Order Number	Description
SVC20EXALTBPLYFEE	EXALT B Daily Rental Fee for Capital
SVC20EXALTBDMNTFEE	EXALT B Monthly Rental Fee for Capital
SVC20EXALTBANNFEE	EXALT B Annual Rental Fee for Capital

DIRECT VISUALISATION SYSTEM



SpyGlass™ DS Capital Equipment

Order Number	Description	Units
M00546650	SpyGlass DS Digital Controller	box 1

SpyGlass™ DS Disposable Equipment



Order Number	Description	Units
M00546610	SpyScope™ DS II Access and Delivery Catheter	box 1
M00546550	SpyGlass™ Retrieval Basket	box 1
M00546560	SpyGlass™ Retrieval Snare	box 1
M00546470	SpyBite™ Max Biopsy Forceps	box 1
SIT-575-15	Hydra Irrigation Tubing	box 15

SpyGlass™ DS Service Codes

SpyGlass DS Digital EverCare POS

Order Number	Description
SRVM005SPYDIG1009	SpyGlass Digital EverCare POS 1 Year
SRVM005SPYDIG1010	SpyGlass Digital EverCare POS 2 Years
SRVM005SPYDIG1011	SpyGlass Digital EverCare POS 3 Years
SRVM005SPYDIG1012	SpyGlass Digital EverCare POS 4 Years

SpyGlass DS Digital EssentialCare POS

Order Number	Description
SRVM005SPYDIG1001	SpyGlass Digital EssentialCare POS 1 Year
SRVM005SPYDIG1002	SpyGlass Digital EssentialCare POS 2 Years
SRVM005SPYDIG1003	SpyGlass Digital EssentialCare POS 3 Years
SRVM005SPYDIG1004	SpyGlass Digital EssentialCare POS 4 Years

SpyGlass DS Digital EverCare Post Sale

Order Number	Description
SRVM005SPYDIG1013	SpyGlass Digital EverCare POST 1 Year
SRVM005SPYDIG1014	SpyGlass Digital EverCare POST 2 Years
SRVM005SPYDIG1015	SpyGlass Digital EverCare POST 3 Years
SRVM005SPYDIG1016	SpyGlass Digital EverCare POST 4 Years

SpyGlass DS Digital EssentialCare Post Sale

Order Number	Description
SRVM005SPYDIG1005	SpyGlass Digital EssentialCare POST 1 Year
SRVM005SPYDIG1006	SpyGlass Digital EssentialCare POST 2 Years
SRVM005SPYDIG1007	SpyGlass Digital EssentialCare POST 3 Years
SRVM005SPYDIG1008	SpyGlass Digital EssentialCare POST 4 Years

DIRECT VISUALISATION SYSTEM (CONTINUED)

SpyGlass™ Discover Digital Catheter



Order Number	Description	Units
M00546780	SpyGlass Discover Digital Catheter	box 1
M00542950	SpyGlass Discover Controller	box 1

SpyGlass™ Discover Service Codes

SpyGlass Discover EssentialCare POS

Order Number	Description
SRVM005SPYDISEC01	SpyGlass Digital EssentialCare POS 1 Year
SRVM005SPYDISEC02	SpyGlass Digital EssentialCare POS 2 Years
SRVM005SPYDISEC03	SpyGlass Digital EssentialCare POS 3 Years
SRVM005SPYDISEC04	SpyGlass Digital EssentialCare POS 4 Years

SpyGlass Discover EssentialCare Post Sale

Order Number	Description
SRVM005SPYDISEC05	SpyGlass Digital EssentialCare POST 1 Year
SRVM005SPYDISEC06	SpyGlass Digital EssentialCare POST 2 Years
SRVM005SPYDISEC07	SpyGlass Digital EssentialCare POST 3 Years
SRVM005SPYDISEC08	SpyGlass Digital EssentialCare POST 4 Years

Autolith™ Touch Biliary EHL System

Order Number	Description	Units
M00546680*	Autolith Touch EHL Generator	box 1
M00546620	1.9Fr, 375cm Biliary EHL Probe	box 1
M00546750	Autolith Touch Extender Cable	box 1
M00546760	Autolith Touch Foot Pedal	box 1

*Contains Extender Cable and Foot Pedal

Autolith™ Service Codes

Autolith EssentialCare POS

Order Number	Description
SRVM005ATLT1001	Autolith EssentialCare POS 1 Year
SRVM005ATLT1002	Autolith EssentialCare POS 2 Years
SRVM005ATLT1003	Autolith EssentialCare POS 3 Years
SRVM005ATLT1004	Autolith EssentialCare POS 4 Years



Autolith EssentialCare PLUS POS

Order Number	Description
SRVM005ATLTECP1001	Autolith EssentialCare PLUS POS 3 Years
SRVM005ATLTECP1002	Autolith EssentialCare PLUS POS 4 Years
SRVM005ATLTECP1003	Autolith EssentialCare PLUS POS 5 Years

Autolith EssentialCare PLUS Post Sale

Order Number	Description
SRVM005ATLTECP1004	Autolith EssentialCare PLUS POST 3 Years
SRVM005ATLTECP1005	Autolith EssentialCare PLUS POST 4 Years
SRVM005ATLTECP1006	Autolith EssentialCare PLUS POST 5 Years

Pulmonary



Biliary
Rx Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics



Contents



Alphabetical Index Product Code Index

IMAGING SYSTEMS

EXALT™ Model B Single-Use Bronchoscope



Order Number	Item	Units
M00542760	EXALT Monitor (mount included)	box 1
M00542780	EXALT Monitor (mount not included)	box 1
M00542770	EXALT Monitor (mount included) in alternative languages	box 1
M00542790	EXALT Monitor (mount not included) in alternative languages	box 1
M00542820	EXALT Monitor mount	box 1
M00542751	EXALT Model B Single-Use Bronchoscope – Large	box 5
M00542731	EXALT Model B Single-Use Bronchoscope – Regular	box 5
M00542711	EXALT Model B Single-Use Bronchoscope – Slim	box 5
8500-098BL9	Power Cord (EU) – Plug Type F	box 1
8520-098BL2	Power Cord (CH) – Plug Type J	box 1
8540-098BL2	Power Cord (DK) – Plug Type K	box 1
8550-098BL3	Power Cord (IT) – Plug Type L	box 1
8560-098BL3	Power Cord (IL) – Plug Type H	box 1
8580-098BL5	Power Cord (ZA) – Plug Type M	box 1
9650-098BL7	Power Cord (UK) – Plug Type G	box 1

EXALT™ Model B Service Codes

EXALT Model B Monitor EssentialCare POS

Order Number	Description
SRVM005EXBMONEC01	EXALT B Monitor EssentialCare POS 1 Year
SRVM005EXBMONEC02	EXALT B Monitor EssentialCare POS 2 Years

EXALT Model B Monitor EssentialCare PLUS POS

Order Number	Description
SRVM005EXBMONECP01	EXALT B Monitor EssCare PLUS POS 1 Year
SRVM005EXBMONECP02	EXALT B Monitor EssCare PLUS POS 2 Years

EXALT Model B Monitor Electric Safety Test

Order Number	Description
SRVM005EXBMONEST1	EXALT B Monitor Elect Safety Test (EST)
SRVM005EXBMONSTK1	EXALT B Monitor Elect Safety Test (STK)

EXALT Model B Monitor EssentialCare Post Sale

Order Number	Description
SRVM005EXBMONEC05	EXALT B Monitor EssentialCare POST 1 Year
SRVM005EXBMONEC06	EXALT B Monitor EssentialCare POST 2 Years

EXALT Model B Monitor EssentialCare PLUS Post Sale

Order Number	Description
SRVM005EXBMONECP05	EXALT B Monitor EssCare PLUS POST 1 Year
SRVM005EXBMONECP06	EXALT B Monitor EssCare PLUS POST 2 Years

EXALT Model B Rental Fees

Order Number	Description
SVC20EXALTBPLYFEE	EXALT B Daily Rental Fee for Capital
SVC20EXALTBDMNTFEE	EXALT B Monthly Rental Fee for Capital
SVC20EXALTBANNFEE	EXALT B Annual Rental Fee for Capital

ENDOSCOPIC BRONCHIAL ULTRASOUND



Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needles Olympus Scope Compatible




Order Number	Description	Minimum Working Channel (mm)	Sheath OD (mm)	Units
M00558220	Expect Pulmonary 22ga – Olympus Scope Compatible	2.0	1.6	box 1
M00558221	Expect Pulmonary 22ga – Olympus Scope Compatible	2.0	1.6	box 5
M00558250	Expect Pulmonary 25ga – Olympus Scope Compatible	2.0	1.4	box 1
M00558251	Expect Pulmonary 25ga – Olympus Scope Compatible	2.0	1.4	box 5
M00558731	Expect Pulmonary Needle Adaptor – Olympus Scope Compatible			box 10

Acquire™ Pulmonary Endobronchial Ultrasound Fine Biopsy Needles Olympus Scope Compatible



Order Number	Description	Needle Size (gauge) / (mm)	Minimum Working Channel (mm)	Sheath Diameter (mm)	Units
M00552350	22ga Acquire Pulmonary – Olympus scope compatible	22/0.72	2.0	1.6	box 1
M00552351	22ga Acquire Pulmonary – Olympus scope compatible	22/0.72	2.0	1.6	box 5
M00552400	25ga Acquire Pulmonary – Olympus scope compatible	25/0.52	2.0	1.4	box 1
M00552401	25ga Acquire Pulmonary – Olympus scope compatible	25/0.52	2.0	1.4	box 5

AIRWAY STENTS



Ultraflex™ Tracheobronchial Covered Stent System Distal Release



Order Number	Description	Stent OD (mm)	Expanded Stent Length (mm)	Covered Length (mm)
M00564820	Ultraflex Tracheal NG Stent System	14	80	65
M00564830	Ultraflex Tracheal NG Stent System	16	40	25
M00564840	Ultraflex Tracheal NG Stent System	16	60	45
M00564850	Ultraflex Tracheal NG Stent System	16	80	65
M00564860	Ultraflex Tracheal NG Stent System	18	40	25
M00564870	Ultraflex Tracheal NG Stent System	18	60	45
M00564880	Ultraflex Tracheal NG Stent System	18	80	65
M00564890	Ultraflex Tracheal NG Stent System	20	40	25
M00564900	Ultraflex Tracheal NG Stent System	20	60	45
M00564910	Ultraflex Tracheal NG Stent System	20	80	65
M00564740	Ultraflex Tracheobronchial Stent System	8	40	25
M00564750	Ultraflex Tracheobronchial Stent System	10	30	15
M00564760	Ultraflex Tracheobronchial Stent System	10	40	25
M00564770	Ultraflex Tracheobronchial Stent System	12	30	15
M00564780	Ultraflex Tracheobronchial Stent System	12	40	25
M00564790	Ultraflex Tracheobronchial Stent System	14	30	15
M00564800	Ultraflex Tracheobronchial Stent System	14	40	25
M00564810	Ultraflex Tracheobronchial Stent System	14	60	45

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Jagwire™ Pulmonary Guidewire, M00515171 / .038in (0.97mm) Amplatz Super Stiff™ Guidewire, M00550090

Additional Specifications for Ultraflex Tracheobronchial Stents: Delivery System Working Length is 95cm, Maximum Catheter Shaft OD is 16F (5.33mm), Compressed Stent on catheter yields Delivery System OD range 15-25F (5.00-8.33mm)



Ultraflex™ Tracheobronchial Uncovered Stent System Proximal Release




Order Number	Description	Stent OD (mm)	Expanded Stent Length (mm)
M00564710	Ultraflex Tracheobronchial Stent System	14	40
M00564720	Ultraflex Tracheobronchial Stent System	14	60
M00564640	Ultraflex Tracheobronchial Stent System	8	20
M00564650	Ultraflex Tracheobronchial Stent System	8	40
M00564660	Ultraflex Tracheobronchial Stent System	10	20
M00564670	Ultraflex Tracheobronchial Stent System	10	40
M00564680	Ultraflex Tracheobronchial Stent System	12	20
M00564690	Ultraflex Tracheobronchial Stent System	12	40
M00564700	Ultraflex Tracheobronchial Stent System	14	20

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Jagwire™ Pulmonary Guidewire, M00515171 / .038in (0.97mm) Amplatz Super Stiff™ Guidewire, M00550090

Additional Specifications for Ultraflex Tracheobronchial Stents: Delivery System Working Length is 95cm, Maximum Catheter Shaft OD is 16F (5.33mm), Compressed Stent on catheter yields Delivery System OD range 15-18F (5.00-6.00mm)

AIRWAY STENTS (CONTINUED)



Ultraflex™ Tracheobronchial Uncovered Stent System Distal Release




Order Number	Description	Stent OD (mm)	Expanded Stent Length (mm)
M00564530	Ultraflex Tracheobronchial Stent System	16	40
M00564540	Ultraflex Tracheobronchial Stent System	16	60
M00564560	Ultraflex Tracheobronchial Stent System	18	40
M00564570	Ultraflex Tracheobronchial Stent System	18	60
M00564590	Ultraflex Tracheobronchial Stent System	20	40
M00564600	Ultraflex Tracheobronchial Stent System	20	60
M00564500	Ultraflex Tracheobronchial Stent System	10	30
M00564510	Ultraflex Tracheobronchial Stent System	12	30
M00564520	Ultraflex Tracheobronchial Stent System	14	30

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Jagwire™ Pulmonary Guidewire, M00515171 / .038in (0.97mm) Amplatz Super Stiff™ Guidewire, M00550090

Additional Specifications for Ultraflex Tracheobronchial Stents: Delivery System Working Length is 95cm, Maximum Catheter Shaft OD is 16F (5.33mm), Compressed Stent on catheter yields Delivery System OD range 15-20F (5.00-6.67mm)

BIOPSY FORCEPS



Radial Jaw™ 4 Pulmonary Biopsy Forceps



Order Number	Description	Jaw OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00515181	Radial Jaw 4 Pulmonary Standard Capacity	1.8	100	2.0	Blue	box 5
M00515182	Radial Jaw 4 Pulmonary Standard Capacity	1.8	100	2.0	Blue	box 20
M00515191	Radial Jaw 4 Pulmonary Standard Capacity w/Needle	1.8	100	2.0	Blue	box 5
M00515192	Radial Jaw 4 Pulmonary Standard Capacity w/Needle	1.8	100	2.0	Blue	box 20
M00515201	Radial Jaw 4 Pulmonary Large Capacity	2.2	100	2.8	Blue	box 5
M00515202	Radial Jaw 4 Pulmonary Large Capacity	2.2	100	2.8	Blue	box 20


MINI-FORCEPS

CoreDx™ Pulmonary Mini-Forceps



Order Number	Description	Jaw OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00515220	Pulmonary Mini-Forceps	1.0	130	1.2	box 1

TRANSBRONCHIAL ASPIRATION NEEDLES



eXcelon™ Transbronchial Aspiration Needles



Order Number	Description	Gauge	NEEDLE OD (mm)	ID (mm)	Needle Length (mm)	Catheter Length (cm)	Sheath OD (mm)	Units
M00564101	eXcelon Transbronchial Aspiration Needle	19	1.07	0.69	15	130	1.8	box 5
M00564111	eXcelon Transbronchial Aspiration Needle	20	0.90	0.58	15	130	1.8	box 5
M00564121	eXcelon Transbronchial Aspiration Needle	21	0.81	0.50	15	130	1.8	box 5


CYTOLOGY BRUSHES

Celebrity™ Cytology Brushes – Bronchoscope



Order Number	Description	Required Working Channel (mm)	Bristle OD (mm)	Sheath Length (cm)	Units
M00516001	Celebrity Cytology Brush	2.0	1.0	140	box 10
M00516011	Celebrity Cytology Brush	2.0	1.5	140	box 10
M00516071	Celebrity Cytology Brush	2.0	1.9	100	box 10
M00516151	Celebrity Cytology Brush	2.0	1.9	150	box 10

AIRWAY RETRIEVAL BASKETS



Zero Tip™ Airway Retrieval Baskets



Order Number	Description	Sheath OD (mm)	Working Length (cm)	Basket Opening (mm)
M00513200	Zero Tip Airway Retrieval Basket	0.8	120	12
M00513210	Zero Tip Airway Retrieval Basket	1.0	120	16

Packaged one per box

DILATION BALLOONS



CRE™ Pulmonary Balloon Dilatation Catheter



Order Number	Description	Balloon Length (cm)	Initial Balloon OD (mm)	Intermediate Balloon OD (mm)	Maximum Balloon OD (mm)	Working Catheter Length (cm)
M00550300	CRE Pulmonary Balloon Dilatation Catheter	5.5	12	13.5	15	110
M00550310	CRE Pulmonary Balloon Dilatation Catheter	5.5	15	16.5	18	110
M00550320	CRE Pulmonary Balloon Dilatation Catheter	5.5	18	19.0	20	110
M00550330	CRE Pulmonary Balloon Dilatation Catheter	3.0	8	9.0	10	110
M00550340	CRE Pulmonary Balloon Dilatation Catheter	3.0	10	11.0	12	110
M00550350	CRE Pulmonary Balloon Dilatation Catheter	3.0	12	13.5	15	110

Packaged one per box


GUIDEWIRES

Jagwire™ Pulmonary Guidewire



Order Number	Description	OD (in) / (mm)	Length (cm)	Units
M00515171	Jagwire Pulmonary Guidewire	.035 / 0.89	180	box 2

INFLATION SYSTEM



Alliance™ II Inflation System



Order Number	Description	Units
M00550620	Alliance II Integrated Inflation/Lithotripsy Device	box 1
M00550601	Alliance II Single-Use Syringe/Gauge Assembly	box 5

Closure



Biliary
RX Biliary™ System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics



Contents



Alphabetical Index Product Code Index

CLIPPING DEVICES

X-Tack™



New / Old Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Units
M00505500 / XTACK-160-H*	X-Tack Endoscopic HeliX Tacking System (Gastric)	> 230	2.8	box 1
M00505510 / XTACK-235-H*	X-Tack Endoscopic HeliX Tacking System (Colonic)	> 230	2.8	box 1

* The UPN change will be effective starting April 18th.

X-Tack™ Accessories



New / Old Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Units
M00505550 / CNH-C01-213-L*	Suture Cinch – Long	≥ 204	2.8	box 1
M00505540 / CNH-G01-000*	OverStitch Suture Cinch	≥ 141	2.8	box 1
ES26165-3	Ensizor Endoscopic Scissors Gastric Length	165	2.8	box 3
ES26235-3	Ensizor Endoscopic Scissors Colonic Length	235	2.8	box 3

* The UPN change will be effective starting April 18th.

ENDOSCOPIC SUTURING

OverStitch™ Systems



New / Old Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Units
M00505400 / ESS-G02-160***	OverStitch Endoscopic Suturing System	*	-	box 1
M00505441 / PLY-G02-020-APL**	OverStitch 2-0 Polypropylene Suture	*	2.8	box 12
M00505540 / CNH-G01-000**	OverStitch Suture Cinch	160	3.2	box 1
M00505470 / THX-165-028***	OverStitch Tissue Helix	*	2.8	box 1
M00505480 / NXT-THXP-130***	OverStitch NXT Tissue Helix Pro	165	2.8	box 1
M00505420 / ESS-G02-NXT***	OverStitch NXT Endoscopic Suture System	*	-	box 1

* Suitable for endoscopes with max 110 cm working length.

** The UPN change will be effective starting April 18th.

*** The UPN change will be completed by the end of 2025.

OverStitch™ Accessories



New / Old Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Units
M00505560 / OVT-027-160*	OverTube Endoscopic Access System	-	-	box 3
ES26165-3	Ensizor Endoscopic Scissors Gastric Length	165	2.8	box 3

* The UPN change will be completed by the end of 2025.

At the time of the UPN change, units will change from box 3 to box 1.

ENDOSCOPIC SUTURING (CONTINUED)

NEW Mantis™



Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Units
M00521420	Mantis	235	2.8	box 1
M00521421	Mantis	235	2.8	box 10
M00521422	Mantis	235	2.8	box 20
M00521423	Mantis	235	2.8	box 40

ENDOLUMINAL VACUUM THERAPY

NEW Eso-SPONGE®



Order Number	Description	Inner Diameter (mm)	Outer Diameter (mm)	Units
5526540	Eso-SPONGE with 15 mm Overtube	15	19	box 1
5526550	Eso-SPONGE with 13 mm Overtube	13	17	box 1

NEW Endo-SPONGE®



Order Number	Description	Units
5526510	Endo-SPONGE	box 10
5526520	Endo-SPONGE	box 5
5526530	Endo-SPONGE	box 1

NEW EVT Vacuum Sources



Order Number	Description	Units
5526604	Redyrob Trans Plus	box 10
5526650	Endo-SPONGE® Pump*	box 1
5526651	Endo-SPONGE® Pump Carrying Bag	box 100
5526653	Disposable Cannister 300 ml with tubing and Y-connector	box 16
19116	MV1 Mobile Programmable Low-Vacuum-System*	box 1
18032	Collection Bottle MV1	box 15
18022	Bacterial Filter MV1	box 20

* The "Do Not Reuse" symbol does not apply to this product.

Haemostasis



Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics



CLIPPING DEVICES



Resolution™ Clips



Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Units
M00522600	Resolution Clip	155	2.8	11	box 1
M00522601	Resolution Clip	155	2.8	11	box 10
M00522602	Resolution Clip	155	2.8	11	box 20
M00522610	Resolution Clip	235	2.8	11	box 1
M00522611	Resolution Clip	235	2.8	11	box 10
M00522612	Resolution Clip	235	2.8	11	box 20

* MR Conditional. Non-clinical testing has demonstrated the Resolution 360™ Clip is Magnetic Resonance (MR) Conditional according to ASTM F2503. Refer to the directions for use for MR conditions of use.



Resolution 360™ Clips



Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Units
M00521230	Resolution 360 Clip	235	2.8	11	box 1
M00521231	Resolution 360 Clip	235	2.8	11	box 10
M00521232†	Resolution 360 Clip	235	2.8	11	box 20
M00521233†	Resolution 360 Clip	235	2.8	11	box 40
M00521240	Resolution 360 Clip	155	2.8	11	box 1
M00521241	Resolution 360 Clip	155	2.8	11	box 10
M00521242	Resolution 360 Clip	155	2.8	11	box 20
M00521243	Resolution 360 Clip	155	2.8	11	box 40

* MR Conditional. Non-clinical testing has demonstrated the Resolution 360™ Clip is Magnetic Resonance (MR) Conditional according to ASTM F2503. Refer to the directions for use for MR conditions of use.


† Please check product availability in your country with your sales representative

CLIPPING DEVICES (CONTINUED)

Resolution 360™ ULTRA Clips



Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Units
M00521400	Resolution 360 ULTRA Clip	235	2.8	17	box 1
M00521401	Resolution 360 ULTRA Clip	235	2.8	17	box 10
M00521402	Resolution 360 ULTRA Clip	235	2.8	17	box 20
M00521403	Resolution 360 ULTRA Clip	235	2.8	17	box 40




Speedband Superview Super 7™ Multiple Band Ligators



Order Number	Description	Minimum Working Channel (mm)	Scope OD (mm)	Bands / Ligating Unit	Units
M00542250	Super 7 Band Ligator	2.8	8.6-11.5	7	box 1
M00542251	Super 7 Band Ligator	2.8	8.6-11.5	7	box 2
M00542253	Super 7 Band Ligator	2.8	8.6-11.5	7	box 4

Irrigation catheter supplied in each kit

BIPOLAR ELECTROHAEMOSTASIS CATHETERS




Injection Gold Probe™ Electrohaemostasis/Injection Catheters



Order Number	Description	OD (F) / (mm)	Gauge	NEEDLE OD (mm)	ID (mm)	Minimum Working Channel (mm)	Working Length (cm)	Units
M00560150	Injection Gold Probe Catheter	7 / 2.33	25	0.51	0.24	2.8	210	box 1
M00560160	Injection Gold Probe Catheter	10 / 3.33	25	0.51	0.24	3.7	210	box 1

Type Generator/Connector: Please refer to Electrosurgical Bipolar Cable Adaptors

BIPOLAR ELECTROHAEMOSTASIS CATHETERS (CONTINUED)




Gold Probe™ Electrohaemostasis Catheters



Order Number	Description	OD (F) / (mm)	Minimum Working Channel (mm)	Working Length (cm)	Units
M00560070	Gold Probe Catheter	7 / 2.33	2.8	300	box 1
M00560100	Gold Probe Catheter	10 / 3.33	3.7	300	box 1
M00560220	Gold Probe Catheter	7 / 2.33	2.8	350	box 1

Type Generator/Connector: Please refer to Electrosurgical Bipolar Cable Adaptors

NEEDLES




Interject™ Sclerotherapy Needles*



Order Number	Description	Sheath Design	Sheath OD (F) / (mm)	Gauge	NEEDLE OD (mm)	ID (mm)	Maximum Needle Extension Length (mm)	Working Length (cm)	Units
M00518151	Interject Sclerotherapy Needle	Contrast	1.8	23	0.64	0.32	4	200	box 5
M00518161	Interject Sclerotherapy Needle	Contrast	1.8	25	0.51	0.24	4	200	box 5
M00518251	Interject Sclerotherapy Needle	Contrast	1.8	23	0.64	0.32	6	200	box 5
M00518261	Interject Sclerotherapy Needle	Contrast	1.8	25	0.51	0.24	6	200	box 5
M00518351	Interject Sclerotherapy Needle	Contrast	2.3	23	0.64	0.32	4	240	box 5
M00518361	Interject Sclerotherapy Needle	Contrast	2.3	25	0.51	0.24	4	240	box 5

* With Star Catheter Technology

NEEDLES (CONTINUED)




Interject™ Clear Sclerotherapy Needles*



Order Number	Description	Sheath Design	Sheath OD (F) / (mm)	Gauge	NEEDLE OD (mm)	ID (mm)	Maximum Needle Extension Length (mm)	Working Length (cm)	Units
M00518111	Interject Sclerotherapy Needle	Clear	1.8	25	0.51	0.24	4	200	box 5
M00518301	Interject Sclerotherapy Needle	Clear	2.3	23	0.64	0.32	4	240	box 5
M00518311	Interject Sclerotherapy Needle	Clear	2.3	25	0.51	0.24	4	240	box 5

* With Star Catheter Technology

ELECTROSURGICAL BIPOLAR CABLE ADAPTOR




Order Number	Description	Fix Pin Spacing (mm)	Generator Compatibility*	Units
M00505200	Bipolar Cable Adaptor	22	ERBE, EMC Dolley, Olympus (ESG100, PSD60)**	box 1
M00505210	Bipolar Cable Adaptor	28	Valleylab®, Lamidey (Duo), Bovie (Aaron 3250™), Olympus (UES10, UES30, Surgemaster), Bowa, Soring	box 1
M00505220	Bipolar Cable Adaptor	28 with Third Pin	ConMed® (System 5000), Bovie (Aaron 2250™), Ethicon®	box 1

Bipolar Cable Adaptors are used to connect Gold Probe™ and Injection Gold Probe™ with the above mentioned generators.

* For generators and their associated cable adaptors not listed above, please contact your Boston Scientific Representative.

** Some variations in faceplates may exist within a generator model. We recommend that you confirm pin spacing prior to ordering.

Tissue Resection



Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics



Tissue Resection

EMR DEVICES FOR UPPER GI



Captivator™ EMR



Order Number	Description	Units
M00561600	Captivator EMR for Standard Gastrosopes	box 1
M00561601	Captivator EMR for Standard Gastrosopes	box 5
M00561610	Captivator EMR for Large Gastrosopes	box 1
M00561611	Captivator EMR for Large Gastrosopes	box 5


ELECTROSURGICAL KNIFE

NEW ORISE™ ProKnife



Order Number	Description	Catheter OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00519351	1.5mm ORISE ProKnife	2.21	230	2.8	box 1
M00519361	2.0mm ORISE ProKnife	2.21	230	2.8	box 1
M00519371	3.0mm ORISE ProKnife	2.21	230	2.8	box 1

SNARES



Rotatable Snares*




Order Number	Description	Stiffness	Loop Width (mm)	Sheath OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00561821	Small Oval	Medium Stiff	13	2.4	240	2.8	box 10
M00561831	Medium Oval	Medium Stiff	20	2.4	240	2.8	box 10

* With Star Catheter Technology.


Tissue Resection

SNARES (CONTINUED)



Sensation™ Short Throw Snares


Order Number	Description	Stiffness	Loop Width (mm)	Sheath OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562691	Small Oval	Flexible	13	2.4	240	2.8	box 5
M00562692	Small Oval	Flexible	13	2.4	240	2.8	box 20
M00562693	Small Oval	Flexible	13	2.4	240	2.8	box 40
M00560321	Small Oval	Medium Stiff	13	2.4	240	2.8	box 5
M00562671	Medium Oval	Flexible	27	2.4	240	2.8	box 5
M00562672	Medium Oval	Flexible	27	2.4	240	2.8	box 20
M00562673	Medium Oval	Flexible	27	2.4	240	2.8	box 40
M00560311	Medium Oval	Medium Stiff	27	2.4	240	2.8	box 5
M00562651	Large Oval	Medium Stiff	30	2.4	240	2.8	box 5
M00562652	Large Oval	Medium Stiff	30	2.4	240	2.8	box 20
M00560331	Medium Crescent	Medium Stiff	27	2.4	240	2.8	box 5



Captiflex™ Single-Use Snares


Order Number	Description	Stiffness	Loop Width (mm)	Sheath OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562471	Extra Small Oval	Flexible	11	2.4	240	2.8	box 10
M00562421	Small Oval	Flexible	13	2.4	240	2.8	box 10
M00562422	Small Oval	Flexible	13	2.4	240	2.8	box 40
M00562401	Medium Oval	Flexible	27	2.4	240	2.8	box 10
M00562402	Medium Oval	Flexible	27	2.4	240	2.8	box 40

SNARES (CONTINUED)



Captivator™ Single-Use Snares

Order Number	Description	Stiffness	Loop Width (mm)	Sheath OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562301	Small Oval	Stiff	13	2.4	240	2.8	box 10
M00562321	Medium Oval	Stiff	27	2.4	240	2.8	box 10
M00562391	Large Oval (thin wire)	Flexible	30	2.4	240	2.8	box 10
M00561311	Large Oval	Medium Stiff	30	2.4	240	2.8	box 10
M00562451	Small Hexagonal	Stiff	13	2.4	240	2.8	box 10
M00562341	Medium Hexagonal	Stiff	27	2.4	240	2.8	box 10
M00562371	Medium Crescent	Stiff	27	2.4	240	2.8	box 10




Captivator™ II Single-Use Snares

Order Number	Description	Stiffness	Loop Width (mm)	Sheath OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00561221	10 mm Rounded	Stiff	10	2.4	240	2.8	box 10
M00561222	10 mm Rounded	Stiff	10	2.4	240	2.8	box 20
M00561223	10 mm Rounded	Stiff	10	2.4	240	2.8	box 40
M00561231	15 mm Rounded	Stiff	15	2.4	240	2.8	box 10
M00561232	15 mm Rounded	Stiff	15	2.4	240	2.8	box 20
M00561233	15 mm Rounded	Stiff	15	2.4	240	2.8	box 40
M00561241	20 mm Rounded	Stiff	20	2.4	240	2.8	box 10
M00561242	20 mm Rounded	Stiff	20	2.4	240	2.8	box 20
M00561243	20 mm Rounded	Stiff	20	2.4	240	2.8	box 40
M00561191	25 mm Rounded	Stiff	25	2.4	240	2.8	box 10
M00561192	25 mm Rounded	Stiff	25	2.4	240	2.8	box 20
M00561193	25 mm Rounded	Stiff	25	2.4	240	2.8	box 40
M00561291	33 mm Rounded	Stiff	33	2.4	240	2.8	box 10

Tissue Resection

SNARES (CONTINUED)




Captivator™ COLD Single-Use Snare



Order Number	Description	Loop Width (mm)	Sheath OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00561101	10mm Rounded-Stiff	10	2.4	240	2.8	box 10
M00561102	10mm Rounded-Stiff	10	2.4	240	2.8	box 20
M00561103	10mm Rounded-Stiff	10	2.4	240	2.8	box 40

HOT BIOPSY FORCEPS



Radial Jaw™ 4 Hot Biopsy Forceps



Order Number	Description	Jaw OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00515031	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	box 5
M00515032	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	box 20
M00515033	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	box 40

ACCESSORIES

Active Cords

Order Number	Description	Units
M00561250	Active Cord, Fits Olympus® Brand Generator	box 1
M00561270	Active Cord, Fits Endostat™, Endostat™ II, Endostat™ III, Bovie®, Valleylab®, Aspen Labs and ERBE® Brand Generators	box 1



Tissue Resection

RETRIEVAL DEVICES

RescueNet™



Order Number	Description	Net Size (mm)	Sheath OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
DGN-538-5	RescueNet Retrieval Device	30	2.5	230	2.8	box 5

Rescue™ Forceps



Order Number	Description	Jaw Opening (mm)	Sheath OD (mm)	Minimum Working Channel (cm)	Working Length (cm)	Units
M00527651	Rescue Rat Tooth Grasping Forceps	7	2.4	2.8	230	box 5
M00527751	Rescue Rat Tooth Grasping Forceps – Pediatric	4.5	1.8	2.0	180	box 5
M00527551	Rescue Alligator Grasping Forceps Long	9	2.4	2.8	230	box 5
M00527451	Rescue Alligator Grasping Forceps	8	2.4	2.8	230	box 5
M00538351	Rescue Rat Tooth / Alligator Grasping Forceps	8	2.4	2.8	230	box 5

POLYP TRAPS


TrapEase™ Polyp Traps



Order Number	Description	Units
STE-298-25	TrapEase Quad Trap	box 25
STE-297-50	TrapEase Single Trap	box 50




Esophageal and Enteral Metal Stents



Esophageal and Enteral Metal Stents

ESOPHAGEAL STENTS



WallFlex™ Fully Covered RMV Esophageal Longer Loop Stent System

Order Number	Description	Stent OD (mm)	Proximal / Distal Flares OD (mm)	Stent Length (mm)	Catheter Diameter (F)	Working Length (cm)
M00516270	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	18	25 / 23	103	18.5	120
M00516280	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	18	25 / 23	123	18.5	120
M00516290	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	18	25 / 23	153	18.5	120
M00516300	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	23	28 / 28	105	18.5	120
M00516310	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	23	28 / 28	125	18.5	120
M00516320	WallFlex Fully Covered RMV Esophageal Longer Loop Stent System	23	28 / 28	155	18.5	120

Packaged one per box

Limited availability in all European countries (except UK and France). Please contact your local Territory Manager.

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621

WallFlex™ Fully Covered RMV Esophageal Stent

Order Number	Description	Stent OD (mm)	Stent Length (mm)	Working Length (cm)
M00516210	WallFlex Fully Covered RMV Esophageal Stent	18	103	78
M00516220	WallFlex Fully Covered RMV Esophageal Stent	18	123	78
M00516230	WallFlex Fully Covered RMV Esophageal Stent	18	153	78
M00516240	WallFlex Fully Covered RMV Esophageal Stent	23	105	78
M00516250	WallFlex Fully Covered RMV Esophageal Stent	23	125	78
M00516260	WallFlex Fully Covered RMV Esophageal Stent	23	155	78

Packaged one per box

Limited availability in all European countries (except UK and France). Please contact your local Territory Manager.

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621



Esophageal and Enteral Metal Stents

ESOPHAGEAL STENTS (CONTINUED)



WallFlex™ Partially Covered Esophageal Stent System



Order Number	Description	Stent OD (mm)	Proximal / Distal Flares OD (mm)	Stent Length (mm)	Covered Length (mm)	Catheter Diameter (F)	Working Length (cm)
M00516900	WallFlex Partially Covered Esophageal Stent System	18	23 / 23	103	71	18.5	120
M00516910	WallFlex Partially Covered Esophageal Stent System	18	23 / 23	123	90	18.5	120
M00516920	WallFlex Partially Covered Esophageal Stent System	18	23 / 23	153	121	18.5	120
M00516930	WallFlex Partially Covered Esophageal Stent System	23	28 / 28	105	73	18.5	120
M00516940	WallFlex Partially Covered Esophageal Stent System	23	28 / 28	125	93	18.5	120
M00516950	WallFlex Partially Covered Esophageal Stent System	23	28 / 28	155	123	18.5	120

Packaged one per box

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621

WallFlex™ Fully Covered Esophageal Stent



Order Number	Description	Stent OD (mm)	Stent Length (mm)	Working Length (cm)
M00516700	WallFlex Fully Covered Esophageal Stent	18	103	78
M00516710	WallFlex Fully Covered Esophageal Stent	18	123	78
M00516720	WallFlex Fully Covered Esophageal Stent	18	153	78
M00516730	WallFlex Fully Covered Esophageal Stent	23	105	78
M00516740	WallFlex Fully Covered Esophageal Stent	23	125	78
M00516750	WallFlex Fully Covered Esophageal Stent	23	155	78


Packaged one per box

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621



Esophageal and Enteral Metal Stents

ESOPHAGEAL STENTS (CONTINUED)



Ultraflex™ Covered Esophageal NG Stent System

Proximal Release




Order Number	Description	Stent OD (mm)	Proximal Flare OD (mm)	Stent Length (mm)	Covered Length (mm)
M00513840	Ultraflex Esophageal NG Stent System	18	23	100	70
M00513850	Ultraflex Esophageal NG Stent System	18	23	120	90
M00513860	Ultraflex Esophageal NG Stent System	18	23	150	120
M00514240	Ultraflex Esophageal NG Stent System	23	28	100	70
M00514250	Ultraflex Esophageal NG Stent System	23	28	120	90

Distal Release

Order Number	Description	Stent OD (mm)	Proximal Flare OD (mm)	Stent Length (mm)	Covered Length (mm)
M00513730	Ultraflex Esophageal NG Stent System	18	23	100	70
M00513740	Ultraflex Esophageal NG Stent System	18	23	120	90
M00513750	Ultraflex Esophageal NG Stent System	18	23	150	120
M00514200	Ultraflex Esophageal NG Stent System	23	28	100	70
M00514210	Ultraflex Esophageal NG Stent System	23	28	120	90

Packaged one per box

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621



Ultraflex™ Uncovered Esophageal NG Stent System

Proximal Release



Order Number	Description	Stent OD (mm)	Proximal Flare OD (mm)	Stent Length (mm)
M00513800	Ultraflex Esophageal NG Stent System	18	23	70
M00513810	Ultraflex Esophageal NG Stent System	18	23	100
M00513820	Ultraflex Esophageal NG Stent System	18	23	120
M00513830	Ultraflex Esophageal NG Stent System	18	23	150
M00514220	Ultraflex Esophageal NG Stent System	23	28	70
M00514230	Ultraflex Esophageal NG Stent System	23	28	120

Distal Release

Order Number	Description	Stent OD (mm)	Proximal Flare OD (mm)	Stent Length (mm)
M00513700	Ultraflex Esophageal NG Stent System	18	23	70
M00513710	Ultraflex Esophageal NG Stent System	18	23	100
M00513720	Ultraflex Esophageal NG Stent System	18	28	150

Packaged one per box

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621



Contents

Alphabetical Index

Product Code Index

Biliary

RX Biliary System™

Biliary Lumen Patency

Endoscopic Ultrasound (EUS)

Guidewires

Biopsy

Dilation

Imaging Systems

Pulmonary

Closure

Haemostasis

Tissue Resection

Esophageal and Enteral Metal Stents

Infection Prevention

Endobariatrics

Esophageal and Enteral Metal Stents

ESOPHAGEAL STENTS (CONTINUED)

Agile™ Fully Covered Esophageal Stent System – Removable



Order Number	Description	Stent OD (mm)	Proximal / Distal Flares OD (mm)	Stent Length (mm)	Delivery System Diameter (Fr)	Working Length (cm)
M00517580	Agile Fully Covered Esophageal Stent System – Removable	14	19	60	10.5	180
M00517590	Agile Fully Covered Esophageal Stent System – Removable	14	19	100	10.5	180
M00517600	Agile Fully Covered Esophageal Stent System – Removable	14	19	120	10.5	180
M00517610	Agile Fully Covered Esophageal Stent System – Removable	14	19	150	10.5	180
M00517620	Agile Fully Covered Esophageal Stent System – Removable	18	23	60	10.5	180
M00517630	Agile Fully Covered Esophageal Stent System – Removable	18	23	100	10.5	180
M00517640	Agile Fully Covered Esophageal Stent System – Removable	18	23	120	10.5	180
M00517650	Agile Fully Covered Esophageal Stent System – Removable	18	23	150	10.5	180

Agile™ Partially Covered Esophageal Stent System



Order Number	Description	Stent OD (mm)	Proximal / Distal Flares OD (mm)	Stent Length (mm)	Delivery System Diameter (Fr)	Working Length (cm)
M00517200	Agile Partially Covered Esophageal Stent System	14	19	60	10.5	180
M00517210	Agile Partially Covered Esophageal Stent System	14	19	100	10.5	180
M00517220	Agile Partially Covered Esophageal Stent System	14	19	120	10.5	180
M00517230	Agile Partially Covered Esophageal Stent System	14	19	150	10.5	180
M00517240	Agile Partially Covered Esophageal Stent System	18	23	60	10.5	180
M00517250	Agile Partially Covered Esophageal Stent System	18	23	100	10.5	180
M00517260	Agile Partially Covered Esophageal Stent System	18	23	120	10.5	180
M00517270	Agile Partially Covered Esophageal Stent System	18	23	150	10.5	180

Agile™ OTW Fully Covered Esophageal Stent System – Removable



Order Number	Description	Stent OD (mm)	Proximal / Distal Flares OD (mm)	Stent Length (mm)	Delivery System Diameter (Fr)	Working Length (cm)
M00517790	Agile OTW Fully Covered Esophageal Stent System – Removable	23	28	60	18.5	78
M00517800	Agile OTW Fully Covered Esophageal Stent System – Removable	23	28	100	18.5	78
M00517810	Agile OTW Fully Covered Esophageal Stent System – Removable	23	28	120	18.5	78
M00517820	Agile OTW Fully Covered Esophageal Stent System – Removable	23	28	150	18.5	78



Agile™ OTW Partially Covered Esophageal Stent System



Order Number	Description	Stent OD (mm)	Proximal / Distal Flares OD (mm)	Stent Length (mm)	Delivery System Diameter (Fr)	Working Length (cm)
M00517710	Agile OTW Partially Covered Esophageal Stent System	23	28	60	18.5	78
M00517720	Agile OTW Partially Covered Esophageal Stent System	23	28	100	18.5	78
M00517730	Agile OTW Partially Covered Esophageal Stent System	23	28	120	18.5	78
M00517740	Agile OTW Partially Covered Esophageal Stent System	23	28	150	18.5	78

Esophageal and Enteral Metal Stents

COLONIC AND DUODENAL STENTS



WallFlex™ Colonic Stent System with Anchor Lock Delivery System



Order Number	Description	Diameter Flare / Body (mm)	Stent Length (cm)	Working Length (cm)	Catheter Diameter (F) / (mm)
M00565100	WallFlex Colonic Stent	27 / 22	6	230	10 / 3.33
M00565110	WallFlex Colonic Stent	27 / 22	9	230	10 / 3.33
M00565120	WallFlex Colonic Stent	27 / 22	12	230	10 / 3.33
M00565130	WallFlex Colonic Stent	27 / 22	6	135	10 / 3.33
M00565140	WallFlex Colonic Stent	27 / 22	9	135	10 / 3.33
M00565150	WallFlex Colonic Stent	27 / 22	12	135	10 / 3.33
M00565040	WallFlex Colonic Stent	30 / 25	6	230	10 / 3.33
M00565050	WallFlex Colonic Stent	30 / 25	9	230	10 / 3.33
M00565060	WallFlex Colonic Stent	30 / 25	12	230	10 / 3.33
M00565070	WallFlex Colonic Stent	30 / 25	6	135	10 / 3.33
M00565080	WallFlex Colonic Stent	30 / 25	9	135	10 / 3.33
M00565090	WallFlex Colonic Stent	30 / 25	12	135	10 / 3.33

Packaged one per box

Recommended Guidewire Systems 230cm: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 450cm: M00556161 Straight Dream Tip™, M00556171 Angled Dream Tip™ or .035 in (0.89mm) Jagwire™ Guidewire Stiff Shaft, 450cm: M00556601 Straight Tip, M00556611 Angled Tip or .035in (0.89mm) Super Stiff Wire, 500cm: H954180011

Systems 135cm: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm: M00556121 Straight Tip, M00556131 Angled Tip or .035in (0.89mm) Jagwire™ Guidewire Stiff Shaft, 260cm: M00557661 Straight Tip

WallFlex™ Colonic Soft Stent System with Anchor Lock Delivery System





Order Number	Description	Diameter Flare / Body (mm)	Stent Length (cm)	Working Length (cm)	Delivery System Diameter/ Working Channel (F/mm)	Guidewire Diameter (inch/mm)
M00555700	WallFlex Colonic Soft Stent	23 / 18	6	230	9 / 3.0	0.035 / 0.89
M00555710	WallFlex Colonic Soft Stent	23 / 18	9	230	9 / 3.0	0.035 / 0.89
M00555720	WallFlex Colonic Soft Stent	23 / 18	12	230	9 / 3.0	0.035 / 0.89
M00555730	WallFlex Colonic Soft Stent	25 / 20	6	230	9 / 3.0	0.035 / 0.89
M00555740	WallFlex Colonic Soft Stent	25 / 20	9	230	9 / 3.0	0.035 / 0.89
M00555750	WallFlex Colonic Soft Stent	25 / 20	12	230	9 / 3.0	0.035 / 0.89
M00555760	WallFlex Colonic Soft Stent	27 / 22	6	230	9 / 3.0	0.035 / 0.89
M00555770	WallFlex Colonic Soft Stent	27 / 22	9	230	9 / 3.0	0.035 / 0.89
M00555780	WallFlex Colonic Soft Stent	27 / 22	12	230	9 / 3.0	0.035 / 0.89



Esophageal and Enteral Metal Stents

COLONIC AND DUODENAL STENTS



WallFlex™ Duodenal Stent System with Anchor Lock Delivery System



Order Number	Description	Diameter Flare / Body (mm)	Stent Length (cm)	Working Length (cm)	Catheter Diameter (F) / (mm)
M00565010	WallFlex Duodenal Stent	27 / 22	6	230	10 / 3.33
M00565020	WallFlex Duodenal Stent	27 / 22	9	230	10 / 3.33
M00565030	WallFlex Duodenal Stent	27 / 22	12	230	10 / 3.33

Packaged one per box

Recommended Guidewires: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 450cm: M00556161 Straight Dream Tip™, M00556171 Angled Dream Tip™ or .035 in (0.89mm) Jagwire™ Guidewire Stiff Shaft, 450cm: M00556601 Straight Tip, M00556611 Angled Tip or .035in (0.89mm) Super Stiff Wire, 500cm: H954180011


WallFlex™ Duodenal Soft Stent System with Anchor Lock Delivery System



Order Number	Description	Diameter Flare / Body (mm)	Stent Length (cm)	Working Length (cm)	Delivery System Diameter/ Working Channel (F/mm)	Guidewire Diameter (inch/mm)
M00555800	WallFlex Duodenal Soft Stent	23 / 18	6	230	9 / 3.0	0.035 / 0.89
M00555810	WallFlex Duodenal Soft Stent	23 / 18	9	230	9 / 3.0	0.035 / 0.89
M00555820	WallFlex Duodenal Soft Stent	23 / 18	12	230	9 / 3.0	0.035 / 0.89
M00555830	WallFlex Duodenal Soft Stent	25 / 20	6	230	9 / 3.0	0.035 / 0.89
M00555840	WallFlex Duodenal Soft Stent	25 / 20	9	230	9 / 3.0	0.035 / 0.89
M00555850	WallFlex Duodenal Soft Stent	25 / 20	12	230	9 / 3.0	0.035 / 0.89
M00555860	WallFlex Duodenal Soft Stent	27 / 22	6	230	9 / 3.0	0.035 / 0.89
M00555870	WallFlex Duodenal Soft Stent	27 / 22	9	230	9 / 3.0	0.035 / 0.89
M00555880	WallFlex Duodenal Soft Stent	27 / 22	12	230	9 / 3.0	0.035 / 0.89



Infection Prevention



Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention
Endobariatrics



Contents



Alphabetical Index Product Code Index

IRRIGATION & INSUFFLATION

Hydra™ Irrigation System

Order Number	Description	Compatibility	Units
SIT-575-15	Hydra Irrigation Tubing	Olympus™ AFU-100, ERBE™ EIP 2, Medivators™ EndoStratus™ & EndoGator™	box 15
SIT-577-15	Hydra Irrigation Tubing	Olympus OFP2	box 15
SIT-361-100	⊗ Hydra Disposable Auxiliary Water Jet Connectors	Olympus 160, 180 & 190 Series	box 100
SCT-609-10	Hydra Water Bottle Cap with CO2	Air/Water Tubing Length: 16" (41 cm) CO2 Tubing Length: 18" (46 cm) Olympus	box 10

CLEANING BRUSHES

HedgeHog™ Cleaning Brushes

Order Number	Description	Units
SBS-227-50	HedgeHog Single-End Channel Brush	box 50
SBD-228-50	HedgeHog Double-End Channel Brush	box 50
SBD-289-50	HedgeHog Dual-End Channel/Valve Brush	box 50
SBD-291-50	HedgeHog Sweeper Brush	box 50
SBB-382-50	HedgeHog Double Bundle	box 50

VALVES

Orca™

Order Number	Description	Units
SUV-617-50	Orca Sterile, Single Use Air/Water and Suction Valves	box 50

OrcaPod™

Order Number	Description	Units
SUV-628-50	OrcaPod Sterile, Single Use Air/Water and Suction Valves with Seal Biopsy Valve	box 50
SUV-629-50	OrcaPod Sterile, Single Use Air/Water and Suction Valves, Seal Biopsy Valve and Hydra Water Jet Connector	box 50

Seal™

Order Number	Description	Units
SBC-501-100	Seal Single Use Biopsy Valve	box 100

Porpoise™

Order Number	Description	Units
M00501891	Porpoise™ Air/Water Cleaning Adapter (For use with Olympus 160, 180 and 190 gastrointestinal endoscopes)	box 50

SCOPE TRANSPORT

CinchPad™ Procedure & Endoscope Transport Pad

Order Number	Description	Units
SPX-017-25	CinchPad Endoscope Transport Pad	box 25

BITE BLOCKS

Blox™ Bite Blocks

Order Number	Description	Units
SBT-546-100	Blox Bite Block 60Fr with Dental Rim	box 100
SBT-114-100	Blox Bite Block 54Fr with Dental Rim	box 100
SBT-321-100	Blox Bite Block 54Fr without Dental Rim	box 100

SPONGES

Koala™ Sponge

Order Number	Description	Units
M00501762	Koala Cleaning Sponge	box 50
SES-240-50	Koala Cleaning Sponge Imbibed w/ Pure	box 50

Compliance Endokits

Order Number	Description	Units
CIK-115-20	OP4 I-KIT 1.1 YellowSponge	box 20
CIK-122-20	OP3 I-KIT 1.1 YellowSponge	box 20

LUBRICATING GELS

EndoGlide™ Lubricating Gels

Order Number	Description	Units
M00501902	EndoGlide+ 2.0 oz Lubricant Gel	box 50
M00501905	EndoGlide+ 1.1 oz Lubricant Gel	box 100

Endobariiatrics



Biliary
RX Biliary System™
Biliary Lumen Patency
Endoscopic Ultrasound (EUS)
Guidewires
Biopsy
Dilation
Imaging Systems
Pulmonary
Closure
Haemostasis
Tissue Resection
Esophageal and Enteral Metal Stents
Infection Prevention



Contents



Alphabetical Index

Product Code Index

ENDOSCOPIC SUTURING

OverStitch™ Systems



New / Old Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Units
M00505400 / ESS-G02-160***	OverStitch Endoscopic Suture System	*	-	box 1
M00505420 / ESS-G02-NXT***	OverStitch NXT Endoscopic Suture System	*	-	box 1
M00505441 / PLY-G02-020-APL**	OverStitch 2-0 Polypropylene Suture	*	2.8	box 12
M00505540 / CNH-G01-000**	OverStitch Suture Cinch	160	3.2	box 1
M00505470 / THX-165-028***	OverStitch Tissue Helix	*	2.8	box 1
M00505480 / NXT-THXP-130***	OverStitch NXT Tissue Helix Pro	165	2.8	box 1

* Suitable for endoscopes with max 110 cm working length.

** The UPN change will be effective starting April 18th.

*** The UPN change will be completed by the end of 2025.

INTRAGASTRIC BALLOON

BIB™ System



New / Old Order Number	Description	Units
M00505680 / B-40800*	BIB 6 Month Balloon	box 1

* The UPN change will be completed by the end of 2025.

Orbera365™



New / Old Order Number	Description	Units
M00505690 / B-50012*	Orbera365 12 Month Balloon	box 1
AS2290718	Endoline Viper	box 5
AS2640718	Endoline Punc	box 5

* The UPN change will be completed by the end of 2025.

OVERSTITCH ACCESSORIES

Overtube™



New / Old Order Number	Description	Units
M00505560 / OVT-027-160*	OverTube Endoscopic Access System	box 3

* The UPN change will be completed by the end of 2025.

Ensizor™




Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Units
ES26165-3	Ensizor Endoscopic Scissors Gastric Length	165	2.8	box 3
ES26235-3	Ensizor Endoscopic Scissors Colonic Length	235	2.8	box 3

INTRAGASTRIC BALLOON ACCESSORIES

Extraction Kit



Order Number	Description	Units
B-40815	Extraction Kit	box 5



**Boston
Scientific**
Advancing science for life™

www.bostonscientific.eu

© 2025 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-2174601-AA

All trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material is not intended for use in France.

Biliary	RX Biliary System™	Biliary Lumen Patency	Endoscopic Ultrasound (EUS)	Guidewires
Biopsy	Dilation	Imaging Systems	Pulmonary	Closure
Haemostasis	Tissue Resection	Esophageal and Enteral Metal Stents	Infection Prevention	Endobariatrics



Contents



Alphabetical
Index

Product Code
Index