

Endoscope Image Processor HDVS-S200D

- ✓ High definition resolution
- ✓ Multifunctional use



Disposable Imaging Catheter



Exclusively Distributed by Boston Scientific

Ordering information

Endoscopic image processing device



Model	HDVS-S200D
Size (mm)	320mm (W) * 231mm (D) * 77mm (H)
Weight (kg)	About 3.0
Output Interface	(DVI*1, SDI*2) (3)
Storage interface	USB*1
Features	White balance, freeze, zoom, LED intensity, image adjustment, photo, video, USB recording, electronic endoscope port

Single-use video cholangioscope




Model / specification	SCV-P-01G	SCV-P-02G
Outer Diameter of Distal End (mm)	3.70mm	2.95mm
Working Channel	2.00mm	1.20mm
Working length (mm)	2200mm	
Bending angle	≥50° in four directions	
Field of view	120°	
Depth of field (mm)	2-50mm	

Note:Indications, contraindications, warnings and instructions for use can be found on the product label or package insert. The products are for reference only and have not yet been approved for sale in some countries.
Name of registrant/manufacturer: Scivita Medical Technology (Suzhou) Co., Ltd.
© 2023 by Boston Scientific or its subsidiaries. All rights reserved.
Please refer to the package insert for contraindications or precautions.

All trademarks are the property of their respective owners. Photos taken by Scivita Medical Technology (Suzhou) Co.,Ltd
CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each Scivita Medical device. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.

CE 0123

Endoscopic image processing device



*Some features may require a USB drive for storage which is not included with device