



# FARAPULSE™

PFA Platform

## Product overview




## FARAPOINT™

Pulsed Field Ablation Catheter

### FARAPOINT Pulsed Field Ablation (PFA) Catheter

| Reference Number                            |                                   |                                                                                                                                               |
|---------------------------------------------|-----------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| M004PFCE81M310                              | FARAPOINT FF Curve PFA Catheter   |                                                                                                                                               |
| Catheter Specifications                     |                                   |                                                                                                                                               |
| Catheter type                               | Focal                             |                                                                                                                                               |
| Energy modality                             | PFA only                          |                                                                                                                                               |
| Shaft diameter                              | 8F                                |                                                                                                                                               |
| Usable length                               | 115 cm                            |                                                                                                                                               |
| Number of electrodes                        | 4                                 |                                                                                                                                               |
| Electrode spacing                           | 1.5 mm – 6.0 mm – 1.5 mm          |                                                                                                                                               |
| Electrode sizing                            | 2.0 mm – 1.5 mm – 1.5 mm – 2.0 mm |                                                                                                                                               |
| Tip diameter                                | 2.5 mm                            |                                                                                                                                               |
| Compatibility                               |                                   |                                                                                                                                               |
| Sheath                                      | 8.5F or larger                    | <i>Boston Scientific compatible:<br/>VersaCross™ Steerable Sheath (8.5F)<br/>SureFlex™ Steerable Guiding Sheath<br/>TorFlex™ Fixed Sheath</i> |
| Accessories                                 |                                   |                                                                                                                                               |
| PFA Rx Gen 2 Connection Cable for FARAPOINT | REF M004PF41M444                  |                                                                                                                                               |


**FARAPOINT is sterile and single-use only**



# FARASTAR™

Pulsed Field Ablation Generator

## FARASTAR Pulsed Field Ablation (PFA) Generator Ordering Information



| Reference Number                                                |                                                     |
|-----------------------------------------------------------------|-----------------------------------------------------|
| M004PF61M401                                                    | FARASTAR Pulsed Field Ablation Generator            |
| M004PF61M407                                                    | FARASTAR Recording System Module                    |
| Shipping Information – FARASTAR Pulsed Field Ablation Generator |                                                     |
| Weight                                                          | 552 lbs (250 kg)                                    |
| Dimensions                                                      | L 120 x H 200 x W 80 cm                             |
| Shipping Information – FARASTAR Recording System Module         |                                                     |
| Weight                                                          | 552 lbs (250 kg)                                    |
| Dimensions                                                      | L 120 x H 70 x W 80 cm                              |
| Technical Specifications                                        |                                                     |
| Input voltage                                                   | 120 V / 60 Hz                                       |
| Output voltage                                                  | Proprietary FARAPULSE biphasic, bipolar waveform    |
| Weight                                                          | 243 lbs (110 kg)                                    |
| Dimensions – FARASTAR Pulsed Field Ablation Generator           | L 79 x W 61 x H 117 cm (160 cm with monitor)        |
| Dimensions – FARASTAR Recording System Module                   | L 41 x W 19 x H 46 cm                               |
| User Interface                                                  | Articulating touchscreen                            |
| Accessories                                                     |                                                     |
| M004PF61M404                                                    | FARASTAR Stimulation Module Cable                   |
| M004PF61M405                                                    | FARASTAR 10-Pin EGM Cable                           |
| M004PF61M409                                                    | FARASTAR Stimulation Module Male Cable              |
| M004PF61M410                                                    | FARASTAR Stimulation Module Female Cable            |
| M004PF61M411                                                    | FARASTAR Stimulation Module Y-Cable Long            |
| M004PF61M412                                                    | FARASTAR Stimulation Module Y-Cable Short           |
| M004PF61M413                                                    | FARASTAR Recording System Module Catheter Pin Cable |
| M004PF61M415                                                    | FARASTAR Recording System Module ECG Trunk Cable    |
| M004PF61M416                                                    | FARASTAR Recording System Module ECG Output Module  |
| M004PF61M417                                                    | FARASTAR Recording System Module EGM Input Module   |
| M004PF61M606                                                    | FARASTAR ABL Cable MSM                              |
| M004PF61M618                                                    | FARASTAR 4-Pin EGM Cable FARAPPOINT                 |

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at [www.IFU-BSCI.com](http://www.IFU-BSCI.com). Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.

EP- 2293601-AA. FARAPULSE™ Pulsed Field Ablation Platform - C € 0459 ; VersaCross Connect - C € 1639.

**Boston  
Scientific**  
Advancing science for life™

[www.bostonscientific.eu](http://www.bostonscientific.eu)

© 2025 Boston Scientific Corporation  
or its affiliates. All rights reserved.