


RX Cytology Brush

Wireguided Cytology Brush

**Boston
Scientific**
Advancing science for life™



The **RX Cytology Brush Wireguided Cytology Brush** is an essential diagnostic tool that is part of a complete portfolio of devices for the treatment of pancreatobiliary diseases.

RX Cytology Brush


Wireguided Cytology Brush



The RX Cytology Brush Wireguided Cytology Brush is a double-lumen wireguided catheter for the collection of cells and the detection of malignant neoplasms in the biliary duct.

DEVICE FEATURES

- Handle has reinforced metal tube for reliable actuation.
- Distal radiopaque marker aids in catheter visualization under fluoroscopy for proper positioning.
- Open channel allows for quick exchange of the guidewire.



Brush lumen may be used to inject contrast dye through the proximal luer.

ORDERING INFORMATION

Order Number	Product Description	Brush Width (mm)	Catheter O.D. (F/mm)	Working Length (cm)	Recommended Guidewire (in/mm)	Units
M00545000	RX Cytology Brush	2.1	8/2.67	200	.035/0.89	Each

All trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.

**Boston
Scientific**

Advancing science for life™

www.bostonscientific.eu

© 2014 Boston Scientific Corporation or its affiliates. All rights reserved.

ENDO-253201-AA
CE 2797