



NRG™
Transseptal Needle



TECHNICAL SPECIFICATIONS & ORDERING INFORMATION

Features	Specifications
Needle curves	C0, C1
Needle lengths	71 cm, 89 cm, 98 cm
Dilator compatibility	0.032"
Needle gauge	Prox: 18 ga, Dist: 21 ga

Choose your needle configuration

Product Code	Needle Length	Curve	Compatible Transseptal Sheaths
M004NRG30010	71 cm	C0	8F or 8.5F fixed curve – 63 cm
M004NRG30030	71 cm	C1	
M004NRG30050	89 cm	C0	8.5F fixed curve – 81 cm
M004NRG30090	98 cm	C0	8.5F steerable curve – 71 cm
M004NRG30110	98 cm	C1	

