



**SureFlex™**  
Steerable Guiding Sheath



## TECHNICAL SPECIFICATIONS & ORDERING INFORMATION

Features	Specifications
French size compatibility	8.5F (2.8 mm)
Sheath usable length	71 cm
Sheath overall length	91 cm
Sheath radiopaque marker	Platinum/iridium band
Dilator usable length	94 cm
Guidewire length	180 cm
Compatible guidewire	0.032" (0.81 mm)
Bi-directional curve angles	90° CCW, 180° CW

Product Code	Distal Curve Diameter
TSK3001	Small (17 mm)
TSK3003	Medium (22 mm)
TSK3005	Large (50 mm)