




VersaCross™ RF Transseptal Solution



TECHNICAL SPECIFICATIONS

VersaCross™ Transseptal Sheath




Features	Specifications
Inner diameter	8.5F
Outer diameter	11F
Sheath usable length	63 cm
Dilator usable length	67 cm
Dilator Curves	D0 (Standard curve) D1 (Large curve)
Guidewire max outer diameter	0.035" (0.81 mm)

VersaCross™ RF Wire

Features	Specifications
Tip configuration	J-tip, Pigtail
RF wire length	180 cm, 230 cm
Wire diameter	0.035"

VersaCross™ Steerable Sheath


Features	Specifications
Inner diameter	8.5F
Outer diameter	12F
Sheath usable length	72 cm
Dilator usable length	95 cm
Sheath curves	Medium (22 mm) Large (50 mm)
Bi-directional curve angles	90° CCW, 180° CW
Guidewire max outer diameter	0.035" (0.81 mm)



ORDERING INFORMATION

VersaCross™ Transseptal Solution includes:


- **VersaCross™ RF Wire** (Pigtail or J-Tip)
- Single-use RFP 100A Connector Cable
- **VersaCross™ Transseptal Sheath**
- Dilator with **TRUform™ Shapeable Technology**
- 0.035" Guidewire



Product Code	Wire Tip	Wire Length	Sheath Length	Sheath Curve	Dilator Curve
VXSK3001	J	180 cm	63 cm	45°	D0
VXSK3011	J	230 cm	63 cm	45°	D0
VXSK3021	Pigtail	180 cm	63 cm	45°	D0
VXSK3031	Pigtail	230 cm	63 cm	45°	D0
VXSK3032	Pigtail	230 cm	63 cm	45°	D1
VXS85-35-63-45-D0	Single Fixed Sheath		63 cm	45°	D0
VXS85-35-63-45-D1	Single Fixed Sheath		63 cm	45°	D1

VersaCross™ Steerable Transseptal Solution includes:

- **VersaCross™ RF Wire** (Pigtail or J-Tip)
- Single-use RFP 100A Connector Cable
- **VersaCross™ Steerable Transseptal Sheath**
- Dilator with **TRUform™ Shapeable Technology**
- 0.035" Guidewire



Product Code	Wire Tip	Wire Length	Sheath Length	Curl	Dilator Curve
VSTK3003	J	180 cm	72 cm	M	D0
VSTK3010	J	230 cm	72 cm	M	D1
VSTK3015	Pigtail	180 cm	72 cm	M	D0
VSTK3019	Pigtail	230 cm	72 cm	L	D0
VSTK3021	Pigtail	230 cm	72 cm	M	D0
VSTK3022	Pigtail	230 cm	72 cm	M	D1
VST85-35-BD-71M-D0	Single Steerable Sheath		72 cm	M	D0
VST85-35-BD-71M-D1	Single Steerable Sheath		72 cm	M	D1

VersaCross RF Transseptal Solution compatible only with RFP-100A RF Puncture Generator.