

*The therapy referred to here may be one of several left heart procedures that require transseptal access. Additional procedural steps and/or devices may be required to deliver the therapy. Before use, consult Instructions for Use for any devices accordingly.

†POLARSHEATH™ Steerable Sheaths as shown in this brochure is sold separately.

‡The VersaCross Connect Transseptal Dilator is for use with a 12F (4.04 mm) ID POLARSHEATH™ Steerable Sheath that is 68 cm in length, specifically, model: M004 CRBS3050.


CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com.

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries.

This material not intended for use in France.



VersaCross Connect™
Access Solution for POLARSHEATH™



**Boston
Scientific**
Advancing science for life™

www.bostonscientific.eu

© 2025 Boston Scientific Corporation
or its affiliates. All rights reserved.

EP-1680505-AD C € 1639

Greater Confidence in Left Heart Access


ZERO EXCHANGE

POLARSHEATH™ Steerable Sheath delivery.
For use with the POLARx™ Cryoablation System.


SIMPLE, SEAMLESS, SMOOTH ZERO EXCHANGE

POLARSHEATH™ Steerable Sheath Delivery[†]


EXCHANGELESS **3-in-1** RF Wire



1 Start
easily track up the SVC



2 Access
with purpose-built RF




3 Deliver
with confidence


INSTANTLY DEPLOY
0.035" WIRE

To deliver **POLARSHEATH™**
Steerable Sheath
with confidence

DELIVER THERAPY ON TARGET



TRUform™
Shapeable Technology



Precise RF
Puncture Technology



VersaCross Connect™

Access Solution for POLARSHEATH™

Personalize your solution

1. Choose your VersaCross™ RF Wire



Length: 180 cm
Diameter: 0.035"

2. Connect with VersaCross Connect™ Transseptal Dilator



Inner diameter: 0.035"
Outer diameter: 12F (4.04 mm)
Useable length: 84 cm

Ordering information

Part. No.	Wire	Wire length	Dilator curve	Dilator length
VXAK3031	Pigtail	180 cm	D1	84 cm
VXAK3035	J-Tip	180 cm	D1	84 cm


OMNIVIZ™

Technology


Reliably locate the **VersaCross Connect™** solution on fluoroscopy and ultrasound. Positional markers indicate position of RF tip within dilator. Connect to DuoMode™ Cable to track and mark RF tip on mapping systems.

FAST TRACK


to your therapy delivery
in a **SINGLE SOLUTION**



VersaCross Connect™
Access Solution†



**Needle-based
workflow**




VersaCross Connect™

Access Solution for POLARSHEATH™


DESIGNED for use with the
POLARSHEATH™
Steerable Sheath†

**3-IN-1
RF WIRE**
for ZERO EXCHANGE
delivery



Smooth transition between
VersaCross Connect™ Dilator
and POLARSHEATH™ Sheath†

Snap-lock designed for
POLARSHEATH™ Sheath†



†POLARSHEATH™ Steerable Sheaths as shown in this brochure is sold separately.