



EMBOLD™ Detachable Coil System
SIMPLE BY DESIGN



FIBERED



PACKING



SOFT



**MULTI-CATHETER
COMPATIBILITY**

Accommodates .021" to
.027" ID microcatheters
for fewer exchanges and
complications.



**KINKLESS
PERFORMANCE**

Nitinol delivery system
to ensure kinkless
performance.



**HANDLE-FREE
DETACHMENT**

Provides complete control,
reducing procedure
steps and streamlining
workflow.



**BEST-IN-CLASS
OCCLUSION POWER**

Highly thrombogenic design
achieves faster occlusion
with fewer coils.

EMBOLD™ Detachable Coil System

| EMBOLD™ Fibered Coils | | | |
|-----------------------|-------------|------------------|----------------|
| Diameter (mm) | Length (cm) | UPN | GTIN |
| 2 | 4 | M001394240020040 | 08714729983521 |
| | 6 | M001394240020060 | 08714729983538 |
| 3 | 6 | M001394240030060 | 08714729983545 |
| | 12 | M001394240030120 | 08714729983552 |
| 4 | 8 | M001394240040080 | 08714729983569 |
| | 15 | M001394240040150 | 08714729983576 |
| 5 | 8 | M001394240050080 | 08714729983583 |
| | 15 | M001394240050150 | 08714729983590 |
| 6 | 10 | M001394240060100 | 08714729983606 |
| | 20 | M001394240060200 | 08714729983613 |
| 8 | 20 | M001394240080200 | 08714729983620 |
| | 20 | M001394240100200 | 08714729983637 |
| 10 | 30 | M001394240100300 | 08714729983644 |
| | 50 | M001394240100500 | 08714729983651 |
| 12 | 20 | M001394240120200 | 08714729983668 |
| | 30 | M001394240120300 | 08714729983675 |
| 14 | 50 | M001394240120500 | 08714729983682 |
| | 20 | M001394240140200 | 08714729983699 |
| 18 | 30 | M001394240140300 | 08714729983705 |
| | 50 | M001394240140500 | 08714729983712 |
| 20 | 50 | M001394240180500 | 08714729983729 |
| 22 | 50 | M001394240200500 | 08714729983736 |
| 22 | 60 | M001394240220600 | 08714729983743 |
| 32 | 60 | M001394240320600 | 00191506010485 |

| EMBOLD™ Soft Coils | | | |
|--------------------|-------------|------------------|----------------|
| Diameter (mm) | Length (cm) | UPN | GTIN |
| 2 | 2 | M001394250020020 | 00191506026363 |
| | 4 | M001394250020040 | 00191506026370 |
| | 6 | M001394250020060 | 00191506026387 |
| 3 | 8 | M001394250020080 | 00191506026394 |
| | 4 | M001394250030040 | 00191506026400 |
| | 8 | M001394250030080 | 00191506026417 |
| 4 | 12 | M001394250030120 | 00191506026424 |
| | 15 | M001394250030150 | 00191506026431 |
| | 8 | M001394250040080 | 00191506026448 |
| 5 | 10 | M001394250040100 | 00191506026455 |
| | 15 | M001394250040150 | 00191506026462 |
| | 30 | M001394250040300 | 00191506026479 |
| 6 | 10 | M001394250050100 | 00191506026486 |
| | 15 | M001394250050150 | 00191506026493 |
| | 20 | M001394250050200 | 00191506026509 |
| 8 | 30 | M001394250050300 | 00191506026516 |
| | 15 | M001394250060150 | 00191506026523 |
| | 20 | M001394250060200 | 00191506026530 |
| 6 | 30 | M001394250060300 | 00191506026547 |
| | 60 | M001394250060600 | 00191506026554 |
| | 30 | M001394250080300 | 00191506026561 |
| 8 | 60 | M001394250080600 | 00191506026578 |

| EMBOLD™ Packing Coils | | | |
|-----------------------|-------------|------------------|----------------|
| Diameter (mm) | Length (cm) | UPN | GTIN |
| N/A | 15 | M001394260000150 | 00191506026585 |
| | 30 | M001394260000300 | 00191506026592 |
| | 45 | M001394260000450 | 00191506026608 |
| | 60 | M001394260000600 | 00191506026615 |



www.bostonscientific.eu

© 2025 Boston Scientific Corporation or its affiliates. All rights reserved.

PI-1669010-AA

CE 0344

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.

All trademarks are the property of their respective owners.