


WALLSTENT™ Endoprosthesis

**PROVEN SOLUTIONS,
DEPENDABLE CHOICE**





WALLSTENT™ Endoprosthesis

The Most Comprehensive Range of Sizes

Diameters from 5 mm to 24 mm* and deployed lengths from 23 mm to 145 mm**

Reconstrainability

Designed for ease of placement, the WALLSTENT Endoprosthesis allows reconstrainability with the stent deployed up to the limit marker band

Visibility

Exclusive HALO™ Technology and platinum core provides excellent fluoroscopic visibility intended to facilitate precise positioning

Compression Resistance

Braided construction and Elgiloy™ Material designed to provide excellent compression resistance

Closed-Cell Design

Intended to provide increased scaffolding for optimal lesion coverage and a smooth inner lumen

Conformability

Excellent adaption to anatomical contours designed to provide exceptional stent-to-wall apposition

Flexibility

Highly flexible stent designed to smoothly cross lesions

*Available sizes vary per indication. See product information table for specific size availability.

**Approximate implanted stent length. Refer to DFU sizing chart for more information.

Product Information for WALLSTENT-Uni™ Endoprosthesis

Multi-Purpose Vascular Stent with Unistep™ Plus Delivery System

Short instrument												
Order Number	Fully Open Stent		Size of WALLSTENT-Uni™ implanted in Vessel (mm)						Delivery System			
	Diam. (mm)	Length (mm)	Vessel Diam.	Stent Length	Vessel Diam.	Stent Length	Vessel Diam.	Stent Length	Catheter Size (F)	Recomm. Introd. Sheath (F)	Recomm. Introd. Sheath (mm)	
M001731000	5	20	4	26	3	30	2	33	6	≤ 2.1	6	2.00
M001731020	5	40	4	54	3	63	2	68	6	≤ 2.1	6	2.00
M001731040	5	55	4	78	3	89	2	97	6	≤ 2.1	6	2.00
M001731060	5	80	4	116	3	134	2	145	6	≤ 2.1	6	2.00
M001731080	6	20	5	25	4	29	3	32	6	≤ 2.1	6	2.00
M001731100	6	45	5	52	4	60	3	65	6	≤ 2.1	6	2.00
M001731120	6	60	5	74	4	85	3	93	6	≤ 2.1	6	2.00
M001731140	6	90	5	111	4	128	3	139	6	≤ 2.1	6	2.00
M001731160	7	20	6	25	5	28	4	30	6	≤ 2.1	6	2.00
M001731180	7	40	6	50	5	57	4	63	6	≤ 2.1	6	2.00
M001731200	7	60	6	72	5	82	4	90	6	≤ 2.1	6	2.00
M001731220	7	90	6	108	5	123	4	134	6	≤ 2.1	6	2.00
M001731240	8	20	7	29	6	36	5	41	6	≤ 2.1	6	2.00
M001731260	8	40	7	49	6	56	5	61	6	≤ 2.1	6	2.00
M001731280	8	60	7	70	6	79	5	86	6	≤ 2.1	6	2.00
M001731300	8	80	7	105	6	119	5	130	6	≤ 2.1	6	2.00
M001731320	10	20	9	27	8	33	7	38	6	≤ 2.1	6	2.00
M001731340	10	42	9	48	8	54	7	58	7	≤ 2.4	7	2.33
M001731360	10	68	9	69	8	77	7	83	7	≤ 2.4	7	2.33
M001731380	10	94	9	103	8	115	7	124	7	≤ 2.4	7	2.33
M001731400	12	20	11	26	10	31	9	36	8	≤ 2.9	9	3.00
M001731420	12	40	11	47	10	51	9	55	8	≤ 2.9	9	3.00
M001731440	12	60	11	66	10	73	9	79	8	≤ 2.9	9	3.00
M001731460	12	90	11	100	10	110	9	119	8	≤ 2.9	9	3.00
M001731480	14	20	13	27	12	33	11	38	9	≤ 3.2	10	3.33
M001731490	14	40	13	46	12	50	11	54	9	≤ 3.2	10	3.33
M001731500	14	60	13	65	12	72	11	77	9	≤ 3.2	10	3.33
M001731510	14	90	13	98	12	107	11	115	9	≤ 3.2	10	3.33
M001731520	16	20	15	23	14	28	13	32	9	≤ 3.2	10	3.33
M001731530	16	40	15	45	14	49	13	52	9	≤ 3.2	10	3.33
M001731540	16	60	15	64	14	70	13	75	9	≤ 3.2	10	3.33
M001731550	16	90	15	97	14	105	13	112	9	≤ 3.2	10	3.33
M001731560	18	40	17	45	16	48	15	51	10	≤ 3.5	11	3.67
M001731570	18	60	17	64	16	69	15	73	10	≤ 3.5	11	3.67
M001731580	18	90	17	95	16	103	15	110	10	≤ 3.5	11	3.67
M001731590	20	40	19	40	18	44	17	47	10	≤ 3.5	11	3.67
M001731600	20	55	19	57	18	63	17	68	10	≤ 3.5	11	3.67
M001731610	20	80	19	86	18	94	17	102	10	≤ 3.5	11	3.67
M001731620	22	35	21	35	20	40	19	44	10	≤ 3.5	11	3.67
M001731630	22	45	21	50	20	57	19	62	10	≤ 3.5	11	3.67
M001731640	22	70	21	75	20	85	19	93	10	≤ 3.5	11	3.67
M001731650	24	35	23	35	22	39	21	43	11	≤ 3.9	12	4.00
M001731660	24	45	23	50	22	56	21	61	11	≤ 3.9	12	4.00
M001731670	24	70	23	75	22	84	21	92	11	≤ 3.9	12	4.00


Total instrument length 100 cm
Usable instrument length 75 cm
Recommended guidewire max. 0.035" (0.89 mm), 180 cm

Product Information for WALLSTENT-Uni™ Endoprosthesis

Multi-Purpose Vascular Stent with Unistep™ Plus Delivery System

Long instrument												
Order Number	Fully Open Stent		Size of WALLSTENT-Uni™ implanted in Vessel (mm)						Delivery System			
	Diam. (mm)	Length (mm)	Vessel Diam.	Stent Length	Vessel Diam.	Stent Length	Vessel Diam.	Stent Length	Catheter Size (F)	Recomm. Introd. Sheath (F)	Recomm. Introd. Sheath (mm)	
M001731010	5	20	4	26	3	30	2	33	6	≤ 2.1	6	2.00
M001731030	5	40	4	54	3	63	2	68	6	≤ 2.1	6	2.00
M001731050	5	55	4	78	3	89	2	97	6	≤ 2.1	6	2.00
M001731070	5	80	4	116	3	134	2	145	6	≤ 2.1	6	2.00
M001731090	6	20	5	25	4	29	3	32	6	≤ 2.1	6	2.00
M001731110	6	45	5	52	4	60	3	65	6	≤ 2.1	6	2.00
M001731130	6	60	5	74	4	85	3	93	6	≤ 2.1	6	2.00
M001731150	6	90	5	111	4	128	3	139	6	≤ 2.1	6	2.00
M001731170	7	20	6	25	5	28	4	30	6	≤ 2.1	6	2.00
M001731190	7	40	6	50	5	57	4	63	6	≤ 2.1	6	2.00
M001731210	7	60	6	72	5	82	4	90	6	≤ 2.1	6	2.00
M001731230	7	90	6	108	5	123	4	134	6	≤ 2.1	6	2.00
M001731250	8	20	7	29	6	36	5	41	6	≤ 2.1	6	2.00
M001731270	8	40	7	49	6	56	5	61	6	≤ 2.1	6	2.00
M001731290	8	60	7	70	6	79	5	86	6	≤ 2.1	6	2.00
M001731310	8	80	7	105	6	119	5	130	6	≤ 2.1	6	2.00
M001731330	10	20	9	27	8	33	7	38	6	≤ 2.1	6	2.00
M001731350	10	42	9	48	8	54	7	58	7	≤ 2.4	7	2.33
M001731370	10	68	9	69	8	77	7	83	7	≤ 2.4	7	2.33
M001731390	10	94	9	103	8	115	7	124	7	≤ 2.4	7	2.33
M001731410	12	20	11	26	10	31	9	36	8	≤ 2.9	9	3.00
M001731430	12	40	11	47	10	51	9	55	8	≤ 2.9	9	3.00
M001731450	12	60	11	66	10	73	9	79	8	≤ 2.9	9	3.00
M001731470	12	90	11	100	10	110	9	119	8	≤ 2.9	9	3.00

Total instrument length 160 cm
Usable instrument length 135 cm
Recommended guidewire max. 0.035" (0.89 mm), 260-300 cm



Boston Scientific
Advancing science for life™

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations.

PI-846001-AA Printed in Germany by medicalvision.

www.bostonscientific.eu

© 2020 Boston Scientific Corporation or its affiliates. All rights reserved.