

FETCH™ 2


Aspiration Catheter

Exceptional Deliverability.[§] Uncompromised Aspiration.[#]


A robust shaft and low crossing profile maximizes deliverability performance and robustness while maintaining excellent aspiration.

FETCH2 Compared to Export™ AP


64% more trackable* Distal End Flexibility Bench Test




48% more flexible† Tip Flexibility Bench Test



40% more kink resistant‡ Kink Resistance Bench Test



FETCH2's distal aspiration lumen cross-sectional area is maximized with a crescent shaped design.



DISTAL
Aspiration Lumen
Guidewire Lumen

[§] Fetch2 compared to Export™ AP. Bench test completed on Fetch2 (n = 5) and Export AP (n = 3). Data on file at BSC.

[#] In vitro testing done at BSC with 6F Export™ and Fetch2 devices. Data on file at BSC.

* Fetch™2 compared to Export™ AP. Bench test completed on Fetch™2 (n = 5) and Export™ AP (n = 3). Data on file at BSC.



† Tip flexibility bench test completed on Fetch2 (n = 5) Export AP (n = 3). Data on file at BSC.

‡ Kink resistance bench test completed on Fetch2 (n = 6) and Export AP (n = 6). Data on file at BSC.

** 1 lbf = 4.448222 N


Testing completed by Boston Scientific. Bench test results may not necessarily be indicative of clinical performance.

FETCH™ 2 Compared to Export™ Advance



Catheter Specifications ^{††}	FETCH2	Export AP	Export Advance
Working Length	135 cm	140 cm	140 cm
Hydrophilic Coated Distal Shaft	25 cm	27 cm	38 cm
Distal Outer Shaft Diameter	0.056" (1.4 mm)	0.068" (1.7 mm)	0.067" (1.7 mm)

FETCH2



[§] Pushability bench test completed on FETCH2 (n = 5) and Export Advance (n = 5).
^{***} Tip Flexibility bench test completed on FETCH2 (n = 5) and Export Advance (n = 5).
^{††} Information taken from FETCH2, Export AP, and Export Advance product labeling.

Testing completed by Boston Scientific. Data on file at BSC.
 Bench test results may not necessarily be indicative of clinical performance.
 All trademarks are the property of their respective owners.

Product Ordering Information

Product Name	UPN
FETCH2	109400-005

Boston Scientific
 Advancing science for life™

www.bostonscientific.eu

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations. This material is not approved for use or distribution in France.

IC-320105-AA June 2015 Printed in Germany by medicalvision.

© 2015 Boston Scientific Corporation or its affiliates. All rights reserved.
 DINCAR2653EA