

Window Indicator Band

Catheter Specifications


System Compatibility	Ultra
Vessel Diameter	≥ 6 mm
Working Length	105 cm
Shaft Diameter	8F (2.67 mm)
Double Marker Band	15 mm
Guidewire Compatibility	0.035" (0.89 mm) over-the-wire
Sheath Compatibility	8F (2.67 mm)
Power Pulse™	Yes

Order Information:

AngioJet™ ZelanteDVT™ Thrombectomy Set Ultra

GTIN (1 per box)	UPN (1 per box)
08714729904731	114610-002

For use only with the AngioJet Ultra System



ANGIOJET™ ZELANTEDVT™ Thrombectomy Catheter

DOMINATE DVTs



¹Presented by Dr. Mark Garcia at CIRSE 2013; Final PEARL Data Aug 2013

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations.

Indications, operating specifications and availability may vary by country. Check with local product representation and country-specific Information For Use for your country. This material is not approved for use or distribution in France.

PI-340601-AA SEP 2015 Printed in Germany by medicalvision.

ZelanteDVT™ provides the power to remove large venous clot burdens

The AngioJet™ ZelanteDVT thrombectomy catheter is specifically designed to treat deep vein thrombosis (DVT) in large-diameter upper and lower peripheral veins. ZelanteDVT is the most powerful thrombectomy catheter in the market-leading AngioJet portfolio. The catheter makes it possible to more efficiently remove large thrombus burden in veins and restore blood flow, giving you the power to dominate DVTs.

Four times the thrombus-removal power of Solent™ Omni and Solent Proxi catheters

Power Pulse™ enabled for infusion of physician-specified fluids, including thrombolytic agents

Single inflow window for torqueable and directional thrombectomy power



The AngioJet™ Thrombectomy System provides the power and flexibility to remove venous thrombus and restore flow in the most challenging of DVT cases. Recent PEARL Registry data showed¹:

38% of DVT treatments completed in less than 6 hours;
75% completed in less than 24 hours

No catheter-directed thrombolysis (CDT) needed in **40%** of venous cases; Lytic delivered by AngioJet (including Power Pulse™) in **86%** of venous cases

87% of AngioJet venous cases were completed in two or fewer sessions

84% of venous-treated patients had freedom from rethrombosis following 365 days

¹AngioJet ZelanteDVT was not studied in the PEARL Registry.

