

Encore™ 40 Inflator



- Inflates **up to 40 ATMs**.
- **25 cc capacity**, complete with pressure gauge.
- Large printed pressure gauge that glows in the dark for **easy reading in dark labs**.
- Quick latch mechanism allows for **rapid deflation**, even from max pressure.
- Clear barrel provides a **precise view of the content and volume** allowing user to perform thorough purging.

Order Number	Description	Units per box
M001394472540	25CC 40ATM Inflation Device	1

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Material not intended for use in France.

PI-901901-AA © 2020 Boston Scientific Corporation or its affiliates. All rights reserved.