


ATHLETIS™ Ultra-High-Pressure Balloon

LOWER COMPLIANCE, HIGHER PERFORMANCE

Dilatation Efficiency is a measure of a PTA balloon's ability to open resistant lesions at the lowest pressure possible.

Athletis achieves virtually the same balloon dilatation using lower pressure and remaining more non-compliant than Conquest 40.




All products 8 mm balloons.
Average measurements taken by Boston Scientific (n = 15). Data on file (TR 92200478).
Bench test results may not necessarily be indicative of clinical performance.

Athletis' proprietary braided balloon design contributes to its high performance and lower compliance compared to Conquest 40.


44% MORE
NON-COMPLIANT

Lower Compliance Compared to Conquest 40 (Bench Test)



*Data from Conquest 40 brochure.
All products 8 mm balloons.

Average measurements taken by Boston Scientific (n = 15). Data on file (TR 92200478).
Bench test results may not necessarily be indicative of clinical performance.





www.bostonscientific.eu

ATHLETIS

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Material not intended for use in France.

Bench Test results may not necessarily be indicative of clinical performance.

Athletis is a registered or unregistered trademark of Boston Scientific Corporation or its affiliates. All other trademarks are property of their respective owners.

© 2020 Boston Scientific Corporation
or its affiliates. All rights reserved.

PI-873404-AA