

MUSTANG™

0.035" Balloon Dilatation Catheter

Product Codes

UPN Number	Catalog Number	Sizes (mm x mm x cm)	Qty	UPN Number	Catalog Number	Sizes (mm x mm x cm)	Qty
H74939171030240	3917103024	3.0 x 20 x 40		H74939171040440	3917104044	4.0 x 40 x 40	
H74939171030270	3917103027	3.0 x 20 x 75		H74939171040470	3917104047	4.0 x 40 x 75	
H74939171030210	3917103021	3.0 x 20 x 135		H74939171040410	3917104041	4.0 x 40 x 135	
H74939171030340	3917103034	3.0 x 30 x 40		H74939171040640	3917104064	4.0 x 60 x 40	
H74939171030370	3917103037	3.0 x 30 x 75		H74939171040670	3917104067	4.0 x 60 x 75	
H74939171030310	3917103031	3.0 x 30 x 135		H74939171040610	3917104061	4.0 x 60 x 135	
H74939171030440	3917103044	3.0 x 40 x 40		H74939171040840	3917104084	4.0 x 80 x 40	
H74939171030470	3917103047	3.0 x 40 x 75		H74939171040870	3917104087	4.0 x 80 x 75	
H74939171030410	3917103041	3.0 x 40 x 135		H74939171040810	3917104081	4.0 x 80 x 135	
H74939171030640	3917103064	3.0 x 60 x 40		H74939171041040	3917104104	4.0 x 100 x 40	
H74939171030670	3917103067	3.0 x 60 x 75		H74939171041070	3917104107	4.0 x 100 x 75	
H74939171030610	3917103061	3.0 x 60 x 135		H74939171041010	3917104101	4.0 x 100 x 135	
H74939171030840	3917103084	3.0 x 80 x 40		H74939171041270	3917104127	4.0 x 120 x 75	
H74939171030870	3917103087	3.0 x 80 x 75		H74939171041210	3917104121	4.0 x 120 x 135	
H74939171030810	3917103081	3.0 x 80 x 135		H74939171041570	3917104157	4.0 x 150 x 75	
H74939171031040	3917103104	3.0 x 100 x 40		H74939171041510	3917104151	4.0 x 150 x 135	
H74939171031070	3917103107	3.0 x 100 x 75		H74939171041870	3917104187	4.0 x 180 x 75	
H74939171031010	3917103101	3.0 x 100 x 135		H74939171041810	3917104181	4.0 x 180 x 135	
H74939171031270	3917103127	3.0 x 120 x 75		H74939171042070	3917104207	4.0 x 200 x 75	
H74939171031210	3917103121	3.0 x 120 x 135		H74939171042010	3917104201	4.0 x 200 x 135	
H74939171031570	3917103157	3.0 x 150 x 75		H74939171050240	3917105024	5.0 x 20 x 40	
H74939171031510	3917103151	3.0 x 150 x 135		H74939171050270	3917105027	5.0 x 20 x 75	
H74939171031870	3917103187	3.0 x 180 x 75		H74939171050210	3917105021	5.0 x 20 x 135	
H74939171031810	3917103181	3.0 x 180 x 135		H74939171050340	3917105034	5.0 x 30 x 40	
H74939171032070	3917103207	3.0 x 200 x 75		H74939171050370	3917105037	5.0 x 30 x 75	
H74939171032010	3917103201	3.0 x 200 x 135		H74939171050310	3917105031	5.0 x 30 x 135	
H74939171040240	3917104024	4.0 x 20 x 40		H74939171050440	3917105044	5.0 x 40 x 40	
H74939171040270	3917104027	4.0 x 20 x 75		H74939171050470	3917105047	5.0 x 40 x 75	
H74939171040210	3917104021	4.0 x 20 x 135		H74939171050410	3917105041	5.0 x 40 x 135	
H74939171040340	3917104034	4.0 x 30 x 40		H74939171050640	3917105064	5.0 x 60 x 40	
H74939171040370	3917104037	4.0 x 30 x 75		H74939171050670	3917105067	5.0 x 60 x 75	
H74939171040310	3917104031	4.0 x 30 x 135		H74939171050610	3917105061	5.0 x 60 x 135	

Excellence on every dimension.

Mustang™ 0.035" Balloon Catheter Product Codes

UPN Number	Catalog Number	Sizes (mm x mm x cm)	Qty	UPN Number	Catalog Number	Sizes (mm x mm x cm)	Qty
H74939171050840	3917105084	5.0 x 80 x 40		H74939171070340	3917107034	7.0 x 30 x 40	
H74939171050870	3917105087	5.0 x 80 x 75		H74939171070370	3917107037	7.0 x 30 x 75	
H74939171050810	3917105081	5.0 x 80 x 135		H74939171070310	3917107031	7.0 x 30 x 135	
H74939171051040	3917105104	5.0 x 100 x 40		H74939171070440	3917107044	7.0 x 40 x 40	
H74939171051070	3917105107	5.0 x 100 x 75		H74939171070470	3917107047	7.0 x 40 x 75	
H74939171051010	3917105101	5.0 x 100 x 135		H74939171070410	3917107041	7.0 x 40 x 135	
H74939171051240	3917105124	5.0 x 120 x 40		H74939171070640	3917107064	7.0 x 60 x 40	
H74939171051270	3917105127	5.0 x 120 x 75		H74939171070670	3917107067	7.0 x 60 x 75	
H74939171051210	3917105121	5.0 x 120 x 135		H74939171070610	3917107061	7.0 x 60 x 135	
H74939171051570	3917105157	5.0 x 150 x 75		H74939171070840	3917107084	7.0 x 80 x 40	
H74939171051510	3917105151	5.0 x 150 x 135		H74939171070870	3917107087	7.0 x 80 x 75	
H74939171051870	3917105187	5.0 x 180 x 75		H74939171070810	3917107081	7.0 x 80 x 135	
H74939171051810	3917105181	5.0 x 180 x 135		H74939171071040	3917107104	7.0 x 100 x 40	
H74939171052070	3917105207	5.0 x 200 x 75		H74939171071070	3917107107	7.0 x 100 x 75	
H74939171052010	3917105201	5.0 x 200 x 135		H74939171071010	3917107101	7.0 x 100 x 135	
H74939171060240	3917106024	6.0 x 20 x 40		H74939171071270	3917107127	7.0 x 120 x 75	
H74939171060270	3917106027	6.0 x 20 x 75		H74939171071210	3917107121	7.0 x 120 x 135	
H74939171060210	3917106021	6.0 x 20 x 135		H74939171071570	3917107157	7.0 x 150 x 75	
H74939171060340	3917106034	6.0 x 30 x 40		H74939171071510	3917107151	7.0 x 150 x 135	
H74939171060370	3917106037	6.0 x 30 x 75		H74939171071870	3917107187	7.0 x 180 x 75	
H74939171060310	3917106031	6.0 x 30 x 135		H74939171071810	3917107181	7.0 x 180 x 135	
H74939171060440	3917106044	6.0 x 40 x 40		H74939171072070	3917107207	7.0 x 200 x 75	
H74939171060470	3917106047	6.0 x 40 x 75		H74939171072010	3917107201	7.0 x 200 x 135	
H74939171060410	3917106041	6.0 x 40 x 135		H74939171080240	3917108024	8.0 x 20 x 40	
H74939171060640	3917106064	6.0 x 60 x 40		H74939171080270	3917108027	8.0 x 20 x 75	
H74939171060670	3917106067	6.0 x 60 x 75		H74939171080210	3917108021	8.0 x 20 x 135	
H74939171060610	3917106061	6.0 x 60 x 135		H74939171080340	3917108034	8.0 x 30 x 40	
H74939171060840	3917106084	6.0 x 80 x 40		H74939171080370	3917108037	8.0 x 30 x 75	
H74939171060870	3917106087	6.0 x 80 x 75		H74939171080310	3917108031	8.0 x 30 x 135	
H74939171060810	3917106081	6.0 x 80 x 135		H74939171080440	3917108044	8.0 x 40 x 40	
H74939171061040	3917106104	6.0 x 100 x 40		H74939171080470	3917108047	8.0 x 40 x 75	
H74939171061070	3917106107	6.0 x 100 x 75		H74939171080410	3917108041	8.0 x 40 x 135	
H74939171061010	3917106101	6.0 x 100 x 135		H74939171080640	3917108064	8.0 x 60 x 40	
H74939171061240	3917106124	6.0 x 120 x 40		H74939171080670	3917108067	8.0 x 60 x 75	
H74939171061270	3917106127	6.0 x 120 x 75		H74939171080610	3917108061	8.0 x 60 x 135	
H74939171061210	3917106121	6.0 x 120 x 135		H74939171080840	3917108084	8.0 x 80 x 40	
H74939171061570	3917106157	6.0 x 150 x 75		H74939171080870	3917108087	8.0 x 80 x 75	
H74939171061510	3917106151	6.0 x 150 x 135		H74939171080810	3917108081	8.0 x 80 x 135	
H74939171061870	3917106187	6.0 x 180 x 75		H74939171081040	3917108104	8.0 x 100 x 40	
H74939171061810	3917106181	6.0 x 180 x 135		H74939171081070	3917108107	8.0 x 100 x 75	
H74939171062070	3917106207	6.0 x 200 x 75		H74939171081010	3917108101	8.0 x 100 x 135	
H74939171062010	3917106201	6.0 x 200 x 135		H74939171081270	3917108127	8.0 x 120 x 75	
H74939171070240	3917107024	7.0 x 20 x 40		H74939171081210	3917108121	8.0 x 120 x 135	
H74939171070270	3917107027	7.0 x 20 x 75					
H74939171070210	3917107021	7.0 x 20 x 135					

UPN Number	Catalog Number	Sizes (mm x mm x cm)	Qty	UPN Number	Catalog Number	Sizes (mm x mm x cm)	Qty
H74939171081570	3917108157	8.0 x 150 x 75		H74939171100340	3917110034	10.0 x 30 x 40	
H74939171081510	3917108151	8.0 x 150 x 135		H74939171100370	3917110037	10.0 x 30 x 75	
H74939171081870	3917108187	8.0 x 180 x 75		H74939171100310	3917110031	10.0 x 30 x 135	
H74939171081810	3917108181	8.0 x 180 x 135		H74939171100440	3917110044	10.0 x 40 x 40	
H74939171082070	3917108207	8.0 x 200 x 75		H74939171100470	3917110047	10.0 x 40 x 75	
H74939171082010	3917108201	8.0 x 200 x 135		H74939171100410	3917110041	10.0 x 40 x 135	
H74939171090240	3917109024	9.0 x 20 x 40		H74939171100640	3917110064	10.0 x 60 x 40	
H74939171090270	3917109027	9.0 x 20 x 75		H74939171100670	3917110067	10.0 x 60 x 75	
H74939171090210	3917109021	9.0 x 20 x 135		H74939171100610	3917110061	10.0 x 60 x 135	
H74939171090340	3917109034	9.0 x 30 x 40		H74939171100840	3917110084	10.0 x 80 x 40	
H74939171090370	3917109037	9.0 x 30 x 75		H74939171100870	3917110087	10.0 x 80 x 75	
H74939171090310	3917109031	9.0 x 30 x 135		H74939171100810	3917110081	10.0 x 80 x 135	
H74939171090440	3917109044	9.0 x 40 x 40		H74939171120240	3917112024	12.0 x 20 x 40	
H74939171090470	3917109047	9.0 x 40 x 75		H74939171120270	3917112027	12.0 x 20 x 75	
H74939171090410	3917109041	9.0 x 40 x 135		H74939171120210	3917112021	12.0 x 20 x 135	
H74939171090640	3917109064	9.0 x 60 x 40		H74939171120340	3917112034	12.0 x 30 x 40	
H74939171090670	3917109067	9.0 x 60 x 75		H74939171120370	3917112037	12.0 x 30 x 75	
H74939171090610	3917109061	9.0 x 60 x 135		H74939171120310	3917112031	12.0 x 30 x 135	
H74939171090840	3917109084	9.0 x 80 x 40		H74939171120440	3917112044	12.0 x 40 x 40	
H74939171090870	3917109087	9.0 x 80 x 75		H74939171120470	3917112047	12.0 x 40 x 75	
H74939171090810	3917109081	9.0 x 80 x 135		H74939171120410	3917112041	12.0 x 40 x 135	
H74939171100240	3917110024	10.0 x 20 x 40		H74939171120640	3917112064	12.0 x 60 x 40	
H74939171100270	3917110027	10.0 x 20 x 75		H74939171120670	3917112067	12.0 x 60 x 75	
H74939171100210	3917110021	10.0 x 20 x 135		H74939171120610	3917112061	12.0 x 60 x 135	
				H74939171120840	3917112084	12.0 x 80 x 40	
				H74939171120870	3917112087	12.0 x 80 x 75	
				H74939171120810	3917112081	12.0 x 80 x 135	

Physician's Signature

**Powerful dilatation combined
with superior deliverability
and a complete size matrix.**

MUSTANG™ BALLOON DILATATION CATHETER

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Directions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

INDICATIONS FOR USE: The Mustang Balloon Dilatation Catheter is indicated for Percutaneous Transluminal Angioplasty (PTA) in the peripheral vasculature, including iliac, femoral, popliteal, tibial, peroneal, subclavian, and renal arteries and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. The Mustang Balloon Dilatation Catheter is also indicated for post-dilatation of balloon expandable and self-expanding stents in the peripheral vasculature. Mustang Balloon Dilatation Catheters with balloons up to 120 mm in length are indicated for the treatment of biliary strictures.

CONTRAINDICATIONS: None Known.

WARNINGS: The Mustang PTA balloon catheter is not for use in the coronary arteries. Any use for procedures other than those indicated in these instructions is not recommended.

PRECAUTIONS: The Mustang PTA Balloon Dilatation Catheter shall only be used by physicians trained in the performance of percutaneous transluminal angioplasty. The Mustang PTA Balloon Dilatation Catheter should be used with caution for procedures involving calcified lesions or synthetic vascular grafts due to the abrasive nature of these inflation sites. The Mustang PTA Balloon Dilatation Catheters are not intended for injection of contrast medium. Precautions to prevent or reduce clotting should be taken when any catheter is used: • Consider systemic anticoagulation • Flush or rinse all products with sterile saline or a similar solution prior to use.

ADVERSE EVENTS: The complications that may result from a balloon dilatation procedure include, but are not limited to: • Allergic reaction (device, contrast medium and medications) • Arteriovenous fistula • Embolization (air, device, plaque, etc.) • Hematoma • Hemorrhage, including bleeding at puncture site • Pseudoaneurysm • Sepsis/infection • Thromboembolic episodes • Vessel injury, e.g. dissection, perforation, rupture • Vessel occlusion • Vessel spasm

Additionally, the following complications may result from a balloon dilatation procedure within the biliary system: • Cholangitis • Hemobilia • Pancreatitis • Duct injury, e.g. dissection, perforation, rupture • Duct occlusion • Duct spasm

Mustang is a registered or unregistered trademark of Boston Scientific Corporation or its affiliates.

**Boston
Scientific**

Advancing science for life™

Peripheral Interventions

300 Boston Scientific Way
Marlborough, MA 01752 USA

www.bostonscientific.com

To order product or for more information
contact customer service at 1.888.272.1001.

© 2014 Boston Scientific Corporation
or its affiliates. All rights reserved.

PI-238002-AA NOV2014