

Expel™ Drainage and Introducer Products

Product Catalog and Ordering Information

C-Codes*

C1729	Drainage Catheter
C2617	Stent, non-coronary, temporary, no delivery system
C1894	AccuStick™ System

* Boston Scientific is not responsible for the correct use of codes on submitted claims; this information does not constitute reimbursement or legal advice.

Conversion Chart

French (F)	mm
6	2
8	2.67
10	3.33
12	4
14	4.67
16	5.33
18	6
24	8
26	8.67

Order Number Description

Catheter Type
Catheter Length
39310-08200
Expel Product Family
French Size
0 = single
K1 = kit

Catheter Type


MPD Singles (J-Tip)	10
MPDL Singles (Loop)	11
MPDL Tight Loops	
MPDL Loop Kits	
MPD Large Diameter	20
Nephrostomy Singles	41
Nephrostomy Kits	
Internal / External	51
Ureteral Stents	61
Ureteral Stents Non-Standard Length	

Boston Scientific
Advancing science for life™

EXPEL™ Drainage and Introducer Products

**EXCEPTIONAL DELIVERABILITY,
FLEXIBILITY AND SECURITY**

PRODUCT CATALOG AND ORDERING INFORMATION



Boston Scientific
Advancing science for life™

Expel, Flexima, Twist Loc, Accustick II and Percufix are registered or unregistered trademarks of Boston Scientific Corporation or its affiliates.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations.

Information contained herein is for distribution outside the U.S. only. Illustrations for information purposes – not indicative of actual size or clinical outcome. This material is not approved for use or distribution in France.

PI-356601-AB Dec 2015 Printed in Germany by medicalvision.

www.bostonscientific.eu

© 2015 Boston Scientific Corporation or its affiliates. All rights reserved. DINPER4650EA

FORGET COMPROMISE. YOU CAN *HAVE IT ALL* WITH BOSTON SCIENTIFIC'S DRAINAGE PORTFOLIO.


Built upon 25 years of percutaneous drainage experience, Boston Scientific's Expel Drainage Catheter is resetting the standard. Experience the unprecedented buckling resistance and flexibility of our new, proprietary Flexithane™ Material and the simplicity of our tamper-resistant Twist-Loc™ Hub.

Flexithane Catheter Material

- Best-in-class buckling resistance and most balanced performance[†]
- Smooth tapered tip protected within the pigtail
- Hydrophilic coating applied to distal portion of catheter
- Highly radiopaque Flexithane Catheter Material

Visual Markers[‡]

- Radiopaque marker band denoting the location of the proximal-most drainage hole
- Printed length markings



Unlocked Hub



Locked Hub

Twist-Loc Hub[§]

- Simple, secure locking mechanism with triple confirmation:
 - 1 Tactile twist
 - 2 Audible "click"
 - 3 Visual color-block indicator
- Smooth, low profile design

[†]Twist-Loc Hub is only available on locking loop catheters

[‡]All Expel catheters and stents are equipped with radiopaque markers

[§]Bench testing performed by Boston Scientific Corporation. Data on file. N=3. Bench test results may not necessarily be indicative of clinical performance.

Expel™ Multi-Purpose Drainage Catheters

Boston Scientific Multi-Purpose Drainage Catheters are designed for drainage of abscesses and fluid collections.

EXPUL MULTI-PURPOSE DRAINAGE CATHETERS* FLEXITHANE™ MATERIAL

ORDER NUMBER	DESCRIPTION Catheter Type / French Size / Length (cm)
--------------	--

Catheters

39310-0620	MPD / 6 / 20
39310-0820	MPD / 8 / 20
39310-1020	MPD / 10 / 20
39310-1225	MPD / 12 / 25
39310-1425	MPD / 14 / 25
39311-0825	MPDLoop / 8 / 25
39311-1025	MPDLoop / 10 / 25
39311-1225	MPDLoop / 12 / 25
39311-1425	MPDLoop / 14 / 25

Kits

39311-0825K1	MPDLoop / 8 / 25
39311-1025K1	MPDLoop / 10 / 25

Each Kit Contains: 21 Gauge x 15 Diagnostic Needle; 0.018" Floppy Tip Guidewire; AccuStick™ II Introducer; 0.035" J Tip Guidewire; Dilators; Expel™ Locking Loop Catheter; Percufix™ Catheter Cuff Kit; Stopcock; 30 cm Connecting Tube.



MPD Loop
(with Locking Pigtail)



MPD
(with Non-locking J-tip)

EXPUL MULTI-PURPOSE DRAINAGE CATHETERS* FLEXITHANE™ MATERIAL

ORDER NUMBER	DESCRIPTION Catheter Type / French Size / Length (cm)
--------------	--

SPECIALTY SIZES

Short Length, Tight Loop Catheters

39311-0615	MPDLoop TL / 6 / 15
39311-0815	MPDLoop TL / 8 / 15
39311-1015	MPDLoop TL / 10 / 15

Long Length Catheters

39311-0845	MPDLoop / 8 / 45
39311-1045	MPDLoop / 10 / 45
39311-1245	MPDLoop / 12 / 45
39311-1445	MPDLoop / 14 / 45

*All catheters packaged with both metal and flexible stiffening cannulae.

Large Capacity Catheters

39320-1630	MPD / 16 / 30
39320-1830	MPD / 18 / 30
39320-2440	MPD / 24 / 40
39320-2640	MPD / 26 / 40

Plastic cannula, metal cannula, trocar (except 24 F and 26 F), plug, end cap, connector tube.



Expel Catheter
Tight Loop Tip



Large Capacity Catheter

Expel Urinary Drainage Catheters

Boston Scientific Urinary Drainage Catheters and stents are designed to provide internal and external drainage of the urinary tract.

EXPUL NEPHROSTOMY CATHETERS* FLEXITHANE™ MATERIAL

ORDER NUMBER	DESCRIPTION Catheter Type / French Size / Length (cm)
--------------	--

Catheters

39341-0825	Neph / 8 / 25
39341-1025	Neph / 10 / 25
39341-1225	Neph / 12 / 25
39341-1425	Neph / 14 / 25

Kits

39341-0825K1	Neph / 8 / 25
39341-1025K1	Neph / 10 / 25
39341-1225K1	Neph / 12 / 25

Each Kit Contains: 21 Gauge x 15 Diagnostic Needle; 0.018" Floppy Tip Guidewire; AccuStick™ II Introducer; 0.0385" J Tip Guidewire; Dilators; Expel™ Locking Loop Catheter; Percufix™ Catheter Cuff Kit; Stopcock; 30 cm Connecting Tube.

*All catheters packaged with both metal and flexible stiffening cannulae.



EXPUL NEPHROURETERAL STENTS* FLEXITHANE™ MATERIAL

ORDER NUMBER	DESCRIPTION Stent Type / French Size / Length (cm)
--------------	---

Internal/External (I/E) Stents

39351-0822	I/E / 8 / 22
39351-0824	I/E / 8 / 24
39351-0826	I/E / 8 / 26
39351-0828	I/E / 8 / 28
39351-1022	I/E / 10 / 22
39351-1024	I/E / 10 / 24
39351-1026	I/E / 10 / 26
39351-1028	I/E / 10 / 28

*All stents packaged with both metal and flexible stiffening cannulae.



EXPUL URETERAL STENTS* FLEXITHANE™ MATERIAL

ORDER NUMBER	DESCRIPTION Stent Type / French Size / Length (cm)
--------------	---


Stents

39361-0620	US / 6 / 20
39361-0622	US / 6 / 22
39361-0624	US / 6 / 24
39361-0626	US / 6 / 26
39361-0628	US / 6 / 28
39361-0820	US / 8 / 20
39361-0822	US / 8 / 22
39361-0824	US / 8 / 24
39361-0826	US / 8 / 26
39361-0828	US / 8 / 28
39361-1020	US / 10 / 20
39361-1022	US / 10 / 22
39361-1024	US / 10 / 24
39361-1026	US / 10 / 26
39361-1028	US / 10 / 28

Short Length Stents

39361-0612	US / 6 / 12
39361-0614	US / 6 / 14
39361-0616	US / 6 / 16
39361-0812	US / 8 / 12
39361-0814	US / 8 / 14
39361-0816	US / 8 / 16
39361-1012	US / 10 / 12
39361-1014	US / 10 / 14
39361-1016	US / 10 / 16

*All stents packaged with flexible stiffening cannula (except 6.3 F) and stabilizer.



Biliary Drainage Catheters

Boston Scientific Biliary Drainage Catheters are designed to facilitate drainage of the biliary tree.

FLEXIMA BILIARY DRAINAGE CATHETERS* WITH RO MARKER AND TEMPTIP HYDROPHILIC TIP

Includes radiopaque marker and also dissolvable tip at distal end of catheter to facilitate drainage flow.

ORDER NUMBER	DESCRIPTION Catheter Type / French Size / Length (cm)
M001282600	Biliary / 8 / 35
M001282610	Biliary / 10 / 35

*All catheters packaged with both metal and flexible stiffening cannulae.



Flexima Biliary Drainage Catheter

FLEXIMA BILIARY DRAINAGE CATHETERS* WITH RO MARKER

Special highly radiopaque marker designed to provide clear view of catheter location under fluorocopy for proper placement in the biliary tract.

ORDER NUMBER	DESCRIPTION Catheter Type / French Size / Length (cm) / Durometer
--------------	--

Catheters

M001272600	Biliary / 8 / 35 / Regular
M001272610	Biliary / 10 / 35 / Regular
M001272620	Biliary / 12 / 35 / Regular
M001272630	Biliary / 14 / 35 / Regular
M001272640	Biliary / 8 / 35 / Firm
M001272650	Biliary / 10 / 35 / Firm
M001272660	Biliary / 10 / 35 / Soft
M001272670	Biliary / 10 / 35 / Soft

*All catheters packaged with both metal and flexible stiffening cannulae.



Flexima Biliary RO Marker

FLEXIMA BILIARY DRAINAGE CATHETERS*

Large drainage holes designed to facilitate drainage of thick viscous bile.

ORDER NUMBER	DESCRIPTION Catheter Type / French Size / Length (cm) / Durometer
--------------	--

Catheters

M001271560	Biliary / 8 / 35 / Regular
M001271570	Biliary / 10 / 35 / Regular
M001271580	Biliary / 12 / 35 / Regular
M001271590	Biliary / 14 / 35 / Regular
M001271540	Biliary / 8 / 35 / Firm
M001271550	Biliary / 10 / 35 / Firm
M001271600	Biliary / 10 / 35 / Soft
M001271620	Biliary / 14 / 35 / Soft

Kits

M001271660	Biliary / 10 / 35 / Regular
M001271630	Biliary / 8 / 35 / Firm
M001271640	Biliary / 10 / 35 / Firm

Each Kit Contains: 21 Gauge x 15 Diagnostic Needle; 0.018" Floppy Tip Guidewire; AccuStick II Introducer; 0.035" J Tip Guidewire; Dilators; Expel LockingLoop Catheter; Percufix Catheter Cuff Kit; Stopcock; 30 cm Connecting Tube.

FLEXIMA BILIARY DRAINAGE CATHETERS* WITH TEMPTIP HYDROPHILIC TIP

Dissolvable tip at distal end of catheter designed to create a large opening to facilitate drainage.

ORDER NUMBER	DESCRIPTION Catheter Type / French Size / Length (cm) / Durometer
--------------	--

Catheters

M001281560	Biliary / 8 / 35 / Regular
M001281570	Biliary / 10 / 35 / Regular
M001281580	Biliary / 12 / 35 / Regular
M001281590	Biliary / 14 / 35 / Regular



Flexima Biliary Hub

Non-Vascular Introducers and Accessories

Boston Scientific Offers a Comprehensive Portfolio of Non-Vascular Accessories to complement our drainage catheter and stent offerings.

ACCUSTICK™ II WITH RO MARKER INTRODUCER SYSTEMS

For gaining initial drainage access and facilitating placement of 0.038" (1.016 mm) working guidewire alongside an 0.018" (0.508 mm) guidewire.

ORDER NUMBER	DESCRIPTION
20-702	AccuStick II with 0.038" J Tip Wire
20-703	AccuStick II with 0.038" Straight Wire
20-704	AccuStick II without 0.038" Wire
20-705	AccuStick II without Guidewires
20-710	AccuStick II with 0.038" Stainless Steel J Tip and 0.018" Nitinol Guidewire
20-108	0.018", 60 cm, 7.5 cm Floppy Tip Guidewire
20-604	21 Gauge 15 cm Diagnostic Needle (1ea.)
39370-21201	21 Gauge 20 cm Diagnostic Needle (10pk)

Each AccuStick II System Contains: 21 Gauge x 15 cm Diagnostic Needle; 0.018" Guidewire; AccuStick II Sheath/Dilator.



CONNECTING TUBE FOR DRAINAGE BAG

Designed to provide a simple, secure connection to standard drainage fluid collection bags.

ORDER NUMBER	DESCRIPTION
22-910	CTube, 30 cm
39380-00151	Large Capacity CTube, 15 cm (5pk)

Each Kit Contains: 30 cm Connecting Tube; Stopcock.




FASCIAL DILATORS

(Accept 0.038" (1.016 mm) Guidewires).

ORDER NUMBER	DESCRIPTION Fascial Dilator / French Size / Length (cm)
48-250	TD / 10 / 20

Each Kit Contains: 30 cm Connecting Tube; Stopcock.



PERCUFIX™ CATHETER CUFF KIT

Designed for secure, comfortable fixation of catheter to skin site after catheter placement.

ORDER NUMBER	DESCRIPTION
20-700	PCC

Each Kit Contains: Percufix Catheter Cuff; Two Cable Ties; 6" x 8" (152.4 x 203.2 mm) Sheet of Transparent Moisture Vapor Permeable Dressing.

