

# Expel™ Drainage and Introducer Products

Product Catalog and Ordering Information

## Conversion Chart

French (F)	mm
6	2
8	2.67
10	3.33
12	4
14	4.67
16	5.33
18	6
24	8
26	8.67

## Order Number Description

Catheter Type
Catheter Length  
**39310-08200**  
Expel Product Family
French Size
0 = single  
K1 = kit

### Catheter Type


MPD Singles (J-Tip)	10
MPDL Singles (Loop)	11
MPDL Tight Loops	
MPDL Loop Kits	20
MPD Large Diameter	
Nephrostomy Singles	41
Nephrostomy Kits	
Internal / External	51
Ureteral Stents	61
Ureteral Stents Non-Standard Length	

**Boston Scientific**  
Advancing science for life™

**EXPEL™** Drainage and Introducer Products

**EXCEPTIONAL DELIVERABILITY,  
FLEXIBILITY AND SECURITY**

**PRODUCT CATALOG AND ORDERING INFORMATION**



\*All catheters and stents are packaged with both metal stiffening cannula and an optional plastic stiffener. Pigtail includes a pigtail straightener (except for 14F). MPDL includes an optional trocar for direct stick puncture.

Expel, Flexima, Twist Loc, Accustick II and Percufix are registered or unregistered trademarks of Boston Scientific Corporation or its affiliates.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations.

Information contained herein is for distribution outside the U.S. only. Illustrations for information purposes – not indicative of actual size or clinical outcome. This material is not approved for use or distribution in France.

PI-356601-AD May 2016 Printed in Germany by medicalvision.

**Boston Scientific**  
Advancing science for life™

[www.bostonscientific.eu](http://www.bostonscientific.eu)

© 2016 Boston Scientific Corporation or its affiliates. All rights reserved. DINPER4650EA



# FORGET COMPROMISE. YOU CAN *HAVE IT ALL* WITH BOSTON SCIENTIFIC'S DRAINAGE PORTFOLIO.


*Built upon 25 years of percutaneous drainage experience, Boston Scientific's Expel Drainage Catheter is resetting the standard. Experience the unprecedented buckling resistance and flexibility of our new, proprietary Flexithane™ Material and the simplicity of our tamper-resistant Twist-Loc™ Hub.*

## ***Flexithane Catheter Material***

- Best-in-class buckling resistance and most balanced performance†
- Smooth tapered tip protected within the pigtail
- Hydrophilic coating applied to distal portion of catheter
- Highly radiopaque Flexithane Catheter Material

## ***Visual Markers‡***

- Radiopaque marker band denoting the location of the proximal-most drainage hole
- Printed length markings



Unlocked Hub



Locked Hub

## ***Twist-Loc Hub§***

- Simple, secure locking mechanism with triple confirmation:
  - 1 Tactile twist
  - 2 Audible "click"
  - 3 Visual color-block indicator
- Smooth, low profile design

† Bench testing performed by Boston Scientific Corporation. Data on file. N=3. Bench test results may not necessarily be indicative of clinical performance.

‡ All Expel catheters and stents are equipped with radiopaque markers.

§ Twist-Loc Hub is only available on locking loop catheters.

# Expel™ Multi-Purpose Drainage Catheters

Boston Scientific Multi-Purpose Drainage Catheters are designed for drainage of abscesses and fluid collections.

## EXPET MULTI-PURPOSE DRAINAGE CATHETERS\*

### FLEXITHANE™ MATERIAL

ORDER NUMBER	DESCRIPTION Catheter Type / French Size / Length (cm)
--------------	--

#### Catheters

H7493931006200	MPD / 6 / 20
H7493931008200	MPD / 8 / 20
H7493931010200	MPD / 10 / 20
H7493931012250	MPD / 12 / 25
H7493931014250	MPD / 14 / 25
H7493931108250	MPDLoop / 8 / 25
H7493931110250	MPDLoop / 10 / 25
H7493931112250	MPDLoop / 12 / 25
H7493931114250	MPDLoop / 14 / 25

#### Kits

H749393110825K1	MPDLoop / 8 / 25
H749393111025K1	MPDLoop / 10 / 25

**Each Kit Contains:** 21 Gauge x 15 Diagnostic Needle; 0.018" Floppy Tip Guidewire; AccuStick™ II Introducer; 0.035" J Tip Guidewire; Dilators; Expel™ Locking Loop Catheter; Percufix™ Catheter Cuff Kit; Stopcock; 30 cm Connecting Tube.



MPD Loop  
(with Locking Pigtail)



MPD  
(with Non-locking J-tip)

## EXPET MULTI-PURPOSE DRAINAGE CATHETERS\*

### FLEXITHANE™ MATERIAL

ORDER NUMBER	DESCRIPTION Catheter Type / French Size / Length (cm)
--------------	--

#### SPECIALTY SIZES

##### Short Length, Tight Loop Catheters

H7493931106150	MPDLoop TL / 6 / 15
H7493931108150	MPDLoop TL / 8 / 15
H7493931110150	MPDLoop TL / 10 / 15

##### Long Length Catheters

H7493931108450	MPDLoop / 8 / 45
H7493931110450	MPDLoop / 10 / 45
H7493931112450	MPDLoop / 12 / 45
H7493931114450	MPDLoop / 14 / 45

##### Large Capacity Catheters<sup>^</sup>

H7493932016300	MPD / 16 / 30
H7493932018300	MPD / 18 / 30
H7493932024400	MPD / 24 / 40
H7493932026400	MPD / 26 / 40

<sup>^</sup> Plastic cannula, metal cannula, trocar (except 24 F and 26 F), plug, end cap, connector tube.



Expel Catheter  
Tight Loop Tip



Large Capacity Catheter

# Expel Urinary Drainage Catheters

Boston Scientific Urinary Drainage Catheters and stents are designed to provide internal and external drainage of the urinary tract.

## EXPET NEPHROSTOMY CATHETERS\*

### FLEXITHANE™ MATERIAL

ORDER NUMBER	DESCRIPTION Catheter Type / French Size / Length (cm)
--------------	--

#### Catheters

H7493934108250	Neph / 8 / 25
H7493934110250	Neph / 10 / 25
H7493934112250	Neph / 12 / 25
H7493934114250	Neph / 14 / 25

#### Kits

H749393410825K1	Neph / 8 / 25
H749393411025K1	Neph / 10 / 25
H749393411225K1	Neph / 12 / 25

**Each Kit Contains:** 21 Gauge x 15 Diagnostic Needle; 0.018" Floppy Tip Guidewire; AccuStick™ II Introducer; 0.035" J Tip Guidewire; Dilators; Expel™ Locking Loop Catheter; Percufix™ Catheter Cuff Kit; Stopcock; 30 cm Connecting Tube.



## EXPET NEPHROURETERAL STENTS\*

### FLEXITHANE™ MATERIAL

ORDER NUMBER	DESCRIPTION Catheter Type / French Size / Length (cm)
--------------	--

#### Internal/External (I/E) Stents

H7493935108220	I/E / 8 / 22
H7493935108240	I/E / 8 / 24
H7493935108260	I/E / 8 / 26
H7493935108280	I/E / 8 / 28
H7493935110220	I/E / 10 / 22
H7493935110240	I/E / 10 / 24
H7493935110260	I/E / 10 / 26
H7493935110280	I/E / 10 / 28



## EXPET URETERAL STENTS\*\*

### FLEXITHANE™ MATERIAL

ORDER NUMBER	DESCRIPTION Catheter Type / French Size / Length (cm)
--------------	--


#### Stents

H7493936108200	US / 6 / 20
H7493936108220	US / 6 / 22
H7493936108240	US / 6 / 24
H7493936108260	US / 6 / 26
H7493936108280	US / 6 / 26
H7493936108200	US / 8 / 20
H7493936108220	US / 8 / 22
H7493936108240	US / 8 / 24
H7493936108260	US / 8 / 26
H7493936108280	US / 8 / 28
H7493936110200	US / 10 / 20
H7493936110220	US / 10 / 22
H7493936110240	US / 10 / 24
H7493936110260	US / 10 / 26
H7493936110280	US / 10 / 28

#### Short Length Stents

H7493936106120	US / 6 / 12
H7493936106140	US / 6 / 14
H7493936106160	US / 6 / 16
H7493936108120	US / 8 / 12
H7493936108140	US / 8 / 14
H7493936108160	US / 8 / 16
H7493936110120	US / 10 / 12
H7493936110140	US / 10 / 14
H7493936110160	US / 10 / 16

\*\* All stents packaged with flexible stiffening cannula (except 6.3 F) and stabilizer.



## Biliary Drainage Catheters

Boston Scientific Biliary Drainage Catheters are designed to facilitate drainage of the biliary tree.

### FLEXIMA BILIARY DRAINAGE CATHETERS\* WITH RO MARKER AND TEMPTIP HYDROPHILIC TIP

Includes radiopaque marker and also dissolvable tip at distal end of catheter to facilitate drainage flow.

ORDER NUMBER	DESCRIPTION Catheter Type / French Size / Length (cm)
<b>Catheters</b>	
M001282600	Biliary / 8 / 35
M001282610	Biliary / 10 / 35



Flexima Biliary Drainage Catheter

### FLEXIMA BILIARY DRAINAGE CATHETERS\* WITH RO MARKER

Special highly radiopaque marker designed to provide clear view of catheter location under fluoroscopy for proper placement in the biliary tract.

ORDER NUMBER	DESCRIPTION Catheter Type / French Size / Length (cm) / Durometer
<b>Catheters</b>	
M001272600	Biliary / 8 / 35 / Regular
M001272610	Biliary / 10 / 35 / Regular
M001272620	Biliary / 12 / 35 / Regular
M001272630	Biliary / 14 / 35 / Regular
M001272640	Biliary / 8 / 35 / Firm
M001272650	Biliary / 10 / 35 / Firm
M001272660	Biliary / 10 / 35 / Soft
M001272670	Biliary / 12 / 35 / Soft



Flexima Biliary RO Marker

### FLEXIMA BILIARY DRAINAGE CATHETERS\*

Large drainage holes designed to facilitate drainage of thick viscous bile.

ORDER NUMBER	DESCRIPTION Catheter Type / French Size / Length (cm) / Durometer
<b>Catheters</b>	
M001271560	Biliary / 8 / 35 / Regular
M001271570	Biliary / 10 / 35 / Regular
M001271580	Biliary / 12 / 35 / Regular
M001271590	Biliary / 14 / 35 / Regular
M001271540	Biliary / 8 / 35 / Firm
M001271550	Biliary / 10 / 35 / Firm
M001271600	Biliary / 10 / 35 / Soft
M001271620	Biliary / 14 / 35 / Soft
<b>Kits</b>	
M001271660	Biliary / 10 / 35 / Regular
M001271630	Biliary / 8 / 35 / Firm
M001271640	Biliary / 10 / 35 / Firm



Flexima Biliary Hub

**Each Kit Contains:** 21 Gauge x 15 Diagnostic Needle; 0.018" Floppy Tip Guidewire; AccuStick II Introducer; 0.035" J Tip Guidewire; Dilators; Expel LockingLoop Catheter; Percufix Catheter Cuff Kit; Stopcock; 30 cm Connecting Tube.

### FLEXIMA BILIARY DRAINAGE CATHETERS\* WITH TEMPTIP HYDROPHILIC TIP

Dissolvable tip at distal end of catheter designed to create a large opening to facilitate drainage.

ORDER NUMBER	DESCRIPTION Catheter Type / French Size / Length (cm) / Durometer
<b>Catheters</b>	
M001281560	Biliary / 8 / 35 / Regular
M001281570	Biliary / 10 / 35 / Regular
M001281580	Biliary / 12 / 35 / Regular
M001281590	Biliary / 14 / 35 / Regular

## Non-Vascular Introducers and Accessories

Boston Scientific Offers a Comprehensive Portfolio of Non-Vascular Accessories to complement our drainage catheter and stent offerings.

### ACCUSTICK™ II WITH RO MARKER INTRODUCER SYSTEMS

For gaining initial drainage access and facilitating placement of 0.038" (1.016 mm) working guidewire alongside an 0.018" (0.508 mm) guidewire.

ORDER NUMBER	DESCRIPTION
M001207020	AccuStick II with 0.038" (1.016 mm) J Tip Wire
M001207030	AccuStick II with 0.038" (1.016 mm) Straight Wire
M001207040	AccuStick II without 0.038" (1.016 mm) Wire
M001207050	AccuStick II without Guidewires
M001207100	AccuStick II with 0.038" (1.016 mm) Stainless Steel J Tip and 0.018" (0.508 mm) Nitinol Guidewire
M001201081	0.018" (0.508 mm), 60 cm, 7.5 cm Floppy Tip Guidewire
M001206041	21 Gauge (0.723 mm) Diagnostic Needle



**Each AccuStick II System Contains:** 21 Gauge x 15 cm Diagnostic Needle; 0.018" Guidewire; AccuStick II Sheath/Dilator.

### PERCUFIX™ CATHETER CUFF KIT

Designed for secure, comfortable fixation of catheter to skin site after catheter placement.

ORDER NUMBER	DESCRIPTION
M001207001	PCC

**Each Kit Contains:** Percufix Catheter Cuff; Two Cable Ties; 6" x 8" (152.4 x 203.2 mm) Sheet of Transparent Moisture Vapor Permeable Dressing.



### CONNECTING TUBE FOR DRAINAGE BAG

Designed to provide a simple, secure connection to standard drainage fluid collection bags.

ORDER NUMBER	DESCRIPTION
M001229101	CTube, 30 cm
H7493937000151	Large Capacity CTube, 15 cm

**Each Kit Contains:** 30 cm Connecting Tube; Stopcock.



### FASCIAL DILATORS

(Accept 0.038" (1.016 mm) Guidewires).

ORDER NUMBER	DESCRIPTION Fascial Dilator / French Size / Length (cm)
M001482501	TD / 10 / 20

**Each Kit Contains:** 30 cm Connecting Tube; Stopcock.

