



## EXPEL™

Drainage Catheter with Twist-Loc™ Hub

# EXCEPTIONAL DELIVERABILITY, FLEXIBILITY AND SECURITY



The top four performance attributes are plotted on each axis with the competitive performance average displayed as the white box. Lines falling outside the white box depict performance better than average.†

\*Compression + Bending = Flexibility


†Bench testing performed by Boston Scientific Corporation. Data on file. N=3. Bench test results may not necessarily be indicative of clinical performance.

## Flexithane™ Catheter Material

- Best-in-class buckling resistance and most balanced performance†
- Smooth tapered tip protected within the pigtail
- Hydrophilic coating applied to distal portion of catheter
- Highly radiopaque Flexithane™ Catheter Material

## Visual Markers‡

- Radiopaque marker band denoting the location of the proximal-most drainage hole
- Printed length markings



Radiopaque Marker Band

## Twist-Loc Hub§

- Simple, secure locking mechanism with triple confirmation:

- 1 Tactile twist
- 2 Audible “click”
- 3 Visual color-block indicator

- Smooth, low profile design



Unlocked Hub



Locked Hub

## Broad Product Portfolio

- A variety of catheters and stents
- Additional lengths
- Large diameter catheters



**EXPEL DRAINAGE CATHETER: DELIVERABLE. FLEXIBLE. SECURE.**

†Bench testing performed by Boston Scientific Corporation. Data on file. N=3. Bench test results may not necessarily be indicative of clinical performance.

‡All Expel catheters and stents are equipped with radiopaque makers

§Twist-Loc Hub is only available on locking loop catheters

All cited trademarks are the property of their respective owners. AccuStick, ESP, Expel, Flexithane, Percufix and Twist-Loc are registered or unregistered trademarks of Boston Scientific Corporation or its affiliates.

**CAUTION:** The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Material not for use or distribution in France.

PI-357501-AA NOV 2015. Printed in The Netherlands by DeBudelse/creativeservices.