

AMS700™
Three-Piece Penile Implant

AMS' MOST POPULAR INFLATABLE PROSTHESIS – THE AMS 700™ PENILE PROSTHESIS PRODUCT LINE

AMS 700™ has an overall
patient and partner
satisfaction rate of **92%**
and **96%** respectively¹

Available with or without preconnected cylinders and pump for either an infrapubic or penoscrotal approach.
Offering a combination of features focused on both surgeon and patient satisfaction.

EXCELLENT GIRTH EXPANSION ▲ EXCEPTIONAL RIGIDITY ▲ NATURAL APPEARING FLACCIDITY

EASE OF IMPLANT EASE OF USE

Conceal™ Low-Profile Reservoir

- ▲ Parylene coated to enhance reservoir durability²
- ▲ Optimizes fit to offer a sleek, flexible solution

Cylinder and Rear Tip Extender (RTE) Design

- ▲ Reduced proximal profile designed to ease placement into proximal corpora even in difficult to dilate cases²
- ▲ Snap-Fit™ RTEs optimize fit and provide a secure, rigid connection²
- ▲ Broad range of RTE lengths provide custom fit for patients (available in 0.5 cm increments, ranging from 0.5-7.5 cm)

Momentary Squeeze (MS) Pump

- ▲ One-touch button designed for easy deflation²
- ▲ Lock-out valve designed to resist auto inflation²

InhibiZone™ Antibiotic Treatment

- ▲ Only antibiotic impregnated inflatable penile prosthesis on the market³
- ▲ Creates a zone of inhibition effective against the bacteria commonly associated with inflatable penile prosthesis infections⁴
- ▲ Provides a significant reduction in rate of revisions due to infection⁴

IMPORTANT SAFETY INFORMATION⁴

- ▲ Not for patients with urogenital infections or active skin infections in region of surgery.
- ▲ InhibiZone™ treated version of this device is not for patients with known allergy or sensitivity to rifampin (rifampicin), or to minocycline, or other tetracyclines.
- ▲ Will make latent natural or spontaneous erections, as well as other interventional treatment options impossible. Penile shortening, curvature or scarring may result.
- ▲ Patients with diabetes, spinal cord injuries, or open sores may have an increased risk of infection associated with a prosthesis.


¹ F. Montorsi, et al. AMS Three-Piece Inflatable Implants for Erectile Dysfunction: A Long-Term Multi-Institution Study in 200 Consecutive Patients: Eur Urol. 2000; 37: 50-55.

² Data on file at AMS

³ FDA/PDP D970012-S065

⁴ AMS 700™ with MS Pump™ Penile Prosthesis Product Line Instructions for Use. American Medical Systems. 2013.

**AMS
700™ MS Series**



AMS 700™ MS LGX INFLATABLE PENILE PROSTHESIS

▲ AMS Exclusive: Cylinders designed to optimize Length and Girth eXpansion with up to 20% elongation possible, depending on patients' anatomy¹

PRODUCT	INHIBIZONE™ ITEM NUMBER
LGX Preconnect MS 12 cm, PS, IZ	72404250
LGX Preconnect MS 15 cm, PS, IZ	72404251
LGX Preconnect MS 18 cm, PS, IZ	72404252
LGX Preconnect MS 21 cm, PS, IZ	72404253
LGX Preconnect MS 12 cm, IP, IZ	72404255
LGX Preconnect MS 15 cm, IP, IZ	72404256
LGX Preconnect MS 18 cm, IP, IZ	72404257
LGX Preconnect MS 21 cm, IP, IZ	72404258

AMS 700™ MS CXR INFLATABLE PENILE PROSTHESIS

▲ All the benefits of an IPP including girth expansion, for patients that require smaller cylinders for those with scarred or sclerotic corpora¹

PRODUCT	INHIBIZONE™ ITEM NUMBER
CXR Preconnect MS, 10 cm PS, IZ	72404260
CXR Preconnect MS, 12 cm PS, IZ	72404261
CXR Preconnect MS, 14 cm PS, IZ	72404262
CXR Preconnect MS, 16 cm PS, IZ	72404263
CXR Preconnect MS, 18 cm PS, IZ	72404264
CXR Preconnect MS, 10 cm IP, IZ	72404265
CXR Preconnect MS, 12 cm IP, IZ	72404266
CXR Preconnect MS, 14 cm IP, IZ	72404267
CXR Preconnect MS, 16 cm IP, IZ	72404268
CXR Preconnect MS, 18 cm IP, IZ	72404269

AMS 700™ MS CX INFLATABLE PENILE PROSTHESIS

▲ Controlled eXpansion cylinders that optimize girth and offers additional length options.

PRODUCT	INHIBIZONE™ ITEM NUMBER
CX Preconnect MS 12 cm, PS, IZ	72404230
CX Preconnect MS 15 cm, PS, IZ	72404231
CX Preconnect MS 18 cm, PS, IZ	72404232
CX Preconnect MS 21 cm, PS, IZ	72404233
CX Preconnect MS 24 cm, PS, IZ	72404234*
CX Preconnect MS 12 cm, IP, IZ	72404235
CX Preconnect MS 15 cm, IP, IZ	72404236
CX Preconnect MS 18 cm, IP, IZ	72404237
CX Preconnect MS 21 cm, IP, IZ	72404238
CX Preconnect MS 24 cm, IP, IZ	72404239*

*Available as special order

RESERVOIRS

PRODUCT	INHIBIZONE™ ITEM NUMBER
Conceal™ Reservoir	720185-01
Reservoir, 65 mL, PC/IZ	72404155
Reservoir, 100 mL, PC/IZ	72404156

*All parts are also available without InhibiZone™ Antibiotic Treatment

Summary of Warnings and Contraindications

The AMS 700™ Series Inflatable Penile Prosthesis is intended for use in the treatment of chronic, organic, male erectile dysfunction (impotence). These devices are contraindicated in patients who have active urogenital infections or active skin infections in the region of surgery or (for the AMS 700 with InhibiZone™) have a known sensitivity or allergy to rifampin, minocycline, or other tetracyclines). Implantation will make latent natural or spontaneous erections, as well as other interventional treatment options, impossible. Men with diabetes, spinal cord injuries or open sores may have an increased risk of infection. Failure to evaluate and treat device erosion may result in infection and loss of tissue. Implantation may result in penile shortening, curvature, or scarring. Possible adverse events include, but are not limited to, urogenital pain (usually associated with healing), urogenital edema, urogenital ecchymosis, urogenital erythema, reservoir encapsulation, patient dissatisfaction, auto-inflation, mechanical malfunction, and impaired urination.

Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions and potential adverse events.



Rx Only

AMS 700™ MS Series, InhibiZone™, Snap-Fit™ and Conceal™ are trademarks of American Medical Systems, Inc.
© 2013 American Medical Systems, Inc. All Rights Reserved. Minnetonka, MN 55343
AMSUS/700-00230/ June 2013 www.AmericanMedicalSystems.com 1-800-328-3881
AU, CA, NZ and U.S. Use Only