

# Ambicor™

## Inflatable Penile Prosthesis



## Two-piece inflatable penile prosthesis

The AMS Ambicor prosthesis has an overall patient and partner **satisfaction rate** of **96%** and **91%**, respectively<sup>1</sup>

### ABOUT THE DEVICE

- High pressure cylinder designed to provide support for erectile dysfunction (ED) patients with ED
- 2-piece inflatable penile prosthesis designed for ease of placement
  - No separate reservoir
  - Pre-filled and pre-connected implant
- Stackable, snap-fit RTEs for customised fit
- Requires few squeezes to inflate cylinders<sup>2</sup>
- High patient and partner satisfaction<sup>1</sup>
- One-step bend deflation
- **Diameter:** 12.5mm, 14mm, 15.5mm
- **Length:** 14cm, 16cm, 18cm, 20cm or 22cm
- **Rear Tip Extenders (RTEs):** add length to proximal cylinder end – RTEs are packaged with the Ambicor device and include one pair each of 0.5 cm, 1 cm, 2 cm and 3 cm


### IMPORTANT SAFETY INFORMATION<sup>3</sup>

- Not for patients with urogenital infections or active skin infections in region of surgery
- Will make latent natural erections, as well as other interventional treatment options impossible
- Penile shortening, curvature or scarring may result
- Patients with diabetes, spinal cord injuries or open sores may have an increased risk of infection associated with a prosthesis

\* An extra set of RTEs is available, separate from the Ambicor Prosthesis packaging, that includes one pair each of 0.5 cm, 1 cm, 2 cm and 3 cm


## DIMENSIONAL SPECIFICATIONS

| Cylinder diameter    | Cylinder length | Product number |
|----------------------|-----------------|----------------|
| 12.5mm               | 14cm            | 72401450       |
| 12.5mm               | 16cm            | 72401451       |
| 12.5mm               | 18cm            | 72401452       |
| 14mm                 | 16cm            | 72401453       |
| 14mm                 | 18cm            | 72401454       |
| 14mm                 | 20cm            | 72401455       |
| 15.5mm               | 18cm            | 72401456       |
| 15.5mm               | 20cm            | 72401457       |
| 15.5mm               | 22cm            | 72401458       |
| Snap-Fit RTE Package |                 | 72402890       |



## AMS AMBICOR INFLATABLE PENILE PROSTHESIS IS HIGHLY RECOMMENDED<sup>1</sup>

n=131



**93%** of patients would recommend to others  
n=131

**90%** of partners would recommend to other couples in similar situations  
n=131

1. Levine LA, Estrada CR, Morgentaler A. Mechanical reliability and safety of, and patient satisfaction with the Ambicor inflatable penile prosthesis: results of a 2 center study. *J Urol.* 2001 Sep;166(3):932-7.

2. AMS Ambicor™ Inflatable Penile Prosthesis Operating Room Manual, American Medical Systems, Inc. 2013.

3. AMS Ambicor™ Inflatable Penile Prosthesis Instructions for Use, American Medical Systems, Inc. 2013.

All cited trademarks are the property of their respective owners.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material is not intended for use in France.

MH-394221-AB. Printed in the The Netherlands by De Budelse/Gosling.