



Ureterscopy Accessories

Leak-resistant accessories designed to make your cases flow smoothly



SAPS™ Single Action Pumping System

A pumping system designed for physicians or assistants to effectively control irrigation and clear the field of view with a simple one-handed press of the syringe.

- One-way valves with automatic refill of the syringe for convenient, immediate access to irrigation fluid

SAPS™ CF Continuous Flow version pictured



UroLok™ II Scope Adaptor

Our UroLok Scope Adaptor prevents backflow and secures slippery, coated devices. Designed to connect with irrigation bags and/or contrast syringe using a standard luer adapter. Additional product length helps keep connections out of the way during procedures.



Gateway™ Advantage Y-Adapter

A scope adaptor designed with a leak-resistant seal and two ports to smoothly and securely introduce multiple instruments into the scope.

- Main port allows passage of devices up to 8F
- Blue luer seal accepts devices up to 3F

Ordering Information

SAPS™ CF Continuous Flow Single Action Pumping System

Order Number	GTIN Number	Description	Unit
M0067201011	08714729045076	Single Action Pumping System CF	bx 5

System Includes: 10mL Vacuum Syringe, Continuous Flow Adjustable Length Tubing System with Dual-Check Valve and Roller Clamp

SAPS™ Single Action Pumping System

Order Number	GTIN Number	Description	Unit
M0067201001	08714729045069	Single Action Pumping System	bx 5

System Includes: 10mL Vacuum Syringe, Check Valve with Adjustable Length Connecting Tubes

UroLok™ II Scope Adaptor

Order Number	GTIN Number	Description	Unit
M0067301401	08714729066569	UroLok II Scope Adaptor	bx 10

Gateway™ Advantage Y-Adapter

Order Number	GTIN Number	Description	Unit
M0067301501	08714729845331	Gateway Advantage Y-Adapter	bx 10

Encore™ 26 Inflation Device

Order Number	GTIN Number	Description	Unit
M0067101131	08714729761815	20mL Encore 26 Inflation Device with Pressure Gauge	bx 5
M0067101140	08714729755814	20mL Encore 26 Inflation Device with Pressure Gauge	ea

Drainage Bag Connectors

Order Number	GTIN Number	Description	Unit
M0067401200	08714729010999	Connection for Female Luer Lock Fitting to Drainage Bag	ea
M0067401220	08714729011002	Connection for 3 through 6F Ureteral Catheters to Drainage Bag	ea

Ureteral Catheter Connectors

Order Number	GTIN Number	Description	Unit
M0067401100	08714729010982	Connection for 3 through 6F Catheters with Luer Lock Fitting	ea
M0067401110	08714729010975	Connection for 7 through 10F Catheters with Luer Lock Fitting	ea

Connecting Tubes for Drainage Bag

Order Number	GTIN Number	Description	Unit
M0067401000	08714729010944	30cm Connecting Tube for 3 through 6F Catheters	ea
M0067401010	08714729012832	30cm Connecting Tube for 7 through 10F Catheters	ea
M0067401020	08714729010951	30cm "Y" Connecting Tube for two 3 through 6F Catheters	ea

Applicable Medicare C-codes: There are no applicable Medicare C-codes for the above products.

To order, contact your sales representative.

**Boston
Scientific**
Advancing science for life™

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. Please check availability with your local sales representative or customer service.

All images are the property of Boston Scientific. All trademarks are the property of their respective owners.

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

© 2018 Boston Scientific Corporation
or its affiliates. All rights reserved.

URO-446715-AB SEP 2018